QUALIRIS QC EXPERIENCE

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QUALIRIS QC EXPERIENCE

Control plasmas for external quality assessment

• Kit containing 6 x 1-ml Vials

(REF 01046)

April 2019

English

1/ INTENDED USE

The **Qualiris QC Experience** kit provides unassayed plasmas intended for external quality assessment with the Qualiris by Stago programme.

The following parameters may be performed with each reagent:

- prothrombin time (PT)
- Owren's Prothrombin Time
- activated partial thromboplastin time (APTT)
- fibrinogen (Clauss' method)
- thrombin time (TT)
- Reptilase® time
- factors II, V, VII, VIII, IX, X, XI, XII, XIII
- antithrombin (AT)
- protein C
- protein S
- plasminogen
- antiplasmin
- von Willebrand factor (VWF:Ag, VWF:RCo).

PT, Owren's PT and APTT are general screening tests for haemostasis. TT and Reptilase® time explore the formation of fibrin. Factor XIII completes the cross-linking of fibrin. Factors VII, XII, XI, IX, VIII, V, X, II and fibrinogen are involved in the coagulation pathway. Antithrombin, protein C and protein S are physiological inhibitors of the coagulation process. Plasminogen and antiplasmin are involved in fibrinolysis. The von Willebrand factor is involved in primary haemostasis and coagulation.

2/ KIT REAGENTS

Kit of 6 vials of lyophilized human plasmas numbered from 1 to 6.

WARNING - POTENTIAL BIOHAZARDOUS MATERIAL

These reagents contain materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, approved methods are used to test the plasma for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

3/ PRECAUTIONS

Store at 2-8 °C. For *in vitro* use only. These reagents are to be used by certified medical laboratory personnel only. The disposal of waste materials must be carried out according to current local regulations. Exercise great care in the handling of these reagents and of patient samples.

4/ REAGENTS PREPARATION

Reconstitute each vial with exactly 1 ml of distilled water. Allow the solution to stand at room temperature (18-25 °C) for 30 minutes. Then, swirl the vial gently before use.

5/ REAGENTS STABILITY

The reagents in intact vials are stable until the expiration date indicated on the box label, when stored at 2-8 °C.

Once reconstituted, reagents remain stable for 4 hours at 20 ± 5 °C.

6/ REAGENTS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

- Reagents for performing tests listed in Section 1.
- Common clinical laboratory equipment and materials (distilled water...).

7/ DIRECTIONS FOR USE

The **Qualiris QC Experience** plasmas are to be used in the same manner as patients' plasmas according to the calendar below. Then report the results as indicated on the Qualiris QC Experience programme.

The S1 vials will be used for the campaign of the first half of the year and the S2 vials will be used for the campaign of the second half of the year.

	First half of the year	Second half of the year
Controls to test	1, 2 and 3	4, 5 and 6