

## CLOTEST\*: RAPID UREASE TEST FOR THE DIAGNOSIS OF H. PYLORI<sup>1</sup>

### ACCURACY


- 98% sensitivity and 97% specificity<sup>2</sup>
- Diagnosis of 75% of H. Pylori infections with no false positives within 20 minutes<sup>3</sup>

### CONVENIENCE

- Easy to read colour changes: from yellow to magenta
- Positive control test separately available:  
Jack Bean Urease Test

### AFFORDABILITY

- Less expensive diagnosis than culture or histology



# CLOTEST\*

## Rapid Urease Test

### Golden Rules

Refer to the CLOTEST\* Instructions for Use in the box for complete instructions.

**Storage:** Standard refrigerator range: 2-8°C

**Shelf Life:** 18 months


**Precautions:** Do not use if the gel is not yellow,  
seal is damaged,  
gel is dehydrated and  
expiration date has passed.

#### Preparation:

- Discontinue the use of antibiotics and bismuth preparations 3 weeks prior to the biopsy
- Put the slide 7-10 minutes at room temperature to obtain a faster result
- Make sure the biopsy sample is completely immersed in the gel
- Put the CLOTEST\* slide on a warming plate (30-40°C) to speed the reaction time (<3h)

#### Optimal sampling:

- The pyloric region (at least 2 cm away from the pylorus)
- The great or lesser curvature
- Avoid sampling from eroded tissue



#### Positive control test: Jack Bean Urease Test

In case the slide is still yellow after 24 hours:

- Insert one CLOTEST\* positive control sample in the gel
- Re-seal the lid over the slide and wait 10 minutes
- Look at the gel for a positive colour change (contact your sales representative in case of no colour change)



PRODUCT CODE	DESCRIPTION	PACKAGING
60480	CLOTEST* Rapid Urease Test	25 each
60407	Jack Bean Urease Control Tablets	50 each

References: 1. Rosaida M.S. et al. Evaluation of a new biopsy urease test for the diagnosis of Helicobacter Pylori infection. European Journal of Gastroenterology and Hepatology, 2004; 16: 195-199. 2. Dye KD, Marshall MJ, Frierson HF, Barrett LJ, Guerrant RL, McCallum RW. Is CLOtest alone adequate to diagnose Campylobacter pyloridis. Am J Gastroenterol 1988; 83: 1032 (abstract). 3. Marshall BJ, Warren JR, Francis GJ, Langton SR, Goodwin CS, Blincow E. Rapid Urease test in the management of Campylobacter pyloridis associated gastritis. Am J Gastroenterol 1987; 82(3): 200-210.

# AVANOS

For more information, please send an email to [customerservice.uk.ie@avanos.com](mailto:customerservice.uk.ie@avanos.com) or visit [www.avanos.co.uk](http://www.avanos.co.uk).



### Description

Item# 60480:

AVANOS CLOtest\* Rapid Urease Test (RUT) accurately and conveniently detects the urease enzyme of *Helicobacter pylori* in gastric mucosal biopsies.

Its use is intended for the presumptive diagnosis of *H. pylori* infection.

### Indication

The CLOtest\* RUT detects 75% of *H. pylori* infections within 20 minutes with no false positives.

By one hour 85% of positive patients will be detected by CLOtest\*.

and at three hours 90% are detected.

Easy to read colour changes.

Single use product.

### Counter indication

Do not reuse this device.

Consult the Instruction for Use.

### Main materials

Urease indicator gel sealed inside a plastic slide.

The gel contains urea, USP (29 mg/mL) phenol red (a pH indicator), buffers and a bacteriostatic agent to prevent the growth of contaminating urease-positive organisms.

Not formulated with Natural rubber latex.

Not formulated with DEHP.



### Packaging

Dispenser of 25 units CLOtest\*in plastic slides.

Box size: 7,94 x 15,0 x 4,44 cm (LxWxD)

Box weight: 181 g (0,400 lbs)

Bar coding: GS1-128 symbology, linear on carton box.

Item #	Single Unit	Carton box
60480	none	350770927459

### Manufacturing

Manufactured in USA.

The manufacturing site's quality is ISO13485 compliant.

### Regulatory information

Product CE marked as per 98/79/EEC Directive relative to In Vitro Diagnostic Medical Device.

Class of the device: I-(IVD).

### Storage

Store in a dry and cool place at 2°- 8°C (36°-48°F) temperature.

Do not use the product if the gel is not yellow, if the seal is damaged

and the gel appears dehydrated, or if the expiration date has passed.

Keep as much as practicably possible in its dispenser.

### Shelf life

18 months, from date of manufacture.

### Sterilisation

These products are non-sterile.

## EC Declaration of Conformity

<b>Name of Products:</b>	Diagnostic Products
<b>Product Codes:</b>	60480 – CLOtest* 60407 – CLOtest* Reagent Urease Type III
<b>GMDN Code:</b>	52787 (Helicobacter pylori urease IVD, kit, chromogenic)
<b>Legal Manufacturer (Place of Issue):</b>	Avanos Medical, Inc 5405 Windward Parkway Alpharetta, Georgia (GA) 30004, USA
<b>EU Authorized Representative:</b>	Avanos Medical Belgium BVBA Leonardo Da Vincilaan 1 1930 Zaventem Belgium
<b>Product Standards:</b>	N/A
<b>Start of CE:</b>	07/10/2019
<b>Conformity Assessment Route:</b>	Annex III
<b>Device Classification, Rules:</b>	General IVD
<b>CE certificate number:</b>	N/A (Self-certification)
<b>Notified Body:</b>	N/A (Self-certification)
<b>Quality System Certificate:</b>	FM 702907

I, the undersigned, hereby declare that the above specified medical devices meet the applicable provisions of European In Vitro Diagnostic Medical Devices Directive 98/79/EC and are in accordance with Annex III of the EC Directive, supported by the Conformity Assessment Procedure, and adhering to the essential requirements in accordance with Annex I of the Directive 98/79/EEC

This declaration is based on the existing Technical Documentation as per Annex III CE marking carried out as per Annex X of the Directive 98/79/EC

All supporting documentation that contains proof of compliance to the aforementioned Directives is retained under the premises of Avanos Medical, Inc This declaration applies to all devices from the signature date forward

**Authorized Signature:**

  
\_\_\_\_\_  
Name Thomas Kozma  
Title Director, Regulatory Affairs  
Avanos Medical, Inc  
Date 17 July 2019