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

PAXMAN SCALP COOLING SYSTEM **INSTRUCTIONS FOR USE**


CE
CE0086

UK ROW ISSUE 1 0718


PAXMAN°
PIONEERS IN SCALP COOLING

COOLING CAP INFORMATION

Patients' heads vary enormously in shape and size. We have developed a number of different sized caps to accommodate these varying sizes.

Each silicone cap has a serial number on the cap which can also be used to identify the size of the cap as the first letter(s) of this serial number will be either S, M, or L depending on the size of the cap. Fig 11.

1.4.3

The caps are made from medical grade soft silicone. They are light in weight but flexible and strong.

1.4.4

The silicone caps are supplied with a neoprene cover. The purpose of this cover is to absorb condensation from the external surface of the cap during use and prevent it from running down the back of patients' necks.

The neoprene cover also acts as an excellent insulator helping with the efficiency of the cooling cap.

The bungee cord and chin strap also ensure a snug, tight fit of the cooling cap. The cap covers are colour coordinated for size via the Paxman logo as shown in Fig 12.


Fig 11




Fig 12

CLEANING OF THE PAXMAN SCALP COOLING SYSTEM AND CAPS / CAP STORAGE

NOTE

Always read the product label of an EPA-registered product label thoroughly before use.

INFORMATION

Do not tumble dry the cap or the cap cover. Allowing the cap cover to dry naturally is preferred. Tumble drying may affect the lifespan of the cap cover.

CAUTION

The silicone cap and cap cover are supplied non-sterile. Neither the cap nor cap cover are suitable for cleaning by use of an autoclave.

CAUTION

It is NOT advisable to store caps hooked from the support arm by the chin strap.

Upon completion of scalp cooling and when the cap has been removed from the patient's head, the following procedure should be used to clean the cap.

Cleaning and Disinfection Procedure:

- Use the infection control procedure of the facility.
- The Paxman Cooling Cap should be thoroughly cleaned and disinfected before and after use.

Cleaning

- Remove cap cover from the silicone cap.
- Wash the internal and external surface of the cap with warm soapy tap water (100-110°F). *J. 23.*
- Rinse the cap with tap water to remove soap residues.
- Dry the cap thoroughly with a lint-free cloth.
- Use a 70% isopropyl alcohol (IPA) wipe to thoroughly wipe down the cooling cap.
- Allow the cooling cap to air dry.

Disinfecting

- Use a 10% bleach/water solution made fresh daily or purchase a bleach cloth/ disinfecting wipe (e.g., Super Sani-Cloth® & Sani-Cloth® HB Germicidal) to thoroughly wipe down the cooling cap.
- Allow the cooling cap to air dry.
- It is a violation of Federal Law to use an EPA registered product in a manner inconsistent with its label and labelling.
- Wipe over the internal and external surface of both the silicone cap and the cap cover (including the chinstrap) with a sanitizing wipe. Any hospital-grade range of disinfectant / cleaning wipes is suitable for use with the cap and cap cover. *J. 23.*
- In the event of the cap cover becoming soiled, the cap cover can be hand or machine washed in water no warmer than 40°C/104°F using standard mild detergent.
- Clean the external surfaces of the Paxman Scalp Cooling System machine using wipes dampened with soapy water of no more than 20°C/68°F or hospital-grade disinfecting wipes, ensuring liquid does not come into contact with the touchscreen display.
- When the units are not in use they should be stored in the plastic box provided.

QUICK RELEASE COUPLINGS

The system is designed for ease of use when connecting and disconnecting caps and coolant lines. The system incorporates simple to operate, quick release couplings. The design of the couplings is such that no spillage of coolant will occur.



Fig 9.

- When connecting the couplings, the male and female connectors should be held with one in each hand, and simply pushed firmly together to connect as shown in Fig 9. A click will be heard when the couplings are successfully connected.


Fig 10.

- When disconnecting the couplings, simply hold the couplings in both hands and press the white release button as shown in Fig 10. The couplings will spring apart. Continue to pull the couplings apart to complete the disconnect procedure.

CAP FITTING TO THE PATIENT

Fig 23.

Fig 24.

Fig 25.


Fig 26.


Fig 27.

Fig 28.

Option 2:

Fitting the cap excluding the cap cover.

Please note that all the standard protocols of patient preparation should be carried out before fitting the cap. It is also important to still utilize the headband or medical gauze with this cap to prevent direct contact between the silicone cap and bare skin.

- Fit the required inner silicone cap to the patients head, line the front of the silicon cap up with the hairline and around the ears Fig 30.
- Turn the outer cap cover inside out Fig 31.
- Place the inside out cap cover on top of the patients head and silicone cap placing your hand on to hold in place Fig 32.
- Roll down the cap cover ensuring all the inner silicone cap is totally covered Fig 33 & 24.
- Fully tighten the chin strap to set the cap in position Fig 25.
- Pull tight on the bungee cord at the back of the cap Fig 26 and lock in place utilising the bungee cord toggle Fig 27.
- Loosen off the chin strap to a more comfortable setting Fig 28.
- Connect the Cap to the Paxman System.

Fig 30.


Fig 31.



Fig 32.

Fig 33.

TECHNICAL DATA AND SPECIFICATIONS

Weight

Refrigeration unit 43kg

Caps between 795g and 900g (inc. coolant and cover)

Dimensions

H 640mm x W 320mm x D 420mm. Extended height 1650mm

Power supply

(100-120v or 220-240v) 50/60Hz single phase 850VA

All electrical circuits earthed and protected with fused system

All fuses are Bussmann quick blow fuses – 10 A at 250 Vac with F-speed

Condition of use

Operational

Temperature: 10-30°C / 50-86°F

Atmospheric pressure: 70-106 kPa

Relative humidity: 30-75%

Transport and Storage

Temperature: -40 to 70°C / -40 to 158°F

Atmospheric pressure: 50-106kPa

Relative humidity: 10-90%

Temperature control

Electronic thermostat with touch screen controller

Coolant flow rate and temperature continuously monitored in line

Alarms

Visual and audible for: high temperature conditions

Loss of coolant flow

Inlet pressure of coolant flow

Pressure into cap - 17psi

Flow rate of coolant 10-50 ml/s

REFRIGERATION - Hermetically sealed unit, CFC-free R134A refrigerant.

COOLANT - Orbis C. See data sheet section 10

Device classification

CE marked in accordance with annex V of the Medical Devices Directive 93/42/EEC for a Class IIa Device. Notified Body No. CE 0086

Safety Standards

BS EN 60601-1: 2006 + IEC60601-1-2006 for medical electrical equipment.

BS EN 60601-1-2: 2007 + IEC60601-1-2006 electromagnetic compatibility

Electrical Classification

Class 1 type BF

Product specifications are subject to continued evaluation and improvements and subject to change without prior notification.

CONTRAINDICATIONS / WARNINGS / SIDE EFFECTS

Side effects

-  Known side effects associated with scalp cooling therapy include:
- Discomfort due to feeling cold.
 - Headache.
 - Light-headedness or dizziness.
 - Nausea.

EQUIPMENT SUPPLIED

Equipment consists of:

Single Paxman Scalp Cooling System



1. Electronically operated refrigeration unit complete with 2m length power cord, 1 set of coolant lines with quick release couplings, covered with insulation and shroud, single adjustable coolant line support arm.
2. 1 x removable cap support,
3. 1 x Accessories box containing 3 standard cooling caps, 1 x small, 1 x medium and 1 x large* complete with insulation covers and quick release couplings.
4. 1 litre of coolant (top up)
5. Manual
6. Quick Reference Guide.
7. Patient information leaflets.

* It is optional to have any combination of caps


Single Paxman Scalp Cooling System

1. Electronically operated refrigeration unit complete with 2m length power cord, 2 sets of coolant lines with quick release couplings, covered with insulation and shroud, single adjustable coolant line support arm.
2. 2 x removable cap support
3. 1 x Accessories box containing 4 standard cooling caps, 2 x small, 1 x medium and 1 x large* complete with insulation covers and quick release couplings.
4. 1 litre of coolant (top up)
5. Manual.
6. Quick Reference Guide.
7. Patient information leaflets.

* It is optional to have any combination of caps

DEVICE DESCRIPTION

Paxman Scalp Cooling System (PSCS) is a compact, mobile refrigeration unit that has been tested to IEC 60601 safety standards for medical devices.

The unit circulates a refrigerated liquid coolant, at low pressure, through a special cooling cap that is positioned on the patient's head. The circulation of the refrigerated coolant through the cap extracts heat from the patient's scalp and maintains it at the optimum temperature for scalp cooling.

The Paxman Scalp Cooling System meets all the requirements for overall design, biocompatibility, and electrical safety. The results of the non-clinical testing confirm the output meets the design inputs and specifications. The following non-clinical testing was performed.

Biocompatibility

The system passed biocompatibility testing in accordance with the requirement for ISO 10993-1 and FDA Guidance. Specifically, Cytotoxicity testing per ISO 10993-5, Sensitization testing per ISO 10993-10, and Intracutaneous Reactivity testing per ISO 10993-10.

Electrical Safety Testing

The system passed electrical safety testing in accordance with requirements for IEC 60601-1 medical electrical equipment.

Electromagnetic Compatibility


The system passed electromagnetic compatibility (EMC) testing to meet requirements for IEC 60601-1-2 medical electrical equipment. There is no known equipment within the intended operational environment of the scalp cooling system that will interfere with, or be affected by the scalp cooling system.

Warning

-  Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Scalp Cooling System, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Electrical Safety

-  • The Paxman Scalp Cooling System is a Class 1 Electrical device and must be connected directly to a Protective Earth System.
-  • Extension Leads or Power Strips should NOT be used.
-  • Grounding reliability can only be achieved when equipment is connected to an equivalent receptacle marked as hospital grade.
-  • Only use an appropriately approved power cord for each country of use.
-  • The Plug is the main isolation device.
-  • Isolate the device by removing the plug before working on the device.

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COOLANT LEVEL

Coolant Dispenser

As far as the user is concerned the coolant level will only require topping up when the screen advises doing so.

When the screen advises the user to top up the coolant, Fig 53, it means that the coolant level has dropped below the high-level sensor.

Adding 1lt of coolant will bring the coolant level back to its optimum level.

Bespoke 1 lt dispensers of coolant are available to ensure that the device cannot be overfilled with coolant. Fig 52.


Fig 52


Fig 53


Fig 54


Fig 55

CAUTION

Use only the Paxman approved coolant and bag dispenser. DO NOT attempt to add coolant to the device in any other way except by the approved bag dispenser.

INFORMATION

The device can be topped-up with coolant when it is switched on or when it is switched off.

Topping Up the Coolant Level

The coolant top-up point at the back of the device, Fig 56, has a non return valve at the end to minimize spillage. This should be unclopped and the valve should be turned to point upwards for top-up Fig 54. This will allow the coolant dispenser to be securely attached.

Push the coolant dispenser connector firmly into the top up point valve until it is held secure.

Coolant will automatically start flowing into the device, Fig 55.

When the dispenser is empty, turn the top up point back to its original position with the valve pointing down whilst leaving the dispenser connected. Allow any residual coolant to flow back into the dispenser before removing it. This will close the valve and seal the device.

Secure the valve in its original position in the clip.

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