

# **cobas<sup>®</sup> 6800/8800 Systems**

# **Qualification Specifications**

Version 4.0

This page intentionally left blank.

## Table of Content

<b>PREFACE .....</b>	<b>5</b>
REVISION HISTORY .....	5
COPYRIGHTS AND TRADEMARKS .....	6
ABOUT THIS DOCUMENT .....	7
SCOPE .....	7
ABBREVIATIONS .....	7
REFERENCES .....	7
<b>IQ SPECIFICATIONS DSS SW 1.2 .....</b>	<b>8</b>
IQ TECHNICAL SPECIFICATIONS VERIFICATION .....	8
IQ: MOVABLE CHECK .....	9
IQ: INSTRUMENT BELT CHECK .....	9
IQ :MODULE CONNECTION CHECKS .....	10
<b>IQ SPECIFICATIONS DSS SW 1.3 AND 1.4 .....</b>	<b>11</b>
IQ TECHNICAL SPECIFICATIONS VERIFICATION .....	11
IQ MOVABLE CHECK .....	12
IQ INSTRUMENT BELT CHECK .....	12
IQ MODULE CONNECTION CHECK .....	13
<b>OQ SPECIFICATIONS DSS SW 1.2 .....</b>	<b>14</b>
OQ TRANSFER MODULE .....	14
<i>OQ Infrastructure</i> .....	14
<i>OQ Teaching and Calibration</i> .....	14
<i>OQ Rack handler</i> .....	14
<i>OQ Sample transfer head</i> .....	15
OQ PROCESSING MODULE A/B .....	15
<i>OQ Infrastructure</i> .....	15
<i>OQ Processing Module Fluidics</i> .....	15
<i>OQ Teachings</i> .....	15
<i>OQ Processing module heads</i> .....	16
OQ ANALYTIC MODULE .....	17
<i>OQ Teaching</i> .....	17
<i>OQ Handling Checks</i> .....	17
<i>OQ Analytic Unit A1/A2/B1/B2</i> .....	17
OQ CONSUMABLE LOADING AREA .....	18
<i>OQ Handling Checks</i> .....	18
<i>OQ Teaching</i> .....	18
<i>OQ Handling</i> .....	18
OQ INSTRUMENT .....	19
<i>OQ Safety interlock</i> .....	19
<i>OQ Handheld barcode reader</i> .....	19
OQ SAMPLE SUPPLY MODULE .....	19
OQ CHECK RUN .....	19
<b>OQ SPECIFICATIONS DSS SW 1.3 AND 1.4 .....</b>	<b>20</b>

Qualification Specifications

OQ PROCESSING MODULE A/B EXTENDED PRIME.....	20
OQ REAGENT DRAWER ADJUSTMENT .....	20
OQ INSTRUMENT TEACH.....	20
OQ INSTRUMENT TIGHTNESS CHECK .....	20
OQ RACK HANDLER CALIBRATION .....	20
OQ REAGENT TRANSFER A/B PIPETTING CHECKS .....	21
OQ ANALYTIC UNIT A1/A2/B1/B2 RAMP CHECK .....	21
OQ ANALYTIC UNIT A1/A2/B1/B2 FUNCTIONAL CHECK.....	22
OQ MAIN HANDLER HANDLING CHECK.....	22
OQ SAFETY INTERLOCK CHECK.....	22
OQ HANDHELD BARCODE READER .....	23
OQ SAMPLE SUPPLY MODULE.....	23
OQ CHECK RUN .....	23
<b>SWQ SW 1.2 SPECIFICATIONS.....</b>	<b>24</b>
SWQ: SAMPLE SUPPLY MODULE .....	24
SWQ: CHECK RUN .....	24
SWQ: POOL COMMUNICATION TEST RUN.....	24
<b>SWQ SW 1.3 SPECIFICATIONS.....</b>	<b>25</b>
SWQ: SAMPLE SUPPLY MODULE .....	25
<b>SWQ SW 1.4 SPECIFICATIONS.....</b>	<b>26</b>
SWQ: SAMPLE SUPPLY MODULE .....	26

## Preface

### Revision History

---

Version	Revision Date	Revision Information
1.0	10-Jun-2014	This is the first issue of this document
1.1	01-Dec-2014	Spelling corrections Changed “instrument” into “Systems”
1.2	27-Aug-2015	Added the Software Qualification Procedure to the reference Added SWQ Specifications section
2.0	18-Aug-2016	SWQ Specifications are differentiated by their related SWQ documents for SW 1.1 and 1.2 Editorial changes
3.0	01-Aug-2017	Added DSS SW version specific sections. Removed numbering for consistency reasons. Removed document versions under <i>References</i> . Combined Module Connection Checks. Reduced set of OQ Checks with DSS SW 1.3.
4.0	See GRIPS	Added SWQ Specifications for for SW 1.4 Removed SWQ Specifications for SW 1.1 Removed Check Run and Pool Communication Test Run from SWQ Specifications for SW 1.3 (Impact Assessment DH-03406.12-047) IQ and OQ Specifications for SW 1.1 removed.

---

Published by Roche Diagnostics International Ltd, a subsidiary of F. Hoffmann-La Roche Ltd, 4070 Basel, Switzerland.

Questions or comments regarding the content of this document can be directed to your Roche representative.

Every effort has been made to ensure that all the information contained in this document is correct at the time of printing.

However, Roche Diagnostics International Ltd reserves the right to make any changes necessary without notice as part of ongoing product development.

## **Copyrights and Trademarks**

The following trademarks are acknowledged:

**cobas** is a trademark of Roche.

© 2020, F. Hoffmann-La Roche Ltd  
All rights reserved

## About this Document

This document summarizes the acceptance criteria and specifications for the **cobas<sup>®</sup>** 6800/8800 Systems Installation Qualification (IQ), Operational Qualification (OQ) and Software Qualification (SWQ) for all supported software versions.

Complementary Information can be found in the IQ- and OQ-Certificate PDF files generated by the DSS Software.

## Scope

These specifications are only valid for the **cobas<sup>®</sup>** 6800/8800 Installation Qualification and Operational Qualification.

## Abbreviations

GRIPS	Global Repository of Information about Products and Services
iSDoc	integrated Service Documentation
SWQ	Software Qualification
IQ	Installation Qualification
OQ	Operational Qualification
Inacc.	Inaccuracy
Impr.	Imprecision
IG	Instrument Gateway Server
IC	Instrument Control
IM	Instrument Manager

## References

The following documents are available on GRIPS and refer to the related specifications:

Document
<b>cobas<sup>®</sup></b> 6800/8800 Systems - Installation Qualification (IQ)
<b>cobas<sup>®</sup></b> 6800/8800 Systems - Operational Qualification (OQ)
<b>cobas<sup>®</sup></b> 6800/8800 Systems Software 1.2 - Software Qualification (SWQ)
<b>cobas<sup>®</sup></b> 6800/8800 Systems Software 1.3 - Software Qualification (SWQ)
<b>cobas<sup>®</sup></b> 6800/8800 Systems Software 1.4 - Software Qualification (SWQ)

## IQ Specifications DSS SW 1.2

### IQ Technical Specifications Verification

Instrument Configuration	Conditions	Specifications
Stand-alone	Environmental Temperature	15°C – 32°C
	Relative Humidity	30% – 80% (25% – 80%; 18.2°C – 32°C) (20% – 80%; 21.7°C – 32°C) (no condensation)
With IG	Environmental Temperature	15°C – 28°C
	Relative Humidity	30% – 80% (25% – 80%; 18.2°C – 28°C) (20% – 80%; 21.7°C – 28°C) (no condensation)
With <b>cobas p 680</b>	Environmental Temperature	15°C – 30°C
	Relative Humidity	30% – 80% (25% – 80%; 18.2°C – 30°C) (20% – 80%; 21.7°C – 30°C) (no condensation)
Line Voltage <b>cobas®</b> 6800/8800 instrument	Phase 1	180VAC – 264VAC
	Phase 2	
	Phase 3	
Line Voltage Sample Supply Module	Plug 1 (Main)	180VAC – 264VAC
	Plug 2 (UPS)	



## Qualification Specifications

### IQ: Movable Check

Test	Test Detail	Specification
Set the IC configuration ID	Movable instrument configuration is set	Pass
Movable Platform	Movable functionalities working correctly	Pass

### IQ: Instrument Belt Check

Module	Unit	Belt	Specification
Transfer Module	Transfer Module Handler	X-Belt	32.0 – 36.0 Hz
		Y-Belt (left/right)	27.5 – 30.5 Hz
	Sample Transfer	Y-Belt (left/right)	30.0 – 33.0 Hz
	Rack Handler	X-Belt	65.5 – 75.5 Hz
		Y-Belt	30.0 – 34.0 Hz
		IN-Belt	37.0 – 45.0 Hz
Processing Module	Processing Module Handler	X-Belt (front/rear)	37.5 – 42.5 Hz
		Y-Belt	20.5 – 23.5 Hz
	Processing Transfer Unit	X-Belt	33.0 – 37.0 Hz
	Reagent Transfer Unit	Y-Belt	20.0 – 23.0 Hz
	Solid Waste Station	Y-Belt	50.5 – 57.0 Hz
Consumable Loading Area (cobas <sup>®</sup> 6800)	Main Handler	X-Belt	22.5 – 25.5 Hz
Consumable Loading Area (cobas <sup>®</sup> 8800)	Main Handler	X-Belt	13.0 – 15.0 Hz

**Qualification Specifications**

Analytic Module (cobas <sup>®</sup> 6800)	Amplification Plate Handler	X-Belt	54.0 – 66.0 Hz
		Y-Belt	40.0 – 50.0 Hz
		Z-Belt	18.0 – 22.0 Hz
Analytic Module (cobas <sup>®</sup> 8800)	Amplification Plate Handler	X-Belt	20.0 – 24.0 Hz
		Y-Belt	40.0 – 50.0 Hz
		Z-Belt	18.0 – 22.0 Hz

**IQ :Module Connection Checks**

Test	Test Detail	Specification
Transfer module connection check	All transfer module units reacting and communicating	Pass
Processing module A connection check	All processing module A units reacting and communicating	Pass
Processing module B connection check (cobas <sup>®</sup> 8800)	All processing module B units reacting and communicating	Pass
Analytic module connection check	All analytic module units reacting and communicating	Pass
Consumable loading Area connection check	All consumable loading area units reacting and communicating	Pass

## IQ Specifications DSS SW 1.3 and 1.4

### IQ Technical Specifications Verification

Instrument Configuration	Conditions	Specifications
Stand-alone	Environmental Temperature	15°C – 32°C
	Relative Humidity	30% – 80% (25% – 80%; 18.2°C – 32°C) (20% – 80%; 21.7°C – 32°C) (no condensation)
With IG	Environmental Temperature	15°C – 28°C
	Relative Humidity	30% – 80% (25% – 80%; 18.2°C – 28°C) (20% – 80%; 21.7°C – 28°C) (no condensation)
With <b>cobas p</b> 680	Environmental Temperature	15°C – 30°C
	Relative Humidity	30% – 80% (25% – 80%; 18.2°C – 30°C) (20% – 80%; 21.7°C – 30°C) (no condensation)
Line Voltage <b>cobas<sup>®</sup></b> 6800/8800 instrument	Phase 1	180VAC – 264VAC
	Phase 2	
	Phase 3	
Line Voltage Sample Supply Module	Plug 1 (Main)	180VAC – 264VAC
	Plug 2 (UPS)	

## Qualification Specifications

### IQ Movable Check

Test	Test Detail	Specification
Set the IC configuration ID	Movable instrument configuration is set	Pass
Movable platform	Movable functionalities working correctly	Pass

### IQ Instrument Belt Check

Module	Unit	Belt	Specification
Transfer Module	Transfer Module Handler	X-Belt	32.0 – 36.0 Hz
		Y-Belt (left/right)	27.0 – 31.0 Hz
	Sample Transfer	Y-Belt (left/right)	30.0 – 33.0 Hz
Processing Module	Processing Module Handler	X-Belt (front/rear)	37.0 – 43.0 Hz
		Y-Belt	21.0 – 23. Hz
	Reagent Transfer Unit	Y-Belt	20.0 – 23.0 Hz
Consumable Loading Area (cobas <sup>®</sup> 6800)	Main Handler	X-Belt	22.0 – 26.0 Hz
Consumable Loading Area (cobas <sup>®</sup> 8800)	Main Handler	X-Belt	12.0 – 16.0 Hz
Analytic Module (cobas <sup>®</sup> 6800)	Amplification Plate Handler	X-Belt	54.0 – 66.0 Hz
		Y-Belt	40.0 – 50.0 Hz
		Z-Belt	18.0 – 22.0 Hz
Analytic Module (cobas <sup>®</sup> 8800)	Amplification Plate Handler	X-Belt	20.0 – 24.0 Hz
		Y-Belt	40.0 – 50.0 Hz
		Z-Belt	18.0 – 22.0 Hz

**IQ Module Connection Check**

Test	Test Detail	Specification
Module Connection Check	The communication to all modules is verified	Pass

## **OQ Specifications DSS SW 1.2**

### **OQ Transfer Module**

#### **OQ Infrastructure**

<b>Test</b>	<b>Test Detail</b>	<b>Specification</b>
Fan speed	Ventilation fan transfer module	$\geq 2000$
Service light	Service light rotating	Pass
Module light	All light colors are visible.	Pass
Global indicator	All light colors are visible.	Pass
Airflow screen	Can move down and up	Pass
Fuses	All fuses are good	Pass
Temperature sensor	Temperature changed	Pass

#### **OQ Teaching and Calibration**

<b>Test</b>	<b>Test Detail</b>	<b>Specification</b>
Transfer module handler teach	Teach the transfer module handler	Pass
Rack handler teach	Teach the rack handler	Pass
Sample pipettors teach	Teach the sample pipettors	Pass
Rack handler calibration	Calibrate the rack handler camera	Pass

#### **OQ Rack handler**

<b>Test</b>	<b>Test Detail</b>	<b>Specification</b>
Rack handler check	All racks are read and moved	Pass

## Qualification Specifications

### OQ Sample transfer head

Test	Test Detail	Specification
Sample transfer head tightness check	Sample transfer head has no leaks	Pass

### OQ Processing Module A/B

#### OQ Infrastructure

Test	Test Detail	Specification
Fan speed	Ventilation fan processing module	$\geq 2000$
Service light	Service light rotating	Pass
Module light	All light colors are visible.	Pass
Global indicator	All light colors are visible.	Pass
Airflow screen	Can move down and up	Pass
Fuses	All fuses are good	Pass
Temperature sensor	Temperature changed	Pass

#### OQ Processing Module Fluidics

Test	Test Detail	Specification
Extended prime	Prime without error and leaks	Pass without any leaks

#### OQ Teachings

Test	Test Detail	Specification
Processing module handler teach	Teach processing module handler	Pass
Reagent transfer head teach	Teach reagent transfer head	Pass
Processing transfer head teach	Teach processing transfer head	Pass
Interim reagent storage	Teach interim reagent storage	Pass

**Qualification Specifications**


---

**OQ Processing module heads**

Test	Test Detail	Specification
Reagent transfer tightness check	Perform tightness check	Pass
Processing transfer 1/2 tightness check	Perform tightness check	Pass
Reagent head balance calibration	The OWS is calibrated with a 10g weight	Calibration is OK
Reagent head pipetting check (Reagent head pipetting A-H)	Initial pipetting 800µl	Inacc: <3%
	Pipetting 50µl (high volume)	Inacc: <3.8% Impr: <1.25%
	Pipetting 10µl (low volume)	Inacc: <7.5% Impr: <3.75%



## OQ Analytic Module

### OQ Teaching

Test	Test Detail	Specification
Amplification plate handler teach	Teach the plate handler	Pass

### OQ Handling Checks

Test	Test Detail	Specification
Analytic module handling check	Move the amplification plate to the different position in the analytic module	Pass

### OQ Analytic Unit A1/A2/B1/B2

Test	Test Detail	Specification
Analytic unit ramp check	Check the cyclor block functionality	Pass
Analytic unit functional check	Channels imprecision	<0.5%
	Channels inhomogeneity	<35%
	Min Excitation power channel 1	>=5170
	Max Excitation power channel 1	<=35159
	Min Excitation power channel 2	>=1650
	Max Excitation power channel 2	<=5552
	Min Excitation power channel 3	>=5830
	Max Excitation power channel 3	<=18426
	Min Excitation power channel 4	>=11770
	Max Excitation power channel 4	<=35637
	Min Excitation power channel 5	>=1980
	Max Excitation power channel 5	<=5866

## **OQ Consumable Loading Area**

### **OQ Handling Checks**

<b>Test</b>	<b>Test Detail</b>	<b>Specification</b>
Reagent drawer adjustment	Drawer adjustment	Pass
Reagent storage check	Checks the reagent and control cassette handling in the reagent storage	Pass

### **OQ Teaching**

<b>Test</b>	<b>Test Detail</b>	<b>Specification</b>
Calibrate main handler	Calibrates the main handler	Pass
Consumable loading area teach	Teach the consumable area on all positions	Pass

### **OQ Handling**

<b>Test</b>	<b>Test Detail</b>	<b>Specification</b>
Main handler handling check	Handling of the different consumables through the instrument	Pass

## **OQ Instrument**

### **OQ Safety interlock**

Test	Test Detail	Specification
Safety interlock check	All interlock unit are functioning correctly	Pass

### **OQ Handheld barcode reader**

Test	Test Detail	Specification
Handheld barcode reader	Read tube and rack barcode	The information read corresponds to the tube and rack barcode

### **OQ Sample Supply Module**

Test	Test Detail	Specification
Sample supply module time synchronization	Synchronize the sample supply module	Time of the sample supply module corresponds to the IM.
Sample supply module check	Check all sensor and mechanical checks	Pass

### **OQ Check Run**

Test	Test Detail	Specification
Check Run	System mechanical, electrical, and communication test.	Run completed without error messages, batch, control, and final sample results communicated throughout the system.

## OQ Specifications DSS SW 1.3 and 1.4

### OQ Processing module A/B extended prime

Test	Test Detail	Specification
Processing module A/B extended prime	Prime without error and leaks	Pass without any leaks

### OQ Reagent Drawer Adjustment

Test	Test Detail	Specification
Reagent Drawer Adjustment	Drawer adjustment	Pass

### OQ Instrument Teach

Test	Test Detail	Specification
Instrument Teach	Teaches all reference points in all modules	Pass

### OQ Instrument Tightness check

Test	Test Detail	Specification
Instrument Tightness check	Checks the tightness of the Sample transfer, the Processing transfer Reagent transfer	Pass

### OQ Rack handler calibration

Test	Test Detail	Specification
Rack handler calibration	Calibrates the rack handler camera	Pass

**Qualification Specifications**


---

**OQ Reagent transfer A/B pipetting checks**

Test	Test Detail	Specification
Reagent head pipetting intro	OWS is positioned	Pass
Reagent head balance calibration	The OWS is calibrated with a 10g weight	Pass
Reagent head A/B pipetting checks (Reagent head pipetting A-H)	Initial pipetting 800µl	Inacc: <3%
	Pipetting 50µl (high volume)	Inacc: <3.8% Impr: <1.25%
	Pipetting 10µl (low volume)	Inacc: <7.5% Impr: <3.75%

**OQ Analytic unit A1/A2/B1/B2 ramp check**

Test	Test Detail	Specification
Analytic unit ramp check	Check the cyclor block functionality	Pass

## Qualification Specifications

### OQ Analytic unit A1/A2/B1/B2 functional check

Test	Test Detail	Specification
Analytic unit functional check	Channels imprecision	<0.5%
	Channels inhomogeneity	<35%
	Min Excitation power channel 1	>=5170
	Max Excitation power channel 1	<=35159
	Min Excitation power channel 2	>=1650
	Max Excitation power channel 2	<=5552
	Min Excitation power channel 3	>=5830
	Max Excitation power channel 3	<=18426
	Min Excitation power channel 4	>=11770
	Max Excitation power channel 4	<=35637
	Min Excitation power channel 5	>=1980
	Max Excitation power channel 5	<=5866

### OQ Main handler handling check

Test	Test Detail	Specification
Main handler handling check	Handling of the different consumables through the instrument	Pass

### OQ Safety interlock check

Test	Test Detail	Specification
Safety interlock check	All interlock unit are functioning correctly	Pass

**Qualification Specifications**

---

**OQ Handheld barcode reader**

Test	Test Detail	Specification
Handheld barcode reader	Read tube and rack barcode	The information read corresponds to the tube and rack barcode

**OQ Sample Supply Module**

Test	Test Detail	Specification
Sample supply module time synchronization	Synchronize the sample supply module	Time of the sample supply module corresponds to the IM.
Sample supply module check	Check all sensor and mechanical checks	Pass

**OQ Check Run**

Test	Test Detail	Specification
Check Run	System mechanical, electrical, and communication test.	Run completed without error messages, batch, control, and final sample results communicated throughout the system.

## SWQ SW 1.2 Specifications

### SWQ: Sample Supply Module

Test	Test Detail	Specification
Sample supply module time synchronization	Synchronize the sample supply module	Time of the sample supply module corresponds to the IM.
Sample supply module check	Check all sensor and mechanical checks	Pass

### SWQ: Check Run

Test	Test Detail	Specification
Check Run	System mechanical, electrical, and communication test.	Run completed without error messages, batch, control, and final sample results communicated throughout the system.

### SWQ: Pool Communication Test Run

Test	Test Detail	Specification
Pool Communication Test Run	Data communication with Instrument Gateway	Pass



## SWQ SW 1.3 Specifications

### SWQ: Sample Supply Module

Test	Test Detail	Specification
Sample supply module time synchronization	Synchronize the sample supply module	Time of the sample supply module corresponds to the IM.
Sample supply module check	Check all sensor and mechanical checks	Pass

# SWQ SW 1.4 Specifications

## SWQ: Sample Supply Module

Test	Test Detail	Specification
Sample supply module time synchronization	Synchronize the sample supply module	Time of the sample supply module corresponds to the IM.