MAGLUMI[®] SARS-CoV-2 S-RBD IgG (CLIA)

INTENDED USE

The kit is an *In Vitro* chemiluminescence immunoassay for the quantitative determination of S-RBD IgG antibodies to SARS-CoV-2 in human serum and plasma using the MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer.

SUMMARY AND EXPLANATION OF THE TEST

SARS-CoV-2, formerly called 2019-nCoV, is a betacoronavirus and cause an acute respiratory disease (coronavirus disease-2019, COVID-19). Coronavirus is an envelope virus with four structural proteins: spike (S) protein, membrane (M) protein, envelope (E) protein and nucleocapsid (N) protein. Among the 4 coronavirus structural proteins, the spike (S) and nucleocapsid (N) proteins are the main immunogens. Specially, S protein is a major protective antigen that elicits highly potent neutralizing antibodies (nAbs) and plays an essential role in viral attachment, fusion, entry, and transmission. S protein comprises an N-terminal S1 subunit responsible for virus-receptor binding and a C terminal S2 subunit responsible for virus-cell membrane fusion. S1 is further divided into an N-terminal domain (NTD) and a receptor-binding domain (RBD), RBD within S1 directly interacts with host receptors. Human angiotensin converting enzyme 2 (hACE2) is the receptor to which SARS-CoV-2 binds to enter host cells.

nAbs are antibodies that are capable of keeping an infectious agent (for instance, virus) from infecting a cell by neutralizing or inhibiting its biological effect. RBD within S1 subunit is the most critical target for SARS-CoV nAbs. Such nAbs can interrupt the interaction of RBD and its receptor ACE2. Thus, SARS-CoV-2 S-RBD IgG antibody level in human serum or plasma correlates with protective immune responses in individuals who have recovered from SARS-CoV-2 infection and also reflect herd immunity at a population level, informing clinical management of patients with past or ongoing COVID-19 infection.

PRINCIPLE OF THE TEST

The SARS-CoV-2 S-RBD IgG (CLIA) assay is an indirect chemiluminescence immunoassay.

The sample, buffer and magnetic microbeads coated with S-RBD recombinant antigen are mixed thoroughly and incubated, forming immune-complexes. After precipitation in a magnetic field, decant the supernatant, and perform a wash cycle. Then add ABEI labeled with anti-human IgG antibody, and incubate to form complexes. After precipitation in a magnetic field, decant the supernatant, and perform another wash cycle. Subsequently, the Starter 1+2 are added to initiate a chemiluminescent reaction. The light signal is measured by a photomultiplier as relative light units (RLUs), which is proportional to the concentration of SARS-CoV-2 S-RBD IgG presented in the sample.

KIT COMPONENTS

Material Provided

| Component | Contents | 100 tests (REF:130219017M) | 50 tests (REF:130619017M) |
|----------------------------|---|-------------------------------|------------------------------|
| Magnetic Microbeads | Magnetic microbeads coated with recombinant SARS-CoV-2 S-RBD antigen, PBS buffer, NaN ₃ (<0.1%). | 2.5 mL | 1.5 mL |
| Calibrator Low | Low concentration of SARS-CoV-2 S-RBD IgG antibodies, PBS buffer, NaN ₃ (<0.1%). | 1.0 mL | 1.0 mL |
| Calibrator High | High concentration of SARS-CoV-2 S-RBD IgG antibodies, PBS buffer, NaN ₃ (<0.1%). | 1.0 mL | 1.0 mL |
| Buffer | NaCl and BSA, NaN ₃ (<0.1%). | 22.5 mL | 12.0 mL |
| ABEI Label | Anti-human IgG antibody labeled with ABEI, Tris-HCl buffer, NaN ₃ (<0.1%). | 22.5 mL | 12.0 mL |
| Diluent | PBS buffer, NaN ₃ (<0.1%). | 5.5 mL | 3.5 mL |
| Control 1 | SARS-CoV-2 S-RBD IgG antibodies, PBS buffer, NaN ₃ (<0.1%). | 1.0 mL | 1.0 mL |
| Control 2 | SARS-CoV-2 S-RBD IgG antibodies, PBS buffer, NaN ₃ (<0.1%). | 1.0 mL | 1.0 mL |

All reagents are provided ready-to-use.

Required But Not Provided

| | | |
|-----------------|-------------------|-----------------------------|
| MAGLUMI Series: | Reaction Module | REF: 630003 |
| | Starter 1+2 | REF: 130299004M; 130299027M |
| | Wash Concentrate | REF: 130299005M |
| | Light Check | REF: 130299006M |
| | Reaction Cup | REF: 130105000101 |
| | Maglumi 600 | REF: 23020018 |
| | Maglumi 800 | REF: 23020003 |
| | Maglumi 1000 | REF: 23020009 |
| | Maglumi 2000 | REF: 23020006 |
| | Maglumi 2000 Plus | REF: 23020007 |
| | Maglumi 4000 | REF: 23020014 |
| | Maglumi 4000 Plus | REF: 23020037 |
| | MAGLUMI X8 | REF: 010101008801 |

Please order all above from Shenzhen New Industries Biomedical Engineering Co., Ltd. (SNIBE) or our authorized representative.

CALIBRATION

Traceability: This method has been standardized against the SNIBE internal reference substance.

Test of assay specific calibrators allows the RLU values to adjust the assigned master curve. Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent Radio Frequency Identification (RFID) CHIP.

Recalibration is recommended if any of the following conditions occurs:

- After each exchange of lots (Reagent or Starter 1+2).
- Every week and/or each time a new reagent kit is used.
- After instrument service is required.
- If controls lie outside the expected range.

QUALITY CONTROL

Follow government regulations or accreditation requirements for quality control frequency.

The controls in kit are only applicable with MAGLUMI system and used matching with the same LOT numbers of corresponding reagents. For each target value and range refer to the label.

To monitor system performance, quality control materials (negative control and positive control) are required. Treat all quality control samples with the same level of care as patient samples. A satisfactory level of performance is achieved when analyte values obtained are within the acceptable Control Range for the system or within your range, as determined by an appropriate internal laboratory quality control scheme. If the quality control results do not fall within the Expected Range or within the laboratory's established Range, measurement of the quality control should be repeated. If the quality control results still do not fall within the range, do not report results and take the following actions:

- Verify that the materials are not expired.
- Verify that required maintenance was performed.
- Verify that the assay was performed according to the instruction for use.
- Rerun the assay with fresh quality control samples.
- If necessary, contact your local technical support provider or distributor for assistance.

SPECIMEN COLLECTION AND PREPARATION

- Human serum or plasma may be used with the SARS-CoV-2 S-RBD IgG (CLIA) assay. For serum samples, tubes without additive/accessory, or tubes containing clot activator or clot activator with gel could be applied for the assay. For plasma samples, the anticoagulants including K2-EDTA, K3-EDTA, Na2-EDTA were tested and found acceptable.
- Please pay attention to the risk of infection during sample collection and preparation.
- Ensure that complete clot formation in specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy, may exhibit increased clotting time.
- If the specimen is centrifuged before a complete clotting, the presence of fibrin may cause erroneous results. Samples must be free of fibrin and other particulate substance.
- Do not use grossly hemolyzed specimens as well as specimens containing particulate matter or exhibiting obvious microbial contamination. Inspect all specimens for bubbles, and remove bubbles before analysis for optimal results.
- Specimens removed from the separator gel, cells or clot may be stored up to 3 days at 2-8°C, or 2 months frozen at -20°C.
- Frozen and thawed for 2 times. Frozen specimens must be mixed thoroughly after thawing by low speed vortexing or by gently inverting.
- For optimal results, specimens should be free of fibrin, red blood cells, or other particulate matter. Such specimens may give inconsistent results and must be transferred to a centrifuge tube and centrifuged at ≥ 10,000RCF (Relative Centrifugal Force) for 10 minutes. Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.
- Before shipping specimens, it is recommended that specimens be removed from the separator, red blood cells or clot. When shipped, specimens should be packaged and labeled in compliance with applicable state, federal and international regulations covering the transport of clinical specimens and infectious substances. Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances. Do not exceed the storage limitations listed above.
- The sample volume required for a single determination is 10 µL.

WARNING AND PRECAUTIONS FOR USERS

- **IVD** For *In Vitro* Diagnostic Use.
- Follow the package insert carefully. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

Safety Precautions

- **CAUTION:** This product requires the handling of human specimens. It is recommended that all human sourced materials be considered potentially infectious and handled in accordance with the 29 CFR 1910.1030 Occupational exposure to bloodborne pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.
- All samples, biological reagents and materials used in the assay should be considered potentially able to transmit infectious agents. They should therefore be disposed of in accordance with the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with prevailing regulatory requirements.
- This product contains Sodium Azide. Dispose of contents and container must be in accordance with all local, regional and national regulations.
- Refer to safety data sheets, which are available on request.

Handling Precautions

- Do not use reagent kits beyond the expiration date.
- Do not interchange reagent components from different reagents or lots.
- Prior to loading the Reagent Kit on the system for the first time, the Reagent Kit requires mixing to re-suspend magnetic microbeads that have settled during shipment. For magnetic microbeads mixing instructions, refer to the Preparation of the Reagent section of this package insert.
- To avoid contamination, wear clean gloves when operating with a reagent kit and sample.
- Over time, residual liquids may dry on the septum surface. These are typically dried salts which have no effect on assay efficacy.
- To avoid evaporation of the liquid in the opened reagent kits in refrigerator, it is recommended that the opened reagent kits to be sealed with reagent seals contained within the packaging. The reagent seals are "single use", and if more seals are needed, please contact SNIBE or our authorized representative.
- For detailed discussion of handling precautions during system operation, refer to the SNIBE service information.

STORAGE AND STABILITY

- Store at 2-8°C. Do not freeze.
- Keep upright for storage to facilitate later proper resuspension of magnetic microbeads.
- Keep away from sunlight.
- The stability study is still on-going, the following data is obtained by referring to similar products:

| Stability of the reagent | |
|--------------------------|----------------------------------|
| unopened at 2-8°C | until the stated expiration date |
| opened at 2-8°C | 6 weeks |
| onboard | 4 weeks |

- To ensure the best kit performance, it is recommended to place opened kits in the refrigerator after the end of the intraday test work.
- The controls in kit can be stable until the expiration date at 2-8°C in an unopened state, and can be stable for 6 weeks when properly stored at 2-8°C after opening, or stored up to 2 months frozen at -20°C or colder. Frozen and thawed for 2 times.

TEST PROCEDURE

Preparation of the Reagent

- Take the reagent kit out of the box and observe the sealing film and other parts of the reagent kit to see if there is any leakage. In case of leakage, please contact your local distributor immediately. And then tear off the kit sealing film carefully.
- Open the reagent area door; Hold the reagent handle to get the RFID label close to the RFID reader (for about 2s); the buzzer will beep; one beep sound indicates successful sensing.
- Keeping the reagent straight insert to the bottom along the blank reagent track.
- Observe whether the reagent information is displayed successfully in the user software interface, otherwise repeat the above steps.
- Resuspension of the magnetic microbeads takes place automatically when the kit is loaded successfully, ensuring the magnetic microbeads are totally resuspended homogenous prior to use.

Assay Calibration

- Click **<Calibration>** or **<Batch Calibration>** button to execute calibration operation; For specific information on ordering calibrations, refer to the Calibration Section of the Operating Instructions.
- Execute recalibration according to the calibration interval required in this package insert.

Quality Control

- In order to avoid manually error in entry of QC information, the provided barcode labels of quality control in the kit can be used attached on the test tubes.
- If users do not use the provided barcode labels for positive and negative controls contained within the packaging, then quality controls should be ordered manually.
- The controls in kit must be edited as "CONTROL" in MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer user software. For specific information on ordering quality controls, refer to the Quality Control Section of the Operating Instructions.

Sample Testing

- Order the samples in the Sample Area of the software and click the **<Start>** button to execute testing. For specific information on ordering patient specimens, refer to the Sample Ordering Section of the Operating Instructions.

To ensure proper test performance, strictly adhere to the operating instructions of MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer.

DILUTION

The high concentration samples can be diluted automatically by analyzers or manually. The recommended dilution is 1:9 with the Diluent in kit. After manual dilution, multiply the result by the dilution factor. After dilution by the analyzer, the analyzer software automatically takes the dilution into account when calculating the sample concentration.

The automatic sample dilution is available after dilution settings are done in the MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer user software. Please refer to the operating instructions of MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer.

LIMITATIONS

- This test is suitable only for investigating single samples, not for pooled samples.
- Bacterial contamination or repeated freeze-thaw cycles may affect the test results.
- Assay results should be utilized in conjunction with other clinical and laboratory methods to assist the clinician in making individual patient diagnostic decisions.
- Assay results should not be used as the sole basis for the diagnosis and exclusion of novel coronavirus pneumonia, but only as a supplement to existing viral nucleic acid detection reagents and imaging features.
- If the SARS-CoV-2 S-RBD IgG results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.
- HAMA antibodies in test samples may cause interference in immunoassays.^{1,2}

RESULTS

Calculation of Results

The analyzer automatically calculates the concentration in each sample by means of a calibration curve which is generated by a 2-point calibration master curve procedure. The results are expressed in AU/mL. For further information please refer to the operating instructions of MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer.

Interpretation of Results

Results of reference range study in China using the SARS-CoV-2 S-RBD IgG assay was as follows:

- Non-reactive: A result less than 1.00 AU/mL (<1.00 AU/mL) is considered to be non-reactive.
- Reactive: A result greater than or equal to 1.00 AU/mL (≥1.00 AU/mL) is considered to be reactive.
- For samples with concentration near the cut-off or positive, follow-up tests should be performed.
- Results may differ between laboratories due to variations in population. It is recommended that each laboratory establish its own expected ranges.

PERFORMANCE CHARACTERISTICS

Precision

Precision for SARS-CoV-2 S-RBD IgG assay was determined as described in the CLSI EP05-A3. 2 controls and 3 human serum pools containing different concentration of analyte were assayed in duplicate at three sites on five days, with 3 runs per day, one lot of reagent for each run and 2 replicates per run. The result is summarized in the following table:

| Sample | Mean Value (AU/mL) | N | Repeatability | | Between-Lot | | Between-Day | | Between-Site | | Reproducibility | |
|--------|--------------------|----|---------------|------|-------------|------|-------------|------|--------------|------|-----------------|-------|
| | | | SD (AU/mL) | %CV | SD (AU/mL) | %CV | SD (AU/mL) | %CV | SD (AU/mL) | %CV | SD (AU/mL) | %CV |
| QC1 | 0.396 | 90 | 0.007 | 1.77 | 0.016 | 4.04 | 0.039 | 9.85 | 0.024 | 6.06 | 0.049 | 12.37 |
| QC2 | 3.916 | 90 | 0.176 | 4.49 | 0.078 | 1.99 | 0.054 | 1.38 | 0.073 | 1.86 | 0.213 | 5.44 |
| S1 | 0.550 | 90 | 0.042 | 7.64 | 0.009 | 1.64 | 0.049 | 8.91 | 0.000 | 0.00 | 0.065 | 11.82 |
| S2 | 2.421 | 90 | 0.099 | 4.09 | 0.018 | 0.74 | 0.066 | 2.73 | 0.167 | 6.90 | 0.206 | 8.51 |
| S3 | 5.109 | 90 | 0.115 | 2.25 | 0.046 | 0.90 | 0.043 | 0.84 | 0.213 | 4.17 | 0.250 | 4.89 |

Limit of Blank (LoB)

The LoB for SARS-CoV-2 S-RBD IgG assay is 0.100 AU/mL.

Limit of Detection (LoD)

The LoD for SARS-CoV-2 S-RBD IgG assay is 0.180 AU/mL.

Linear Range

0.180-100 AU/mL (defined by the Limit of Detection and the maximum of the master curve).

Interference

Two samples (one negative sample, one positive sample) were spiked with potential endogenous interference and exogenous interference. The results of the interferences are listed in the following table:

| Interference | No interference up to | Interference | No interference up to | Interference | No interference up to |
|-----------------------------|-----------------------|-------------------------|-----------------------|---------------------------|-----------------------|
| Bilirubin | 40 mg/dL | Levofloxacin | 1.776 mg/dL | Mometasone | 2.5 mg/dL |
| Triglycerides | 1000 mg/dL | Azithromycin | 1.201 mg/dL | Budesonide | 3.2 mg/dL |
| Hemoglobin | 2000 mg/dL | Ribavirin | 90 mg/dL | Mucin | 260 mg/dL |
| HAMA | 30 ng/mL | Meropenem | 80.15 mg/dL | Zanamivir | 1.2 mg/dL |
| RF | 1500 IU/mL | Tobramycin | 2.4 mg/dL | Peramivir | 60 mg/dL |
| ANA | 400 AU/mL | Oseltamivir | 1.0 mg/dL | Lopinavir | 48 mg/dL |
| Anti-Mitochondrial | 1:64(titer) | Oxymetazoline | 2.5 mg/dL | Ritonavir | 120 mg/dL |
| Total IgG | 1600 mg/dL | Sodium chloride | 45 mg/dL | Arbidol | 36 mg/dL |
| Total IgM | 280 mg/dL | Beclomethasone | 2.5 mg/dL | Flunisolide | 2.5 mg/dL |
| Interferon α | 1500 U/mL | Dexamethasone | 18 mg/dL | Histamine dihydrochloride | 4.5 mg/dL |
| Phenylephrine hydrochloride | 1.0 mg/dL | Triamcinolone acetonide | 5.5 mg/dL | Biotin | 5.0 mg/dL |
| Fluticasone propionate | 2.5 mg/dL | Ceftriaxone sodium | 81.03 mg/dL | - | - |

Cross-Reactivity

The cross-reactivity study for the SARS-CoV-2 S-RBD IgG (CLIA) assay was designed to evaluate potential cross reactants. The results are listed in the following table:

| Category | N of samples | Reactive | Category | N of samples | Reactive |
|---|--------------|----------|-----------------------------------|--------------|----------|
| Human Coronavirus antibodies (HKU1, OC43, NL63, 229E) | 24 | 0 | Measles virus antibodies | 8 | 0 |
| Influenza A virus antibodies | 30 | 0 | CMV antibodies | 9 | 0 |
| Influenza B virus antibodies | 13 | 0 | Rotavirus antibodies | 6 | 0 |
| Respiratory syncytial virus antibodies | 7 | 0 | Norovirus antibodies | 6 | 0 |
| Rhinovirus antibodies | 18 | 0 | Mumps virus antibodies | 6 | 0 |
| Adenovirus antibodies | 42 | 0 | Varicella zoster virus antibodies | 7 | 0 |

| | | | | | |
|------------------------|----|---|------------------------------|----|---|
| Enterovirus antibodies | 24 | 0 | M.Pneumoniae antibodies | 7 | 0 |
| EB virus antibodies | 16 | 0 | Human immunodeficiency virus | 10 | 0 |

High-Dose HOOK

No high-dose hook effect was seen when for SARS-CoV-2 S-RBD IgG concentrations up to 3500 AU/mL

Clinical Sensitivity

The clinical sensitivity of the SARS-CoV-2 S-RBD IgG assay was determined in China by testing 351 samples confirmed COVID-19 infected specimens.

| Days Post Onset of Symptoms | N of samples | Reactive | Sensitivity | 95% CI |
|-----------------------------|--------------|----------|-------------|--------------|
| 0-7 | 55 | 41 | 74.5% | 63.0%-86.1% |
| 8-14 | 94 | 92 | 97.9% | 95.0%-100.0% |
| ≥ 15 | 202 | 202 | 100.0% | 99.9%-100.0% |

The positive rate of S-RBD IgG antibodies may be affected by the infection period of the test subject (when blood sampling) in different studies.

Clinical Specificity

The clinical specificity of the SARS-CoV-2 S-RBD IgG assay was determined in China by testing 229 samples collected from suspected COVID-19 patients with PCR confirmation of negative result.

| N of samples | Non-reactive | Specificity | 95% CI |
|--------------|--------------|-------------|--------------|
| 229 | 228 | 99.6% | 98.7%-100.0% |

REFERENCES

1. Primus F J, Kelley E A, Hansen H J, et al. "Sandwich"-type immunoassay of carcinoembryonic antigen in patients receiving murine monoclonal antibodies for diagnosis and therapy. *Clinical Chemistry*, 1988, 34(2):261-264.
2. Schroff R W, Foon K A, Beatty S M, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Research*, 1985, 45(2): 879-885.

Shenzhen New Industries Biomedical Engineering Co., Ltd.

No.23, Jinxiu East Road, Pingshan District, 518122 Shenzhen, P.R. China

Tel: +86-755-21536601

Fax: +86-755-28292740

Shanghai International Holding Corp. GmbH (Europe)

Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

Fax: +49-40-255726

SYMBOLS EXPLANATIONS

| | | | |
|--|---------------------------------------|--|---|
| | Consult instructions for use | | Manufacturer |
| | Temperature limit (Store at 2-8°C) | | Use-by date |
| | Contains sufficient for <n> tests | | Keep away from sunlight |
| | This way up | | Authorized representative in the European Community |
| | In vitro diagnostic medical device | | Kit components |
| | Catalogue number | | Batch code |

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

FOR PROFESSIONAL USE ONLY

Store at room temperature (15-30°C)

CAUTION: COMPLETELY READ THE
INSTRUCTIONS BEFORE PROCEEDING

SYMBOLS EXPLANATIONS

Authorised representative in the
European Community

Manufacturer

Consult instructions for use

Kit components

In vitro diagnostic medical device

Batch code

Catalogue number

Use-by date

Temperature limit
(Store at room temperature (15-30°C))

Contains sufficient for

Keep away from sunlight

1. Intended Use

Reaction Modules required for the implementation of the MAGLUMI assays on the MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer (Including Maglumi 600, Maglumi 800, Maglumi 1000, Maglumi 1000 Plus, Maglumi 2000, Maglumi 2000 Plus, Maglumi 4000 and Maglumi 4000 Plus).

2. Material provided

6x64 reaction modules (with 6 cavities each).

For additional details such as proper test performance of MAGLUMI assays, please refer to the instructions for the assay to be performed.

3. Handling

Loading:

- Place a maximum of 8 reaction modules on the transport belt as illustrated below.
- Make sure that the transport belt is not moving.
- The correct position of the reaction modules is defined by a notch in the module and the guiding rail on the loading station.

Modules are for single use only

Front View

↑
Positioning notch

↓ ↓
Transport belt to stacker Guiding rail for
Reaction Modules

4. Storage

- Unopened at room temperature (15-30°C) in the original box.
- Do not expose to sunlight or humidity.
- After use, reseal the box to protect the modules from contamination.
- Do not use after expiry date (two years after the production date) on the label.

5. Warning and precautions for users

- Observe the normal precautions required for handling all laboratory reagents.
- Disposal of all waste material should be in compliance with local guidelines.

MAGLUMI[®] Starter 1+2

1. INTENDED USE

Starter reagent kit required to generate the chemiluminometric light signal necessary for the implementation of MAGLUMI assays on the MAGLUMI Fully-auto chemiluminescence immunoassay (CLIA) analyzer.

| Catalog Number | Specification | Indicates Analyzers |
|----------------|--|--|
| 130299004M | Starter 1: 230 mL × 1 Starter 2: 230 mL × 1 | Maglumi 600, Maglumi 800, Maglumi 1000, Maglumi 1000 Plus, Maglumi 2000, Maglumi 2000 Plus, Maglumi 4000 and Maglumi 4000 Plus |
| 130299027M | Starter 1: 1.5 L × 1 Starter 2: 1.5 L × 1 | MAGLUMI X8 |

2. MATERIAL PROVIDED

| Contents | |
|--|-------------------------------------|
| Starter 1 | Catalyst in 1.5% NaOH |
| Starter 2 | 0.18% H ₂ O ₂ |
| All reagents are provided ready-to-use in optimized concentrations for the MAGLUMI CLIA analyzers. | |

For additional details such as proper test performance of MAGLUMI assays, please refer to the instructions for the assay to be performed.

3. STORAGE AND STABILITY

- The starter reagents are to be stored unopened in their original vials.
- Store at 2-30 °C until the expiry date on the label (12 months after the production date).
- On board stability: 4 weeks.
- Keep away from direct sunlight.

4. HANDLING

Whenever the starter reagents are replaced, follow the procedure described in the User Manual.

- Place Starter 1 and Starter 2 into the loading position in "1" and "2". To avoid any mix-up, the tubings for the starter reagents are marked. Correct positioning is absolutely essential.
- Do not shake violently to avoid bubble formation.
- Always close starter cover after loading reagents.
- Never leave starter reagent vials open.
- Do not exchange the cover of starter reagent vials 1 and 2.

5. WARNING AND PRECAUTIONS FOR USERS

- For *In Vitro* Diagnostic Use.
- Observe the normal precautions required for handling all laboratory reagents.
- Disposal of all waste material should be in compliance with local guidelines.
- The starter reagents include 1.5% NaOH and 0.18% H₂O₂. If splashes of the NaOH solution or the H₂O₂ solution get into the eye, immediately and thoroughly flush with water or a suitable buffer solution. If necessary, a physician should be consulted.
- Safety data sheets are available on request.

Pursuant to EC Regulation 1272/2008 (CLP) hazardous reagents are classified and labeled as follow:

| | |
|---|--|
| REAGENTS: | Starter 1 |
| CLASSIFICATION: | Skin Irritant, Category 2 |
| SIGNAL WORD: | Warning |
| SYMBOLS / PICTOGRAMS: |  GHS07 Exclamation mark |
| HAZARD STATEMENTS: | H315 Causes skin irritation. |
| PRECAUTIONARY STATEMENTS: | P264: Wash hands thoroughly after handling. P280: Wear protective gloves/protective clothing. P362+P364: Take off contaminated clothing and wash before reuse. P302+P352: IF ON SKIN: Wash with plenty of soap and water. |
| CONTAINS: (only substances prescribed pursuant to Article 18 of EC Regulation 1272/2008). | Sodium hydroxide; caustic soda (Starter 1) |

Shenzhen New Industries Biomedical Engineering Co., Ltd.
No.23, Jinxiu East Road, Pingshan District, 518122 Shenzhen, P.R. China
Tel: +86-755-21536601 Fax: +86-755-28292740

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
Tel: +49-40-2513175 Fax: +49-40-255726

SYMBOLS EXPLANATIONS

| | | | |
|---|---|---|---|
|  | Consult instructions for use |  | Manufacturer |
|  | Temperature limit (Store at 2-30°C) |  | Use-by date |
|  | This way up |  | Keep away from sunlight |
|  | <i>In vitro</i> diagnostic medical device |  | Authorised representative in the European Community |
|  | Catalogue number |  | Kit components |
|  | Batch code | | |

MAGLUMI™ Wash Concentrate

1. INTENDED USE

MAGLUMI Wash Concentrate is intended to be diluted for preparation of System Liquid which is required for washing the magnetic microbeads of the MAGLUMI Fully-auto chemiluminescence immunoassay (CLIA) analyzer (Including Maglumi 600, Maglumi 800, Maglumi 1000, Maglumi 1000 Plus, Maglumi 2000, Maglumi 2000 Plus, Maglumi 4000, Maglumi 4000 Plus and MAGLUMI X8).

2. MATERIAL PROVIDED

MAGLUMI Wash Concentrate (REF: 130299005M) :
1 × 0.714 L

Each bottle contains 14 times concentrated Tris-HCl buffer solution for 750 assays.

For additional details such as proper test performance of MAGLUMI assays, please refer to the instructions for the assay to be performed.

* Assay number may differ because of working condition.

3. PREPARATION

- Add one bottle of MAGLUMI Wash Concentrate into a 10 liter canister.
- Fill the canister with purified water slowly up to the 10 L mark. Use only purified water up to standard ASTM D1193-06(2011) Type IV Grade C.
- Mix the liquid gently.
- **Avoid foam formation!**
- Pouring the remainders or continuous reloading of the ready-to-use System Liquid in the canister is possible, provided that the canister is rinsed once weekly.
- Avoid foam formation during pouring.

4. STORAGE AND SHELF-LIFE

- Storage of Wash Concentrate: At 2-30 °C until the expiry date (12 months after the production date) on the label.
- Shelf-life of ready-to-use System Liquid: 4 weeks at 2-30 °C.
- **Keep away from direct light.**

5. RINSING OF CANISTER

- The canister need to use 1% sodium hypochlorite solution for rinsing (recommended once a week).
- Discard residual liquid before cleaning the canister.
- Add 2L 1% sodium hypochlorite solution into the canister.
- Close the canister.
- Shake the canister to and fro.
- Empty the canister and rinse it more than three times with purified water.
- Use only purified water up to standard ASTM D1193-06(2011) Type IV Grade C.

6. WARNING AND PRECAUTIONS FOR USERS

- FOR PROFESSIONAL USE ONLY
- Observe the normal precautions required for handling all laboratory reagents.
- Handle all washing waste as potentially biohazardous material.
- Disposal of all waste material should be in compliance with local guidelines.
- Avoid contamination of skin and mucosa.
- Safety data sheets are available on request.

Shenzhen New Industries Biomedical Engineering Co., Ltd.

No.23, Jinxu East Road, Pingshan District, 518122 Shenzhen, P.R. China

Tel: +86-755-21536601

Fax: +86-755-28292740

Shanghai International Holding Corp. GmbH (Europe)

Eiffestrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175

Fax: +49-40-255726

SYMBOLS EXPLANATIONS

| | | | |
|--|---|--|---|
| | Consult instructions for use | | Manufacturer |
| | Temperature limit (Store at 2-30 °C) | | Use-by date |
| | This was up | | Keep away from sunlight |
| | <i>In vitro</i> diagnostic medical device | | Authorized representative in the European Community |
| | Catalogue number | | Kit components |
| | Batch code | | |

MAGLUMI[®] System Tubing Cleaning Solution

1. INTENDED USE

The MAGLUMI System Tubing Cleaning Solution is used to clean the pipettor and the washer needles. It is useful to minimize reagent carryover by reducing protein precipitation in the tubing, and improve routine maintenance of the MAGLUMI system. The kit can only be used with MAGLUMI Fully-auto chemiluminescence immunoassay (CLIA) analyzer (Including Maglumi 600, Maglumi 800, Maglumi 1000, Maglumi 1000 Plus, Maglumi 2000, Maglumi 2000 Plus, Maglumi 4000, Maglumi 4000 Plus and MAGLUMI X8), with tubing cleaning program.

2. MATERIAL PROVIDED

System tubing cleaning solution(REF: 130299007M):

500 mL/bottle, containing sodium hypochlorite.

One empty reagent integral with 7 vials:

Only the 3 blank vials near the handle of the kit are for tubing cleaning solution, max volume 25mL each.

3. MAINTENANCE FREQUENCY

MAGLUMI System Tubing Cleaning Solution is used in the monthly maintenance procedure, as the last step in the normal maintenance session.

The solution is sufficient for 12 times of maintenance. Maintenance with MAGLUMI System Tubing Cleaning Solution cannot replace routine maintenance of the MAGLUMI system.

Some MAGLUMI kits may require more frequent maintenance sessions with clean solution; in addition, if the workload is very heavy for the analyzer, the maintenance also can be carried weekly.

4. HOW TO USE

- Insert the empty reagent integral filled with MAGLUMI System Tubing Cleaning Solution integrally in the reagent area track 1.
- Run the tubing cleaning program in the software.
 - For one washing cycle: Transfer the solution from MAGLUMI System Tubing Cleaning Solution bottle into the 3 blank reagent integral vials near the kit's handle. The reagent integral is added with 13mL liquid each time.
 - For two washing cycles: Transfer the solution from MAGLUMI System Tubing Cleaning Solution bottle into the 3 blank reagent integral vials near the kit's handle. The reagent integral is added with 25mL liquid each time.
- The whole procedure needs about 40 minutes.

5. STORAGE

- Upon receipt, store the kit at 2-30 °C.
- Keep away from the sunlight. The shelf life of solution is 12 months.

6. WARNING AND PRECAUTIONS FOR USERS

- For *In Vitro* Diagnostic Use.
- Before using the MAGLUMI System Tubing Cleaning Solution, please read the operation instructions carefully.
- Do not shaking violently to avoid bubble formation.
- Handle the waste liquid as potentially hazardous biological material. They must be considered potentially able to transmit infectious agents and should therefore be disposed of in accordance with the prevailing regulations and guidelines of the agencies holding jurisdiction competent authorities.
- If splash contacts your eyes, rinse your eyes immediately; in case of swallowing by accident, drink plenty of water, and seek the medical care immediately.

Pursuant to EC Regulation 1272/2008 (CLP) hazardous reagents are classified and labeled as follow:

| | |
|-----------------------|---|
| REAGENTS: | System tubing cleaning solution |
| CLASSIFICATION: | Skin Irritant, Category 2 |
| SIGNAL WORD: | Warning |
| SYMBOLS / PICTOGRAMS: |  GHS07 Exclamation mark |
| HAZARD STATEMENTS: | H315 Causes skin irritation. |

| | |
|--|---|
| PRECAUTIONARY STATEMENTS: | <p>P264: Wash hands thoroughly after handling.</p> <p>P280: Wear protective gloves/protective clothing/eye protection/face protection.</p> <p>P302+P352: IF ON SKIN: Wash with plenty of soap and water.</p> <p>P362: Take off contaminated clothing.</p> |
| CONTAINS:(only substances prescribed pursuant to Article 18 of EC Regulation 1272/2008). | Sodium hypochlorite, solution 2.0%-4.7% Cl active. |

Shenzhen New Industries Biomedical Engineering Co., Ltd.

No.23, Jinxu East Road, Pingshan District, 518122 Shenzhen, P.R. China
 Tel: +86-755-21536601 Fax: +86-755-28292740

Shanghai International Holding Corp. GmbH (Europe)

Eiffestrasse 80, 20537 Hamburg, Germany
 Tel: +49-40-2513175 Fax: +49-40-255726

SYMBOLS EXPLANATIONS

| | | | |
|---|---|---|---|
|  | Consult instructions for use |  | Manufacturer |
|  | Temperature limit (Store at 2-30 °C) |  | Use-by date |
|  | This way up |  | Keep away from sunlight |
|  | <i>In vitro</i> diagnostic medical device |  | Authorised representative in the European Community |
|  | Catalogue number |  | Kit components |
|  | Batch code | | |

MAGLUMI[®] Light Check

Intended Use

Light Check reagent required to check the validity of Starter 1+2 as well as the functioning of the measuring and pipetting units. A Light Check has to be carried out on the MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer (Including Maglumi 600, Maglumi 800, Maglumi 1000, Maglumi 1000 Plus, Maglumi 2000, Maglumi 2000 Plus, Maglumi 4000, Maglumi 4000 Plus and MAGLUMI X8) either once daily prior to starting the first measurement or whenever a new lot of starter reagents is used. This control program checks the proper functioning of the instrument and starter reagents to avoid data error due to instrument defects or mislocated or expired starter reagents.

Material provided

| | | |
|--|---|------------|
| Package (REF: 130299006M) | 5 vials per box. Each vial contains 2 mL Light Check liquid. | |
| Component | ABEI (N-(4-Aminobutyl)-N-ethylisoluminol) | 0.77 mg/mL |
| | BSA | 10 mg/mL |
| The Light Check reagent is provided in liquid. | | |

For additional details such as proper test performance of MAGLUMI series assays, please refer to the instructions for the assay to be performed.

Preparation of Reagents

- Carefully open the vial containing the liquid material.
- Prior to use, turn upside down slowly and let it stand for 5 minutes.
- Foam **must** be avoided and cleaned before use.

Storage and Stability

- Unopened in the original vial at 2-8°C until the expiry date (12 months after the production date) on the label.
- **After use, immediately close the Light Check vial and store it in a refrigerator at 2-8°C.**
- Keep away from direct sunlight.
- Opened Light Check reagent can be steady for 28 days when stored at 2-8°C.

Handling

- Place the Light Check Vial into the matched patient rack; make the barcode label of the vial face the barcode reader of the sample area (refer to operation instructions of MAGLUMI series Fully-auto chemiluminescence immunoassay analyzer for detailed information).
- Seal the vial and store the rest of the Light Check at 2-8°C immediately after been used.

Target values for Light Check *

- RLU: 450,000-650,000.
- (* Target values may differ due to analyzer setting; contact SNIBE authorized technician for further information)

Warning and precautions for users

- FOR PROFESSIONAL USE ONLY.
- Observe the normal precautions required for handling all laboratory reagents.
- Disposal of all waste material should be in accordance with local guidelines.

Shenzhen New Industries Biomedical Engineering Co., Ltd.
No.23, Jinxiu East Road, Pingshan District, 518122 Shenzhen, P.R. China
Tel: +86-755-21536601 Fax: +86-755-28292740

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany
Tel: +49-40-2513175 Fax: +49-40-255726

SYMBOLS EXPLANATIONS

| | | | |
|--|---|--|---|
| | Consult instructions for use | | Manufacturer |
| | Temperature limit (Store at 2-8°C) | | Use-by date |
| | This way up | | Keep away from sunlight |
| | <i>In vitro</i> diagnostic medical device | | Authorised representative in the European Community |
| | Catalogue number | | Kit components |
| | Batch code | | |

Reagent Note

International
Technical
Service
Department

| | | | |
|----------------|--|--------------------|------------------------------|
| Assay | MAGLUMI® SARS-CoV-2 S-RBD IgG (CLIA) MAGLUMI® SARS-CoV-2 Neutralizing Antibody (CLIA) | | |
| Subject | WHO Standard (20/136) Unit Conversion | | |
| REF No. | RN21040201 | Issued Date | April 2 nd , 2021 |

| | |
|----------------------------------|---|
| Immediate action required | |
| Recommendation | |
| For your information | √ |

Dear customers

We have calibrated the SARS-CoV-2 S-RBD and SARS-CoV-2 Neutralizing antibody kits with the WHO standard (20/136).

By analyzing and calculating the test results of WHO standard (20/136), the following unit conversion relationships are obtained.

1. For SARS-CoV-2 S-RBD IgG reagent, 1 AU/mL is equivalent to 4.33 BAU/mL.

2. For SARS-CoV-2 Neutralizing Antibody reagent, 1 µg/mL is equivalent to 405 IU/mL. (It is based on the conversion of 0.3 µg/mL equivalent to 121.5 IU/mL with the SARS-CoV-2 Neutralizing Antibody assay. by testing 90 individuals confirmed SARS-CoV-2 VNT₅₀ ≥ 20, all values were ≥ 0.300 µg/mL.)

Snibe International Technical Service Department

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