## Naudojimo instrukcija

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8-

## Declaration of Conformity

### **gke Steri Record®** Chemical Process Indicators type 1 according to EN ISO 11140-1 to monitor sterilization processes

Self-adhesive documentation labels with indicator to be used in hand labeller:

Art. No.	Description	Rolls/pack	Adhesiveness	Sterilization process
230-871 - 874	C-S-L-1-SA-G,-B,-R,-Y	12	single	Steam ✓
240-871 - 874	C-S-L-1-DA-G,-B,-R,-Y		double	
240-883	C-S-L-1-DA-R	4	double	Ethylene oxide Formaldehyde Hydrogen peroxide
232-875	C-E-L-1-SA-P	12		
242-875	C-E-L-1-DA-P	12		
242-885		2		
243-874	C-F-L-1-DA-Y	12		
243-884		2		
244-871	C-V-L-1-DA-G	12		
244-881		2		
244-873	C-V-L-1-DA-R	12		
244-883		2		

Double self-adhesive process indicators for label printers:

Art. No.	Description	Amount/roll	Application	Sterilization process
211-377, -366	C-S-L-1-38x39-SC-DA	700	To print single labels	Steam
211-342, -442, -282, -482	C-S-L-1-60x40-SC-DA	800		
211-344, -345	C-S-L-1-60x40-SU-DA	800		
211-378	C-S-L-1-60x72-SM-DA	400		
211-381		450		
211-349, -449	C-S-L-1-80x40-SS-DA	800		
211-373	C-S-L-1-80x50-W-SU-DA	750		
211-382	C-S-L-1-80x50-WF-SA	1.000		
211-367, -369, -370, -469	C-S-L-1-SL-DDOC	800		
211-375		2.000		
211-390, -394, -490, -494	C-S-L-1-70-DA	70/80 m	Endless roll to be cut in the printer	

Process indicators for containers to be used in steam sterilization processes:

Art. No.	Description	Amount	Application
211-341	C-S-L-1-CL-40x60-SA	1.000/roll	Container label
211-362	C-S-L-1-CS	1.008/pack	Container seal

Notice: On all **gke** packages, an additional letter code has been added to the 6-digit article number. The additional letter code refers to the language and/or customized version. It is only added on the outside label, the inside of the pack is identical to the article numbers and the above table. All articles with the same 6-digit number have the same specifications.

**They conform with the European and International Standard**

## EN ISO 11140-1 Type 1 ✓

We hereby declare under our sole responsibility, that all above mentioned process indicators have been tested in the **gke** application laboratory using a resistometer according to EN ISO 18472 (certificate available on request) and are in conformity with the European and International Standard EN ISO 11140-1 type 1. The chemical indicator strips have an expiry date of 3 years after production. The expiry date is printed on the package.

The indicators and the adhesive do not contain Latex, DEHP or PVC and the adhesive does not damage the packaging materials when used under normal process conditions. All above products manufactured by **gke**-GmbH are produced in accordance with the requirements of the Medical Device Regulation (MDR) 2017/745.

This document certifies that the above performance criteria and the **gke** test requirements for quality control are met. The continuous quality is guaranteed by our quality management system according to EN ISO 13485\*.

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# Declaration of Conformity

## **gke Steri Record®** Indicators to monitor sterilization processes

are available from **gke** as the following product variations according to EN ISO 11140-1:

Art. no.	Product code	Pack size	Sterilization process	Indicator type according EN ISO 11140-1			
<b>BDS introduction packs</b>							
211-110	C-S-BDS-EU-HPCD-KIT	1 + 100	Steam	Type 2			
211-150	C-S-BDS-EUH-RCPCD-KIT						
211-120	C-S-BDS-EU-RCPCD-KIT						
211-130	C-S-BDS-USA-RCPCD-KIT						
<b>Indicator refill pack for all gke BDS-Tests</b>							
211-111	C-S-BDS-SV1	100					
211-112		250					
211-115		500					
<b>Batch monitoring introduction packs</b>							
211-260	C-S-PM-HL-HPCD-KIT	1 + 100			EO	n.a.	
211-263	C-S-PM-HL-RCPCD-KIT						
211-264	C-S-PM-HL-OCPCD-KIT						
211-253	C-S-PM-SHL-RCPCD-KIT						
211-254	C-S-PM-SHL-OCPCD-KIT						
<b>Indicator refill packs for all gke batch monitoring systems</b>							
211-252	C-S-PM-SV1	250	LTFSF	n.a.			
211-255		500					
212-202	C-E-PM	250					H <sub>2</sub> O <sub>2</sub>
213-203	C-F-PM						
214-203	C-V-PM						
<b>Package monitoring multi-variable indicators</b>							
211-403	C-S-P-4-SV1	3200			Steam	Type 4	
212-360	C-E-P-4-SV3				EO		
214-221	C-V-P	400			H <sub>2</sub> O <sub>2</sub>	n.a.	
214-222		800					
214-223		3200					
<b>Package monitoring integrating indicators</b>							
211-224	C-S-P-5-SV1	400	Steam	Type 5			
211-225		800					
211-226		3200					
211-230	C-S-P-5-78x48-SA-SV1	1000					
211-235		500					
211-220	C-S-P-5-58x35-SV1	1000			EO	n.a.	
211-222		200					
212-206	C-E-P-5-58x35-W-SV3	200			LTFSF		
213-206	C-F-P-5-58x35-W-SV4						
214-206	C-V-P-58x35-W						
<b>Package monitoring emulating indicators</b>							
211-241	C-S-P-6-SV1	2000	Steam	Type 6			
211-242		500					
211-243		250					
211-238	C-S-P-6-SV2	2000					
211-239		500					
211-240		250					

**Notice:** On all **gke**-packages, an additional letter code has been added to the 6-digit article number. The additional letter code refers to the language and/or customized version. It is only added on the outside label, the inside of the pack is identical to the article numbers in the above table. All articles with the same 6-digit number have the same specifications.

We hereby declare under our sole responsibility that above mentioned integrating indicators have been tested in the **gke** application laboratory and are in conformity with the European and ISO standard EN ISO 11140-1 for steam, ethylene oxide, formaldehyde and hydrogen peroxide sterilization processes. The corresponding indicator types are available in the above table. The indicators listed above are tested in a resistometer according EN ISO 18472.

According to ISO 11140-1, 5.9 the indicators do not release any particles or toxic substances in sufficient quantities to cause a health hazard. The colour of the indicators remains constant after the sterilization process has passed successfully and does not fade back to the colour before sterilization for at least 5 years under the condition the indicator is not stored in contact with oxidation agents like H<sub>2</sub>O<sub>2</sub> or other chemicals.

The above mentioned products manufactured by **gke**-GmbH are produced in accordance with the requirements of the Medical Device Directive 93/42/EEC and the Medical Device Regulation 2002:618.

This document certifies that the above performance criteria and the **gke** requirements for quality control are met. The continuous quality is guaranteed by our quality management system according to EN ISO 13485\*.

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751-019 EN V15

## Declaration of Conformity

### gke Steri-Record® Batch Monitoring Systems (BMS) including Refill Packs for Steam Sterilization Processes

Art.-No. *	gke Product Name	Content		Product Description of the Process Challenge Device (PCD)	Indicator according EN ISO 11140-1 type 2 using the following reference for equivalence	
		PCD	Indicator strips			
211-253	C-S-PM-SHL-RCPCD-KIT	1	100	Compact PCD (C-PCD) round version	7 kg cotton pack with biological indicators according EN ISO 11138-3	
200-020	PM-SHL-RCPCD	1	-			
211-254	C-S-PM-SHL-OCPCD-KIT	1	100	Compact PCD (C- PCD) oval version		
200-024	PM-SHL-OCPCD	1	-			
211-260	C-S-PM-HL-HPCD-KIT	1	100	Helix PCD	PCD according EN 867-5**	
200-150	PM-HPCD-2-150	1	-			
211-263	C-S-PM-HL-RCPCD-KIT	1	100	Compact CD, (C- PCD) round version		
200-021	PM-HL-RCPCD	1	-			
211-264	C-S-PM-HL-OCPCD-KIT	1	100	Compact PCD (C- PCD) oval version		
200-026	PM-HL-OCPCD	1	-			
200-029	PM-HDH-RCPCD	1	-	Compact-PCD, round version, red		higher requirements than EN 867-5**
200-030	PM-VHDH-RCPCD	1	-	Compact-PCD, round version, brown		much higher requirements than EN 867-5**
211-252	C-S-PM-SV1	-	250	Indicator refill packs for all Batch Monitoring Systems (BMS) above SV = 134°C; 3 min and 121°C; 15 min	Indicator strips for above mentioned test devices	
211-255		-	500			
211-211	C-S-PM-SV2	-	100			
211-212		-	250			
211-215		-	500			

\* Notice: On all gke-packages, an additional letter code has been added to the 6-digit article number. The additional letter code refers to the language and/or customized version. It is only added on the outside label, the inside of the pack is identical to the article numbers in the above table. All articles with the same 6-digit number have the same specifications.  
\*\* Test reports are available on request.

All gke Steri-Record®-batch monitoring systems (BMS) are indicator systems according to EN ISO 11140-1 type 2, consisting of a process challenge device (PCD) with a capsule hosting a chemical indicator inside and tested in a resistometer according EN ISO 18472.

Two PCD designs are available:

1. The Compact-PCD®s consist of an external plastic case with an internal stainless-steel tube, available in round and oval design. Both versions meet the same specification. gke guarantees a life span of more than 10.000 cycles under the condition the PCD is used according to the directions for use.
2. The Helix PCDs consist of a 1,5 m, 2 mm inner diameter PTFE-tube with a stainless-steel capsule at one end to host the indicator.

The combination of a PCD and indicator (= indicator system type 2) is used to monitor the temperature and time integral as well as the efficacy of the penetration characteristics of steam in a sterilization process.

During the sterilization process the main physical parameters, pressure and temperature, can be monitored and are usually recorded by the sterilizers. The gke Steri-Record®-batch monitoring systems in addition monitor air removal, potential leaks and the presence of non-condensable gases in steam to assure the total penetration of steam into packs and into hollow devices and therefore sterility at the worst cases inside the process. The PCDs have been calibrated, using the minimum sterilization process conditions of 134 °C, 3 minutes or 121°C, 15 minutes are achieved, all bars of the chemical indicator in the PCD change colour from yellow to black. If some bars of the chemical indicator remain yellow or if the color changes only to yellow-brown after a longer sterilization period, non-condensable gases are present inside of the PCD with the consequence of a potential malfunction in part of the process. The sterilization times may be increased up to 10 minutes at 134 °C or 30 minutes at 121 °C.

For 134°C, 18 min (prion program) we recommend to use the specific indicator strip (refill pack art. no. 211-211; -212; -215, see above).

The PCDs can be used in small and large sterilizers for solid instruments, porous loads and complex minimal invasive instrumentation. All BMS can be used for routine batch monitoring in steam sterilization processes. During process validation the information should be provided which BMS should be used for routine monitoring.

Above gke Steri-Record® batch monitoring systems can't be used for sterilization processes which achieve the air removal with single vacuum, gravity displacement, steam flashing or overpressure air removal cycles. Depending on the load configuration and packaging an individual validation for those processes is necessary to select a PCD accordingly for batch monitoring purposes.

The test devices and the indicators do not release any particles or toxic substances. The colour of the indicators remains constant after sterilization and does not fade back to the colour before sterilization.

The test results are only valid if original gke Steri-Record® indicator strips with original gke Steri-Record® test devices (PCD's) are used and the instructions for use are applied.

This document certifies that the above performance criteria and the gke test requirements for quality assurance are met. The continuous quality is of our products is assured by our quality management system according to EN ISO 13485\*.

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751-006 EN V08

# Colour reference chart

## for all *gke* Batch Monitoring Systems (BMS) (SV = 134°C; 3 min or 121°C; 15 min) for steam sterilization processes

### PASS RESULTS

The batch can be released if all 4 colour segments of the processed indicator are black as shown below, and if all other critical process parameters such as temperature, pressure and sterilization time are in accordance with the cycle reference.


all colour segments are black

→ sufficient steam penetration

### FAIL RESULTS


The examples shown below demonstrate insufficient steam penetration which may be due to various reasons (see note below).

If one of the results shown below is found after a sterilization process, the respective batch must not be released and requires re-processing.


colour segments are black at one end and yellow-brownish at the other end

→ insufficient steam penetration and air removal

all colour segments are brownish, but not black

→ temperature achieved, but no steam penetration and air removal

all colour segments are yellow

→ insufficient temperature and no steam penetration and air removal

#### Note

A failure of the batch monitoring test is not a conclusive proof that the fault is necessarily associated to the sterilizer itself (insufficient air removal or leakages), but may well be due to external causes such as unpurged steam pipes at the start of the process (if no prior Bowie-Dick Simulation Test was run), non-condensable gases in the external steam supply or changing temperature of the cooling water. Therefore batch monitoring tests made in subsequent cycles may show different results.

720-002 EDM V02.05/2020

+EGKE720002EDM0/SS129925

19-

---

# ***gke Steri-Record*<sup>®</sup> Indicator System according to EN ISO 11140-1 Type 2 to monitor sterilization processes**

---


Fig. 1: Different sensitive process challenge devices, round section

## **Application**

In the sterilizer standards EN 285 and EN 13060 type tests are described, e.g. BD-Test or Helix Test according to EN 867-5 (new EN ISO 11140-6). They are only testing the standard requirements of the sterilizers without any relation to the load.

They are quite often used as so-called batch monitoring systems and sold as various test systems on the market to monitor the load. However, these tests have no relation to the load. Therefore, it cannot be automatically concluded that difficult loads can be sterilized successfully. A monitoring system for the load should not just pass a sterilizer but must represent the load.

This adoption can be achieved in two ways:

### **1. Process Monitoring System (PMS):**

The sterility of the load is assured if the sterilization process has been validated with a worst-case load. The continuous monitoring of the load shall be carried out with a test system which monitors the maximum performance of the sterilizer. Modern sterilizers provide much better performance than the required type tests mentioned above. For this purpose **gke** offers a variety of different PCDs (process challenge devices) with increasing sensitivities to monitor the process at highest level. These test devices are called process monitoring systems (PMS).

### **2. Batch Monitoring System (BMS):**

If the characteristics of the load to be sterilized are known, a PCD can be designed and validated to be more difficult to sterilize than the defined load using the method according to DIN 58921 (medical device simulator). **gke** has already tested several typical loads with the method according to DIN 58921 and offers special PCDs for Tattoo, Dental and Ophthalmic loads. For these test systems separate data sheets with detailed information are available.

## **Product Description**

### **1. Process Challenge Device:**

The **gke Steri-Record**<sup>®</sup> Compact-PCD<sup>®</sup> consists of an external plastic case with an internal stainless steel tube and a capsule holding the indicator.

The four colour-coded PCDs are available in two different construction versions:

- round version with stainless steel bracket. This PCD can be put vertically on a loading rack or placed horizontally on the bracket.
- oval version (green and orange): This PCD can be placed horizontally on the flat side of the PCD, suitable for small sterilizers (2.5 cm height).

All PCDs are re-usable and can be used for a considerable number of sterilization cycles.

## 2. Indicator Strip:

The indicators have a size of 40 x 6 mm and fit into the PTFE holder of the capsule. This size is identical to standard biological indicators according to EN ISO 11138 series and fit also in those PCDs and can be used alternatively, but shall not be opened after sterilization and aseptically transferred in a microbiological laboratory. Chemical indicators for all sterilization processes can be ordered from **gke**. For details see our product catalogue.

## Performance Characteristics

The **gke** PCDs are type 2 indicator systems according to EN ISO 11140-1, consisting of a „specific test load“ (Process Challenge Device = PCD) and „indicator“ (indicator strip). They have been validated from a laboratory accredited according to the standard EN ISO 17025 (test reports on request).

**gke** offers a variety of PCDs with different levels of difficulty. Using a test with the highest level which will still pass in a validated process the sterilizer is always tested at its limits. Failures that are not detected with standard methods will now be visible and the maximum possible safety to test sterilization processes can be guaranteed. If a process has been validated with a defined load and a suitable test system, this method secures the sterilization at the most difficult locations. This is not possible by only recording the physical data with documentation, e.g. with a printer.

The identification of a suitable PMS to monitor sterilization has to be verified during process validation depending on the load. The **gke** representatives support in the selection of the test system depending on the load.


Fig. 2: PCDs, round and oval construction

## 1. PMS Compact-PCD<sup>®</sup>, colour: green

This PMS Compact-PCD<sup>®</sup>, colour: green, represents a challenging test concerning air removal and steam penetration that can be used for routine monitoring of simple hollow devices, tubes, porous loads and solid instruments. However, this PMS is less demanding regarding air removal and steam penetration than the hollow load test described in EN 867-5.

## 2. PMS Compact-PCD<sup>®</sup>, colour: orange

The PMS Compact-PCD<sup>®</sup>, colour: orange, simulates the hollow load test according to EN 867-5. This type test is included in the standard for large sterilizers (EN 285) as well as in the standard for small class B sterilizers (EN 13060).

The standard EN 867-5 describes a test cycle with 9 injections of steam to 950 mbar. In this test cycle the hollow load test will be successfully penetrated, if a vacuum of 450 mbar will be reached in between each steam pulse.

## 3. HDH-PMS Compact-PCD<sup>®</sup>, colour: red

The HDH-PMS Compact-PCD<sup>®</sup>, colour: red, (High Demand Hollow) is used for routine monitoring of complex instruments or long tubes, whereas the hollow load test according to EN 867-5 is not anymore sufficient.

To successfully penetrate the **gke**-HDH-PMS in the same test cycle a vacuum of at least 210 mbar is necessary.

## 4. VHDH-PMS Compact-PCD<sup>®</sup>, colour: brown

The VHDH-PMS Compact-PCD<sup>®</sup>, colour: brown, (Very High Demand Hollow) makes increasing demands on the air removal performance. The test cycle according to EN 867-5 with 9 injections of steam requires vacuum cycles of at least 160 mbar to successfully penetrate the VHDH-PMS.

These data are only valid under laboratory conditions using a reproducible process in a test sterilizer according to EN ISO 18472. All other variables that affect the air removal and steam penetration are strictly defined in such a program and must remain unchanged in all test cycles.


Fig. 3: VHDH-PMS Compact-PCD<sup>®</sup>, colour brown, round

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