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
Stereodylator®

and accessories

STEREODYNATOR® with accessories and optional device cart

Treatment time	max. 60 minutes
Display	15" TFT LCD Touchscreen
Power supply	110-240V, 50/60Hz
Power consumption	max 200 VA
Class of protection	I
Class of equipment	Ila according to MDD
Dimensions	350 x 290 x 350 mm (W x H x D)
Weight	max. 13 kg without accessories
Optional suction mode	Venturi principle for up to 2 channels and 4 suction cups
Optional ultrasound mode	1/3 MHz
Warranty period	24 months
	Subject to technical changes without prior notice.

The following versions of STEREODYNATOR® are available:

article	Part N°
STEREODYNATOR®	025-0-1000
STEREODYNATOR® with module for Ultrasound Therapy, incl. 1 Ultrasound probe 5 cm ² , incl. 1 Contact gel 250 ml	025-0-1000-U
STEREODYNATOR® incl. Suction application module	025-0-1000-V
STEREODYNATOR® with module for Ultrasound Therapy, incl. 1 Ultrasound probe 5 cm ² , incl. 1 contact gel 250 ml and incl. Suction application module	025-0-1000-UV
Device cart for STEREODYNATOR®	026-0-3000 + 026-0-3020

STEREODYNATOR®

This comfort therapy unit is the new benchmark in Electro Therapy. The unit works with a three-dimensional interference current generated through the superposition of three middle-frequency currents. It can be used as an Electro Therapy device, ultrasound therapeutic device and for the combination of electrical stimulus current with ultrasound.

The combination therapy combines the analgesic and hyperemic effect of the electrical stimulus and the mechanical and thermic effects of ultrasound.

Characteristics:

- three-dimensional stimulation effect
- multi-site stimulation effect intensity dynamics
- dynamic behaviour of stimulation site coupled with endogenous/exogenous stimulation shifting

Application:

Practical experiences made clear that the combination of both physiotherapeutic effects are particularly suitable for diagnostics and the treatment of pain in the musculoskeletal system and the peripheral nerves. Endogenous heat, micro-massage and ionic shift bring about the extraordinary therapeutic successes, even after a short period of treatment. Also with circulatory disturbances and rheumatic-infl ammatory processes this therapy is a striking addition to the HF-Heat Therapy or to drug therapy.

Advantageous is the little fatigue of the muscle due to bipolar currents and the influence of the vegetative systems with diadynamic currents.

Application in treatments:

- muscle toning and detoning
- galvanization and iontophoresis
- treatment of paralyses with complete or partial muscular degeneration
- treatment of atrophy due to lack of movement and weak muscles
- pain, muscle spasm, functional diseases of the musculoskeletal system
- pelvic floor stimulation
- treatment of urinary and fecal incontinence

Application in diagnostics:

- MF-test according to Dr. Lange
- recording of I/T-curves (shown numerically and graphically on the display)
- extensive non-invasive electro diagnostics of peripheral paralyses
- qualitative and quantitative determination of faradic excitability
- determination of rheobase, chronaxy and accommodation neurodiagnostic examinations with galvano-palpation

Optional:

- integrated module for Ultrasound Therapy (1/3 MHz) with titanium head as well as for the combination of Electro- and Ultrasound Therapy (optional)
- integrated module for Venturi suction application
- device cart – in design and function matching with the new Stereodylator®

ONE TOUCH TO START

Der STEREODYNATOR®- the new sophisticated Electro Therapy device with „that certain extra“: special quick start of therapy via list of favorites or patient's file
See for more information: www.stereodylator.gbo-med.de


Electrotherapy with stereodynamic current

The use of stereodynamic interference is an important and valuable extension to conventional stimulation-current therapy. A stereodynamic interferential current results from a superposition of three medium frequency currents, each flowing in a different direction. Selecting suitable stimulation frequencies leads to a three-dimensional (stereodynamic) stimulation field with endogenous dynamic behavior over time.

Site of interference

Multi-positional stimulation effects

-  Interferential current with equal X and Z parts, minimal Y contribution
-  Interferential current with equal X, Y and Z parts

Effects

When stereodynamic interferential current is applied, there is essentially no sensation of current flow. Even sensitive patients can be treated with higher stimulation current intensities. As a result, the three-dimensional effect can reach even deep-lying tissue regions.

The properties of stereodynamic current are complex: reactive hyperemia is only one of the well-known effects. Through the summation of a large number of pulses, tensed musculature is detonicized, with stimulations subduing of pain – even in the region of refelectory changes.

The motorically effective stimulation frequencies cause a motion „through“ the muscles, while avoiding fatigue of individual muscle fibers. Active kinethotherapy following muscular over-exercition is also promoted.

Traumatic edemas are rapidly absorbed, owing to an improved permeability of the cellular membranes. A slight massaging effect is observed.

Allevation of pain, better blood circulation with an increased supply of oxygen, and accelerated elimination of metabolic waste products are among the principal benefits of this method. Tensed muscles are also detonicized by means of an effective kineto-therapeutical procedure with a minimum of sensory strain of the patient.

The range of indications includes disorders of the skeletal, supportive and locomotory system, circulatory disorders, neuralgia and neuritis.

- Improvement of the metabolic process
- Effect of restoring function to functionally damaged internal organs
- Reduction in swelling with effusion and edemas
- Improved resorption in the case of accidental injuries
- Vegetative and physisic stimulation

Table of currents

STEREODYNATOR® offers you a wide range of different currents for electro therapy. For details see at the user's manual or look at the online menu on the touchscreen of the device.

The following currents are available within the therapy programs of STEREODYNATOR®:

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Interference currents (2 channel):

- Interference
- Dipol Vector
- Isoplanar Vector

Medium frequency currents:

- Russian stimulation (acc. to Kots)
- Bodyshaping
- Cellulitis
- Fatburning
- Lymphdrainage
- Mobilisation
- Trauma
- Sinus modulated, Electro massage
- Gym burning mode
- Gym training
- Gym current

Other forms of currents:

- Galvanic current

Impuls currents:

- High voltage current
- Impulsed galvanisation IG 30
- Impulsed galvanisation IG 50
- Faradic current F1
- FM Frequency modulation 7-14 Hz
- Micro current biphased
- Micro current monophased
- Ultra stimulation current acc. to Traeubert
- TENS biphased
- TENS monophased
- TENS Burst biphased
- TENS Burst monophased

Stereodynamic Interference currents: (3 channel)

- Electro massage
- Lipolysis
- Cellulite
- Muscle pump
- Sedat
- Myomol
- Edema
- Reha
- Pain
- Shaking massage
- Sedation-Modulation
- Vegetat. Stim I-III
- Universal

Stereodynamic currents: (3 channel)




- Gym current
- Gym burning mode
- Gym training

Diagnostic programs:

- Monophased(T/R) single pulse
- Biphased (T/R) single pulse
- I/T-Kurve
- MF-Check acc. to Dr. Lange



Accessories

	Part-No.	Specification	Use	Application	Remarks
	025-0-0100	Electrode handpiece with finger-tip switch, 2 m cable	For neuro diagnosis		
	025-0-0101	Metal electrode for Electrode handpiece	For handpiece with finger-tip switch		
	025-0-0102	Sponge for Electrode handpiece	For handpiece with finger-tip switch		
	45-38-880 EH725	Conductive rubber electrode with cable, 3 x 4 cm², 4 mm connector	With sponge pockets 45-39-102 EH725	With rubber- or elastic-strap	
	45-38-898 EH725	Conductive rubber electrode with Cable 6 x 8 cm², 4 mm connector	With sponge pockets 45-39-110 EH725	With rubber- or elastic-strap	
	45-38-492 EH725	Conductive rubber electrode with cable, 8 x 12 cm², 4 mm connector	With sponge pockets 45-38-203 EH725	With rubber- or elastic-strap	
	45-38-229 EH725	Tin sheet metal 18 x 28 cm², (3 pcs. DIN A 4) with cutting template	For making tin-plate electrodes	With electrode connector and rubber strap	

Accessories

	Part-No.	Specification	Use	Application	Remarks
	45-38-237 EH725	Electrode connector	For plate electrodes and tin sheet metal	With rubber- or elastic-strap	Banana jack 4 mm
	45-02-078 EH023	Connecting terminal	For: plate electrodes and tin sheet metal		Banana jack 4 mm
	45-38-476 EH725	Intermediate layer of viscose sponge, 30 x 20 x 1 cm² (5 pcs.)	For: plate electrodes and tin sheet metal	With rubber- or elastic-strap	Moisten well
	45-38-468 EH725	Tin plate electrode, 9 cm x 12 cm with fixed cable and 4 mm jack	With viscose- or fleece-intermediate layer, fleece- or viscose sponge pockets	With rubber- or elastic-strap	
	45-38-203 EH725	Sponge pockets large, 10 pcs. For electrodes 8 x 12 cm	For plate electrodes or conductive rubber - electrodes	With rubber- or elastic-strap	Moisten well; can also be used for iontophoresis
	45-39-110 EH725	Sponge pockets medium, 10 pcs. For electrodes 6 x 8 cm	For conductive rubber - electrodes	With rubber- or elastic-strap	Moisten well; can also be used for iontophoresis
	45-39-102 EH725	Sponge pockets small, 10 pcs. For electrodes 4 x 6 cm	For conductive rubber - electrodes	With rubber- or elastic-strap	Moisten well; can also be used for iontophoresis

Accessories

Part-No.	Specification	Use	Application	Remarks
22-73-821 Q1408	Viscose sponge pockets, 10 pcs.	For plate electrodes or conductive rubber electrodes (large)	With rubber- or elastic-strap	Molsten well; can also be used for iontophoresis
011-0-0033	Elastic strap with velcro-fastener, 70 cm long, 7,5 cm wide 2 pcs.	For plate-, rubber-, DUO- and Star-electrodes		
45-38-971 EH725	Rubber strap 1 pcs., 1 row punched, 42 cm long, 2,7 cm wide	For plate-, rubber-, DUO- and Star-electrodes		
45-38-963 EH725	Rubber strap 1 pcs., 1 row punched, 135 cm long, 2,7 cm wide	For plate-, rubber-, DUO- and Star-electrodes		
45-38-955 EH725	Rubber strap 1 pcs., 5 row punched, 135 cm long, 7 cm wide	For plate-, rubber-, DUO- and Star-electrodes		
22-74-801 Q1320	Fastening button 1 pcs.	For all rubber straps		
011-0-0024	Bow electrode with 2 x electrode cup 4,2 cm 2 x electrode cup 1,6 cm 2 x sponge 6 cm 2 x sponge 2,2 cm			

Accessories

Part-No.	Specification	Use	Application	Remarks
22-73-605 Q1409	Electrode sponge (Viscose), diameter 6 cm 10 pcs.	For electrode cups large, inner diameter 4,2 cm		Bow electrode only
22-73-613 Q1409	Electrode sponge (Viscose), diameter 2,2 cm 10 pcs.	For electrode cups small, inner diameter 1,6 cm		Bow electrode only
45-38-435 EH725	Patient cable, 2 pole, length 2,5 m	All stimulators of the new series		
45-38-500 EH722	Suction electrode set, single pole, standard, diameter 50 mm	For new series with vacuum option		
45-38-633 EH722	Electrode hose blue	For new series with vacuum option		
45-38-641 EH722	Electrode hose red	For new series with vacuum option		
020-453826	Suction electrode small, diameter 35 mm		Self adhesive	Venturi-electrode, injector system, no moisture sucked in

Accessories

Part-No.	Specification	Use	Application	Remarks
45-38-526 EH722	Suction electrode large, diameter 70 mm, set consisting of: Suction head 1 pole, suction cup, electrode plate, 12 intermediate layers.	For new series with vacuum option	Self adhesive	Venturi-electrode, injector system, no moisture sucked in
011-1-0147	Suction electrode diameter 50 mm, set consisting of Suction head 1 pole, suction cup, electrode plate, 12 intermediate layers	For new series with vacuum option		Can be used as replacement for Series 7XX
011-0-0015	Suction electrode set, consists of 2 x 011-1-0147	For new series with vacuum option		Replacement for 45-00-039EH008
45-38-517 EH722	Electrode plate standard, diameter 40 mm	For suction electrode 50 mm		
45-38-525 EH722	Electrode plate large, diameter 60 mm	For suction electrode 70 mm		
45-38-575 EH722	Suction cup standard, diameter 50 mm	For suction electrode standard 50 mm		
45-38-583 EH722	Suction cup large, diameter 70 mm	For suction electrode 70 mm		

Accessories

Part-No.	Specification	Use	Application	Remarks
45-38-559 EH722	Suction head 1-pole	For all venturi suction electrodes		
020-453827	Electrode inlay small, diameter 25 mm, 30 x 4 mm	For suction electrodes small	Insert in suction electrode	Moisten well, replace frequently
45-38-930 EH722	Electrode inlay standard, diameter 52 mm, 12 pcs.	For Suction electrodes standard	Insert in suction electrode	Moisten well, replace frequently
45-38-948 EH722	Electrode inlay large, diameter 72 mm, 12 pcs.	For Suction electrode large	Insert in suction electrode	Moisten well, replace frequently
011-0-0025	Y-cable	For connecting two electrodes to one pole		
017-0-0077	Vaginal electrode With 2 mm connector			
017-0-0076	Anal electrode With 2 mm connector			

Accessories

Part-No.	Specification	Use	Application	Remarks
45-39-151 EH725	Therapy set 1			Patient cable, 2 tin plate electrodes, 2 rubber straps 42 cm, 2 rubber straps 135 cm, 4 fastening buttons, 10 sponge pockets
45-39-169 EH725	Therapy set 2			Patient cable, 2 rubber electrodes 120 x 60, 2 elastic straps with velcro fastener, 10 sponge pockets
45-39-177 EH725	Diagnosis set			Patient cable, 3 tin sheets 18 x 28 cm, 2 rubber straps 42 cm, 2 rubber straps 135 cm, 4 fastening buttons, 6 intermediate fleece layers
026-0-3000	UNICAR® Basic rack with large device tray incl. Drawer			
026-0-3020	Cable and electrode tray for Stereodylator			






Accessories

Part-No.	Specification	Use	Application	Remarks
45-38-856 EH722	STEREO- Suction electrode set, 3 circuits			Contains Suction- electrode I, STEREO Suction- electrode II, Paper sheets (100 pcs.)
45-38-534 EH722	STEREO- suction head I with cable and connector			
45-38-542 EH722	STEREO- suction head II with cable and connector			
45-38-674 EH722	Multi circuit suction cup	For all DUO- and STEREO- suction heads		
45-38-600 EH722	Electrode button	For all DUO- and STEREO- suction heads		
71-21-841 EH060	Paper sheets (100 pcs.)	Consumables for DUO and STEREO suction electrodes		Molsten well: Replace frequently




Accessories

	Part-No.	Specification	Use	Application	Remarks
	45-38-716 EH725	Star electrode pair standard, diameter 130 mm		With rubber- or elastic-strap	
	45-38-708 EH725	Star electrode pair small, diameter 80 mm		With rubber- or elastic-strap	
	45-38-724 EH725	Star electrode pair large, diameter 180 mm		With rubber- or elastic-strap	
	45-39-260 EH719	Standard accessory kit			Electrode set, 10 pockets, 2 rubber straps small 135 cm, 2 rubber straps large 135 cm, 4 fixing buttons
	49-07-309 E3020	Electrode pocket standard for star electrode, 10 pcs.			Moisten well
	49-07-317 E3020	Electrode pocket small for star electrode, 10 pcs.			Moisten well
	49-07-325 E3020	Electrode pocket large for star electrode, 10 pcs.			Moisten well


Accessories (Ultrasound)

	Part-No.	Specification	Use	Application	Remarks
	023-0-0141	Ultrasound probe with cable and plug, emitting surface 5 cm ²			
	023-0-0140	Ultrasound probe with cable and plug, emitting surface 2,5 cm ²			
	025-0-0141	Ultrasound probe with cable and plug, emitting surface 5 cm ²			Alternative part for 023-0-0141 used with another probe holder
	025-0-0140	Ultrasound probe with cable and plug, emitting surface 2,5 cm ²			Alternative part for 023-0-0140 used with another probe holder
	45-39-128 EH725	Ultrasonic contact gel 250 ml			

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
Interference Current Stimulation Device

Stereodynator

Instructions for use

gbo Medizintechnik AG has taken care in the preparation of this manual, but does not assume any liability, expressed or implied, of any kind nor does it assume any responsibility for errors or omissions.

All rights reserved. No part of this manual may be reproduced, in any form or by any means (electronic, mechanical, or otherwise) without the prior written permission of gbo Medizintechnik AG.

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Table of Contents

1 INTRODUCTION	5
1.1 Intended Use.....	5
1.2 Note concerning the operating personnel.....	5
1.3 Admission.....	5
1.4 Device description.....	5
2 OPERATING CONCEPT	7
2.1 Control Panel.....	7
2.1.1 Touch screen with direct help system.....	7
2.1.2 Frequently used therapies.....	8
2.1.3 Direct help menu.....	9
2.2 Intensity regulator.....	10
3 START OF OPERATION	11
3.1 Transport and Assembly.....	11
3.2 Connecting and Switch-On.....	11
3.2.1 Device Settings.....	12
3.3 Placing out of operation.....	13
4 THERAPY	13
5 ELECTRODES	14
5.1 Electrode Positioning.....	14
5.1.1 Monophase Electrode Techniques.....	15
5.1.2 Biphasic Electrode Techniques.....	15
5.2 Single-Pole Electrodes.....	16
5.3 Three-pole electrodes.....	17
5.4 Suction Electrodes.....	18
5.5 Electrode Handpiece.....	19
6 ULTRASOUND THERAPY MODULE	20
7 MAINTENANCE	21
7.1 Safety Controls.....	21
7.2 Cleaning, Disinfection and Care.....	22
7.2.1 Cleaning of the Device.....	22
7.2.2 Cleaning of the Electrodes.....	22
7.2.3 Disinfecting the Electrodes.....	22
7.2.4 Cleaning of the Suction Electrodes.....	22
7.2.5 Disinfecting of the Suction Electrodes.....	23
7.2.6 Cleaning of the elastic straps.....	23
7.2.7 Cleaning of the Electrode Sponges and Electrode Pockets.....	23
8 WARNINGS AND SAFETY PRECAUTIONS	24
8.1 Contraindications for use of ultrasound.....	25
8.2 Contraindications for use of current stimulation.....	27
9 EXPLANATION OF THE SIGNS USED	28

10 TECHNICAL DATA	29
Current types.....	30
11 ACCESSORIES	33
12 TROUBLESHOOTING	34
12.1 Further errors.....	34
13 APPENDIX	35
14 INDEX	38

Warnings and safety precautions

Warning!

Warnings which have to be observed by all means!

Caution!

Observe the instructions for use!

Note!

Information that will facilitate your work.

1 Introduction

1.1 Intended Use

Therapy with low frequency, medium frequency and interference currents as well as ultrasound therapy incl. combination therapy current stimulation - Ultrasound.

1.2 Note concerning the operating personnel

The operation of the device must be performed only by trained medical professionals.

1.3 Admission

We recommend that the operator should be trained when using the unit.

1.4 Device description

The device with interference current named **Stereodynator** is a microprocessor controlled electrical stimulus device for Electro Therapy. The wide range of usage qualifies this medical device for use in physiotherapeutic departments of clinics as well as in modern well-equipped private medical practices.

All accessories of the previous series can also be used with the new generation in the same manner. A suction application aid for vacuum therapy as well as a module for ultrasound therapy with continuous and impulse ultrasound waves are optionally available. The control elements for suction application aid and ultrasound have already been integrated in the software and can be activated by upgrading the device.

A stereodynamic interference current is generated through the superposition of three middle-frequency currents flowing in different directions. The stimulating lower frequencies are generated through interference of two phase shifted circuits in the area of superposition. The additional third circuit generates, as opposed to classical interference methods, on one side a slow change in intensity and on the other side a rhythmic shifting of the interference field. These dynamics of the stimulation location and the intensity decrease the habituation effect and therefore improve the therapeutic effect. The **Stereodynator** uses the three-dimensional interference method. The special characteristics are as follows:

- Local stimulation effect.
- Multi-site stimulation effect.
- Intensity dynamics.
- Dynamic behavior of stimulation site coupled with endogenous/exogenous stimulation shifting.

Programs improve the handling: The **Stereodynator** is equipped with several programs which automatically set the treatment frequency and the treatment time.

The stimulating effect can alternatively be administered endogenously or exogenously. During endogenous application, the interference is generated through superposition of the electrical circuits in the body. This allows an intense effect. During exogenous application, the interference is created in the medical device. The stimulating effect takes place directly underneath the electrodes on the body surface.

In addition to the three-circuit interference current, the **Stereodynator** offers a complete selection of single circuit currents for all known therapeutic procedures.

The **Stereodynator** is particularly suitable for:

- Pain therapy with three-dimensional interference current.
- Muscle toning and muscle detoning.
- Galvanization and iontophoresis.
- Pelvic floor stimulation.
- Treatment of urinary and fecal incontinence.
- Treatment of paralyses with complete or partial muscular degeneration.
- Treatment of atrophies due to inactivity or weakened muscles after longer periods of inactivity.
- Electrical stimulation therapy without electrolytic side effects and only slight muscle fatigue.
- Treatment of pain, muscle spasms, functional diseases of the locomotor system, such as sports injuries, peripheral circulatory disturbances, influencing of the vegetative system with diadynamic currents, ultrastimulation current, microstimulation current, TENS- and TENS Burst currents, high-voltage currents and middle-frequency currents.

In **medical diagnosis** it is suitable for:

- Qualitative and quantitative determination of faradic excitability.
- Determination of the rheobase, chronaxy and accommodability.
- Recording of stimulus intensity/stimulus time characteristics (I/T curves) with defined, measurable and reproducible rectangular and triangular pulses (exponential current).
- MF test according to Dr. Lange.
- Extensive non-invasive electro-diagnostics of peripheral paralyses.
- Neurodiagnostic examinations with galvano-palpations.

2 Operating Concept

The **Stereodynator** is equipped with a touch screen display, one on/off switch, four connecting plugs for patient cables, two connecting plugs for ultrasound probes and the intensity regulator.

The **Stereodynator** is being operated via the control panel of the large touch screen. The user will obtain explanations regarding the device's operation and functionality of individual elements via a direct help system during operation. Consequently, operation is explained in broad outline only.

2.1 Control Panel


Figure 1: Front view

2.1.1 Touch screen with direct help system

The touch screen should be operated with marked pressure of the fingertip. Pointed objects are unsuited for operation and can destroy the touch screen.

Every button is clearly labeled. Through the optics, the user can differentiate operating buttons from non-operating buttons.

In the case of doubt the user can obtain further information via the direct help system.

In order to adjust numerical data the arrow keys may either be pressed repeatedly or pressed constantly for automatic increase/decrease of the value.

The scrollbars are used like with a PC: press the scrollbar and then pull into the desired direction.

2.1.2 Frequently used therapies

All frequently used functions are directly accessible on the front screen. To this purpose, the seven most used therapies and the three therapies applied last are listed on the right side. From here, they can directly be selected. An additional symbol indicates whether a certain therapy is being used more often (upwards arrow) and also names the type of therapy (electro therapy, stereo-, combination- or ultrasound therapy).


Figure 2 View of Touch screen display.

During daily operation, the user can thus recall typically used therapies with only one touch. Additional therapies can either be recalled according to a patient's name, i.e. via log book, or after indications or name of therapy.

2.1.3 Direct help menu

The user of the **Stereodynoter** is guided on the screen. The direct help menu gives information about the selected buttons.

To open the direct help menu, first touch the „?“-button on the screen and then the button about which you desire information. A new window opens where the information is shown. In order to close this window you need to touch the “X”-button.

Available buttons are indicated in light grey. Unavailable buttons are indicated in dark grey. If you push an unavailable button, information about the use of this button is shown on the screen automatically.


Figure 3: Direct help menu as shown on the display

2.2 Intensity regulator

During therapy the output voltage is set with the **Intensity regulator**. It is furnished in the form of a rotor pulse generator. The current increases by turning clockwise and it is reduced by turning counter-clockwise. The current intensity is shown on the display in mA. In order to reach a certain intensity the user may have to turn the regulator clockwise for more than one rotation if necessary.

3 Start of Operation

3.1 Transport and Assembly

After unpacking check all components for external integrity.

The **Stereodylator** is a portable unit. To place the unit, each flat surface is appropriate. It must not be placed on the floor. Keep a wall distance of at least 20 cm. The device should neither be placed in front of radiators nor should it be covered by pillows or blankets while in operation. The lower side of the unit must be a plane area. The device is not made for outdoor operation.

For maximum usability there is a special device cart available (see chapter 11, accessories) which matches in design and function with the new **Stereodylator**. This device cart is perfect for the placement of the unit and various accessories. With the device cart the **Stereodylator** is mobile and always ready for use. For therapy it can easily be moved to the patient.

The **Stereodylator** corresponds to the regulations DIN/VDE 0750, EN 60601. It is a device of protection class I. Within the scope of the Medical Device Directive (MDD) the current stimulation device belongs to class IIa.

Warning!

Hints for the clinical environment

- The unit is not designed to be operated in places with the inherent risk of explosions. If it is used in dangerous areas of anaesthesia departments, the possibility of an explosion cannot be excluded.

If the patient and/or the patient cable is directly exposed to a radiator of a medical device for high frequency heat therapy, damage of the device or danger to the patient cannot be excluded. As a rule, a clearance distance of 3 m is sufficient.

3.2 Connecting and Switch-On

The current stimulation device has been adjusted for the connection to a supply voltage of 100 to 240 V. It is not necessary to switch over the voltage – the device adjusts automatically to the right voltage.

Irrespective of the adjusted supply voltage, the device is appropriate for mains frequencies of 50 to 60 Hz.

Connect the **Stereodylator** with the mains cable to a socket with protective ground. The protective ground must work correct

The device is switched on by the main switch on the back of the device.

Warning!

To avoid the risk of electric shock, this unit must be connected only to a supply with protective conductor!

Warning!

It is only allowed to connect USB-Sticks for data backup and update!
Never connect a PC or another USB-device to this connector!

It is not allowed to use the two holder for Ultrasound probes as carrier handle!


Figure 4: Back and side view of the device

3.2.1 Device Settings

You reach the device settings menu by pressing the button „settings“ (see Fig. 2). In the device settings menu all settings can be performed. With the direct help function you are able to retrieve all information corresponding to the special menu item.

3.3 Placing out of operation

In order to disconnect the device just disconnect it from the mains power supply. No other measures are to be taken.

4 Therapy

The **Stereodylator** guides the user through the application of the therapy. With the direct help menu all buttons and their functions are explained directly on the screen.

On the right side of the screen therapy programs can be selected directly from this list of favourites or with the sub-menus "More Therapies" or "Indications".

First, choose the channel you want to use and then you can select the desired current type with the program lists. The direct help menu will give you information about the following steps of the therapy.


Figure 5: Touchscreen with chosen therapy program (indication)

5 Electrodes

Basically all standard electrodes are connectable in electrotherapy. Plate electrodes and rubber electrodes are appropriate for large-area treatments.

The electrode handpiece with fingertip switch and connecting cable is particularly suitable to **make a diagnosis**. With the switch single pulses can be started. With the handpiece also small areas can be treated. Additionally you need an electrode set.

The **therapy** can be executed without any problem with suction electrodes or plate electrodes. For each indication of the chosen therapy program you find a picture on the screen of the **Stereodylator** which shows how to apply the electrodes.

5.1 Electrode Positioning

The **positioning of the patient** and/or muscle to be treated is very important. The therapy shall always be carried out in a comfortable and relaxed position. The joints shall be placed in an angular position, in a mid-position, from which both the flexor as well as the extensor muscles can be stimulated.

In general the preparation and application of the electrodes is carried out as follows:

1. Always put a wet electrode intermediate layer or electrode pocket between the electrode and the skin in order to avoid current sensations on the skin. The layer must always be about 1 cm larger than the electrode itself.
2. Put the electrode completely into the electrode pocket, so that the pocket protrudes on the open side with about 1 cm.
3. After multiple use the electrode intermediate layer or electrode pocket can slightly shrink. If the distance of about 1 cm is no longer given, the intermediate layer or pocket must be replaced.
4. Prior to the initial use, the macroporous intermediate layers and pockets made of viscose sponge must be thoroughly washed out under running water to remove manufacturing residues.
5. To moisten the intermediate layer or pocket put them into normal tap water, or better into 1% salt solution. The intermediate layers or pockets should be soaked, then removed and dropped; do **not** wring them out.
6. If necessary, flatten deformed electrodes, insert them into the pocket and/or put them on the intermediate layer. Place them on the area to be treated and fit them to the given body surface. Use only thick and well moistened sponges.
7. Apply the holding strap in such a manner that the electrode fits perfectly.
8. The sponge pockets or intermediate layers have to be cleaned after the expiration of treatment, (see chapter 7.2.4).

**Note!**

- The above mentioned applies for the application of suction electrodes accordingly!

**Warning!**

- Do not apply the electrodes on injured skin. Even minor abrasions can cause a burning to the patient. Herewith the intensity of current by the patient can be erroneously judged. If a treatment cannot be avoided, apply zinc ointment to cover.
- It is recommended not to exceed the current density of 2 mA/cm² in all electrode surfaces. In case of galvanization the limiting value is 0.2 mA/cm² and in case of iontophoresis it is 0.1 mA/cm².
- Do not use disposable (self-adhesive) electrodes in case of currents with galvanic portion, there is the danger of acid burns!!

5.1.1 Monophase Electrode Techniques

In case of monophase techniques, direct and indirect stimulation are distinguished. Generally the cathode is used as a stimulating electrode. In case of direct stimulation it is applied on the nervous stimulant point (all muscles react) and in case of indirect stimulation on the muscle stimulant point (only the stimulated muscle reacts). Apply the electrolytically inactive electrode (mostly a large anode) distally.

For differentiated therapeutic treatments and easy electro-diagnostics it is recommended to use the electrode handpiece as an electrolytically active electrode.

5.1.2 Biphas Electrode Techniques

A characteristic of this method is that both (mostly equally large) electrodes are applied on the muscle to be stimulated or on the functionally related muscle group. The cathode is to be placed distally.

5.2 Single-Pole Electrodes

The use of single-pole electrodes is applicable for all therapy currents. For large-area treatment we recommend the use of plate electrodes, for all other applications the use of self-cut electrodes. The last-mentioned need an electrode connector for current supply. Both electrode types have to be connected to the device with a patient cable.

The electrode size depends on the area that should be passed by the current. Under small electrodes the current flow will be more concentrated and more localized than when applying large electrodes. The application area depends on the clinical picture.

**Warning!**

As far as current types are concerned for which there exists the threat of cauterization (currents with galvanic portion), the maximum recommended current density is 0.5 mA/cm² electrode surface. In case of galvanization the limit value is 0.2 mA/cm² and in case of iontophoresis 0.1 mA/cm².

Electrode	Surface	Max. Current / Galv / Iontophoresis
Conductive rubber electrode 3 cm x 4 cm	12 cm ²	24 mA / 2.4 mA / not useful/ 1.2 mA
Conductive rubber electrode 6 cm x 8 cm	48 cm ²	96 mA / 9.6 mA / not useful/ 4.8 mA
Conductive rubber electrode 8 cm x 12 cm	96 cm ²	192 mA / 19.2 mA / 9.6 mA
Plate electrode 9 cm x 12 cm	108 cm ²	216 mA / 21.6 mA / 10.8 mA

Table 1: Maximum Current of conductive rubber electrodes and plate electrodes

**Note!**

- A current density of more than 2 mA/cm² on all electrode surfaces requires a higher attention of the operator.

5.3 Three-pole electrodes

The triple-pole electrodes can be used for stereodynamic interference currents.

The conductive rubber electrodes of the three circuits are placed in two flexible three-pole star electrode sets (I and II) and are attached to a common patient cable. There are three sizes available (small, standard and large). Electrode pockets are necessary for the application and fixing material for the bandage (see chapter 11 *Accessories*).

Each stereodynamic electrode set is placed into a well moistened pocket so that the black conductive rubber electrodes face the humidity carrier.

Electrode position:

- A)** If the electrodes are placed opposite to each other (**transverse current flow, transregional**) the cable connections must point at the opposite direction.
- B)** If the electrodes are placed on one side (**longitudinal current flow**), their cables must point at the same direction.

Only if these electrode positions are followed, the flow paths of current of the three circuits can transpose and optimal stereodynamic stimulation currents are generated.

In special cases (e.g. depending on the application in the case of shoulder-arm syndrome) it is possible to deviate from this general rule.

You can also use triple-pole suction electrodes if your device is equipped with a suction application aid. This also applies to the application with suction electrodes.


Figure 6:
Three-pole electrodes

5.4 Suction Electrodes

Instead of the above-mentioned electrodes, suction electrodes can be used for faster application. Depending on the device configuration, the suction application aid is either optional or standard (part no. 025-0-1000-V or 025-0-1000-UV).

Warning!

- For current types for which there exists the threat of cauterization (currents with galvanic portion), the maximum recommended current density is 0.5 mA/cm² electrode surface. In case of galvanization the limit value is 0.2 mA/cm².
- In case of iontophoresis do not use suction electrodes.

Suction electrode	Surface	Max. current value / Galv
Electrode plate, small (Sponge = Ø 36 mm)	10 cm ²	20 mA / <i>not useful</i> / 2 mA
Electrode plate, standard (Sponge = Ø 52 mm)	21 cm ²	42 mA / <i>not useful</i> / 4,2 mA
Electrode plate, large (Sponge = Ø 72 mm)	40 cm ²	80 mA / 8 mA

Table 2: Maximum current of the suction electrodes

Caution!

- Set the vacuum intensity only as high that the vacuum electrodes keep in a good way.
A much too high vacuum intensity may cause a hematoma.

5.5 Electrode Handpiece

The electrode handpiece is made for simple electro-diagnostics. The electrode handpiece (electrolytically active electrode) is firmly connected to the patient cable. With a second wire the plate electrode is connected (electrolytically inactive electrode). In addition an electrode set is required consisting of metal shaft, electrode pockets and rubber rings.

The following procedure is advisable:

1. Connect the electrode handpiece.
2. The electrode handpiece will be recognized by the device.
3. Place the inactive electrode to the patient
4. Select diagnostic current
5. Increase Intensity and start single pulse

Note!

The fingertip electrode only works with currents of the diagnostic range. The push button of the handpiece electrode works as start button for a single pulse.

6 Ultrasound Therapy Module

Depending on the device configuration, the ultrasound therapy module is either an optional or a standard device (part no. 025-0-1000-U or 025-0-1000-UV). With the Ultrasound probe either Ultrasound therapies as also combination therapies with Ultrasound and current can be done.

Note!

- The ultrasound frequency can be changed to 1 or 3 MHz.
- The frequency for pulsed mode is 100 Hz.
- The selectable parameters are continuous or pulsed mode 1:1, 1:2, 1:5, 1:10 or 1:20.
- The coupling status of the probe is shown by an LED on the probe housing and also as a graphic illustration on the display. The LED will be off in case of insufficient coupling. The LED will be on in case of good contact.
- As a standard setting the ultrasound power is shown in W/cm^2 on the display.

Warning!

- Do not forget the contact gel!
- The ultrasound probe is to be treated with care. External influences such as a mechanical shock or impact can alter the characteristics of the ultrasound probe. We recommend carrying out a visual examination at least once a year to check for fissures that would allow liquids to soak in, as well as to make sure that the cables and connectors are flawless.

	Power for			
	5 cm ² ultrasound probe		2.5 cm ² ultrasound probe	
1 MHz	0.5 to 15 W	0.1 to 3 W/cm ²	0.1 to 7.5 W	0.1 to 3 W/cm ²
3 MHz	0.1 to 7.5 W	0.1 to 1.5 W/cm ²	0.1 to 3.75 W	0.1 to 1.5 W/cm ²

Table 3: Setting options of the ultrasound power, increments are shown in bold.

Note!

Unless the ultrasound probe is coupled to the patient, no therapy time will pass. When the probe is coupled again, the therapy time will resume. Hence, the ultrasound probe can be removed during the treatment.

Note!

The Ultrasound probe can be exchanged. After switching on the unit or respectively after plugging the Ultrasound probe each probe will be calibrated individually.

7 Maintenance

Functionality, reliability and safety characteristics of the **Stereodynator** are guaranteed only upon handling the device in accordance with the operating instructions. Safety control, maintenance work, repair work and modifications must only be carried out by the manufacturer or by service agents authorized by him. In case of a failure, parts which influence the safety of the device must only be replaced by original spare parts of the manufacturer. The electric installation must correspond to the requirements in accordance with VDE/IEC. **The device does not contain any parts which need maintenance work done by the user.**

Warning!

No parts of the device are allowed to be serviced while a patient is connected.

7.1 Safety Controls

The device is subject to the provisions of the Medical Device Directive. The safety controls have to be carried out on the basis of this directive. Thereby, the operator regulation has to be observed in particular.

Irrespective of the legal rules or beyond the scope of the Medical Device Directive, it is recommended to have the device checked by the manufacturer or by a service agency authorized by him at 12- months intervals.

The check shall consist of at least the following criteria:

- Electrical safety check in accordance with the test plan of the manufacturer.
- Check of the device in respect of external integrity.
- Check of all display and operating elements in respect of damages.
- Check of all inscriptions in respect of legibility.

7.2 Cleaning, Disinfection and Care

7.2.1 Cleaning of the Device

For cleaning and disinfection of the current stimulation device and its accessories (except electrode fleece, electrode pockets or felts, (see chapter 7.2.3) there should not be used any agents containing higher portions of phenol derivatives, alcohol, compounds of chlorine or peracetic acid. It is recommended to use disinfectants on aldehyde basis. Never use an abrasive. Fingerprints on the touchscreen should be cleaned with a dry cloth or some isopropanol.

Warning!

Prior to cleaning or disinfection unplug the mains plug out of the socket!

The current stimulation device is suited for wiping disinfection. Make sure that no liquids soak the device. Under no circumstances the plug or socket must get wet. Do not sprinkle the device for cleaning or disinfection. The device is not suited for heat sterilization or for sterilization with gases.

7.2.2 Cleaning of the Electrodes

You can clean the electrodes after a treatment by using some warm water (appr. 40 °C) and a soap suds. Dry the electrodes with a cleaning tissue or let them air-dry.

7.2.3 Disinfecting the Electrodes

The rubber electrodes can be disinfected by spray/wipe disinfection. Spray onto the black side of the electrodes, leave to soak and wipe off with a lint free cloth. (Please follow the recommendations of the manufacturer.)

We recommend Bacillol-AF for the disinfection of the electrodes.

Note!

The conductivity of the electrodes is optimized by the use of a certain amount of graphite. Consequently, when using and disinfecting the electrodes, black color may come off.

7.2.4 Cleaning of the Suction Electrodes

You can clean the electrodes after a treatment by using some warm water (appr. 40 °C) and a soap suds. Dry the electrodes with a cleaning tissue or let them air-dry.

It is possible to disassemble the single-pole suction electrodes to clean the silicon ring and the metal plate separately. In case of strong stains you may clean the metal plate with a sponge for cleaning of stainless steel.

The silicon ring of the three-pole suction electrodes can be removed and cleaned separately. In case of strong stains the surface of the three plates may be cleaned with a sponge for cleaning of stainless steel.

7.2.5 Disinfecting of the Suction Electrodes

The suction electrodes can be disinfected by spray/wipe disinfection. Spray onto the electrodes, leave to soak and wipe off with a lint free cloth. (Please follow the recommendations of the manufacturer.)

We recommend Bacillol-AF for the disinfection of the electrodes.

7.2.6 Cleaning of the elastic straps

The elastic straps can be washed in a washing machine. Please observe the washing instructions of the instruction label or on the package.

7.2.7 Cleaning of the Electrode Sponges and Electrode Pockets

The Electrode sponges and pockets can be washed in a washing machine with 95 °C, or be disinfected with disinfectant solution. We recommend Dismozon[®] pur (chemical agent: Magnesium monoperoxyphthalat Hexahydrat).

Note!

- After disinfection Electrode sponges and pockets have to be wringed out closely at least 6 times under warm running water!

8 Warnings and Safety Precautions

Warning!

- For patients with implanted electronic device electrical stimulation treatment is to be carried out only after having checked any risks.
- Jewellery and eyeglasses have to be taken off during the treatment.
- Turn off cellular phones and radiophones or place them in a distance of 3 m from the device.
- Active implants (cardiac pacemakers, defibrillators...) are extremely vulnerable. Here the therapy should only be carried out under continuous pulse and ECG control.
- If the patient and/or the patient cable is in the direct range of a high-frequency, short-wave or micro-wave therapeutic device, a damage to the device or an injury to the patient cannot be excluded. Please keep a clearance distance of 3 m.
- A simultaneous connection of the patient to a high-frequency surgery device can lead to burns under the electrical stimulus electrodes.
- It is recommended not to exceed the current density of 2 mA/cm² on all electrode surfaces.
- For current types with the threat of cauterization (all currents with galvanic portion), the maximum recommended current density is 0.5 mA/cm² electrode surface. For galvanization the limit value is 0.2 mA/cm² and for iontophoresis 0.1 mA/cm².
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- Do not use disposable (self-adhesive) electrodes with currents with galvanic portion since there is the threat of cauterization!!
- The unit is not designed to be operated in places with the inherent risk of explosions. If it is used in dangerous areas of anesthesia departments, the possibility of an explosion cannot be excluded.
- In case of short circuit the lithium battery, the battery can explode. Only service staff should change the battery.
- In case of any visible failure contact gbo Medizintechnik AG or one of the service agencies authorized by gbo Medizintechnik AG immediately.

**Note!**

The current stimulation treatment should not be done thru or over the head. Not directly on the eyes, not covering the mouth, not at the front side of throat (especially the carotid sinus), or with electrodes on breast and upper back, or with electrode surfaces crossing the heart..

8.1 Contraindications for use of ultrasound

An exact diagnosis must be made before the therapy is started. The diagnosis determines the specific treatment with ultrasound.

Regardless of the low dosage and the relatively wide field of application, distinct contraindications have become evident. The "harmfulness" frequently discussed in the literature is of secondary importance, because no organic damage can be caused with the stated low dose, provided that the static exposure to ultrasound is avoided.

**Warning!****Ultrasound is contraindicated in case of:**

- Changes to the skin, particularly with infectious diseases and birthmarks
- Tumorous diseases of all stages
- Feverish conditions
- Poor general condition and general atrophy
- Active tuberculosis, regardless of stage and localization
- Acute inflammations
- Stomach ulcers
- Following recent Thorium-X treatment, X-ray depth treatment
- Diabetes mellitus
- Pregnancy
- Vascular diseases of the extremities (thrombophlebitis, thrombosis, varicosis)
- Disorders of coagulation of blood
- Acute joint rheumatism.

**Warning!**

In addition to this group of diseases one should aim to exclude certain organs from a direct treatment. **Do not expose the following organs to direct treatment:**

- Eyes, brain and spinal marrow
- Laminectomy-related spinal incisions
- Anaesthetised areas
- Heart and lungs
- **No exposure** to heart segments with functional heart complaints.
- **No exposure** to the epiphysis zones of children.

**Note!**

- This list is not exhaustive. In individual cases the attending doctor should decide on contraindications and criteria for the treatment.
- In case of all visible failures contact gbo Medizintechnik AG or one of the service agencies authorized by gbo Medizintechnik AG immediately.

8.2 Contraindications for use of current stimulation

Despite the low dosage is definite contraindications have emerged.

	<p>Warning!</p> <p>Stimulation should not be applied in the following syndromes:</p> <ul style="list-style-type: none"> • Open wounds in the treatment area • Malignant tumors in the treatment area • Feverish conditions <p>Stimulation current to be applied to the following symptoms only after consultation with the doctor:</p> <ul style="list-style-type: none"> • Patients with implanted electronic devices • Patients with metallic prostheses (only biphasic currents should be used, preferably medium-frequency current used) • Pregnancy
--	---

	<p>Additional precautions!</p> <ul style="list-style-type: none"> • In patients with diagnosed heart problems • In treatment on areas of the skin which lack normal sensation assets • For muscle stimulation after recently running operations in which muscle contraction may disrupt the healing process
--	---

	<p>Note!</p> <ul style="list-style-type: none"> • This list is not exhaustive. In individual cases the attending doctor should decide on contraindications and criteria for the treatment. • In case of all visible failures contact gbo Medizintechnik AG or one of the service agencies authorized by gbo Medizintechnik AG immediately.
--	---

9 Explanation of the Signs used

CE – conformity sign

Operating instructions!

Follow operating instructions!

Application part ungrounded, protection degree Type BF.

This product complies with WEEE Directive 2002/96/EG (waste electrical and electronic equipment). Separate collection for electrical and electronic equipment. The waste removal at the end of the service life will be done by the manufacturer.

Ultrasound probe connection

Patient cable connection

Vac I ⊕

Suction electrode I connection

Vac II ⊖

Suction electrode II connection

10 Technical Data

Mains voltage and frequency:	100 – 240 V, 50-60 Hz	
Power consumption*:	max. 200VA	
Mains fuse :	T5A H250V	
Mode of operation:	Continuous operation	
MDD device class:	IIa	
Protection class	I	
Type of applied part:	BF	
Protection class of housing:	IP X0	
Dimensions*:	max. 36 cm x 28.5 cm x 37.5 cm (W x H x D)	
Weight*:	as base device: as –V variant: 13.2 kg without accessories as –UV variant: 14 kg without accessories	
Color:	Aluminium nature anodized and grey RAL 7016	
Display:	15" TFT LCD Touchscreen	
Battery:	CR 2032	
environmental conditions:	operation of the device:	temperature range +10 °C ... +35 °C relative humidity of air 30 ... 75 %
	transport and storage:	temperature range +5 °C ... +50 °C relative humidity of air < 90 %, non-condensing
Current stim. module:	max. 125 mA _{peak} ± 5%	
Output current:	max. 200 V ± 5%	
Output voltage:	See chapter 0	
Current types:	Carrier Frequency	2kHz – 10kHz ± 1% Sinus wave
	Low frequency	0.1Hz – 200Hz ± 1%
	Accuracy	± 10%
Applied part	Contact surfaces at the electrodes.	
Therap. ultrasound module:	3W/cm ² at 1 MHz, 1.5 W/cm ² at 3.3 MHz	
Maximum Ultrasound Intensity:	1 MHz: +0% / -20% all other frequencies: ± 20%	
Accuracy of the output:	5 cm ² and 2.5 cm ² ± 20%	
Effective radiating area:	CW, 5%, 10%, 20%, 30% and 50%	
Pulse ratio:	100 Hz ± 5%	
Pulse sound frequency:	1 MHz ± 10% / 3.3 MHz ± 10%	
Ultrasonic frequency:	IP X7	
Protection class of probe:	43°C	
Max surface temperature:	Titanium surface of the ultrasonic head	
Applied part		
Suction application module:	180 mbar ± 20% in maximum position	
Suction:	0 - 30 pulses per minute	
Suction massage frequency:		

*) dependent on the device type and/or whether the device is equipped with therapeutic ultrasound modules and suction application modules.

By request of technical personnel gbo Medizintechnik can offer spare part lists and circuit diagrams.

The IEC connector is used for all-pole disconnection from the mains.

The exchange of components by service personnel is described in the service manual.

gbo Medizintechnik AG reserves the right to modify design and specifications without prior notice.

Current types

The **Stereodynator** offers a large variety of current types for therapy. For details see the direct help menu on the display of the device. The following current types are available as therapeutic programs with the **Stereodynator**:

Current	Graphics	Pulse duration T in ms	Pause duration R in ms	D.C. portion in %	Frequency in Hz	Max. output Voltage V _{peak}	Max. output current mA _{peak}
Diadynamic currents							
CP (Courtes périodes)				47.25	50 / 100	200	80
CP-id (MF-part lowered)				45	50 / 100	200	80
DF (Diphase fixe)				63	100	200	80
LP (Longues périodes)				46	50/100	200	80
MF (Monophasé fixe)				31.5	50	200	80
MM (Monophasé modulé)				20	100	200	80
RS (Rythme syncope)				15.25	50	200	80
Pulsed currents							
Galvanic current				100	0	70	70
High voltage		0.02 – 0.06	5 ms – 1 s	0	1 – 200	200	80
Impulse galvanization IG 30		30.0	50.0	18	12.5	200	80
Impulse galvanization IG 50		50.0	70.0	0.3	8.3	200	80
Frequency modulation FM / 7 – 14 Hz		1	70.9 – 142.4	0.4	7 - 14	200	80
Micro current biphasic		2,5 - 166	2,5 - 166	0	3 – 200	200	80

Current	Graphics	Pulse duration T in ms	Pause duration R in ms	D.C. portion in %	Frequency in Hz	Max. output Voltage Vpeak	Max. output current mApeak
Mikro current monophasic		2,5 - 166	2,5 - 166	50	3 - 200	200	80
Faradic current		1,0	19,0	2,6	50	200	80
Ultrastimulation current		2	5	28,6	143	200	80
TENS biphasic		0,02 - 0,2	5 - 1000	0	1 - 200	200	80
TENS monophasic		0,02 - 0,2	5 - 1000	max 3,8	1 - 200	200	80
TENS Burst biphasic		0,02 - 0,2	5 - 1000	0	1 - 200	200	80
TENS Burst monophasic		0,02 - 0,2	5 - 1000	max 3,8	1 - 200	200	80
Middle frequency current 8 kHz							
Russian Stimulation (Kofs)		1 - 1000	1 - 1000	0	0,5 - 200	200	80
Gym Current			1 - 10	0	0 - 200	200	80
MF sinus modulated				0	0,1 - 200	200	80
Interference currents							
Interference				0	0,1 - 200	200	80
Dipol Vector				0	0,1 - 200	200	80
Isoplanar Vector				0	0,1 - 200	200	80
Stereodynamic Interference currents							
Stereodyn. Interference				0	0,1 - 200	200	80

Current	Graphics	Pulse duration T in ms	Pause duration R in ms	D.C. portion in %	Frequency in Hz	Max. output Voltage Vpeak	Max. output current mApeak
3 Channel currents							
Gym Burning Mode					0,1 - 200	200	80
Gym Current			1 - 10		0,1 - 200	200	80
Diagnoses							
Monophasic (T/R) single pulse mode		0,2 - 2000				200	80
Biphasic (T/R) single pulse mode		0,2 - 2000				200	80
I/T-Curve		0,2 - 1000				200	80
Rheobase/Chronaxe		0,2 - 1000				200	40/80
MF-Test acc. to Dr. Lange						200	80

Impedance: 50 - 5000 Ohm

Carrier frequency: 2 - 10 kHz sinusoidal

11 Accessories

For the Stereodylator we offer a large variety of accessories.

specification	part number
Standard accessory kit for Stereodylator® Electrode set, 10 pockets, 2 rubber straps small, 2 rubber straps large, 4 fixing buttons	45-38-260EH719
Therapy set 2 for Stereodylator® Patient cable, 2 rubber electrodes, 2 elastic straps with velcro fastener, 10 sponge pockets	45-39-169EH725
Suction electrode set single pole standard for Stereodylator® diameter 50 mm	45-38-500EH722
STEREO-suction electrode set 3 circuits for Stereodylator® STEREO-suction electrode I, STEREO-suction electrode II, paper sheets (100 pcs)	45-38-856EH722
External air filter	008-1-0442

Detailed information can be found in the brochure "Stereodylator and accessories". This brochure is available as a pdf-file by download on our website www.gbo-med.de.

The following versions of the Stereodylator are available:

specification	part number
Stereodylator®	025-0-1000
Stereodylator® incl. ultrasonic therapy module , incl. 1 ultrasonic transducer 5 cm ² , incl. 1 contact gel, 250 ml	025-0-1000-U
Stereodylator® incl. suction application aid	025-0-1000-V
Stereodylator® incl. ultrasonic therapy module , incl. 1 ultrasonic transducer 5 cm ² , incl. 1 contact gel, 250 ml and incl. suction application aid	025-0-1000-UV
Device cart for Stereodylator®	026-0-1000



Figure 7:
Stereodylator
with
accessories and
device cart

Note!

Use only original accessories
to ensure the safe operation of
the device.

12 Troubleshooting

Each problem that occurs during operation of the unit will be shown in a message window on top of the display and also signaled by an acoustic tone. Most of the problems can be solved by the instructions displayed.

In general:

1. The malfunction is shown on the display.
2. An acoustic signal is heard.
3. Follow the instructions on the display.

Suggestions:

- Turn off and turn on the unit.
- A full selftest of all device functions is performed and the unit is reinitialized. If the unit does not come up to normal state, please call your local service.

12.1 Further errors

symptom	cause and action
Device cannot be switched on, no display shown.	Please check the main plugs and sockets. If necessary; contact your service agents or the manufacturer.
No acoustic signal is heard. (End of therapy...)	Please check the settings of the acoustic tones in the menu. The volume must be greater than 0.
The following message appears: The desired suction intensity cannot be reached. Please clean or replace filter.	Clean or replace the external air filter on the back side of the device.

Please contact your service agent or the manufacturer if the problems cannot be solved by the measures mentioned above.

Please note that the unit must be placed on a plane horizontal surface. The device should neither be placed in front of radiators nor should it be covered with pillows or blankets while in operation. Do not cover the ventilation slots on the bottom of the unit either.

13 Appendix

Comments According to the Medical Device Directive

The Stereodylator is a mains operated current stimulation device of protection class I.

The device is in accordance with the EC Medical Device Directive (93/42/EEC) and therefore carries the CE-sign with the number of the "notified body for medical devices". The respective graphical symbol is placed on the type plate.

According to the Medical Device Directive, Stereodylator is a device of class IIa.

The manufacturer is only responsible for the security, operational reliability and functionality of the device if:

- the device is used in accordance with the instructions for use;
- the electrical installation of the location where the device will be used corresponds to the respective current requirements of electrical safety;
- the device is not used in hazardous environments and humid locations;
- the mountings, add ons, internal adjustments, modifications or repairs are realized only by personnel authorized by the manufacturer;
- the operator regulation of this EC-directive is observed within the scope of the Medical Device Directive.

You may obtain technical support by the manufacturer or the dealers or service authorized by the manufacturer. The manufacturer projects a product life of 10 years.

Stereodylator is an electronic device. Disposal has to be done according to regulations for electronic devices. Consumables have to be disposed as residual waste.

On request, the manufacturer will provide you with further technical descriptions for all serviceable parts of the device, such as circuit diagrams, spare part lists and adjustment instructions as far as these are of use for the qualified technical staff of the user.

Comments on electromagnetic compatibility (EMC)

Medical, electrical devices are subject to special precautions concerning the EMC. They must be installed and operated according to the EMC-advice given in the accompanying documents. In particular medical, electrical devices may be influenced by portable and mobile RF-communication devices.

The manufacturer guarantees the conformity of the unit with the EMC-requirements only when using accessories which are listed in the EC declaration of conformity. The usage of other accessories may cause an increased emission of electromagnetic disturbances or may lead to a reduced electromagnetic immunity.

The unit must not be arranged physically close to other devices or stacked with them. If such an order is necessary nevertheless, the unit must be observed in order to check it for the intentional operation.

You find more EMC-comments in the chapter "Warnings and Safety Precautions" of this manual as well as in the Technical Information on the next two pages.

In accordance with the EMC-regulations for medical products we are obliged by law to provide the following information.

Guidance and manufacturer's declaration — electromagnetic emissions

The equipment is intended for the use in the electromagnetic environment specified below. The customer or the user of the equipment should ensure the compliance with this stipulation		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions, CISPR 11	Group 1	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions, CISPR 11	Class B	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions, IEC 61000-3-2 (*)	Class A	
Voltage fluctuation/flicker emissions, IEC 61000-3-3 (*)	complies	

(*) Note: For devices with a power consumption between 75 W and 1000 W only.

Guidance and manufacturer's declaration — electromagnetic immunity

The equipment is intended for the use in the electromagnetic environment specified below. The customer or the user of the equipment should ensure the compliance with this stipulation.			
Immunity test	IEC 60601- test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD), IEC 61000-4-2	±6 kV contact	±6 kV contact	Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
	±8 kV air	±8 kV air	
Electrical fast transient/burst, IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
	±1 kV for input/output lines	±1 kV for input/output lines	
Surge, IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
	±2 kV common mode	±2 kV common mode	
Voltage dips, short interruptions and voltage variations on power supply input lines, IEC 61000-4-11	<5% U _n for ½ cycle (>95% dip)	<5% U _n for ½ cycle (>95% dip)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment requires continued operation during power mains interruptions, it is recommended that the equipment be powered from an uninterruptible power supply or a battery.
	40% U _n for 5 cycles 60% dip)	40% U _n for 5 cycles 60% dip)	
	70% U _n for 25 cycles 30% dip)	70% U _n for 25 cycles 30% dip)	
	<95% U _n for 5 s (>5% dip)	<95% U _n for 5 s (>5% dip)	
Power frequency (50/60 Hz) magnetic field, IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U_n is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration — electromagnetic immunity

The equipment is intended for the use in the electromagnetic environment specified below. The customer or the user of the equipment should ensure the compliance with this stipulation.

Immunity test	IEC 60601- test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF, IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{air}	d=1.2√P
Radiated RF, IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3V/m	d=1.2√P for 80 MHz to 800 MHz d=2.3√P for 800 MHz to 2.5 GHz
			Where P is the maximum output power rating of the transmitter in watts according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Interference may occur in the vicinity of equipment marked with the following symbol: 

Recommended separation distances to portable and mobile RF communication equipment

The equipment is intended to be operated in an electromagnetic environment where radiated RF interference is controlled. The user can help to avoid interferences by meeting the minimum separation distances between portable and mobile RF communication equipment (transmitters) according to the maximum output power of the communication equipment.

Rated power of the transmitter (W)	Separation distance according to the transmission frequency (m)		
	150 kHz to 80 MHz d=1.2√P	80 MHz to 800 MHz d=1.2√P	800 MHz to 2.5 GHz d=2.3√P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

14 Index

<p>A</p> <p>Accessories 33</p> <p>B</p> <p>biphase electrode technics 15</p> <p>C</p> <p>Cleaning 22 Connecting 11 Contraindications for use of current stimulation 27 Contraindications for use of ultrasound 25 Control Panel 7 Current types 30</p> <p>D</p> <p>direct help 7 Direct help menu 9 disinfection 22</p> <p>E</p> <p>electrode handpiece 19 Electrode Positioning 14 Electrodes 14, 22, 23</p> <p>I</p> <p>Intended Use 5 Intensity regulator 10</p>	<p>M</p> <p>Maintenance 21 monophase electrode technics 15</p> <p>O</p> <p>Operating Concept 7</p> <p>P</p> <p>precautions 27</p> <p>S</p> <p>safety controls 21 Signs 28 single-pole electrodes 16 suction electrodes 18 Switch-On 11</p> <p>T</p> <p>Technical Data 29 Therapy 13 Three-pole electrodes 17 Touch screen 7 Transport and Assembly 11 Troubleshooting 34</p> <p>U</p> <p>ultrasound module 20</p> <p>W</p> <p>warnings 11, 24</p>
---	--

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