QuikRead CRP QuikRead U-ALB QuikRead FOB QuikRead Strep A

QuikRead[®] CRP

Test before prescribing

- **RELIABLE:** as accurate as a laboratory test
- **QUICK:** result available exactly when needed
- **EASY TO USE:** simple to operate – even by non-laboratory personnel

ORION
DIAGNOSTICA

QuikRead[®] CRP

QuikRead CRP is a simple test for quantitative measurement of C-reactive protein (CRP).

CRP is normally present at very low concentrations in the blood of healthy people. CRP concentrations are markedly increased in bacterial infections, whereas viral infections mostly induce a very modest elevation or none at all. In addition to clinical signs, measurement of CRP can therefore be used as a basic tool to distinguish bacterial infections from viral infections. Monitoring CRP levels also provides an objective means of assessing treatment response, as CRP levels fall rapidly as a result of effective therapy.

When the test is performed near the patient, the CRP result will be available during patient consultation and can effectively guide antibiotic use. The QuikRead system – consisting of a small instrument and a ready-to-use kit – is especially designed for use in primary health care settings being instantly ready for operation and giving immediate test results.

Features and benefits of QuikRead CRP

RELIABLE: as accurate as a laboratory test

- quantitative result
- reproducible CRP result within the 8–160 mg/l range

QUICK: result available exactly when needed

- measuring time is less than 1 minute


EASE OF USE: simple to operate – even by non-laboratory personnel

- performed on a finger prick blood sample (alternatively serum or plasma sample)
- no pipetting of reagents, built-in calibration, instrument robust and simple to operate

For more information about CRP you will find out in “The value of C-reactive protein in primary health care”.

Product	Cat. No.
QuikRead 101 Instrument	06040
• package including adapter, accuracy dispenser, manual and dust cover	
• can be operated on 8 x AA (1.5) batteries or an 100–240 AC power supply	
QuikRead CRP Kit, Prefilled cuvettes, 50 tests (with capillaries and plungers)	134191
(without capillaries and plungers)	134192
QuikRead CRP Kit, 50 tests (with capillaries and plungers)	67961
(without capillaries and plungers)	68789
QuikRead Dispenser, 1 ml	05560
QuikRead CRP Control, 1 ml	68296
• recommended for quality control	
QuikRead Capillaries, 50 pcs	67962
QuikRead Plungers, 50 pcs	67966

QuikRead is a registered trademark of Orion Diagnostica Oy.

Orion Diagnostica Oy, P.O. Box 83, FI-02101 Espoo, Finland
Tel. +358 10 4261, Fax +358 10 426 2794, www.oriondiagnostica.com, www.quikread.com

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Intended use
Orion Diagnostica QuikRead CRP Control is intended for quality control of CRP assays by the QuikRead® Instrument or QuikRead® 101 Instrument or QuikRead go® Instrument. For In vitro diagnostic use.

Table with 2 columns: QuikRead QuikRead 101, QuikRead go. It details the use of the control in the same way as unknown whole blood specimen would be used.

Reagents
The control is human-based material. It should be handled in the same way as patient sample.
Preparation
QuikRead CRP Control is ready-to-use. Allow the material to stabilize for 30 minutes at room temperature (18, 25°C). Mix well. Do not shake.
Storage
Stable until the expiry date when stored at 2...8°C. Opened vial. Stable for 2 months at 2...8°C. Do not use any vial that has visible evidence of microbial growth.

Expected values
The individual result should be within the range indicated on the control vial label. Each laboratory should establish its own mean and precision data.

Warnings and precautions
Reagents contain <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
Do not use the product after the expiry date marked on the outer package.
Do not exceed the stability periods for opened reagents.
QuikRead CRP Control contains material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HBsAg). As a precaution, handle QuikRead CRP Control as if capable of transmitting infectious agents.

Disposal
Dispose of contents according to national and local law.
Materials of the components
Paper: Instructions for use
Cardboard: Outer box
Glass: Control vials
Plastic: Vial caps
Rubber: Stopper
When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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