pRESET Thrombectomy Device

phenox

INSTRUCTIONS FOR USE

2	DEUTSCH GEBRAUCHSANWEISUNG	20	DANSK BRUGSANVISNING	38	БЪЛГАРСКИ ИНСТРУКЦИИ ЗА УПОТРЕБА
4	ENGLISH INSTRUCTIONS FOR USE	22	NORSK BRUKSANVISNING	40	EESTI KASUTUSJUHEND
6	FRANÇAIS MODE D'EMPLOI	24	POLSKI INSTRUKCJA UŻYTKOWANIA	42	HRVATSKI UPUTE ZA UPORABU
8	ITALIANO ISTRUZIONI PER L'USO	26	SLOVENŠČINA NAVODILA ZA UPORABO	44	SLOVENČINA NÁVOD NA POUŽITIE
10	ESPAÑOL INSTRUCCIONES DE USO	28	ČEŠTINA NÁVOD K POUŽITÍ	46	LIETUVIŲ K. NAUDOJIMO INSTRUKCIJA
12	PORTUGUÊS INSTRUÇÕES DE UTILIZAÇÃO	30	ROMÂNĂ INSTRUCȚIUNI DE UTILIZARE	48	TÜRKÇE KULLANIM TALİMATI
14	NEDERLANDS GEBRUIKSAANWIJZING	32	MAGYAR HASZNÁLATI UTASÍTÁS	50	РУССКИЙ ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
16	SVENSKA BRUKSANVISNING	34	ΕΛΛΗΝΙΚΑ ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	52	اللغة العربية إرشادات الاستخدام
18	SUOMI KÄYTTÖOHJEET	36	LATVISKI LIETOŠANAS PAMĀCĪBA		

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Package content

1 x **pRESET** Thrombectomy Device
(of version **pRESET STANDARD**, **pRESET LITE** and **pRESET LUX**)


Device description

The **pRESET** (STANDARD and LITE) consists of a self-expanding Nitinol structure ①, carries one x-ray visible marker on its proximal ② and two on its distal end ③ and is firmly attached to an insertion wire ④.

The **pRESET LUX** consists of a self-expanding radiopaque Nitinol structure ①, carries one x-ray visible marker each on its proximal ② and its distal end ③ and is firmly attached to an insertion wire ④.

The device is stored in compressed form in an introducer sheath (not displayed).

The **pRESET** (all versions) is introduced into the target vessel through a suitable microcatheter and deployed inside the thrombus. After complete deployment the slow withdrawal of the instrument occurs under continuous aspiration via the guiding catheter or aspiration catheter.

Intended Purpose

The intended purpose of **pRESET** Thrombectomy Devices is acute treatment of ischemic stroke by removal of thrombi from intracranial vessels e.g., the internal carotid artery, middle cerebral artery, or basilar artery or from vessels situated distally of them, with a temporary contact duration less than 60 minutes.

Indications

The **pRESET** Thrombectomy Device is designed for mechanical clot retrieval from intracranial arteries as acute ischemic stroke treatment

- for patients who are ineligible for intravenous thrombolysis or
- for patients who failed thrombolysis therapy and
- as a supplement treatment of an initiated thrombolysis therapy.

Contraindications

The use of **pRESET** is contraindicated in patients with

- preexisting hemorrhage in the brain tissue supplied by the target vessel and
- established infarction in the brain tissue supplied by the target vessel.

Compatibility

The following **pRESET** devices and microcatheters are compatible and are recommended for use in minimum vessel diameters:

pRESET Device	Insertion Wire Length	Microcatheter	Min. Vessel Ø	Relaxed Ø x Usable Length
PRE-LT-3-20	2 m	ID 0.0165"/0.017": - ECHELON 10/14 - EXCELSIOR SL 10 - HEADWAY 17	1,5 mm	3 x 20 mm
PRE-LT-4-20	2 m		1,5 mm	4 x 20 mm
PRE-LUX-4-20	2 m	ID 0.021": - TREVO PRO 18 - PROWLER SELECT PLUS - VIA 21 - PHENOM 21	2 mm	4 x 20 mm

PRE-4-20	2 m	ID 0.021": - TREVO PRO 18 - PROWLER SELECT PLUS - REBAR 18 - VIA 21 - PHENOM 21 - RAPIDTRANSIT	2 mm	4 x 20 mm
PRES-5-40	2 m		2 mm	5 x 40 mm
PRE-6-30	2 m		3 mm	6 x 30 mm
PRES-6-50	2 m	ID 0,021": - TREVO PRO 18 - PROWLER SELECT PLUS - REBAR 18 - VIA 21 - PHENOM 21	3 mm	6 x 50 mm

Recommended Procedure

Preparation of procedure and patient

- Record and document any available information about the date, time and nature of the onset of clinical stroke symptoms.
- Record and document data on past and recent medical history, including medication.
- As far as possible, inform the patient and document the patient's consent to the planned intervention, pointing out the possible complications and potential consequences (see Complications). In cases where patients cannot give consent themselves, their relatives should, as far as possible, be asked if they know what the presumed patient's will would be. Otherwise, in case of emergency, the rules of emergency care for incapacitated patients apply subject to different institutional or national requirements.
- A native CT scan of the brain is required to rule out intracranial hemorrhage or large hypodense ischemic brain lesions. In that case the use of **pRESET** is contraindicated.
- Advise the anesthesiologist to maintain appropriate levels of systemic blood pressure. Hypotension is to be strictly avoided.**
Induce general anesthesia with neuromuscular relaxation from the very beginning of the angiographic diagnosis and eventual endovascular treatment. Complete hemodynamic monitoring is recommended, but must not delay the endovascular procedure.
- Prepare both groins, introduce 8F sheath preferably in the right femoral artery.
- For vessel diagnostic conduct angiography of the complete anterior and posterior circulation. Take magnified views of the occluded vessel(s), and if necessary oblique views.
- Define the occluded vessel(s) as target(s) for further therapy.
- Use the largest possible guiding catheter (VA: 6F, ICA: 8F). In the case of significant elongation or stenosis of the access vessel, consider the implantation of a suitable stent. In the case of stenting the administration of suitable inhibitors of platelet aggregation is to be followed.
According to user reports, an 8F guiding catheter with balloon occlusion and a flexible 5 to 6F catheter for aspiration to get closer to the thrombus site for aspiration can be useful for thrombectomy.
- It is important that the device is introduced in sufficiently large target vessels only.** Heed the minimum vessel diameter and the catheter compatibility.

Introduction of microcatheter

- Never catheterize against any resistance!**
Introduce a suitable microcatheter (see Compatibility table) into the target vessel, with careful passage beyond the thrombus. The application of the so-called "road map" technique is recommended.
- Withdraw the microguidewire from microcatheter. The careful injection of a small amount (e.g., 0.2 ml) of contrast medium may be used to verify the proper position of the tip of the microcatheter beyond the thrombus and in the vessel lumen.


Preparation and introduction of pRESET

- In case of major resistance during insertion of pRESET, remove the device and eventually the microcatheter and catheterize again.**
With a rotating hemostatic valve (RHV) and slow but constant saline flushing connected to the microcatheter, **pRESET** will be introduced into the microcatheter. During this process, the RHV will be opened, the introducer sheath with the device inside will pass the valve, the valve will gently be closed and the introducer sheath will be flushed back until water exits from the proximal end of the introducer sheath.
After the introducer sheath with the device inside is flushed, the introducer sheath will be advanced until close contact with the hub of the microcatheter. The introducer sheath will be firmly kept in this position; **pRESET** is advanced from the introducer sheath into the microcatheter by gentle pushing of the insertion wire. The process of device advancement will be continued until about sixty (60) cm of the insertion wire are pushed forward. Then the introducer sheath will be removed.
The device is pushed further forward until the marker ("FLUOROSAFE MARKER": distance to distal device tip approx. 125 cm in compressed condition) of the insertion wire reaches the entry of the hemostatic valve. This procedure does not need fluoroscopy because the marker ("FLUOROSAFE MARKER") identifies the position to which the system can be advanced without the device tip leaving the microcatheter.
- Never push pRESET beyond the distal tip of the microcatheter! To do so may result in dissection or perforation of the concerning artery.**
Under continuous fluoroscopy the device will be advanced until the tip of the microcatheter will be reached by the distal marker of the device.

Deployment of pRESET

- Deploy the device completely by gentle and slow withdrawal of the microcatheter until the microcatheter tip is located sufficiently proximal. This is achieved if the distal microcatheter marker is placed proximal to the proximal x-ray visible **pRESET** marker.

During this process, use continuous fluoroscopy to keep the device in its original position. The device is ideally deployed in the thrombus. Is the **pRESET** longer than the thrombus, the instrument should protrude distally. The illustration right hand depicts the radiographic image of a **pRESET STANDARD & LITE** (top) and **pRESET LUX** (down).

16. Wait five (5) minutes before starting the withdrawal. The saline flush will be removed from the guiding catheter RHV and will be replaced by a fifty (50) ml syringe with Luer Lock. Two (2) additional syringes of the same size will be prepared and available. Take care for firmly closed valves to avoid aspiration of air.

17. If applicable, inflate the balloon of the guide catheter acc. to its IFU.

Withdrawal of **pRESET** under aspiration

18. **Some resistance may be observed at this point; the more distally the device has been introduced, the more resistance might be encountered. Always pull gently, never abrupt or in steps!**

During continuous aspiration on the designated catheter, the unit microcatheter + **pRESET** will **SLOWLY** be pulled back.

For the efficacy of the device in the posterior circulation, the condition of the vertebral artery contralateral to the vessel catheterized with the guiding catheter is important: If this artery is not supplying the basilar artery, no further action is required. If this artery is fully patent, aspiration on the guiding catheter will only affect the basilar artery when the non-catheterized basilar artery is temporarily blocked. This may be achieved preferably with a shape adaptive balloon catheter or by probing the A. vertebralis with a diagnostic catheter having an occluding effect.

19. Follow under fluoroscopy as the device is travelling through the vessel under aspiration. As soon as the **pRESET** has entered the aspirating catheter, withdrawal can be accelerated. Pull the device in front of the RHV. Disconnect the RHV and the **pRESET** from the aspirating catheter. The aspiration started slowly with the pulling back of the device, was increased from start to finish and has yielded one hundred to one hundred and fifty (100-150) ml of blood up to this point. Aspirate for safety reasons another fifty (50) ml of blood. As needed close the aspirating catheter, e.g., using a syringe. If applicable, deflate the balloon of the guiding catheter acc. to its IFU

20. Gently inject approx. ten (10) ml of contrast medium into the guiding catheter to see if the target vessel was recanalized by use of **pRESET**.
 - a) If the target vessel is patent, decide if further vessels are occluded and need recanalization. If this is not the case, withdraw the guiding catheter, wait for ten to fifteen (10-15) minutes and re-inject the concerning supra-aortic vessel to confirm permanent patency.
 - b) If the target vessel is *not* patent, decide upon further action, which may include a repeated **pRESET** procedure, use of another mechanical device or fibrinolysis.

Reuse of **pRESET** in same patient

21. In case of a repeated or an additional procedure an intact device can be reused in the same patient up to two (2) times.
 - 1) Removal of **pRESET** from the RHV: Open the hemostatic valve (RHV) of the aspirating catheter **COMPLETELY** before pulling out the device.
 - 2) Cleaning: Gently rinse **pRESET** using sterile water or saline solution; blood or thrombus must not adhere to the device.
 - 3) Inspection: Check the device to exclude any damage.
 - 4) Preparation for new procedure: Push the proximal end of the insertion wire into the corresponding introducer sheath. Then move the introducer sheath up to the device. Then pull the device **SLOWLY** and completely back in the introducer sheath. If the device cannot be pulled into the sheath it has to be discarded.
 - 5) The reapplication of the **pRESET** follows steps 10 to 20 of this instructions for use.

Precautions

- Do not use **pRESET** for other purposes than the intended one.
- Do not deploy the device in an artery smaller than the corresponding minimum vessel diameter (Ø) according to the table given in the Compatibility table to avoid vessel damage.
- Use the device only in microcatheters listed in the Compatibility table. Microcatheter or other catheters with larger inner diameters may lead to a stressful passage for **pRESET** device and to a disconnection of the device or parts of it.
- Apply the device only under continuous fluoroscopic monitoring.
- Copious flushing of the introducer sheath is required to remove potentially trapped air.
- Pay attention during device passage from introducer sheath into microcatheter that the introducer sheath does not move back proximally. If this happens the instrument and/or the insertion wire may be kinked and/or weakened and therefore might be disconnected.
- Pushing the device beyond the tip of the microcatheter may dissect or perforate the target vessel.
- Pulling or pushing on the insertion wire with inadequate force and torsion of the insertion wire may disconnect **pRESET** from the insertion wire. In such cases the retrieval by a snare device is recommended.
- In case of excessive retraction resistance the reintroduction of the microcatheter for splinting of the **pRESET** insertion wire or removal of the device is recommended. If the microcatheter cannot be pushed completely over the device remove both, device and microcatheter as a unit.
- If vasospasm is suspected in the affected vessel region, all necessary counteractive measures (e.g. medication treatment) should be taken before **pRESET** treatment.
- **pRESET** can be used up to three (3) times in the same patient. The contact with the patient of all applications must in sum be less than 60 minutes.
- Between the applications the device must be cleansed of blood and thrombus and inspected for damages. If the device or the introducer sheath is damaged or the device cannot be pulled into the introducer sheath it has to be discarded!
- **pRESET** is a delicate device and demands cautious treatment. Do not advance or withdraw the device against resistance. Do not rotate the insertion wire. If necessary, remove **pRESET** via microcatheter.
- The device must only be used by specialised and trained physicians.
- The deployed **pRESET** must not be retracted through implanted vascular implants (e.g., stents or flow diverters)!

General Information

- Keep away from heat. Store in a cool, dry place.
- Use before date of expiry, as sterility is not guaranteed otherwise..
- Do not use damaged packages to ensure product integrity and sterility.
- The device must be checked for damage before use. Do not use deformed or damaged devices, as function cannot be assumed otherwise.
- The device is for single use only. The device must not be resterilized or recycled for other patients because it cannot be reliably cleaned by the user.
- **pRESET**, retrieved clots, blood and if applicable its packaging must be disposed of appropriately in marked containers.

Complications

The following and further complications may occur during or despite the use of the **pRESET**:

- Air embolism, embolism in distal vessels, vessel occlusion, thrombosis and cerebral ischemia
- Perforation, rupture, dissection and other arterial lesions
- Vasospasm, formation of a pseudoaneurysm, intracranial hemorrhage
- Allergic reaction, infection
- Space-occupying infarction, neurologic deficit including all of the consequences of a stroke
- Persistent vegetative state, death

Symbols and their meanings

Symbols given on label:

	Caution
	Consult instructions for use
	Do not reuse
	Do not resterilize
	Batch code
	Sterilized using ethylene oxide
	Do not use if package is damaged
	Non-pyrogenic
	Use-by date
	Catalogue number
	Content
	Keep dry. Keep away from sunlight.
	This product was brought onto the market in accordance with Directive 93/42/EEC concerning medical devices.
	Manufacturer
	Material: NiTi (Nickel Titanium alloy)
	Latex-free
	Phthalate-free
	Medical device
FLUOROSAFE MARKER	Marker (black) on insertion wire
RADIOPAQUE	Radiopaque device structure

Limits of Liability

phenox GmbH shall not be liable for damage caused by use other than that intended or reuse of the product.

phenox and **pRESET** are registered trademarks of phenox GmbH in the Federal Republic of Germany and other states.

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phenox GmbH | Lise-Meitner-Allee 31 | D-44801 Bochum, Germany