ACL/PCL Fixation Systems

BIOSURE® HA and BIOSURE PK Interference Screws

BIOSURE Interference Screw

The BIOSURE family of interference screws provides an all-purpose range of screws and instrumentation. BIOSURE Interference Screws incorporate the latest proven technologies to help you deliver exceptional graft fixation. Available in 24 sizes, including 6 mm screws for double bundle repairs and reverse thread screws for right knees.

Absorbable

Our absorbable BIOSURE HA interference screw combines poly-L-lactic acid (PLLA) with hydroxyapatite (HA), a calcium phosphate similar to the mineral compound of bone. Once absorbed, it has the fracture resistance of human corticle bone. Unlike other materials where HA has been derived from animal or coralline sources, our material is completely synthetic and does not exhibit natural variability or immunological hazards.

Non-absorbable

Our non-absorbable BIOSURE PK interference screw incorporates PEEK-OPTIMA® polymer from Invibio®. With a modulus similar to cortical bone, PEEK-OPTIMA material exhibits an ideal combination of strength, stiffness and toughness, together with biocompatibility. BIOSURE PK interference screws offer the revisability of an absorbable screw and inherent radiolucent characteristics from a material that has gained confidence and acceptance in the medical community for 10 years.

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
BIORCI®-HA Screws

BIORCI-HA Screw System Implants

The BIORCI-HA Screw System has proven effectiveness for graft fixation and site repair. BIORCI-HA screws use standard BIORCI instruments and offer 23 screw sizes from which to choose.

Reference #	Description
7207686	7/7 x 20 mm BIORCI-HA
7207674	7/7 x 25 mm BIORCI-HA
7207677	7/7 x 30 mm BIORCI-HA
7207675	7/8 x 25 mm BIORCI-HA
7209023	7/8 x 25 mm BIORCI-HA reverse thread
7207676	7/8 x 30 mm BIORCI-HA
7209012	8 x 20 mm BIORCI-HA
7207678	8 x 25 mm BIORCI-HA
7209021	8 x 25 mm BIORCI-HA reverse thread
7207679	8 x 30 mm BIORCI-HA
7207680	8 x 35 mm BIORCI-HA
7209013	9 x 20 mm BIORCI-HA
7207681	9 x 25 mm BIORCI-HA
7209020	9 x 25 mm BIORCI-HA reverse thread
7207682	9 x 30 mm BIORCI-HA
7207683	9 x 35 mm BIORCI-HA
7209016	10 x 25 mm BIORCI-HA
7209019	10 x 25 mm BIORCI-HA reverse thread
7209017	10 x 30 mm BIORCI-HA
7209018	10 x 35 mm BIORCI-HA
7209014	11 x 30 mm BIORCI-HA
7209022	11 x 35 mm BIORCI-HA
7209015	12 x 35 mm BIORCI-HA
7207636	1.5 mm x 9" (23 cm) Guide Wire, non-sterile (6 per package)
014859	1.5 mm x 12" (30 cm) Guide Wire, sterile (1 per package)

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
Ordering Information


From bioabsorbable screws to powered instrumentation, Smith & Nephew offers everything you need for consistently strong cruciate reconstructions.

BIORCI Screw System Implants

1.5 mm cannulation, sterile, 1 per box

7207554	Screw, BIORCI 7 mm (7 mm Head) x 25 mm	7207564	Screw, BIORCI 9 x 30 mm
7207555	Screw, BIORCI 7 mm (7 mm Head) x 30 mm	7207565	Screw, BIORCI 9 x 35 mm
7207557	Screw, BIORCI 7 mm (8 mm Head) x 25 mm	7207567	Screw, BIORCI 10 x 30 mm
7207558	Screw, BIORCI 7mm (8 mm Head) x 30 mm	7207568	Screw, BIORCI 10 x 35 mm
7207560	Screw, BIORCI 8 x 25 mm	7207607	Screw, BIORCI 7mm (8 mm Head) x 25 mm, Reverse Thread
7207561	Screw, BIORCI 8 x 30 mm	7207608	Screw, BIORCI 8 x 25 mm, Reverse Thread
7207562	Screw, BIORCI 8 x 35 mm		
7207563	Screw, BIORCI 9 x 25 mm		

BIORCI Absorbable Screws offer optimal strength and selection for soft tissue and B-T-B reconstructions.

Pullout Strength Chart comparing BIORCI strength of fixation vs. competitors. BIORCI bioabsorbable screw is 28% more effective at maintaining fixation than the closest competitor.**

Comparison chart showing a 58% increase in tibial pullout strength utilizing longer length screws.**

BIORCI-HA Screws

1.5 mm cannulation, sterile, 1 per box

7207686	Screw, BIORCI-HA 7/7 x 20 mm	7209018	Screw, BIORCI-HA 10 x 35 mm
7207674	Screw, BIORCI-HA 7/7 x 25 mm	7209014	Screw, BIORCI-HA 11 x 30 mm
7207677	Screw, BIORCI-HA 7/7 x 30 mm	7209022	Screw, BIORCI-HA 11 x 35 mm
7207675	Screw, BIORCI-HA 7/8 x 25 mm	7209015	Screw, BIORCI-HA 12 x 35 mm
7207676	Screw, BIORCI-HA 7/8 x 30 mm	7209023	Screw, BIORCI-HA 7/8 x 25 mm Reverse Thread
7209012	Screw, BIORCI-HA 8 x 20 mm	7209021	Screw, BIORCI-HA 8 x 25 mm Reverse Thread
7207678	Screw, BIORCI-HA 8 x 25 mm	7209020	Screw, BIORCI-HA 9 x 25 mm Reverse Thread
7207679	Screw, BIORCI-HA 8 x 30 mm	7209019	Screw, BIORCI-HA 10 x 25 mm Reverse Thread
7207680	Screw, BIORCI-HA 8 x 35 mm	7207636	1.5 mm x 9" (23 cm) Guide W
7209013	Screw, BIORCI-HA 9 x 20 mm	014859	1.5 mm x 12" (30 cm) Guide
7207681	Screw, BIORCI-HA 9 x 25 mm	3134033	1.5 mm x 17" (30 cm) Guide
7207682	Screw, BIORCI-HA 9 x 30 mm		
7207683	Screw, BIORCI-HA 9 x 35 mm		
7209016	Screw, BIORCI-HA 10 x 25 mm		
7209017	Screw, BIORCI-HA 10 x 30 mm		

BIORCI-HA Screw System Instrumentation

7208661	7 mm BIORCI Tap	7207359	BIORCI Starter
7208662	8 mm BIORCI Tap	7207369	BIORCI Tray
7208663	9 mm BIORCI Tap	7207766	Wingless Compaction Router
7208664	10 mm BIORCI Tap	3134186	Single Flute Endoscopic Drill Bit
7208656	BIORCI Torx® Driver		

Kopija tiktai

2015-11-30
Vadybininkas
Evaldas Ročys

¹ C.M.Agrawal, M. Fan, C. Zhu and K.A. Athanasiou: A New Technique to Control the pH in the Vicinity of Bioabsorbable Implants.

* In vitro data on file with Smith & Nephew. ** Data on file with Smith & Nephew.

Complete BIORCI Screw Systems can be tailored to include guides, manual and powered instrumentation. Ask your local sales representative for assistance.

Endoscopy
Smith & Nephew, Inc.
Andover MA 01810
USA

www.smith-nephew.com
T +1 978 749 1000
F +1 978 749 1108
US Customer Service +1 800 343 5717

BIORCI is a trademark of Smith & Nephew, Inc. Reg US Pat & TM Off.
Torx is a registered trademark of Textron, Inc.
©2001, 2002, 2007 Smith & Nephew, Inc. All rights reserved. Printed in USA.

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2015-11-30

Kopija tikra

Developed by Smith & Nephew, BIORCI® Bioabsorbable Interference Screws are among the strongest bioabsorbable screws available for arthroscopists performing soft tissue or bone-tendon-bone (B-T-B) reconstructions. With one of the highest pull-out strengths and the greatest resistance to breakage among bioabsorbable screws, BIORCI® Screws are designed to hold the graft in place for secure anatomic aperture reconstruction. BIORCI-HA® Screws go a step further, combining the proven strength and integrity of the BIORCI Screw System with hydroxylapatite (HA).

BIORCI® and BIORCI®-HA

Enhanced Bioabsorbable Screws

BIORCI and BIORCI-HA bioabsorbable screws are a complete line of bioabsorbable implants allowing for multiple approaches of repair and the flexibility to accommodate a variety of anatomical placements.

Longer thread engagement provides greater pullout strength within the tibial tunnel.

BIORCI bioabsorbable screw is the only system that includes reverse-threaded screws to keep the soft tissue grafts in optimal position.

Unique and patented Torx driver-to-screw interface disperses force throughout the screw during insertion.

Smith & Nephew is the developer of the first soft tissue fixation screw. We offer the broadest choice in technique-specific screws today, with quality and reliability second to none.

A safe, proven alternative to metal screws, BIORCI-HA screws are strong when needed, then absorb discreetly.

Research indicates the pH levels surrounding the absorption of polymers is reduced with the combined absorption of hydroxylapatite!


BIORCI-HA screws offer 20 mm and 25 mm lengths for use with bone-tendon-bone grafts.

The 20 mm and 25 mm length screws are sized for optimal bone-tendon-bone graft fixation. Soft threads minimize gra laceration, and range of diameters ens a tight fit in all tunnels.

Instrumentation

The bone-tendon-bone approach utilizes an endoscopic cannulated drill bit for femoral tunnel preparation and a cannulated drill bit for drilling the tibial tunnel.

Unique Torx® interface helps minimize breakage


Laser markings ensure that implants are fully seated

Vadybininkas
Egidius Ročys

Kopija tikra

15-11-30

7

BIORCI screws are the only implants that offer a unique oversized head incorporated in their femoral design.

BIORCI screws are optimized for surgeons utilizing the anteromedial approach.

The enlarged head provides extra contact area at the anatomic site. The head also reduces any stress risers caused by "point loading" that may occur on other full-threaded screws.

Aperture fixation using the BIORCI screws results in a short graft construct between the femoral and tibial tunnels.

Instrumentation

To optimize the preparation of the femoral tunnel during drilling, select RCI winged routers. These patented routers compact the tunnel while drilling and provide a countersink for the enlarged head.

Maximize your length of graft/tunnel engagement for greater pullout strength.

BIORCI tibial-specific screws lead the industry in addressing tibial fixation issues.

These longer length screws provide more soft threads to ensure maximum length interference between the graft and tibial tunnel.

Designed in conjunction with leading surgeons to maximize the fixation along the tibial tunnel.

Instrumentation

The tibial tunnel can be prepared with either the RCI tibial compaction routers or a standard cannulated drill bit.

BIORCI screws with a rounded head and soft threads for maximum protection of the grafts.

The latest addition to the RCI screw line.

Femoral screws with a matched rounded head. These screws incorporate the same Torx® interface as all BIORCI screws to minimize breakage when your approach is through the tibial tunnel.

The 1.5 mm guide wire is the largest in the industry to assist you in accurate parallel placement of the screws.

Instrumentation

The transtibial approach utilizes an endoscopic cannulated drill bit for femoral tunnel preparation and a cannulated drill bit for drilling the tibial tunnel.

Reverse Thread – because your patients do not have two left knees.

Use the patented Reverse Thread BIORCI bioabsorbable screw for your right knees.

To produce mirror images on right and left limbs, the final femoral graft position within the tunnel must be the same.

Standard (clockwise) screw thread orientation produces clockwise rotation of the graft in left knees, resulting in a posterior final graft position. However, standard (clockwise) screw thread orientation results in an anterior final femoral graft position within the femoral tunnel in right knees.


Therefore, to achieve a posterior final graft position in a right knee, to mirror the rotation of the left knee, a reverse (counterclockwise) threaded screw may be selected.

Instrumentation

To optimize the preparation of the femoral tunnel during drilling, select the RCI winged routers.


RCI FEMORAL ROUTER

WINGLESS COMPACTION ROUTER

RCI TIBIAL ROUTER

CANNULATED DRILL BIT

ENDOSCOPIC CANNULATED DRILL BIT

CANNULATED DRILL BIT

Vadybininkas
Evaldas Ročus


To maximize product efficiency it is recommended that the surgical site is prepared with the BIORCI Starter or Tap prior to insertion of the absorbable implant.

ACL/PCL Fixation Systems

ENDOBUTTON® Family of Fixation Devices

ENDOBUTTON CL Ultra Fixation Device

Delivering improved strength and stiffness, the ENDOBUTTON CL ULTRA Fixation device offers the strongest soft tissue femoral fixation currently available. Ideal for anatomic or traditional ACL reconstruction techniques. Preloaded with ULTRABRAID® #5 lead suture (white) and DURABRAID® #5 flipping suture (green) for added procedure efficiency.

Reference #	Description
72203331	10 mm ENDOBUTTON CL ULTRA
72200146	15 mm ENDOBUTTON CL ULTRA
72200147	20 mm ENDOBUTTON CL ULTRA
72200148	25 mm ENDOBUTTON CL ULTRA
72200149	30 mm ENDOBUTTON CL ULTRA
72200150	35 mm ENDOBUTTON CL ULTRA
72200151	40 mm ENDOBUTTON CL ULTRA
72200152	45 mm ENDOBUTTON CL ULTRA
72200153	50 mm ENDOBUTTON CL ULTRA
72200154	55 mm ENDOBUTTON CL ULTRA
72200155	60 mm ENDOBUTTON CL ULTRA

ENDOBUTTON DIRECT Fixation Device

The titanium ENDOBUTTON device maximizes the amount of graft in the tunnel, an ideal solution for accommodating shorter femoral tunnels that result from more anatomic placement.

Reference #	Description
72201206	5 mm ENDOBUTTON DIRECT
72202980	5 mm ENDOBUTTON DIRECT with suture
72201207	6 mm ENDOBUTTON DIRECT
72202981	6 mm ENDOBUTTON DIRECT with suture
72201208	7 mm ENDOBUTTON DIRECT
72202982	7 mm ENDOBUTTON DIRECT with suture
72201209	8 mm ENDOBUTTON DIRECT
72202983	8 mm ENDOBUTTON DIRECT with suture
72201210	9 mm ENDOBUTTON DIRECT
72202984	9 mm ENDOBUTTON DIRECT with suture
72201211	10 mm ENDOBUTTON DIRECT
72202985	10 mm ENDOBUTTON DIRECT with suture
72201212	ENDOBUTTON DIRECT Holder
72201600	ENDOBUTTON DIRECT Scale

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10

ACL/PCL Fixation Systems


ENDOBUTTON® Family of Fixation Devices

ENDOBUTTON CL BTB Fixation Device
 ENDOBUTTON CL BTB provides dependable fixation for bone-tendon-bone ACL reconstruction procedures.

L.

Reference #	Description
72203330	15 mm ENDOBUTTON CL BTB
7210080	20 mm ENDOBUTTON CL BTB
7210081	25 mm ENDOBUTTON CL BTB
7210082	30 mm ENDOBUTTON CL BTB
7210083	35 mm ENDOBUTTON CL BTB
7210084	40 mm ENDOBUTTON CL BTB
7210085	45 mm ENDOBUTTON CL BTB
7210086	50 mm ENDOBUTTON CL BTB
7210087	55 mm ENDOBUTTON CL BTB
7210088	60 mm ENDOBUTTON CL BTB

g
B
A

XTENDOBUTTON Fixation Device
 The Smith & Nephew XTENDOBUTTON Fixation Device, used in conjunction with any Smith & Nephew ENDOBUTTON device, allows a 6 to 10 mm area of the femoral tunnel to be covered. In addition to providing an effective solution for cortical blowouts, this device can accommodate a variety of femoral tunnel drilling techniques including inside-out, outside-in and medial portal.

Reference #	Description
72200134	XTENDOBUTTON Fixation Device

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ACL/PCL Fixation Systems

ENDOBUTTON® Family of Fixation Devices


CL Device

Continuous loop of suture eliminates the need for knot tying. Ideal for a post or screw-and-washer technique.


Reference #	Description
72200839	25 mm CL Device
72200840	30 mm CL Device
72200841	35 mm CL Device
72200842	40 mm CL Device
72200843	45 mm CL Device
72200844	50 mm CL Device

ENDOBUTTON Fixation Device

This titanium button is the original endoscopic femoral soft tissue fixation device, and where the ENDOBUTTON revolution began. It can be used for single incision, soft tissue or bone-tendon-bone fixation.

Reference #	Description
013186	ENDOBUTTON, 4.0 mm x 12 mm

TRUKOR® Depth Gauge

This TRUKOR Depth Gauge easily estimates the native dimensions of the ACL for tunnel sizing.

Reference #	Description
72200948	TRUKOR Depth Gauge

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
Kopija tikra

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14.

Rotator Cuff Soft Tissue Fixation – Metallic Anchors

TWINFIX® Ti 5.0 and TWINFIX Ti 6.5 Suture Anchors

TWINFIX Ti Suture Anchors

The unique thread design of the TWINFIX Ti anchors provides extraordinary pullout strength and easy insertion in all bone conditions, and provides superior solutions to the orthopaedic surgeon.

To offer the surgeon flexibility, depending on personal preference and patient profile, TWINFIX anchors are available with DURABRAID® and ULTRABRAID® suture. ULTRABRAID is constructed from UHMW (ultra-high molecular weight polyethylene) and possesses a unique and proprietary braid configuration. This provides ULTRABRAID with its extraordinary strength properties while contributing to its superior handling qualities.

TWINFIX Ti 5.0 and TWINFIX Ti 6.5 with DURABRAID Suture

Reference #	Description
72200766	5.0 mm TWINFIX Ti Suture Anchor with two #2 preloaded DURABRAID sutures, sterile
7210300	5.0 mm TWINFIX Ti Suture Anchor with two #2 preloaded DURABRAID sutures with needles, sterile
72200769	6.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded DURABRAID sutures, sterile
7210214	6.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded DURABRAID sutures with needles, sterile

TWINFIX Ti 5.0 and TWINFIX Ti 6.5 with ULTRABRAID Suture

Reference #	Description
72200755	5.0 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID sutures*, sterile
7210708	5.0 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID sutures* with needles, sterile
72200758	6.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID sutures*, sterile
7210709	6.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID sutures* with needles, sterile
72200687	5.0 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID sutures*, sterile, Suture Anchor Pack (3 per box)
72200759	5.0 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID sutures*, sterile, no needles (5 per box)

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Kopija tikra

16

Instability Soft Tissue Fixation – Metallic Anchors

TWINFIX® Ti 3.5 and TWINFIX Ti 2.8 Suture Anchors

3.

TWINFIX Ti 3.5 and TWINFIX Ti 2.8 Suture Anchor

The unique thread design of the TWINFIX Ti 3.5 anchor and the TWINFIX Ti 2.8 anchor provides extraordinary pullout strength and easy insertion in all bone conditions, and provides superior solutions to the orthopaedic surgeon. Multiple suture options can save costs by reducing the number of anchors per procedure. The multiple point fixation also provides a stronger repair at the reattachment site.

TWINFIX Ti 3.5 with DURABRAID® suture

Reference #	Description
72200763*	3.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded DURABRAID sutures, sterile
7210302	3.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded DURABRAID sutures with needles, sterile

TWINFIX Ti 2.8 with DURABRAID Suture

Reference #	Description
72200762*	2.8 mm TWINFIX Ti Suture Anchor with one #2 preloaded DURABRAID suture, sterile
72200761	2.8 mm TWINFIX Ti Suture Anchor with one #1 preloaded DURABRAID suture, sterile

TWINFIX Ti 3.5 with ULTRABRAID® Suture

Reference #	Description
72200752	3.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID sutures, sterile
7210707	3.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID sutures with needles†, sterile
72200686†	3.5 mm TWINFIX Ti Suture Anchor with two #2 preloaded ULTRABRAID (3 per box) sutures, sterile, Suture Anchor Pack
72200760†	3.5 mm TWINFIX Ti Suture Anchor with one #2 preloaded ULTRABRAID sutures, sterile, no needles (5 per box)

TWINFIX Ti 2.8 with ULTRABRAID Suture

Reference #	Description
72200750	2.8 mm TWINFIX Ti Suture Anchor with one #2 preloaded ULTRABRAID suture, sterile
72200685†	2.8 mm TWINFIX Ti Suture Anchor with one #2 preloaded ULTRABRAID (3 per box) suture, sterile, Suture Anchor Pack
72200749†	2.8 mm TWINFIX Ti Suture Anchor with one #2 preloaded ULTRABRAID suture, sterile, no needles (5 per box)

TWINFIX Ti 2.8 and 3.5 Instruments

Reference #	Description
7209502	2.5 mm Drill for TWINFIX Ti 3.5 mm and 5.0 mm anchor
7209500	6" Drill Guide for Glenoid Placement – 6" x 3.5 mm Spiked Tip for TWINFIX Ti 2.8 mm and 3.5 mm anchor
7209506	1.8 mm Drill for TWINFIX Ti 2.8 mm anchor
7209503	Awl for TWINFIX Ti 2.8 mm anchor

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
FAST-FIX® (all-inside)

FAST-FIX 360

The FAST-FIX 360 meniscal repair system offers exceptional fixation strength, easier implant deployment, smaller insertion points that minimize disruption to the meniscus, a built in depth penetration limiter and a stiffer needle shaft for enhanced control. All this to help optimize the meniscus repair.

A.


Reference #	Description
72202467	FAST-FIX 360 straight
72202468	FAST-FIX 360 curved
72202469	FAST-FIX 360 reverse curved
72202674	FAST-FIX 360 Straight Knot Pusher/Cutter and Slotted Cannula sets
72202675	FAST-FIX 360 Curved Knot Pusher/Cutter and Slotted Cannula sets

ULTRA FAST-FIX

Reference #	Description
72201490	ULTRA FAST-FIX, straight
72201493	ULTRA FAST-FIX, AB, straight
72201491	ULTRA FAST-FIX, curved
72201494	ULTRA FAST-FIX AB, curved
72201492	ULTRA FAST-FIX, reverse curved
72201495	ULTRA FAST-FIX AB, reverse curved
72201537	Knot Pusher/Suture Cutter, straight

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

THE TECHNIQUE

www.ostecca.com

The ULTRA FAST-FIX[®] Meniscal Repair System is a unique all-inside implant system, offering the fixation strength of an open-vertical mattress stitch without the invasive surgical procedure normally required for suture-based repairs.

When the original FAST-FIX[®] Meniscal Repair System was introduced in 2001, it set the benchmark for minimally-invasive, all-inside repairs. Thanks to its preloaded implants, pre-tied sliding knot, and innovative pusher/cutter device, this innovative system lets you deploy two implants vertically or horizontally on either side of the meniscus, tighten the suture and trim the excess.

Building on its proven clinical success, the all-new ULTRA FAST-FIX system adds easier knot sliding and stronger suture. The result is a faster, more secure meniscal repair system that will help maximize the chances of a successful meniscus tear recovery.

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
2015-11-30

Kopija tikra


21

Slide the pre-tied knot, tighten the suture, and trim the excess.

4 Use the reverse-curved needle to repair tears on the inferior surface (tibial side).

5 Completed repair with vertical and horizontal mattress stitches.

4.

The only product of its kind.

Fast and easy

Unlike conventional suture-based repair systems, the ULTRA FAST-FIX system is an implant system with a pre-tied, self-sliding knot that eliminates the need for intra-articular knot tying.

Strong and clinically proven

Provides a strong, reproducible and reliable meniscal repair with biomechanical properties equal to that of the open, vertical mattress stitch suture technique.

Contains no hard device heads

Minimizes trauma to articular cartilage

Easy knot sliding with ULTRABRAID[®] suture

Offers advantages over traditional polyester suture, including higher knot-breaking strength, increased lubricity and a stronger resistance to fraying.

PEEK-OPTIMA[®] high strength non-absorbable, polymer implants

PEEK-OPTIMA polymer lends confidence that the implant will resist breakage when used with high-strength ULTRABRAID suture.


PLLA absorbable implants

Bioabsorbable version of ULTRA FAST-FIX.

Curved and reverse-curved needles

Curved needles promote safer and easier access to a multitude of tear sites. The reverse-curved needle is designed for repairing tears on the inferior surface. Because the needle's point is on the opposite side of the curve, it can safely enter the inferior area without skiving the meniscus or the tibial plateau.

A Stronger and Stiffer Repair than Competitive Devices

After 1000 cycles between 5N and 20N. All values are means +/- standard errors. Bovine meniscus model.

An evaluation of the mechanical properties of the ULTRA FAST-FIX Meniscal Repair System showed a significant increase in ultimate strength and stiffness, and a decrease in elongation compared to ETHIBOND[®], RAPIDLOC[®] and the conventional, clinically successful FAST-FIX device[®]. This in vitro mechanical evaluation of the repair devices employed cyclic and load-to-failure mechanical tests on a repaired bovine meniscus.

Vadykavikas
I. Zontop, T. Eggers, A. Musahl, V. Weimann, A. Petersen W.
Cyclic testing of flexible all-side meniscus suture anchors.
AM J Sports Med. 2005; 33:1-7.

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23.

Instability Soft Tissue Fixation

BIORAPTOR® CURVED and OSTEORAPTOR® CURVED Suture Anchors and Guide System

The BIORAPTOR Curved and OSTEORAPTOR Curved Suture Anchors and Guide System are designed specifically to manage the surgical challenges of arthroscopic instability repair.

Optimal Anchor Placement

Curved drill guides with a 25° distal bend improve the trajectory of drilling and optimize anchor placement in the glenoid*

Increased Access

- Anchors can be placed confidently in the most inferior aspect of the glenoid without the need to compromise the subscapularis tendon†
- The Curved Guide System removes the need to create additional working portals to access the repair site of a SLAP lesion†


Suture Anchor Flexibility

Achieve secure fixation from a choice of materials and sizes when using the BIORAPTOR 2.3 mm PK or OSTEORAPTOR 2.3 mm PLLA-HA or OSTEORAPTOR 2.9 mm PLLA-HA, preloaded on flexible inserters with high-strength ULTRABRAID® Suture


Reference #	Description
72203280	2.3 mm BIORAPTOR Curved PK Suture Anchor, with one ULTRABRAID (#2) Suture (COBRAID-black)
72203281	2.3 mm BIORAPTOR Curved PK Suture Anchor, with one ULTRABRAID (#2) Suture (COBRAID-blue)
72203290	2.3 mm OSTEORAPTOR Curved Suture Anchor, with one ULTRABRAID (#2) Suture (COBRAID-blue)
72203291	2.3 mm OSTEORAPTOR Curved Suture Anchor, with one ULTRABRAID (#2) Suture (COBRAID-black)
72203424	2.9 mm OSTEORAPTOR Curved Suture Anchor, with one ULTRABRAID (#2) Suture (blue)
72203423	2.9 mm OSTEORAPTOR Curved Suture Anchor, with two ULTRABRAID (#2) Sutures (blue, COBRAID-blue)

Instruments

72203156	Curved Drill Guide for 2.3 mm Curved Suture Anchors – Crown Tip
72203298	Curved Drill Guide for 2.3 mm Curved Suture Anchors – Fishmouth Tip
72203459	Curved Drill Guide for 2.3 mm Curved Suture Anchors – Reverse Fishmouth Tip
72203159	Flexible Obturator for 2.3 mm Curved Suture Anchor System – Blunt
72203460	Flexible Obturator for 2.3 mm Curved Suture Anchor System – Trocar Tip
72203160	Flexible Drill for 2.3 mm Curved Suture Anchor System
72203372	Curved Drill Guide for 2.9 mm Curved Suture Anchors – Crown Tip
72203373	Curved Drill Guide for 2.9 mm Curved Suture Anchors – Fishmouth Tip
72203484	Curved Drill Guide for 2.9 mm Curved Suture Anchors – Reverse Fishmouth Tip
72203374	Flexible Obturator for 2.9 mm Curved Suture Anchor System – Blunt
72203485	Flexible Obturator for 2.9 mm Curved Suture Anchor System – Trocar Tip
72203375	Flexible Drill for 2.9 mm Curved Suture Anchor System

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HEALICOIL® PK Suture Anchors (non-absorbable)

HEALICOIL PK Suture Anchors

The unique open-architecture design of the HEALICOIL PK Suture Anchor utilizes significantly less material than standard solid-core anchors, dramatically minimizing the volume of material implanted in the shoulder. The HEALICOIL PK anchor inserts easily into all densities of humeral bone. Torque is distributed along the length of the anchor during insertion and improved anchor thread engagement in the bone delivers superior pullout strength.

The HEALICOIL PK Suture Anchors can be used for primary or double row repairs and in the medial row of a suture bridging technique.

Non-absorbable

Our non-absorbable PK suture anchors incorporate PEEK-OPTIMA® polymer from Invibio®. PEEK-OPTIMA material provides an ideal combination of strength, biocompatibility, revisability and inherent radiolucent characteristics.

Reference #	Description
72203378	4.5 mm HEALICOIL PK Suture Anchor with two ULTRABRAID® (#2) Sutures (blue, blue-COBRAID), sterile
72203379	5.5 mm HEALICOIL PK Suture Anchor with Two ULTRABRAID (#2) Sutures (blue, blue-COBRAID), sterile
72203380	5.5 mm HEALICOIL PK Suture Anchor with Three ULTRABRAID (#2) Sutures (blue, blue-COBRAID, COBRAID-black), sterile

Instruments

72202621	3.8 mm Tapered Awl, disposable (general hole preparation for normal bone quality)
72201915	3.8 mm Tapered Awl, reusable (general hole preparation for normal bone quality)
72203335	4.5 mm HEALICOIL Awl Dilator, reusable (hole preparation for normal and hard bone quality)
72203336	5.5 mm HEALICOIL Awl Dilator, reusable (hole preparation for normal and hard bone quality)
72203482	3.5 mm HEALICOIL Spade Tip Drill for use with 4.5 mm HEALICOIL Suture Anchor
72203483	4.5 mm HEALICOIL Spade Tip Drill for use with 5.5 mm HEALICOIL Suture Anchor

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

Rotator Cuff Soft Tissue Fixation – Non-absorbable Polymer Anchors

FOOTPRINT PK Suture Anchor

FOOTPRINT PK Suture Anchor

The FOOTPRINT PK Suture Anchor consists of a non-absorbable polymer (PEEK-OPTIMA® polymer from Invibio®). The sizes include a 5.5 mm anchor, and a back-up size of 6.5 mm. The anchor is composed of two parts, the "outer shell" and the "inner plug." The design allows the anchor to be tapped in. Then, independent of implantation, the inner plug is screwed down, securing the medial anchor suture limbs that have been captured in the nose of the anchor.

The FOOTPRINT PK Suture Anchor offers a rotator cuff repair technique that allows the surgeon to precisely apply the desired amount of suture tension that is appropriate and customized for each patient. The suture tension can be adjusted even after the anchor is implanted, unlike any competitive product on the market.

Reference #	Description
72201697	5.5 mm FOOTPRINT PK Suture Anchor
72201825	5.5 mm FOOTPRINT PK Suture Anchor with disposable awl
72201698	6.5 mm FOOTPRINT PK Suture Anchor
72201826	6.5 mm FOOTPRINT PK Suture Anchor with disposable awl

Instruments

72201450	3.8 mm Straight Awl for use with FOOTPRINT PK, TWINFIX® PK FT and KINSA® RC Anchors, reusable
72201915	3.8 mm Tapered Awl, reusable (gold awl)
72200867	Mallet with nylon ends
72200879	Replacement end caps for mallet (set of 2)
72201707	Spade Drill, 4.5 mm for use in hard bone with the FOOTPRINT PK 5.5 mm anchor and the TWINFIX FT 5.5 mm Anchor
72201708	Spade Drill, 5.5 mm for use in hard bone with the FOOTPRINT PK 6.5 mm anchor and the TWINFIX FT 6.5 mm Anchor
72201899	Disposable Tension Tuner

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FOOTPRINT Ultra PK Suture Anchor

FOOTPRINT Ultra PK Suture Anchor

Available in two sizes, 4.5 mm and 5.5 mm, the FOOTPRINT Ultra PK anchor is intended for use in the lateral row in double row suture bridging techniques.

- Uniquely designed, inner plug adjustability means more **precise** fixation and suture **tensioning** of final repair
- **Reduced suture slippage** means less gap formation and rotator cuff displacement*
- Exclusive tap-in anchor design with offset barbed geometry threads means **high pullout strength** and **easy insertion***

Reference #	Description
72202901	4.5 mm FOOTPRINT Ultra Suture Anchor
72202902	5.5 mm FOOTPRINT Ultra Suture Anchor

Instruments

72202621	3.8 mm Tapered Awl-Disposable
72201915	3.8 mm Tapered Awl-Reusable
72202116	3.5 mm Spade Tip Drill for use with the 4.5 mm anchor hole preparation device for hard bone
72201707	4.5 mm Spade Tip Drill for use with the 5.5 mm anchor hole preparation device for hard bone
72201450	3.8 mm Straight Awl-Reusable
72202986	3.8 mm Straight Awl-Disposable

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2015 -11- 30

Kopija tikra

* Data on file at Smith & Nephew, ITR-3679, 2008

⚙️ Special, Made to Order products are available to meet the Individual needs of surgeons. Please appreciate that longer lead time for delivery should be expected. For more information, please contact your sales representative or Customer Service.

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
8,

SAPHYRE Ablation Probes

SAPHYRE Ablation Probes deliver powerful, consistent RF energy for vigorous soft tissue ablation and hemostasis. The jewel-cut electrodes deliver focused energy to target tissue with reduced current output requirements, generating less heat than competitive ablation probes. In addition, the SAPHYRE Bipolar Ablation Probe Coolback Shaft reduces the risk of secondary tissue heating at the return electrode site.

Reference #	Description	Shaft Color
7209686	SAPHYRE Ablation Probe, 90°	Pewter
7209685	SAPHYRE Ablation Probe, 60°	Pewter
7209684	SAPHYRE Ablation Probe, 90° high profile	Pewter

SAPHYRE II Ablation Probes

The SAPHYRE II Bipolar Ablation Probes have integrated suction for enhanced visualization in the sterile field.

Reference #	Description	Shaft Color
7210111	SAPHYRE II Ablation Probe with suction, 90°	Pewter
7210113	SAPHYRE II Ablation Probe with suction, 60°	Pewter
7210112	SAPHYRE II Ablation Probe with suction, 90° high profile	Pewter
7210491	SAPHYRE II Ablation Probe with suction, 40°	Pewter

2018-11-30
Kopija tikra

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NOTE: For use with the VULCAN RF System

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VULCAN[®] RF Arthroscopy System

DYONICS[®] GLIDER[®] Monopolar Articular Cartilage Probe

Designed exclusively for treatment of chondromalacia, indicated specifically for chondroplasty. The DYONICS GLIDER monopolar probe debrides and smooths diseased cartilage along the contoured surfaces of the joint while limiting the depth of chondrocyte death to 150 microns or less.†


L

Articular Cartilage Probe

The GLIDER probe is an early-intervention tool for treatment of articular cartilage disease. Through controlled application of RF energy on Grade II and III chondromalacia, the GLIDER probe is designed to maximize preservation of healthy cartilage tissue while restoring a smooth gliding surface.

Creates a smooth cartilage surface

In vitro studies have demonstrated the DYONICS GLIDER probe creates a smooth cartilage surface with a superficial layer of cell death.

Reference #	Description	Shaft Color
7210438*	DYONICS GLIDER Articular Cartilage Probe (requires 3.71 software)	Black

Optional	Description
7210977	Reuseable slotted cannula

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2015 -11- 30

Kopija tikra

[Handwritten signature]

* Item shown
† Data on file at Smith & Nephew, 2008.
NOTE: For use with the VULCAN RF System

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