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### Matavimo procesas

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## Automatic QC

### Measuring process

The solution pack contains three levels of QC solution. The analyzer is designed to run each level once every 24 hours. However, it is possible to set up a schedule to run QC more often, if required, as described in section *QC Schedule* in chapter 1.

The QC solutions that come from pouches in the solution pack enter the sample path through the inlet as a normal blood sample. The only difference is the position of the inlet that remains in the closed position.

Stage	Description																		
1.	When an automatic QC is scheduled to be run, it will postpone measurements etc., unless the analyzer is busy measuring a blood sample. In this case, the scheduled QC will be run after the analyzer has completed the measurement.																		
2.	<p>The QC measuring procedure begins:</p> <ul style="list-style-type: none"> <li>For QC level B only: A measurement of high oxygen is performed on a gas from a pouch that is aspirated before the QC solution.</li> <li>Measurement of the QC solution is performed as soon as it is positioned in the measuring chambers.</li> </ul>																		
3.	The result is saved in the Quality Control log.																		
4.	The result is compared with the defined control range, measuring range and statistics range.																		
5.	<p>The absence of any markings next to a parameter indicates that a parameter was measured without any fault.</p> <table> <tr> <th>Marking</th><th>Explanation</th></tr> <tr> <td>?</td><td>Error in the previous calibration, or analyzer malfunction.</td></tr> <tr> <td>W</td><td>A violated Westgard Rule.</td></tr> <tr> <td>R</td><td>A violated RiLiBÄK rule.</td></tr> <tr> <td>↑ ↓</td><td>Parameter value is outside the control range, but inside the statistics range.  Only the values within the statistics range are considered accepted and are included in the QC statistics.</td></tr> <tr> <td>↑↑ ↓↓</td><td>Parameter value is outside the statistics range and is not included in the statistics.</td></tr> <tr> <td>↑↑↑ ↓↓↓</td><td>Parameter value is outside the range of indication. Measurement is not included in the statistics.</td></tr> <tr> <td>*</td><td>Parameter values with user-defined corrections – see section <i>Parameters and input setup</i>, chapter 1 for details</td></tr> <tr> <td>.....</td><td>Parameter value could not be calculated, most likely due to a system error or malfunction. These values will for the most part be accompanied by a "?". To obtain a possible explanation, press <b>Message</b>.</td></tr> </table>	Marking	Explanation	?	Error in the previous calibration, or analyzer malfunction.	W	A violated Westgard Rule.	R	A violated RiLiBÄK rule.	↑ ↓	Parameter value is outside the control range, but inside the statistics range.  Only the values within the statistics range are considered accepted and are included in the QC statistics.	↑↑ ↓↓	Parameter value is outside the statistics range and is not included in the statistics.	↑↑↑ ↓↓↓	Parameter value is outside the range of indication. Measurement is not included in the statistics.	*	Parameter values with user-defined corrections – see section <i>Parameters and input setup</i> , chapter 1 for details	.....	Parameter value could not be calculated, most likely due to a system error or malfunction. These values will for the most part be accompanied by a "?". To obtain a possible explanation, press <b>Message</b> .
Marking	Explanation																		
?	Error in the previous calibration, or analyzer malfunction.																		
W	A violated Westgard Rule.																		
R	A violated RiLiBÄK rule.																		
↑ ↓	Parameter value is outside the control range, but inside the statistics range.  Only the values within the statistics range are considered accepted and are included in the QC statistics.																		
↑↑ ↓↓	Parameter value is outside the statistics range and is not included in the statistics.																		
↑↑↑ ↓↓↓	Parameter value is outside the range of indication. Measurement is not included in the statistics.																		
*	Parameter values with user-defined corrections – see section <i>Parameters and input setup</i> , chapter 1 for details																		
.....	Parameter value could not be calculated, most likely due to a system error or malfunction. These values will for the most part be accompanied by a "?". To obtain a possible explanation, press <b>Message</b> .																		
6.	After the measurement is complete, it is followed by a rinse. For further information on the rinse process see page 3-6.																		

## Manual QC samples

### Measuring process

The following table describes the analytical process of a manual QC measurement using the manual QC option.

Stage	Description
1.	The analyzer is ready to accept a QC sample. <ul style="list-style-type: none"><li>• "Ready" message is displayed</li><li>• Traffic light is GREEN or YELLOW.</li></ul> The desired parameters are available
2.	At the <b>Ready</b> screen, the user lifts the inlet handle to the syringe position. Press <b>Ampoule – QC</b> .
3.	The adapter is pressed against the inlet gasket and the inlet probe extends into the sample, which is automatically aspirated. <b>NOTE:</b> It is mandatory to use the adapter to minimize the risk of possible glass pieces from the ampoule getting into the system of the analyzer.
4.	When the aspiration is finished, close the inlet. The screen is now ready to accept QC information.
5.	After the measurement is complete, it is followed by a rinse. For further information on the rinse process see page 3-6.