safeCLINITUBES

989-890. 2012041.

en: Sensitive to light. Store in a dark place
cs: Citlivý ke světlu. Ukládejte v temnotě
da: Følsom over for lys. Opbevares mørkt
de: Lichtempfindlich, im Dunkeln aufbewahren
el: Ευαίσθητο στο φώς. Αποθήκευση σε σκοτεινό μέρος
es: Sensible a la luz. Almacenar en lugar oscuro
et: Tundlik valgusele. Säilitada pimedas
fi: Valolle herkkä. Säilytettävä pimeässä
fr: Sensible à la lumière. Conserver à l'obscurité
hu: Fényre érzékeny Sötét helyen tárolja
it: Fotosensibile. Conservare al buio
ja: 感光性ですので暗所で保存してください。
lt: Jautrus šviesai. Laikyti tamsioje vietoje
lv: Jūtīgs pret gaismu. Uzglabāt tumšā vietā
nl: Lichtgevoelig, niet blootstellen aan direct zonlicht. Opslag in een donkere ruimte
no: Lysfølsom, lagres mørkt
pl: Wrażliwe na światło
pt: Sensível a luz. Armazenar num local escuro
ro: Sensibil la lumină. A se păstra în spații întunecoase
ru: Светочувствителен. Хранить в темном месте
sk: Citlivé na svetlo. Skladujte na tmavom mieste
sr: Osetljivo na svetlost. Čuvati na tamnom mestu
sv: Ljuskänslig. Förvara mörkt
zh: 对光敏感。请避光保存。

en	INSTRUCTIONS	It	ISTRUKCIJOS
cs	INSTRUKCE	Iv	INSTRUKCIJAS
da	BRUGSANVISNING	nl	GEBRUIKSAANWIJZING
de	GEBRAUCHSANLEITUNG	no	INSTRUKSJON
el	ΟΔΗΓΙΕΣ	pl	INSTRUKCJE
es	INSTRUCCIONES	pt	INSTRUÇÕES
et	KASUTUSJUHISED	ro	INSTRUCȚIUNI
fi	OHJEET	ru	ИНСТРУКЦИИ
fr	MODE D'EMPLOI	sk	NÁVOD
hu	UTASÍTÁSOK	sr	INSTRUKCIJA
it	ISTRUZIONI	sv	BRUKSANVISNING
ja	インストラクション	zh	指示

CLINITUBES

en: Product Information/cs: Informace o výrobku/da: Produktinformation/de: Produktinformationen/el: Πληροφορίες προϊόντος/es: Información del producto/et: Toote kirjeldus/fi: Tuoteseloste/fr: Informations sur les produits/hu: Termék információ/it: Informazione sul prodotto/ja: 製品情報/lt: Informacija apie gamini/lv: Informācija par produktu/nl: Productinformatie/no: Produktinformasjon/pl: Informacje o produkcie/pt: Informação sobre o produto/ro: Informații despre produs/ru: Перечень продуктов/sk: Informácia o produkte/sr: Informacije o proizvodu/sv: Produkt-information/zh: 产品信息

en: Ref (glass) cs: Ref (sklo) da: Ref (glas) de: Ref (Glas) el: Αναφ. (γυάλινο) es: Ref (vidrio) et: Ref (klaas) fi: Ref (lasi) fr: Réf. (verre) hu: Cikkszám (üveg) It: Rif. (vetro) ja: Ref. (ガラス) lt: Medžiaga (stiklas) lv: Ref. (stikls) nl: Ref (glass) no: Ref (glass) pl: Ref (szkło) pt: Ref (vidro) ro: Ref (sticlă) ru: Код (стекл.) sk: Ref (sklo) sv: Ref (staklo) zh: 参考 (玻璃)	Ref (plastic) Ref (plast) Ref (Kunststoff) Αναφ. (πλαστικό) Ref (plástico) Ref (plastik) Ref (muovi) Réf. (plastique) Cikkszám (műanyag) Rif. (plastica) Ref. (プラスチック) Medžiaga (plastikas) Ref. (plastmasa) Ref (kunststof) Ref (plastik) Ref (plastik) Ref (plástico) Ref (plastic) Код (пласт.) Ref (umelá hmota) Ref (plastika) Ref (plast)	Name Označení Navn Name Όνομα Modelo Nimi Nimi Type Megnevezés Nome 型式 Pavadinimas Nosaukums Type Navn Nazwa Modelo Denumire Наименование Назov Ime Typ Typ 製品名	Volume Objem Volumen Volumen Όγκος Volumen Maht Tilavuus Volume Úrtartalom Volume Volume 容量 Töris Tilpums Volume Volum Objetość Volume Volum Объём Objem Zapremina Volym 容量	Capillaries per kit Počet kapilár v sadě Kapillærør pr. kit Kapillaren pro kit Αναφ. (πλαστικό) Capillares por kit Kapillaaride kogus Kapillaarien määrä Capillaires par kit Kapillárisok száma a csomagban Capillari per kit キットあたりのキャピラリー数 Kapiliarų kiekis komplekte Kapilari komplektā Capillaren per buisje Kapillærør per kit Liczba kapilar w zestawie Capillares por kit Capilare per trusă Кон-во капилляров в наборе Počet kapilár na súpravu Kapilara po priboru Kapillärer per förpackning 毛細管/盒	Heparin type Druh heparinu Heparintype Heparintyp Τύπος ηπαρίνης Tipo de heparina Heparini mark Heparinityppi Type d'héparine Heparin típus Tipo di eparina Heparino típus Heparina tips Type heparine Type heparin Rodzaj heparyny Tipo de heparina Tip heparinā Тип гепарина Typ heparinu Typ heparina Typ av heparin 肝素类型	Heparin IU/mL blood Heparin - IU/mL krve Heparin IU/mL blod Heparin IU/mL Blut Ηπαρίνη IU/mL αίματος/ Heparina IU/mL sangre Heparini kogus IU/mL vere kohta Heparinia IU/mL verta Héparine IU/mL de Sang Heparin IU/mL vér Eparina IU/mL sangue Heparino IU/mL kraujo Heparins IU/ml asins Heparine I.E./ml blood/ Type heparin Ilość heparyny IU/mL krwi Heparina IU/ml. sangue Heparinā IU/ml sāgne Гепарин ME/мл крови Heparin IU/ml krv Heparin IU/mL krvi Heparin IU/mL blod 肝素IU/mL 血液	Mixing wires & capillary caps Michaci drátky & čepičky kapilár Blandestifter & kapillærhætter Mischstäbchen & Kapillärverschlusskappen Σύρματα ανάμιξης και καλύμματα τριχοειδούς αγγείου Hierrecillos mezcladores y tapones de capilares Magnetsegajad ja otsakorgid Sekoituslangat & tulpat Agitateurs et bouchons de capillaire Keverőtűskék és lezáró kupakok Ancorette di miscelazione e tappini per capillari Mikszingfűvay és kanyilari-kyapp Magnetsegajad & otsakorgid Magnetsegajad & otsakorgid Masiymo stypeliai ir kapilario kamsteliai Mengijertjes & afsluitdopjes Blandestav & Kapillærkorker Szyfty mieszażące i zatyczki kapilar Ferrinhos de mistura e tampões de capilar Fire de amestecare şi capace capilare Проволочки-смесители и колпачки для капилляров Miešadje drótky & kapilárne čiarocky Zicice za mesanje & zatvaraci za kapilare Blandstift & kapillärhåttor 混合丝 & 毛細管帽	Magnet Magnet Magnet Magnet Μαγνήτης Imanes Magnet Magneetti Aimant Măgnes Magnete Магнетас Magnēts Magneet Magnet Imās Magnet Магнит Magnet Magnet Magnet 磁鉄	Recommended for use with the following analyzers: Doporučeno pro uživatele přístroje řady: Anbefalet til brug sammen med følgende apparater: Empfohlen zur Verwendung an folgenden Analysatoren: Συστήνεται για χρήση στους κάτωθι αναλυτές: Se recomienda usar con los siguientes analizadores: Soovitatakse kasutada järgmistele analüsaatoritele: Suositellaan käytettäväksi seuraavien analysaattoreiden kanssa: Utilisation recommandée avec les analyseurs suivants: Felhasználása a következő analizátorokhoz javasolt: Raccomandato per l'uso con i seguenti analizzatori: 次の装置での使用をお勧めします。 Rekomenduojaama naudoti tokiems analizatoriams: Rekomendēts izmantošanai uz sekojošiem analizatoriem: Aanbevolen voor gebruik met de volgende analysers: Anbefalt for bruk med følgende analysatorer: Polecany do stosowania z następującymi analizatorami: Recomendado para utilização nos seguintes analisadores: Recomandat pentru utilizare cu următoarele analizoare: Рекомендуем использовать со следующими анализаторами: Doporučené pre použitie s týmito analyzátorami: Preporučuje se upotreba uz sledeće analizatore: Rekomenderad användning på följande analysatorer: 推荐在下列分析仪上使用:
905-661 942-869		D956G-70-35x1 D957G-70-35x5	35 µL 35 µL	75 375	* *	70 70	√ √	√ -	ABL5, ABL500 Series, ABL600 Series, ABL700 Series, ABL800 BASIC ***, ABL800 FLEX, ABL8x7 FLEX, OSM2
942-885 942-873	942-889	D941G-240-55x5 D941G-NA-240-55x10 D941P-NA-240-55x3	55 µL 55 µL 55 µL	375 750 750	** ** **	240 240 240	√ - -	- - -	ABL700 Series, ABL800 BASIC ***, ABL800 FLEX, ABL8x7 FLEX, OSM2
	942-888	D957P-70-70x1	70 µL	250	*	70	√	-	ABL77 Series, ABL90 FLEX,
942-875 942-876	942-890 942-891	D941G-240-85x5 D941P-240-85x1 D941G-NA-240-85x10 D941P-NA-240-85x3	85 µL 85 µL 85 µL 85 µL	375 250 750 750	** ** ** **	240 240 240 240	√ - - -	- - - -	ABL300/330, ABL500 Series, ABL 600 Series, ABL77 Series
905-663 942-878	942-892	D956G-70-100x1 D957G-70-100x5 D957P-70-100x1	100 µL 100 µL 100 µL	75 375 250	* * *	70 70 70	√ √ √	√ - -	ABL330, ABL5, ABL50, ABL700 Series, ABL800 BASIC ***, ABL800 FLEX, ABL90 FLEX
905-664 942-880	942-893	D956G-70-125x1 D957G-70-125x5 D957P-70-125x1	125 µL 125 µL 125 µL	50 250 250	* * *	70 70 70	√ √ √	√ - -	ABL30, ABL505/555, ABL600 Series, EML100/105, ABL77 Series, ABL8x7 FLEX
942-882		D941G-80-140x5	140 µL	250	**	80	√	-	ABL30
905-666 942-884		D956G-70-210x1 D957G-70-210x5	210 µL 210 µL	50 250	* *	70 70	√ √	√ -	ABL600 Series, ABL700 Series, ABL800 BASIC ***, ABL800 FLEX, ABL8x7 FLEX, OSM3

Radiometer Medical ApS
Åkandervej 21
2700 Brønshøj
Denmark
www.radiometer.com

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, AQUIRE, PICO, CLINITUBES and QUALICHECK are trademarks of or used under license by Radiometer Medical ApS.

© 2012 Radiometer Medical ApS. All Rights Reserved.
989-890. 2012041.

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Intended use: Heparinized plastic capillaries intended for blood sampling.

WARNING – Risk of infection
Blood sampling should be performed by authorized personnel only. Always handle blood samples and glass capillaries with care. Avoid direct contact with the sample by using certified rubber gloves.

Sampling Procedures Using Mixing Wires and Capillary Caps

- Select the puncture site and increase the circulation in the region.
- Loosely mount a capillary cap on one end of the capillary.
- Insert a mixing wire in the capillary and allow it to slide to the same end as the loosely mounted capillary cap.
- Apply antiseptic to the sampling site. Then puncture the skin so that drops of blood rapidly form. Fill the capillary with blood from the middle of the blood drop to prevent air from getting into the capillary.
- Position the capillary cap on the capillary cap on the other end of the capillary.
- Move the mixing wire 20 times along the full length of the capillary using the magnet to mix the blood and the anticoagulant together.
- Analyze the sample within 10 minutes. Note that for samples with a $pO_2 > 80$ mmHg (10.7 kPa) a positive bias from 1 to 9 % is observed. The bias is dependent on the pO_2 , pH and the temperature.
- Before the sample is analyzed, mix it again. Do not remove the mixing wire to before the sample is aspirated for measurement. Slide the mixing wire to the end of the capillary opposite to that from which the blood is to be aspirated.
- Remove both capillary caps.
- Aspirate the sample into the analyzer or system.

Sampling Procedure Without Using Mixing Wires and Capillary Caps
Please refer to points 1, 4 and 10 in the list of procedures above.

WARNING – Risk of infection
Discard the waste as infectious waste ¹⁾.
¹⁾ CLSI document GPSS-A2, Clinical Laboratory Waste Management

WARNING – Risk of infection and erroneous results
Reuse of single-use devices may lead to infection of patients and erroneous results.

For further instructions please refer to the relevant Operator's Manual for your analyzer or system.

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