NycoCardTM

CRP

EN	5-8	NL	29-32	HU	53-56
DE	9-12	EL	33-36	RO	57-60
FR	13-16	PL	37-40	NO	61-64
IT	17-20	LT	41-44	SV	65-68
ES	21-24	CS	45-48	DA	69-72
PT	25-28	RU	49-52	FI	73-76

REF 1116078

1116120 Rev. B 2016/01

SYMBOLS/ SYMBOLE/ SYMBOLES/ SIMBOLI/ SÍMBOLOS/ SÍMBOLO/ SYMBOLEN/
 ΣΥΜΒΟΛΑ/ SYMBOLE/ SIMBOLIAI/ SYMBOLY/ СИМВОЛЫ/ JELMAGYARÁZAT/
 SIMBOLURI/ SYMBOLER/ SYMBOLER/ SYMBOLER/ SYMBOLIT

Conformity to the European directive 98/79/EC
 Konformität mit der europäischen Richtlinie 98/79/EC
 Conformité avec la directive européenne 98/79/CE
 Conformità con la direttiva europea 98/79/CE
 Cumplimiento de la directiva europea 98/79/CE
 Conformidade com a Directiva Europeia 98/79/EC
 Voldoet aan de Europese richtlijn 98/79/EG
 Συμμόρφωση με την οδηγία 98/79/ΕΚ για προϊόντα IVD
 Zgodność z dyrektywą IVD 98/79/EC
 Atitiktis Europos direktyvai 98/79/EB
 Ve shodě s evropskou směrnicí 98/79/ES
 Соответствие с европейской директивой 98/79/EC
 Megfelel a 98/79/EK európai irányelvnek
 Conformitate cu directiva europeană 98/79/EC
 Överensstämelse med EU-direktiv 98/79/EC
 Överensstämmer med EU-direktiv 98/79/EC
 Överensstämelse med EU-direktiv 98/79/EC
 Yhdenmukainen Euroopan direktiivin 98/79/EC kanssa

IVD

In Vitro Diagnostic Medical Device
 Medizinisches Gerät zur *In-vitro*-Diagnostik
 Dispositif médical de diagnostic *in vitro*
 Dispositivo medicale per la diagnostica *in vitro*
 Dispositivo médico de diagnóstico *in vitro*
 Dispositivo médico de diagnóstico *in vitro*
 Medisch hulpmiddel voor *in-vitro* diagnostiek
In vitro διαγνωστικό ιατροτεχνολογικό προϊόν
 Urządzenie medyczne do diagnostyki *in vitro*
In vitro diagnostinis medicinos prietaisai
 Diagnostické zdravotnícké prostriedky *in vitro*
 Медицинское приспособление для диагностики *In Vitro*
In vitro diagnosztikai orvostechnikai eszköz
In vitro dispozitiv medical pentru diagnostic
 Medisinsk utstyr til *in vitro* diagnostikk
In vitro diagnostisk medicinsk utrustning
In vitro diagnostisk medicinsk udstyr
In vitro diagnostinen lääkinnällinen laite

2°C

8°C

2...8°C

2...25°C

Storage temperature (upper and lower limits)
 Lagerungstemperatur (Ober- und Untergrenzen)
 Température de conservation (limites supérieure et inférieure)
 Temperatura di conservazione (limiti superiore e inferiore)
 Temperatura de almacenamiento (límites superior e inferior)
 Temperatura de conservação (limites superior e inferior)
 Opslag temperatuur (boven- en ondergrenzen)
 Θερμοκρασία αποθήκευσης (υψηλά και χαμηλά όρια)
 Temperatura przechowywania (zakres temperatur)
 Laikymo temperatūra (viršutinė ir apatinė riba)
 Skladujte při teplotě (horní a dolní hranice)
 Температура хранения (верхняя и нижняя границы)
 Tárolási hőmérséklet (alsó és felső határérték)
 Temperatura de depozitare (limita inferioară și superioară)
 Oppbevaringstemperatur (øvre og nedre grenser)
 Förvaringstemperatur (övre och nedre gränser)
 Oppbevaringstemperatur (øvre og nedre grænser)
 Säilytyslämpötila (ylimmä ja alimmat rajat)

REF

Catalogue number
Bestellnummer
Numéro de référence du catalogue
Numero di catalogo
Número de catálogo
Número do catálogo
Catalogus nummer
Αριθμός καταλόγου
Numer katalogowy
Katalogo numeris
Katalogové číslo
Номер по каталогу
Katalógusszám
Catalogul numărul
Katalognummer
Artikelnummer
Katalognummer
Luettelonumero

LOT

Lot number
Chargen-Nr.
Numéro de lot
Numero di lotto
Número de lote
Número de lote
Lot nummer
Αριθμός παρτίδας
Numer partii
Partijos numeris
Číslo šarže
номер партии
Gyártási szám
Numarul de lot
Lotnummer
Lotnummer
Lot-number
Eränumero

Expiry date (year-month)
Verfallsdatum (Jahr-Monat)
Date de péremption (année-mois)
Data di scadenza (anno-mese)
Fecha de caducidad (año-mes)
Data de validade (ano-mês)
Niet gebruiken na (jaar-maand)
Ημερομηνία λήξης (έτος-μήνας)
Data ważności (rok-miesiąc)
Tinkamumo laikas (metai-mėnuo)
Datum expirace (rok-měsíc)
Срок годности (год-месяц)
Lejáratí idő (év-hónap)
Data expirării (an-lună)
Utløpsdato (år-måned)
Utgångsdatum (år-månad)
Udløbsdato (år-måned)
Viimeinen käyttöpäivämäärä
(vuosi-kuukausi)

Consult the instruction for use
Bitte die Gebrauchsanleitung beachten
Consulter le mode d'emploi
Consultare le istruzioni per l'uso
Consulte las instrucciones de uso
Consulte as instruções de utilização
Raadpleeg de gebruiksaanwijzing
Συμβουλευτείτε τις οδηγίες χρήσης
Sprawdzić w instrukcji użytkowania
Skaitykite naudojimo instrukciją
Viz návod k použití.
См. инструкцию по использованию
További információ a használati útmutatóban talál
Consultați instrucțiunile de utilizare
Se bruksanvisningen
Se bruksanvisningen
Se brugervejledningen
Katso lisätietoja käyttöohjeesta

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PRODUCT DESCRIPTION

Intended use

NycoCard™ CRP is an *in vitro* diagnostic test for the quantitative determination of C-reactive protein (CRP) in human whole blood and in human serum and plasma. The measurement of CRP provides information for the detection and evaluation of infection, tissue injury, inflammatory disorders and associated diseases.

Summary and explanation of the test

C-reactive protein (CRP) is one of the cytokine induced acute-phase proteins, the levels of which rise during a general, unspecific response to infections and non-infectious inflammatory processes^{3,5}. In healthy persons the serum or plasma CRP levels are below 5 mg/L¹. This threshold is often exceeded within four to eight hours after an acute inflammatory event, with CRP values reaching approximately 20 to 500 mg/L². As elevated CRP levels are always associated with pathological changes, the CRP assay provides information for the diagnosis, therapy, and monitoring of inflammatory diseases²⁻⁵.

Test principle

NycoCard™ CRP is a solid phase, sandwich-format, immunometric assay. In the test well of the device there is a membrane coated with immobilised CRP-specific monoclonal antibodies. A diluted sample is applied to the test device. When the sample flows through the membrane, the C-reactive proteins are captured by the antibodies. CRP trapped on the membrane will then bind the gold-antibody conjugate added, in a sandwich-type reaction. Unbound conjugate is removed from the membrane by the washing solution. A paper layer underneath the membrane absorbs excess liquid. In the presence of a pathological level of CRP in the sample, the membrane appears red-brown with colour intensity proportional to the CRP concentration of the sample. The colour intensity is measured quantitatively with the NycoCard™ READER II.

Kit contents: 48 Tests

TD/Test Device	1 x 48 units
Plastic device containing a membrane coated with monoclonal anti-CRP antibodies.	
R1/Dilution Liquid	1 x 48 x 0.4 mL
Borate buffer (pH 9.0) and detergents.	
R2/Conjugate	1 x 3.5 mL
Solution containing monoclonal anti-CRP antibodies labeled with ultra-small gold particles.	
R3/Washing Solution	1 x 3.0 mL
Phosphate buffered NaCl solution (pH 7.4) and detergents.	

Material required

(not supplied with the kit)

- Pipette (50 µL) and pipette tips.
- Capillary tube holder.
- NycoCard™ READER II.
- Capillary tubes 5 µL.
- Alere Afinion™ CRP Control.

Warnings and precautions

- For *in vitro* diagnostic use.
- Do not mix components from different kit containers.
- Do not use test kits after the expiry date.
- The reagents contain sodium azide in concentrations < 0.1% as a preservative. Sodium azide is a toxic agent. Avoid contact with eyes and skin. Flush abundantly with water upon disposal if reagents are spilled.
- R1/Dilution Liquid contains detergents that are irritating to eyes and skin. Flush with water if spilled or in contact with skin.

TEST CHARACTERISTICS

Analytical specificity

Monoclonal antibodies specific to human CRP are used in the test. No other human blood components are found to cross react with CRP in the NycoCard™ CRP system.

Standardisation

NycoCard™ CRP is calibrated against ERM®-DA/IFCC reference preparation.

Measuring range

NycoCard™ READER II, SN ≥ 54170:

Whole blood samples: 8-200 mg/L.

Serum or plasma samples: 5-120 mg/L.

Measuring interval: 1 mg/L.

Reference range

< 5 mg/L¹

Precision

A precision study was run according to CLSI guideline EP5-A2. Within-device (total) precision was determined for a plasma sample and Alere Afinion™ CRP Control levels C I and C II for 20 days,

each sample assayed in duplicate, twice a day. An EDTA-whole blood sample was assayed for 4 days, each sample assayed in 20 replicates per day. The data are summarized in the table below.

Sample	N (days)	Mean CRP (mg/L)	Within-device precision (total) CV (%)
Serum	20	97	5.2
Control C I	20	22	6.4
Control C II	20	69	6.6
EDTA-whole blood	4	51	3.7

Limitations of the test

Heparin, citrate or EDTA, used as anti-coagulant, does not interfere with the test. Increased levels of bilirubin, lipids or rheumatoid factor (RF) have none or negligible effect on the test results. Very high leukocytes levels ($> 40 \times 10^9/L$) in the sample might cause falsely elevated results (see "Trouble shooting"). Corrections should be made for hematocrit deviating from 40% (see "Effect of hematocrit").

STABILITY AND STORAGE

Unopened kits

The expiry dates apply to kit components stored in sealed, original containers at 2-8°C. Avoid direct sunlight or exposure to temperatures above 25°C. Do not freeze.

Opened kits

Avoid direct sunlight or exposure to temperatures above 25°C.

TD/Test Device: Stable until expiry date at 2-8°C or 6 weeks at 15-25°C (refrigerate at the end of the working day). The sealing foil should be removed immediately before use.

R1/Dilution Liquid: Stable until expiry date, refrigerated or at room temperature (2-25°C).

R2/Conjugate: Stable until expiry date at 2-8°C or 6 weeks at 15-25°C (refrigerate at the end of the working day).

R3/Washing Solution: Stable until expiry date at 2-8°C or 6 weeks at 15-25°C (refrigerate at the end of the working day).

Samples

Whole blood (with anticoagulant): Stable for 3 days at 2-8°C.

Serum or plasma: Stable for 3 days at 2-8°C. For prolonged storage, keep samples frozen (-20°C).

Sample diluted with R1: Analyse diluted sample within 1 hour (diluted EDTA blood or plasma within 15 minutes).

TEST PROCEDURE

Important procedural notes!

- Bring test tube(s) with R1/Dilution Liquid to room temperature (15-25°C) before use.
- R2/Conjugate, R3/Washing Solution and TD/Test Device can be used cold or equilibrated to room temperature.
- Mark the R1 tube and the test device with suitable patient or control identification.
- Always use a new pipette tip.
- Do not touch the test membrane with the pipette tip.
- The procedure steps should be performed successively without any interruptions.
- Caps should be screwed tightly on after use.

Sample material

Capillary blood, venous blood with or without anticoagulant (heparin, citrate or EDTA), serum and plasma (heparin, citrate or EDTA) can be used.

Quality control

The Alere Afinion™ CRP Control should be used to confirm the efficacy of the reagents and correct performance of the test. Test this control according to the same procedure as for a patient sample.

Note! Use the "CRP Serum/Plasma" menu for reading the result with Nycocard™ READER II. The measured value should be within the acceptance limits stated on the vial label.

Test procedure (illustrations on page 78)

1 Dilute sample

Fill a 5 µL capillary with patient sample or Alere Afinion™ CRP Control, and drop the capillary into the tube with R1/Dilution Liquid. Close the tube and mix well for 10 seconds.

Note! Avoid air bubbles in the capillary and excess sample on the outside of the capillary.

2 Apply sample

Apply 50 µL diluted sample or diluted Alere Afinion™ CRP Control to the TD/Test Device. Allow the sample to soak into the membrane (approx. 30 seconds).

Note! Avoid air bubbles on the membrane. Do not touch the membrane with the pipette tip.

3 Apply R2/Conjugate

Apply one drop R2/Conjugate to the TD/Test Device. Allow the reagent to soak into the membrane (approx. 30 seconds).

Note! The droplet bottle should be held vertically, about 1 cm above the membrane.

4 Apply R3/Washing Solution

Apply one drop R3/Washing Solution to the TD/Test Device. Allow the reagent to soak into the membrane (approx. 20 seconds).

Note! The droplet bottle should be held vertically, about 1 cm above the membrane.

5 Read the result

Read the result within 5 minutes using the Nycocard™ READER II. Follow the READER II user instruction manual.

Note! Use the "CRP Whole Blood" menu for reading whole blood samples, and the "CRP Serum/Plasma" menu for reading serum or plasma samples and the Alere Afinion™ CRP Control.

INTERPRETATION OF RESULTS

Interpret the NycoCard™ CRP test results with careful consideration to the patient's medical history, clinical examinations and other laboratory results. If the test result is questionable or if clinical signs and symptoms appear inconsistent with the test result, re-test the sample or confirm the result using another method. Analyse control materials frequently to verify the performance of the NycoCard™ READER II test system.

Effect of hematocrit

The "CRP Whole Blood" menu of the NycoCard™ READER II is calibrated to read the CRP serum concentration of a blood sample with a hematocrit of 40%.

If the actual hematocrit value deviates from 40%, the result should be corrected by multiplying with the respective factor in the table:

Hct%	Factor	Hct%	Factor
20-29	0.8	56-58	1.4
30-36	0.9	59-61	1.5
37-42	1.0	62-63	1.6
43-47	1.1	64-65	1.7
48-51	1.2	66-67	1.8
52-55	1.3	68-69	1.9

Reference range,
Hct: Women: 35-44%. Men: 39-48%.

TROUBLE SHOOTING

Possible cause of unexpected low result	Corrective action
Incorrect sample volume.	Retest the sample. Make sure that the capillary is filled completely. Avoid air bubbles.
Incorrect volume of diluted sample applied.	Check the pipette (change the Mini-Pet once a year). Make sure that the pipette tips fit well on the pipette. Retest the sample.
Too long delay (> 2 minutes) between steps 2 and 3.	Retest the sample. Steps 2 to 4 should be run without any interruptions.
Possible cause of unexpected high result	Corrective action
Very high leukocytes count (> 40x10 ⁹ /L) might give turbid sample dilutions. The time before the sample is absorbed into the test device will increase due to a partially blocked membrane.	Clarify the diluted sample by centrifugation. Retest the sample. Alternatively, use serum or plasma.
Incorrect sample volume.	Retest the sample. Avoid excess sample on the outside of the capillary.
Incorrect volume of diluted sample applied.	Check the pipette (change the Mini-Pet once a year). Make sure that the pipette tips fit well on the pipette. Retest the sample.
A reagent has been applied to the TD/Test Device before the preceding reagent has soaked completely into the membrane.	Retest the sample. Make sure that each reagent soaks completely into the membrane before applying the next reagent.
Too long delay (> 2 minutes) between steps 3 and 4.	Retest the sample. Steps 2 to 4 should be run without any interruptions.

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© 2016 Alere. All rights reserved. The Alere Logo and Alere are trademarks of the Alere group of companies.