Isotonac·3, Isotonac·4, Hemolynac·3, Hemolynac·3N, Hemolynac·5, Hemolynac·310, Hemolynac·510, Cleanac·3, Cleanac, Cleanac·710

PRODUKTO IDENTIFIKACIJA

Pavadinimas: Isotonac·3, Isotonac·4, Hemolynac·3, Hemolynac·3N, Hemolynac·5, Hemolynac·310, Hemolynac·510, Cleanac·3, Cleanac, Cleanac·710

PASKIRTIS

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BIBLIOGRAFIA

1. Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition, CLSI Document H26-A2, Vol. 30, No. 14, 2010.
2. Health Care Financing Administration. Clinical Laboratory Improvement Amendments of 1998; Final Rule. Federal Register, Vol. 57, [42 CFR § 493.1217, 493.1218], 1992.

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NIHON KOHDEN FIRENZE S.r.l

Via Torta, 72/74 - 50019 Sesto Fiorentino-Firenze, Italy
Tel +39 055 3045 1 - Fax +39 055 308548

Isotonac-3, Isotonac-4, Hemolynac-3, Hemolynac-3N, Hemolynac-5, Hemolynac-310, Hemolynac-510 Cleanac-3, Cleanac, Cleanac-710

ENGLISH

PRODUCT IDENTIFICATION

name: Isotonac-3, Isotonac-4, Hemolynac-3, Hemolynac-3N, Hemolynac-5, Hemolynac-310, Hemolynac-510, Cleanac-3, Cleanac, Cleanac-710

INTENDED USE

For in vitro diagnostic use only. Reagents to be used on NIHON KOHDEN analyzers in combination as indicated in Table 1.

Isotonac-4, Hemolynac-310, Hemolynac-510 and Cleanac-710 are only intended for human blood. In regard to the intended use, please refer to Table 2.

STORAGE AND SHELF LIFE AFTER FIRST OPENING

Store between 1 and 30 °C. Check reagent integrity and expiry date, before installing the device. Once opened, the reagent must be used within the expiry date indicated in Table 4. When not in use, use the original cap to close the container except Hemolynac-310, Hemolynac-510 and Cleanac-710. Do not force.

COMPOSITION

Active ingredients

Ref. Table 3

WARNING AND PRECAUTIONS

Labelling information

Ref. Table 5.

Disposal information

Unused product and emptied containers should be disposed in accordance with local environmental regulations. Instrument wastes, as potential biohazardous, must be handled with care and disposed in accordance with local environmental regulations.

SPECIMEN

The devices have been designed to be used with whole blood samples. Refer to Operator's Manual for sample collection and storage conditions.

PROCEDURE

Refer to Operator's Manual for reagents installation and priming, and for instrument operation. When installing a new batch of diluent or lysing reagent, check the background according to Operator's Manual. Refer to Operator's Manual for working temperature.

Use each reagent in its own packaging (do not transfer). Follow GLP recommendations. If the instrument connections were not compatible, work according to Good Laboratory Practice (GLP), in order to avoid any kind of contamination, either from external environment to the device, or from the device to the external environment.

METHODOLOGY AND CALCULATION OF ANALYTICAL RESULTS

Refer to Operator's Manual for information on methodology and calculation of analytical results.

INTERNAL QUALITY CONTROL

Background verification

Before running testing samples check background count according to Operator's Manual, and verify that all the values for each parameter are within Operator's Manual limits. If not, follow Operator's Manual prescriptions. In case of problems attributed to one of the reagents, install a new unit of reagent and verify the background again. Do not run testing samples until the background is within the specified limits.

Calibration and calibration verification procedures

Instrument calibration should be verified in the following conditions:

1. after first installation
 2. at least every six months (or according to manufacturer's requirements for calibration frequency, if more stringent)
 3. after a major repair or a major change (replacement of critical parts that may influence test performance, software upgrade, complete change of reagents)
 4. when quality control data reflect a significant change of baseline accuracy (drift), or an out-of-control condition
 5. in case of replacement of the lysing reagent Hemolynac-3 with the lysing reagent Hemolynac-3N and vice versa
- Calibration procedures should be preceded by a preventive maintenance (verification of baseline accuracy and precision: ref. CLSI document H26-A2). Calibration and calibration verification must be performed using appropriate calibration material (fresh whole blood samples which have been analyzed by reference methods, or appropriate calibrator), following Operator's Manual prescriptions. When calibration is required, verify instrument calibration immediately after calibrating, according to Operator's Manual instructions (refer also to CLSI document H26-A2).

Daily Quality Control Procedures

In order to assure the accuracy and reliability of test results, it is recommended to perform quality control procedures, according to Operator's Manual instructions, using appropriate control material (retained patient specimens or stabilized control blood), immediately before, during, or immediately after each analysis run (intended as an analysis sequence). It is recommended to run two different levels of stabilized control blood (refer to CLSI document H26-A2).

BIBLIOGRAPHY

1. Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition, CLSI Document H26-A2, Vol. 30, No. 14, 2010.
2. Health Care Financing Administration. Clinical Laboratory Improvement Amendments of 1998; Final Rule. Federal Register, Vol. 57, [42 CFR § 493.1217, 493.1218], 1992.

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First Edition: 16 Mar 2004

Printed: 2004 03 30

Hematology Control for WBC 5 Part Differential MEK-5DN/5DH/5DL

The CE mark is a protected conformity mark of the European Community. The products herewith comply with the requirements of the IVD Directive 98/79/EEC.

General

The MEK-SD hematology control for WBC 5 part differential is a tri-level control designed to monitor the performance of Nihon Kohden hematology analyzers. Refer to the assay tables for specific instrument models.

This hematology control is an in vitro diagnostic reagent composed of human erythrocytes, mammalian leukocytes and platelets in a plasma-like fluid with preservatives. It is a stable material that provides a means of checking the accuracy and precision of hematology analyzers. It should be handled in the same manner as a blood sample.

To achieve optimum performance and safety, thoroughly read this instruction and the hematology analyzer operator's manual before using this hematology control.

WARNING

- **POTENTIAL BIOHAZARDOUS MATERIAL.** Handle as if product were capable of transmitting disease. Human blood components used in preparation of this product have been tested by FDA approved methods and found non reactive for the presence of hepatitis B surface antigen (HBsAg), and antibody to hepatitis C virus (HCV), and human immunodeficiency virus (HIV-1/HIV-2). However, no known test method can offer complete assurance that specimens of human origin will not transmit infectious disease. When handling or disposing of vials, follow precautions used for patient specimens as described in CDC and FDA recommendations and OSHA blood borne pathogen rules.
- Wear rubber gloves to protect yourself from infection.

CAUTION

- Do not swallow the hematology control. It may be infectious.
- Dispose of the hematology control by following your local laws for disposing of infectious medical waste.

NOTE

- Use the hematology control on Nihon Kohden hematology analyzers which use light scatter technique.
- Do not use the hematology control after the expiration date. Otherwise the measurement data may not be accurate.

- Handle the hematology control in the same manner as a blood sample.
- Use the hematology control at temperature from 18 to 30°C. Otherwise the measurement data may not be accurate.
- Mix the hematology control thoroughly before use.
- Store the hematology control at temperatures between 2 and 8°C. Do not freeze the reagent. If the reagent is frozen or is stored at temperature outside the range, the measurement data may not be accurate.
- Platelet cannot be measured by precipitation or centrifuge.
- The hematology control cannot be used for samples whose WBC is differentiated manually.
- Follow the instructions in the operator's manual of the hematology analyzer for handling the samples.
- Opened hematology control should be used within 14 days.

Instructions for Use

1. Take the hematology control out of the refrigerator and allow it to warm at temperatures of 18 to 30°C for 15 minutes before mixing.
2. Mix the hematology control by holding it horizontally between the palms of the hands. Roll it back and forth for 20 to 30 seconds and occasionally turn it upside down. Mix the hematology control until the red cells are completely suspended. Mix thoroughly but do not shake it.
3. Measure the hematology control as instructed in the hematology analyzer operator's manual. Immediately before measuring, gently turn the hematology control upside down 8 to 10 times.
4. Return the hematology control to the refrigerator within 30 minutes of use.

NOTE

- Incomplete mixing of the hematology control affects the measurement result.
- Dark red supernatant fluid, discoloration of the product or unacceptable results may indicate deterioration. Do not use the hematology control if deterioration is suspected.

Expected Results

Refer to the assay sheet.

Storage

Store at 2 to 8°C.

Expiration Date

The expiration date is shown on the container label.

For this user's guide in your local language and assay sheet, contact your Nihon Kohden distributor.

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