3.1

LISS/Coombs

ID-Card

English

B004014 02.13

Indirect and direct antiglobulin test

Product-Identification: 50531

INTRODUCTION

Polyspecific anti-human globulin (AHG) reagents are used for routine alloantibody detection and identification, compatibility tests and the direct antiglobulin test (DAT).

The most important function of the polyspecific AHG reagent is to detect the presence of IgG. The importance of anticomplement in the AHG reagent is debatable since antibodies detectable only by their ability to bind complement are rather rare. However, anti-C3d activity is important for the DAT in the investigation of autoimmune hemolytic anemia (AIHA) [4]. A positive DAT generally indicates that the red cells are coated *in vivo* with immunoglobulin and/or complement.

The microtubes of the ID-Card "LISS/Coombs" contain polyspecific AHG, to be used for antibody screening, antibody identification, crossmatch and the DAT. For the indirect antiglobulin test (IAT), labour intensive washing procedures are eliminated, due to the fact that the red cell suspension is added to the microtube before the plasma/serum, creating a barrier over the gel suspension, thus avoiding neutralization of the AHG by serum IgG proteins.

The anti-human globulin IgG used in the ID-Card "LISS/Coombs" is not heavy chain specific and thus may also be capable of reacting with the Kappa (k) and Lambda (l) light chains of IgA and IgM molecules.

The ID-Card "LISS/Coombs" is suitable for the DAT, for the compatibility test, for antibody screening and identification with "ID-DiaCell" and "ID-DiaPanel".


BIO-RAD

REAGENTS**IVD**

ID-Card "LISS/Coombs" with 6 microtubes containing polyspecific AHG (rabbit anti-IgG and monoclonal anti-C3d, cell line C139-9) within the gel matrix.

Preservative: < 0.1% NaN₃.

laikymo temperatūra

*Do not store near any heat, air conditioning sources or ventilation outlets.
Stability: see expiry date on label.*

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspension.
 - ID-DiaCell, ID-DiaPanel: Test cell reagents.
- (see related package insert)

FURTHER MATERIALS REQUIRED

- ID-Dispenser
- ID-Pipetor
- ID-Tips (pipetor tips)
- Suspension Tubes
- ID-Working table
- ID-Incubator 37 °C
- ID-Centrifuge 6, 12 or 24

SAMPLE MATERIAL

For optimal results, the determination should be performed using a freshly drawn sample, or in accordance with local laboratory procedures for sample acceptance criteria. Preferably, blood samples should be drawn into citrate, EDTA or CPD-A anticoagulant. Samples drawn into plain tubes (no anticoagulant) may also be used.

When the use of serum instead of plasma is required, the serum must be well cleared, by centrifugation at 1500 g for 10 minutes, before use avoid fibrin residues, which may interfere with the reaction pattern.

PREPARATION OF BLOOD SAMPLE

a) Red cell suspension for DAT or autocontrol

Prepare a 0.8% red cell suspension in ID-Diluent 2 as follows:

Allow the diluent to reach room temperature before use.

1. Dispense 1.0 ml of ID-Diluent 2 into a clean tube.
2. Add 10 µl of packed red cells, mix gently.

The cell suspension may be used immediately.

b) Red cell suspension for crossmatch procedures

Prepare a 0.8% red cell suspension in ID-Diluent 2 as above. Where whole blood directly from a segment of the blood bag is used, add 20 µl of blood to the 1.0 ml of ID-Diluent 2.

c) Plasma or serum for indirect antiglobulin test (IAT) procedures

Where samples are not for immediate testing they should be stored at 2–8 °C after separation (see also under "Sample material") for a maximum of 48 hours, thereafter at -20 °C or in accordance with local/national policies/guidelines.

TEST PROCEDURES

Do not use ID-Cards which show signs of drying, have bubbles, damaged seals, drops of gel or supernatant in the upper part of the microtubes or on the underside of the aluminium foil.

I. Direct antiglobulin test (DAT)

1. Identify the appropriate microtubes of the ID-Card "LISS/Coombs" with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.

The logo for BIO-RAD, consisting of the words "BIO-RAD" in white, uppercase letters inside a green rounded rectangular border.

3. Pipette 50 µl of the red cell suspension to the appropriate microtube.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

II. Antibody screening (IAT)

Use ready-to-use test cell reagents "ID-DiaCell".

Allow the test cell reagents and samples to reach room temperature before use.

1. Identify the appropriate microtubes of the ID-Card "LISS/Coombs" with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each test cell reagent to the appropriate microtubes (marked with the corresponding test cell).
4. When an autocontrol is to be included, pipette 50 µl of the sample's own red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's or donor's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

III. Antibody Identification (IAT)

Use the ready-to-use test cell reagent "ID-DiaPanel".

Allow the test cell reagents and samples to reach room temperature before use.

1. Identify two ID-Cards "LISS/Coombs" with the patient's or donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of each "ID-DiaPanel" test cell to the appropriate microtube (marked 1 to 11).
4. Pipette 50 µl of the sample's own red cell suspension to the 12th microtube (autocontrol).

5. Add 25 µl of the patient's or donor's plasma or serum to all 12 microtubes.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

IV. Compatibility test

1. Identify the appropriate microtubes of the ID-Card "LISS/Coombs" with recipient's and donor's name or number.
2. Remove the aluminium foil from as many microtubes as required by holding the ID-Card in the upright position.
3. Pipette 50 µl of the donor red cell suspensions to the appropriate microtubes.
4. For the autocontrol, pipette 50 µl of the patient's own red cell suspension to the appropriate microtube.
5. Add 25 µl of the patient's plasma or serum to each microtube.
6. Incubate the ID-Card for 15 minutes at 37 °C in the ID-Incubator.
7. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
8. Read and record the results.

INTERPRETATION OF THE RESULTS

A) Principle

Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.

Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for:

I. Direct antiglobulin test (DAT)

- A negative reaction indicates absence of detectable IgG antibodies or C3d complement component on the red cells.
- A positive reaction (± to ++++) indicates that the patient's red cells are sensitized (red cells coated with IgG antibodies and/or C3d).

The logo for BIO-RAD, consisting of the words "BIO-RAD" in white, uppercase letters inside a green rounded rectangular border.

II. Antibody screening

- A negative reaction indicates the absence of detectable irregular antibodies in the patient's or donor's serum or plasma.
- A positive reaction indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents "ID-DiaCell I-II" or "ID-DiaCell I-II-III" corresponds to the lot number indicated on the antigen table.
- Following the reaction pattern and the antigen configuration, the type of antibody present may be indicated. Perform the usual further tests to identify the antibody.
- A positive reaction with one or more test cells and a negative autocontrol suggest the presence of a specific antibody (see Remarks).
- A positive reaction with all test cells and a positive autocontrol may be due to non-specific reactions.
- A positive reaction with all test cells and a positive autocontrol but with one or more test cells showing a stronger positive reaction than the autocontrol, the patient sample should be submitted for further testing, to investigate the possibility of an underlying allo-antibody.

III. Antibody Identification

- A positive reaction indicates the presence of irregular antibodies. Enter the reactions obtained on the antigen table. Verify that the lot number of the test cell reagents "ID-DiaPanel" corresponds to the lot number indicated on the antigen table.
 - Following the reaction pattern and the antigen configuration, the type of antibody present can, in most cases, be identified (autocontrol must be negative).
 - A positive reaction with all "ID-DiaPanel" test cells and a negative autocontrol may be due to non-specific reactions or may indicate the presence of an alloantibody directed against a high frequency antigen.
 - A positive reaction with all "ID-DiaPanel" test cells and a positive autocontrol may be due to non-specific reactions.
-

- A positive reaction with all "ID-DiaPanel" test cells and the autocontrol but with one or more test cells showing stronger reactions than the autocontrol, may indicate an underlying allo-antibody and further investigation should be undertaken.

IV. Compatibility test

- A negative reaction indicates compatibility of the donor blood with the recipient.
- A positive reaction indicates incompatibility of the donor blood with the recipient, due to presence of antibodies directed against antigens on the donor red cells. Further investigation to identify the antibody specificity should be performed.

REMARKS

As for all procedures covered by GLP rules, the sensitivity of the above procedures should be validated using antibodies of known potency. The Bio-Rad Anti-D reference reagent provides the means to regularly perform controls for all antibody detection procedures (see related package insert).

LIMITATIONS

- a) ID-Cards which show air bubbles or gel drops in the upper part of the microtubes and/or the seal, must be centrifuged before use.
- b) Certain drugs are known to cause positive reactions in anti-human globulin procedures.
- c) Some pathological conditions are also reported as causing positive reactions in anti-human globulin procedures.
- d) Cells that have become polyagglutinable, due to cryptantigen exposure e.g. T antigen, either *in vivo* or *in vitro* may react with all human sera. Further investigation of such reactions is required.
- e) Bacterial or other contamination of materials used can cause false positive or false negative results.

- f) Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
- g) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- h) Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
- i) The use of other test cell reagents than "ID-DiaCell" or "ID-DiaPanel" may modify the reaction patterns.
- j) Too heavy or too weak red cell suspensions can cause aberrant results.

PERFORMANCE CHARACTERISTICS

Specificity/sensitivity

Performance evaluation studies have been realized with ID-Cards "LISS/Coombs". Results obtained showed that performance characteristics of each application of the card are in conformity with the application field and comparable to other approved products.

Reproducibility

Intra-assay reproducibility (repeatability) and inter-assay reproducibility of the ID-Cards "LISS/Coombs" have been evaluated internally. Neither false positive nor false negative results were observed. Differences between reactions in positive samples were less than one reaction strength.

BIBLIOGRAPHY

1. Coombs, R.R.A., Mourant, A.E. and Race, R.R.: Lancet 1945; 2: 15.
2. Coombs, R.R.A., Mourant, A.E. and Race, R.R.: Brit. J. Exp. Path 1945; 26: 255.
3. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.
4. Technical Manual; Roback, J.D. (ed.); 17th ed. 2012; American Association of Blood Banks.

PRODUCTS

ID-Card "LISS/Coombs"

4 x 12.	REF 004014
24 x 12.	REF 004017
60 x 12.	REF 004016
112 x 12.	REF 004015

These products are guaranteed to perform as described on the label and in the instruction sheet. The manufacturer declines all responsibility arising out of the use or sale of these products in any way or for any purpose other than those described therein.

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Name:	Function:	Issued in:	Date:	Signature
Galéa, Diane	Site Quality Management Representative	Cressier FR	14.10.2019	

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English

B001254 06.13

A, B, DVI-/A, B, DVI-**monoclonal****Confirmation of blood groups ABO and RhD****Product-Identification: 50053****INTRODUCTION**

Compatibility testing must include the examination of previous transfusion records for the recipient's serological history [1]. If previous testing was done, results of current testing must be compared with the interpretation of the previously determined blood group. Concurrence between both test results gives some assurance that there have been no identification errors and that the tests have been performed and interpreted correctly.

Current recommendations for RhD typing suggest that for transfusion recipients and antenatal patients, anti-D reagents should not detect the DVI phenotype. The D antigen consists of many epitopes, and in the 9 epitope model, DVI red cells lack all but three. This means that individuals possessing the DVI phenotype may produce an anti-D to the missing epitopes after immunisation by fetal or transfused RhD positive red cells. To ensure that appropriate therapeutic measures are instigated, DVI patients' red cells should be assigned 'Rh negative' status. Conversely, donor bloods should be tested with anti-D that does detect DVI and assigned 'Rh positive' status, to avoid the unit being transfused to an RhD negative or partial-D patient.

The ABD confirmation card can be used for the ABO/RhD blood group control of patients.

REAGENTS

ID-Card "DiaClon ABD-Confirmation for Patients" contains monoclonal anti-A [cell line LM 297 / 628 (LA-2)], anti-B [cell line LM 306 / 686 (LB-2)] and anti-D [cell lines TH-28, RUM-1, LDM1] within the gel matrix.
Preservative: < 0.1% Na₂S₂O₃.

Caution: All reagents should be treated as potentially infectious.

Do not store near any heat, air conditioning sources or ventilation outlets.

Stability: see expiry date on label.

laikymo temperatūra

ADDITIONAL REAGENTS REQUIRED

- ID-Diluent 2: modified LISS for red cell suspensions.

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PREPARATION OF BLOOD SAMPLE

Prepare a 5% red cell suspension in ID-Diluent 2 as follows:
Allow the diluent to reach room temperature before use.

1. Dispense 0.5 ml of ID-Diluent 2 into a clean tube.
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.

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CONTROLS

Known positive and negative samples should be included in accordance with the relevant guidelines of quality assurance.

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1. Identify the ID-Card with the unique patient or donor number / details as appropriate.
2. Remove the aluminium foil from as many microtubes as required by holding the ID card in the upright position.
3. Add 10 or 12.5 µl of the red cell suspension of the first sample to the first 3 microtubes of the ID-Card.
4. Centrifuge the ID-Card for 10 minutes in the ID-Centrifuge.
5. Read and record the results.

Note: If only one patient is to be tested, remove the aluminium foil from the first 3 microtubes. The remaining microtubes may be used at a later date for a subsequent sample.

INTERPRETATION OF THE RESULTS

A) Principle [2]

- Positive: Agglutinated cells forming a red line on the surface of the gel or agglutinates dispersed in the gel.
 Negative: Compact button of cells on the bottom of the microtube.

B) Reactions for blood groups ABO

Anti-A	Anti-B	Blood group
+++ to ++++	negative	A
negative	+++ to ++++	B
+++ to ++++	+++ to ++++	AB
negative	negative	O

Weaker reactions than +++ may indicate A or B subgroups. For correct interpretation, a complete grouping test should be performed (anti-A, anti-B, anti-AB). In the presence of weak or very weakly expressed antigens the reaction can be negative. The anti-B of monoclonal origin does not react with the acquired B antigen.

C) Reactions for RhD

+++ to ++++	± to ++*	negative
RhD positive	RhD weak positive	RhD negative

* ±, trace or weak reactions should be subject to further investigations to distinguish between weak and partial D types as appropriate for the category of sample being tested.

The anti-D test sera were selected so as **not to react** with DVI variants.

Note that most guidelines do not recommend further testing for weak or partial-D in patients.

LIMITATIONS

- ID-Cards which show air bubbles in the gel or drops in the upper part of the microtubes and / or the seal, must be centrifuged before use.
- ABD-confirmation does not replace the complete ABO/Rh determination. It should only be used as a confirmation of previously determined blood groups where, according to general rules, ABO-typing is confirmed with reverse grouping and ABO/Rh-typing is validated by the presence of a negative control test.
- Bacterial or other contamination of materials used can cause false positive or false negative results.
- Fibrin residues in the red cell suspension may trap non-agglutinated cells presenting a fine pink line on top of the gel while most of the cells are on the bottom of the microtube after centrifugation.
- Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- Use of suspension solutions other than ID-Diluent 2 may modify the reactions.
- Too heavy or too weak red cell suspensions can cause aberrant results.

BIBLIOGRAPHY

- Holland P.V. ed.; Standards for Bloodbanks and Transfusion Service, 13th ed., Arlington,VA: American Association of Blood Banks, 1984;
- Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reactions. Transfusion 1990; 30: 109-113.

PRODUCTS

ID-Card "DiaClon ABD-Confirmation for Patients"

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Changes to the version 07.11 are shaded grey.

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Name:	Function:	Issued in:	Date:	Signature
Galéa Diane	Site Quality Management Representative	Cressier FR	29.03.2017	

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English 3.3

B009290 03.13

Modified LISS for red cell suspensions

Product-Identification: 05761

INTRODUCTION

Low ionic strength solution (LISS) increases the rate of antibody association and thus enhances antigen/antibody reactions.

1D-Diluent 2[®] is a modified low ionic strength solution made for the ID-System, for preparing 5% red cell suspensions for blood grouping as well as 0.8% red cell suspensions for crossmatching, autocontrol, direct anti-human globulin test, blood grouping of newborns and test cells prepared in the laboratory.


REAGENTS

1D-Diluent 2[®]: modified LISS for Red Cell Suspension, in 100 and 500 ml bottles or racks of 60 x 700 µl for the IH-Analyzers.

Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

Once opened and if handled in accordance with good laboratory practice (GLP) principles and stored correctly, the 100 and 500 ml bottles can be used until the date of expiry or for up to 6 months afterwards, whichever is the sooner.

Stability: see expiry date on label.

Laikymo temperatūra

FURTHER MATERIALS REQUIRED (FOR 100 AND 500 ML BOTTLES)

- ID-Dispenser
- Pipette 25 and 50 µl
- ID-Tips (pipetor tips)
- Suspension Tubes
- ID-Working Table

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PREPARATION OF BLOOD SAMPLE (FOR 100 AND 500 ML BOTTLES)

A) For blood group determinations

Prepare a 5% red cell suspension of the patient's red cells in 1D-Diluent 2[®] as follows:

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1. Dispense 0.5 ml of 1D-Diluent 2[®] into a clean tube
2. Add 50 µl of whole blood or 25 µl of packed cells, mix gently.

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B) For crossmatching, direct anti-human globulin test, autocontrol in antibody screening and other tests as specified in the relevant box inserts

Prepare a 0.8% red cell suspension in 1D-Diluent 2[®] as follows:

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1. Dispense 1.0 ml of 1D-Diluent 2[®] into a clean tube.
2. Add 10 µl of packed red cells, mix gently.

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C) For preparation of test cell reagents in the laboratory

1. Wash the red cells with isotonic saline solution or with 1D-Diluent 2[®] three times or more, until the supernatant is clear.
2. Decant the supernatant and resuspend the packed cells to a 0.8% suspension with 1D-Diluent 2[®].
3. Under aseptic conditions place the suspension into sterile glass vials or tubes, close with sterile caps.

Red cell reagents thus prepared and stored at 2–8 °C are stable for one day.

Important: consult the related instruction sheet of the ID-Cards for precise working procedures.

PREPARATION OF BLOOD SAMPLE (RACKS WITH 60 X 700 µL)

Please refer to the user manual of the IH-Analyzers.

LIMITATIONS

- a) Bacterial or other contamination of materials used can cause false positive or false negative results.
- b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.
- c) "ID-Diluent 2" has been formulated in order to reduce the incidence of detection of non-significant red cell bound complement components.
- d) Some antibodies in the Kell system may react more weakly in LISS techniques.
- e) LISS solutions may affect autoantibody activity and therefore may cause problems with certain samples.

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Changes to the version 08.12 are shaded grey.

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0123

Declaration of conformity

MANUFACTURER: **DiaMed GmbH**
ADDRESS: **Pra Rond 23**
1785 Cressier FR
Switzerland
Phone: +41(0)26 674 51 11

PRODUCT NAME: **ID-Diluent 2**
Id-n°: **05761**
REF: **009260 / 009280 / 009290**

We hereby declare that the above mentioned product meets the provisions of the following Directives:

APPLICABLE DIRECTIVE: Directive 98/79/EC of the European Parliament and of the Council of 27 October 1998 on in vitro Diagnostic medical devices

CLASSIFICATION: Others

CONFORMITY ROUTE: Annex III

GMDN Code: 52718

Generic Device Group Term: Low ionic strength saline (LISS) immunohaematology reagent IVD

NOTIFIED BODY: /

Name:	Function:	Issued in:	Date:	Signature
Galéa Diane	Site Quality Management Representative	Cressier FR	01.02.2017	

DiaMed GmbH
Pra Rond 23
1785 Cressier FR / Šveicarija
Telefonas: +41 (0)26 674 51 11
Faksas: +41 (0)26 674 54 45

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Vaidas Vilmantas
MB „Beikeris“. I. k. 304539005

Test Cell Reagents for the ID-System

English 3.4

B004350 02.15

Test cell reagents for antibody screening: ID-DiaCell I-II, ID-DiaCell I-II-III, ID-DiaCell IP-IIP-IIIP (papainized), ID-DiaCell I-II-III Asia, ID-DiaCell Pool; for antibody identification: ID-Panel, ID-Panel-P; special antigens: ID-Di^a (Diego) Positive, ID-I Negative Cell

INTRODUCTION

The reliability of antibody detection is largely dependent on the availability of test cells with appropriate antigens and on the sensitivity of the test methods used.

The requirements for antigen configuration are stringent: it must allow the safe detection of all clinically significant antibodies. For the Rh system, MNSs, Duffy and Kidd, the antigens must be in homozygous form. The Lewis antigens must be present, as should the rare antigen Kp^a.

It is generally considered most effective to perform screening tests by both anti-human globulin (AHG) and enzyme test procedures. Due to higher sensitivity of the indirect antiglobulin test (IAT) with procedures such as the ID-System, some scientists in various countries have formed the opinion that the enzyme test has become somewhat less important.

However, enzyme techniques are useful when increased sensitivity in antibody screening is desired or where more than one antibody may be present. They enhance the reactions of certain antibodies, notably in the Rh, Kell and Kidd systems, whereas antibodies to enzyme-sensitive antigens may not be detected, notably in the Duffy and MNS Systems.

The test cell reagents are specially designed for the ID-System.

REAGENTS

All test cell reagents are of human origin, in a buffered suspension medium at 0.8% (\pm 0.1%).
Preservatives: the antibiotics trimethoprim and sulfamethoxazole.

For antibody screening, single donors, blood group O:

ID-DiaCell I-II	$R_1^w R_1 + R_2 R_2$ for IAT and NaCl test
ID-DiaCell I-II-III	$R_1^w R_1 + R_2 R_2 + rr$ for IAT and NaCl test
ID-DiaCell IP-IIP-IIIP	papainized, for enzyme technique
ID-DiaCell Pool	$R_1 R_1 + R_2 R_2$ (2 pooled cells for donor screening)
ID-DiaCell I-II-III Asia	$R_1 R_1 + R_2 R_2 + cell$ of the GP.MUR phenotype, for IAT and NaCl test

For antibody identification, single donors, blood group O:

ID-DiaPanel	11 test cells for IAT and NaCl test
ID-DiaPanel-P	11 test cells papainized, for enzyme technique

Special antigens:

These cells are for use as an adjunct with other routine antibody screening cell sets.

ID-Di^a (Diego) positive
ID-I negative cell

Shipment on standing order, every 4 weeks.

Caution: The source materials from which these products were manufactured, were found non-reactive for HBsAg, HCV and HIV (1+2) when tested with licensed reagents. However, no known test method can assure that infectious agents are absent. Products from human blood should be considered potentially infectious.

laikymo temperatūra

ADDITIONAL REAGENTS REQUIRED

- ID-Card "LISS Coombs + Enzyme Test" 3 microtubes with polyspecific anti-human globulin AHG) serum and 3 microtubes containing neutral gel (Id-no: 50581).
- ID-Card "LISS/Coombs" 6 microtubes with polyspecific AHG serum (Id-no: 50531).
- ID-Card "Coombs Anti-IgG" 6 microtubes containing rabbit anti-IgG (Id-no: 50540).
- ID-Card "NaCl, Enzyme test and cold agglutinins" 6 microtubes containing neutral gel (Id-no: 50520).
- ID-Card "Reverse Grouping with Antibody Screening" 3 microtubes containing neutral gel and 3 microtubes containing polyspecific anti-human globulin (AHG) serum (Id-no: 50510).
- ID-Diluent 2: modified LISS for red cell suspensions (Id-no: 05761).

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CONTROLS

Controls should be included in accordance with the relevant guidelines of quality assurance.

USE OF THE ID-TEST CELL REAGENTS

- All test cell reagents are for use with the ID-Cards of the ID-System only.
- **Strictly follow the test procedures as described in the specific package inserts of the ID-Cards to be used.**
- Always gently resuspend the red cells, by inverting the vial several times before use and also before placing the vials into a pipetting automate.
- **Make sure that the test cells are at room temperature (18–25 °C) when in use. Atlikimo/laikymo temperatūra**
- During the working procedures, check that the test cell reagents remain in suspension. If there is settling of the cells, resuspend again.
- **For the ID-System, precise pipetting is of importance.** Use the ID-Pipetors for serial pipetting.
- Avoid contamination of the test cell reagents.
- When recording the reactions, ensure that the lot number of the antigen table corresponds with the lot number of the reagent vials.
- After use, close the vials and replace them in the refrigerator.

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Note: Test cell reagents for antibody screening containing pooled cells can show a double cell population appearance depending on the antibody present. This is considered to be a positive result.

B) Reactions

see package inserts of the corresponding ID-Cards

LIMITATIONS

- a) Bacterial or other contamination of materials used can cause false positive or false negative results.
- b) Strict adherence to the procedures and recommended equipment is essential. The equipment should be checked regularly according to GLP procedures.

BIBLIOGRAPHY

1. Technical Manual of the American Association of Blood Banks, 13th edition, 1999.
2. Lapierre, Y., Rigal, D., Adam, J. et al.: The gel test; A new way to detect red cell antigen-antibody reaction. Transfusion 1990;30:109-113.

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Changes to the version 02.13 are shaded gray.

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LITERATŪRA

1. Technical Manual of the American Association of Blood Banks, 13th edition, 1999.
2. Lapiere, Y., Rigal, D., Adam, J. Et al.: The gel test; A new way to detect red cell antigen-antibody reaction. Transfusion 1990;30:109-113.

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CONFORMITY ROUTE: Annex IV

GMDN Code: 52544

Generic Device Group Term: Antibody screening reagent red blood cell IVD, kit, agglutination

NOTIFIED BODY: TÜV Product Service GmbH
Ridlerstrasse 65-80339 München (Germany)
CE-N° 0123

Name:	Function:	Issued in:	Date:	Signature
Galéa Diane	Site Quality Management Representative	Cressier FR	31.01.2017	

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Gamintojas: **DiaMed GmbH**
Patvirtinimo agentūra: TUV Product service GmbH