DIRECTIONS FOR USE  
ESTABLISHMENT LABS S.A.  
STERILE SILICONE BREAST IMPLANTS  
MOTIVA IMPLANT MATRIX®

Directions For Use  
Sterile Silicone Breast Implants  
Motiva Implant Matrix®  
Establishment Labs S.A.

INTRODUCTION

The aim of this Product Insert is to provide an overview of essential information about Establishment Labs S.A. Sterile Silicone Breast Implants Motiva Implant Matrix®, including device description, indications for use, contraindications, warnings, precautions, relevant topics that must be discussed with the patient, adverse events, other reported conditions, returned goods policy, product evaluation, warranty and medical device reporting.


DEVICE DESCRIPTION

Establishment Labs S.A. Sterile Silicone Breast Implants Motiva Implant Matrix® are mammary augmentation/ reconstruction devices constructed of successive cross-linked layers of silicone elastomer and a low diffusion barrier shell technology that provide these implants their elasticity and integrity. All implants are composed of the above-described shell, a patch, and silicone gel fill. The shell is filled with ProgressiveGel™, ProgressiveGel™ PLUS or ProgressiveGel™ Ultima™, an Establishment Labs S.A. silicone gel proprietary formula. All raw materials are supplied by an FDA approved silicone materials formula in the USA.

The following are the reference ranges of Establishment Labs S.A. Sterile Silicone Breast Implants Motiva Implant Matrix®:

Motiva Implant Matrix® - Round								
Base (cm)	MINI		DEMI		FULL		CORSE	
	P(cm)	V (cc)	P(cm)	V (cc)	P (cm)	V (cc)	P (cm)	V (cc)
8.5	2.2	105	3.1	135	3.5	145	4.0	180
9	2.3	125	3.3	155	3.7	175	4.2	210
9.5	2.4	140	3.4	180	3.9	205	4.5	240
9.75	2.4	150	3.4	190	4.0	220	4.6	260
10	2.5	160	3.5	205	4.1	235	4.8	280
10.25	2.5	170	3.5	215	4.2	255	4.9	300
10.5	2.6	185	3.6	230	4.3	275	5.1	325
10.75	2.6	205	3.7	245	4.4	295	5.2	350
11	2.7	220	3.8	265	4.5	315	5.4	380
11.25	2.7	230	3.8	285	4.6	335	5.5	410
11.5	2.8	245	3.9	300	4.7	355	5.7	440
11.75	2.8	260	3.9	320	4.8	375	5.8	475
12	2.9	275	4.0	340	4.9	400	6.0	510
12.25	2.9	290	4.0	360	5.0	425	6.1	550
12.5	3.0	310	4.1	380	5.1	450	6.3	590
13	3.1	360	4.3	425	5.3	500	6.6	650
13.5	3.2	400	4.4	475	5.5	550	6.9	725
14	3.3	430	4.5	525	5.7	625	7.2	825
14.5	3.4	475	4.6	575	5.9	700	7.5	925
15	3.5	525	4.8	625	6.1	775	7.8	1050

P=Projection  
V=Volume

INDICATIONS

Establishment Labs S.A. Sterile Silicone Breast Implants Motiva Implant Matrix® are indicated for the following procedures in female patients:

- **Breast augmentation for women of at least 18 years of age**, including previous augmentation to increase the breast size and revision surgery to correct or improve the result of a previous breast augmentation surgery.
- **Breast Reconstruction.** Breast reconstruction, including previous reconstruction to replace breast tissue that has been removed due to cancer or trauma, or that has failed to develop properly due to a severe breast anomaly, as well as revision surgery to correct or improve the results of a previous breast reconstruction surgery.

CONTRAINDICATIONS

Breast augmentation with implants is contraindicated in:

- Women with existing carcinoma of the breast, without mastectomy.
- Women with advanced fibrocystic disease considered premalignant, without accompanying subcutaneous mastectomy.
- Women with active infections.
- Women who are currently pregnant or nursing.
- Women with any disease, including uncontrolled diabetes, which is clinically known to impact wound-healing ability.
- Women who show tissue characteristics clinically incompatible with mammoplasty, such as tissue damage resulting from radiation, inadequate tissue, compromised vascularity or ulceration.
- Women with any condition – or treatment – determined by the surgeon to constitute an unjustifiable surgical risk (e.g., unstable cardiovascular disease, coagulopathies, chronic pulmonary problems, etc).

WARNINGS

Care during surgical insertion and subsequent procedures:

- Do not allow sharp instruments, such as scalpels or needles, to come in contact with the device during the implantation or other surgical procedures. Patients should be instructed to inform other treating physicians to also observe this warning.
- Do not immerse the implant in iodine solution. If iodine solution is used in the pocket, make sure that it is rinsed thoroughly with deionized water so that no residual solution remains in the pocket.
- Do not allow the implant to come in contact with cauterization devices.
- Do not alter the implant or attempt to repair or insert a damaged implant.
- Ensure that excessive force is not applied to a very small area of the shell during insertion of the device through the incision. Instead, apply force over as large an area of the implant as possible during insertion.
- The incision should be of appropriate length to accommodate the volume and profile of the implant. This will reduce the potential for creating excessive stress to the implant when inserting it. Forcing implants through too small an opening may result in local weakening of the breast implant shell, potentially leading to shell damage and possible implant rupture.
- Periareolar and axillary incision sites may make the insertion more difficult,

increasing the risk of damage to the implant. Periareolar incision may considerably reduce the possibility of future breast-feeding.

- Do not use the periumbilical approach to place the implant.
- Avoid creating wrinkles or folds in the device during the insertion. It is recommended to run a finger around the implant before closing to make sure the implant is flat.
- Do not place more than one implant per breast pocket.
- Do not treat capsular contracture by closed capsulotomy or forceful external compression, which will likely produce implant damage, rupture, folds, and/or hematoma.
- Procedures such as open capsulotomy, breast pocket revision, hematoma/seroma aspiration, biopsy, and lumpectomy might result in damage to the implant shell, so they must be carefully performed. Care should be taken when re-positioning the implant during subsequent procedures to avoid contamination of the implant. Use of excessive force during any subsequent procedure can contribute to local weakening of the breast implant shell, potentially leading to shell damage and possible implant rupture.
- Do not re-use or re-sterilize any product that has been previously implanted. Breast implants are intended for single use only.
- Do not use microwave diathermy in patients with breast implants, as it has been associated with tissue necrosis, skin erosion, and implant extrusion.

## PRECAUTIONS

### 1. Specific Populations

Safety and effectiveness of breast augmentation surgery have not been established for the following populations and/or conditions:

- Patients with autoimmune diseases (e.g., lupus, scleroderma).
- Patients whose immune system is compromised (e.g., currently receiving immunosuppressive therapy such as steroids).
- Patients with conditions or medications that may interfere with wound healing ability (e.g., poorly controlled diabetes, or corticosteroid therapy) or blood clotting (e.g., concomitant Warfarin therapy).
- Patients with reduced blood supply to breast or overlying tissue.
- Patients undergoing radiation therapy.
- Women with ptotic breasts where nipples fall below the inframammary fold, without concurrent mastopexia.
- Previous repeated contour correction failures.
- Patients with clinical diagnosis of depression or other mental health disorders, including BDD (body dysmorphic disorder) and eating disorders. The patient should be advised to discuss any history of mental health disorders with her surgeon prior to surgery. Patients with a diagnosis of depression, or other mental health disorders, should wait until stabilization of these conditions prior to undergoing breast implantation surgery.
- There may be other patients with complicated medical histories, who are judged to present risk factors that might interfere with the safety and effectiveness of the breast implantation. As with any surgical procedure, the patient's medical history should be carefully reviewed to ensure that she is an appropriate candidate for breast implant surgery.

### 2. Surgical Precautions

**Preliminary product examination**— Immediately before insertion, examine the device by gently manipulating it while carefully checking for rupture, leakage sites or particulate contamination.

**Surgical technique and implant selection**— There are several surgical techniques that can be used to perform the implantation of a silicone gel-

filled breast implant. Therefore, the surgeon is advised to use his/her clinical judgment in choosing the procedure that is best for the patient, consistent with this product insert. After setting realistic aesthetic goals that assure mutual understanding between doctor and patient, the surgeon must choose from current and accepted surgical techniques to minimize the incidence of adverse reactions and achieve the best results.

The implant size should be consistent with the patient's chest wall dimensions, including base width measurements, characteristics of the tissue and projection of the implant.

Textured implants, larger implants, subglandular placement, and an insufficient amount of tissue available to cover the implant may cause implants to be more palpable.

Implants of larger sizes may increase the risk of complications such as extrusion, hematoma, infection, palpable implant folds, and visible skin wrinkling.

**CAUTION:** Only surgeons with qualified training and certified by the corresponding national medical board of their countries should use this product. The use of this product by unqualified practitioners may result in extremely poor aesthetic outcomes and serious adverse effects.

## RELEVANT TOPICS THAT MUST BE DISCUSSED WITH THE PATIENT

### Patient Counseling Information

This document and the patient information booklet must be carefully reviewed before counseling a patient about Establishment Labs S.A. Sterile Silicone Breast Implants Motiva Implant Matrix® and breast augmentation surgery. Doctors must thoroughly read and understand the contents of this document and make sure that any questions or concerns have been resolved prior to proceeding with use of the device. Breast implantation is an elective procedure, and the patient must understand its potential risks and benefits in order to make an informed decision. For this reason, the patient should be instructed to read the document called **“Motiva Implant Matrix®: Information for the Patient”**. The doctor must discuss with the patient warnings, contraindications, precautions, important factors to consider, complications, and all other aspects of the document. The physician should inform the patient about the potential complications and that medical management of serious complications may include additional surgery and explantation.

### Informed Consent

Each patient should be given the Establishment Labs S.A. **“Motiva Implant Matrix®: Information for the Patient”** during her surgical consultation. It is the surgeons' responsibility to ensure this happens and it is a requirement for the use of the device. The patient must be given enough time to read and completely understand the information regarding the risks, benefits and recommendations associated with silicone gel-filled breast implant surgery.

In order to document a successful informed decision process, the patient, a witness and the surgeon should sign the **“Informed Consent Document”**, which will be part of the patient's medical file.

Some of the relevant topics patients need to be aware of when considering the use of silicone gel-filled breast implants are:

**Rupture**— Breast implants can rupture when the shell develops a tear or hole. Rupture can occur at any time after implantation, but it is more likely to happen when the implant has been in place for a long time. Rupture of a silicone gel-filled breast implant is most often silent (the patient does not experience any symptoms and there are no physical signs of changes with the implant) rather than symptomatic. Therefore, patients should be advised

to have regular MRIs over their lifetime to screen for silent rupture even if they are not having any apparent problems. The first MRI should be performed at 3 years postoperatively, then regularly at 2-year intervals, and submitted to the treating surgeon. Patients should be provided with a list of radiology centers with experience with breast implant MRI films to scan for signs of rupture. The importance of these MRI evaluations should be emphasized. If rupture is noted on an MRI, the patient should be strongly encouraged to have her implant removed.

**Mammography**– Patients should be advised to have routine mammography exams performed according to their surgeon's recommendations. The importance of these exams should be emphasized. Patients should be instructed to inform their examiners about the presence, type, and placement of their implants, and to request a diagnostic mammography, rather than a screening mammography. Breast implants may complicate the interpretation of mammographic images by obscuring underlying breast tissue and/or by compressing overlying tissue. Accredited mammography centers, technicians with experience in examining patients with breast implants, and the use of displacement techniques are needed to adequately visualize breast tissue in the implanted breast. The current recommendations for preoperative/screening mammograms are no different for women with breast implants than for those women without implants. Pre and post-surgical mammographies may be performed to determine a baseline for routine future studies in augmentation patients.

**Explantation**– Implants are not lifetime devices, and there is a possibility that patients will undergo implant removal(s), with or without replacement, over the course of their life. When implants are explanted without replacement, changes to the patient's breasts may be irreversible. Complication rates are higher following revision surgery (removal with replacement).

**Reoperation**– Rupture, unacceptable cosmetic outcomes (dimpling, wrinkling, and other potentially permanent cosmetic changes of the breast) and other complications may require additional surgeries to the patients' breast. Patients should be advised that their risk of future complications increases with revision surgery as compared to primary augmentation or reconstruction surgery. For example, the risk of severe capsular contracture doubles for both augmentation and reconstruction patients with implant replacement compared to first time implantation. There is a risk of an accidental compromise of implant shell integrity during reoperation, potentially leading to product failure.

**Infection**– Signs of acute infection reported in association with breast implants include edema, erythema, tenderness, pain, and fever. As with other invasive surgeries, Toxic Shock Syndrome (TSS), a life-threatening condition, has been reported in rare instances following breast implant surgery. Symptoms of TSS occur suddenly and can include high fever (102°F, 38.8°C or higher), vomiting, diarrhea, sunburn-like rash, red eyes, dizziness, lightheadedness, muscle aches and drops in blood pressure, which may cause fainting. Patients should immediately contact their physician for diagnosis and treatment if any of these symptoms occur.

**Breast Examination Techniques**– Patients should perform breast self-examinations monthly and be shown how to distinguish the implant from their breast tissue. The patient should not manipulate or squeeze the implant excessively. The patient should be told that the presence of lumps, persistent pain, swelling, hardening, or change in the implant shape could suggest

symptomatic rupture of the implant. If the patient has any of these signs, she should be advised to report them, and possibly have an MRI evaluation to screen for rupture.

**Lactation**– Breast implant surgery may interfere with the ability to successfully breast feed, either by reducing or eliminating milk production. Particularly, periareolar incision may considerably reduce the possibility of breast-feeding.

**Avoiding Damage During Treatment**– Patients should inform other treating physicians of the presence of implants to minimize the risk of damage to the implants.

**Topical Medications**– The patient should consult a physician or a pharmacist before the use of topical medicines (e.g. steroids) in the breast area.

**Trauma**– The patient should consult your surgeon or a physician if you suspects any complications, in particular in the case of trauma or compression caused, for example, by extreme massaging of the breast region, by some sport activities or by using seat belts.

**Smoking**– Smoking may interfere with the healing process.

**Radiation to the Breast**– Establishment Labs S.A. has not tested the *in vivo* effects of radiation therapy in patients who have breast implants. The literature suggests that radiation therapy may increase the likelihood of capsular contracture, necrosis, and implant extrusion.

**Insurance Coverage**– Patients should check with their insurance company regarding coverage issues before undergoing surgery.

**Mental Health and Elective Surgery**– It is important that all patients seeking an elective procedure such as breast augmentation have realistic expectations that focus on improvement rather than perfection. Ask the patient to openly discuss, prior to surgery, any history of depression or other mental health disorders.

#### **Postoperative Care:**

The patient should be advised that she will likely feel tired and sore for several days following the operation, and that her breasts may remain swollen and sensitive to physical contact for a month or longer. She also could experience a feeling of tightness in the breast area as the skin adjusts to the new breast size. The patient should avoid any strenuous activities for at least a couple of weeks but should be able to return to work within a few days. Breast massage may also be recommended as appropriate.

#### **Life expectancy of the breast implant:**

Breast implants are not lifetime devices. However, the life expectancy of a silicone breast implant cannot be precisely estimated, as there are many factors beyond the manufacturer's control that can affect the longevity of a device. The time period varies from woman to woman. Some women could need replacement surgery few years after the augmentation procedure while others can have their implants intact for 10-20 years. Therefore, the life expectancy of the implant cannot be guaranteed.

#### **ADVERSE EVENTS**

Because breast implant surgery is more often performed using general anesthesia, it is associated with the same risks as other invasive surgical procedures.

After breast implant surgery, patients might experience swelling, hardness, discomfort, itching, bruising, twinges and pain over the first few weeks. Potential adverse events that may occur with silicone gel-filled breast implant surgery include:

### Capsular Contracture

Normally, capsules of collagen fibers form as an immune response around a foreign body, such as a breast implant, tending to isolate it. Capsular contracture occurs when the capsule tightens and squeezes the implant. This can cause the implant to turn rigid (from slightly firm to quite hard) and the firmest ones can cause varying degrees of discomfort, pain and palpability. In addition to the firmness, capsular contracture can result in a deformed breast, visible surface wrinkling and/or displacement of the implant. Detection of breast cancer by mammography may also be more difficult. Capsular contracture may be more common following infection, hematoma, and seroma, and the chance of it happening may increase over time. Capsular contracture occurs more commonly in patients undergoing revision surgery than in patients undergoing primary implantation surgery. Capsular contracture is a risk factor for implant rupture, and it is the most common reason for reoperation in augmentation and reconstruction patients.

Capsular contracture is graded into 4 levels depending on its severity. Baker Grade I: the breast is normally soft and looks natural; Baker Grade II: the breast is a little firm but looks normal; Baker Grade III: the breast is firm and looks abnormal; Baker Grade IV: the breast is hard, painful, and looks abnormal. Patients should also be advised that additional surgery might be needed in cases where pain and/or firmness are severe (Baker Grades III or IV) and that capsular contracture may happen again after additional surgeries.

Correction of capsular contracture may require surgical removal or release of the capsule, or removal and possible replacement of the implant itself. Closed capsulotomy (external manipulation of the capsule in order to "pop" the tissue capsule and open it up) used to be a common procedure for treating capsular contracture, but most manufacturers, including Establishment Labs S.A., contraindicate it because it can cause implant rupture.

### Rupture

Breast implants can potentially remain intact for decades in the body, but all such devices will fail at some point.

Breast implants rupture when the shell develops a tear or hole. Rupture can occur at any time after implantation, but it is more likely to occur the longer the implant is in place. The following may cause implants to rupture: damage by surgical instruments, implant stress and weakening during implantation, age and design of the implant, submuscular rather than subglandular location, occurrence of post-operative hematomas or seromas, folding or wrinkling of the implant shell, excessive force to the chest (e.g., during closed capsulotomy, which is contraindicated), trauma, compression during mammographic imaging, and severe capsular contracture.

Silicone gel-filled implant ruptures are most often silent. (MRI examination is currently the best screen method for silent rupture.) This means that most of the time neither doctor nor patient will know if the implant has a tear or hole in the shell. This is why a first MRI is recommended after 3 years, and then at regular intervals every 2 years thereafter, to screen for ruptures. Sometimes there are symptoms associated with gel implant rupture, such as lumps surrounding the implant or in the armpit, change or loss of size or shape of the breast or implant, pain, tingling, swelling, numbness, burning, or hardening of the breast.

When MRI findings of rupture are found, or if there are signs or symptoms of rupture, the implant should be removed, with or without replacing it. If the patient

develops symptoms that suggest implant rupture, she should be told to have an MRI evaluation to screen for rupture.

If rupture occurs, silicone gel may either remain within the scar tissue capsule surrounding the implant (intracapsular rupture), move outside the capsule (extracapsular rupture), or move beyond the breast (migrated gel). There is also a possibility that rupture may progress from intracapsular to extracapsular and beyond.

Below is a summary of information related to the health consequences of implant rupture, which **have not been fully established**, in women who had a variety of implant models from different manufacturers.

- Local breast complications that were associated with rupture in the literature include breast firmness, a change in breast shape or size, and breast pain. These symptoms are not specific to rupture and may be also experienced by women who have capsular contracture.
- There have been rare reports of gel migration to nearby tissues such as the chest wall, armpit, or abdominal wall, and to more distant locations down the arm or into the groin. This has led to nerve damage, granuloma formation and/or breakdown of tissues in direct contact with the gel in a few cases. There have been reports of silicone presence in the liver of patients with silicone breast implants. Movement of silicone gel material to lymph nodes in the armpit also has been reported, even in women without evidence of rupture, leading to lymphadenopathy.
- Concerns have been raised over whether ruptured implants are associated with the development of connective tissue or rheumatic diseases and/or symptoms such as fatigue and fibromyalgia. A number of epidemiology studies have evaluated large populations of women with breast implants from a variety of manufacturers and implant models. These studies do not support an association of breast implants and rheumatic disease.

### Pain

Most women undergoing augmentation or reconstruction with a mammary implant will experience some post-operative breast and/or chest pain. While this pain normally recedes in most women as they heal after surgery, it can become a chronic problem in other women.

Hematoma, migration, infection, implants that are too large, or capsular contracture can cause chronic pain. Sudden, severe pain may be associated with implant rupture. The surgeon should instruct the patient to immediately report if there is significant pain or if pain persists.

### Changes in Nipple and Breast Sensation

Breast surgery can result in an increased/decreased breast and/or nipple sensitivity. Typically the sensation is lost after complete mastectomy where the nipple itself is removed, and can be severely lessened by partial mastectomy. The range of changes varies from intense sensitivity to no feeling in the nipple or breast following surgery. While some of these changes can be temporary, they can also be permanent, and may affect the patient's sexual response or ability to nurse.

### Infection

Infection can occur with any surgery or implant. Most infections resulting from surgery appear within a few days to weeks after the operation. However, infection is possible at any time after surgery. In addition, breast and nipple piercing procedures may increase the possibility of infection. Infections in tissue with an implant present are harder to treat than infections in tissue without an implant. If an infection does not respond to antibiotics, the implant may have to be removed, and another implant may be placed after the infection is resolved.

As with other surgical procedures, toxic shock syndrome in rare instances has been noted in women after breast implant surgery. This is a life-threatening condition and its symptoms include sudden fever, vomiting, diarrhea, fainting, dizziness, and/or sunburn-like rash. Patients should be instructed to contact their doctor immediately for diagnosis and treatment if they have these symptoms.

### **Hematoma/Seroma**

Hematoma is a collection of blood within the space around the implant, and a seroma is a build-up of fluid around the implant. Having a hematoma and/or seroma following surgery may result in infection and/or capsular contracture later on. Symptoms from a hematoma or seroma may include swelling, pain, and bruising. If a hematoma or seroma occurs, it will usually be soon after surgery. However, they can also occur at any time after injury to the breast. While the body absorbs small hematomas and seromas, some will require surgery, typically involving draining, and potentially placing a surgical drain in the wound temporarily for proper healing. A small scar can result from surgical draining. Implant rupture also can occur from surgical draining if there is damage to the implant during the procedure.

### **Breast-feeding**

Although most women with breast implants who attempt nursing have successfully breast-fed their babies, it is not known if there are increased risks for a woman with breast implants or if the children of women with breast implants are more likely to have health problems. At this time, it is not known if it is possible for a small amount of silicone to pass from the breast implant silicone shell into breast milk during breast-feeding, or what the potential consequences might be. A periareolar surgical approach may further increase the chance of breast-feeding difficulties. However, the American Academy of Pediatrics has stated that there is no reason why a woman with implants should refrain from nursing.

### **Calcification**

Calcium deposits can form in scar tissue surrounding the implant and may cause pain and firmness, and be visible on a mammography. These deposits must be identified as different from calcium deposits that are a sign of breast cancer. Additional surgery may be necessary to remove and examine calcifications. Calcium deposits also occur in women who undergo breast reduction procedures, in patients who have had hematoma formation, and even in the breasts of women who have not undergone any breast surgery. The occurrence of calcium deposits significantly increases with age.

### **Delayed Wound Healing**

Some patients may experience a prolonged wound healing time. Smoking may interfere with the healing process. Delayed wound healing may increase the risk of infection, extrusion, and necrosis. Wound healing times may vary depending on the type of surgery or incision.

### **Implant Extrusion**

Lack of adequate tissue coverage, local trauma or infection may result in exposure and extrusion of the implant. This has been reported with the use of steroid drugs or after radiation therapy of breast tissue. If tissue breakdown occurs and the implant becomes exposed, implant removal may be necessary, which may result in additional scarring and/or loss of breast tissue.

### **Necrosis**

Necrosis is the formation of dead tissue around the implant. This may prevent wound healing and require surgical correction and/or implant removal. Permanent scar deformity may occur following necrosis. Factors associated with necrosis include infection, use of steroids in the surgical pocket, smoking,

chemotherapy/radiation, and excessive heat or cold therapy.

### **Granulomas**

These are benign lumps that can form when body cells surround foreign material such as silicone. Like any lump, it should be further evaluated to rule out a malignancy.

### **Breast Tissue Atrophy/Chest Wall Deformity**

The pressure of the breast implant may cause the breast tissue to thin and shrink (with increased implant visibility and palpability), potentially leading to chest wall deformity. This can occur while implants are still in place or following implant removal without replacement. Either of these conditions may result in additional surgeries and/or unacceptable dimpling/creasing of the breast.

### **Lymphadenopathy**

Literature reports associate lymphadenopathy with both intact and ruptured silicone breast implants. One study reported that armpit lymph nodes from women with both intact and ruptured silicone gel implants had abnormal tissue reactions, granulomas, and the presence of silicone. These reports occurred in cases of women who had implants from a variety of manufacturers and implant models.

### **Unsatisfactory Results**

Unsatisfactory results such as wrinkling, asymmetry, implant displacement/migration, incorrect size, implant palpability/visibility, scar deformity, and/or hypertrophic scarring, may occur. Some of these results may cause discomfort. Pre-existing asymmetry may not be entirely correctable by implant surgery. Revision surgery could be indicated to increase patient satisfaction, but this involves additional considerations and risks. Careful preoperative planning and surgical technique can minimize but not always prevent unsatisfactory results.

### **Other Reported Conditions**

There have been reports in the literature of other conditions in women with silicone breast implants. Many of these conditions have been studied to evaluate their potential association with breast implants. However, no causal relationship has been established between breast implants and the conditions listed below.

#### **Connective Tissue Disease (CTD)**

Since the early 1990s, nearly a dozen comprehensive systemic reviews have been commissioned by government health ministries in several countries to examine the alleged links between silicone gel breast implants and systemic diseases. A clear consensus has emerged from these independent scientific reviews that there is no clear evidence of a causal link between the implantation of silicone breast implants and connective tissue disease.

#### **Cancer**

Breast cancer reports in the medical literature reveal that patients with breast implants are not at a greater risk than those without breast implants for developing breast cancer. Some reports have suggested that breast implants may interfere with or delay breast cancer detection by mammography and/or biopsy; however, other reports in the medical literature indicate that breast implants do not significantly delay breast cancer detection or adversely affect cancer survival prognosis in implanted women. Some studies even suggest lower rates of breast cancer in women with breast implants.

#### **Neurological Disease, Signs, and Symptoms**

Some women with breast implants have experienced neurological disturbances (e.g., visual symptoms or alterations in sensation, muscle strength, walking, balance, thinking or memory) or diseases (e.g., multiple sclerosis) and they believe those symptoms are related to their implants. However, there is no

evidence in published literature of a causal relationship between breast implants and neurological disease.

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Small quantities of silicone may diffuse through the elastomer envelope of silicone gel-filled implants. The detection of small quantities of silicone in the periprosthetic capsule, axillary lymph nodes and other distal regions in patients with apparently intact gel-filled implants has been reported in the literature. Some studies on long-term implants have suggested that gel-bleed may contribute to the development of capsular contracture and lymphadenopathy. On the other hand, evidence against gel-bleed being a significant contributing factor to capsular contracture and other local complications is provided by the fact that there are similar or lower complication rates for silicone gel-filled breast implants than for saline-filled breast implants.

### Gel Fracture

Gel fracture may occur with cohesive silicone as a result of handling during surgery or alternatively due to the development of capsular contracture and may result in device distortion. This can lead to both patient and surgeon dissatisfaction with the aesthetic outcome of surgery and may require a second procedure.

### Interference with Mammography

Breast implants (especially in subglandular placement) may complicate the interpretation of mammographic images by obscuring the underlying breast tissue and/or by compressing overlying tissue. Despite of the fact that the presence of breast implants lessens tissue compression range during mammography, a number of studies looking at breast cancers in women with implants have found no significant difference in stage of disease at time of diagnosis, and prognosis appears to be similar in implanted and non implanted patients. Accredited mammography centers, technicians with experience in imaging patients with breast implants, and use of displacement techniques are needed to adequately visualize breast tissue in the implanted breast. Anterior breast tissue is best visualized with displacement views and posterior breast tissue with compression views. The decrease in visible area of 35% with compression views is improved to 25% with displacement views. The current recommendations for preoperative/screening mammograms are no different for women with breast implants than for those women without implants.

### Anaplastic Large Cell Lymphoma (ALCL) Lymphomas, including anaplastic large T-cell lymphoma (ALCL)

According to European safety information and from the U.S. Food and Drug Administration (FDA), as well as scientific literature, a possible association was identified between the breast implants and the rare development of anaplastic large cell lymphoma (ALCL), a type of non- Hodgkin's lymphoma. Women with breast implants may have a very small, but increased risk of developing anaplastic large cell lymphoma in an area adjacent to the implant. This specific entity is included in the WHO 2016 classification under the terminology "BIA-ALCL".

The U.S. Food and Drug Administration in 2011 released a safety communication stating that women with breast implants may have a very small, but increased risk of developing anaplastic large cell lymphoma, or ALCL, in the scar tissue or fluid adjacent to a breast implant. ALCL is not breast cancer; it is a rare type of cancer involving the cells of the immune system.

Cases of ALCL in the breast have also been reported in women without breast implants. When occurring in association with breast implants, patients were diagnosed following symptoms such as pain, lumps, swelling, fluid accumulation, or asymmetry. In the cases reported, ALCL was most commonly diagnosed years after implant placement and usually located immediately surrounding the breast implant. ALCL associated with breast implants is distinct and separate from primary lymphoma of the breast, which has a worse prognosis and is a majority B-cell phenotype (Domchek, 2002)<sup>1</sup>. In contrast, ALCL is always a T-cell lymphoma, and in epidemiological studies occurs at a higher rate in women with breast implants compared to the general population (DeJong, 2008)<sup>2</sup>.

Confirmed cases of ALCL may require surgery and chemotherapy as part of an individualized treatment plan by appropriate specialists. Women without symptoms do not require screening or implant removal and should continue with their normal breast health routine. The U.S. Food and Drug Administration and other regulatory bodies have confirmed that breast implants have a reasonable assurance of safety and efficacy.

## INSTRUCTIONS FOR USE

### Single Use

This product is intended to be used only in one patient for a single procedure. DO NOT reuse explanted implants. To reuse a single-use device could expose patients and staff to risks which outweigh the perceived benefits of using such devices. This product is not intended to be reprocessed in any way and/or used again, not even on the same patient. The reuse of single-use devices can affect their safety, performance and effectiveness, exposing patients and staff to unnecessary risk, such as infection, the inability to assure the proper cleaning and decontamination, the presence of residues of cleaning agents, reaction to endotoxins, exposure to other bio-hazards and/ or device failure. This practice may also have legal implications that vary according to each jurisdiction<sup>3</sup>.

### Product Traceability

The product traceability stickers, provided with each device and located within the internal product packaging, provide product-specific information and should be attached to the patient's chart for identification purposes. Stickers are also available for the Device Tracking Form, the Patient ID Card and the hospital files, if applicable. You should encourage your patient to participate in the Establishment Labs S.A. device-tracking program, entering their implant(s) information at [www.motivaimplants.com](http://www.motivaimplants.com).

This will help ensure that Establishment Labs S.A. has a record of each patient's contact information so that they can be contacted in the event of a recall or other problems with the implants of which they should be made aware.

### Sterile Product

Each sterile silicone breast implant is supplied in a sealed, double sterile barrier primary package. Use standard procedures to maintain sterility during transfer of the breast implant to the sterile field. Remove the breast implant and accessories from their packages in an aseptic environment, using talc-free gloved hands.

1. Domchek SM, Hecht JL, Fleming MD, Pinkus GS, Canellos GP. Lymphomas of the breast: primary and secondary involvement. *Cancer*. 2002 Jan 1;94(1):6-13. <http://onlinelibrary.wiley.com/doi/10.1002/cncr.10163/pdf>

2. de Jong D, Vasmel WL, de Boer JP, Verhave G, Barbé E, Casparie MK, van Leeuwen FE. Anaplastic large-cell lymphoma in women with breast implants. *JAMA*. 2008 Nov 5;300(17):2030-5. <https://www.ncbi.nlm.nih.gov/pubmed/18984890>

Sterility of the implant is maintained only if the thermoform packages, including the package seals, are intact.

**DO NOT** use the product if the thermoform packages or seals have been damaged.

**DO NOT** re-sterilize the product.

Avoid prolonged exposure to extreme storage conditions. We recommend keeping these devices at room temperature, at atmospheric pressure, in dry conditions and away from direct sunlight.

**DO NOT** implant any device that may appear to have particulate contamination, damage, or loss of shell integrity. A sterile back-up implant must be readily available at the time of surgery.

**DO NOT** implant any device that may appear to have leaks or scratches.

### How to Open Sterile Product Package


**DO NOT** expose the breast implant to talc, sponges, towels, or other contaminants.

1. A non-sterile team member should open the outer package.
2. Remove the inner package and invert it over the sterile field, allowing the sealed inner thermoform package to gently slide into the field.
3. Use the pull-tab to open the lid of the inner thermoform package.
4. Retrieve the breast implant and examine it for any particulate contamination, damage, or loss of shell integrity. If satisfactory, return the breast implant to the inner thermoform tray. At this point, you might slightly rinse the implant with a small amount of saline to remove the static and cover the tray while the lid until implantation to prevent contact with airborne and surgical field particulate contaminants.

### Surgical Technique and Implant Selection

There are several surgical techniques that can be used to perform the implantation of silicone gel-filled breast implants. Therefore, the surgeon is advised to use his/her clinical judgment in choosing the procedure that is best for the patient, consistent with this product insert. After setting realistic aesthetic goals that assure mutual understanding between doctor and patient, the surgeon must choose from current and accepted surgical techniques to minimize the incidence of adverse reactions and achieve the best results.

The surgeon must carefully assess the implant size, projection and surface, as well as the incision placement, pocket dissection and implant placement criteria according to the patients' anatomy and desired aesthetic outcomes.


### Incision Site Selection

- The periareolar incision is usually more concealed, but may considerably reduce the possibility of future breast-feeding as compared to other incision

sites. A periareolar incision may be associated with a higher risk of changes in nipple sensation.

- The inframammary incision is generally less concealed than the periareolar, but it is associated with less breast-feeding difficulties.
- The axillary incision is the least concealed of all incision sites.
- The periumbilical approach should not be used with Sterile Silicone Breast Implants Motiva Implant Matrix® for a number of reasons, including potential damage to the implant shell.

### Implant Placement Selection

- The possible benefits of submuscular placement are that it may result in less palpable implants, less likelihood of capsular contracture, and easier mammographies. This placement may be preferable if the patient has thin or weakened breast tissue. However, submuscular placement is associated with a longer surgical procedure, a more prolonged recovery period and more pain. Also, it may make it more difficult to perform some reoperation procedures.
- Subglandular placement may make surgery and recovery shorter, be less painful, and provide easier access for reoperation than the submuscular placement. However, this placement may result in more palpable implants, greater risk of capsular contracture, ptosis and increased difficulty in imaging the breast with mammography.
- Dual plane placement has been associated by some authors with the benefits of submuscular placement with the advantages of a faster recovery and less pain and postoperative discomfort.

### During the Surgical Procedure:

- It is advisable to have more than one size of breast implant in the operating room at the time of surgery to allow for flexibility in determining the appropriate size to be used.
- A backup implant should also be available.
- Insufficient pocket dissection increases the risk of rupture and incorrect positioning of the implant. Make sure that the incision is large enough to facilitate insertion without excessive manipulation and handling of the device, and to avoid damage. A well-defined, dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.
- DO NOT** use lubricants during placement since they increase the risk of pocket contamination and may also affect the tissue-capsule interface.
- DO NOT** damage the breast implant with sharp surgical instruments such as needles and scalpels, blunt instruments such as clamps and forceps, or by

- over-handling and manipulation during introduction into the surgical pocket.
- **DO NOT** use excessive force during breast implant placement.
- **DO NOT** manipulate the implant for radial expansion, compression or dissection of the pocket.
- **DO NOT** use more than one implant per breast pocket.

### **Maintaining Hemostasis/Avoiding Fluid Accumulation**

The risk of postoperative hematoma and seroma may be reduced by carefully handling the hemostasis during surgery, and possibly also by postoperative use of a closed drainage system. Persistent or excessive bleeding must be controlled before implantation.

Any postoperative evacuation of hematoma or seroma must be conducted with care to avoid contamination of or damage to the breast implant.

### **Instructions and precautions for removal**

Among the most common reasons for breast implant removal are complications such as capsular contracture, suspected rupture and implant malposition, as well as the patient's desire to change the implant size or shape. The surgeon is advised to use his/her clinical judgment in choosing from current and accepted breast implants removal and replacement surgical techniques to minimize the incidence of adverse reactions and achieve the best results for the patient.

## **SPECIFIC INSTRUCTIONS FOR USE**

### **APPLICABLE TO BREAST IMPLANTS CONTAINING A MICROTRANSPONDER**

#### **Description**

Breast implants containing a microtransponder include a miniaturized, implantable, radio frequency identification device (RFID), which is placed in the breast implant filler material. Scanners to scan and read the information in the microtransponders are purchased separately. The microtransponder is a passive device that contains an electronic circuit that is activated externally by a low-power electromagnetic field, emitted by a handheld battery powered scanner. The microtransponder is used to store an electronic unique device identification number ("Qid™"). The ESN number is used by patient-approved physicians and other health professionals, to access a database that will provide the implanted device specific information and other information provided by the patient.

#### **Indications**

The microtransponder is indicated for use as a miniature, implantable microchip that is inserted into the breast implant. The microtransponder in the breast implant provides the patient a unique identification number that may be used to access a database containing the breast implant information (serial and lot numbers; reference number; volume, size and projection, model, surface type, manufacturing date, etc.) and any other information that may be supplied by the patient or the physician in relation to the surgical procedure or the patient condition.

#### **Contraindications**

Breast implants containing a microtransponder should not be used in patients known to have allergies or sensitivity to the composition of the microtransponder (USP Type III glass).

#### **Precautions**

MRI Compatible. Patients with the breast implants containing a microtransponder may safely undergo MRI diagnostics in up to 3 Tesla cylindrical systems. See

section below: INSTRUCTIONS FOR PATIENTS UNDERGOING MRI for detailed instructions.

### **Instructions For Patients Undergoing MRI**

The patient should be monitored continuously throughout the MRI procedure using visual and audio means (e.g. intercom system). Instruct the patient to alert the MRI system operator of any unusual sensations or problems so that, if necessary, the MRI system operator can immediately terminate the procedure. Provide the patient with a means to alert the MRI system of any unusual sensations or problems. Do not perform the MRI if the patient is sedated, anesthetized, confused, or otherwise unable to communicate with the MRI system operator. Establishment Labs S.A. Sterile Silicone Breast Implants Motiva Implant Matrix® are MRI conditional.

Patients should be advised to have regular MRIs over their lifetime to screen for silent rupture even if they are not having any apparent problems. As mentioned before, it is advised to have the first MRI 3 years postoperatively, then regularly at 2-year intervals. The patient implanted with Sterile Silicone Breast Implants can undergo MRI scan under the following conditions:

- Static magnetic field of 1.5-Tesla and 3 -Tesla only.
- Maximum spatial gradient magnetic field of 4.000-gauss/cm (40-T/m).
- Maximum MR system reported, whole body average specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Control Operating Mode.
- Under the scan defined conditions, the Motiva Implant Matrix® Breast Implant with Qid™ is expected to produce a maximum temperature rise of 3.0°C after 15 minutes of continuous scanning (i.e., per pulse sequence).

In non-clinical testing, the magnetically induced displacement force and magnetically induced torque were tested and no clinically significant displacement or torque was detected. It is important to note that Motiva Implants® with Q Inside Safety Technology™ (Qid™) contain a ferrite component in the microtransponder that provides unique device identification data through an external reader. This microtransponder may create an imaging void during breast implant MRI (known as artifact effect) that can block visualization of a small area near the Qid™. In non-clinical testing, the image artifact caused by Sterile Silicone Breast Implants Motiva Implant Matrix® breast implants with Qid™ extends approximately 15 mm from this implant when imaged using a gradient echo pulse sequence and a 3-Tesla MR system.

In selected cases, the following imaging techniques can be used as an alternative to better visualize the region affected by the implant Qid™:

- Ultrasound.
- Tomosynthesis.
- Digital compression mammogram.
- Subtraction contrast Mammography.
- Scintimammography.

### **Additional Instructions For Use**

Additional Directions for Use for breast implants which include a microtransponder:

- Verify the microtransponder in the implant before opening the sterile barriers, using the corresponding scanner, if available.
- Re-verify the microtransponder in the implant after implantation using the corresponding scanner, if available.

### **Cautions**

If the breast area is subsequently subjected to physical trauma as a result of

accident or injury, the patient shall consult her physician to ensure that the microtransponder is properly functioning. If for any reason the microtransponder would stop being scannable by the appropriate scanner, this sole situation will not impair the breast implant to continue fulfilling its function appropriately and does not constitute a complication.

## SPECIFIC INSTRUCTIONS FOR USE

### APPLICABLE TO BREAST IMPLANTS CONTAINING RADIOPAQUE ORIENTATION MARKS

#### Description

Breast implants containing radiopaque orientation mark include one or several marks, such as dots, lines or tees made out of radiopaque material to identified potential rotation of the device after implantation during MRI procedure. Radiopaque orientation marks are available in some Motiva Implant Matrix® product lines and are specifically a tool for the physician.

The physician is responsible for deciding whether the result of the breast implant rotation detected during MRI procedure after implantation, if any, is greatly enough to justify reaccommodate the breast implant, an explantation of the device or a replacement.

#### Indications

The radiopaque orientation marks are indicated for detection of breast implant rotation during MRI procedure. The physician is responsible for deciding whether the result of the breast implant rotation detected during MRI procedure after implantation, if any, is greatly enough to justify reaccommodate the breast implant, an explantation of the device or a replacement.

#### Contraindications

No contraindications are available at the moment.

#### Precautions

MRI Compatible.

Radiopaque marks may not be seen during other equivalent procedures. Radiopaque marks may be confused with tissue nodules. Always verify the specific type of mark in the device before analyzing the MRI results to avoid confusion and misleading decisions.

#### Caution

Radiopaque marks may be confused with tissue nodules. Always verify the specific type of mark in the device before analyzing the MRI results to avoid confusion and misleading decisions.

### DEVICE TRACKING

Silicone gel-filled breast implants are subject to Device Tracking. Compliance with this requirement is mandatory. This means that it is required to report to Establishment Labs S.A., either directly or through a representative, the lot and serial numbers of the device(s) implanted in a patient, the date of her surgery, her ID number, her personal contact information, and information relating to the surgeon's practice. This information will be recorded on a Device Tracking Form supplied by Establishment Labs S.A. with each silicone gel-filled breast implant. Establishment Labs S.A. strongly recommends that all patients receiving silicone gel-filled breast implants participate in the Establishment Labs S.A. device-tracking program, entering their implant information at [www.motivainplants.com/#implantRegistration](http://www.motivainplants.com/#implantRegistration). This will help ensure that Establishment Labs S.A.

has a record of each patient's contact information so that they can be contacted in the event of a recall or other problems with the implants of which they should be made aware.

### PROCEDURE RECORDING AND DEVICE ID CARD

Each breast implant is supplied with five Patient Record Labels showing the reference number, lot number, serial number, side (left or right), and volume of the implant. Patient Record Labels are located on the internal product packaging attached to the main label. To complete the Patient ID Card, adhere one Patient Record Label for each implant on the back of each patient's ID Card. Another label should be affixed to the patient's chart. A third label should be attached to the practitioner's records and the fourth label is provided for hospital records when applicable. A spare label is also provided. If a Patient Record Label is unavailable, the appropriate information may be copied by hand from the device label.

#### PATIENT ID CARD

Every patient must have a record of her surgical procedure in case of future consultations or additional surgeries. Each implant is provided with a Patient ID card, which must be given to the patient for personal reference. Apart from the information stated on the Record Labels that should be affixed to the back of the card, the Patient ID card includes the patient's name, position of the implant (submuscular, subglandular, dual plane, other), date of implantation and name of the treating surgeon.

### INFORMATION ON EXPECTED LIFETIME

In practice, it is not possible to predict accurately the actual lifetime of an individual implant. It is well understood that several factors are beyond the control of the manufacture. These factors might have a significant effect on the lifetime of an individual device. The factors include the actual implantation procedure, the anatomy and state of health of the patient, the behavior and activities (for example sporting activities), as well as predictable and unpredictable external mechanical influences.

### REPORTING AND DEVICE RETRIEVAL EFFORTS

In the event of an explantation, the reason for explantation should be reported on the Establishment Labs S.A. Incident Report Form and the explanted device returned to the local Establishment Labs S.A. representative. In case there is no local representative available, report directly at Establishment Labs S.A. Zona Franca Coyoil. Distrito San José de Alajuela, Costa Rica; Phone: +506 2434-2400 / Fax: +506 2434-2450 or [customerservice@establishmentlabs.com](mailto:customerservice@establishmentlabs.com).

The explanted device must be decontaminated and properly packaged before returning it according to the Establishment Labs S.A. Explanted Implants Return Protocol. In case local health codes do not allow for the implant to be returned, please contact us directly at [customerservice@establishmentlabs.com](mailto:customerservice@establishmentlabs.com) for specific instructions.

### PRODUCT EVALUATION

Establishment Labs S.A. requires that any complications resulting from the use of this device be brought to the immediate attention of the company through the Motiva Implant Matrix® Incident Report Form addressed to Establishment Labs S.A. Zona Franca Coyoil. Distrito San José de Alajuela, Costa Rica; Phone: +506 2434-2400 / Fax: +506 2434-2450 or [customerservice@establishmentlabs.com](mailto:customerservice@establishmentlabs.com).

## RETURNED GOODS POLICY

Product returns should be handled through the local Establishment Labs S.A. representative. In case there is no local representative available, report directly at Establishment Labs S.A. Zona Franca Coyoil. Distrito San José de Alajuela, Costa Rica; Phone: +506 2434-2400 / Fax: +506 2434-2450 or customerservice@establishmentlabs.com.

All package seals must be intact for goods to be eligible for return. Returned products may be subject to a restocking charge. For more information, please contact the local Establishment Labs S.A. representative.

## ESTABLISHMENT LABS S.A. ALWAYS CONFIDENT WARRANTY® LIMITED WARRANTY, LIMITATION OF LIABILITY, AND DISCLAIMER OF OTHER WARRANTIES

The complete terms, conditions and limitations of the Establishment Labs S.A. Always Confident Warranty® can be reviewed in the website [www.motivaimplants.com](http://www.motivaimplants.com) or can be provided by the local Establishment Labs S.A. representative. No warranty or program of Establishment Labs S.A. covers any costs, fees or expenses related to any medical treatment and/or the surgical replacement of the implants. Establishment Labs S.A. warrants that this product is free of manufacture defects at the time of its shipment. Establishment Labs S.A. shall not be responsible for any incidental or consequential loss, damage or expenses directly or indirectly arising from the use of this product. Establishment Labs S.A.'s sole responsibility in the event that Establishment Labs S.A. determines the product was defective when shipped by Establishment Labs S.A., shall be the replacement of the product. Establishment Labs S.A. assumes no further liability. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law, or otherwise, including, but not limited to, any implied warranties of merchantability, suitability for use or performance.

## MANUFACTURER

**Establishment Labs S.A.:** Coyoil Free Zone & Business Park Building 4th Street, Building B-15, Alajuela, Costa Rica. Phone: +506 2434-2400 Fax: +506 2434-2450

[customerservice@establishmentlabs.com](mailto:customerservice@establishmentlabs.com)

[www.establishmentlabs.com](http://www.establishmentlabs.com)

[www.motivaimplants.com](http://www.motivaimplants.com)

## EUROPEAN REPRESENTATIVE

**Emergo Europe:** Prinsessegracht 20, 2514 AP The Hague, The Netherlands.

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## GARANTIA ALWAYS CONFIDENT DE ESTABLISHMENT LABS S.A

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## FABRICANTE

**Establishment Labs S.A.:** Coyoil Free Zone & Business Park Building 4th Street, Building B-15, Alajuela, Costa Rica. Telefone: +506 2434-2400 Fax: +506 2434-2450  
[customerservice@establishmentlabs.com](mailto:customerservice@establishmentlabs.com)  
[www.establishmentlabs.com](http://www.establishmentlabs.com)  
[www.motivaimplants.com](http://www.motivaimplants.com)

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Establishment Labs S.A.: Coyol Free Zone & Business Park Building  
4th Street, Building B-15, Alajuela, Costa Rica, 전화: +506 2434-  
2400 팩스: +506 2434-2450  
customerservice@establishmentlabs.com  
www.establishmentlabs.com  
www.motivaimplants.com

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Emergo Europe: Prinsessegracht 20, 2514 AP The Hague, The  
Netherlands.

EDC Motiva BVBA

Oude Bosuilbaan 107 Deurne 2100, Belgium.  
전화: +32 3460-1133 팩스: +32 3460-1132

## MOTIVA IMPLANTS FEEDBACK FORM

This survey provides valuable customer input that helps us offer the best possible service. **Thank you** for completing and returning the Motiva Implants Feedback Form by mail to Establishment Labs S.A. Coyoil Free Zone & Business Park Building 4th Street, Building B-15, Alajuela, Costa Rica, or by fax to +506 2434-2450. You can also fill in the form on-line at [www.motivainplants.com](http://www.motivainplants.com)

Name: .....	<b>Implant Information</b>
Gender: ..... Age: .....	Reference Number: .....
Address: .....	Serial Number: .....
City: .....	Product: .....
Country: .....	Lot Number: .....

Please rate your satisfaction level with each of the following statements:

- |                        |                               |
|------------------------|-------------------------------|
| 1 = very satisfied     | 4 = somewhat dissatisfied     |
| 2 = somewhat satisfied | 5 = very dissatisfied         |
| 3 = neutral            | Leave blank if not applicable |

<b>PACKAGING</b>	1	2	3	4	5
1. How satisfied are you with the presentation of the external product packaging?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Do you find the packaging easy to open?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfied are you with the presentation of the internal packaging (thermoforms)?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Do you find the device easy to retrieve from its packaging?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. Please suggest any packaging improvements: .....					

<b>PATIENT ID CARD</b>	1	2	3	4	5
1. How satisfied are you with the presentation of the Patient ID Card?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Do you think the Patient ID Card is easy to find inside the box?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfied are you with the information collected on the Patient ID Card?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Please suggest any Patient ID Card improvements: .....					

<b>PRODUCT INSERT</b>	1	2	3	4	5
1. How satisfied are you with the presentation of the Product Insert?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. How satisfied are you with the information provided in the Product Insert?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. Please suggest any Product Insert improvements: .....					

<b>LABELING</b>	1	2	3	4	5
1. Rate your satisfaction with the information provided on the external packaging labels.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. Rate your satisfaction with the information provided on the internal packaging labels.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. Please suggest any Product Labeling improvements: .....					
<b>PRODUCT</b>	1	2	3	4	5
1. Rate your satisfaction with the product range available inside Motiva Implant Matrix.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
2. How satisfying is the appearance of the implant?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
3. How satisfying is the feel of the implant?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
4. Rate your satisfaction with the implant ease of insertion into the surgical pocket.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
5. How satisfying is the quality of the product?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
6. Rate your overall satisfaction with the product.	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
7. Please rate Motiva Implants to similar products on the market. <input type="radio"/> Much better <input type="radio"/> Somewhat better <input type="radio"/> Average <input type="radio"/> Somewhat worse <input type="radio"/> Much worse If your answer is "somewhat worse", or "much worse", please explain: .....					
8. Would you use Motiva again in the future? <input type="radio"/> Definitely <input type="radio"/> Probably <input type="radio"/> Probably not <input type="radio"/> Definitely not <input type="radio"/> Not sure If your answer is "probably not", "definitely not" or "not sure", please explain: .....					
9. Would you recommend Motiva breast implants to a colleague? <input type="radio"/> Definitely <input type="radio"/> Probably <input type="radio"/> Probably not <input type="radio"/> Definitely not <input type="radio"/> Not sure If your answer is "probably not", "definitely not" or "not sure", please explain: .....					
<b>SUPPORT</b>	1	2	3	4	5
How satisfied are you with delivery of the products?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
How satisfied are you with the support provided by your distributor?	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
Please suggest how Establishment Labs can improve your customer experience: .....					

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Copyright 2016 Establishment Labs S.A.