



CCM-iBIS-SG

Network Connection Instruction

# Local Network: Connecting iBIS, NAS, and PC as local network system



# Port Allocation : Local Network Connection



1. Connect iBIS to Port No. 1 of the NAS

2. Connect Phototune PC to Port No.2 of the NAS, and set Image Data Folder in Phototune Config Unity to the IP address of Port No. 2

<CCM-iBIS-SL/SG>  
Device Name : CCM-NEXT-01  
IP : 192.168.40.11

<NAS>  
Device Name : CCM-NAS-01  
LAN Port No. 1\_IP : 192.168.40.51  
LAN Port No. 2\_IP : 192.168.40.61

<Phototune PC>  
Device Name : Not specified  
IP : 192.168.40.xx (can be any number)

# Clinical Network : Connecting the system to other networks

## NAS & Switching Hub

Use designated LAN cable(3m)  
KB-T7PK-03NV(CAT7\_ScTP)

KB-T7PK-03NV(CAT7\_ScTP)

< HUB >

\*Optional Item  
EHC-F16MN-HW

< External Device >  
Electrical Health Record etc.

Must be grounded

## iBIS & NAS

Use designated LAN cable(3m)  
KB-T7PK-03NV(CAT7\_ScTP)

KB-T7PK-03NV(CAT7\_ScTP)

## Switching Hub & PC

Use designated LAN cable(3m)  
KB-T7PK-03NV(CAT7\_ScTP)

Data Transfer

Data Transfer

< CCM-iBIS-SL/SG >

Device Name : CCM-NEXT-01

IP : 192.168.40.11

Must be grounded

< NAS >

Device Name : CCM-NAS-01

LAN Port No. 1\_IP : 192.168.40.51

LAN Port No. 2\_IP : 192.168.40.61

Must be grounded

< Phototune PC >

Device Name : Not specified

IP : 192.168.40.xx

(xx can be any number)

# Port Allocation : Clinical Network Connection

2. Connect Switching Hub to Port No. 2 of the NAS

1. Connect iBIS to Port No. 1 of the NAS

< HUB >

< External Device >  
Electrical Health Record etc.

4. Connect Phototune PC to the Switching Hub, and set Image Data Folder in Phototune Config Unity to the IP address of Port No. 2

3. Change IP address of NAS Port No.2 according to the clinical network IP if necessary

< CCM-iBIS-SL/SG >

Device Name : CCM-NEXT-01

IP : 192.168.40.11

< NAS >

Device Name : CCM-NAS-01

LAN Port No. 1\_IP : 192.168.40.51

LAN Port No. 2\_IP : 192.168.40.61

< Phototune PC >

Device Name : Not specified

IP needs to be the same work group as that of clinical network