Cardiac

AFIAS Tn-I Plus

INTENDED USE

AFIAS Tn-I Plus is a fluorescence Immunoassay (FIA) for the quantitative determination of cardiac troponin-I (Tn-I) in human whole blood/serum/plasma. It is useful as an aid in management and monitoring of acute myocardial infarction (AMI).

For *in vitro* diagnostic use only.

INTRODUCTION

Cardiac troponins are currently the most sensitive and specific biochemical markers of myocardial necrosis. There are three types of troponin in heart muscle fibers. Those are troponin-C, -I, and -T. Together they contribute to make cardiac muscle fibers contract. The clinical measurement of serum Tn-I has become an important tool in the diagnosis of acute myocardial infarction. Serum Tn-I is a more reliable than creatine kinase as a prognostic marker in people with ischemic chest pain. National and international scientific organizations have suggested the use of troponins, Tn-I and Tn-T, when implementing new diagnostic strategies in patients with acute coronary syndrome.

PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibodies in buffer bind to antigens in sample, forming antigen-antibody complexes and migrate onto nitrocellulose matrix to be captured by the other immobilized antibodies on the strip.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by instrument for AFIAS tests to show Tn-I concentration in the sample.

COMPONENTS

AFIAS Tn-I Plus consists of a 'cartridges'.

- Each cartridge packaged in an aluminum pouch has three components including cartridge part, detector part and diluent part.
- Cartridge part contains the membrane called a test strip, which has streptavidin at the test line, and chicken IgY at the control line.
- Detector part contains anti Tn-I-fluorescence conjugate, anti Tn-I-biotin conjugate, anti chicken IgY-fluorescence conjugate, bovine serum albumin (BSA), sucrose, bromophenol blue, and sodium azide as a preservative.
- Diluent part contains tween 20, CA-630, and sodium azide as a preservative.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow the instructions and procedures described in this 'Instruction for use'.
- Use fresh samples only and avoid direct sunlight.
- It is possible to use frozen samples. Please refer to "SAMPLE COLLECTION AND PROCESSING."

- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange test components between different lots or use test components after the expiration date, either of which might yield incorrect test result(s).
- The cartridge should remain sealed in its aluminum pouch until just before use. Do not use the cartridge if the pouch found damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local requirement. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- Allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used cartridges and pipette tips should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
- No Biotin interference was observed in **AFIAS Tn-I Plus** when biotin concentration in the sample was below 5 ng/mL. If a patient has been taking biotin at dosage of more than 0.03 mg a day, it is recommended to test again 24 hours after discontinuation of biotin intake.
- AFIAS Tn-I Plus** will provide accurate and reliable results subject to the below conditions.
 - AFIAS Tn-I Plus** should be used only in conjunction with the instrument for AFIAS tests.
 - Have to use recommended anticoagulant sample.

Recommended anticoagulant
Sodium heparin, Lithium heparin, Sodium citrate

STORAGE AND STABILITY

Storage condition			
Component	Storage Temperature	Shelf life	Note
Cartridge	2 - 8 °C	20 months	Unopened
		1 month	Resealed

- Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal.

LIMITATION OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause the false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors,

degradation of the test components/reagents or presence of interfering substances in the test samples.

- Any clinical diagnosis based on the test result should be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

MATERIALS SUPPLIED

REF SMFP-35

Components of **AFIAS Tn-I Plus**

Cartridge Box Contains	
- Cartridge	24
- Pipette tip (Zipper bag)	24
- ID chip	1
- Instruction for use	1
- Spare cartridge zipperbag	1

MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from **AFIAS Tn-I Plus**.

Please contact our sales division for more information.

AFIAS-6	REF FPRR020
AFIAS-1	REF FPRR019
Boditech Tn-I Plus Control	REF CFPO-212
Boditech Tn-I Plus Calibrator	REF CFPO-213

SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS Tn-I Plus** is human whole blood/serum/plasma.

- It is recommended to test the sample within 24 hours after collection.
- The serum or plasma should be separated from the clot by centrifugation within 3 hours after the collection of whole blood. If longer storage is required, e.g. if the test could not be performed within 24 hours, serum or plasma should be immediately frozen below -20 °C. The freezing storage of sample up to 3 months does not affect the quality of results.
- However, the whole blood sample should not be kept in a freezer in any case.
- Once the sample was frozen, it should be thawed one time and only for test, because repeated freezing and thawing can result in the changed test values.

TEST SETUP

- Check the components of the **AFIAS Tn-I Plus** as described below: Cartridges, Pipette tips, ID chip, Cartridge spare zipper bag and Instruction for use
- If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for AFIAS tests
- Empty the tip box.
- Insert the ID chip into the "ID chip port".
(Please refer to the instrument for AFIAS tests 'Operation Manual' for complete information and operating instructions.

TEST PROCEDURE

- Take 100 µL of sample with a pipette and dispense it into the sample well in the cartridge.
- Insert the cartridge into the cartridge holder
- Insert a tip into the tip hole of the cartridge.
- Tap the 'START' icon on the screen.
- The test result will be displayed on the screen after 12 minutes.

※ Note: Refer to the instrument for AFIAS tests Operation Manual to select sample type.

INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculates the test result automatically and displays Tn-I concentration of the test sample in terms of ng/mL.
- The working range of the **AFIAS Tn-I Plus** is 0.01-15.00 ng/mL.
- Expected Values**
 - In studies performed with the **AFIAS Tn-I Plus** assay involving 70 healthy volunteers in Korea, the upper reference limit (99th percentile) for Tn-I was 0.04 ng/mL. The lowest concentration with a CV less than or equal to 10 % with the **AFIAS Tn-I Plus** assay was 0.04 ng/mL.
 - Due to the release kinetics of Tn-I, a result below the decision limit within the first hours of the onset of symptoms does not rule out myocardial infarction with certainty. If myocardial infarction is still suspected, repeat the test at appropriate intervals.
 - A cut-off of 0.3 ng/mL Tn-I is recommended for diagnosis of AMI, as this yields optimal performance of 91 % of sensitivity and 92.1 % of specificity. However, laboratories should establish their own diagnostic cut-off concentration based on the clinical practice at their perspective institutions.

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are not provided with AFIAS Tn-I Plus. For more information regarding obtaining the control materials, contact Boditech Med Inc.'s Sales Division for assistance. (Please refer to the instruction for use of control material.)

PERFORMANCE CHARACTERISTICS

Analytical sensitivity	
Limit of Blank (LoB):	0.004 ng/mL
Limit of Detection (LoD):	0.01 ng/mL
Limit of Quantitation (LoQ):	0.03 ng/mL

■ Analytical specificity

- Cross-reactivity

Biomolecules such as below the ones in the table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS Tn-I Plus** test results did not show any significant cross-reactivity with these biomolecules.

Cross-reactivity material	Conc. (ng/mL)
CK-MB	60
NT-proBNP	1,000
Myoglobin	1,000
D-Dimer	1,000

- Interference

Interference materials such as below the ones in the table were added to the test sample(s) the same as the below concentrations. **AFIAS Tn-I Plus** test results did not show any significant interference with these materials except for EDTA.

Interference material	Conc.
Bilirubin	350 µmol/L
Cholesterol	13 mmol/L
D-glucose	1,000 mg/dL
Hemoglobin	2 g/L
L-Ascorbic acid	350 µmol/L
Triglyceride mixture	500 mg/dL
Heparin	3000 U/L
Sodium citrate	2 mg/mL
EDTA	3.4 µmol/L
Biotin	5 ng/mL

■ Precision

- Intra-assay (Between lot)

One person tested three different lots of **AFIAS Tn-I Plus**, ten times at each concentration of the control standard.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.94	0.95	0.92
7.50	7.65	7.32	7.51
Conc. (ng/mL)	Mean	SD	CV (%)
0.23	0.23	0.013	5.5
0.94	0.93	0.057	6.1
7.50	7.50	0.444	5.9

- Inter-assay (Between person)

Three different persons tested **AFIAS Tn-I Plus**; five times at each concentration of the control standard.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.95	0.96	0.94
7.50	7.75	7.50	7.48
Conc. (ng/mL)	Mean	SD	CV (%)
0.23	0.23	0.013	5.4
0.94	0.95	0.056	5.9
7.50	7.58	0.400	5.3



■ Accuracy

The accuracy was confirmed by 3 different lots testing ten times each different concentration.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.95	0.91	0.93
7.50	7.56	7.35	7.37
Conc. (ng/mL)	Mean	CV (%)	Recovery(%)
0.23	0.23	5.8	101.0
0.94	0.93	5.3	99.0
7.50	7.43	5.7	99.0

■ Comparability

Using Access2 (Beckman Coulter) as a comparison reference device for **AFIAS Tn-I Plus**, 100 clinical samples were independently tested for its Tn-I concentration following each instrument's procedure. Results of both the test methods were analyzed and their comparability was investigated with linear regression and coefficient of correlation (R).

REFERENCES

1. Mauro Panteghini, Franca Pacani, Kiang-Teck J.Yeo, Fred S. Apple, Robert H. Christenson. Francesco Dati, Johannes Mair, Jan Ravkilde, and Alan H.B. We. Evaluation of Imprecision for Cardiac Troponin Assays at Low-Range Concentrations. 2004;50:2:327-332.
2. Alan McNeil, PhD, FRACP, FRCPA. The Trouble with Troponin. Heart, Lung and Circulation 2007;16:S13-S16.
3. David M. Bunk and Micahel J. Welch. Characterization of a New Certified Reference Material for Human Cardiac Troponin I. Clinical Chemistry 2002;52:2:212-219
4. Jaffe AS, Ravkilde J, Roberts R, Naslund U, Apple FS, Galvani M, Katus H. It's time for a change to a troponin standard. Circulation 2000;102:1216-1220.
5. Jillan R. Tate, David Heathcote, Gus Koerbin, Gary Thean,

David Andriske, Jone Bonar, Janice Gill. The harmonization of cardiac troponin I measurement is independent of sample time collection but is dependent on the source of calibrator. Clinica Chimica Acta 324:2002:13-23.

6. Ohman EM, Armstrong PW, Christenson RH, et al. Cardiac troponin T levels for risk stratification in acute myocardial ischemia. N Engl J Med 1996;335:1333-41.
7. Antman EM, Tanasijevic MJ, Thompson B, et al. Cardiac specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996;335:1342-9.

Note: Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For technical assistance; please contact:

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
AFIAS Tn-I Plus
Cat. No. : SMFP-35

Classification:

Others (Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route:

Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

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I. General Information

- A. Product Name/Catalogue Number
: AFIAS Tn-I Plus / SMFP-35
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

II. Hazard Identification

GHS classification of the product and national or regional

	Cartridge
GHS classification	Not applicable
IMDG Code	Not applicable

GHS label elements including precautionary

	Cartridge
Symbols	No data available
Signal word	No data available
Hazard statements	No data available

Precautionary statements

	Cartridge
Prevention	No data available
Reaction	No data available
Storage	No data available
Disposal	No data available

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS433
	Rev. No.	04
	Rev. Date	2020. 09. 14

Other hazards

NEPA	Cartridge	
	Test strip	Housing
	Health: 1	
	Flammability: 1	
	Physical hazard: 0	

III. Composition/Information on Ingredients

[Cartridge]

	Component	Synonyms	Classification No.	Amount (%)
Test strip	Test line : Glutaraldehyde-Streptavidin	No data available	No data available	0.0000191
	Control line : Chicken IgY	No data available	No data available	0.0000097
Housing	Poly(methyl methacrylate)	PMMA	9011-14-7	55 ~ 70
	Acrylonitrile butadiene styrene	ABS	9003-56-9	28 ~ 45

IV. First Aid Measures

- In case of eye contact
 - : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
 - : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
 - : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring.
- If swallowed
 - : If irritation or signs of toxicity occur, seek medical attention.
- * The following symptoms may occur: Irreversible eye damage.
- * Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- * This material should be considered as being potentially infectious.

V. Fire Fighting Measures

- Recommended extinguishing media
 - : Suitable extinguishing agents including CO₂, WATER SPRAY or regular form.
- Specific hazard from the chemical
 - : Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS433
	Rev. No.	04
	Rev. Date	2020. 09. 14

- Special measure for fire-fighters

: Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.

- Protective Equipment and precautions for Fire fighters

: Chemical resisted protective equipment for fire-fighter.

VI. Accidental Release Measures

- Personal precautions

* If not danger, stop to leak

* Do not contact with exposed material without protective equipment

* Note the avoid conditions and materials

* Do not breathe dust / fume / gas / mist / vapours / spray

- Environmental precautions

: No data available

- Methods and materials for purification and cleaning up

: After absorbed exposed materials, clean up with cleaner and water

VII. Handling and Storage

- Precaution for safe handling

* Avoid skin friction.

* Note the avoid conditions and materials from fire or flame.

* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask

- Method for safe handling

* Note the avoid conditions and materials.

* Store container tightly closed in a well-ventilated area.

* Storage temperature: low temperature

VIII. Exposure controls and Personal protection

- Engineering management

: Good insulation and ventilation should be sufficient to control airborne levels.

- Personal protective equipment

Respiratory protection	Chemical resisted respiratory protection
Hand protection	Chemical resisted protective gloves
Eye protection	Chemical resisted eye protection
Skin and body protection	Chemical resisted protective clothes

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS433
	Rev. No.	04
	Rev. Date	2020. 09. 14

- Standard of exposure

Chemical	No data available
Biological	No data available

IX. Physical and Chemical properties

Cartridge		
	Test strip	Housing
State	No data available	Natural color or pigment processed color
Odor	Odorless	Little odorless
Odor threshold value	No data available	No data available
pH	No data available	Natural
Melting/ Freezing point	No data available	180 ~ 200°C/ No data available
Initial boiling point or boiling range	No data available	No data available
Flash point	No flash occurred under 90 °C (Closed-cup)	No data available
Evaporating rate	No data available	No data available
Lower explosion limit	No data available	No data available
Upper explosion limit	No data available	No data available
Vapor pressure	No data available	No data available
Water solubility	No data available	Not soluble in water, soluble in acetone
Density	No data available	1.11
Vapor density	No data available	No data available
Specific gravity	No data available (water=1)	No data available (water=1)
N-octan/water partition coefficient	No data available	No data available
Auto-ignition temperature	No data available	No data available
Decomposition temperature	No data available	No data available
Viscosity	No data available	No data available

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS433
	Rev. No.	04
	Rev. Date	2020. 09. 14

Molecular weight	No data available	50,000 ~ 250,000
------------------	-------------------	------------------

X. Stability and Reactivity

	Cartridge	
	Test strip	Housing
Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful	- Stable at high temperature and pressure
Conditions to avoid	Heat, high temperature	Heat, high temperature above 200°C
Materials to avoid	No data available	Oxidizing agent: fire or explosive risk
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency conditions	

XI. Toxicological Information

Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available
Aspiration hazard	No data available

XII. Ecological Information

Toxicity	No data available
----------	-------------------

MATERIAL SAFETY DATA SHEET	Document No.	BT-MSDS433
	Rev. No.	04
	Rev. Date	2020. 09. 14

Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	No data available

XIII. Disposal Considerations

Cartridge	Test strip	Not available
	Housing	Consider the notice on Waste management act

XIV. Transport Information

IMDG Code	No data available
DRF	No data available
RID	No data available
ADR	No data available
AND	No data available
IATA	Not regulated as dangerous goods
- Not a dangerous good under DOT, IATA, ADR, IMDG or RID	

XV. Regulatory Information

Korea Industrial Safety and Health Act	No data available
Korea Hazardous Materials Safety Control Act	No data available
Korea Toxic Chemical Control Act	No data available
Korea Wastes Control Act	No data available
Other internal and foreign acts	No data available

XVI. Other Notes

Limitations: The information and recommendations set forth in this MSDS are believed to be correct as of this date.

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