Cardiac

AFIAS  
NT-proBNP

INTENDED USE

**AFIAS NT-proBNP** is a fluorescence Immunoassay (FIA) for the quantitative determination of NT-proBNP in human whole blood/serum/plasma. It is useful as an aid in the diagnosis of persons suspected of having congestive heart failure.  
For *in vitro* diagnostic use only.

INTRODUCTION

N-terminal pro-brain natriuretic peptide (NT-proBNP) is produced predominately by the cardiac ventricular myocytes<sup>[1]</sup> and is released in response to myocardial stress and filling pressure<sup>[2]</sup> and is involved in maintaining intravascular volume homeostasis<sup>[3,4]</sup>. After stimulation of heart muscle cells, the natriuretic peptides are produced as prohormones (proBNP) and this is cleaved into two fragments which are secreted into the bloodstream as the 32 amino acids active BNP and the N-terminal fragment of 76 amino acids designated as NT-proBNP. NT-proBNP immunoassays are widely used and are now considered to be a useful marker and have a high degree of diagnostic accuracy in clinical practice and cardiovascular research as a diagnostic tool for the occurrence and severity of heart failure (HF) <sup>[5,6,7]</sup>. Therefore NT-proBNP measurements in human blood are helpful not only for the cardiac disease diagnosis but also for evaluation of patients with suspected HF and assessment of severity of the disease.

PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibodies in buffer bind to antigens in the sample, forming antigen-antibody complexes and migrate onto nitrocellulose matrix to be captured by the other immobilized antibodies on the test strip.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by instrument for AFIAS tests to show NT-proBNP concentration in the sample.

COMPONENTS

- AFIAS NT-proBNP** consists of a ‘cartridges’.
- Each sealed aluminum pouch contains two cartridges.
  - Each cartridge packaged in an aluminum pouch has three components including a cartridge part, detector part and detector diluent part.
  - Cartridge part contains the membrane called a test strip, which has streptavidin at the test line, and chicken IgY at the control line.
  - Detector part contains anti-NT-proBNP fluorescence conjugate, anti-chicken IgY-fluorescence conjugate, anti-NT-proBNP-biotin conjugate, bovine serum albumin (BSA), sucrose, bromophenol blue, tween 20 and sodium azide as a preservative in phosphate buffered saline (PBS).
  - Detector diluent part contains mouse IgG, bovine serum albumin (BSA), tween 20, CA-630, NaCl and sodium azide as a preservative in phosphate buffered saline (PBS).

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this ‘Instruction for use’.
- Use fresh samples only and avoid direct sunlight.
- It is possible to use frozen samples. Please refer to “SAMPLE COLLECTION AND PROCESSING.”
- Do not reuse cartridges and pipette tips. A cartridge should be used for testing one sample only.
- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange test components between different lots or use test components after the expiration date, either of which might yield incorrect test result(s).
- The cartridge should remain sealed in its aluminum pouch until just before use. Do not use the cartridge, if the pouch found damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local requirements. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- Allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used cartridges and pipette tips should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
- No Biotin interference was observed in **AFIAS NT-proBNP** when biotin concentration in the sample was below 5 ng/mL. If a patient has been taking biotin at dosage of more than 0.03 mg a day, it is recommended to test again 24 hours after discontinuation of biotin intake.
- AFIAS NT-proBNP** will provide accurate and reliable results subject to the below conditions.
  - AFIAS NT-proBNP** should be used only in conjunction with the instrument for AFIAS tests.
  - Have to use recommended anticoagulant sample.

Recommended anticoagulant
K <sub>2</sub> EDTA, K <sub>3</sub> EDTA, Sodium heparin, Lithium heparin

STORAGE AND STABILITY

Storage condition			
Component	Storage Temperature	Shelf life	Note
Cartridge	2 - 30 °C	20 months	Unopened
		1 month	Resealed

- Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal.

LIMITATION OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.

- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is the most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause the false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result should be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

MATERIALS SUPPLIED

REF SMFP-36	
Components of <b>AFIAS NT-proBNP</b>	
Cartridge Box Contains	
- Cartridge	24
- Pipette tip (Zipper bag)	24
- ID chip	1
- Instruction for use	1
- Spare cartridge zipperbag	1

MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from **AFIAS NT-proBNP**.

Please contact our sales division for more information.	
AFIAS-6	REF FPRR020
AFIAS-1	REF FPRR019
Boditech NT-proBNP Control	REF CFPO-245
Boditech NT-proBNP Calibrator	REF CFPO-271

SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS NT-proBNP** is human whole blood/serum/plasma.

- It is recommended to test the sample within 24 hours after collection.
- The serum or plasma should be separated from the clot by centrifugation within 3 hours after the collection of whole blood. If longer storage is required, e.g. if the test could not be performed within 24 hours, serum or plasma should be immediately frozen below -20 °C. The freezing storage of sample up to 3 months does not affect the quality of results.
- However, the whole blood sample should not be kept in a freezer in any case.
- As a repeated freeze-thaw cycle may affect the test result, do not refreeze previously frozen samples.

TEST SETUP

- Check the components of the **AFIAS NT-proBNP** as described below: Cartridges, Pipette tips, ID chip, Spare cartridge zipper bag and Instruction for use.
- If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for AFIAS tests

- Empty the tip box.
- Insert the ID chip into the “ID chip port”.  
(Please refer to the instrument for AFIAS tests ‘Operation Manual’ for complete information and operating instructions.

TEST PROCEDURE

- 1) Take 100 µL of sample with a pipette and dispense it into the sample well in the cartridge.  
2) Insert the cartridge into the cartridge holder  
3) Insert a tip into the tip hole of the cartridge.  
4) Tap the ‘START’ icon on the screen.  
5) The test result will be displayed on the screen after 12 minutes.

※ Note: Refer to the instrument for AFIAS tests Operation Manual to select sample type.

INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculates the test result automatically and displays NT-proBNP concentration of the test sample in terms of pg/mL.
- The working range of the AFIAS NT-proBNP is 10-30,000 pg/mL.
- The cut-off of the AFIAS NT-proBNP is 300pg/mL.

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with AFIAS NT-proBNP. For more information regarding obtaining the control materials, contact Boditech Med Inc.’s Sales Division for assistance.  
(Please refer to the instruction for use of control material.)

PERFORMANCE CHARACTERISTICS

- Analytical sensitivity**  
Limit of Blank (LoB): 5.94 pg/mL  
Limit of Detection (LoD): 8.03 pg/mL  
Limit of Quantitation (LoQ): 34.0 pg/mL
- Analytical specificity**
  - Cross-reactivity  
Biomolecules listed in the following table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS NT-proBNP** test results did not show any significant cross-reactivity with these biomolecules.

Cross-reactivity material	Conc. [ng/mL]
CK-MB	60
Troponin complex	1,000
Myoglobin	1,000
D-Dimer	1,000

- Interference

Interference materials listed in the following table were added to the test sample(s) the same as the below concentrations listed below. **AFIAS NT-proBNP** test results did not show any significant interference with these materials.

Interference material	Conc.
Bilirubin	350 µmol/L
Cholesterol	13 mmol/L
D-glucose	1,000 mg/dL
Hemoglobin	2 g/L
L-Ascorbic acid	350 µmol/L
Triglyceride mixture	500 mg/dL
Heparin	3000 U/L
EDTA	3.4 µmol/L
Biotin	5 ng/mL

■ Precision

- Between lot

One person tested three different lots of **AFIAS NT-proBNP**, ten times at each concentration of the control standard.

- Between person

Three different persons tested **AFIAS NT-proBNP**; ten times at each concentration of the control standard.

- Between day

One person tested **AFIAS NT-proBNP** during five days; five times at each concentration of the control standard.

- Between site

One person tested **AFIAS NT-proBNP** at three different sites; five times at each concentration of the control standard.

Conc. [pg/mL]	Between lot		Between person	
	AVG	CV (%)	AVG	CV (%)
234	235	6.0	233	6.6
1875	1890	6.0	1861	5.8
15000	15037	5.0	15105	6.2
Conc. [pg/mL]	Between day		Between site	
	AVG	CV (%)	AVG	CV (%)
234	233	5.9	230	6.5
1875	1835	4.8	1901	6.0
15000	14946	5.9	14997	5.6

■ Accuracy


The accuracy was confirmed by testing 3 different lots, ten times at each concentration of the control standard.

Mean of 3 lots			
Conc. [pg/mL]	Lot 1	Lot 2	Lot 3
234	231	234	236
1875	1947	1864	1984
15000	15406	14798	14903
Conc. [pg/mL]	Mean	CV (%)	Recovery (%)
234	234	6.3	99.9
1875	1932	5.5	103.0
15000	15036	6.0	100.2

■ Comparability

NT-proBNP concentrations of 100 clinical samples were quantified independently with **AFIAS NT-proBNP** and **Comparator A** as per prescribed test procedures. Test results were compared, and their comparability was investigated with linear regression and coefficient of correlations (R). Linear regression and coefficient of

correlation between the two tests were  $y=0.914x-0.438$ ,  $R = 0.98$ , respectively.

REFERENCES

1. A, Puschendorf B, Mair J. Cardiac natriuretic peptides: new laboratory parameters in heart failure patients. Clin Lab 2001; 47: 265-67
2. Maeda K, Tsutamoto T, Wada A, Hisanaga T, Kinoshita M. Plasma brain natriuretic peptide as a biochemical marker of high left ventricular end-diastolic pressure in patients with symptomatic left ventricular dysfunction. Am Heart J. 1998 May; 135(5 Pt 1):825-32.
3. Pfister R, Schneider CA. Natriuretic peptides BNP and NT-pro-BNP: established laboratory markers in clinical practice or just perspectives? Clin Chim Acta 2004; 349: 25-38.
4. Cowie M.R., Struthers A.D., Wood D.A., Coats A.S., Thompson S.G., PooleWilson P.A., et al. Value of natriuretic peptides in assessment of patients with possible new heart failure in primary care. Lancet. 1997 Nov 8;350(9088):1349-53.
5. Hobbs F.D., Davis R.C., Roalfe A.K., Hare R., Davies M.K., Kenkre J.E. Reliability of N-terminal pro-brain natriuretic peptide assay in diagnosis of heart failure: cohort study in representative and high risk community populations. BMJ. 2002 Jun 22;324(7352):1498.
6. Hogenhuis J, Voors AA, Jaarsma T, Hoes AW, Hillege HL, Kragten JA, van Veldhuisen DJ. Anaemia and renal dysfunction are independently associated with BNP and NT-proBNP levels in patients with heart failure. Eur J Heart Fail. 2007 Aug;9(8):787-94. Epub 2007 May 25.
7. Ewald B, Ewald D, Thakkestian A, Attia J. Meta-analysis of B type natriuretic peptide and N-terminal pro B natriuretic peptide in the diagnosis of clinical heart failure and population screening for left ventricular systolic dysfunction. Intern Med J 2008;38: 101-13.

**Note:** Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For technical assistance; please contact:

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3. Pfister R, Schneider CA. Natriuretic peptides BNP and NTpro- BNP: established laboratory markers in clinical practice or just perspectives? Clin Chim Acta 2004; 349: 25-38.

4. Cowie M.R., Struthers A.D., Wood D.A., Coats A.S., Thompson S.G., PooleWilson P.A., et al. Value of natriuretic peptides in assessment of patients with possible new heart failure in primary care. Lancet. 1997 Nov 8;350(9088):1349-53.

5. Hobbs F.D., Davis R.C., Roalfė A.K., Hare R., Davies M.K., Kenkre J.E. Reliability of N-terminal pro-brain natriuretic peptide assay in diagnosis of heart failure: cohort study in representative and high risk community populations. BMJ. 2002 Jun 22;324(7352):1498.

6. Hogenhuis J, Voors AA, Jaarsma T, Hoes AW, Hillege HL, Kragten JA, van Veldhuisen DJ. Anaemia and renal dysfunction are independently associated with BNP and NT-proBNP levels in patients with heart failure. Eur J Heart Fail. 2007 Aug;9(8):787-94. Epub 2007 May 25.

7. Ewald B, Ewald D, Thakkestian A, Attia J.Meta-analysis of B type natriuretic peptide and N-terminal pro B natriuretic peptide in the diagnosis of clinical heart failure and population screening for left ventricular systolic dysfunction. Intern Med J 2008;38: 101-13.

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
AFIAS NT-proBNP  
Cat. No. : SMFP-36

Classification:

Others (Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route:

Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer.

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## I. General Information

- A. Product Name/Catalogue Number  
: AFIAS NT-proBNP / SMFP-36
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

## II. Hazard Identification

GHS classification of the product and national or regional

	Cartridge
GHS classification	Not applicable
IMDG Code	Not applicable

GHS label elements including precautionary

	Cartridge
Symbols	No data available
Signal word	No data available
Hazard statements	No data available

Precautionary statements

	Cartridge
Prevention	No data available
Reaction	No data available
Storage	No data available
Disposal	No data available

Other hazards

	Cartridge
	Housing
	Health: 1
NEPA	Flammability: 1
	Physical hazard: 0

## MATERIAL SAFETY DATA SHEET

Document No.	BT-MSDS428
Rev. No.	05
Rev. Date	2021. 04. 15

### III. Composition/Information on Ingredients

[Cartridge]

	Component	Synonyms	Classification No.	Amount (%)
Housing	Poly(methyl methacrylate)	PMMA	9011-14-7	57 ~ 67
	Acrylonitrile butadiene styrene	ABS	9003-56-9	33 ~ 43

### IV. First Aid Measures

- In case of eye contact
  - : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
  - : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
  - : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring.
- If swallowed
  - : If irritation or signs of toxicity occur, seek medical attention.
- \* The following symptoms may occur: Irreversible eye damage.
- \* Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- \* This material should be considered as being potentially infectious.

### V. Fire Fighting Measures

- Recommended extinguishing media
  - : Suitable extinguishing agents including CO<sub>2</sub>, WATER SPRAY or regular form.
- Specific hazard from the chemical
  - : Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters
  - : Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters
  - : Chemical resisted protective equipment for fire-fighter.

### VI. Accidental Release Measures

- Personal precautions
  - \* If not danger, stop to leak
  - \* Do not contact with exposed material without protective equipment
  - \* Note the avoid conditions and materials
  - \* Do not breathe dust / fume / gas / mist / vapours / spray

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS428
	Rev. No.	05
	Rev. Date	2021. 04. 15

- Environmental precautions  
: No data available
- Methods and materials for purification and cleaning up  
: After absorbed exposed materials, clean up with cleaner and water

## VII. Handling and Storage

- Precaution for safe handling
  - \* Avoid skin friction.
  - \* Note the avoid conditions and materials from fire or flame.
- \* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
  - \* Note the avoid conditions and materials.
  - \* Store container tightly closed in a well-ventilated area.
- \* Storage temperature: low temperature

## VIII. Exposure controls and Personal protection

- Engineering management  
: Good insulation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

Respiratory protection	Chemical resisted respiratory protection
Hand protection	Chemical resisted protective gloves
Eye protection	Chemical resisted eye protection
Skin and body protection	Chemical resisted protective clothes

- Standard of exposure

Chemical	No data available
Biological	No data available

## IX. Physical and Chemical properties

	Cartridge
	Housing
State	Natural color or pigment processed color
Odor	Little odorless
Odor threshold value	No data available
pH	Natural
Melting/	180 ~ 200°C/

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS428
	Rev. No.	05
	Rev. Date	2021. 04. 15

Freezing point	No data available
Initial boiling point or boiling range	No data available
Flash point	No data available
Evaporating rate	No data available
Lower explosion limit	No data available
Upper explosion limit	No data available
Vapor pressure	No data available
Water solubility	Not soluble in water, soluble in acetone
Density	1.11
Vapor density	No data available
Specific gravity	No data available (water=1)
N-octan/water partition coefficient	No data available
Auto-ignition temperature	No data available
Decomposition temperature	No data available
Viscosity	No data available
Molecular weight	50,000 ~ 250,000

## X. Stability and Reactivity

	Cartridge
	Housing
Chemical stability	Stable at high temperature and pressure
Conditions to avoid	Heat, high temperature above 200°C
Materials to avoid	Oxidizing agent: fire or explosive risk
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency

## XI. Toxicological Information

Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS428
	Rev. No.	05
	Rev. Date	2021. 04. 15

Specific target organ toxicity (single exposure)	No data available
Specific target organ toxicity (repeated exposure)	No data available
Aspiration hazard	No data available

## XII. Ecological Information

Toxicity	No data available
Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	No data available

## XIII. Disposal Considerations

Cartridge	Housing	Consider the notice on Waste management act
-----------	---------	---

## XIV. Transport Information

IMDG Code	No data available
DRF	No data available
RID	No data available
ADR	No data available
AND	No data available
IATA	Not regulated as dangerous goods

Not a dangerous good under DOT, IATA, ADR, IMDG, or RID.

## XV. Regulatory Information

Korea Industrial Safety and Health Act	No data available
Korea Hazardous Materials Safety Control Act	No data available
Korea Toxic Chemical Control Act	No data available
Korea Wastes Control Act	No data available
Other internal and foreign acts	No data available

## XVI. Other Notes

Limitations: The information and recommendations set forth in this MSDS are believed to be correct as of this date.

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