# Boditech PCT Plus Control

## INTENDED USE

Boditech PCT Plus Control is intended for the quality control of PCT Plus Assay Kit manufactured by Boditech Med Inc.  
For *in vitro* diagnostic use only.

## COMPONENTS

Boditech PCT Plus Control consists of 'Boditech PCT Plus Control Level 1', 'Boditech PCT Plus Control Level 2', 'Instructions for use' and 'Control value & Barcode sheet'.

- Boditech PCT Plus Control is provided in lyophilized form.
- The control contains Recombinant Procalcitonin and Sodium azide in Horse serum.
- The control materials are contained in vials and the vial are further packaged in a box.

## SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech PCT Plus Control should not be used past the expiration date.
- Boditech PCT Plus Control is solely designed for the quality control of PCT Plus assay kits manufactured by Boditech Med Inc.
- All waste materials should be disposed of in accordance with the requirements of your local waste management authorities.

## STORAGE AND STABILITY

- Storage and stability condition of Boditech PCT Plus Control.

	Unopened	Opened (after reconstitution)	
Temperature	+2 to +8 °C	+2 to +8 °C	-20 to -80 °C
Expiration date	Until expiration date on the label.	1 day	7 days

- Close the opened vial tightly after use.
- After use, any residual substance should not be returned to the original vial.
- Bacterial contamination of reconstituted Boditech PCT Plus Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

## PROCEDURE

Boditech PCT Plus Control is supplied in lyophilized form.

1. Carefully reconstitute each vial of lyophilized control material with exactly 1 mL of sterilized distilled water.
2. Close the bottle and allow it to stand for 30 minutes before use. Ensure the contents are completely dissolved by swirling the gently. (To avoid formation of foam, do not shake the vial.)

Please refer to the instructions of use of the test cartridges for detailed test procedure.

In the event of damage to the package, contact the **Boditech Med Inc.'s Technical Services**.

## MATERIALS SUPPLIED

### REF CFPO-225

Boditech PCT Plus Control Box (2 vials)

- |  |   |
|--|---|
| - Boditech PCT Plus Control Level 1 (1 mL) | 1 |
| - Boditech PCT Plus Control Level 2 (1 mL) | 1 |
| - Instructions for use                     | 1 |
| - Control value & Barcode sheet            | 1 |

## INTERPRETATION OF THE RESULT

The test result of the 'Boditech PCT Plus Control' should be consistent with the expected result of control value sheet. If the test results fall outside the expected result, repeat the test.

※ Possible causes for erroneous test result

- Errors in testing process
- Incorrect storage condition of Boditech PCT Plus Control
- Use of expired or contaminated Boditech PCT Plus Control.
- Faulty Boditech's PCT Plus assay kits.
- Faulty Boditech's instruments.

## QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.

**Note:** Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For Technical Assistance

**Boditech Med Inc.'s Technical Services** at

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European Representative: OBELIS S.A  
Bd. Général Wahis 53,  
1030 Brussels,  
Belgium

Product: Boditech PCT Plus Control  
Cat. No. : CFPO-225

Classification: Others (Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer.

Standards applied: ISO 15223-1:2016, EN ISO 13485:2012, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2012,  
EN ISO 17511:2003, EN ISO 18113-1:2011, EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, February 7, 2019

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## I. General Information

- A. Product Name/Catalogue Number  
: Boditech PCT Plus Control / CFPO-225
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

## II. Hazard Identification

GHS classification of the product and national or regional

GHS classification	Chronic aquatic toxicity (Category 3), H412
IMDG Code	Not applicable

GHS label elements including precautionary

Symbols	No data available
Signal word	No data available
Hazard statements	H412 Harmful to aquatic life with long lasting effects.

Precautionary statements

Prevention	No data available
Reaction	No data available
Storage	No data available
Disposal	No data available

Other hazards

NEPA	No data available
------	-------------------

## MATERIAL SAFETY DATA SHEET

Document No.	BT-MSDSC51
Rev. No.	00
Rev. Date	2018. 08. 08

### III. Composition/Information on Ingredients

Component	CAS Number	Weight %
Horse serum	Not applicable	99.8
Sodium Azide	26628-22-8	0.1
Recombinant Procalcitonin	No data available	0.1

We recommend handling all chemicals with caution.

### IV. First Aid Measures

- In case of eye contact
  - : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
  - : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
  - : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring
- If swallowed
  - : If irritation or signs of toxicity occur, seek medical attention.
- \* The following symptoms may occur: Irreversible eye damage.
- \* Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- \* This material should be considered as being potentially infectious.

### V. Fire Fighting Measures

- Recommended extinguishing media
  - : Suitable extinguishing agents including CO<sub>2</sub>, WATER SPRAY or regular form.
- Specific hazard from the chemical
  - : Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters
  - : Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters
  - : Chemical resisted protective equipment for fire-fighter.

### VI. Accidental Release Measures

- Personal precautions
  - \* If not danger, stop to leak
  - \* Do not contact with exposed material without protective equipment
  - \* Note the avoid conditions and materials

## MATERIAL SAFETY DATA SHEET

Document No.	BT-MSDSC51
Rev. No.	00
Rev. Date	2018. 08. 08

- \* Do not breathe dust / fume / gas / mist / vapours / spray
- Environmental precautions
  - : No data available
- Methods and materials for purification and cleaning up
  - : After absorbed exposed materials, clean up with cleaner and water

### VII. Handling and Storage

- Precaution for safe handling
  - \* Avoid skin friction.
  - \* Note the avoid conditions and materials from fire or flame.
- \* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
  - \* Note the avoid conditions and materials.
  - \* Store container tightly closed in a well-ventilated area.
- \* Storage temperature: low temperature

### VIII. Exposure controls and Personal protection

- Engineering management
  - : Good insulation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

Respiratory protection	Chemical resisted respiratory protection
Hand protection	Chemical resisted protective gloves
Eye protection	Chemical resisted eye protection
Skin and body protection	Chemical resisted protective clothes

- Standard of exposure

Chemical	No data available
Biological	No data available

### IX. Physical and Chemical properties

State	Liquid at 20 °C
Odor	No data available
Odor threshold value	No data available
pH	7.83
Melting/Freezing point	No data available
Initial boiling point or boiling range	No data available
Flash point	No flash occurred under 90 °C (Closed-cup)

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDSC51
	Rev. No.	00
	Rev. Date	2018. 08. 08

Evaporating rate	No data available
Lower explosion limit	No data available
Upper explosion limit	No data available
Vapor pressure	No data available
Water solubility	No data available
Density	No data available
Vapor density	No data available
Specific gravity	1.01 (water=1)
N-octan/water partition coefficient	No data available
Auto-ignition temperature	No data available
Decomposition temperature	No data available
Viscosity	No data available
Molecular weight	No data available

## X. Stability and Reactivity

Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful
Conditions to avoid	Heat, high temperature
Materials to avoid	No data available
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency

## XI. Toxicological Information

	Detection buffer
Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDSC51
	Rev. No.	00
	Rev. Date	2018. 08. 08

Reproductive toxicity	No data available
Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available
Aspiration hazard	No data available

## XII. Ecological Information

Toxicity	No data available
Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	Harmful to aquatic life with long lasting effects.

## XIII. Disposal Considerations

Not available

## XIV. Transport Information

IMDG Code	No data available
DRF	No data available
RID	No data available
ADR	No data available
AND	No data available
IATA	Not regulated as dangerous goods

Not a dangerous good under DOT, IATA, ADR, IMDG, or RID.

## XV. Regulatory Information

Korea Industrial Safety and Health Act	No data available
Korea Hazardous Materials Safety Control Act	No data available
Korea Toxic Chemical Control Act	No data available
Korea Wastes Control Act	No data available
Other internal and foreign acts	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDSC51
	Rev. No.	00
	Rev. Date	2018. 08. 08

#### XVI. Other Notes

Limitations: The information and recommendations set forth in this MSDS are believed to be correct as of this date.

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