Boditech Tn-I Plus Control

INTENDED USE

Boditech Tn-I Plus Control is intended for *in vitro* diagnostic use in the quality control of Tn-I Assay Kit. **For *in vitro* diagnostic use only.**

SUMMARY AND PRINCIPLE

The use of Boditech Tn-I Plus Control may be considered as an objective assessment of the precision of Tn-I Assay Kits and is an integral part of Good Laboratory Practices. Boditech Tn-I Plus Control is provided in lyophilized form.

COMPONENTS

Boditech Tn-I Plus Control consists of 'Boditech Tn-I Plus Control level 1', 'Boditech Tn-I Plus Control level 2', 'Instruction for Use' and 'Control value & Barcode Sheet'.

- The control contains Human Cardiac Troponin I/C Complex Antigen and horse serum
- Each control vial packed in a box.

SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech Tn-I Plus Control should not be used past the expiration date.
- Boditech Tn-I Plus Control is solely designed to be provided instrument-specific calibration curves of Boditech Readers and Tn-I Assay Kits.
- Human source materials from which Boditech Tn-I Plus Control is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

STORAGE AND STABILITY

- Storage and stability condition of Boditech Tn-I Plus Control.

| Unopened | |
|-------------|-------------------------------------|
| Temperature | Expiration date |
| +2 to +8 °C | Until expiration date on the label. |

| Opened (After reconstitution) | |
|-------------------------------|-----------------|
| Temperature | Expiration date |
| +2 to +8 °C | 1 day |
| -20 to -80 °C | 7 days |

- Close the opened Boditech Tn-I Plus Control bottle tightly after use.
- Once the Boditech Tn-I Plus Control was frozen, it should be used ONE TIME ONLY for test, because repeated freezing and thawing can result in the change of test values.
- After use, any residual product should NOT BE RETURNED to the original vial.

- Bacterial contamination of reconstituted Boditech Tn-I Plus Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

INSTRUCTIONS FOR USE

Boditech Tn-I Plus Control is supplied in lyophilized form.

1. Carefully reconstitute each vial of lyophilized with exactly 1 mL of sterilized distilled water.
2. Close the bottle and allow to stand for 30 minutes before use. Ensure contents are completely dissolved by swirling gently. Avoid formation of foam. Do not shake.

Please refer to package inserts of the test cartridges for detailed test procedure.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

In the event of damage to the package, contact the **Boditech Med Inc.'s Technical Services.**

MATERIALS SUPPLIED

| REF | CFPO-212 | 2.6 t.s. |
|-----|---|----------|
| | Boditech Tn-I Plus Control Box (2 vials) | |
| | Boditech Tn-I Plus Control level 1 (1 mL) | 1 |
| | Boditech Tn-I Plus Control level 2 (1 mL) | 1 |
| | Instruction for Use | 1 |
| | Control value & Barcode Sheet | 1 |

QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.

For Technical Assistance

Boditech Med Inc.'s Technical Services at

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
OBELIS S.A
Bd. Général Wahis 53,
1030 Brussels,
Belgium

Product:

Boditech Tn-I Plus Control
Cat. No. : CFPO-212

Classification:

Others (Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route:

Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

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I. General Information

- A. Product Name/Catalogue Number
: Boditech Tn-I Plus Control / CFPO-212
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
 - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer
: Boditech Med Inc.
- D. Address
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.
: +82-33-243-1400

II. Hazard Identification

GHS classification of the product and national or regional

| | |
|--------------------|----------------|
| GHS classification | Not applicable |
| IMDG Code | Not applicable |

GHS label elements including precautionary

| | |
|-------------------|-------------------|
| Symbols | No data available |
| Signal word | No data available |
| Hazard statements | No data available |

Precautionary statements

| | |
|------------|-------------------|
| Prevention | No data available |
| Reaction | No data available |
| Storage | No data available |
| Disposal | No data available |

Other hazards

| | |
|------|-------------------|
| NEPA | No data available |
|------|-------------------|

MATERIAL SAFETY DATA SHEET

| | |
|--------------|--------------|
| Document No. | BT-MSDSC42 |
| Rev. No. | 00 |
| Rev. Date | 2018. 01. 26 |

III. Composition/Information on Ingredients

| Component | CAS Number | Weight % |
|----------------------|------------|----------|
| Bovine Serum Albumin | N/A | 99.9% |
| Sodium Azide | 26628-22-8 | 0.1% |

We recommend handling all chemicals with caution.

IV. First Aid Measures

- In case of eye contact
 - : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
 - : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
 - : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring
- If swallowed
 - : If irritation or signs of toxicity occur, seek medical attention.
- * The following symptoms may occur: Irreversible eye damage.
- * Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- * This material should be considered as being potentially infectious.

V. Fire Fighting Measures

- Recommended extinguishing media
 - : Suitable extinguishing agents including CO₂, WATER SPRAY or regular form.
- Specific hazard from the chemical
 - : Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters
 - : Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters
 - : Chemical resisted protective equipment for fire-fighter.

VI. Accidental Release Measures

- Personal precautions
 - * If not danger, stop to leak
 - * Do not contact with exposed material without protective equipment
 - * Note the avoid conditions and materials
 - * Do not breathe dust / fume / gas / mist / vapours / spray

MATERIAL SAFETY DATA SHEET

| | |
|--------------|--------------|
| Document No. | BT-MSDSC42 |
| Rev. No. | 00 |
| Rev. Date | 2018. 01. 26 |

- Environmental precautions
: No data available
- Methods and materials for purification and cleaning up
: After absorbed exposed materials, clean up with cleaner and water

VII. Handling and Storage

- Precaution for safe handling
 - * Avoid skin friction.
 - * Note the avoid conditions and materials from fire or flame.
- * Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
 - * Note the avoid conditions and materials.
 - * Store container tightly closed in a well-ventilated area.
- * Storage temperature: low temperature

VIII. Exposure controls and Personal protection

- Engineering management
: Good insulation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

| | |
|--------------------------|--|
| Respiratory protection | Chemical resisted respiratory protection |
| Hand protection | Chemical resisted protective gloves |
| Eye protection | Chemical resisted eye protection |
| Skin and body protection | Chemical resisted protective clothes |

- Standard of exposure

| | |
|------------|-------------------|
| Chemical | No data available |
| Biological | No data available |

IX. Physical and Chemical properties

| | |
|--|--|
| State | Liquid at 20 °C |
| Odor | No data available |
| Odor threshold value | No data available |
| pH | 7.83 |
| Melting/Freezing point | No data available |
| Initial boiling point or boiling range | No data available |
| Flash point | No flash occurred under 90 °C (Closed-cup) |
| Evaporating rate | No data available |

| | | |
|-----------------------------------|--------------|--------------|
| MATERIAL SAFETY DATA SHEET | Document No. | BT-MSDSC42 |
| | Rev. No. | 00 |
| | Rev. Date | 2018. 01. 26 |

| | |
|-------------------------------------|-------------------|
| Lower explosion limit | No data available |
| Upper explosion limit | No data available |
| Vapor pressure | No data available |
| Water solubility | No data available |
| Density | No data available |
| Vapor density | No data available |
| Specific gravity | 1.01 (water=1) |
| N-octan/water partition coefficient | No data available |
| Auto-ignition temperature | No data available |
| Decomposition temperature | No data available |
| Viscosity | No data available |
| Molecular weight | No data available |

X. Stability and Reactivity

| | |
|---|--|
| Chemical stability | - Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful |
| Conditions to avoid | Heat, high temperature |
| Materials to avoid | No data available |
| Hazardous decomposition products formed under fire conditions | Irritation or toxic gases and fume may be occurred in fire-emergency |

XI. Toxicological Information

| | |
|--|-------------------|
| | Detection buffer |
| Information on the likely routes of exposure | No data available |
| Acute toxicity | No data available |
| Skin corrosion/Irritation | No data available |
| Serious eye damage/Eye irritation | No data available |
| Respiratory/Skin sensitization | No data available |
| Germ cell mutagenicity | No data available |
| Carcinogenicity | No data available |
| Reproductive toxicity | No data available |

| | | |
|-----------------------------------|--------------|--------------|
| MATERIAL SAFETY DATA SHEET | Document No. | BT-MSDSC42 |
| | Rev. No. | 00 |
| | Rev. Date | 2018. 01. 26 |

| | |
|---|-------------------|
| Specific target organ toxicity(single exposure) | No data available |
| Specific target organ toxicity(repeated exposure) | No data available |
| Aspiration hazard | No data available |

XII. Ecological Information

| | |
|-------------------------------|-------------------|
| | Detection buffer |
| Toxicity | No data available |
| Persistence and degradability | No data available |
| Bio-accumulative potential | No data available |
| Mobility in soil | No data available |
| Other adverse effects | No data available |

XIII. Disposal Considerations

Not available

XIV. Transport Information

| | |
|-----------|----------------------------------|
| IMDG Code | No data available |
| DRF | No data available |
| RID | No data available |
| ADR | No data available |
| AND | No data available |
| IATA | Not regulated as dangerous goods |

Not a dangerous good under DOT, IATA, ADR, IMDG, or RID.

XV. Regulatory Information

| | |
|--|-------------------|
| Korea Industrial Safety and Health Act | No data available |
| Korea Hazardous Materials Safety Control Act | No data available |
| Korea Toxic Chemical Control Act | No data available |
| Korea Wastes Control Act | No data available |
| Other internal and foreign acts | No data available |

| | | |
|-----------------------------------|--------------|--------------|
| MATERIAL SAFETY DATA SHEET | Document No. | BT-MSDSC42 |
| | Rev. No. | 00 |
| | Rev. Date | 2018. 01. 26 |

XVI. Other Notes

Limitations: The information and recommendations set forth in this MSDS are believed to be correct as of this date.