

Space^{plus} Perfusor®

Instructions for use



CE 0123

Table of contents

1	About these instructions for use	5	8	Display screens.....	22
1.1	Purpose of these instructions for use	5	8.1	General status displays	22
1.2	Warnings.....	6	8.2	Home menu.....	22
1.3	Abbreviations.....	6	8.3	Menus	22
2	Symbols	7	8.3.1	Drug data menu.....	23
2.1	Symbols on the product and packaging.....	7	8.3.2	Patient data menu	23
2.2	Status indicators on the display and housing.....	8	8.3.3	Infusion data menu	23
3	Intended purpose	9	8.4	Run screen.....	23
4	Safety instructions.....	10	8.4.1	Syringe run screen.....	24
4.1	General use.....	10	9	Main functions.....	25
4.2	Inspection on delivery	10	9.1	Switching on, switching off and standby	25
4.3	Before and during use.....	10	9.2	Programming an infusion	25
4.4	Cleaning	11	9.3	Starting an infusion	25
4.5	Protecting patients.....	11	9.4	Pausing/ending an infusion	25
4.6	Alarm volume and patient call	12	9.5	Locking/unlocking the display	26
4.7	Setting up the device.....	12	9.6	Administering a bolus.....	26
4.8	Device environment.....	13	9.6.1	The bolus menu.....	26
4.9	Using the syringe correctly.....	13	9.6.2	Programming a bolus.....	26
4.10	Transporting the device	14	9.6.3	Reusing a programmed bolus	27
4.11	Avoiding damage to the device.....	14	9.6.4	Administering a manual bolus	27
4.12	Accessories	14	9.6.5	Setting the standard bolus rate.....	27
4.13	Maintenance, servicing and spare parts	15	9.6.6	The bolus run screen	27
4.14	Software and updates	15	9.7	Changing the syringe.....	27
4.15	Safety standards.....	15	9.7.1	Opening the claws manually	28
5	Device overview	16	9.8	Priming the line	28
5.1	Front.....	16	9.9	Balance Et info	29
5.2	Back	18	9.9.1	Infused Volume.....	29
6	Ports and connections.....	19	9.9.2	Infusion status	29
6.1	Attaching/detaching the pump clamp.....	19	9.9.3	Therapy history	30
6.2	Connecting/disconnecting pumps.....	19	9.10	Keep vein open (KVO mode).....	30
6.3	Connecting cables.....	20	9.11	Changing the device settings.....	31
6.3.1	Power cable.....	20	9.11.1	Activating or deactivating the code lock	31
6.3.2	Accessory cable for bolus button and patient call.....	20	9.11.2	Selecting pressure level.....	32
7	Operating the device.....	21	10	The drug library.....	33
7.1	Categories.....	21	10.1	Accessing the drug library	33
7.2	Entering values	21	10.2	Programming an infusion with the library	33
7.3	Deleting values.....	21	10.3	The limits.....	35
7.4	Disabled controls.....	21	10.4	The review screen.....	36
			10.5	Changing the drug data during an infusion ...	36

Table of contents

11	Description of the infusion profiles.....	37	24	Electromagnetic compatibility	59
11.1	Rate/volume/time	37	24.1	Electromagnetic interference emissions.....	61
11.1.1	Overview of infusion parameters.....	37	24.2	Electromagnetic immunity.....	62
11.1.2	Changing the infusion parameters.....	37	24.3	Recommended safe distances.....	66
11.2	Dose calculation	37	25	Instructions for use for accessories.....	67
12	Hints & alarms	39	25.1	12 V interface cable (871923112)	67
12.1	Hints	39	25.2	Staff call interface cable (8718031)	67
12.2	Service hints	39	26	Ordering data.....	69
12.3	Alarm display.....	39	26.1	Infusion pump	69
12.4	Alarm priorities.....	39	26.2	Interface cable	69
12.5	Alarm types.....	40	26.3	Recommended accessories	69
12.5.1	Notifications.....	40	26.4	Power cables.....	69
12.5.2	Reminders.....	41	26.5	Disposables.....	70
12.5.3	Pre-alarms.....	41	26.5.1	Syringes.....	70
12.5.4	Operating alarms.....	41	Index		72
12.5.5	Device alarms	43			
13	Software & updates	44			
13.1	Updating the drug library.....	44			
13.2	Activating additional updates	44			
14	Cleaning & disinfection	44			
15	Battery mode	46			
15.1	Notes for optimal battery operation	46			
16	Decommissioning	46			
17	Warranty	46			
18	Maintenance and repair	47			
19	Disposal	47			
20	Technical Safety Check (TSC)/service.....	47			
20.1	Electrical Safety Inspection.....	47			
21	Start-up and trumpet curves.....	48			
21.1	Significance for clinical practice	48			
21.2	Typical start-up and trumpet curves.....	49			
22	Time to alarm	51			
22.1	Omnifix® 50 mL.....	51			
22.2	Compatible syringes	52			
23	Technical data	53			
23.1	Pump.....	53			
23.2	WiFi interface.....	57			

About these instructions for use

1 About these instructions for use

1.1 Purpose of these instructions for use

These instructions for use are part of the device and describe how to use the device safely and correctly.

- Read these instructions for use before using the pump!
- Keep the instructions for use near the pump!
- Please read and follow the accompanying documents!

Main functions

Therefore, in KVO mode, a very low, non-therapeutic rate is delivered.

If the KVO mode is activated in the device settings, the pump asks you whether you would like to start KVO after the programmed infusion time or the programmed volume.

The KVO rate depends on the rate of infusion and can vary depending on the configuration.

Menu	Meaning
Language	Select language
Date & time	Set the date and time

9.11 Changing the device settings

You can change the device settings via the [Home menu](#).

Which settings you can change depends on the configuration of the pump and may vary.

Menu	Meaning
	Activate or deactivate code lock
	Activate or deactivate KVO
e	Set the standard bolus rate
Plevel	Set the cut-off pressure (from)
Audio volume	Set audio volume f
Display ss	Set the brightness f

1;2;3 p.d. 2.8; 1.8

10.2 Programming an infusion with the library

Once you have accessed the drug library, the pump will guide you step-by-step through the following settings:

- Care unit
- Pump location
- Drug
- Patient profile
- Concentration of the drug
- Infusion profile



patient by ensuring

1;2;3 p.d. 1.12; 1.15



12.4 Alarm priorities

The pump distinguishes between two types of alarm priorities.

High priority alarms are visually indicated by a **red** alarm window and the flashing of the red **LED**.

Lower-priority alarms will be marked with a **yellow** alarm window and the steady illumination of the **yellow** LED.

In addition, the acoustic signals of the two priorities are different from one other.

Depending on the configuration, the alarm can be forwarded via the staff call.

The different alarm types and priorities are described below:

Alarm type	Alarm colour (LED & alarm window)
Notification	Yellow
Reminders	Yellow /grey
Pre-alarm	Yellow
Operating alarm	Red
Device alarm	Red

Hints & alarms

1;2;3 p.d. 1.12

12.5 Alarm types

12.5.1 Notifications

Some notifications act like alarms because they have to be given special attention.

A notification provides information about an event on the pump. This may be directly related to delivery or to a basic function.

If a notification is generated during delivery, the infusion will continue.

Display notification	Meaning
----------------------	---------

Battery is

Display notification	Meaning
----------------------	---------

Display touched for g

The display has been touched for more than 15 seconds without interruption.

Pump too hot	The pump has detected that the temperature is too high.
--------------	---

Pump too cold	
---------------	--

- Put the pump in stand-

Hints & alarms

Display notification	Meaning
KVO finished	KVO time has elapsed. <ul style="list-style-type: none"> – End the infusion or continue the infusion.
Syringe empty	No volume is left in the syringe. Due to the different tolerances of syringes from different manufacturers, some fluid may remain in the syringe. Restarting causes the syringe to be completely emptied until a pressure alarm is triggered. <ul style="list-style-type: none"> – Change the syringe or end the infusion.
Pressure too high	There is an occlusion in the infusion line. The set pressure level has been exceeded. Due to the different tolerances of syringes from third-party manufacturers, a pressure alarm can be triggered by high frictional forces on the syringe. The pump automatically reduces the post occlusion bolus. <ul style="list-style-type: none"> – Check that the line is not bent or damaged, that all connections are open and that all the filters are unblocked.

Display notification	Meaning
Code entry required	Someone has stopped the pump without entering the code or has tried to switch it off while the code lock was active. <ul style="list-style-type: none"> – Enter the correct code.
Syringe holder open	The syringe holder was opened while an infusion was in progress. <ul style="list-style-type: none"> – Close the syringe holder.
Syringe not fixed correctly	The syringe is not fixed correctly. <ul style="list-style-type: none"> – Replace the syringe.
Pump too hot	The pump has detected that the battery temperature is too high. <ul style="list-style-type: none"> – Lower the ambient temperature
Pump too cold	The pump has detected that the battery temperature is too low. <ul style="list-style-type: none"> – Increase the ambient temperature.
Battery empty	The battery is discharged. The battery alarm is signalled for 3 min. Then the pump will automatically turn off. <ul style="list-style-type: none"> – Connect the pump to the mains power.

Technical data

23 Technical data

Atitikimas direktyvai
2017/745

23.1 Pump

The technical data applies to all operation locations listed in the intended purpose including home use as well as emergency and transport situations (road ambulances, fixed-wing and rotary-wing air ambulances).

Parameter	Value
Type of device	Infusion syringe pump
Product classification	according to Directive 93/42 EEC and Regulation (EU) 2017/745:
Protection against liquids and particles	IP44
Power supply	
- Mains power supply	
- 12V supply	
Power consumption	
- Mains power supply	
- 12V supply	
Internal battery	
- Battery life	
- Recharging time	

Index

A

Abbreviations 6
Accessories 67
Accessory cable 20

B

Back 18
Balance & info 29
Battery mode 46, 53
Battery operation and maintenance 46
Bolus, administration 26

C

Cleaning 11
Code lock, activating/deactivating 31
Compatible syringes 52
Connections 19

D

Decommissioning 46
Delivery accuracy 55
Display, locking/unlocking 26
Disposal 47
Drug library 33

E

Electromagnetic compatibility 59
Electromagnetic immunity 62
Electromagnetic interference emissions 61

F

Front 16

H

Hints & alarms 39
Home menu 22

I

Infusion profiles 37
Infusion, programming 25

Intended purpose 9

Interface cable, 12 V CP 67, 69
Interface cable, staff call CP 67, 69
Interfaces 54

K

KVO Mode 30

M

Maintenance 15
Menus 22

O

Operating conditions 54
Operation 21
Ordering data 69

R

Recommended safe distances 66
Repair 47
Run screen 23

S

Safety instructions 10
Safety standards 15
Service 47
Servicing 47
Setting-up 12
Software & updates 44
Staff call 54
Start-up curves 48
Status displays 8
Storage conditions 54
Symbols 7
Symbols on the product
and packaging 7
Syringe, changing 27

T

Technical data 53, 68

Temperature 54

Time to alarm 51

Transporting 14

Trumpet curves 48

W

Warnings 6

Warranty 46

Manufacturer:
B. Braun Melsungen AG
34209 Melsungen
Germany
Tel +49(0) 56 61 71-0
www.bbraun.com

Sales:
B. Braun Melsungen AG
Hospital Care Division
34209 Melsungen
Germany
Tel +49(0) 56 61 71-0
Fax: +49(0) 56 61 71-20 44
www.bbraun.de



39012107
2021-01-29
Date of last revision: January 2021