Discovery XR656 HD

Digital Radiographic System powered by Helix™

Key Features

- Helix 2.0 image processing
- AutoRad workflow automation
- VolumeRAD with Metal Artifact Reduction
- Dual Energy
- Auto Image Paste with Auto Spine
- FlashPad HD high resolution detectors
- Live video streaming
- Bariatric elevating table
- iCenter asset management
- Stand-alone Helix Workstation

Discovery XR656 HD, GE's advanced digital radiographic system, powered by Helix™ 2.0 advanced image processing helps you achieve clinical excellence with ease and efficiency.

Designed to handle 2D and 3D exams via standard and advanced radiographic applications, various system configurations are available to align to your department's radiographic requirements. The main configurations include:

- Table and wall stand (standard or extended arm)
- Wall stand only system (standard or extended arm)
- Table only system
- OTS only system (without Table and Wall stand)

Table of contents

Introduction.....	3
Advanced Applications	5
FlashPad HD Detectors.....	7
Helix™ 2.0 Advanced Image Processing.....	8
Workflow Features	9
OTS User Interface	11
Optional Workflow Features	11
Stand-Alone Helix™ Workstation†	12
Acquisition Workstation	12
DICOM Connectivity	13
Network Connectivity, Remote Service and Security	14
Overhead Tube Suspension.....	15
Elevating Table.....	16
Digital Tilting-Arm Wall Stand.....	17
Generator.....	18
X-ray Tube.....	18
Dose Monitoring and DAP	19
Automatic Collimator	19
Accessories†	20
Environmental Conditions	21
Installation Power.....	21
Typical Room Layout.....	22
Compliance to Standards.....	23
Warranty	23
About GE Healthcare	24

Introduction

Outstanding image quality & diagnostic confidence at low dose. Effortless precision.
Powered by Helix™ advanced image processing.


FlashPad HD high resolution wireless digital detectors that capture extraordinary detail at low dose with a DQE of 75% at 0 lp/mm and 100 micron, 5.0 lp/mm resolution

Helix™ 2.0 advanced image processing delivers extraordinary anatomical detail and consistent performance despite variations in patient size, exposure techniques, collimation, and presence of metal implants

Live streaming camera⁺ Live video of the patient allows technologists to stay connected to their patients, monitor patient safety, and potentially reduce rejects from patient motion or incorrect orientation.


AutoRAD⁺ workflow automation suite featuring Auto Positioning, Auto Tracking, Auto Protocol Assist, Auto Field-of-View, Auto QA, Quick Charge, Quick Share, and Quick Connect.

VolumeRAD with Metal Artifact Reduction⁺


Digital Tomosynthesis provides multi-level image slices provides similar data to CT at very low doses. Metal Artifact Reduction (MAR) reduces undershoot and multi-plane streaks improving image quality and clinical confidence while allowing more accurate placement of screws, cages, or plate fusions.

Dual Energy Subtraction⁺ Beyond image processing, a dual energy chest exam is performed by acquiring two images at different energy levels separated by less than 160 milliseconds. Three images are generated, leading to the detection of abnormalities that may have been obscured in a conventional radiograph. - - -

Auto Image Paste⁺ enhanced with Auto Spine for seamless long bone and scoliosis exams offering outstanding precision and speed at the table and wall stand.

1.1.30

Improve operational efficiency with powerful analytics and advanced service technologies that never sleep. Discovery XR656 HD is compatible with:

GE's X-ray Quality Application featuring Repeat Reject Analytics⁺ provides web-based dashboards to manage quality assurance, uncover the root cause of rejected X-ray exams, plan targeted training and help reduce unnecessary radiation dose.

iCenter asset management⁺ software platform designed to help you optimize utilization of your X-ray equipment and balance workload using the full power of healthcare data analytics

InSite™ remote service support. Remote diagnostics and troubleshooting for fast resolutions, often without a field engineer visit.

The standalone **Helix workstation⁺** is an additional workstation that provides access to the patient worklist, image reprocessing and QC functions on a separate workstation to free up the system for subsequent patient acquisitions. Compatibility to reprocess images from all FlashPad HD-based systems.

A full range of accessories are also available, such as: mobile stretchers, weight bearing exam tools, patient table accessories and grids.

Advanced Applications

VolumeRAD[†]


VolumeRAD Digital Tomosynthesis brings 3-dimensional imaging into your radiography room, helping resolve inconclusive x-rays "here and now" with minimal added radiation. This cost-effective technology improves patient experience by providing additional information that could assist to make care decisions faster.

Metal Artifact Reduction (MAR) for VolumeRAD reduces undershoot and multi-plane streaks outside the area of interest improving image quality and clinical confidence while allowing the radiologist to more accurately determine placement of screws, cage, or plate fusions.

VolumeRAD Advanced Application	
Minimum FOV	10 cm x 10 cm
Number of projections	25 - 60
Maximum sweep angle	+/- 20 degrees
Sweeps available	Table horizontal, wall stand cross-lateral, wall stand vertical, wall stand horizontal
Acquisition resolution	200 um
Thoracic dose	0.09 mSv
Acquisition time	15 Sec
Reconstruction time w/o MAR	60 Sec
Reconstruction time w/ MAR	80 Sec

VolumeRAD has been shown to:

- Detect wrist fractures with 95% sensitivity and occult wrist fractures with 42%¹ sensitivity compared to standard X-ray²
- Decrease the need for additional MRI/CT exams by approximately 15% for suspected hip fractures²
- Reduce the need for an additional MR or CT exam by up to 50% if coupled with an additional clinical exam for hip fractures²
- Lower the patient dose by 95% as compared to a CT exam of a wrist fracture²
- Reduce the need for additional CT imaging by up to 80% when used as an adjunct for inconclusive chest exams²

¹ Based on the meta-analysis of three peer-reviewed scientific publications (on file). The average sensitivity was 42%.


² Based on the meta-analysis of peer-reviewed scientific publications (on file).

Dual Energy Subtraction†

Beyond image processing, the attenuation of bone and soft tissue to different energy x-rays are directly measured. With one acquisition containing two exposures at different energy levels, the standard, bone, and soft-tissue images are created. This allows the radiologist to get additional information on calcified structures and foreign objects.

This technology is typically applied to chest and abdomen anatomies.

Dual Energy Advanced Application	
Time between Images	160 msec
Low-Energy image (typical)	60-80 kVp
High-Energy image (typical)	110-150 kVp
Dose of chest image relative to a single-energy image	120%

Auto Image Paste†

Acquire multiple images in one fast, seamless, highly automated exam to generate images that are larger than the detector size.


Auto Image Paste has been enhanced with **AutoSpine**, an intelligent algorithm that follows the contour of the spine for vertical equalization enabling a natural balance of brightness & contrast along the patient body in lateral spine exams

Auto Image Pasting Advanced Application		
	Table	Wall Stand
Maximum coverage	110 cm	150 cm
Number of images	2 to 4	2 to 5
Average acquisition time (3 images, 90 cm coverage)	≤10 sec	
Time from acquisition start to final pasted image	≤22 sec	

1.1.28

Auto Image Paste includes imaging of the spine for scoliosis evaluation and imaging of the legs for orthopedic evaluations. A patient stand comfortable during acquisition at the wall stand.

Auto Image Pasting Supported Anatomies		
	Table	Wall Stand
Spine AP	✓	✓
Spine PA	✓	✓
Spine Lateral	✓	✓
Leg AP	✓	✓
Leg PA		✓

FlashPad HD Detectors

Digital image acquisition supports fast and efficient exam procedures, eliminating time spent handling film and cassettes, as well as reliability issues inherent in cassette tray systems, thus helping to reduce overall exam times and improve patient satisfaction.

At the core of Discovery XR656 HD is GE's **FlashPad HD wireless flat panel detectors** that provide outstanding image quality at low dose for general radiography applications.

	HD 4343 (17"x17")	HD 3543 (14"x17")	HD 2530 (10"x12")	
Detector technology	Single panel (non-tiled) amorphous silicon with a Cesium Iodide (CSI) scintillator			1.1.59
Dimensions	ISO 4090 cassette size 460 x 460 x 15.5 mm	ISO 4090 cassette size 384 x 460 x 15.5 mm	ISO 4090 cassette size 282 x 332 x 15.5 mm	
Pixel matrix	4288 x 4288	3524 x 4288	2508 x 3004	
Weight (w/ Battery)	3.8 kg (8 lbs)	3.2 kg (7 lbs)	1.8 kg (4 lbs)	
Battery operation (@ 60 exposures/hr)	3.5 hours	4 hours	2 hours	

FlashPad HD combines decades of expertise in designing digital detectors to provide fine detail with low-dose operation.

Detector imaging performance			
Pixel size	100 μ m	T4	1.1.62
Limiting resolution	5 lp/mm		
DQE (typical, RQA5)	75% @ 0lp/mm 60% @ 1lp/mm 40% @ 3lp/mm		
MTF (typical, RQA5)	70% @ 1lp/mm 40% @ 2lp/mm 15% @ 4lp/mm		
Image depth	16 bit		1.1.63
Dyn range (RQA5)	6 uR –9 mR		
Load capacity (max)	150 kg distributed weight, 100 kg concentrated weight (45mm diam area)		
Ingress protection	IPX4 (splashing water)		1.1.61

The integrated Quality Assurance Procedure (QAP) enables easy image quality checks by the customer. A phantom, optimized for Digital Image Quality testing, is included with the system. Changes are highlighted and can be corrected before they become a problem.

1.1.58 **Quick Connect** adaptive wireless technology enables automatic channel switching to improve image transfer and avoid wireless interference with other equipment on the hospital network

Quick Connect Communication	
Wireless interface	802.11n, 5Ghz
Ethernet interface	1000 Mbps
Tether available	yes

QuickShare allows registered detectors to work across multiple compatible GE systems with no additional configuration required. Pairing enables registered detectors to connect wirelessly within seconds.

1.1.65 **Quick Charge** allows in-bin charging whereby the FlashPad HD detectors charge when they are in the table or wall stand housing. This gives you peace of mind knowing that the detectors will be ready when you need them.

Quick Charge	
Battery time to full charge	2.5 hrs
Battery indicator and charger	yes

Helix™ 2.0 Advanced Image Processing

Helix 2.0 advanced image processing algorithms harness the full high-resolution power of FlashPad HD detectors to deliver exceptional image quality despite challenging exams conditions.

Helix 2.0 algorithms are designed to deliver outstanding resolution, excellent edge presentation, consistency, and noise handling. The algorithms incorporate the following capabilities:

Helix 2.0 Advanced Image processing	
Auto Shuttering	Automatically detects collimator edges and adjusts to the selected field of view: • ACED (Automatic Collimator Edge Detection) An algorithm that provides automated collimator masking of the image. • ICED (Intelligent Collimator Edge Detection) An intelligent algorithm that relies solely on image information to locate collimation edges present in an x-ray image.
Raw Radiation Rejection	Identify raw radiation pixels and improve post processing image display.
1.1.72 Grid Line Reduction	Suppresses grid lines on the image without impacting anatomical details.
Detail Preserving Noise Reduction	Suppresses the mottle noise in denser areas of the anatomy while preserving detail in the rest of the image. The algorithm takes into account tissue penetration and dose reaching the receptor while preserving the detailed structures.
Localized Contrast Enhancement	Improves local contrast in the lung and other areas without affecting the overall look of the image.
T2 Auto-grid	Software that can be used in lieu of a physical anti-scatter grid to improve image contrast in general radiographic images by reducing the effects of scatter radiation.
Multi-resolution processing	Improved edge presentation, exquisite bone detail, differentiation of soft tissues and visualization of metal.
Tissue Equalization	Enhance contrast in thick and thin regions of the anatomy to correct over penetrated and under-penetrated areas within the image.
AI Smart Windowing	Delivers the correct display brightness and contrast without needing window level and width adjustments, Provides consistent brightness and contrast across variations in exposure technique. The smart windowing uses an artificial intelligence algorithm to determine the correct settings.
Multiple customizable looks	Four pre-set image processing selections “looks” are optimized for each anatomical view with the ability to define multiple custom looks for each anatomical view/patient size combination.

Workflow Features

The Worklist is the starting point for patient set up and selecting procedures for acquisition. The Worklist provides an automated method of obtaining exam and protocol information for a patient directly from a DICOM worklist server.

Worklist Key Features

- Column-customized worklist
- Auto-refresh
- Emergency patient imaging
- Patient edit/auto-folders
- Barcode reader option for patient data entry and selection

Auto-Positioning protocol link enables pre-determined receptor and corresponding auto-positioning protocol per view.

Auto-Positioning Key Features

- Seven predefined positions
- Ten user-customizable positions
- OTS tube angle, longitudinal, lateral, rotations and vertical positioning
- Table receptor longitudinal positioning
- Wall stand vertical positioning and tilting

1.1.6, 1.1.70 The Acquisition screen is where the exam is set up and exposure details are adjusted. Set of default adult and pediatric protocols allows quick selection of the appropriate techniques for common procedures/exams with the ability to define unlimited number of custom protocols.

Acquisition Key Features

- Column-customized worklist
- Auto-refresh
- Emergency patient imaging
- Patient edit/auto-folding
- Barcode reader option for patient data entry and selection

Image annotations includes both system and user annotations.

Annotations Key Features

System data with configurable font size and on/off control:

- Patient identification
- Exam date & time 1.1.74
- Exam identification
- Image identification
- Dose and DAP values
- Hospital/Institution name
- X-Ray parameters: kVp, mA, ms, mAs, and DEI values
- RRA classification & reason
- Anatomy information
- Processing information
- User measurements: Size and angle for line, ellipse and Cobb
- Image size and zoom
- Operator entered annotations

User image annotations:

- Preset and custom annotations
- Line, ellipse and Cobb tools
- L/R markers

Image management provides fast access to the image and exam database for case reviews and file management.

Image Management Key Features

- Copy patient images into a second patient entry patient entry
- De-identify makes an anonymous copy of the exam removing all patient identifying information
- PACS Test Images to test the quality of images sent to PACS
- TG18 test images that are used to calibrate the display monitor.
- Image exchange via Read/Write (write once, multiple access) CD/DVD-ROM
- DICOM viewer provided on the archive CD/DVD

The Quick Toolbar provides quick access to commonly used functions while reviewing images. Adapt the Quick Toolbar to include your most frequent operations.

Quick Toolbar Functions	
Pan image	Moves the image within the viewing area
Image magnifying glass	3x magnification of a small part of the image
Brightness / Contrast	Adjust the brightness and contrast via the pointer
User annotations	Individual tools to add and erase custom annotations
Add RL marker	Places a Right or Left marker on the image
Manual shutter auto-reprocess	Adjust the shutter and reprocess the image
Quick Enhance	Applies another predefined image processing look
Patient orientation	Change the patient orientation of an image.
Rotation/Flip	Rotation and flipping of the image
Add manual shutter	Adjust the automatically applied shutter
Line/Ellipse/Cobb	Measurement tools of the size and angle of features
Invert	Reverses light and dark areas of the image
Film manager	Printing of multiple images in a series
Manual print	Printing of the currently selected image
Restore image	Returns the image to its original state

The toolbar includes **Quick Enhance**, a one-touch function that can reprocess images with a different custom look with no additional dose to the patient and no additional clicks for the user. Customize Quick Enhance by anatomy for multiple uses including instrument check, implant visualization, line placement.

The Image Display Tools contains the controls to flip, rotate, adjust brightness, adjust contrast, invert, and apply windowing to the selected image.

Image Display & Preview Key Features

- Window width and level
- Gray scale invert
- Interpolated zoom with roam
- Image flips (horizontal, vertical) with automatic indicator
- Image Rotate with 90° increments
- Free rotation of full 360°
- Image orientation management
- L/R markers
- Free text annotation
- Manual shuttering

Optional Workflow Features

Auto Protocol Assist[†]

System will automatically transition directly to the Acquire screen when the protocol code downloaded from the HIS/ RIS (automatically performed with worklist refresh) matches the exam code contained in the protocol database. This tool eliminates the user steps required to select patient exam types and initiate an exam.

Repeat/Reject Analysis[†]

An automated quality assurance tool that allows for repeat or reject images to be captured and categorized by technologist

Reports can be exported in DVD, CD or USB media for ease of use. Discovery XR656 HD is also compatible with GE's Xray Quality Application featuring Repeat Reject Analytics.

X-ray Quality Application is a web-based solution that can connect to multiple compatible radiography systems and help you identify root causes of rejects, enhance training, drive efficiency and help reduce dose.

Auto Field of View[†]

Auto Field of View enables the user to pre-define the collimation size on an individual view basis.

1.1.30 Live Streaming Camera[†]

A camera situated on the tube provides live video streaming of the patient imaging area onto the operator console to monitor the patient status, movement and orientation before an X-ray is taken.

Video Stream Specifications

Image resolution	1140 x 1140
Frame rate	30 fps

OTS User Interface

The OTS User Interface allows you to make receptor, FOV, kV and mAs selections without returning to the acquisition workstation. This interface also provides the functions to move the OTS in multiple directions.

1.1.27

OTS User Interface	
Display type	LCD touch interface
Orientation	Automated horizontal and vertical display for table and wall stand applications
Technique adjust	kVp, mAs
Receptor selection	Table, wall stand, table-top or cassette
Patient information	Patient name and date of birth for in-room verification
Collimator control	Field-of-view selection
Information display	Exam inhibit notification area
Position display	SID, Tube Angle, Column Rotation, Detector Orientation
Position control	Lock and detent control
Auto Field of View (FOV) default selections	43 cm x 43 cm (17 in x 17 in) 35 cm x 43 cm (14 in x 17 in) 43 cm x 35 cm (17 in x 14 in) 24 cm x 30 cm (9 in x 12 in) 30 cm x 24 cm (12 in x 9 in) 18 cm x 24 cm (7 in x 9 in) 24 cm x 18 cm (9 in x 7 in)

1.1.67 Acquisition Workstation

The Acquisition Workstation is the primary interface to the network and provides image post-processing capabilities. The System Controller Module provides single point control, directing and coordinating overall system operation, while monitoring all system modules automatically through software.

Workstation Specifications (nominal)

1.1.77 Monitor	24 inch (61 cm) LCD Color Monitors (1920 x 1200 pixels) - Non-touch ≥250 cd/m2 calibrated brightness
CPU	Intel Xeon processor, 3.5 GHz, 4 cores
Hard disk storage	1 TB Enterprise-Class SATA
Image storage	≥ 17000 images Programmable auto delete function
RAM	32 GB
GPU	NVIDIA
Image processing times	Fast preview images: ≤1 second (docked) or ≤2 seconds (wireless or Tether) Final Conditioned Image including Auto-Shuttering: 1.1.67, 1.1.69 ≤6 seconds (docked) or ≤8.5 seconds (wireless or tether) Volume Rad reconstruction: 60 seconds
Expose to expose cycle time	≤5 seconds @ 70% HU
Time to boot the system after normal shutdown	≤180 seconds
System reset time	≤290 seconds

Stand-Alone Helix™ Workstation†

Supplemental workstation allows an additional technologist or radiologist to edit images and reprocess without system impact, keeping the system available for patient exams. Provides a central location to perform occasional quality control by physicists and technologists for exams needing reprocessing.

- Accepts images from any FlashPad HD-based system, including both Fixed and Mobile systems.
- Push any GE HD image to the Helix workstation, perform reprocessing, and then push to PACS
- Powerful GPU processor enables fast reprocessing and reconstruction
- Identical software and hardware specifications as the Acquisition Workstation

DICOM Connectivity

Images may be transmitted manually or automatically through the DICOM interface to printers, archival devices, servers, or review workstations. System Access and Authorization Control support HIPAA Compliance.

Refer to the DICOM Conformance Statement for complete definition of supported DICOM connectivity services.

DICOM 3.0 Services

1.1.75	Modality Worklist (SCU)	Interface with HIS/RIS with programmable auto refresh
	MPPS (SCU)	Feedback the status of exams to the HIS/RIS
1.1.75	Storage (SCU)	Manual and auto send image (DX or CR IOD) to multiple PACS
1.1.75	Storage commitment (SCU)	Send commitment state
	Query/Retrieve (SCU)	Query/Retrieve images from PACS
	Query/Retrieve (SCP)	Provide Query/Retrieve service instance to other system
	Media exchange	CD/DVD DICOM image export and import.
1.1.75	Grayscale print	Manual and Auto print with print layout options at the console
	Verification services	C-Echo as SCU and SCP
1.1.75	Dose Structure Report (optional)	Send Dose values for each study to an archiving system.
	Audit trail	Collect usage logs for security auditing with DICOM conformance format.

IHE Integration Profiles

Scheduled workflow	Acquisition Modality: Patient Based Worklist Query/Broad Worklist Query Assisted Acquisition Protocol Setting
Patient Information Reconciliation	Acquisition Modality
REM (optional)	Radiation Exposure Monitoring for Dose structure report
TLS [†]	Transmission security

Several popular printers have been validated for connectivity and image quality. Recommendations are available from your sales representative. Minimum printer requirements: 10 and 12-bit printers.

Printed images are not intended for diagnostic use unless produced with a printer capable of at least 1,000 gradations of gray scale (or at least 10 bits)

Non-DICOM laser cameras will require an upgrade to DICOM connectivity.

Printing Features

- Orthopedic magnification/Print
- Multi-format printing: 1x1, 2x1, 1x2, 2x2
- Print preview function
- Auto portrait and landscape print

Network Connectivity, Remote Service and Security

IPv4 and IPv6 are supported and configurable.

GE Healthnet Services can provide physical network connectivity solutions – Layer 1 and 2 Ethernet (IEEE 802.3) interoperability – and include network components and physical installation.

InSite Remote Service – IIP (Integrated Insite Platform) remote system tool that supports remote communications between customer and GE. Problems may be diagnosed, resolution expedited, or problem fixed remotely with IIP without the need for a Field Service Engineer be onsite

IT Security Pack[†]

IT Security Pack Features	
Anti-Virus and Intrusion Detection	Enabled. Anti-virus software monitors files (read or execute) and report errors to user.
Anti-Virus scanning	Real-time or on-demand
Virus Signature Update	Updates via local USB or site ePO server
FIPS Compliance	Complies with FIPS 140-2 validated tools and libraries for data encryption or hashing
TLS	Encrypted data transition using Transport Layer Security

1.1.17

Overhead Tube Suspension

The Overhead Tube Suspension (OTS) system with motorized movement delivers excellent levels of operational support designed for efficient operation and precise positioning.

Overhead Tube Suspension – Features

Motion control	Manual and Auto-positioning
Auto-positioning	5 axis servo motion: longitudinal, lateral, vertical, tube angulation and column rotation
Auto tracking	Vertical tracking to align with wall stand detector Vertical tracking to maintain table SID
Override mode	Two-position key switch toggles between automatic and override manual control for complex or emergency positioning
Auto-detents	Assisting the user with locating and securing detents
Cable management	Cable chain or cable drapes depending on room layout
Remote control	Infrared remote controller (optional) for auto-positioning and FOV selection

1.1.25,
1.1.26

1.1.29

Auto-Positioning enables the user to select a predefined system position and automatically move the equipment by simply holding the “Auto-Positioning” buttons at the acquisition console or with the remote control. This feature is designed to help reduce user fatigue and increase the productivity of the operator.

Auto-Positioning will incorporate angulation of the tube, longitudinal, lateral, rotational and vertical positioning of OTS, table detector longitudinal positioning, wall stand detector vertical and tilting positioning

Auto-Positioning comes with more than 9 default positions at the table, wall stand and park position at various SIDs and vertical and horizontal orientations. The user can create additional defined positions as needed.

Simultaneous movement of the OTS might be limited in some situations due to safety considerations during auto-positioning

The OTS has available movement on 5 axes for use in advanced applications. The suspension provides convenient movement and accurate positioning of the equipment.

Overhead Tube Suspension - Specifications

Longitudinal travel range	2.3 m (for a 3.4 m rail) 4.7m (for a 5.8 m rail)	1.1.18
Lateral travel range	1.0 m (for 2 m bridge) 2.0 m (for 3 m bridge) 3.3 m (for 4 m bridge)	1.1.19
Vertical travel range	180 cm. Dual cable safety system and high precision telescoping column	1.1.20
Tube angulation range	+135/-180 degrees. Continuous rotation with locking at any position. Fixed detents at 0, +90, and -90 when selected	1.1.22
Column rotation range	+/-135 degrees. Continuous rotation with locking at any position	1.1.21
Longitudinal Speed	Up to 30 cm/s	
Lateral speed	Up to 25 cm/s	1.1.18, 1.1.19
Vertical speed	Up to 25 cm/s	1.1.20, 1.1.21 1.1.22
Tube angulation speed	Up to 25 degree/s	
Column rotation speed	Up to 25 degree/s	

Elevating Table

A bariatric X-ray table enables you to serve the X-ray imaging needs of patients of all sizes and mobility levels.

Elevating Table - Main Features		
Motion control	Motor driven elevation	1.1.41
	Manual 8-way floating tabletop	1.1.45
	Servo motion for detector longitudinal	
1.1.26 Auto tracking	Automated detector longitudinal tracking to align with tube (including tube longitudinal and angulation)	
Maximum patient weight	400 kg (882 lbs) static non-elevating	
	320 kg (705 lbs.)	1.1.42
	dynamic elevating with load centered	
Patient safety	Double-press safety foot pedal Two safety switches to disable motion during patient transfer Two emergency stop buttons Table anti-collision safety devices Optical anti-finger pinch sensors	
Detector housing	Supports FlashPad HD 3543 and 4343 with portrait and landscape rotation for FlashPad HD 3543	1.1.56 1.1.44
AEC support	3 cell ion chamber	1.1.5
Foot pedal	Front foot pedal (default) Rear foot pedal (optional)	
Available accessories	Tabletop detector holder Patient hand grips Compression bands	

Elevating Table - Specifications		
Tabletop material	Carbon fiber composite	
Tabletop inherent filtration	≤0.7 mm Al equivalent @ 100 kVp	
Tabletop size	93 cm x 240 cm (37 in x 94 in)	1.1.39, 1.1.40
Tabletop longitudinal travel	68 cm	
Tabletop transversal travel	28 cm	
Detector longitudinal travel range	80 cm	
Max. patient coverage	Long: 188 cm (74 in) Lat: 68 cm (26.8 in)	
Elevating range	50 cm – 85 cm (19.7 in – 33.5 in)	
Elevation time from minimum to maximum height	≤18 seconds	
Foot pedal ingress protection	Liquid proof IP36	
Footprint (incl foot pedals)	100 cm x 71 cm (39.4 in x 28 in)	

Available Table Grids†	
100 cm table grid	Focal range: 90 – 120 cm Vertical orientation Aspect ratio: 12:1 Line density: 70 lp/cm
110 cm table grid	Focal range: 95 – 130 cm Vertical orientation Aspect ratio: 12:1 Line density: 70 lp/cm
120 cm table grid	Focal range: 102 – 146 cm Vertical orientation Aspect ratio: 13:1 Line density: 70 lp/cm

Digital Tilting-Arm Wall Stand

The digital tilting-arm wall stand is designed for use with the wireless digital detector, ion chamber and **removable nonreciprocating grid**. 1.1.49

Tilting Wall Stand - Features

Detector housing motion control	Motorized vertical and tilting movement.	1.1.47, 1.1.52
Motion safety	Manual vertical and tilting movement compatible. Vertical electromagnetic braking and tilting anti-collision safety device	
Auto tracking	OTS auto-tracking of the detector	
Foot switch	Vertical motion control	
Remote control	Infrared remote controller (optional)	
Information display panel	Indicates the exposure hold, grid type, tilting angle, and system ready status	
Detector housing	Support FlashPad HD 3543 and 4343 , with portrait and landscape rotation for FlashPad HD 3543	1.1.56
Extended Arm Option	34 cm (13.4 in) longer than standard wallstand to accommodate stretcher work	
AEC support	Std arm: 3 cell ion chamber Ext arm: 4 cell ion chamber The standard right, left and center ion chambers are available when the detector is vertical. The 4th ion chamber is for when the detector is horizontal to support supine chest and abdominal imaging when used with a stretcher	1.1.5
Accessories	Integrated hand grips and lateral support bar for patient comfort and stabilization	1.1.50, 1.1.51

Tilting Wall Stand - Specifications

Floor to center of detector (min)	28.5 cm (11.2 in)	1.1.48
Vertical travel range	150 cm (59 in)	
Patient coverage (max)	192 cm (75.6 in)	
Tilt range	-20° to 90°	1.1.52
Tilt detents	-20°, 0°, 90°	

Available Wall Stand Grids†

100 cm wall stand grid	Focal range: 90 – 118 cm Vertical orientation Aspect ratio: 13:1 Line density: 70 lp/cm
120 cm wall stand grid	Focal range: 102-146 cm Vertical orientation Aspect ratio: 13:1 Line density: 70 lp/cm
130 cm wall stand grid	Focal range: 90 – 190 cm Vertical or Horizontal Aspect ratio: 10:1 Line density: 70 lp/cm
180 cm wall stand grid	Focal range: 145 – 245 cm Vertical orientation Aspect ratio: 13:1 Line density: 70 lp/cm

Generator

The digital radiographic imaging system is available with a high frequency generator.

X-ray Generator	
Nominal Output	50 kW, 65kW or 80 kW 1.1.2
Tube voltage range	40 to 150 kV 1.1.3
Tube current range	50kW: 10 to 630 mA 65 kW: 10 to 800 mA 80kW: 10 to 1000 mA For small focal spot: 10-400mA
Exposure time range	2 to 2000 ms
Current time product range	0.25 to 630mAs 1.1.4 For Tube large focal spot: 0.63-630 mAs For Tube small focus Spot: 0.25-500 mAs
AEC max. backup	512mAs and/or 2000ms
AEC Nominal Irradiation Shortest Time (NIST)	2ms

X-ray Tube

The Maxiray 100 is a high speed, 4-inch (100mm) diameter rotating anode tube unit 1.1.9 and inserts for high-energy radiographic procedures. The shockproof housing is constructed of aluminum and lined with lead to minimize leakage radiation.

X-ray Tube	
Focal spot	0.6 /1.3 (IEC60366, nom) 1.1.10, 1.1.11
Target angle	12.5°
Anode rotation speed	10,000 rpm (150 Hz to 180 Hz)
Anode power (nom)	Small focal 32 kW 1.1.12, 1.1.13 Large focal 96kW
Max. exposure voltage	150 kV 1.1.14
Inherent filtration @ 150 kVp / 75 kVp	Insert: 0.8 mm Al equiv Housing: 0.3 mm Al equiv
Leakage radiation	≤50mR/h @150 kV, 4 mA
Heat storage capacity	Anode: 260 kJ (350 kHU) 1.1.15 Housing: 1,110 kJ (1,500 kHU)
Maximum heat dissipation rate	Anode: 75 kHU/min (925 watts) Tube Unit: 60 kHU/min (740 watts) with blower operating.
Cooling and protection	The tube unit cooling consists of fan, thermal switch & pressure switch
Weight	29.5kg

Automatic Collimator

The multi-leaf collimator allows you to adjust the radiation field size to the anatomy. The collimator can be used in either the manual or automatic mode.

Automatic Collimator		
Inherent filtration @ 70 kVp	2.0 mm Al (collimator only) 2.7 mm Al (collimator and tube)	
Light field	LED lamp	
Light on-time	Automatically controlled light ON time (configurable 0 to 90 sec)	
Light brightness	≥200 Lux	
Centering indication	Shadow crosshair and longitudinal centering laser line	
Copper filters	None, 0.1 mm, 0.2 mm, 0.3 mm	1.1.36, 1.1.37
Filtration selection	Manual or automatic through protocol set up	
Rotation	≥ +/-90°	1.1.35
Collimation control	Manual and motorized, preset through protocol	1.1.34
Auto FOV response time	≤1 sec for complete closure/opening of rectangular blades	
Display	Digital readout of FOV at SID	
SID measurement tape	Available for convenient and precise SID measurement of the tabletop exam	

Dose Monitoring and DAP

1.1.53
1.1.54

Detector Exposure Indicator (DEI) is a tool for tracking patient over/under-exposure by estimating radiation exposure behind the patient and is a relative measure of exposure to the detector.

Exposure Index (EI) is proportional to detector exposure assuming that the x-ray technique used is the same as that of the calibration technique.

Deviation Index (DI) estimates the deviation of actual detector exposure from target detector exposure.

Dose Area Product (DAP) is automatically annotated onto the digital image for the exposure and is displayed on the acquisition screen post exposure.

Because the DAP is a calculation method, any additional filtration used in the beam, aside from that provided with the system, will introduce an error in the reported dose. It is recommended that additional filtration not be used when dose reporting is enabled.

Entrance Dose (unit: mGy) is an estimate of entrance dose (air-kerma) at a distance in front of the wallstand cover or above the tabletop, depending on which receptor was used for acquisition.

The **Dose Reporting Tool** allows a user to export dose relevant data within a specific time frame. The dose data can also be included as an annotation on the image.

Accessories†

Compression band		Hand grips	
1.1.43 Lateral detector holder		Mobile detector Holder	
Desktop battery charger		Detector charging bin	
Clip-on grid (10x12, 14x17)		Detector Portable Handle (with or without grid for 14x17)	
Foot stool for image pasting barrier		UPS 1.1.78	
Stretchers		1.1.28 Patient barrier for image pasting	

Installation Power

Primary source is required for all installations. Demand includes power for the entire digital radiographic imaging system.

Primary Source Input	
Input voltages	Three phases with/without neutral. 380, 400, 420, 440, 460 and 480 VAC $\pm$ 10%.
Input frequency	50 Hz or 60 Hz
Input current	170A (momentary) 4.5A (continuous)
Input power	112kVA (momentary), 2.2kVA (continuous)
Energy consumption:	Current input: 2A, power input 1.0 kVA when system in standby.

Environmental Conditions


Environmental Conditions	Operating	Non-Operating
Altitude	-30 m to +3,000 m relative to sea level	-30m to 3,000m relative to sea level and support non-pressurized air transport
Temperature	15°C to 32°C	-20°C to 60°C (excluding detector) -5°C to 50°C (Detector)
Humidity	20% to 75% RH, non-condensing	10% to 85% RH, non-condensing
Atmosphere Pressure	106 kPa to 70 kPa	106 kPa to 70 kPa
Audible Noise (1 meter from system)	$\leq$ 60 dBA during motion	NA

Typical Room Layout

Recommended room height: 2.9m (114")

For minimum room requirements refer to the Pre-Installation Manual.

1	System Cabinet	6	Tether Interface Box †
2	Longitudinal Station Rails	7	Access Point
3	Detector + Table	8	Magic PC + Monitors
4	OTS and Bridge	9	Grid Holder
5	Detector + Wallstand	10	Cable Rail /Bracket

Compliance to Standards

The Discovery XR656 HD digital radiographic imaging system is designed to meet applicable performance standards for diagnostic X-ray equipment enunciated by the U.S. Department of Health and Human Services pursuant to the Radiation Control for Health and Safety Act. In addition, the system complies with UL, IEC requirements.

Warranty

The published company warranty in effect on date of shipment shall apply. Right reserved to make changes.

†Denotes optional features

About GE Healthcare

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter – great people and technologies taking on tough challenges. From medical imaging, software and IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

GE Healthcare, Europe

Headquarters Buc, France
+33 800 90 87 19

GE Healthcare, Middle East and Africa

Istanbul, Turkey
+90 212 36 62 900

GE Healthcare, North America

Milwaukee, USA
+1 866 281 7545

GE Healthcare, Latin America

Sao Paulo, Brazil
+55 800 122 345

GE Healthcare, Asia Pacific

Tokyo, Japan
+81 42 585 5111

GE Healthcare, ASEAN

Singapore +65 6291 8528

GE Healthcare, China

Beijing, China +86 800 810 8188

GE Healthcare, India

Bangalore, India
+91 800 209 9003

General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

©2019 General Electric Company - All rights reserved.

GE, GE monogram and imagination at work are trademarks of General Electric Company. www.gehealthcare.com

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Edit View 1 of 2

Make desired changes and press Next.

View Name : ABDOMEN_AP

Position View After : Position First

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ABDOMEN / ABD_SERIES / CHEST_PA Standard Acquisition			kV	mA	mAs (Fixed only)	msec (p = Back-up Time)	Speed	SLD	Ion Chamber (AEC only)	Focal Spot	Cu Filter	Grid	Asymm. Collim.	Auto Position
 Small Adult		Fixed	110	200	1	5.0	400	180	N/A	Small	no	Out	Off	Park
 Med Adult		Fixed	120	320	2	6.25	400	180	N/A	Large	no	Out	Off	Park
 Large Adult		Fixed	120	500	3.2	6.4	400	180	N/A	Large	no	Out	Off	Park
 Small Ped		Fixed	65	100	0.8	8.0	800	180	N/A	Small	no	Out	Off	Park
 Med Ped		Fixed	70	100	0.8	8.0	640	180	N/A	Small	no	Out	Off	Park
 Large Ped		Fixed	75	100	1	10.0	500	180	N/A	Small	no	Out	Off	Park
 Small Adult		AEC	110	200	N/A	5.0	400	180	LR	Small	no	180	Off	WS-V...
*  Med Adult		AEC	120	320	N/A	6.25	400	180	LR	Large	no	180	Off	WS-V...
 Large Adult		AEC	120	500	N/A	6.4	400	180	LR	Large	no	180	Off	WS-V...
 Large Ped		AEC	75	100	N/A	10.0	500	180	LR	Small	no	180	Off	WS-V...

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NEW

Available as ...

Premium Line - The Full Range
Comfort Line - The Economical Alternative

Frontal protection in its pure form in the RA650 includes full-body protection over the entire front section with the required 60% of the girth.

- Suitable for a wide range of applications with only frontal radiation exposure
- Ideal fit, maximum comfort for the wearer
- Possible to slip the apron off
- Optional "without protection" in the back, for reduced weight of the apron
- Standard equipment includes shoulder pads and preparation for adaptable embroidery
- Type testing and certification in accordance with (EU) 2016/425, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016

1.1.79

Sizes*, recommended for:

XS	Chest girth ≤ 84 cm	Waist/hip girth ≤ 91 cm
S	Chest girth approx. 85 - 95 cm	Waist/hip girth approx. 92 - 102 cm
M	Chest girth approx. 96 - 106 cm	Waist/hip girth approx. 103 - 112 cm
L	Chest girth approx. 107 - 117 cm	Waist/hip girth approx. 113 - 123 cm
XL	Chest girth approx. 118 - 128 cm	Waist/hip girth approx. 124 - 134 cm
XXL	Chest girth approx. 129 - 139 cm	Waist/hip girth approx. 135 - 145 cm
3XL	Chest girth approx. 140 - 150 cm	Waist/hip girth approx. 146 - 156 cm

* The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.

FRONT PROTECTION APRON

ERGONOMIC RA650

Fast, simple and comfortable

The apron can be slid off by simply opening the hook-and-loop fasteners in the front area, even through the sterile surgical gown.

Crossbands

Comfortable, wide crossbands allow the apron to be slid off easily by opening the hook-and-loop fasteners in the front area, through the sterile surgical gown as well.

Closure system

MAVIG uses only top quality, low-wear hook-and-loop fasteners.

Apsauga/švino ekvivalentas

Bešvinė

Protection/lead equivalent


NovaLite, **Leadfree** or Standard lead

Front, optional:
0.25 mmPb, 0.35 mmPb
or **0.50 mmPb**
Back, optional:
0.25 mmPb
or without protection

Outer material

Colours/designs (see front page)

Optional:
ComforTex® HPMF,
ComforTex® HPMF Hybrid,
ComforTex® HPMF Plus
(Nautilus design only)

Sizes

from XS to 3XL
(see front page)

Skirt length / Body height:
90 cm / 150 - 160 cm
100 cm / 160 - 170 cm
110 cm / 170 - 180 cm
120 cm / 175 - 190 cm
130 cm / 185 - 195 cm
140 cm / 195 - 205 cm

Optional

On request:
Sleeves
Breast pocket
Embroidery
Made-to-measure protection

Accessories/replacement

Bolero
Thyroid/sternum protection

RA614

Thyroid Gland Protection

1.1.79

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NEW

Available as ...
Premium Line
Comfort Line

Neck widths, depending on model:

Hook-and-loop fastener (RA614KL)

Standard size fits neck sizes of 32 - 45 cm.

Magnetic closure (RA614MA)

Standard size fits neck sizes of 32 - 45 cm.

Interchangeable hook-and-loop fastener (RA614AK)

Standard size fits neck sizes of 37 - 50 cm.

Titan

Blue Admiral

Regatta

Curacao

Lagoon

Mandarin

Orchidee

Baccara

"Nautilus"

(Comfort® - spine plus only)

THYROID AND STERNUM PROTECTION

RA614

The thyroid is a radiation-sensitive organ of the human body and therefore in particular need of protection.

Wearing an appropriate shield should therefore be a matter of course.

Our thyroid and sternum protector variants are characterised by a first-class cut and the best fit, taking into account all aspects of radiation protection.

The press-stud connection allows secure adjustment of the thyroid and sternum protector on MAVIG X-ray protective clothing.

Model Variants RA614

The thyroid and sternum protector with press-stud connection to the vest/coat/apron is suitable for all model series.

As standard, it is manufactured with a firmly sewn-on hook-and-loop fastener. The other closure types, magnetic closure and interchangeable hook-and-loop fastener, are available as options.

RA614KL Thyroid/Sternum Shield with hook & burr closure

RA614AK Thyroid/Sternum Shield similar to RA614KL with an exchangeable hook & burr closure

RA614MA Thyroid/Sternum Shield similar to RA614KL with magnetic closure ^{*1}

Snap Fasteners

for a secure attachment of the thyroid and sternum protection to MAVIG protective clothing

Closure type – via magnet

The thyroid and sternum protection is also available with a practical magnetic closure (variant RA614MA).

Accessories: Hygienic cover


To meet the hygienic requirements in everyday life, we recommend using our washable hygienic cover in the colour Titan made of the special ComforTex® MF material.

Apsauga/ švino
ekvivalentas

Bešvinė

Protection/lead equivalent

Novalite, **Leadfree**
or Standard lead

As standard:
0.50 mmPb

1.1.79

Outer material

Colours/designs
(see front page)

Optional:
ComforTex® HPMF,
ComforTex® HPMF Hybrid,
ComforTex® HPMF Plus
(Nautilus design only)

Sizes/Neck widths

Uniform
Depending on the model

Optional

(pages 30/31)

On request:
Embroidery
Made-to-measure protection

Accessories/replacement

Hygienic cover
(RA614HF)

Hook-and-loop fastener
(R614AKV,
adjustment via
push buttons)

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This EU declaration of conformity is issued under the sole responsibility of MAVIG GmbH who confirms that the products mentioned below meet the provision of the Regulation (EU) 2016/425 on personal protective equipment.

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The notified body ANCCP Certification Agency srl., Via dello Struggino, 6 – 57121 Livorno (LI) – Italy performed the EU type-examination and issued the EU type-examination certificate PPE-1203-18319-Cert.

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The PPE is subject to the scheduled periodic surveillance procedure set out in Annex VIII (Module D) of Regulation (EU) 2016/425 under the supervision of the following Notified Body:

ANCCP Certification Agency srl., Via dello Struggino, 6 – 57121 Livorno (LI) – Italy

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Munich, 2021-07-08

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Company Address:
Stahlgruberring 5
81829 Munich
Germany

Post Box:
Postfach 820362
81803 Munich
Germany

Phone:
+49 (0) 89 / 420 96 - 0

Fax:
+49 (0) 89 / 420 96 - 200

e-Mail:
info@mavig.com

Web:
www.mavig.com

Bank Accounts:
HypoVereinsbank Munich
IBAN: DE80 7002 0270
1620 2253 00
SWIFT-BIC: HYVEDEMMXXX

Postbank Munich
IBAN: DE12 7001 0080
0071 3748 07
SWIFT-BIC: PBNKDEFFXXX

Trade-Register Munich
HRB-No. 76 403

Managing Directors:
Christian Stoian
Alexia Stoian-Lepère

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Example for the formation of the article code:

RA614	5	22	KL	1	1
	①	②	③	④	⑤

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④ Außenmaterial – Hülle/ Outer Material - Cover	
	Code
ComforTex HPMF	1
ComforTex HPMF Hybrid	2
ComforTex HPMF Plus	3

⑤ Strahlenschutzmaterial – Innen/ X-Ray Protective Material - Inlay	
	Code
NovaLite 0,50 mmPb	1
Standard-Blei/ <i>Standard Lead</i>	2
Leadfree	3

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Trade-Register Munich
HRB-No. 76 403

Managing Directors:
Christian Stoian
Alexia Stoian-Lepère

VAT-ID No.:
DE 811268845

Buc, September 13, 2022

STATEMENT

1.1.80

We, **GE Medical Systems Société en Commandite Simple** with a registered seat at 283 Rue de la Minière, 78530 Buc, France, with commercial name **GE Healthcare**, established and reputable manufacturer of radiology equipment and Authorized European representative of **GE HUALUN MEDICAL SYSTEMS CO. Ltd** No1, Yong Chang North Road, Beijing Economic Technological Development Zone, BEIJING 100176 CHINA, the manufacturer of **Definium 656 HD** (*Definium 656 HD is a commercial configuration of the Discovery XR656 HD*) hereby confirm that means required for quality assurance and calibration (including integrated software, phantoms) are part of the scope of delivery of **Definium 656 HD** (*Definium 656 HD is a commercial configuration of the Discovery XR656 HD*).

Seyda KANE

Product Marketing Manager EMEA Mature Markets

On behalf and for GE Medical Systems SCS,

Buc, May 26th 2021

TECHNICAL STATEMENT

To whom it may concern,

We, **GE Medical Systems Société en Commandite Simple** with a registered seat at 283 Rue de la Minière, 78530 Buc, France, with commercial name GE Healthcare, established and reputable manufacturer of radiology equipment and Authorized European representative of **GE HUALUN MEDICAL SYSTEMS CO. Ltd** No1, Yong Chang North Road, Beijing Economic Technological Development Zone, BEIJING 100176 CHINA, the manufacturer of **Definium 656HD** (*Definium 656 HD is a commercial configuration of the Discovery XR656 HD*) hereby confirm following parameters are available in Lithuania area.

- DQE of FlashPad HD detectors is 75 % @0.05 lp/mm in accordance with IEC 62220 standard 1.1.64

On behalf and for **GE Medical Systems SCS**,

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Véronique Soltani
Tender & Offer Leader – Shared Services Europe
Date of signature : May 28, 2021

General specifications	Product Length/Depth	330 mm
	Product Height	180 mm
	Product Width	133 mm
	Product Weight	10.46 kg
	Compliances	CB CE Marked
Product specifications	Certifications	CB IEC/EN 62040-2 CE EAC IEC/EN 62040-1
	Catalog Notes	EMEA & GLOBAL
	Secondary frequency - max	60 Hz
	Wattage	1200 W
	Form factor	Tower
	Includes network card	No
	Input connection	C14
	Input voltage - min	162 V
	Output voltage - max	280 V
	Secondary frequency - min	50 Hz
	Software compatibility	Eaton UPS Companion (enables safe system shutdown, energy usage metering and configuration of UPS settings)
	Battery type	Sealed, lead-acid
	Efficiency	95
	Special features	Automatic Voltage Regulation (AVR) Saves space with a compact tower design HID-compliant USB port, for automatic integration with major operating systems (Windows/Mac OS/Linux) Data line protection to ensure that the internet line (including xDSL) is protected against surges
	Construction type	Free standing model
	Output voltage - min	170 V
	Color	Black
	Extended battery capability	
	Output voltage range	170-280 V
	Topology	Line-interactive
	Surge protection	Yes
	Runtime at half load	5 min
	Input nominal voltage	230V default (220/230/240V)
	Input frequency range	45-63 Hz
	Input voltage - max	290 V
	Voltage type	AC
	Phase (output)	1

	Battery management	Automatic battery test, deep-discharge protection, cold-start capable
	Auto shutdown function	
	Number of Outputs C13	6
	Operating temperature - min	0 °C
1.1.78	VA rating	2000 VA
	Potential free switch contact	
	Receptacle	(6) C13
	Phase (input)	
	Package contents	(2) IEC-IEC output cablesUser manual
	Noise level	<45 dB at 1 meter
	Feed type	
	Runtime graph	View runtime graph
	User interface	LEDs
	Battery quantity	2
	Operating temperature - max	40 °C
	Voltage	230V
	Temperature range	0° to 40°C (32° to 104°F)
	Ethernet interface	
	Communication	USB port (HID compliant)
	Relative humidity	0-90% non-condensing
	Output apparent power	
	Rack mounting kit	
	Output frequency	50/60 Hz
	Runtime at full load	30 sec
	Type of interface	USB
	Primary frequency - min	45 Hz
	Battery rating	12 V / 9 Ah
	Type	UPS
	Output nominal voltage	
	Expansion slots	
	Input voltage range	162-290 V
	Nominal frequency	
	Primary frequency - max	63 Hz
	Altitude	2000 m

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
DECLARATION OF CONFORMITY

Following the provisions of the medical devices directive 93/42/EEC and of the directive 2014/53/EU)

We

Manufacturer

GE HUALUN MEDICAL SYSTEMS CO. Ltd
No.1, Yong Chang North Road, Beijing Economic
Technological Development Zone
BEIJING 100176 CHINA

EU Authorized Representative

GE Medical Systems SCS
283 rue de la Minière
78530 BUC, France

Manufacturing site

GE HUALUN MEDICAL SYSTEMS CO. Ltd
No.1, Yong Chang North Road, Beijing Economic
Technological Development Zone
BEIJING 100176 CHINA

GE Medical Systems, LLC
3000 North Grandview Blvd.
Waukesha, WI. 53188 USA

Declare under our sole responsibility that the class IIb device:

Discovery XR656 HD

Stationary basic diagnostic x-ray system, digital

Ref. : 61361656\61360656 Discovery XR656 HD Digital X-Ray System or see addendum

GMDN Code: **37645** UDI-DI code: **00840682138598**

Classification rule (93/42/EC Annex IX) : 10

To which this declaration relates, is in conformity with the requirements of the medical devices directive 93/42/EEC, the requirements of the Radio Equipment and Radio Equipment Directive (**2014/53/EU**).

This conformity is based on the following elements:

- For the directive 93/42/EEC (MDD)
 - Technical Documentation/DHF Ref./ réf: **DOC2066386** (Technical File: Discovery XR656 HD)/**DOC1798359** (DHF) of the product to which this declaration relates.
 - EC certificate: approval of full quality assurance system (Annex II of the medical devices directive 93/42 EEC) delivered by G-MED (Notified Body 0459) Certificate N° **7846**.
 - Harmonized standards applied on the product to which this declaration relates

EN 60601-1: 2006/AC: 2010	EN 60601-1-2: 2015	EN 60601-1-3: 2008/AC: 2010
EN 60601-1-6: 2010	EN 60601-2-28: 2010	EN 60601-2-54: 2009
EN 62304: 2006/AC: 2008	EN 62366: 2008	
EN ISO 10993-1: 2009/AC: 2010	EN ISO 10993-5: 2009	EN ISO 10993-10:2013
EN ISO 10993-18:2009	EN ISO 15223-1: 2016	EN 1041: 2008
- For the directive 2014/53/EU Radio Equipment Directive (RED)
 - Technical Documentation Ref./ réf: **DOC2176588** of the product to which this declaration relates

EN 62479: 2010	EN 301 489-1 V2.2.0	EN 301 489-33 V2.2.0
EN 301 489-17 V3.2.0	EN 301 893 V2.1.1	EN 302 065-1 V2.1.1
EN 50566: 2013+AC:2014	EN 62311: 2008	EN 62209-2: 2010
EN 60601-1: 2006/AC:2010	EN 60950-1:2006 +A11:2009+A1:2010+A12:2011+A2:2013	

Beijing, China 3 January 2019

Kenny Ma
Regulatory Affairs Manager.

ADDENDUM TO THE DECLARATION OF CONFORMITY DOC2066380
<Discovery XR656 HD> – Accessories and Components

Catalog Number (GEHL make)	Sub-Systems and Major Components
S2100BS	XR656 HD BASE SYS CN
S2100OTS	XR656 HD OTS CN
S39252JB	INBOARD BRIDGE 2/3/4 METE
S39252JD	CABLE DRAPE
S39252JF	D8000 DRIVE BELT SELECT
S39252JE	INBOARD STATIONARY RAIL
S2100SWS	G3 STD WALLSTAND CN
S2100EWS	G3 EXT WALLSTAND CN
S2000CSWS	WS Cable Select
S2100TBL	G3 Table CN
S2000CSTBL	Table Cable Select
S2000APUS	AP - US
S2000APRW	AP - RW
S2000APJP	AP - JP
S3000DO	Grid 8:1 10in x 12in
S3000DQ	Grid 8:1 14x17in
S3000DP	Grid 6:1 10in x 12in
S3000DR	Grid 6:1 14x17in
S2000AC	G3 100CM TABLE GRID
S2001AC	G3 110CM TABLE GRID
S2002AC	G3 120CM TABLE GRID
S2003AC	G3 100CM WS GRID
S2004AC	G3 WS 120cm Grid
S2005AC	G3 130CM FOCUS WS GRID
S2006AC	G3 130CM CROSS TABLE GRID
S2007AC	G3 180CM WS GRID
S2009AC	G3 GRID HOLDER
S1200RZ	BAR CODE READER OPTION
S1202RZ	BAR CODE READER W STOOL
S3921KA	STRETCHER TABLE
S39212UP	UPS
S39222FP	TABLE REAR FOOT PEDAL
S39182JW	Compression Band
S39222RR	Hand Grips
S2000BIN	G3 DETECTOR BIN
S3000DL	Detector Battery 10x12in
S3000DN	Detector Battery 14x17in
S3000DT	Battery Charger 14x17in
S3000DS	Battery Charger 10x12in
S2000TTDH	G3 TBL DET HOLDER
S2000TIBCS	G3 SHORT TETHER

E06731BL	MOBILE LATERAL DETECTOR HOLDER V2
S39213TP	267mm Footstool for Patient Barrier
S39212TP	410mm Footstool for Patient Barrier
S2000WA	e-IFU Tech Pub Collector
S1200WA	FRENCH KEYBOARD
S1200WB	DEFINIUM 8000 GERMAN KEYBOARD
S1200WC	DEFINIUM 8000 SPANISH KEYBOARD
S1200WD	DEFINIUM 8000 ITALIAN KEYBOARD
S1200WE	DEFINIUM 8000 PORTUGUESE KEYBOARD
S1200WF	DEFINIUM 8000 SWEDISH KEYBOARD
S1200WK	DEFINIUM 8000 DANISH KEYBOARD
S1200WL	DEFINIUM 8000 DUTCH KEYBOARD
S1200WM	DEFINIUM 8000 NORWEGIAN KEYBOARD
S1202WR	Russia Keyboard
S39212MF	TABLE AUTO IP PERMANENT
S39212SM	WS- TABLE IP PERMANENT
S39212MK	50KW CONFIGURATION PERMAN
S39212MH	65KW CONFIGURATION PERMAN
S39212MJ	80KW CONFIGURATION PERMAN
S39212MT	RRA PERMANENT
S39212MW	AUTO FOV PERMANENT
S39212MY	APA PERMANENT
S39212NF	PATIENT ORI-TAG PERMANENT
S39212NN	RRA FLEXTRIAL
S39212NP	AUTO FOV FLEXTRIAL
S39212NR	APA FLEXTRIAL
S1201LX	Operation Console RCIM2
S39212SJ	GLOBAL D8K GLOBAL LABEL
S39212NG	TABLE AUTO IP FLEXTRIAL
S39212RS	INFRA RED REMOTE
S1204PLA	WS AUTO IP SW
S1204PLB	WS PASTING BARRIER
S2000TIBCL	G3 LONG TETHER
S2000TIBCM	G3 MEDIUM TETHER
S3000DK	Flashpad HD, 10x12
S3000DM	Flashpad HD, 14x17
S3000DV	Flashpad HD, 10x12, PKI
S3000DW	Flashpad HD, 14x17, PKI
S3001DK	Int Flashpad HD 10x12
S3001DM	Int Flashpad HD 14x17
S3001DV	Int Flashpad HD 10x12 Vrx
S3001DW	Int Flashpad HD 14x17 Vrx
S3001DX	2x Int FlashPadHD 1417 Vx
S3001DY	2x Int FlashPad HD 14x17

S2021CFR	21CFR Label for US
E4502ST	25 KAIC X-Ray Main Disconnect Panel 80 Amp, 480 V 208 V
E1500RD	Adventure Series - Bear Brook for XR650 XR656
E1500RR	Adventure Series - Beachside Adventure for XR650 XR656
E1500RL	Adventure Series - Playful Park Adventure for XR650 XR656
E1500CT	Adventure Series - Deep Blue Sea XR Room
S3000GK	GE Grip Flashpad HD 10x12
S3000GM	GE Grip Flashpad HD 14x17
S21DDBS	XR656 HD BS DoD CN
S3003DD	IT Security Pack Option
S3003DL	Dose Structured Rprt Opt
Catalog Number (WSO make)	Sub-Systems and Major Components
S3100BS	XR656 HD BASE SYS US
S3100OTS	XR656 HD OTS US
S2100ME	INBOARD BRIDGE 2/3/4 METER
S2100MG	CABLE DRAPE
S2100LL	D8000 LONG DRIVE BELT SELECT 2/3M OR 4M SELECT
S2100JC	INBOARD STATIONARY RAIL
S3100SWS	G3 STD WALLSTAND US
S3100EWS	G3 EXT WALLSTAND US
S2000CSWS	WS Cable Select
S3100TBL	G3 TABLE US
S2000CSTBL	Table Cable Select
S2000APUS	AP - US
S2000APRW	AP - RW
S2000APJP	AP - JP
S3000DK	Flashpad HD, 10x12
S3000DM	Flashpad HD, 14x17
S3000DO	Grid 8:1 10in x 12in
S3000DQ	Grid 8:1 14x17in
S3000DP	Grid 6:1 10in x 12in
S3000DR	Grid 6:1 14x17in
S2000AC	G3 100CM TABLE GRID
S2001AC	G3 110CM TABLE GRID
S2002AC	G3 120CM TABLE GRID
S2003AC	G3 100CM WS GRID
S2004AC	G3 WS 120cm Grid
S2005AC	G3 130CM FOCUS WS GRID
S2006AC	G3 130CM CROSS TABLE GRID
S2007AC	G3 180CM WS GRID
S2009AC	G3 GRID HOLDER
S1200RZ	BAR CODE READER OPTION
S1202RZ	BAR CODE READER W STOOL

S1700JM	DEFINIUM 5000 FIXED TABLE CARBON FIBER TABLE TOP
S39212UP	UPS
S39222FP	TABLE REAR FOOT PEDAL
S39182JW	Compression Band
S39222RR	Hand Grips
S2000BIN	G3 DETECTOR BIN
S3000DL	Detector Battery 10x12in
S3000DN	Detector Battery 14x17in
S3000DT	Battery Charger 14x17in
S3000DS	Battery Charger 10x12in
S2000TTDH	G3 TBL DET HOLDER
S2000TIBCS	G3 SHORT TETHER
E06731BL	MOBILE LATERAL DETECTOR HOLDER V2
S39213TP	267mm Footstool for Patient Barrier
S39212TP	410mm Footstool for Patient Barrier
S2000WA	e-IFU Tech Pub Collector
S1200WA	FRENCH KEYBOARD
S1200WB	DEFINIUM 8000 GERMAN KEYBOARD
S1200WC	DEFINIUM 8000 SPANISH KEYBOARD
S1200WD	DEFINIUM 8000 ITALIAN KEYBOARD
S1200WE	DEFINIUM 8000 PORTUGUESE KEYBOARD
S1200WF	DEFINIUM 8000 SWEDISH KEYBOARD
S1200WK	DEFINIUM 8000 DANISH KEYBOARD
S1200WL	DEFINIUM 8000 DUTCH KEYBOARD
S1200WM	DEFINIUM 8000 NORWEGIAN KEYBOARD
S1202WR	Russia Keyboard
S1201JZ	TABLE AUTO IP PERMANENT
S1201JY	WS AUTO IP PERMANENT
S1201KC	50KW CONFIGURATION PERMAN
S1201KA	65KW CONFIGURATION PERMAN
S1201KB	80KW CONFIGURATION PERMAN
S1201KK	RRA PERMANENT
S1201KL	AUTO FOV PERMANENT
S1201KM	APA PERMANENT
S1201KR	PATIENT ORI-TAG PERMANENT
S1201LA	RRA FLEXTRIAL
S1201LB	AUTO FOV FLEXTRIAL
S1201LC	APA FLEXTRIAL
S1201LX	Operation Console RCIM2
S1200K	GLOBAL 1 D8K INTERNATIONAL LABELS
S1201KS	TABLE AUTO IP FLEXTRIAL
S1200TZ	INFRARED REMOTE
S2000TIBCM	G3 MEDIUM TETHER

S2000TIBCL	G3 LONG TETHER
S3001DK	Int Flashpad HD 10x12
S3001DM	Int Flashpad HD 14x17
S3001DV	Int Flashpad HD 10x12 Vrx
S3001DW	Int Flashpad HD 14x17 Vrx
S3001DX	2x Int FlashPadHD 1417 Vx
S3001DY	2x Int FlashPad HD 14x17
S3000DV	Flashpad HD, 10x12, PKI
S3000DW	Flashpad HD, 14x17, PKI
S2100SLBR	XR656 INMETRO label
S2100SLCN	XR656 CFDA label
S1203JY	WS AUTO IP SW
S1204PLB	WS PASTING BARRIER
E4502ST	25 KAIC X-Ray Main Disconnect Panel 80 Amp, 480 V 208 V
E1500RD	Adventure Series - Bear Brook for XR650 XR656
E1500RR	Adventure Series - Beachside Adventure for XR650 XR656
E1500RL	Adventure Series - Playful Park Adventure for XR650 XR656
E1500CT	Adventure Series - Deep Blue Sea XR Room
S3000GK	GE Grip Flashpad HD 10x12
S3000GM	GE Grip Flashpad HD 14x17
S31DDBS	XR656 HD BS DoD US
S3003DD	IT Security Pack Option
S3003DL	Dose Structured Rprt Opt

GE HUALUN MEDICAL SYSTEMS CO. Ltd has verified the mutual compatibility of the accessories in combination with Discovery XR656 HD and included relevant information to users with Discovery XR656 HD instructions for use. This activity was subject to appropriate methods of internal control and inspection.