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Produkto apžvalga

Reumatologijos diagnostika

ANA Detect
ANAscreen
Anti-alpha-Fodrin IgA, IgG
Anti-C1q
Anti-CCPhs (high sensitive)
Anti-Centromere B
Anti-dsDNA IgA, IgG, IgM, Screen
Anti-Histone
Anti-Jo-1
Anti-MCV
Anti-Nucleosome
Anti-Rib-P
Anti-RNP / Sm
Anti-RNP-70
Anti-Scl-70
Anti-Sm
Anti-SS-A 52
Anti-SS-A 60
Anti-SS-A
Anti-SS-B
Anti-ssDNA
ENAscreen
Rheumatoid Factor IgA, IgG, IgM, Screen

Trombozės diagnostika

Anti-Annexin V IgG, IgM
Anti-beta-2-Glycoprotein I IgA, IgG, IgM, Screen
Anti-Cardiolipin IgA, IgG, IgM, Screen
Anti-Phosphatidic Acid IgG, IgM
Anti-Phosphatidyl Inositol IgG, IgM
Anti-Phosphatidyl Serine IgG, IgM
Anti-Phospholipid Screen IgG, IgM
Anti-Prothrombin IgA, IgG, IgM, Screen

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ANCAscreen hs (high sensitive)
Anti-BPI
Anti-Cathepsin G
Anti-Elastase
Anti-GBM
Anti-Lactoferrin
Anti-Lysozyme
Anti-MPO
Anti-PR3
Anti-PR3 hs (high sensitive)

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AMA-M2
Anti-Gliadin IgA, IgG, Screen
Anti-gp210
Anti-Intrinsic Factor
Anti-Parietal Cell
Anti-LC1
Anti-LKM-1
Anti-Parietal Cell
Anti-SLA
Anti-SP100
Anti-Tissue-Transglutaminase IgA, IgG, Screen
ASCA IgA, IgG
Calprotectin
Lactoferrin

Skydliaukės diagnostika

Anti-TG
Anti-TPO

Diabeto diagnostika

Anti-Insulin

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Serume

Anti-B.pertussis toxin
Anti-Borrelia
Anti-Chlamydia pneumoniae
Anti-Chlamydia trachomatis
Anti-EBV (VCA)
Anti-EBV (ZEBRA) IgM
Anti-Helicobacter pylori
Anti-HSV-1
Anti-HSV-1/2
Anti-HSV-2
Anti-Measles Virus
Anti-Mumps Virus
Anti-Mycoplasma pneumoniae
Anti-Parvovirus
Anti-VZV
Anti Yersinia

Likvore

Anti-Borrelia likvore
Anti-EBV (VCA) likvore
Anti-HSV-1/2 likvore
Anti-Measles Virus likvore
Anti-Rubella Virus likvore
Anti-VZV likvore

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GENERAL THEORY OF OPERATION

The Alegria analyser different modules are described in the following sections ([list of modules](#)).

ID Module

The ID module consists of two barcode readers, one internal and one external. The external is optional for the end user to use.

7. The module is responsible for linking patient data with a sample strip by reading the bar code attached on the strip (CODE 128) with an external bar code reader. The module also handles the identification and positioning of the test strips in the instrument with an internal bar code reader (fixed).

The internal barcode reader is mounted to the chassis in a fixed position. The external barcode reader is mounted to a holder on the right side of the instrument. It is connected to a COM port on the back of the instrument and can be removed for customers that do not wish to use it.

Dispensing Module

The Dispensing module should dispense sample and the liquids that are already pre-filled in the strip (conjugate, substrate, cut-off, sample buffer) into other wells in the strip. It moves in radial direction compared to the turning plate and in vertical direction.

The Dispensing module consists of a precision pump and two separate needles, one for use of dispensing conjugate and one for use of dispensing the other liquids (programmed by process developer with the Alegria Developer SW).

The separate needle for conjugate is to ensure that there is no risk for cross contamination of conjugate with the other liquids.

The precision pump is set to dispense volumes in a range of 10-200 µl and smallest step is 1µl, (The volumes to dispense are programmed by process developer with the Alegria Developer SW).

Each needle can move independently of the other down into a well and up above the strip so the carriage can reposition itself over the next well. Radially the distance between the two needles is fixed so when one needle has to be moved in this direction the other one will follow.

The whole needle feed tubing is filled with hydraulic fluid (distilled water) to achieve as high accuracy as possible. All tubing's and needles for the dispense system can be washed with either distilled water or wash buffer.

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ORGENTEC Diagnostika GmbH