Serascan Diana 2

Serascan Diana 2P

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Serascan Diana 2

Serascan Diana 2P

Read carefully before using the product. For use in "in vitro" diagnosis only.

INTENDED USE

Screening of unexpected antibodies, in gel technique.

INTRODUCTION

The screening of unexpected antibodies has the aim of detecting the clinically significant antibodies present in the patient's sample. The autocontrol will indicate, in a positive screening of unexpected antibodies, whether it is due to the presence of an autoantibody, an alloantibody or both.

The Serascan Diana 2 and Serascan Diana 2P reagent red blood cells present the most significant antigenic determinants for most of the blood group systems.

PRINCIPLE

An antibody reacts specifically with the antigen which stimulated its production¹. According to this, an antibody can be identified depending on its pattern of reactivity when confronted with a panel of reagent red blood cells with a known antigenic configuration. The principle of the test is based on the gel technique described by Y. Lapiere² for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample.

COMPOSITION

Each vial of Serascan Diana 2 (I, II) contains 10 ml of human red blood cells of O blood group, in a 0,8% suspension in a buffered solution with preservatives.

The Serascan Diana 2 reagent is manufactured from a single donor for each vial. It is provided in stoppered glass vials with a white-coloured bulb.

Each vial of Serascan Diana 2P (II) contains 10 ml of papainised human red blood cells of O blood group, in a 0,8% suspension in a buffered solution with preservatives.

The Serascan Diana 2P reagent is manufactured from a single donor for each vial. It is provided in stoppered glass vials with a red-coloured bulb.

The antigenic configuration of the Serascan Diana 2 and Serascan Diana 2P reagent red blood cells is given in the enclosed table (specific for each product lot).

Reagent ready to use. Do not use the reagent if cloudiness and/or large-scale haemolysis are detected.

STABILITY

Serascan Diana 2 and Serascan Diana 2P are stable up to the expiry date stated on the label, if stored at 2-8 °C in the position indicated on the outer packaging. Do not freeze.

Once the vial has been opened, handle correctly, avoiding contamination of the content, and once the product has been used, keep it at the indicated storage temperature.

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- Diagnostic Grifols, S.A. gel cards for the screening of unexpected antibodies.

SAMPLES

Recently extracted blood samples, collected with the usual blood bank anticoagulants (obtaining plasma) or without anticoagulants (obtaining serum).

Do not use haemolysed, cloudy or contaminated samples, or with clot presence.

The procedure for extracting, collecting and handling the blood must be performed by qualified technician personnel according to current standards and directives³, and following the instructions of the manufacturer of the material used for collecting the sample.

- Screening of unexpected antibodies: use serum or plasma. If necessary, samples stored at 2-8 °C may be used up to 48 hours after their extraction, or frozen samples (-20 °C to -80 °C) up to the time of performing the test.

TEST METHOD

Serascan Diana 2 and Serascan Diana 2P can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual. Allow the samples and reagents to reach room temperature (18 – 25 °C).

Inspect the state of the reagent before use (see section LIMITATIONS "In relation to the product"). Before use, homogenise the red blood cells of the Serascan Diana 2 and Serascan Diana 2P vials, by gently inversion.

Manual method


Serascan Diana 2 and Serascan Diana 2P reagent red blood cells may be used with any of the usual methods of screening of unexpected antibodies with Diagnostic Grifols, S.A. gel cards. Follow the instructions for use of the Diagnostic Grifols, S.A. gel cards used.

English

An autocontrol should always be performed in parallel with each test, confronting the patient's red blood cells with his/her own serum or plasma.

RESULTS

Result reading

Negative:	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
	+/-	Scarcely small-sized agglutinations in the lower half of the column.
	1+	Some small-sized agglutinations in the column.
	2+	Small or medium-sized agglutinations throughout the column.
Positive:	3+	Upper band of medium-sized agglutinations, in the upper half of the column.
	4+	Band of agglutinated red blood cells in the upper part of the column.
	DP	Double Population (double band of red blood cells, at the bottom and in the upper part of the column).
	H	Haemolysis (pinkish supernatant and/or gel rosado).

Stability of the results: It is recommended an immediate reading of the results after centrifuging the cards.

If necessary, a delayed reading can be made up until 24 hours after processing the cards, if kept in vertical position, refrigerated (2-8 °C) and sealed with parafilm or similar material, to avoid evaporation of the supernatant.

Interpretation of the results

To know the antigenic configuration of the reagent red blood cells, consult the Serascan Diana 2/ Serascan Diana 2P table enclosed, which is specific for each product lot. Compare the pattern of reactivity obtained with the antigenic profile of the reagent red blood cells used.

Reagent red blood cells	Sample red blood cells (autocontrol)	Interpretation
0	0	Absence of allo- and autoantibodies
+	0	Presence of alloantibodies
0	+	Presence of autoantibodies
+	+	Presence of autoantibodies or autoantibodies + alloantibodies

++ = positive
0 = negative

If necessary, consult the instructions for use of the Diagnostic Grifols, S.A. gel cards used. The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.

Notes:

- If the autocontrol is positive, investigate the presence of cold antibodies, "rouleaux" and perform a Direct Coombs test.

QUALITY CONTROL

- It is recommended to include known positive and negative controls in each series of tests. In the screening of unexpected antibodies, it is recommended to include weak positive controls.
- If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used should be made.

PERFORMANCE CHARACTERISTICS

Although there is no described procedure or technique capable of detecting all the possible unexpected antibodies present in a sample population with absolute certainty, performance evaluation studies carried out in different hospitals with population samples, back up that Serascan Diana 2 and Serascan Diana 2P own performance characteristics conforming to the intended use of the product. These studies⁴ have included tests of screening and identification of unexpected antibodies, being the results obtained similar or better than those obtained with other established products of equivalent intended use.

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LIMITATIONS

In relation to the sample:

Do not use haemolysed, cloudy or contaminated blood samples, or with clot presence.

Sample: Serum/Plasma

1. If plasma is used, complement-dependent haemolytic reactions may not be detected.

2. If poor anti-coagulated plasma or partially coagulated serum is used, fibrin residues may trap non-agglutinated red blood cells at the top of the gel, appearing as a pinkish or reddish layer, but the negative reaction can be interpreted as such.

It is recommended to recolt the serum during 10 minutes at 37 °C, centrifuge and repeat the test!

3. The formation of "rouleaux", caused by an excess of proteins in the serum, the presence of abnormal proteins, drugs, plasma expanders, etc., may cause false positive reactions'.

4. If there are antibodies against high-incidence antigens or multiple antibodies in the sample, all the reagent red blood cells may be agglutinated.

In relation to the product:

Before use, inspect the state of the reagent:

1. Do not use the reagent if cloudiness and/or large-scale haemolysis are observed. These may be caused by microbiological contamination or improper handling/storage.

2. Close the vials correctly after use, to prevent contamination. Ensure that the stoppers of the reagent red blood cell vials have not been switched.

3. If, due to improper transport or storage, the vial is broken and/or there is a loss of content, do not use the product.

4. Improper storage or handling may cause a loss of reactivity of the red blood cells.

5. Do not use the product beyond the expiry date.

6. The Serascan Diana 2 and Serascan Diana 2P reagent red blood cells do not detect the presence of anti-A or anti-B antibodies.

7. Low-incidence antigens may not be represented in the Serascan Diana 2 and Serascan Diana 2P reagent red blood cells, so negative reactions with them do not always indicate absence of an antibody in the sample under study.

8. The Serascan Diana 2P (papainised) reagent red blood cells present a reduced or absent reactivity of some antigens, which appear shaded on the table enclosed with the product.

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- Serascan Diana 2 and Serascan Diana 2P are of human origin and are manufactured using a non-reactive material for the HBS antigen nor for anti-HIV and anti-HCV antibodies, when tested with authorised reagents. Nevertheless, there is no known procedure to ensure that products of human origin do not transmit hepatitis or AIDS. Human blood products and samples should be handled as if they were potentially capable of transmitting infectious diseases.

- The product must only be used by qualified personnel.

- Once used, the product must be disposed of in special containers for biological waste.

- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

BIBLIOGRAPHY

1. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

2. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.

3. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.

4. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard - 4th edition.

5. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline - 2nd edition.

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LITERATŪRA

1. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.
2. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
3. Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
4. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard - 4th edition.
5. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline - 2nd edition.
6. Diana System-Reference Studies, Diagnostic Grifols, 2001.

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INTENDED USE

Determination of the antigens of the ABO and Rh (D) system and determination of the reverse ABO group, in gel technique.

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The ABO system is defined by the presence or absence of the A and/or B antigens in the red blood cells and by the presence of antibodies in the serum corresponding to the antigen or antigens missing in the red blood cells.

The determination of Rh (D) is defined by the presence or absence of the D (RH1) antigen in the red blood cells.

The anti-A, anti-B, anti-AB, anti-D and anti-D' reagents are used to perform the ABO and Rh blood group typing, being complemented with the reverse group test (determination of the reverse ABO group).

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Each column contains polymerised dextran microspheres in buffered medium which act as a filter. The dextrans are mixed with a reagent that contains specific antibodies or a buffer.

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REAGENTS

Each microtube of the DG Gel ABO/Rh (2D) card contains polymerised dextrans in buffered medium with preservatives and mixed with different reagents. The different microtubes are identified by the front label of the card:

- microtube A: monoclonal anti-A (mixture of IgM antibodies of murine origin, clones 16243 G2 and 16247 E6).
- microtube B: monoclonal anti-B (IgM antibodies of murine origin, clone 9621 A8).
- microtube AB: monoclonal anti-AB (mixture of IgM antibodies of murine origin, clones 16245 F11 D8, 16247 E6 and 7821 D9).
- microtube D: monoclonal anti-D (IgM antibodies of human origin, clone P3x61).
- microtube D': monoclonal anti-D (mixture of IgG and IgM antibodies of human origin, clones P3x290, P3x35, P3x61 and P3x21223 B10).
This anti-D monoclonal reagent detects weak D and partial variants of the D antigen, including the D^{vi} variant.
- microtube Ctl.: buffered solution without antibodies (control microtube).
- microtubes N: buffered solution without antibodies (reverse ABO group test).

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STABILITY

DG Gel ABO/Rh (2D) is stable up to the expiry date stated on the label, with the seal intact and stored at 2-8 °C in the position indicated on the outer packaging. Do not freeze.

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- 10 µl, 50 µl and 1 ml automatic pipettes.
- Disposable pipette tips.
- Glass tubes.
- DG Gel Sol.
- Centrifuge for DG Gel cards.
- Reagent red blood cells for reverse group test (A₁/B) from Diagnostic Grifols, S.A.

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- Determination of the antigens of the ABO/Rh system: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction.
Red blood cells from bags collected in ACD, CPD, CPDA or SAG-Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2-8 °C.
If red blood cells from the bag segment are used, it is recommended to wash them with physiological saline solution before preparing the suspension. Do not use if clots or haemolysis are observed.
- Determination of the reverse ABO group: use serum or plasma. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction; or frozen samples (from -20 °C to -80 °C), up to the time of performing the test.

TEST METHOD

DG Gel ABO/Rh (2D) can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

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Manual method:

1. Determination of the antigens of the ABO/Rh system (microtubes A/B/AB/D/D'/Ctl.):
 - Prepare a 5% red blood cell suspension in DG Gel Sol (50 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Ensure the resuspension of the red blood cells before use.
 - Add 10 µl of 5% red blood cell suspension into each of the microtubes indicated.
2. Determination of the reverse ABO group (microtubes N):
 - Homogenise the vials of A₁/B reagent red blood cells.
 - Dispense 50 µl of A₁ reagent red blood cells into microtube N_{A1} and 50 µl of B reagent red blood cells into microtube N_B.
 - Add 50 µl of serum or plasma.
3. Centrifuge in centrifuge for DG Gel cards.
4. Read the results.

RESULTS

Result reading:

Negative:	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
	+/-	Scarce small-sized agglutinations in the lower half of the column.
	1+	Some small-sized agglutinations in the column.
	2+	Small or medium-sized agglutinations throughout the column.
Positive:	3+	Upper band of medium-sized agglutinations in the upper half of the column.
	4+	Band of agglutinated red blood cells in the upper part of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

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Notes:

1. The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.
2. The microtube Ctl. must be negative. If it is positive, invalidate the test.
Repeat the determination after washing the red blood cells with physiological saline solution and preparing a new suspension of the washed red blood cells. If the microtube Ctl. of the repetition is negative, the results of the test can be interpreted; if it is positive, invalidate the test.
3. In the event of obtaining a blood-serum discrepancy, it must be investigated before releasing the result.
4. ABO/Rh system: in the event of +/- to 3+ reactions, presence of a weak antigen must be investigated.

5. Negative reactions with the anti-D reagent must be verified using other reagents and techniques which may detect different variants of the D antigen.
6. In the event of obtaining discrepant results in microtubes D and D', it must be interpreted as a weak or partial D antigen. It is recommended to analyse the expression of this antigen.
7. The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.
8. Occasionally, a red blood cell retention in the incubation chamber may occur with positive 4+ samples, without interfering in result reading.

QUALITY CONTROL

1. It is recommendable to include known positive and negative controls in each series of tests.
2. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

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ABO/Rh system

The diagnostic sensitivity and specificity of the antibodies present in the DG Gel ABO/Rh (2D) card for the determination of the antigens of the ABO and Rh systems, have been studied³ in a representative number of positive and negative samples.

Antibody	No. of samples	Sensitivity ^(a)	Specificity ^(b)
Anti-A	3012	100%	100%
Anti-B	3012	100%	100%
Anti-AB	3012	100%	100%
Anti-D	3012 ^(c)	99.80%	100%
Anti-D'	3012 ^(c)	99.96%	100%

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(c) These samples include 48 samples with weak expression of the D antigen. If only these samples are considered, the anti-D reagent of microtube D has a sensitivity of 89.58% and that of the microtube D', 97.91%. An anti-D reagent is valid for use in transfusion medicine if it presents a minimum sensitivity value of 70% with a panel of Rh D positive samples of weak expression³.

Reverse ABO group

The DG Gel ABO/Rh (2D) card owns performance characteristics suitable for the determination of the reverse group, backed up by a study⁷ where the obtained results were similar to those obtained with other established products of equivalent intended use.

Precision:

The precision of the reagents present in the DG Gel ABO/Rh (2D) card has been determined in a study⁸, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

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Sample: red blood cells

1. Red blood cells from individuals with A or B variants may present a weak expression of the antigens. Antigen expression may be weakened in red blood cells of persons with leukaemia or other malignant diseases⁹.
2. Abnormal concentrations of serum proteins, the presence of macromolecular solutions in the serum/plasma or the presence of Wharton's jelly in cord blood samples may cause the non-specific agglutination of the red blood cells. It is recommended to wash the red blood cells before performing the test⁹.
3. Transfused patients or those subjected to bone marrow transplant may present images of double population⁹.
4. Patients with high-potency cold antibodies may coat the red blood cells completely, causing spontaneous agglutination⁹.

Sample: serum/plasma

1. If plasma is used, complement-dependent haemolytic reactions may not be detected.
2. If poor anti-coagulated plasma or partially coagulated serum is used, fibrin residues may trap nonagglutinated red blood cells at the top of the gel, appearing as a pinkish or reddish layer, but the negative reaction can be interpreted as such. It is recommended to recollect the serum during 10 minutes at 37 °C, centrifuge and repeat the test⁹.
3. Patients with low or inexistent levels of isoagglutinins: newborns up to the age of 4-6 months, elderly persons, patients with immunodeficiency or with very diluted antibodies due to plasma exchange procedures⁹.

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1. Do not use the card if microbiological contamination or change in microtube colour is observed.
2. If due to improper transport or storage, dispersed drops at the top of the microtube are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
3. If trapped bubbles in the gel, any microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.

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- The reagents of the DG Gel ABO/Rh (2D) card of human monoclonal origin are manufactured using a non-reactive material for the HB, antigen nor for anti-HIV and anti-HCV antibodies, when tested with authorised reagents. Nevertheless, there is no known procedure to ensure that products of human origin will not transmit hepatitis or AIDS. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood cell suspension of a concentration other than that indicated may cause false positive or negative results.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

BIBLIOGRAPHY

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
5. Internal study Diagnostic Grifols, S.A. (DG-R1009/32-1), May 2003.
6. Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine, 5:171-184, 1995.
7. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
8. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
9. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

PRESENTATION

210338 DG Gel ABO/Rh (2D) 50 cards

Profile: A/B/AB/D/D'/Ctl./N_{A1}/N_B

Revision date: June 2010

This document is available in several languages. The translations have been made from the master document in Spanish. In the event of doubts or discrepancies, the wording in the master document in Spanish shall take precedence

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LITERATŪRA

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
5. Internal study Diagnostic Grifols, S.A. (DG-RI009/32-1), May 2003.
6. Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine, 5:171- 184, 1995.
7. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
8. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
9. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999

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Determination of the Rh (D), in gel technique.

INTRODUCTION

The clinical importance of the Rh blood group stems from the fact that the D antigen is highly immunogenic. This antigen is the most important of the group of antigens that forms this system and it is made up of a mosaic of epitopes¹.

The determination of Rh (D) is defined by the presence or absence of the D antigen (RH1) in the red blood cells. Individuals with D variants lack a part of the antigen.

The anti-D reagent is used to perform the Rh (D) group typing.

PRINCIPLE

The principle of the test is based on the gel technique described by Y. Lapiere² for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample. The DG Gel card is a plastic support consisting of 8 microtubes. Each microtube consists of a column and a dispensing/incubation chamber.

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The microtubes without antibodies are used as controls.

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COMPOSITION

Each microtube of the DG Gel Anti-D card contains polymerised dextran in buffered medium, with preservatives and mixed with different reagents. The different microtubes are identified by the front label of the card:

- microtubes **D**: monoclonal anti-D (mixture of IgG and IgM antibodies of human origin, clones P3x290, P3x35, P3x61 and P3x21223 B10).

This anti-D monoclonal reagent detects weak D and partial variants of the D antigen, including the Dⁱ variant.

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STABILITY

DG Gel Anti-D is stable up to the expiry date stated on the label, with the seal intact and stored at 2-8 °C in the position indicated on the outer packaging. Do not freeze.

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The procedure of extracting, collecting and handling the blood must be performed by qualified technical personnel according to current standards and directives^{3,4} and following the instructions of the manufacturer of the material used for collecting the sample.

• Determination of the Rh: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction.

Red blood cells from bags collected in ACD, CPD, CPDA or SAG-Mantol can also be used until the expiry date indicated on the label of the bag, if stored at 2-8 °C.

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GRIFOLS

QUALITY CONTROL

- It is recommended to include known positive and negative controls in each series of tests .
- If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

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The diagnostic sensitivity y specificity of the antibodies present in the DG Gel Anti-D card to determine the antigens of the Rh system, have been studied¹ in a representative number of positive and negative samples.

Antibody	No. of samples	Sensitivity ^(a)	Specificity ^(a)
anti-D	3012 ^(a)	99,96%	100%

- (a) Sensitivity: (number of true positive results / (number of true positive results + number of false negative results)) x 100
(b) Specificity: (number of true negative results / (number of true negative results + number of false positive results)) x 100
(c) These samples include 48 samples with weak expression of the D antigen. If only these samples are considered, the sensitivity of the anti-D reagent is 97.91%. An anti-D reagent is valid for use in transfusion medicine if it presents a minimum sensitivity value of 70% with a panel of Rh D positive samples of weak expression¹.

Precision:

The precision of the reagents present in the DG Gel Anti-D card has been determined in a study², including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and the differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

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Sample: red blood cells

- Abnormal concentrations of serum proteins, the presence of macromolecular solutions in the serum/plasma or the presence of Wharton's jelly in cord blood samples may cause the non-specific agglutination of the red blood cells. It is recommended to wash the red blood cells before performing the test¹.
- Transfused patients or those subjected to bone marrow transplant may present images of double population¹.
- Patients with high-potency cold antibodies may coat the red blood cells completely, causing spontaneous agglutination¹.

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- Do not use the card if microbiological contamination or change in microtube colour is observed.
- If due to improper transport or storage, dispersed drops at the top of the microtube are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
- If trapped bubbles, any microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.

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- The reagents of the DG Gel Anti-D card of human monoclonal origin are manufactured using a non-reactive material for the HBs antigen nor for anti-HIV and anti-HCV antibodies, when tested it with authorised reagents. Nevertheless, there is no known procedure to ensure that products of human origin will not transmit hepatitis or AIDS. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood suspension of a concentration other than that indicated may cause false positive or negative results.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

	“In vitro” diagnostic medical device
	Batch code
	Use by
	Temperature limitation
	Consult Instructions for Use
	Catalogue number
	Cards

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- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
- NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline - 2nd edition, 1999.
- Εσωτερική μελέτη Diagnostic Grifols, S.A. (DG-R1009/32-1), Μαΐου 2003.
- Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors.Transfusion Medicine, 5: 171-184, 1995.
- Εσωτερική μελέτη Diagnostic Grifols, S.A. (DR-103/03), Ιουλίου 2003.
- Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

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- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
- NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
- Estudio interno Diagnostic Grifols, S.A. (DG-R1009/32-1), Mayo 2003.
- Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors.Transfusion Medicine, 5: 171-184, 1995.
- Estudio interno Diagnostic Grifols, S.A. (DR-103/03), Julio 2003.
- Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

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LITERATŪRA

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
5. Internal study Diagnostic Grifols, S.A. (DG-RI009/32-1), May 2003.
6. Jones J, Scott ML, Voak D. Monoclonal anti-D specificity and Rh D structure: criteria for selection of monoclonal anti-D reagents for routine typing of patients and donors. Transfusion Medicine, 5: 171-184, 1995.
7. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
8. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

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UAB Diamedica

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Tel. 8 5 279 0080

Read carefully before using the product. For use in “in vitro” diagnosis only

INTENDED USE

Indirect Coombs and Direct Coombs tests, in gel technique.
The Indirect Coombs tests include: screening and identification of unexpected antibodies, crossmatch tests, autocontrol and red blood cell typing.

INTRODUCTION

Coombs introduced the antiglobulin test in clinical medicine in 1945¹, and this is why it is known as the Coombs test. This test is based on the use of a polyspecific anti-human globulin (AHG) that allows the detection of red blood cells coated with immunoglobulins or complement fractions. The Indirect Coombs test allows the detection of red blood cell antibodies present in the patient's serum or plasma, by “in vitro” sensitisation of red blood cells.
The Direct Coombs test allows the detection of red blood cells sensitised “in vivo” by immunoglobulins or complement fractions.
The screening of unexpected antibodies has the aim of detecting the clinically significant antibodies present in the patient's sample. In a positive screening of unexpected antibodies, the autocontrol will indicate whether it is due to the presence of an autoantibody, an alloantibody or both.
Crossmatch tests:

- major crossmatch test: the donor's red blood cells confronted with the recipient's serum or plasma will show the presence or absence of unexpected antibodies in the recipient's blood that are specific against the antigens of the donor's red blood cells.
- minor crossmatch test: the recipient's red blood cells confronted with the donor's plasma or serum will show the presence or absence of unexpected antibodies in the donor's blood that are specific against the antigens of the recipient's red blood cells.

PRINCIPLE

The principle of the test is based on the gel technique described by Y. Lapiere² for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample.
The DG Gel card is a plastic support consisting of 8 microtubes. Each microtube consists of a column and a dispensing/incubation chamber.
Each column contains polymerised dextran microspheres in buffered medium which act as a filter. The dextrans are mixed with a reagent that contains anti-human globulin.
The microtubes that contain anti-human globulin act by agglutinating the red blood cells sensitised, “in vivo” or “in vitro”, by IgG antibodies or complement fractions.
During centrifugation and depending on their size, the red blood cell agglutinations are trapped in the surface or along the gel column. The non-agglutinated red blood cells descend to the bottom of the microtube.

COMPOSITION

Each microtube of the DG Gel Coombs card contains polymerised dextrans in buffered medium with preservatives and mixed with anti-human globulin. The microtubes are identified by the front label of the card:

- microtubes **AHG**: Coombs, buffered low ionic strength solution (LISS) with polyspecific anti-human globulin. Mixture of rabbit polyclonal anti-IgG and monoclonal anti-C3d antibodies, IgM antibodies of murine origin, clone 12011 D10

Reagent ready to use. Use the microtubes immediately once the seal has been opened.

STABILITY

DG Gel Coombs is stable up to the expiry date stated on the label, with the seal intact and stored at 2-25 °C in the position indicated on the outer packaging. Do not freeze.

MATERIAL REQUIRED BUT NOT PROVIDED

- 10 µl, 25 µl, 50 µl and 1 ml automatic pipettes.
- Disposable pipette tips.
- Glass tubes.
- DG Gel Sol.
- Incubator for DG Gel cards.
- Centrifuge for DG Gel cards.
- Reagent red blood cells for screening/identification of unexpected antibodies from Diagnostic Grifols, S.A.
- Haemoclassification sera.

SAMPLES

Recently extracted blood samples collected with the usual blood bank anticoagulants (obtaining plasma) or without anticoagulants (obtaining serum).
Do not use haemolysed, cloudy or contaminated samples or with clot presence.
The procedure for extracting, collecting and handling the blood must be performed by qualified technician personnel according to current standards and directives^{3,4} and following the instructions of the manufacturer of the material used for collecting the sample.

- Screening/identification of unexpected antibodies, crossmatch tests and autocontrol: use serum or plasma. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction; or frozen samples (-20 °C to -80 °C), up to the time of performing the test.
- Direct Coombs, crossmatch tests, autocontrol and red blood cell typing: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction.

Red blood cells from bags collected in ACD, CPD, CPDA or SAG-Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2-8 °C.
If red blood cells from the bag segment are used, it is recommended to wash them with physiological saline solution before preparing the suspension. Do not use if clots or haemolysis are observed.

TEST METHOD

DG Gel Coombs can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.
Allow the samples and reagents to reach room temperature (18-25 °C).
Inspect the state of the cards before use (see section LIMITATIONS "In relation to the product").
Identify the cards and samples to be used.
Each card is individually identified by a bar code. Unless a bar code reader is used, identify them manually.
If manual method is used: carefully peel off the metal film that covers the microtubes to prevent cross-contaminations among them and carefully dispense the red blood cell suspension, avoiding the pipette tip to come into contact with the wall or the content of the microtubes.

Manual method:

1. Indirect Coombs:
 - 1.1. Screening/identification of unexpected antibodies:
 - Homogenise the vials of reagent red blood cells for screening/identification of unexpected antibodies.
 - Dispense 50 µl of reagent red blood cells into the corresponding microtubes.
 - Add 25 µl of the patient's serum or plasma.
 - Incubate 15 minutes at 37 °C.
 - 1.2. Crossmatch tests (CT):
 - Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Use the donor's red blood cells for the major CT and the recipient's red blood cells for the minor CT. Ensure the resuspension of the red blood cells before use.
 - Dispense 50 µl of the donor's 1% red blood cell suspension (major CT) or 50 µl of the recipient's 1% red blood cell suspension (minor CT) into the corresponding microtube.
 - Add 25 µl of the recipient's serum or plasma (major CT) or 25 µl of the donor's serum or plasma (minor CT).
 - Incubate 15 minutes at 37 °C.
 - 1.3. Autocontrol:
 - Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol) for the patient. Ensure the resuspension of the red blood cells before use.
 - Dispense 50 µl of the patient's 1% red blood cell suspension into the corresponding microtube.
 - Add 25 µl of the patient's serum or plasma.
 - Incubate 15 minutes at 37 °C.
 - 1.4. Red blood cell typing:
 - Follow the instructions for use of the haemoclassification serum used.
2. Direct Coombs:
 - Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol) for the patient. Ensure the resuspension of the red blood cells before use.
 - Add 50 µl of the patient's 1% red blood cell suspension into the corresponding microtube.
3. Centrifuge in centrifuge for DG Gel cards.
4. Read the results.

RESULTS

Result reading

Negative:	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column
	+/-	Scarce small-sized agglutinations in the lower half of the column
	1+	Some small-sized agglutinations in the column
	2+	Small or medium-sized agglutinations throughout the column
	3+	Upper band of medium-sized agglutinations in the upper half of the column
Positive:	4+	Band of agglutinated red blood cells in the upper part of the column
	DP	Double Population (double band of red blood cells, at the bottom and in the upper part of the column)

Stability of the results: it is recommended an immediate reading of the results after centrifuging the cards. Do not leave processed cards in horizontal position.
If necessary, a delayed reading can be made up until 24 hours after processing the cards if kept in vertical position, refrigerated (2-8 °C) and sealed with parafilm or similar material to avoid evaporation of the supernatant.
Interpretation of the results:
Indirect Coombs and Direct Coombs tests: tests determined by the result obtained in microtubes AHG.
Red blood cell typing: follow the instructions for use of the haemoclassification serum used.

Notes:

- 1. The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.
- 2. It is recommended to check/investigate any discrepant results obtained.
- 3. The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.
- 4. Direct Coombs: in newborn's samples, a negative result indicates the absence of detectable antibodies, and a positive result that the red blood cells are sensitised (coated with intrauterine antibodies).

QUALITY CONTROL

- 1. It is recommendable to include known positive and negative controls in each series of tests. It is recommended to include weak positive controls for screening/identification of unexpected antibodies and crossmatch tests, and positive controls with heterozygote antigen expression in the red blood cell typing.
- 2. If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

PERFORMANCE CHARACTERISTICS

Although there is no described procedure or technique capable of detecting all the possible unexpected antibodies present in a sample population with absolute certainty, performance evaluation studies^{5,6} carried out with DG Gel Coombs card back up that it owns performance characteristics conforming to the intended use of the product. These studies have included Indirect Coombs and Direct Coombs tests, being the results obtained comparable to those obtained with other established products of equivalent intended use. The precision of the reagents present in the DG Gel Coombs card has been determined in a study⁷, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

LIMITATIONS

In relation to the sample:

Do not use haemolysed, cloudy or contaminated samples or with clot presence.

Sample: red blood cells

- 1. The presence of some drugs, dextran solutions or remains of silicone gel from the extraction tube in the sample, can induce a positive result in the Direct Coombs⁸.

Sample: serum/plasma

- 1. If plasma is used, complement-dependent haemolytic reactions may not be detected.
- 2. If poor anti-coagulated plasma or partially coagulated serum is used, fibrin residues may trap non-agglutinated red blood cells at the top of the gel, appearing as a pinkish or reddish layer, but the negative reaction can be interpreted as such. It is recommended to reclothe the serum during 10 minutes at 37 °C, centrifuge and repeat the test⁹.
- 3. The presence of high concentrations of immunoglobulins or other serum proteins in the sample can neutralise the polyspecific anti-human globulin, even after numerous washes⁸.

In relation to the product:

Before use, inspect the microtubes of the cards:

- 1. Do not use the card if microbiological contamination or change in microtube colour is observed.
- 2. If due to improper transport or storage, dispersed drops at the top of the microtube are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
- 3. If trapped bubbles, any microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.
- 4. In the 24-hour delayed reading of processed cards with weak positive samples, a loss in the agglutination intensity may be observed.

In relation to the method:

- 1. The Indirect Coombs test at 37 °C in gel or glass spheres techniques has been reported to show a level of sensitivity lower than the obtained with the tube technique, in the detection of weak agglutination reactions of the ABO system⁹. To ensure the ABO compatibility between the recipient's and donor's blood, a serological (saline at room temperature with immediate centrifuging) or computer-based (electronic crossmatch) compatibility confirmation, is recommended.
- 2. Screening and identification of unexpected antibodies/crossmatch tests: the test is acceptable if the volumes of serum or plasma are increased from 25 µl to 50 µl. This variation in the concentration of antibodies brings the antigen/antibody ratio down and can improve the detection of antibodies at very low concentration⁸.

WARNINGS AND PRECAUTIONS

The use of "in vitro" diagnostic reagents for professional use requires taking into account the following indications:

- The reagents of the DG Gel Coombs card are not of human origin. HBV, HCV or HIV contamination is excluded. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood cell suspension of a concentration other than that indicated may cause false positive or negative results.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use red blood cells treated with enzymes in microtubes with anti-human globulin (AHG).
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

BIBLIOGRAPHY

- 1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- 2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- 3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
- 4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
- 5. Internal study Diagnostic Grifols, S.A. (DG-RI009/32-3), October 2003.
- 6. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
- 7. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
- 8. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.
- 9. Phillips P et al. An explanation and the clinical significance of the failure of microcolumn tests to detect weak ABO and other antibodies. Transfusion Medicine, 7: 47-53, 1997.

PRESENTATION

210342 DG Gel Coombs 50 cards Profile: 8x(AHG)

Revision date: June 2010

This document is available in several languages. The translations have been made from the master document in Spanish. In the event of doubts or discrepancies, the wording in the master document in Spanish shall take precedence.

	"In vitro" diagnostic medical device
	Batch code
	Use by
	Temperature limitation
	Consult Instructions for Use
	Catalogue number
	Cards

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LITERATŪRA

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
5. Internal study Diagnostic Grifols, S.A. (DG-RI009/32-3), October 2003.
6. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
7. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
8. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.
9. Phillips P et al. An explanation and the clinical significance of the failure of microcolumn tests to detect weak ABO and other antibodies. Transfusion Medicine, 7: 47-53, 1997..

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

INTENDED USE

Physical medium for saline and enzymatic tests, in gel technique.

The saline and enzymatic technique tests include: screening and identification of unexpected antibodies, crossmatch tests, autocontrol, red blood cell typing and determination of the reverse ABO group.

INTRODUCTION

The reaction between the red blood cell antigens and their corresponding antibodies is easily detected in an agglutination reaction¹.

The agglutination reaction in saline technique occurs when the antigens contact the antibodies in a saline medium.

The agglutination reaction in the enzymatic technique implies the addition of proteolytic enzymes (papain, trypsin, etc.) to the saline medium to improve the reaction between antigens and antibodies.

The screening of unexpected antibodies has the aim of detecting the clinically significant antibodies present in the patient's sample. In a positive screening of unexpected antibodies, the autocontrol will indicate whether it is due to the presence of an autoantibody, an alloantibody or both.

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The ABO system was the first human blood group system discovered, by Landsteiner in 1900¹, and is still the most important in transfusion practice.

The ABO system is defined by the presence or absence of the A and/or B antigens in the red blood cells (forward ABO group) and by the presence of antibodies in the serum (reverse ABO group), corresponding to the antigen or antigens missing in the red blood cells.

PRINCIPLE

The principle of the test is based on the gel technique described by Y. Lapiere² for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample.

The DG Gel card is a plastic support consisting of 8 microtubes. Each microtube consists of a column and a dispensing/incubation chamber.

Each column contains polymerised dextran microspheres in buffered medium which act as a filter.

The microtubes without antibodies are used in techniques in which the antibodies react directly with the red blood cells in the incubation chamber.

During centrifugation and depending on their size, the red blood cell agglutinations are trapped in the surface or along the gel column. The non-agglutinated red blood cells descend to the bottom of the microtube.

COMPOSITION

Each microtube of the DG Gel Neutral card contains polymerised dextrans in buffered medium and preservatives. The microtubes are identified by the front label of the card:

- microtubes **N**: buffered solution without antibodies (neutral microtubes)

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STABILITY

DG Gel Neutral is stable up to the expiry date stated on the label, with the seal intact and stored at 2-25 °C in the position indicated on the outer packaging. Do not freeze.

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- 10 µl, 25 µl, 50 µl and 1 ml automatic pipettes.
- Disposable pipette tips.
- Glass tubes.
- DG Gel Sol.
- Incubator for DG Gel cards
- Centrifuge for DG Gel cards
- Reagent red blood cells for screening/identification of unexpected antibodies from Diagnostic Grifols, S.A.
- Papain reagent.
- Haemoclassification sera.
- Reagent red blood cells for reverse group test from Diagnostic Grifols, S.A.

SAMPLES

Recently extracted blood samples collected with the usual blood bank anticoagulants (obtaining plasma) or without anticoagulants (obtaining serum).

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The procedure of extracting, collecting and handling the blood must be performed by qualified technician personnel according to current standards and directives^{3,4} and following the instructions of the manufacturer of the material used for collecting the sample.

- Screening/identification of unexpected antibodies, crossmatch tests, autocontrol and determination of the reverse ABO group: use serum or plasma. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction; or frozen samples (from -20 °C to -80 °C), up to the time of performing the test.
- Crossmatch tests, autocontrol and red blood cell typing: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction.

Red blood cells collected in ACD, CPD, CPDA or SAG-Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2-8 °C.

If red blood cells from the bag segment are used, it is recommended to wash them with physiological saline solution before preparing the suspension. Do not use if clots or haemolysis are observed.

TEST METHOD

DG Gel Neutral can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

Allow the samples and reagents to reach room temperature (18-25 °C).

Inspect the state of the cards before use (see section LIMITATIONS "In relation to the product").

Identify the cards and samples to be used.

Each card is individually identified by a bar code. Unless a bar code reader is used, identify them manually.

If manual method is used: carefully peel off the metal film that covers the microtubes to prevent cross-contaminations among them and carefully dispense the red blood cell suspension, avoiding the pipette tip to come into contact with the wall or the content of the microtubes.

Manual method

1. Screening/identification of unexpected antibodies:

1.1. Saline technique:

- Homogenise the vials of reagent red blood cells for screening/identification of unexpected antibodies
- Dispense 50 µl of reagent red blood cells into the corresponding microtubes.
- Add 25 µl of the patient's serum or plasma.
- Incubate at the chosen temperature (37 °C, 18-25 °C or 2-8 °C) during 15 minutes.

1.2. Enzymatic technique:

- Homogenise the vials of reagent red blood cells (papainised) for screening/identification of unexpected antibodies.
- Dispense 50 µl of reagent red blood cells into the corresponding microtubes.
- Add 25 µl of the patient's serum or plasma.
- Incubate at 37 °C during 15 minutes.

2. Crossmatch tests (CT):

- Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Use the donor's red blood cells for the major CT and the recipient's red blood cells for the minor CT. Ensure the resuspension of the red blood cells before use.

2.1. Saline technique:

- Dispense 50 µl of the donor's 1% red blood cell suspension (major CT) or 50 µl of the recipient's 1% red blood cell suspension (minor CT) into the corresponding microtube.
- Add 25 µl of the recipient's serum or plasma (major CT) or 25 µl of the donor's serum or plasma (minor CT).
- Incubate at the chosen temperature (37 °C, 18-25 °C or 2-8 °C) during 15 minutes.

2.2. Enzymatic technique:

Follow the instructions for use of Papain reagent.

3. Autocontrol:

- Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol) for the patient. Ensure the resuspension of the red blood cells before use.

3.1. Saline technique:

- Dispense 50 µl of the patient's 1% red blood cell suspension into the corresponding microtube.
- Add 25 µl of the patient's serum or plasma.
- Incubate at the chosen temperature (37 °C, 18-25 °C or 2-8 °C) during 15 minutes.

3.2. Enzymatic technique:

Follow the instructions for use of Papain reagent.

4. Red blood cell typing:

Follow the instructions for use of the haemoclassification serum used.

5. Determination of the reverse ABO group.

- Homogenise the vials of reagent red blood cells for reverse group test.
 - Dispense 50 µl of reagent red blood cells into the corresponding microtubes.
 - Add 50 µl of serum or plasma.
6. Centrifuge in centrifuge for DG Gel cards
7. Read the results.
- Note: if the saline technique is being done at 2-8 °C (cold agglutinins), it is recommended to previously cool cards and reagents during 30 minutes before performing the test.

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Stability of the results: it is recommended an immediate reading of the results after centrifuging the cards. Do not leave processed cards in horizontal position.

If necessary, a delayed reading can be made up until 24 hours after processing the cards if kept in vertical position, refrigerated (2-8 °C) and sealed with parafilm or similar material to avoid evaporation of the supernatant.

Interpretation of the results

Screening/identification of unexpected antibodies, crossmatch tests and autocontrol: tests determined by the result obtained in microtubes N.

Red blood cell typing: follow the instructions for use of the haemoclassification serum used.

Reverse group: follow the instructions for use of the reagent red blood cells used for reverse group.

- Notes:
- The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.
 - It is recommendable to check/investigate any discrepant results obtained (specially the reverse group obtained with the previous ABO group).
 - The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.
 - Occasionally, red blood cell retention in the incubation chamber may occur with positive 4+ samples, without interfering in result reading.

QUALITY CONTROL

- It is recommendable to include known positive and negative controls in each series of tests.
It is recommended to include weak positive controls for screening/identification of unexpected antibodies and crossmatch tests, and positive controls with heterozygote antigen expression in red blood cell typing.
- If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

PERFORMANCE CHARACTERISTICS

The DG Gel Neutral card owns performance characteristics suitable for the intended use of the product, backed up by a study⁵ which has included saline and enzymatic tests, where results obtained were similar to those obtained with other established products of equivalent intended use.

The precision of the DG Gel Neutral card has been determined in a study⁶, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

LIMITATIONS

In relation to the sample:

- Do not use haemolysed, cloudy or contaminated samples or with clot presence.
- Sample: serum/plasma
- If plasma is used, complement-dependent haemolytic reactions may not be detected
 - If poor anti-coagulated plasma or partially coagulated serum is used, fibrin residues may trap non-agglutinated red blood cells at the top of the gel, appearing as a pinkish or reddish layer, but the negative reaction can be interpreted as such.
- It is recommended to replot the serum at 37 °C during 10 minutes, centrifuge and repeat the test⁷.

- Patients with low or inexistent levels of isoagglutinins: newborns up to the age of 4-6 months, elderly persons, patients with immunodeficiency or very diluted antibodies due to plasma exchange procedures⁸.

In relation to the product:

- Before use, inspect the microtubes of the cards:
- Do not use the card if microbiological contamination or change in microtube colour is observed.
 - If due to improper transport or storage, dispersed drops at the top of the microtube, are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
 - If trapped bubbles, any microtubes without supernatant, a decrease in the gel volume or a cracked gel are observed, do not use the card.
 - In the 24-hour delayed reading of processed cards with weak positive samples, a loss in the agglutination intensity may be observed.

In relation to the method:

- Screening and identification of unexpected antibodies/crossmatch tests: the test is acceptable if the volumes of serum or plasma are increased from 25 µl to 50 µl. This variation in the concentration of antibodies brings the antigen/antibody ratio down and can improve the detection of antibodies at very low concentrations⁹.

WARNINGS AND PRECAUTIONS

- The use of “in vitro” diagnostic reagents for professional use requires taking into account the following indications:
- The reagents of the DG Gel Neutral card are not of human origin. HBV, HCV or HIV contamination is excluded. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
 - The product must only be used by qualified personnel.
 - A red blood cell suspension of a concentration other than that indicated may cause false positive or negative results.
 - The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
 - The addition of volumes other than those indicated in the method may modify the reaction.
 - Do not use the card beyond the expiry date.
 - Once used, the product must be disposed of in special containers for biological waste.
 - If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

BIBLIOGRAPHY

- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
- NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
- Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
- Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
- Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

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LITERATŪRA


1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
5. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
6. Internal study Diagnostic Grifols, S.A. (DR-103/03), July 2003.
7. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

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INTENDED USE

Evaluation of positive Direct Coombs samples in gel technique.
It allows to differentiate red blood cells sensitised “in vivo” by IgG type immunoglobulins or the complement C3d fraction.

INTRODUCTION

Coombs introduced the antiglobulin test in clinical medicine in 1945¹, and this is why it is known as the Coombs test. This test is based on the use of a polyspecific anti-human globulin (AHG) that allows the detection of red blood cells coated with immunoglobulins or complement fractions. The Direct Coombs test, with the anti-IgG and anti-C3d reagents, allows to detect and differentiate red blood cells sensitised “in vivo” by IgG type immunoglobulins or the complement C3d fraction.

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Each column contains polymerised dextran microspheres in buffered medium which act as a filter. The dextrans are mixed with a reagent that contains specific antibodies, anti-human globulin or a buffer.

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The microtubes that contain anti-human globulin act by agglutinating the sensitised red blood cells, “in vivo” or “in vitro”, with IgG antibodies or complement fractions.

The microtubes without antibodies are used for controls.

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COMPOSITION

Each microtube of the DG Gel DC Scan card contains polymerised dextrans in buffered medium with preservatives and mixed with different reagents. The different microtubes are identified by the front label of the card:

- microtubes **AHG**: Coombs, buffered low ionic strength solution (LISS) with polyspecific anti-human globulin. Mixture of rabbit polyclonal anti-IgG and monoclonal anti-C3d antibodies, IgM antibodies of murine origin, clone 12011 D10
- microtubes **IgG**: (rabbit polyclonal) anti-IgG in low ionic strength solution
- microtubes **C3d**: monoclonal anti-C3d (IgM antibodies of murine origin, clone 12011 D10) in low ionic strength solution
- microtubes **Ctl.**: buffered solution without antibodies (control microtubes)

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STABILITY

DG Gel DC Scan is stable up to the expiry date stated on the label, with the seal intact and stored at 2-8 °C in the position indicated on the outer packaging. Do not freeze.

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SAMPLES

Recently extracted blood samples collected with the usual blood bank anticoagulants.

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- Direct Coombs: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction. Red blood cells from bags collected in ACD, CPD, CPDA or SAG-Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2-8 °C.

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TEST METHOD

DG Gel DC Scan can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

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Manual method:

1. Direct Coombs (microtubes AHG/IgG/C3d/Ctl.):
 - Prepare a 1% red blood cell suspension in DG Gel Sol (10 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Ensure the resuspension of the red blood cells before use.
 - Add 50 µl of 1% red blood cell suspension to each of the microtubes indicated.
2. Centrifuge in centrifuge for DG Gel cards.
3. Read the results.

RESULTS

Result reading:

Negative:	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
	+/-	Scarce small-sized agglutinations in the lower half of the column.
	1+	Some small-sized agglutinations in the column.
	2+	Small or medium-sized agglutinations throughout the column.
	3+	Upper band of medium-sized agglutinations in the upper half of the column.
Positive:	4+	Band of agglutinated red blood cells in the upper part of the column.
	DP	Double Population (double band of red blood cells, at the bottom and in the upper part of the column).

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Interpretation of the results:

Direct Coombs test: test determined by results obtained in microtubes AHG, IgG and C3d.

Notes:

1. The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.
2. The microtube Ctl. must be negative. If it is positive, invalidate the test.
Repeat the determination after washing the red blood cells with physiological saline solution and preparing a new suspension of the washed red blood cells. If the microtube Ctl. of the repetition is negative, the results of the test can be interpreted; if it is positive, invalidate the test.
3. The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.
4. The positive results obtained in microtube AHG must be contrasted with those obtained in microtubes IgG and C3d. It is recommended to check/investigate any discrepant results obtained.
5. Direct Coombs: in newborn's samples, a negative result indicates the absence of detectable antibodies and a positive result, that the red blood cells are sensitised (coated with intrauterine antibodies).

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PERFORMANCE CHARACTERISTICS

Although there is no described procedure or technique capable of detecting all the possible unexpected antibodies present in a sample population with absolute certainty, a performance evaluation study⁵ carried out with DG Gel DC Scan card back up that it owns performance characteristics conforming to the intended use of the product. This study has included Direct Coombs tests, being the results obtained comparable to those obtained with other established products of equivalent intended use.

The precision of the reagents present in the DG Gel DC Scan card has been determined in a study⁶, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

LIMITATIONS

In relation to the sample:

Do not use haemolysed, cloudy or contaminated samples or with clot presesnce.

Sample: red blood cells

- 1. The presence of some drugs, dextran solutions or remains of silicone gel from the extraction tube in the sample can induce a positive result in the Direct Coombs⁷.

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- 1. Do not use the card if microbiological contamination or change in microtube colour is observed.
- 2. If due to improper transport or storage, dispersed drops at the top of the microtube are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
- 3. If trapped bubbles in the gel, a microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.

WARNINGS AND PRECAUTIONS

The use of “in vitro” diagnostic reagents for professional use requires taking into account the following indications:

- The reagents of the DG Gel DC Scan card are not of human origin. HBV, HCV or HIV contamination is excluded. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood suspension of a concentration other than that indicated may cause false positives or negatives.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

BIBLIOGRAPHY

- 1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- 2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- 3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 4th edition, 1998.
- 4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
- 5. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
- 6. Internal study Diagnostic Grifols, S.A. (DR-103/03), July 2003.
- 7. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

PRESENTATION

210345 DG Gel DC Scan 25 cards Profile: 2x(AHG/IgG/C3d/Ctl.)

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IVD

“In vitro” diagnostic medical device

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LITERATŪRA

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition Blackwell Scientific Publications, Oxford, 1993.
2. Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A4: Procedures for the collection of diagnostic blood specimens by venipuncture Approved Standard, 4th edition, 1998.
4. NCCLS H18-A2: Procedures for the handling and processing of blood specimens; Approved Guideline, 2nd edition, 1999.
5. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
6. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
7. Technical Manual, 13th edition, American Association of Blood Banks, Bethesda, 1999.

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INTENDED USE

Confirmation of the blood groups of the ABO and Rh (D) system, in gel technique.

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The ABO system is defined by the presence or absence of the A (ABO1) and/or B (ABO2) antigens in the red blood cells and by the presence of antibodies in the serum corresponding to the antigen or antigens missing in the red blood cells.

The determination of Rh (D) is defined by the presence or absence of the D antigen (RH1) in the red blood cells.

The anti-A, anti-B and anti-D^W- reagents are used to perform the ABO and Rh blood group typing.

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The microtubes without antibodies are used in techniques in which the antibodies react directly with the red blood cells in the incubation chamber and for controls

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REAGENTS

Each microtube of the DG Gel Confirm P card contains polymerised dextrans in buffered medium, with preservatives and mixed with different reagents. The different microtubes are identified by the front label of the card:

- microtubes A: monoclonal anti-A (IgM antibodies of murine origin, clone Birma-1).
- microtubes B: monoclonal anti-B (IgM antibodies of murine origin, LB-2).
- microtubes D^W–: monoclonal anti-D (IgM antibodies of human origin, clone MS-201).
- microtubes Ctl.: buffered solution without antibodies (control microtubes).

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STABILITY

DG Gel Confirm P is stable up to the expiry date stated on the label, with the seal intact and stored at 2–25 °C in the position indicated on the outer packaging. Do not freeze.

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- 10 µl, 50 µl and 1 ml automatic pipettes.
- Disposable pipette tips.
- Glass tubes.
- DG Gel Sol.
- Centrifuge for DG Gel cards.

SAMPLES

Blood samples recently collected with the routinely used anticoagulants (EDTA, citrate and heparine). Do not use haemolysed, cloudy or contaminated samples or with clot presence.

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- Determination of the antigens of the ABO/Rh system: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2–8 °C can be used up to 48 hours after their extraction.

Red blood cells from bags collected in CPD, CPDA or SAG–Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2–8 °C.

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TEST METHOD

DG Gel Confirm P can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

Allow the samples and reagents to reach room temperature (18–25 °C).

Inspect the state of the cards before use (see section LIMITATIONS “In relation to the product”).

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Manual method:

1. Confirmation of the ABO/Rh groups (microtubes A/B/D^W-/Ctl.):
 - Prepare a 5% red blood cell suspension in DG Gel Sol (50 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Ensure the resuspension of the red blood cells before use.
 - Add 10 µl of 5% red blood cell suspension to each of the microtubes indicated.
2. Centrifuge in centrifuge for DG Gel cards.
3. Read the results.

RESULTS

Result Reading:

Negative:	–	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column
	+ / –	Scarce small-sized agglutinations in the lower half of the column.
	1+	Some small-sized agglutinations in the column.
	2+	Small or medium-sized agglutinations throughout the column.
	3+	Upper band of medium-sized agglutinations in the upper half of the column.
Positive:	4+	Band of agglutinated red blood cells in the upper part of the column.
	DP	Double Population (double band of red blood cells, at the bottom and in the upper part of the column).

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If necessary, a delayed reading can be made up until 24 hours after processing the cards if kept in vertical position, refrigerated (2–8 °C) and sealed with parafilm or similar material to avoid evaporation of the supernatant.

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+ = positive
0 = negative

1. The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.
2. The microtube Ctl. must be negative. If it is positive, invalidate the test.
Repeat the determination after washing the red blood cells with physiological saline solution and preparing a new suspension of the washed red blood cells. If the microtube Ctl. of the repetition is negative, the results of the test can be interpreted; if it is positive, invalidate the test.
3. ABO/Rh system: in the event of +/- to 3+ reactions, presence of a weak antigen must be investigated.
4. Negative reactions with the anti-D^{VI}- reagent must be verified, using other reagents and techniques which may detect different variants of the D antigen.
5. The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.
6. Occasionally, a red blood cell retention in the incubation chamber may occur with positive 4+ samples, without interfering in result reading.

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ABO/Rh system

The diagnostic sensitivity and specificity of the antibodies present in the DG Gel Confirm P card^{5,6}, for the determination of the antigens of the ABO and Rh systems, have been studied in a representative number of positive and negative samples.

Antibody	No. of Samples	Sensitivity ^(a)	Specificity ^(b)
Anti-A	3041 ⁵	100%	100%
Anti-B	3046 ⁵	100%	100%
Anti-D ^{VI-}	1329 ^{(c),6}	99,7%	100%

(a) Sensitivity: $(\text{number of true positive results} / (\text{number of true positive results} + \text{number of false negative results})) \times 100$

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(c) These samples include 49 samples with weak expression of the D antigen. If only these samples are considered, the sensitivity of the anti-D^{III} reagent is 93.9%.

Precision:

The precision^{7,8} of the reagents present in the DG Gel Confirm P card has been determined, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

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Sample: red blood cells

1. Red blood cells from individuals with A or B variants may present a weak expression of the antigens. Antigen expression may be weakened in red blood cells of people with leukaemia or other malignant diseases⁹.
2. Abnormal concentrations of serum proteins, the presence of macromolecular solutions in the serum/plasma or the presence of Wharton's jelly in cord blood samples may cause the non-specific agglutination of the red blood cells. It is recommended to wash the red blood cells before performing the test⁹.
3. Transfused patients or those subjected to bone marrow transplant may present images of double population⁹.
4. Patients with high-potency autoantibodies may coat the red blood cells completely, causing spontaneous agglutination⁹.

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1. Do not use the card if microbiological contamination or change in microtube colour is observed.
2. If due to improper transport or storage, dispersed drops at the top of the microtube are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
3. If trapped bubbles in the gel, any microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.
4. The anti-A reagent may not react with weak variants of A_x phenotype.

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- The reagents of the DG Gel Confirm P card of human monoclonal origin are manufactured using a non-reactive material for the HB_s antigen nor for anti-HIV and anti-HCV antibodies, when tested with authorised reagents. Nevertheless, there is no known procedure to ensure that products of human origin will not transmit hepatitis or AIDS. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood suspension of a concentration other than that indicated may cause false positive or negative results.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorised distributor in your country.

BIBLIOGRAPHY

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapiere Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
4. NCCLS H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.
5. External evaluation Diagnostic Grifols, S.A. (REGD-0002998)M Noviembre 2007.
6. External evaluation and internal study Diagnostic Grifols, S.A. (REGD-0012469), January 2013
7. Internal study Diagnostic Grifols, S.A. (DR-103/03), Julio 2003.
8. Internal study Diagnostic Grifols, S.A. (DG-RI013/07-3), Mayo 2007.
9. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005

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LITERATŪRA

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
4. NCCLS H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.
5. External evaluation Diagnostic Grifols, S.A. (REGD-0002998)M Noviembre 2007.
6. External evaluation and internal study Diagnostic Grifols, S.A. (REGD-0012469), January 2013
7. Internal study Diagnostic Grifols, S.A. (DR-I03/03), Julio 2003.
8. Internal study Diagnostic Grifols, S.A. (DG-RI013/07-3), Mayo 2007.
9. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005.

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INTENDED USE

DG Gel Sol is a reagent for preparing red blood cell suspensions used in gel techniques (DG Gel).

PRINCIPLE

A reduction in the ionic strength of the medium of the red blood cell suspension enhances significantly the association of the antibody with the antigen. DG Gel Sol is a low ionic strength solution that facilitates the attachment of the antibody to the red blood cells by lowering the density of the cation cloud around them, thus improving sensitization.

COMPOSITION

Each DG Gel Sol kit contains:

DG Gel Sol:

- 2 vials of 100 ml, liquid. Each vial contains low ionic strength buffered solution, being its main components glycine 1.37% and glucose 0.85%.

Preparation of the reagents

DG Gel Sol:

Ready to use.

Stability once opened: at 2-8 °C, it has the same stability as when unopened, provided that it is handled and stored correctly, avoiding contamination of the contents and prolonged exposure to temperatures that are not suitable for storage.

Do not freeze.

STABILITY

All unopened reagents are stable up to the expiry date stated on the label if stored at 2-8 °C.

METHOD

DG Gel Sol can be used manually or with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

Manual method (DG Gel):

Follow the instructions for use of the DG Gel card to be used, particularly the points that indicate the use of DG Gel Sol.

LIMITATIONS

Limitations in relation to the reagent:

Before use, inspect the reagent.

1. Do not use if the container is deteriorated and there is a loss of content.
2. Do not use if microbial or fungal contamination is observed.

WARNINGS AND PRECAUTIONS

The use of “in vitro” diagnosis reagents for professional use requires taking into account the following indications:

- The product must only be used by qualified personnel.
- A red blood cell suspension of a concentration other than that indicated may cause false positives or negatives.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Follow the recommended procedures strictly, periodically checking the instruments used.
- Do not use the product beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have doubts or need further information on the use of this product, consult the authorised distributor in your country.

PRESENTATION

210354 DG Gel Sol 2x100ml

Revision date: April 2004

This document is available in several languages. The translations have been made from the master document in English. In the event of doubts or discrepancies, the wording in the master document in English shall take precedence.

IVD	In vitro diagnostic medical device
LOT	Batch code
	Use by
	Temperature limitation
	Consult Instructions for Use
REF	Catalogue number

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INTENDED USE

Determination of the antigens of the ABO and Rh (D) systems and determination of the reverse ABO group test, in gel technique.

INTRODUCTION

The ABO system was the first human blood group system discovered, by Landsteiner in 1900¹, and is still the most important in transfusion practice. The ABO system is defined by the presence or absence of the A (ABO1) and/or B (ABO2) antigens in the red blood cells and by the presence of antibodies in the serum, corresponding to the antigen or antigens missing in the red blood cells.

The clinical importance of the Rh blood group stems from the fact that the D (RH1) antigen is highly immunogenic. The determination of Rh (D) is defined by the presence or absence of the D (RH1) antigen in the red blood cells.

The anti-A, anti-B and anti-D^W reagents are used to perform the ABO and Rh blood group typing, being complemented with the reverse group test (determination of the reverse ABO group).

PRINCIPLE OF THE TEST

The principle of the test is based on the gel technique described by Y. Lapierre² for detecting red blood cell agglutination reactions. The agglutination occurs when the red blood cell antigens contact the corresponding antibodies, present in the reagent or in the serum or plasma sample.

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Each column contains polymerised dextran microspheres in buffered medium which act as a filter. The dextran are mixed with a reagent that contains specific antibodies or a buffer.

The microtubes that contain the specific antibodies into the gel solution act as a reaction medium and the red blood cells agglutinate in contact with the antibodies.

The microtubes without antibodies are used in techniques in which the antibodies react directly with the red blood cells in the incubation chamber and for controls.

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REAGENTS

Each microtube of the DG Gel ABO/Rh (CR) card contains polymerised dextran in buffered medium with preservatives and mixed with different reagents. The different microtubes are identified by the front label of the card:

- microtube **A**: monoclonal anti-A (IgM antibodies of murine origin, clone DAM-1).
- microtube **B**: monoclonal anti-B (IgM antibodies of murine origin, clone 9621 A8).
- microtube **₍₁₎D^W**: monoclonal anti-D (IgM antibodies of human origin, clone P3x61).
- microtube **₍₂₎D^W**: monoclonal anti-D (IgM antibodies of human origin, clone MS-201).
- microtube **Ctl.**: buffered solution without antibodies (control microtube).
- microtubes **N**: buffered solution without antibodies (reverse ABO group test).

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STABILITY

DG Gel ABO/Rh (CR) cards are stable up to the expiry date stated on the label, with the seal intact and stored at 2-25 °C in the position indicated on the outer packaging. Do not freeze.

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- 10 µl, 50 µl and 1 ml automatic pipettes.
- Disposable pipette tips.
- Glass tubes.
- DG Gel Sol.
- Centrifuge for DG Gel cards.
- Reagent red blood cells for reverse A₁/A₂/B/O group test from Diagnostic Grifols, S.A.

SAMPLES

Blood samples recently collected with the routinely used anticoagulants (EDTA, citrate and heparin, for plasma collection) or without anticoagulants (for serum collection).

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The procedure for extracting, collecting and handling the blood must be performed by qualified technician personnel according to current standards and directives^{3,4} and following the instructions of the manufacturer of the material used for collecting the sample.

- Determination of the antigens of the ABO/Rh system: use the red blood cells collected with anticoagulants. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction. Red blood cells from bags collected in CPD, CPDA or SAG-Mannitol can also be used until the expiry date indicated on the label of the bag, if stored at 2-8 °C. If red blood cells from the bag segment are used, it is recommended to wash them with physiological saline solution before preparing the suspension. Do not use if clots or haemolysis are observed.
- Determination of the reverse ABO group: use serum or plasma. If necessary, samples stored at 2-8 °C can be used up to 48 hours after their extraction; or frozen samples (from -20 °C to -80 °C), up to the time of performing the test.

TEST METHOD

DG Gel ABO/Rh (CR) can be used both in a manual method and with semi-automatic or automatic instruments. For the automatic system, see the instrument user manual.

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Inspect the state of the cards before use (see section on LIMITATIONS "In relation to the product").

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Each card is individually identified by a bar code. Unless a bar code reader is used, identify them manually. If manual method is used: carefully peel off the metal film that covers the microtubes to prevent cross-contaminations among them and carefully dispense the red blood cell suspension, avoiding the pipette tip to come into contact with the wall or the content of the microtubes.

Manual method:

1. Determination of the antigens of the ABO/Rh system (microtubes A/B/₍₁₎D^W/₍₂₎D^W/Ctl.):
 - Prepare a 5% red blood cell suspension in DG Gel Sol (50 µl of red blood cell sediment or concentrate in 1 ml of DG Gel Sol). Ensure the resuspension of the red blood cells before use.
 - Add 10 µl of 5% red blood cell suspension to each of the microtubes indicated.
2. Determination of the reverse ABO group (microtubes N):
 - Homogenise the reagent red blood cell vials A₁/A₂/B/O that are going to be used.
 - Dispense 50 µl of the reagent red blood cells A₁, A₂, B or O, depending on your choice, into the three microtubes N.
 - Add 50 µl of serum or plasma.
3. Centrifuge in centrifuge for DG Gel cards.
4. Read the results.

RESULTS

Result reading:

Negative:	-	Band of red blood cells at the bottom of the column and no visible agglutinations in the rest of the column.
Positive:	+/-	Scarce small-sized agglutinations in the lower half of the column.
	1+	Some small-sized agglutinations in the column.
	2+	Small or medium-sized agglutinations throughout the column.
	3+	Upper band of medium-sized agglutinations in the upper half of the column.
	4+	Band of agglutinated red blood cells in the upper part of the column.
DP		Double Population (double band of red blood cells, at the bottom and in the upper part of the column).

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Interpretation of the results:

ABO System

In case of using reagent red blood cells for reverse A₁, B and O group tests, the results should be read as follows:

Forward ABO group			Reverse ABO group			ABO group
Microtube A	Microtube B	Microtube Ctl.	Microtube N+ reagent red blood cells A ₁	Microtube N+ reagent red blood cells B	Microtube N+ reagent red blood cells O	
0	0	0	+	+	0	0
+	0	0	0	+	0	A
0	+	0	+	0	0	B
+	+	0	0	0	0	AB

In case of using reagent red blood cells for reverse A₁, A₂ and B group tests, the results should be read as follows:

Forward ABO group			Reverse ABO group			ABO group
Microtube A	Microtube B	Microtube Ctl.	Microtube N+ reagent red blood cells A ₁	Microtube N+ reagent red blood cells A ₂	Microtube N+ reagent red blood cells B	
0	0	0	+	+	+	0
+	0	0	0	0	+	A
0	+	0	+	+	0	B
+	+	0	0	0	0	AB
+	0	0	+	0	+	A ₂ with anti-A ₁
+	+	0	+	0	0	A ₂ B with anti-A ₁

Rh system (D antigen)

Microtube ₍₁₎ D ^W -	Microtube ₍₂₎ D ^W -	Microtube Ctl.	Interpretation
+	+	0	D positive
0	0	0	D negative
0	+	0	D weak or partial
+	0	0	

+ = positive
0 = negative

Notes:

- The results by themselves alone are not a diagnosis. They must be evaluated together with the patient's clinical information and other data.
- In the forward ABO group test the microtube Ctl. must be negative. If it is positive, invalidate the test. Repeat the test after washing the red blood cells with physiological saline solution and preparing a new suspension of the washed red blood cells. If the microtube Ctl. of the repetition is negative, the results of the test can be interpreted; if it is positive, invalidate the test.
- Some individuals with forward group B or O with low levels of anti-A may not react with the red blood cells reagents for reverse A₂ group test.
- Some individuals with A₂ or A₂B subgroups may present antibodies anti-A₁, where a blood-serum discrepancy can be seen.
- It is recommended to check/investigate any discrepant results obtained, especially in the event of obtaining a blood-serum discrepancy, before releasing the result.
- ABO/Rh system: in the event of +/- to 3+ reactions, presence of a weak antigen must be investigated.
- Negative reactions with the anti-D^W- reagent must be verified using other reagents and techniques, which may detect different variants of the D antigen.
- In the event of obtaining discrepant results in microtubes ₍₁₎D^W- and ₍₂₎D^W-, it must be interpreted as a weak or partial D antigen. It is recommended to analyse the expression of this antigen.
- The observation of complete or partial haemolysis (pinkish supernatant and/or gel column) in microtubes must be interpreted as a positive result, after verifying that it is not due to a problem of extraction and/or handling of the sample.
- Occasionally, red blood cell retention in the incubation chamber may occur with positive 4+ samples, without interfering in result reading.

QUALITY CONTROL

- It is recommendable to include known positive and negative controls in each series of tests.
- If an unexpected control value is obtained, a complete verification of the instrument, reagents and material used must be made.

PERFORMANCE CHARACTERISTICS

Diagnostic sensitivity and specificity:

ABO/Rh system:

The diagnostic sensitivity and specificity of the antibodies present in the DG Gel ABO/Rh (CR) cards^{5,6}, for the determination of the antigens of the ABO and Rh systems, have been studied in a representative number of positive and negative samples.

Antibody	No. of Samples	Sensitivity ^(a)	Specificity ^(b)
Anti-A	3041 ⁵	100%	100%
Anti-B	3046 ⁵	100%	100%
Anti-D (microtube ₍₁₎ D ^W -)	3054 ^{6,5}	99.6%	100%
Anti-D (microtube ₍₂₎ D ^W -)	1329 ^(6,6)	99.7%	100%

(a) Sensitivity: (number of true positive results / (number of true positive results + number of false negative results)) x 100

(b) Specificity: (number of true negative results / (number of true negative results + number of false positive results)) x 100

(c) These samples include 51 samples with weak expression of the D antigen. If only these samples are considered, the sensitivity of the anti-D (microtube ₍₁₎D^W-) reagent is 80.4%.

(d) These samples include 49 samples with weak expression of the D antigen. If only these samples are considered, the sensitivity of the anti-D (microtube ₍₂₎D^W-) reagent is 93.9%.

Reverse ABO group

The DG Gel ABO/Rh (CR) card owns performance characteristics suitable for the determination of the reverse group, backed up by a study⁷ where results obtained were similar to those obtained with other established products of equivalent intended use.

Precision:

The precision of the reagents present in the DG Gel ABO/Rh (CR) card has been determined^{8,9}, including repeatability, inter-lot reproducibility and intra-laboratory reproducibility tests. No false positive or false negative results were obtained, and differences between agglutination intensities in positive samples were 1 agglutination grade or less in all assays.

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Sample: red blood cells

- Red blood cells from individuals with A or B variants may present a weak expression of the antigens. Antigen expression may be weakened in red blood cells of people with leukaemia or other malignant diseases¹⁰.

- Abnormal concentrations of serum proteins, the presence of macromolecular solutions in the serum/plasma or the presence of Wharton's jelly in cord blood samples may cause the non-specific agglutination of the red blood cells. It is recommended to wash the red blood cells before performing the test¹⁰.
- Transfused patients or those subjected to bone marrow transplant may present images of double population¹⁰.
- Patients with high-potency autoantibodies may coat the red blood cells completely, causing spontaneous agglutination¹⁰.

Sample: serum/plasma

- If plasma is used, complement-dependent haemolytic reactions may not be detected.
- If poor anti-coagulated plasma or partially coagulated serum is used, fibrin residues may trap non-agglutinated red blood cells at the top of the gel, appearing as a pinkish or reddish layer, but the negative reaction can be interpreted as such.
- It is recommended to re clot the serum at 37 °C during 10 minutes, centrifuge and repeat the test¹⁰.
- Newborns up to the age of 4-6 months, elderly persons, patients with immunodeficiency or very diluted antibodies due to plasma exchange procedures¹⁰ may have low or inexistent levels of insoagglutinins.

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- Do not use the card if microbiological contamination or change in microtube colour is observed.
- If due to improper transport or storage, dispersed drops at the top of the microtube, are observed, it is recommended to centrifuge the cards before use. If drops do not descend, do not use the card.
- If trapped bubbles, any microtube without supernatant, a decrease in the gel volume or cracked gel are observed, do not use the card.

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- The reagents of the DG Gel ABO/Rh (CR) card of human monoclonal origin are manufactured using a nonreactive material for the HB_s antigen nor for anti-HIV and anti-HCV antibodies, when tested with authorised reagents. Nevertheless, there is no known procedure to ensure that products of human origin will not transmit hepatitis or AIDS. Human blood products and samples must be handled as if they were potentially capable of transmitting infectious diseases.
- The product must only be used by qualified personnel.
- A red blood suspension of a concentration other than that indicated may cause false positive or negative results.
- The use of diluents other than DG Gel Sol for the red blood cell suspension may modify the reaction.
- The addition of volumes other than those indicated in the method may modify the reaction.
- Do not use the card beyond the expiry date.
- Once used, the product must be disposed of in special containers for biological waste.
- If you have any doubts or need further information on the use of this product, consult the authorized distributor in your country.

BIBLIOGRAPHY

- Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
- Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
- NCCLS H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
- NCCLS H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.
- External evaluation Diagnostic Grifols, S.A. (REGD-0002998), November 2007.
- External evaluation and internal study Diagnostic Grifols, S.A. (REGD-0012469), January 2013.
- Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.
- Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
- Internal Study Diagnostic Grifols, S.A. (DG-RI013/07-3), Mayo 2007.
- Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005.

PRESENTATION

210378 DG Gel ABO/Rh (CR) 50 cards Profile:A/B/₍₁₎D^W-/₍₂₎D^W-/Ctl./N/N/N

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LITERATŪRA

1. Mollison PL, Engelfriet CP, Contreras M. Blood Transfusion in Clinical Medicine, 9th edition, Blackwell Scientific Publications, Oxford, 1993.
2. Lapierre Y, et al. The gel test: a new way to detect red cells antigen-antibody reactions. Transfusion, 30: 109-113, 1990.
3. NCCLS H3-A6: Procedures for the collection of diagnostic blood specimens by venipuncture; Approved Standard, 6th edition, 2007.
4. NCCLS H18-A3: Procedures for the handling and processing of blood specimens; Approved Guideline, 3rd edition, 2004.
5. External evaluation Diagnostic Grifols, S.A. (REGD-0002998), November 2007.
6. External evaluation and internal study Diagnostic Grifols, S.A. (REGD-0012469), January 2013.
7. Internal study Diagnostic Grifols, S.A. (AT-10/03), October 2003.

8. Internal study Diagnostic Grifols, S.A. (DR-I03/03), July 2003.
9. Internal Study Diagnostic Grifols, S.A. (DG-RI013/07-3), Mayo 2007.
10. Technical Manual, 15th edition, American Association of Blood Banks, Bethesda, 2005.

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Extended IV Control

English

Quality control for use with DG Gel technique, MDmulticard or conventional tube technique

SUMMARY AND PRINCIPLE

Blood Bank Quality Assurance includes a check of all reagents routinely employed in immunohematologic procedure on the day of use to ensure that the reagents are functioning properly and to verify techniques.

Extended IV Control is intended to allow regular control of materials, work procedures and instrument procedures for (i) the determination of the ABO, Rhesus and K antigens, (ii) the determination of the appropriate ABO blood group antibodies, (iii) the detection of irregular antibodies in both antiglobulin and enzyme technique and (iv) ABO compatibility testing. The product is suitable for use with DG Gel technique, MDmulticard and conventional tube technique.

REAGENT

Warning: Do not use for grouping donors or recipients. Use only in the laboratory for quality assurance of blood bank reagents.

Extended IV Control suspensions have been prepared from red blood cells collected from blood donors. Each tube may be prepared from a single donor or a pool of donors. The suspensions are composed of human erythrocytes (haematocrit $15 \pm 2\%$) in a buffered isotonic nutrient medium with preservatives (0.03% (w/v) neomycin and 0.05% (w/v) chloramphenicol) and added antibodies as appropriate which contain $<0.01\%$ (w/v) sodium azide as a preservative. Frozen/thawed red blood cells may have been used in this product.

Each Extended IV Control kit contains 4 tubes of 6 ml:

Tube	ABO	Rhesus	Other Typing	Antibodies
QCX1	A ₂ B	CcD.Ee	K pos	-
QCX2	A	CC ^o D.ee	Fy ^a neg	Anti-B, Anti-Fy ^a
QCX3	B	ccD.EE	-	Anti-A
QCX4	O	ccddee	K neg	Anti-A, Anti-B, Anti-D (~ 0.05 IU/ml)

Antibody	Clones used	Origin
Anti-A monoclonal	Birma-1	Murine
Anti-B monoclonal	LB2	Murine
Anti-D polyclonal	-	Human
Anti-Fy ^a monoclonal	P3TIM	Murine/human

Direct antiglobulin tests are negative on all red blood cell suspensions of Extended IV Control.

For professional use only. Ready-to-use.

Caution: All blood products should be treated as potentially infectious. All human source material used in the manufacture of this product were tested and found negative for HBsAg and antibodies to HCV, HIV-1/ HIV-2. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents, appropriate care should be taken in the use of this product. The absence of murine viruses has not been established. Once used, the product must be disposed in special containers for biological waste.

Warning: This product contains sodium azide. Sodium azide may react with lead and copper plumbing to form highly explosive metal azides. If discarded into sink, flush with a large volume of water to prevent azide build-up.

STABILITY

Do not expose to extreme temperatures. Store at 2-8°C when not in use. Turbidity may indicate reagent deterioration or microbial contamination; discard if present. Do not use beyond expiration date. If stored appropriately at 2-8°C the product is stable after the first opening until the indicated expiration date.

It is normal to observe a slight change of colour in the supernatant.

PROCEDURE

Reagent provided

Extended IV Control, 4x6 ml, cat.no. 213286

Materials required but not provided

Consider the products and equipment related to the appropriate technique intended to be used.

The reagents must be brought to room temperature (18-25°C) prior testing.

PROCEDURE OUTLINE

Process the Extended IV Control in the same way as patients' samples.

Manual method, DG Gel technique

Use the suspensions of Extended IV Control as control samples with DG Gel reagents of Diagnostic Grifols S.A. and/or 0.8% Reagent Red Blood Cells of Medion Grifols Diagnostics AG. or Diagnostic Grifols S.A.. Follow the instructions for use of the DG Gel card and 0.8% Reagent Red Blood Cells employed.

Manual method, MDmulticard

Use the suspensions of Extended IV Control as control samples with MDmulticard reagents of Medion Grifols Diagnostics AG. Follow the instructions for use of the MDmulticard reagent employed.

Manual method, conventional tube technique

Use the suspensions of Extended IV Control as control samples with tube technique reagents of Medion Grifols Diagnostics AG. Follow the instructions for use of the tube technique reagent and the Reagent Red Blood Cells employed.

Automated method

The suspensions of Extended IV Control can be used in Diagnostic Grifols S.A. instruments for DG Gel. Follow the user's manual of the automatic instrument employed. When working with automated systems, it is recommended to perform a quality control after decontamination processes. Extended IV Control is valid for this purpose.

QUALITY CONTROL

To ensure proper centrifugation, each individual centrifuge should be calibrated for the specific test procedure being performed.

RESULTS

Interpretation

If the test results differ from the awaited results listed in section "REAGENT", check immediately the working procedure, test method, automatic instrument/equipment and other related materials used. Repeat the test after having applied the appropriate corrections.

LIMITATIONS OF PROCEDURE

- 1. Improper techniques may invalidate the results obtained with Extended IV Control.
- 2. False positive or false negative results can occur due to contamination of test material, improper reaction temperature and/or time, improper centrifugation, improper storage conditions of material and omission of test reagents.
- 3. Do not use if evident changes of colour and/or obvious haemolysis is observed.
- 4. If used with Diagnostic Grifols S.A. automatic instruments for DG Gel, the expected ABO interpretation for tube QCX1 is AB instead of A₂B.
- 5. Cells may lose antigen strength during the shelf life of the product giving weaker positive reactions than fresh cells.
- 6. Some interaction between the Rh negative (D) red blood cells and the Anti-D reagent (clones: P3x290, P3x35, P3x61 and P3x21223B10) may be observed as a result of reagents' inherent properties.
- 7. Take under consideration all limitation statements in the instructions for use of the used DG Gel cards, the 0.8% Reagent Red Blood Cells, the conventional tube reagents, the MDmulticard reagents and the automatic instruments.

SPECIFIC PERFORMANCE CHARACTERISTICS

The performance of the product was validated according to the limited field of application (refer to intended use, section “SUMMARY AND PRINCIPLE”). The requirements of the Common Technical Specifications (CTS) 2009/886/EC are fulfilled as far as applicable.

According to the repeatability and reproducibility testing, the total number of tests (n) is displayed below and the sensitivity as well as the specificity was calculated for each technique for which an agreement of 100% for sensitivity and specificity was obtained.

Technique	Antigen Parameters: Control for Antibody Reagents									
	n determination									
	A	B	AB	D	C	E	c	e	C"	K
Gel	320	320	160	480	320	320	240	240	80	80
MDmulticard	160	160	80	320	160	160	160	160	80	80
Tube	80	80	80	80	80	80	80	80	80	40

Technique	Antibody Parameters: Control for Reverse Grouping Cells			
	n determination			
	Anti-A	Anti-B	Anti-AB	No Isoagglutinin
Gel	80	80	80	80
Tube	80	80	80	80

Technique	Antibody Parameters: Control for Antibody Screening Cells and Identification Cells		
	n determination		
	Anti-D	Anti-Fy ^a	No Antibody
Gel	120	120	240
Tube	60	60	120

Sensitivity: probability to obtain a positive result in presence of a positive sample.
Specificity: probability to obtain a negative result in presence of a negative sample.

WARRANTY

This product is warranted to perform as described in its labelling and in the product literature and Medion Grifols Diagnostics AG disclaims any implied warranty of merchantability or fitness for any other purpose, and in no event shall Medion Grifols Diagnostics AG be liable for any consequential damages arising out of the aforesaid express warranty.

IVD

CE0123

Used symbols

IVD

"In Vitro" Diagnostic Medical Device

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Serigrup Diana A₁/B Serigrup Diana 4 Serigrup Diana A₂

Reagent Red Blood Cells for reverse grouping Instructions for use. For *in vitro* diagnostic use.

INTENDED USE

Serigrup Diana A₁/B, Serigrup Diana 4 and Serigrup Diana A₂ Reagent Red Blood Cells are used for the determination of the reverse ABO group in gel technique.

SUMMARY AND EXPLANATION

The ABO system is defined by the presence or absence of the A and/or B antigens in the red blood cells and by the presence of antibodies in the serum, corresponding to the antigen or antigens missing in the red blood cells.

The determination of the serum or reverse ABO group, is based on the presence or absence of anti-A and/or anti-B isagglutinins in the patient's serum, using red blood cells from a known blood group: A₁, A₂, B and O. Anti-A and anti-B isagglutinins are natural antibodies produced by individuals lacking the respective A and B antigens.

Some individuals with subgroups A₁ or A₂ may present anti-A antibodies, which interfere with the procedure for determining the reverse ABO group using A₁ Reagent Red Blood Cells. The Serigrup Diana A₂ reagent is used to investigate such discrepancies.

PRINCIPLE OF THE TEST

Anti-A and anti-B from human serum bind to red blood cells possessing the corresponding antigens, resulting in direct agglutination. Serigrup Diana Reagent Red Blood Cells are utilized in the gel technique to detect antibodies to human blood groups A and B.

REAGENTS

Observable indications

Inspect the conditions of the Reagent Red Blood Cells before use (see Warnings and Precautions). Do not use the product, if you detect clausiness and/or large-scale hemolysis, alterations or changes in color, microbiological contamination or other artifacts. Do not use the product, if the package, vials or caps are damaged or if contents have leaked.

Reagent Red Blood Cells

All Serigrup Diana Reagent Red Blood Cells consist of a 0.8% suspension ready-to-use of human red blood cells of A₁, B, A₂ and O blood group produced by a single donor for each vial. The cells come in a buffered solution with preservatives (neomycin 0.010% (w/v) and chloramphenicol 0.017% (w/v)). Each glass vial of Serigrup Diana contains 10 mL of human red blood cells. Vials come with a blue cap for Serigrup Diana A₁, a black cap for Serigrup Diana A₂, a yellow cap for Serigrup Diana B and a white cap for Serigrup Diana O.

Material required but not provided

Please refer to the Instructions for Use of Grifols gel cards.

STORAGE AND STABILITY

- Do not use beyond the expiration date.
- Store upright (as indicated by the two arrows on the outer packaging) at 2 - 8 °C.
- Do not freeze.

Notes:

- If you identify incorrect temperature conditions during storage or shipment do not use the product.
- To avoid contamination, close the vials with their caps when they are not being used. Ensure that the caps of the Reagent Red Blood Cells vials have not been swapped.
- Once a vial has been used, keep it at the indicated storage temperature.

WARNINGS AND PRECAUTIONS

- The results do not constitute a clinical diagnosis by themselves. Evaluate the results together with the patient's clinical information and other data.
- The product should only be used by qualified personnel.
- Any modifications of the test procedures described in these Instructions for Use require validation by the user.
- The use of volumes other than those indicated in the method, may modify the reaction and lead to incorrect test results, i.e. false-positive or false-negative results.
- Use disposable plastic pipette tips in order to avoid cross-contamination with samples and other reagents.
- Do not mix residues from different vials.
- Serigrup Diana Reagent Red Blood Cells are of human origin and are manufactured using materials that have been tested and found non-reactive for the HBs antigen, and for anti-HIV and anti-HCV antibodies. However, there is no known procedure to ensure that products of human origin will not transmit infectious agents. Human blood products and samples should be handled as if they were potentially capable of transmitting infectious diseases.

English

- Once used, dispose the vial and samples in containers for biological waste, according to local, state and national regulations.
- If you have any questions or need further information on the use of this product, consult the authorized distributor in your country.

SPECIMEN COLLECTION AND PREPARATION

For the determination of the ABO reverse group use serum or plasma. For optimum test results, samples should be stored at 2 - 8 °C for no longer than 72 hours prior to testing. Frozen samples stored up to 5 years at -20 °C or colder may be used after thawing. The collection, separation and handling of the blood should be performed by qualified technical personnel according to current standards¹, and following the instructions of the manufacturer of the material used for collecting the sample. If necessary, samples tested with Serigrup Diana A₁ and Serigrup Diana B, stored at 2 - 8 °C, can be used up to 7 days after collection.

PROCEDURE

- Allow the Reagent Red Blood Cells and the samples to reach room temperature (18 - 25 °C).
- Inspect the condition of the reagents before use.
- Homogenize the red blood cells of Serigrup Diana gently by inversion.
- Follow the instructions for use of the Grifols gel cards used.**

QUALITY CONTROL

Include positive and negative controls with testing on each day of use. The inclusion of weak positive controls is recommended. If an unexpected control result is obtained, an inspection of the instrument, reagents and materials used should be carried out.

RESULTS

Interpretation of results

Determine the reverse grouping in accordance with the following table:

Microtube N Serigrup Diana A ₁	Microtube N Serigrup Diana B	Microtube N Serigrup Diana B	Microtube N Serigrup Diana O	Reverse ABO group
+	+	+	0	O
0	0	+	0	A
+	+	0	0	B
0	0	0	0	AB
+	0	+	0	A ₂ with anti-A ₁
+	0	0	0	A,B with anti-A ₁

+ = positive
0 = negative

Notes:

- If a discrepancy is detected between the forward ABO blood grouping and the reverse grouping, this should be investigated before releasing the result.
- Serigrup Diana A₁ should always be used in combination with Serigrup Diana A₂/B.
- Use the O group cells for Serigrup Diana 4 for control purposes only. Do not use for screening of unexpected antibodies.
- The observation of complete or partial hemolysis (pinkish supernatant and/or gel column) in microtubes should be interpreted as a positive result, after verifying that it is not due to a sample quality issue (caused by incorrect extraction and/or handling) or poor condition of the red blood cells.

LIMITATIONS OF THE PROCEDURE

- Grossly hemolyzed, cloudy or contaminated samples or samples with presence of a clot, may cause false positive or false negative results.
- If plasma is used, complement-dependent hemolytic reactions may not be detected.
- If poorly anti-coagulated plasma or incompletely clotted serum is used, fibrin residues may trap non-agglutinated red blood cells at the top of the gel, appearing as a pinkish or reddish layer. Although the results could be correctly interpreted, in a negative reaction the false appearance of a Double Population could lead to a misinterpretation. In case of incompletely clotted serum samples, it is recommended to re-clot the serum and repeat the test¹.
- The formation of "rouleaux", caused by an excess of proteins in the serum, the presence of abnormal proteins, drugs, plasma expanders, etc., may cause false positive reactions¹.

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