C O N F I D E N T I A L

11, 12 pozicijos

SoFilmPlus[®]


PRODUCT LAUNCH INFORMATION


Somed

1. Product concept

1.1. Appearance 11 poz

12 poz

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Construction

- 1) Polyurethane carrier
- 2) Adhesive with silver ions
- 3) Release liner white with ANTIMICROBIAL print


1.2. Surgical Site Infections

Surgical-site infections are an important health care problem and may result in significantly more post-surgical complications, lengthier hospital stays and increased mortality.


Incision films are a surgeon's best allies for reducing microbial proliferation and minimising the incidence of trans- and post-operative infections, which have always been a serious problem for the patient's health and an important cost-driver for the health care industry.

1.9. Sales Points and Major properties

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

You have a sterile area, the only place where microorganisms can still be introduced is from the skin in the form of sweat. the film keeps the area sterile between the moment of disinfection and incision and that is especially the moment that it is necessary, so that if you cut you do not cut a micro-organism into the skin and that remains sterile outside the wound during the operation.

The moment the incision film comes into contact with sweat, a bridge forms with the zeolite structure and the bacteria that are introduced there. The silver is bound to the zeolite as a positive ion and comes into contact with water (is exchanged with another positive ion) so that the silver that is now in water can start with the termination of bacteria

By osmosis (such as two cups of water a cold and a warm connect with each other than over time are equal), it is possible for the silver ion to enter through the cell membrane. Outside, silver concentration is higher than inside, and a balance is needed. It is well known that silver has a killing effect on bacteria and once it is inside the process can begin.

2. Product specifications

2.1. Product description and application

SoFilm and SoFilmPlus are sterile incision films.

SoFilm and SoFilmPlus are intended to be applied on the patient's skin for covering and protecting during surgical interventions. The incision is made through the film. SoFilmPlus is enhanced with active silver ions for a continuous antimicrobial effect.

SoFilm and SoFilmPlus are consisting of:

- a transparent polyurethane film (1);
- single adhesive coated (2);
- covered with a 120 grams release liner (3) and
- a 5cm wide transparent finger lift.

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SoFilm and SoFilmPlus are individually wrapped, pouched and sterilised by EO.

SoFilm

UV cured, acrylic based, pressure- and temperature-sensitive hotmelt adhesive.

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MEDIBERG print on the release liner

SoFilmPlus

UV cured, acrylic based, pressure- and temperature-sensitive hotmelt adhesive with antimicrobial effect

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ANTIMICROBIAL print on the release liner

2.4. Antimicrobial efficacy

Disinfectant on the human skin


As part of the infection-prevention protocol, disinfectant is applied to the patient's skin in the surgical area before the intervention. However, bacteria spontaneously recolonise the disinfected area of the skin 30 to 60 minutes after disinfection.

sample	before desinfection	† 0 minutes	† 60 minutes
1) Cutasept F (72% Isopropanol + Benzalconium chloride)			
2) AddedPharma Hand disinfectant (Chlorhexidingluconate 0.5% + Alcohol 70%)			

Test laboratory: Quality Labs BT GmbH, 29-Sept-2016

Type of test: Antimicrobial Efficacy Testing of Disinfectants on Material surfaces

A study was conducted on the test product SoFilmPlus in order to determine its antibacterial effectiveness compared to regular incision film.


Test laboratory: Eurofins Biolab s.r.l., 01-Feb-2016

Type of test: ASTM E-2180-07: Determining the Activity of Incorporated Antimicrobial Agent(s) In Polymeric or Hydrophobic Materials

2.5. Dimensions

SoFilm and SoFilmPlus are available in the following dimensions.

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12 poz

2.7. Packaging

Primary packaging (pouch, sterile barrier)

Film dimension	Pouch dimensions	Content per pouch	Weight per pouch
25cm x 20cm	32cm x 10cm	1 piece	17 grams
50cm x 35cm	45cm x 10cm	1 piece	44 grams
50cm x 45cm	55cm x 10cm	1 piece	56 grams
65cm x 45cm	55cm x 10cm	1 piece	68 grams
65cm x 60cm	75cm x 10cm	1 piece	89 grams
95cm x 65cm	75cm x 10cm	1 piece	132 grams
65cm x 85cm	99cm x 10cm	1 piece	129 grams

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EC DECLARATION OF CONFORMITY FULL QUALITY ASSURANCE SYSTEM

according to Annex II excluding (4) of EC Directive
 93/42/EEC and subsequent modifications and integrations.
*(transposed in Italy by the D.Lgs. n. 46 issued on
 24.02.1997 and subsequent modifications and integrations)*

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*The Istituto Superiore di Sanità,
 Notified Body 0373, certifies that
 the total quality assurance system
 enforced by*

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*is in compliance with the applicable
 requirements of Council Directive 93/42/EEC
 and subsequent modifications and integrations.*

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of which this technical sheet is an integral part, refers only to the following product(s) that are subject to surveillance:

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