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CE<sub>0123</sub>

  
**EDAN**

	elite V8	$\geq 2.5$ h	ECG/TEMP module connected, sidestream CO <sub>2</sub> connected, recording at interval of 10 minutes, brightness set to "1")
Charge Time	elite V5/ elite V6	$\leq 10$ h	The monitor is on or in standby mode, 100% charge.
		$\leq 9$ h	The monitor is on or in standby mode, 90% charge.
		$\leq 6$ h	The monitor is off, 100% charge.
		$\leq 5.4$ h	The monitor is off, 90% charge.
	elite V8	$\leq 6$ h	The monitor is on or in standby mode, 100% charge.
		$\leq 5.4$ h	The monitor is on or in standby mode, 90% charge.
Alarm	Low battery alarm is provided.		

## A.6 Display

Model	Display
elite V5	Display screen: 12.1-inch color TFT screen, touch screen is configurable Resolution: 800 × 600 A maximum of 13 waveforms can be displayed on the same screen.
elite V6	Display screen: 15-inch color TFT screen, touch screen is configurable Resolution: 1024 × 768 A maximum of 13 waveforms can be displayed on the same screen.
elite V8	Display screen: 17-inch color TFT screen, touch screen is configurable Resolution: 1280 × 1024 A maximum of 15 waveforms can be displayed on the same screen.


## A.7 Indicators

Power-On LED	Green
AC Power LED	Green
Battery LED	Yellow/green
Physiological Alarm LED	Red/yellow
Technical Alarm LED	Red/yellow/blue


## Chapter 2 Overview

### 2.1 Main Unit

#### Front View

1	Technical alarm indicator	When a technical alarm occurs, the indicator lights on or flashes with different frequencies and colors reflecting the alarm level.
2	Physiological alarm indicator	When a physiological alarm occurs, the indicator lights on or flashes with different frequencies and colors reflecting the alarm level.

### Patient Data at Your Fingertips

The elite V6 Modular Patient Monitor combines the world's leading monitoring technologies in a single platform that seamlessly connects bedside patients to the central monitoring station.

- LAN connection to EDAN MFM-CMS central monitoring station via cable network or built-in Wi-Fi
- Remote login capability via web access from PC, tablet or smartphone
- HL7/XML connectivity with the hospital information system

## elite V6 Modules

MODULE	DESCRIPTION	BENEFITS
M20	<ul style="list-style-type: none"> <li>30-lead ECG, RESP, NIBP, SpO2, 2-LMP</li> <li>Optional: 12-lead ECG, Nebulizer SpO2, Suspend</li> <li>12-lead ECG, Nebulizer SpO2, Suspend</li> <li>RESP, 2-AP, Respiration Substream CO2</li> </ul>	<ul style="list-style-type: none"> <li>Allows seamless monitoring between bedside and transport monitoring</li> </ul>
XM	<ul style="list-style-type: none"> <li>30-lead ECG, RESP, NIBP, SpO2, 2-TEMP</li> <li>Optional: 12-lead ECG, Nebulizer SpO2, Suspend</li> <li>RESP, 2-AP</li> </ul>	<ul style="list-style-type: none"> <li>A multi-parameter module</li> </ul>
VIBP	<ul style="list-style-type: none"> <li>Minimum 8-AP</li> <li>Optional: 12-lead ECG, 2-AP, 1-CP</li> </ul>	<ul style="list-style-type: none"> <li>IBP/Swamping</li> <li>IBP/CO</li> <li>CPP supported with ICP monitoring</li> </ul>
VCO	<ul style="list-style-type: none"> <li>Thermistor/Carotid Catheter Output</li> </ul>	<ul style="list-style-type: none"> <li>Gold standard carotid output methodology</li> </ul>
VICG	<ul style="list-style-type: none"> <li>Invasive Cardiology</li> </ul>	<ul style="list-style-type: none"> <li>Non-invasive carotid output methodology</li> <li>Non-invasive monitoring with real time measurements</li> </ul>
VCO2	<ul style="list-style-type: none"> <li>Respiration Mainstream CO2</li> <li>Respiration Substream CO2</li> <li>EDAN G7 Substream CO2</li> </ul>	<ul style="list-style-type: none"> <li>Low sampling rate at 50/min</li> <li>No water cap</li> <li>Supports water cap design</li> <li>ICARETM algorithm</li> </ul>
VIMV	<ul style="list-style-type: none"> <li>Respiration Respiration Mechanics</li> </ul>	<ul style="list-style-type: none"> <li>Helps detect pulmonary diseases and reduce ventilator associated pneumonia</li> <li>Provides on the support of respiratory disorders</li> </ul>
VAG	<ul style="list-style-type: none"> <li>Mainstream AG</li> <li>Minimum Substream AG/O2</li> <li>Digital Substream AG</li> </ul>	<ul style="list-style-type: none"> <li>Major maintenance technology</li> <li>No water cap</li> <li>Low sampling rate at 30/min</li> <li>Paramagnetic O2 sensor</li> <li>Water cap design</li> </ul>
VBS	<ul style="list-style-type: none"> <li>Bispectral Index</li> <li>BR-EEG</li> </ul>	<ul style="list-style-type: none"> <li>Open access to trending during surgery</li> <li>Reduction of waking</li> <li>Reduces anesthetic dose to speed up recovery and reduce time spent in PACU</li> </ul>

### About Edan


Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
- Ultrasound imaging
- In-Vitro Diagnostics
- Patient Monitoring
- Point-of-Care Testing
- Veterinary
- OB/GYN

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.

**EDAN/USA**

4204 Ashford Drive  
Suite B, San Diego, CA 92117  
1-800-457-4577

### Always Connected

A self-contained transport monitor, the M20, provides seamless data connectivity from bedside monitoring to state-of-the-art transport monitoring. On plug-and-play, it switches role easily between a multi-parameter module and a monitor, connecting data together among different applications, and delivers uninterrupted monitoring during patient transfer.

To keep track of a full-time patient data, the M20 not only carries data between each elite V monitor, but also interfaces with MFM-CMS central station via Wi-Fi.

# iM20

## Patient Monitor

5"

Serving as both multi-parameter module of elite V series monitors and transport monitor, iM20 provides a comprehensive and seamless monitoring solution both inside and outside hospital, connecting ambulatory cares to intensive cares.

### Parameter Configurations

- iSLAP™ ECG algorithm optimized for arrhythmia detection, tachycardia detection, and HR measurement
- SLMR™ 12-lead diagnosis algorithm with 208 kinds of ECG findings in support of world ECG configurations
- iMAP™ SpO2 algorithm with outstanding motion resistance and low perfusion regulatory performance
- Perfusion index (PI) measurement as a blood perfusion indicator
- Optional Nihilor SpO2 with SubSensitivity™ technology
- iCUP™ NIBP algorithm optimized for cardiac patients, hypertensive patients, and neonatal patients
- Optional Suntech NIBP with SunTech Advantage™ FIM technology
- Optional 2-channel invasive blood pressure with waveform averaging and pulse detection
- Optional EIM module with Respiration Lark™ technology made
- 2-channel temperature monitoring

### Built-In Temporary Memory

- 150 hours Trend Review
- 1200 sets NIBP Measurement Review
- 200 sets Alarm Review
- 200 sets Arrhythmia Review
- 50 sets 12-lead Diagnostic Review

No. 1 16.2 16.3



2.1

1.1

2.3

2.4

2.2

### Ambulatory Designs

Designed as a robust transport monitor, iM20 is a solid companion in the harshest conditions.

- With IP44 water and dust protection, iM20 is outstanding for field medical practice, from regular emergency care to rainy weather, to disaster control with dusty or rocky fields.
- Designed with 1 million shock resistance, anti-torqued hits and impacts, shock not interrupt the monitoring, nor does the battery ride on rough road or busy area.

Interfaces and others			
USB Port	1 (OTG, USB1.0/2.0 protocol)		
ECG			
Lead Mode	3-Lead: I, II, III		
	5-Lead: I, II, III, aVR, aVL, aVF, V		
	12-Lead: I, II, III, aVR, aVL, aVF, V1, V2, V3, V4, V5, V6		
Lead naming style	AHA, IEC		
Display Sensitivity (Gain Selection)	1.25 mm/mV (×0.125), 2.5 mm/mV (×0.25), 5 mm/mV (×0.5), 10 mm/mV (×1), 20 mm/mV (×2), 40 mm/mV (×4), AUTO gain		
Sweep speed	6.25 mm/s, 12.5 mm/s, 25 mm/s		
Bandwidth (-3 dB)	Diagnosis: 0.05 Hz to 150 Hz		
	Monitor: 0.5 Hz to 40 Hz		
	Surgery: 1 Hz to 20 Hz		
CMRR	Diagnosis: >95 dB		
	Monitor: >105 dB		
	Surgery: >105 dB		
Notch	In diagnosis, monitor and surgery modes: 50 Hz/60 Hz (Notch filter can be turned on or off manually)		
Recovery time after defibrillation	<5 s		
ESU Protection	Cut mode: 300 W		
	Coagulation mode: 100 W		
	Restore time: ≤10 s		
3.3 Pace pulse detecting lead	one among I, II, III, AVR, AVL, AVF, V1, V2, V3, V4, V5, V6		
Heart Rate			
Range	Adult: 15 bpm to 300 bpm		
	Ped: 15 bpm to 350 bpm		
Accuracy	±1% or ±1 bpm, whichever is greater		
Resolution	1 bpm		
PVC			
Range	Adult: 0 to 300 PVCs/ min		
	Ped/Neo: 0 to 350 PVCs/ min		
Resolution	1 PVCs/min		
ST value			
Range	-2.0 mV to +2.0 mV		
Accuracy	±0.02 mV or 10% (-0.8 mV to +0.8 mV), whichever is greater. Beyond this range: not specified.		
Resolution	0.01 mV		
Arrhythmia analyses			
3.4 ASYSTOLE	Sustain VT	V-Fib/V-Tach	Extreme Tachy
	Extreme Brady	V-Tach	Vent Brady
			Tachy

4

Brady	Wide QRS Tachy	Non-Sustain VT	Afib
Vent Rhythm	Acc. Vent Rhythm	Pause	Pauses/min High
PVCs High	R on T	PVC Bigeminy	PVC Trigeminy
Pacer not Pacing	Pacer not Capture	Missed Beat	VEB
PVC	Couplet	Run PVCs	Multiform PVCs
IPVC	Irr Rhythm	PAC Bigeminy	PAC Trigeminy
Low Voltage(Limb)			

**12-lead ECG synchronization analysis**

Average parameters of heart beat

Heart rate (bpm)

Time limit of P wave (ms)

PR interval (ms)

QRS interval (ms)

QT/QTc (ms)

P-QRS-T AXIS

**RESP**

Method Impedance between RA-LL, RA-LA

Measurement lead Options are lead I and II

Measuring Range  
 Adult 0 rpm to 120 rpm  
 Ped/Neo 0 rpm to 150 rpm *4.2*

Resolution 1 rpm

Accuracy  
 Adult 6 rpm to 120 rpm: ±2 rpm  
 0 rpm to 5 rpm: not specified  
 Ped/Neo 6 rpm to 150 rpm: ±2 rpm  
 0 rpm to 5 rpm: not specified

Gain Selection ×0.25, ×0.5, ×1, ×2, ×3, ×4, ×5

Sweep 6.25 mm/s, 12.5 mm/s, 25.0 mm/s

Apnea Alarm Time 10 s, 15 s, 20 s, 25 s, 30 s, 35 s, 40 s *4.2*

**NIBP**

**EDAN Module**

Method Oscillometry

Mode Manual, Auto, Continuous

Measuring Interval in Auto Mode 1/2/3/4/5/10/15/30/60/90/120/180/240/360/480 min *6.2*

Continuous 5 min, interval is 5 s

Measuring Type SYS, DIA, MAP, PR

Measuring Range  
 Adult Mode  
 SYS: 40 mmHg to 270 mmHg  
 DIA: 10 mmHg to 215 mmHg  
 MAP: 20 mmHg to 235 mmHg *6.1*

Message	Cause	Alarm level
<b>P1 MAP Low</b>	P1 MAP measuring value is below lower alarm limit.	User-selectable
<b>P2 SYS High</b>	P2 SYS measuring value is above upper alarm limit.	User-selectable
<b>P2 SYS Low</b>	P2 SYS measuring value is below lower alarm limit.	User-selectable
<b>P2 DIA High</b>	P2 DIA measuring value is above upper alarm limit.	User-selectable
<b>P2 DIA Low</b>	P2 DIA measuring value is below lower alarm limit.	User-selectable
<b>P2 MAP High</b>	P2 MAP measuring value is above upper alarm limit.	User-selectable
<b>P2 MAP Low</b>	P2 MAP measuring value is below lower alarm limit.	User-selectable
<b>CO<sub>2</sub></b>		
<b>EtCO<sub>2</sub> High</b>	EtCO <sub>2</sub> measuring value is above upper alarm limit.	User-selectable
<b>EtCO<sub>2</sub> Low</b>	EtCO <sub>2</sub> measuring value is below lower alarm limit.	User-selectable
<b>FiCO<sub>2</sub> High</b>	FiCO <sub>2</sub> measuring value is above alarm limits.	User-selectable
<b>CO<sub>2</sub> APNEA</b>	In the set apnea alarm delay time interval, no RESP can be detected using CO <sub>2</sub> module.	High
<b>AwRR High</b>	AwRR measuring value is above upper alarm limit.	User-selectable
<b>AwRR Low</b>	AwRR measuring value is below lower alarm limit.	User-selectable
<b>CI High</b>	CI measuring value is above upper alarm limit.	User-selectable
<b>CI Low</b>	CI measuring value is below lower alarm limit.	User-selectable

## 6.2 Technical Alarm Information

### NOTE:

The ECG alarm information listed in the below table describes the electrode names in America. For the corresponding electrode names in Europe, please refer to the Section *Installing Electrodes*.

Message	Cause	Alarm Level	Action Taken
<b>ECG</b>			
<b>ECG Lead Off</b>	1) The drive electrode or more than one ECG limb electrode falls off the skin;	Low	Make sure that all electrodes, leads and patient cables are properly connected.
	2) ECG cables fall off the monitor.		
<b>ECG LL Lead Off</b>	ECG electrode LL falls off	Low	Make sure that all electrodes, leads and patient cables are properly connected.
	the skin or the ECG cable LL falls off the monitor.		

Message	Cause	Alarm Level	Action Taken
<b>ECG LA Lead Off</b>	ECG electrode LA falls off the skin or the ECG cable LA falls off the monitor.	Low	Make sure that all electrodes, leads and patient cables are properly connected.
<b>ECG RA Lead Off</b>	ECG electrode RA falls off the skin or the ECG cable RA falls off the monitor.	Low	
<b>ECG RL Lead Off</b>	When electrode type is AUTO, ECG electrode RL falls off the skin or the ECG cable RL falls off the monitor, 5/6/10 electrodes switches to 3 electrodes;	Low	
<b>ECG V Lead Off</b>	ECG electrode V falls off the skin or the ECG cable V falls off the monitor.	Low	
<b>ECG V1 Lead Off</b>	ECG electrode V1 falls off the skin or the ECG cable V1 falls off.	Low	Make sure that all electrodes, leads and patient cables are properly connected.
<b>ECG V2 Lead Off</b>	ECG electrode V2 falls off the skin or the ECG cable V2 falls off.	Low	Make sure that all electrodes, leads and patient cables are properly connected.
<b>ECG V3 Lead Off</b>	ECG electrode V3 falls off the skin or the ECG cable V3 falls off.	Low	Make sure that all electrodes, leads and patient cables are properly connected.
<b>ECG V4 Lead Off</b>	ECG electrode V4 falls off the skin or the ECG cable V4 falls off.	Low	Make sure that all electrodes, leads and patient cables are properly connected.
<b>ECG V5 Lead Off</b>	ECG electrode V5 falls off the skin or the ECG cable V5 falls off.	Low	Make sure that all electrodes, leads and patient cables are properly connected.
<b>ECG V6 Lead Off</b>	ECG electrode V6 falls off the skin or the ECG cable V6 falls off.	Low	Make sure that all electrodes, leads and patient cables are properly connected.

	Adult/Pediatric (Sports Mode)	120 s
	Neonate	75 s
Overpressure protection	Adult/Pediatric	< 300 mmHg
	Neonate	< 150 mmHg
<b>PR</b>		
Measuring range	30 bpm to 220 bpm	
Accuracy	±3 bpm or ±2%, whichever is greater	
<b>SpO<sub>2</sub></b>		
<b>EDAN Module</b>		
Measuring Range	0% to 100% 5.1	
Resolution	1%	
Data update period	1 s	
Accuracy	Adult/Pediatric	±2% (70% to 100% SpO <sub>2</sub> ) Undefined (0% to 69% SpO <sub>2</sub> )
	Neonatal	±3% (70% to 100% SpO <sub>2</sub> ) Undefined (0% to 69% SpO <sub>2</sub> )
<b>PI (Perfusion Index)</b>		
Measuring Range	0-10	
Resolution	1	
<b>Pulse Rate</b>		
Measuring Range	25 bpm to 300 bpm 5.2	
Resolution	1 bpm	
Accuracy	±2 bpm	
<b>Nellcor Module</b>		
Measuring Range	1% to 100%	
Resolution	1%	
Data Update Period	1 s	
Accuracy	DS-100A, OXI-A/N(Adult), D-YS (Adult and Pediatric), OXI-P/I (Pediatric)	±3% (70% to 100% SpO <sub>2</sub> )
	MAX-A, MAX-AL, MAX-N, MAX-P, MAX-I, MAX-FAST (Adult and Pediatric)	±2% (70% ~ 100% SpO <sub>2</sub> )
	MAX-A, MAX-AL, MAX-N, MAX-P, MAX-I, MAX-FAST (Adult and Pediatric)	±3% (60% ~ 80% SpO <sub>2</sub> )
	If sensor is used for neonate as recommended, the accuracy will be larger than adult by ±1.	
<b>Pulse Rate</b>		
Measuring Range	20 bpm to 300 bpm	
Resolution	1 bpm	
Accuracy	±3 bpm (20 bpm to 250 bpm)	
<b>TEMP</b>		

- an arterial oxygen simulator (also an electronic pulse simulator).
- 3 During monitoring, if the monitor's reading differs significantly from the patient's physiological condition, it indicates that the signal may be disturbed, resulting in an inaccurate reading. In this case, the artifact can disguise as a similar reading, causing the monitor to fail to send an alarm. In order to ensure reliable monitoring, it is necessary to regularly check whether the sensor is wearing properly and the signal quality is good.

## 11.6 SpO<sub>2</sub> Alarm Delay

There is a delay between a physiological event at the measurement site and the corresponding alarm at the monitor. This delay has two components:

1. The time between the occurrence of the physiological event and when this event is represented by the displayed numerical values. This delay depends on the algorithmic processing time and the sensitivity. The lower the sensitivity configured, the longer the time needed until the numerical values reflect the physiological event.
2. The time between the displayed numerical values exceeding an alarm limit and the alarm indication on the monitor. This delay is the combination of the configured alarm delay time plus the general system delay time.

## 11.7 Perfusion Index (PI)\*

\* Only applicable to the EDAN SpO<sub>2</sub> module

PI is a numeric value indicating perfusion level. It reflects the perfusion level at the monitoring site.

As the measurement of SpO<sub>2</sub> is based on the pulsation caused by the blood flow through the vessel, PI is in relation to the strength of the pulse. Also, you can use PI as a signal quality indicator for the measurement of SpO<sub>2</sub>.

PI is indicated by a value ranging from 0 to 10. The bigger the value is, the better the perfusion and the signal quality will be. The perfusion level and the signal quality are at their maximum when the value reaches 10. When PI is below 2, it indicates the low perfusion and the poor signal quality at the monitoring site; you need to reposition the sensor or find a better site.

The PI value will be displayed in the SpO<sub>2</sub> parameter area.

## 11.8 Measuring SpO<sub>2</sub> and NIBP Simultaneously

53 While measuring SpO<sub>2</sub> and NIBP on the same limb simultaneously, the user can set **NIBP Simul** to **On** in **SpO<sub>2</sub> Setup** menu to lock the SpO<sub>2</sub> alarm status until the NIBP measurement ends. If **NIBP Simul** is set to **Off**, low perfusion caused by NIBP measurement may lead to inaccurate SpO<sub>2</sub> readings and therefore cause false physiological alarms.

0.3

	Pediatric Mode	SYS: 40 mmHg to 230 mmHg DIA: 10 mmHg to 180 mmHg MAP: 20 mmHg to 195 mmHg	
	Neonatal Mode	SYS: 40 mmHg to 135 mmHg DIA: 10 mmHg to 100 mmHg MAP: 20 mmHg to 110 mmHg	
Cuff Pressure Measuring Range	0 mmHg to 300 mmHg		
Pressure Resolution	1 mmHg		
Maximum Mean Error	±5 mmHg		
Maximum Standard Deviation	8 mmHg		
Maximum Measuring Period	Adult/ Pediatric	120 s	
	Neonatal	90 s	
Typical Measuring Period	20 s to 35 s (depend on HR/motion disturbance)		
Overpressure Protection	Adult	297 mmHg±3 mmHg	
	Pediatric	245 mmHg±3 mmHg	
	Neonatal	147 mmHg±3 mmHg	
<b>PR</b>			
Measuring range	40 bpm to 240 bpm		
Accuracy	±3 bpm or 3.5%, whichever is greater		
<b>SunTech Module</b>			
Method	Oscillometric		
Mode	Manual, Auto, Continuous		
Measuring Interval in AUTO Mode	1/2/3/4/5/10/15/30/60/90/120/240 min		
Measuring Type	SYS, DIA, MAP, PR		
Measuring Range	Adult Mode	SYS: 40 mmHg ~ 260 mmHg DIA: 20 mmHg ~ 200 mmHg MAP: 26 mmHg ~ 220 mmHg	
	Pediatric Mode	SYS: 40 mmHg ~ 230 mmHg DIA: 20 mmHg ~ 160 mmHg MAP: 26 mmHg ~ 183 mmHg	
	Neonatal Mode	SYS: 40 mmHg ~ 130 mmHg DIA: 20 mmHg ~ 100 mmHg MAP: 26 mmHg ~ 110 mmHg	
Pressure Resolution	1 mmHg		
Maximum mean error	±5 mmHg		
Maximum standard deviation	8 mmHg		
Maximum measuring period	Adult/Pediatric	130 s	

12

### 13.11 Cleaning Mode

24 The cleaning mode can remove the dust and foreign matters in the air valve to ensure the accuracy of NIBP measurement. To start the cleaning mode, please select **User Maintain > NIBP Maintain > Cleaning Mode**, the monitor displays: **Be sure the cuff is disconnected from monitor**, after confirmation and clicking **Start Cleaning** button, cleaning mode starts. The cleaning mode lasts three minutes. In this mode, the monitor displays **Cleaning in progress**, the remaining time of cleaning mode and cuff value are also displayed. When the counting down finishes, the monitor exits cleaning mode automatically, if the user needs to exit the cleaning mode in advance, please click **Stop** button.

When the air pressure is abnormal, the monitor will automatically turn off the cleaning mode and display the prompt message: **Cleaning failed**.

**NOTE:**

Cleaning mode is only available when the patient type is adult.

### 13.12 Assisting Venipuncture

The user can use the NIBP cuff to cause a pressure close to diastolic pressure, so as to block the venous blood vessel and therefore help venipuncture. To assist venipuncture:

1. Select **NIBP Setup > Venipuncture**;
2. Select the appropriate **Cuff Pressure** according to the patient type;
3. Select **Start**, the monitor displays: **Venipuncture Starting**.
4. Wait until the monitor prompts **In venipuncture process**. If an abnormal alarm occurs before it, no follow-up operation can be carried out. Restart the procedure after checking if necessary;
5. Puncture vein and draw blood sample;
6. Select **Stop** to deflate the cuff. If you do not deflate the cuff, the cuff automatically deflates when the venipuncture time expires (170 seconds for adult and pediatric patient, 85 seconds for neonatal patient).

During venipuncture, pay attention to the cuff pressure and the countdown displayed in the NIBP numerics area. When the remaining time is 30 seconds, the monitor issues a reminder tone and the countdown displays in red, prompting the user that the venipuncture time is to expire.

**NOTE:**

- 1 Only when the monitor exits **Venipuncture** menu, the user can do other operations.
- 2 When the monitor is in DEMO mode, continuous measurement process, manual measurement process, sequence measurement process or auto measurement process, Assisting Venipuncture function is not available.

Channel	2 <i>7.2</i>	
Sensor type	YSI-10K and YSI-2.252K	
Technique	Thermal resistance	
Position	Skin, Oral, Rectum	
Unit	°C , °F	
Measuring Range	0°C to 50°C (32 °F to 122 °F) <i>7.1</i>	
Resolution	0.1°C (0.1 °F)	
Accuracy	Accuracy (not including sensor): ±0.1°C	
	Sensor accuracy: ≤ ±0.2°C	
Transient Response Time	≤30 s	
<b>IBP</b>		
Channel	2	
Technique	Direct invasive measurement	
Measuring range	Art	0 mmHg to +300 mmHg
	PA	-6 mmHg to +120mmHg
	CVP/RAP/LAP/ICP	-10 mmHg to +40 mmHg
	P1/P2	-50 mmHg to +300 mmHg
Resolution	1 mmHg	
Accuracy	±2% or ±1 mmHg, whichever is greater (not including sensor) ICP: 0 mmHg to 40 mmHg: ±2% or ±1 mmHg, whichever is greater; -10 mmHg to -1 mmHg: undefined	
Unit	kPa, mmHg, cmH2O	
<b>PR</b>		
Measuring Range	20 bpm to 300 bpm	
Resolution	1 bpm	
Accuracy	30 bpm to 300 bpm: ±2 bpm or ±2%, whichever is greater; 20 bpm to 29 bpm: undefined	
<b>CO<sub>2</sub> (EFM Module)</b>		
<b>Respironics Sidestream Module</b>		
Applicable Patient Type	Adult, Pediatric and Neonatal	
Method	Infra-red Absorption Technique	
Measure Parameters	EtCO <sub>2</sub> , FiCO <sub>2</sub> , AwRR	
Unit	mmHg/ %/ kPa	
Measuring Range	EtCO <sub>2</sub>	0 mmHg to 150 mmHg
	FiCO <sub>2</sub>	3 mmHg to 50 mmHg
	AwRR	2 rpm to 150 rpm
Resolution	EtCO <sub>2</sub>	1 mmHg
	FiCO <sub>2</sub>	1 mmHg
	AwRR	1 rpm
Accuracy	EtCO <sub>2</sub>	±2 mmHg, 0 mmHg to 40 mmHg

It takes 5 minutes for the temperature measurement to stabilize.

## 14.6 Selecting a Temperature for Monitoring

Select the temperature label according to the measurement site. The label is a unique identifier for each type of temperature.

To select the label,

1. Click the TEMP parameter area to enter **TEMP Setup** menu.
2. Select the appropriate label from the list for **T1** and **T2**.

Label	Description
<b>Tskin</b>	Skin temperature
<b>Trect</b>	Rectal temperature
<b>Tcore</b>	Core temperature

### NOTE

Tcore is only available when **TEMP Sensor** is **YSI-2.252K**.

## 14.7 Calculating Temp Difference

73 The monitor can calculate and display the difference between two temperature values by subtracting the second value from the first. The difference is labeled TD.

## iM20 Patient Monitor Specification

### Physical Specifications

iM20	Dimension	185 mm (W) × 116 mm (H) × 85.3 mm (D)	
	Weight	< 1.5 kg (Including battery, without accessories)	
EFM Module	Dimension	207 mm (W) × 116 mm (H) × 93.4 mm (D)	
	Weight	< 0.58 kg (Including sidestream CO <sub>2</sub> module.)	

### Power Supply

Voltage	DC 11.1 V - 19.8 V		
Current	1.27 A – 2.3 A		

### Power Adapter

AC Power Adapter	Input: 100 V - 240 Vac, 40VA Output: 15 V±5%dc, 24 VA		
DC Power Adapter	Input: 10 V - 16 Vdc Output: 15 V±5%dc		

### Battery

Battery Type	Rechargeable lithium-ion battery		
Capacitance	2400 mAh		
Operating Time	5.5 hrs		
Fast Charging Time	2.5 hrs	The monitor is off, 100% charge	
	2.3 hrs	The monitor is off, 90% charge	

### Display

Display screen	5 inch color TFT, touch screen		
Resolution	800×480		
Wave	A maximum of 13 waveforms (with 12-lead ECG function)		

### Data Storage

Internal Temporary Memory	Trend graph/trend table review	1 h, at 1 s resolution 150 hrs, at 1 min resolution
	Alarm/Monitoring Event data	Up to 200 sets
	NIBP Measurement Review	1200 sets
	Arrhythmia events	Up to 200 sets
	12-lead Diagnosis Review	Up to 50 sets
Non-volatile Memory (internal or external storage device)	A single piece of patient data maximally contains the following information:	
	Trend graph and trend table	240 hours, at 1 min resolution
	NIBP measurement review	1200 sets
	Alarm review	200 sets
	Arrhythmia event	200 sets
	12-lead diagnosis review	50 sets
Waveforms	48 hours	

### Wi-Fi

IEEE	802.11b/g/n
Frequency Band	2.4 GHz ISM band

160

## Chapter 19 Calculation and Titration Table

The monitor provides calculation function and titration table. Calculations are patient data that are not directly measured but calculated by the monitor.

The monitor can perform drug calculation, hemodynamic calculation, oxygenation calculation, ventilation calculation and renal function calculation.

### NOTE:

- 1 The drug calculation function acts only as a calculator. The patient weights in Drug Dose menu and in Patient Information menu are independent of each other. Therefore changing the Weight in Drug Dose menu will not change the weight in the Patient Information menu.
- 2 The calculation results are for reference only and the calculation significance must be determined by the physician.

---

### **WARNING**

---

The correctness of the input parameters and the suitability of the calculated results should be carefully verified. EDAN is not liable for any consequences arising from input or operation errors.

---

## 19.1 Drug Calculation

### 19.1.1 Calculation Procedures

1. The drug calculation window is displayed by selecting **Menu > Common Function > Calculation > Drug Dose**.
2. Select the right pull-down box of the **Drug** option and select the required drug name among the 15 drugs which are listed as follows. And the drug name of **Drug A, Drug B, Drug C, Drug D** and **Drug E** can be defined by the user.
  - Drug A, Drug B, Drug C, Drug D and Drug E
  - Aminophylline
  - Dobutamine
  - Dopamine
  - Epinephrine
  - Heparin
  - Isuprel
  - Lidocaine
  - Nipride
  - Nitroglycerin
  - Pitocin

## 8.5 Viewing Trend Screen

To view the trend screen, the user can press the shortcut key  on the shortcut widget screen or select **Menu > Display Setup > View Selection > TrendScreen**.

Select short trend to open **Short Trend Setup** menu, the user can set:

1. **Parameter**.
2. **Interval**: set the interval to **30 min**, **1 h** and **2 h**.

## 8.6 Viewing OxyCRG Screen

To view the OxyCRG screen, the user can press the shortcut key  on the shortcut widget screen or select **Menu > Display Setup > View Selection > OxyCRG**. This interface is always used in NICU because the SpO<sub>2</sub>, HR and RESP of the neonate are different from those of adults. OxyCRG is in the bottom half part of wave area; it consists of HR trend, SpO<sub>2</sub> trend and RR trend or compressed respiration waveform.

Select OxyCRG waveform to open **OxyCRG Setup** menu, you can set:

1. **Interval**: set the interval to **1 min**, **2 min** and **4 min**.
2. **Parameter**: to select **RESP** or **RR**.
3. **OxyCRG Review**: user can review the 24 hours OxyCRG parameters including HR, SpO<sub>2</sub>, RR. Clicking  or  to left or right move the screen for viewing OxyCRG. Click **Exit OxyCRG Review** to exit the interface.

## 8.7 Viewing Large Font Screen

To open the large font screen, please refer to the following steps:

1. Select the shortcut key  on the shortcut widget screen or.
2. Select **Menu > Display Setup > View Selection > Large Font** to open this interface.

You can view any available parameter by selecting the parameter from the pull-down list on each section.


## Chapter 18 Review

The monitor provides 150-hour trend data of all parameters, storage of 1200 NIBP measurement results, 200 alarm events, 200 arrhythmia events, and 50 sets of 12-lead analysis results. This chapter gives detailed instruction for review of all data.

### 18.1 Trend Graph Review

To review the trend graph, please press the **Trend Graph** key  on the shortcut widget screen or select **Menu > Review > Trend Graph**.

In the trend graph, the y-axis stands for the measurement value and x-axis stands for the time. With the exception of NIBP, other trends are displayed as continuous curves.

- 1 Trend curve area
- 2 Trend data: displays measurement values at the cursor indicated time.
- 3 Cursor
- 4 Cursor time

In the trend graph review window:

- Select **Parameter** and you can choose the required parameters to be displayed in the trend graph.
- To display a different parameter's trend, you can either:
  - ◆ Select ▼ beside the parameter name and choose the desired parameter from the pop-up list (as shown in red circle above).
  - ◆ Press the symbols ▲ and ▼ to switch parameters in batch.
- Select **Zoom** to adjust the trend scale. Once the trend scale on the trend graph review interface is adjusted, the trend scale of the corresponding parameter in **TrendScreen** of the main interface will also change.

## Chapter 6 Alarm Information

### 6.1 Physiological Alarm Information

#### WARNING

The physiological alarms including Asystole, Sustain VT, RESP APNEA, SpO<sub>2</sub> No Pulse, SpO<sub>2</sub> Desat, CO<sub>2</sub> APNEA, AG FiO<sub>2</sub> Low, and AG APNEA cannot be turned off.

Message	Cause	Alarm Level
<b>HR High</b>	HR measuring value is above the upper alarm limit.	User-selectable
<b>HR Low</b>	HR measuring value is below the lower alarm limit.	User-selectable
<b>ST-X High</b>	ST measuring value is above the upper alarm limit. (X stands for I, II, III, aVR, aVL, aVF, V, V1, V2, V3, V4, V5 or V6)	User-selectable
<b>ST-X Low</b>	ST measuring value is below the lower alarm limit. (X stands for I, II, III, aVR, aVL, aVF, V, V1, V2, V3, V4, V5 or V6)	User-selectable
<b>Asystole</b>	No QRS is detected for 4 consecutive seconds	High
<b>V-Fib/V-Tach</b>	4 consecutive seconds' fibrillation wave occurs; Or 5 consecutive ventricular beats, and ventricular HR $\geq 100$ bpm.	High
<b>Run PVCs</b>	$3 \leq$ the number of consecutive PVCs $< 5$	User-selectable
<b>Couplet</b>	2 consecutive PVCs	User-selectable
<b>PVC Bigeminy</b>	A dominant rhythm of N, V, N, V (N = supraventricular beat, V = ventricular beat) was detected.	User-selectable
<b>PVC Trigeminy</b>	A dominant rhythm of N, N, V, N, N, V	User-selectable
<b>R on T</b>	A type of single PVC under the condition that HR $< 100$ , R-R interval is less than 1/3 the average interval, followed by a compensating pause of 1.25X the average R-R interval (the next R wave advances onto the previous T wave).	User-selectable
<b>PVC</b>	Single PVC detected in normal heartbeats, and the number of consecutive single PVC $\geq 4$ within 30 s.	User-selectable
<b>Tachy</b>	Adult: RR interval for 5 consecutive QRS complex $\leq 0.5$ s. Pediatric/neonatal: RR interval for 5 consecutive QRS complex $\leq 0.375$ s.	User-selectable

## Chapter 3 Basic Operation

This user manual is based on the maximum configuration and therefore your monitor may not have all of the functions and options described in the manual. Also, illustrations in this manual serve as examples only and do not necessarily reflect the setup on your monitor. The content displayed on your monitor depends on the way it has been tailored for your hospital.

You may frequently use the follow functions:


- ♦ ECG monitoring (Refer to Chapter *Monitoring ECG* for more information.)
- ♦ SpO<sub>2</sub> monitoring (Refer to Chapter *Monitoring SpO<sub>2</sub>* for more information.)
- ♦ PR monitoring (Refer to Chapter *Monitoring PR* for more information.)
- ♦ NIBP monitoring (Refer to Chapter *Monitoring NIBP* for more information.)
- ♦ Alarm (Refer to Chapter *Alarms* for more information.)

### 3.1 System Components

#### 3.1.1 Main Unit


Front View

*elite V5/ elite V6*


Rear View

elite V5/ elite V6

22

1	Speaker	For alarm tones, pulse tones and so forth.
2	Equipotential grounding terminal	If the monitor is used together with other devices, connect this terminal to eliminate potential ground differences between devices.
3	Power cable safety latch	Used to prevent the power cable from loosing or falling. Place the latch on the power cable and press it down firmly to ensure that it secures the power cable.
4	Power supply interface	for connecting AC power cable.
5	SD card slot	Used to mount an SD memory card.
6	VGA output	It enables the VGA video output.
7	USB interfaces	They support USB1.0/2.0 output.
8	Multifunctional port: Nurse call port/analog output/defibrillator synchronization	If users select it as nurse call, it connects the monitor to the hospital's nurse call system. Alarms indications are alerted through the nurse call system if configured to do so. If users select it as analog output, the monitor outputs the waveform through the port.
9	Network interface	It connects the monitor to the central monitoring system (MFM-CMS) or gateway via standard network cable, which enables MFM-CMS or gateway to achieve bidirectional communication with the monitor.
10	Extended video interface	It connects a secondary display, which extends the display capability of your monitor.
11	PAM connector	It connects the Parameter Amplifier Mainframe to the monitor.
12	RS232 interface	Connect it to communicate with other devices.
13	Heat sink	<u>Adopt No-fan design, which is dust-free, low noise and low consumption.</u>
14	Anti-thief lock	

1	Alarm mute indicator	When the audible alarm is mute, the indicator is in red.
2	Physiological alarm indicator	When a physiological alarm occurs, the indicator lights on or flashes with different frequencies and colors reflecting the alarm level.
3	Technical alarm indicator	When a technical alarm occurs, the indicator lights on or flashes with different frequencies and colors reflecting the alarm level.

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