**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## AIA-PACK® BNP CONTROL SET

### INTENDED USE

The AIA-PACK BNP CONTROL SET is intended for IN VITRO DIAGNOSTIC USE ONLY, by healthcare professionals, for performing quality control procedures with the ST AIA-PACK BNP assay and the CL AIA-PACK BNP assay.

### SUMMARY AND EXPLANATION

The AIA-PACK BNP CONTROL SET contains buffered bovine serum albumin with the assigned levels of brain natriuretic peptide (BNP). The BNP controls shall be run at the beginning of each day on which BNP assays are scheduled.

### CONTENTS

Cat. No.	0025428
2 x 1 mL	AIA-PACK BNP CONTROL LEVEL 1 Buffered bovine serum albumin containing BNP (Lyophilised). See the vial label for the assigned concentration range.
2 x 1 mL	AIA-PACK BNP CONTROL LEVEL 2 Buffered bovine serum albumin containing BNP (Lyophilised). See the vial label for the assigned concentration range.

### WARNINGS AND PRECAUTIONS

1. The AIA-PACK BNP CONTROL SET is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local Tosoh sales representative.
3. Although material derived from human origin is not used for these controls, it is recommended that this product be handled with the same precautions as used for patient samples.
4. Do not use beyond the expiry date.
5. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

Bring the controls to 18°-25°C for use. Using volumetric pipettes, reconstitute the lyophilised controls accurately to the volume of 1 mL with CAP Class I water or the clinical laboratory reagent water (formerly **NCCLS Type I**) defined by **CLSI GP40-A4-AMD guideline**. Allow the lyophilised material to dissolve fully.

### STORAGE AND STABILITY

- Always store the AIA-PACK BNP CONTROL SET in an upright position at 2°-8°C when not in use.
- When stored unopened and refrigerated at 2°-8°C, the control set will remain stable until the expiry date on the label.
- The controls should be used within 7 days of opening or reconstituting, provided the vials are kept tightly sealed and refrigerated at 2°-8°C.

### PROCEDURE

Refer to the Tosoh AIA® System Operators' manual or the Tosoh AIA-CL Analyzer Operator's Manual for additional procedural instructions regarding quality control.

1. Load the appropriate amount of test cups on the instrument.
2. Add the appropriate amount of each control to sample cups. (Refer to the instrument worksheet for the sample volume.)
3. Print a work list and place the sample cups in the position indicated.
4. Select START. Verify that the positions on the work list match the starting position on the screen.

### ASSIGNMENT OF VALUES

The AIA-PACK BNP CONTROL SET contains assigned concentration range of BNP. The assigned range is determined on a lot-by-lot basis and is designed to provide target control levels of BNP. Since the assay values are dependent upon assay procedures as well as several other factors, each laboratory should establish its own range for the assay procedure being monitored.

### LIMITATIONS

The AIA-PACK BNP CONTROL SET is designed solely for use with ST AIA-PACK BNP and CL AIA-PACK BNP assay procedures.

"AIA" and "AIA-PACK" are the registered trademarks of Tosoh Corporation in the European Union, etc.


### TOSOH CORPORATION

3-8-2, Shiba-Koen, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO, BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

Net Volume (after reconstitution for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## ST AIA-PACK cTnI 3rd-Gen

### INTENDED USE

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ST AIA-PACK cTnI 3rd-Gen is designed for in vitro diagnostic use only for the quantitative measurement of cardiac troponin I (cTnI) in human serum, heparinised plasma or EDTA plasma on Tosoh AIA System Analysers. Cardiac troponin I measurement are used as a diagnostic aid for acute myocardial infarction (AMI).

### SUMMARY AND EXPLANATION

Troponin exists as a complex of three distinct proteins: troponin T, troponin C and troponin I<sup>1, 2</sup>. The troponins regulate the calcium-modulated interaction of actin and myosin in striated (skeletal and cardiac) muscle, with troponin I acting as the actinomyosin ATPase inhibitor.

Three isoforms of troponin I have been identified, which are products of different genes, each with unique amino-acid sequences<sup>3</sup>. Cardiac troponin I (cTnI) is found exclusively in cardiac tissue<sup>4</sup>. Though the three forms of troponin I share some homologous structure, cTnI has a sequence of 31 amino-acid residues at the amino terminus of the molecule which renders cTnI distinctly antigenically different from the other isoforms of troponin I, as well as troponin C and troponin T<sup>5</sup>.

Biochemical markers of myocardial injury are considered the "gold standard" for the diagnosis of AMI and are particularly important in non-diagnostic ECG patients<sup>6, 7</sup>. Clinical investigations have shown high cTnI concentrations are extremely useful in aiding the diagnosis of cardiac specific injury in patients with acute myocardial infarction, myocardial contusion<sup>8</sup>, or unstable angina<sup>9</sup>. Following an AMI, cTnI is released into the bloodstream from the damaged heart muscle and, within 4 - 8 hours, levels of cTnI rise above levels established for non-AMI specimens. Peak concentrations are generally reached 12 - 18 hours after AMI's onset and cTnI concentrations remain high for 5 - 10 days<sup>10-12</sup>.

Troponin I shows a temporal rise and fall in AMI patients similar to that seen with CK-MB mass. High levels of cTnI are detectable at a later stage than CK-MB on the cardiac injury. Updated guidelines<sup>13, 14</sup> by the American College of Cardiology and American Heart Association indicate that cTnI is more useful than CK-MB in detecting myocardial injury, especially in the presence of skeletal muscle damage and lesions (rhabdomyolysis, polytraumatism), because cTnI has high cardiac specificity and diagnostic efficiency of myocardial injury<sup>3, 9, 15, 16</sup>.

### PRINCIPLE OF THE ASSAY

The ST AIA-PACK cTnI 3rd-Gen is a two-site immuno-enzymometric assay, performed entirely in ST AIA-PACK cTnI 3rd-Gen test cups. Cardiac troponin I present in the test sample is bound with monoclonal antibody immobilised on magnetic beads and enzyme-labelled monoclonal antibody in the test cups. The magnetic beads are washed to remove any non-bound enzyme-labelled monoclonal antibodies and then incubated with a fluorogenic substrate, 4-methylumbelliferyl phosphate (4MUP). The amount of enzyme-labelled monoclonal antibody that binds to the beads is in direct proportion to the cTnI concentration in the test sample. A standard curve is constructed, and unknown sample concentrations are calculated using this curve. The ST AIA-PACK cTnI 3rd-Gen replaces the ST AIA-PACK cTnI 2nd-Gen product and has been labelled as third-generation in order to avoid any confusion.

### MATERIAL SUPPLIED

Cat. No.	0025215
ST AIA-PACK cTnI 3rd-Gen	5 trays x 20 test cups

Plastic test cups containing twelve lyophilised magnetic beads coated with anti-cTnI mouse monoclonal antibody and 100 µl anti-cTnI mouse monoclonal antibody, conjugated to bovine alkaline phosphatase, with sodium azide as a preservative.

### MATERIALS REQUIRED BUT NOT SUPPLIED

The following materials are required to perform cardiac troponin I analysis using the ST AIA-PACK cTnI 3rd-Gen (Cat. No. 0025215) on Tosoh AIA System Analysers. The following are available separately from Tosoh.

Materials	Cat. No.
AIA Nex•IA/AIA-21	0018539
AIA Nex•IA/AIA-21 LA	0018540
AIA-1800 ST	0019836
AIA-1800 LA	0019837
AIA-2000 ST	0022100
AIA-2000 LA	0022101
AIA-600 II	0019014
AIA-600 II BCR	0019328
AIA-900	0022930
AIA-360	0019945
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK SUBSTRATE REAGENT II / AIA-PACK SUBSTRATE RECONSTITUENT II	
ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET	0025315
ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (1)	0 ng/mL
ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (2)	0.2 ng/mL (approx.)
ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (3)	1.0 ng/mL (approx.)
ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (4)	5.0 ng/mL (approx.)
ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (5)	20.0 ng/mL (approx.)
ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (6)	60.0 ng/mL (approx.)
ST AIA-PACK cTnI 3rd-Gen SAMPLE DILUTING SOLUTION	0025515
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
SAMPLE CUPS	0018581
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970
AIA-PACK SAMPLE TREATMENT CUP	0020971
Additional Requirements for AIA Nex•IA/AIA-21 only:	
PIPETTE TIPS	0018552
PRELOADED PIPETTE TIPS	0018583
Additional Requirements for AIA-600 II, AIA-900, AIA-1800 and AIA-2000 only:	
PIPETTE TIPS	0019215
TIP RACK	0019216
PRELOADED PIPETTE TIPS	0022103

Only materials obtained from Tosoh should be used. Materials obtained elsewhere should not be substituted, since assay performance is characterised based strictly on Tosoh materials.

### WARNINGS AND PRECAUTIONS

1. The ST AIA-PACK cTnI 3rd-Gen is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the aluminium pouch for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. Test cups from different lots or different assays shall not be mixed within a tray.

4. The ST AIA-PACK cTnI 3rd-Gen contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always rinse with large volumes of water to prevent any azide accumulation.
5. The material derived from human blood used in the preparation of only ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET has been tested by FDA-approved methods, and found negative for the presence of HBsAg and antibodies to HIV-1 and HCV. Since no testing method can give complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions as those used for patient samples.
6. Do not use beyond the expiry date.
7. The ST AIA-PACK cTnI 3rd-Gen has been designed so that the high dose "hook effect" will not be a problem for the vast majority of samples. The "hook effect" phenomenon may occur at cTnI concentrations of more than 20,000 ng/mL.
8. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.
9. After opening, the vial of ST AIA-PACK cTnI 3rd-Gen SAMPLE DILUTING SOLUTION should be kept tightly sealed with a clean rubber cap. Sealing with dirty material may cause deterioration of the reagent.
10. The remaining sample diluting solution after use should not be mixed with another vial but be discarded to avoid contamination.
11. Serum, dust, metal, or microorganism contamination may cause degradation of reconstituted substrate solution. Store in a clean environment, away from direct sunlight and ultraviolet light.
12. TOSOH recommends that a new pouch of the test cups should be used for calibration.

## STORAGE AND STABILITY

All unopened materials are stable until the expiry date on the label when stored at the specified temperature.

<b>Materials</b>	<b>Cat. No.</b>
2°–8°C:	
ST AIA-PACK cTnI 3rd-Gen	0025215
ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET	0025315
ST AIA-PACK cTnI 3rd-Gen SAMPLE DILUTING SOLUTION	0025515
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
1°–30°C:	
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970
AIA-PACK SAMPLE TREATMENT CUP	0020971

- After opening the aluminium pouch, ST AIA-PACK cTnI 3rd-Gen test cups can be left on-board of the TOSOH AIA System Analyzers (18-25°C) for a maximum of 10 days (10 x 24 hours). When stored overnight at 2-8°C, the test cups can be used for up to 30 days (30 cycles of 8 hours on board and 16 hours in the refrigerator). Once the aluminium pouch is opened, the test cups must be used within 30 days.
- ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET must be kept tightly sealed and refrigerated at 2-8°C. After opening or reconstituting, the calibrators should be used within 1 day.
- After opening, ST AIA-PACK cTnI 3rd-Gen SAMPLE DILUTING SOLUTION can be left on-board of the TOSOH AIA System Analyzers (18-25 C) for a maximum of 2 days (2 x 24 hours). When stored over night at 2-8°C, the sample diluting solution can be used for up to 6 days (6 cycles of 8 hours on board and 16 hours in the refrigerator). The sample diluting solution should not be used beyond 90 days after opening, even if it is sealed and stored in the refrigerator.
- Reconstituted substrate solution will remain stable for 3 days at 18°-25°C or 30 days at 2°-8°C.

- Working diluent and wash solutions are stable for 30 days at 18°-25°C.
- Reagents should not be used if they appear cloudy or discoloured.

## SPECIMEN COLLECTION AND HANDLING

- Serum, heparinised plasma or EDTA plasma is required for the assay. Citrated plasma SHOULD NOT BE USED.
- When using serum, a venous blood sample is collected aseptically without additives. Store at 18°-25°C until a clot has formed (usually 15-45 minutes), then centrifuge to obtain the serum specimen for assay.
- When using heparinised plasma or EDTA plasma, a venous blood sample is collected aseptically with the specified additive. The blood sample should be collected in accordance with the specimen collection-tube manufacturer's instructions, as well as these package insert instructions, in order to avoid false low results<sup>2, 17)</sup>. Centrifuge and separate plasma from the packed cells as quickly as possible.
- Inadequate centrifugation or the presence of fibrin or particulate matter in the sample may cause an erroneous result.
- Samples containing inhibitors of alkaline phosphatase may cause erroneous results.
- Inspect all samples for air bubbles and foaming. Remove any air bubbles before assay.
- Specimen types should not be used interchangeably during serial monitoring of an individual patient. Measured concentrations may vary slightly between sample types in certain patients.
- Samples may be stored at 2°-8°C for up to 24 hours before analysis. If the analysis cannot be performed within 24 hours, the sample should be stored frozen at -20°C or below for up to 60 days.
- Repeated freeze-thaw cycles should be avoided. Turbid serum samples or samples containing particulate matter should be centrifuged before testing. Before assaying, slowly bring frozen samples to 18°-25°C and mix gently.
- The sample required for analysis is 50 µl.

## PROCEDURE

For the AIA NexIA/AIA-21, AIA-600 II, AIA-900, AIA-1800, AIA-2000 and AIA-360, please refer to their Operator's Manual for detailed instructions.

### I. Reagent Preparation

#### A. Substrate Solution

Bring all reagents to 18°-25°C before preparing the working reagent. Add the entire contents of the AIA-PACK SUBSTRATE RECONSTITUENT II (100 ml) to the AIA-PACK SUBSTRATE REAGENT II (Lyophilised), mix thoroughly to dissolve solid material.

#### B. Wash Solution

Add the entire contents of the AIA-PACK WASH CONCENTRATE (100 ml) to approximately 2.0 l of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by **CLSI GP40-A4-AMD guideline, mix well, and adjust the final volume to 2.5 l.**

#### C. Diluent

Add the entire contents of the AIA-PACK DILUENT CONCENTRATE (100 ml) to approximately 4.0 l of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by **CLSI GP40-A4-AMD guideline, mix well and adjust the final volume to 5.0 l.**

### II. Calibration Procedure

#### A. Calibration Curve


The calibrators for use with the ST AIA-PACK cTnI 3rd-Gen are prepared gravimetrically and compared against internal reference standards.

The calibration curve for the ST AIA-PACK cTnI 3rd-Gen will remain stable for up to 90 days. Calibration stability is monitored by quality-control performance and is dependent on proper

reagent handling and Tosoh AIA System maintenance in accordance with the manufacturer's instructions. Recalibration may be necessary more frequently if controls are out of the established range for this assay or if certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of wash probe or detector lamp adjustment or change).

For further information regarding instrument operation, consult the TOSOH AIA System Operator's Manual.

A sample calibration curve from the AIA-1800 follows and shows the algorithm used to calculate results.

#### B. Calibration Procedure

1. Refer to the appropriate Tosoh AIA System Operators' Manual for procedural instructions.
2. Verify that both the calibrator lot and concentration numbers have been correctly entered on the software.
3. The ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (1) is provided ready for use.
4. The ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (2)-(6) are lyophilised. All levels should be reconstituted with 1.0 ml CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline.
5. Tosoh recommends that all calibrators be run in triplicate.

#### C. Calibration Acceptability Criteria

1. The mean rate for the ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (1) should be  $\leq 0.2$  nmol/(L\*s).
2. Since there is a direct relationship between concentration and rate, the rate should rise as concentrations increase.
3. Repeat values should be within a 10% range.

#### D. Calibration Review and Acceptance

1. Review the calibration curve carefully, applying the criteria listed above.
2. Edit the calibration if necessary, then accept the calibration.

For further information regarding calibration, consult the Tosoh AIA System Operators' Manual.

### III. Quality-Control Procedure

#### A. Commercially-Available Controls

Commercially-available controls shall be run at least once per day. It is recommended that at least two (2) levels of controls, normal and abnormal, be used. Laboratory policy for this particular assay specifies the following:

Control Material: \_\_\_\_\_  
 Frequency: \_\_\_\_\_

Lot number of control material, acceptable limits, and remedial action to be taken if controls fail to meet laboratory criteria, will be found in a separate quality-control document kept by the laboratory.

#### B. Quality-Control Procedure

1. Assay quality-control specimens as instructed in the specific Operators' Manual for your analyser. In addition, refer to the Tosoh AIA System Operators' Manual for detailed instructions on defining and editing files.
2. Quality-control material to be run with this assay is specified by individual laboratory policy.

### IV. Specimen Processing

#### A. Preparation

Following the specific instructions in the Operators' Manual for the analyser, place samples on the instrument appropriately. Bar-coded primary tubes, as well as sample cups, can be run on the AIA Nex•IA/AIA-21, AIA-600 II, AIA-900, AIA-1800, AIA-2000 and AIA-360.

#### B. Assay Procedure

1. Provide a sufficient quantity of ST AIA-PACK cTnI 3rd-Gen test cups for the number of samples to be run.
2. Load patient samples as instructed in the Operators' Manual and proceed with analysis.

Note: The AIA Nex•IA/AIA-21, AIA-600 II, AIA-900, AIA-2000 and AIA-1800 will require AIA-PACK SAMPLE TREATMENT CUP if on-board dilutions are used.

### PROCEDURAL NOTES

1. Lyophilised substrate must be completely dissolved.
2. Ligand assays performed by the TOSOH AIA System Analyzers require that the laboratory use water designated by the CAP as Class I or by the CLSI as the clinical laboratory reagent water. Water should be tested at least once per month and should be free of particulate matter including bacteria. The pH of the water should also be routinely tested. For further information, consult the CLSI document GP40-A4-AMD, Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline-Fourth Edition.
3. If a specimen cTnI concentration is found to be greater than the upper limit of the assay range, 50 ng/mL, the specimen should be diluted with the ST AIA-PACK cTnI 3rd-Gen SAMPLE DILUTING SOLUTION and reassayed in accordance with the Assay Procedure. The recommended dilution for specimens containing more than 50 ng/mL, is 10-fold dilution. It is desirable to dilute the specimen so that the specimen reads between 0.02 and 50 ng/mL. The dilution factor should be entered on the software. For further information on the dilution of specimens, refer to the Tosoh AIA System Operators' Manual.
4. The Tosoh AIA System Analyser can store two different calibration curves for each analyte at one time. Thus, up to two different lots of ST AIA-PACK cTnI 3rd-Gen test cups can be used in the same run.
5. If the assay specifications for this test are not ready on the system software, the specifications must be entered under test code 088.

### CALCULATION OF RESULTS

The Tosoh AIA System Analysers perform all sample and reagent handling operations automatically. The Tosoh AIA System Analysers read the rate of fluorescence produced by the reaction and automatically convert the rate to cTnI concentration in ng/mL.

For samples requiring dilution, the AIA Nex•IA/AIA-21, AIA-600 II, AIA-900, AIA-2000 and AIA-1800 will automatically perform dilutions and calculate results if the dilution factors are entered on the software. For detailed information regarding programming dilutions, consult the appropriate TOSOH AIA System Operators' Manual.

## EVALUATION OF RESULTS

### Quality Control

In order to monitor and evaluate the precision of the analytical performance, it is recommended that commercially-available control samples shall be assayed in accordance with local regulations.

Minimum recommendations for the frequency of running internal control material are:

After calibration, three levels of the internal control are run in order to accept the calibration curve.

The three levels of controls are repeated if certain service procedures are performed (e.g. temperature adjustments, sampling-mechanism changes, maintenance of the wash probe or detector-lamp adjustments or changes).

After daily maintenance, at least two levels of the control shall be run in order to verify overall performance of the Tosoh AIA System Analysers.

If one or more control value(s) exceed the acceptable range, it will be necessary to investigate the validity of the calibration curve before reporting patient results.

Standard laboratory procedures should be followed in accordance with the strict regulatory agency under which the laboratory operates.

### LIMITATIONS

- For diagnostic purposes, results obtained from this assay should be used in conjunction with other data (e.g. symptoms, results of other tests, clinical impressions, therapy, etc.).
- Physicians should be made aware that the published literature <sup>18)</sup> indicates possible interference in troponin I assays - i.e. both false-positive and false-negative results. Studies on patients diagnosed with myocarditis and dilated cardiomyopathy indicate that disease-specific cardiac auto-antibodies may be present in the serum and plasma of these and other patients. Auto-antibodies are generated against cardiac proteins, which could potentially cause false-negative results if the auto-antibodies produced interfere with the epitopes used in a particular manufacturer's assay.
- Patient samples may also contain human anti-mouse antibody (HAMA) due to either natural antibody production or antibodies produced in response to therapy regimens. HAMA may cause either false-positive or false-negative results in assays using mouse monoclonal antibodies. Other heterophile antibodies are also known to interfere in assays of this type. ST AIA-PACK cTnI 3rd-Gen has been designed to minimise the effects of heterophile antibodies, but interference from high titers cannot be ruled out.
- Results from this or any other in vitro diagnostic procedures that do not correlate with the clinical presentation of the patient should be interpreted with extreme caution.
- The use of a single cTnI result on a patient should be discouraged. As with other cardiac markers, a pattern of rise and fall over time should be observed to aid accurate interpretation of the results. Users are reminded that rises in cTnI may occur due to causes other than myocardial infarction, and the absence of a cTnI value over a certain level does not rule out AMI or serious cardiac injury. Assays from different manufacturers may yield different results. Reactivity of the epitopes used in the assay varies with the form of cTnI present in a specific specimen. Also, assay calibration between manufacturers is not standardised.
- Using ST AIA-PACK cTnI 3rd-Gen, the highest measurable concentration of cardiac troponin I in specimens without dilution is 50 ng/mL, and the lowest measurable concentration in specimens is 0.02 ng/mL (assay sensitivity).
- The exact linearity of the ST AIA-PACK cTnI 3rd-Gen depends on the particular calibrator lot in use. Although the approximate value of the highest calibrator is 60.0 ng/mL, the exact concentration may be slightly different. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 50 ng/mL.
- Although haemolysis has an insignificant effect on the assay, haemolysed samples may indicate mistreatment of a specimen before assaying, and the results should be interpreted with caution.

- Lipaemia has an insignificant effect on the assay, except in the case of gross lipaemia, where spatial interference may occur.
- Samples containing fibrin may exhibit either falsely-raised or falsely-decreased results. Fibrin must be eliminated from the sample before assaying begins.
- Specimens from patients taking medicines and/or medical treatment may produce misleading results.
- Samples from patients under Asfotase Alfa (Genetical Recombination) treatment may cause falsely elevated/decreased results.
- For a fuller understanding of the limitations of this procedure, please refer to the SPECIMEN COLLECTION AND HANDLING, WARNINGS AND PRECAUTIONS, STORAGE AND STABILITY and PROCEDURAL NOTES sections of this insert sheet.

### EXPECTED VALUES

Each laboratory should set a reference interval corresponding to the characteristics of the population being tested. As with all diagnostic procedures, clinical results must be interpreted having due regard to concomitant medications being administered to the patient. Cardiac troponin I results should be used in conjunction with the patient's overall clinical picture. The protocol used for specimen-collection intervals will have an effect on cut-off values, since measurable cardiac troponin I may not be detectable until 4 - 6 hours after the onset of symptoms in an AMI. Interpretation should be made in the light of the time-interval between symptoms and specimen-collection.

#### I. Reference Ranges

The interval given here was set in heparinised plasma samples from 343 non-hospitalised, apparently healthy, male and female individuals varying in age and from Asian origin. European studies showed that these reference values are also applicable for a European population.

Number of Samples (n)	343
Reference Interval	< 0.06 ng/mL
99th percentile	0.04 ng/mL

99procentilê = 0.04 ng/ml  
iki 10% CV ties 0.035 ng/ml

#### II. Functional Sensitivity

Based on the imprecision profile, functional sensitivity was 0.035 ng/mL and 0.01 ng/mL at 10% CV and 20% CV respectively.

### CONVERSION FACTORS

Cardiac troponin I (cTnI) concentrations in this application are in units of ng/mL. The obtained values in ng/mL must be multiplied by 1.0 when converting to values in µg/L.

### PERFORMANCE CHARACTERISTICS

#### ACCURACY

- Recovery: Three serum samples and three heparinised plasma samples were spiked with three different levels of cTnI and assayed before and after spiking.

Sample	Initial	cTnI	Expected	Measured	Percent
	Value	Added	Value	Value	Recovery
	(ng/mL)	(ng/mL)	(ng/mL)	(ng/mL)	(%)
Serum A1	0.0	20.1	20.1	19.0	94.7
	0.0	10.0	10.0	9.71	96.9
	0.0	5.01	5.01	4.66	92.9
Serum B1	0.0	20.1	20.1	21.9	109.1
	0.0	10.0	10.0	9.32	92.9
	0.0	5.01	5.01	4.70	93.7
Serum C1	0.0	20.1	20.1	19.6	97.6
	0.0	10.0	10.0	9.30	92.8
	0.0	5.01	5.01	4.68	93.4

Sample	Initial	cTnI	Expected	Measured	Percent
	Value (ng/mL)	Added (ng/mL)	Value (ng/mL)	Value (ng/mL)	Recovery (%)
Plasma A1	0.0	20.1	20.1	21.0	104.5
	0.0	10.0	10.0	9.97	99.4
	0.0	5.01	5.01	4.68	93.3
Plasma B1	0.0	20.1	20.1	21.9	109.0
	0.0	10.0	10.0	9.50	94.8
	0.0	5.01	5.01	4.64	92.5
Plasma C1	0.0	20.1	20.1	21.7	108.5
	0.0	10.0	10.0	9.20	91.7
	0.0	5.01	5.01	4.74	94.6

- b. Dilution: Three serum samples and three heparinised plasma samples containing high concentrations of cTnI were serially diluted with ST AIA-PACK cTnI 3rd-Gen SAMPLE DILUTING SOLUTION and assayed.

Sample	Dilution Factor	Expected	Measured	Percent
		Value (ng/mL)	Value (ng/mL)	Recovery (%)
Serum A2	none	32.7	32.7	100.0
	7.5/10	24.5	24.0	97.8
	5.0/10	16.3	15.8	96.6
	2.5/10	8.17	7.92	96.9
	1.0/10	3.27	3.08	94.4
Serum B2	none	32.5	32.5	100.0
	7.5/10	24.3	24.7	101.4
	5.0/10	16.2	16.3	100.5
	2.5/10	8.12	7.67	94.5
	1.0/10	3.25	3.28	101.0
Serum C2	none	14.3	14.3	100.0
	7.5/10	10.7	11.0	102.8
	5.0/10	7.14	7.30	102.1
	2.5/10	3.57	3.70	103.5
	1.0/10	1.43	1.61	112.7
Plasma A2	none	34.7	34.7	100.0
	7.5/10	26.0	25.8	99.2
	5.0/10	17.4	16.9	97.3
	2.5/10	8.68	8.72	100.4
	1.0/10	3.47	3.69	106.3
Plasma B2	none	31.1	31.1	100.0
	7.5/10	23.3	23.4	100.3
	5.0/10	15.5	15.5	99.5
	2.5/10	7.77	7.61	98.0
	1.0/10	3.11	3.04	97.8
Plasma C2	none	14.1	14.1	100.0
	7.5/10	10.6	11.0	103.9
	5.0/10	7.06	7.18	101.6
	2.5/10	3.53	3.71	104.9
	1.0/10	1.41	1.27	89.9

- c. Linearity: The linearity for ST AIA-PACK cTnI 3rd-Gen was determined, based on the guidance from CLSI Protocol EP6-A. The linearity was measured on the AIA-2000 Instrument and has been demonstrated to be linear from 0.02 to 50 ng/mL, within +/- 10% difference in this interval.

#### PRECISION

- a. Within-run precision was determined using six controls in a total of 20 runs. Within each run, one set of duplicates per control was assayed. The mean of each duplicate was used to obtain the pooled standard deviation (SD), which was then used to calculate the coefficient of variation (CV).

Sample	Mean (ng/mL)	Pooled SD (ng/mL)	CV (%)
Serum A3	0.436	0.0104	2.4
Serum B3	4.01	0.0554	1.4
Serum C3	32.3	0.420	1.3
Plasma A3	0.422	0.0103	2.5
Plasma B3	3.27	0.0484	1.5
Plasma C3	31.4	0.318	1.0

- b. Total precision was determined by the duplicate assay of six controls in 20 separate runs. The means of each run were used to calculate the Pooled standard deviation (SD) and coefficient of variation (CV).

Sample	Mean (ng/mL)	Pooled SD (ng/mL)	CV (%)
Serum A3	0.436	0.0136	3.1
Serum B3	4.01	0.104	2.6
Serum C3	32.3	0.730	2.3
Plasma A3	0.422	0.0146	3.5
Plasma B3	3.27	0.0938	2.9
Plasma C3	31.4	0.624	2.0

#### CORRELATION

- a. The correlation (1) between a commercially-available cTnI test (x) and ST AIA-PACK cTnI 3rd-Gen (y), and the correlation (2) between ST AIA-PACK cTnI 2nd-Gen (x) and ST AIA-PACK cTnI 3rd-Gen (y) were carried out using 97 and 95 patient specimens respectively.

	(1)	(2)
Slope	0.989	0.939
y-Intercept	-0.194	-0.213
Correlation Coefficient	0.998	0.993
Number of Samples	97	95

- b. The correlation (1) between heparinised plasma (x) and serum (y), and the correlation (2) between heparinised plasma (x) and EDTA plasma (y), on ST AIA-PACK cTnI 3rd-Gen were carried out using 96 and 50 patient specimens respectively.

	(1)	(2)
Slope	1.079	0.988
y-Intercept	-0.005	-0.043
Correlation Coefficient	0.997	0.996
Number of Samples	96	50

#### SPECIFICITY

The following substances were tested for cross-reactivity. Cross-reactivity (%) is the percentage of the compound that will be identified as cTnI. If those compounds are present in the specimen at the same concentration as cTnI, the final result will be increased by those percentages.

Compound	Concentration	Cross-reactivity (%)
Cardiac troponin T	500-2,000 ng/mL	0.05
Cardiac troponin C	500-2,000 ng/mL	N.D.*
Skeletal troponin I	500-2,000 ng/mL	N.D.*
*N.D.: Not Detected		

#### SENSITIVITY

- (1) Minimum Detectable Concentration: The minimum detectable concentration (MDC) of cardiac troponin I is estimated to be 0.02 ng/mL. MDC is defined as the concentration of cTnI that corresponds to the rate of fluorescence, which is two standard deviations from the mean rate of fluorescence of 20 repeat determinations of ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (1).
- (2) Limit of Detection: The limit of detection for ST AIA-PACK cTnI 3rd-Gen was determined based on the guidance of CLSI Protocol EP17-A. A blank sample was measured in 60 replicates. Six low level samples were measured in 10 replicates each. As a result, the limit of detection for ST AIA-PACK cTnI 3rd-Gen is approximately 0.008 ng/mL.

#### INTERFERENCE

Interference is defined, for the purposes of this study, with recovery outside of 10% of the known concentration of the specimen after the following substances are added to human specimens.

- Haemoglobin (up to 440 mg/dl), free bilirubin (up to 18 mg/dl) and conjugated bilirubin (up to 18 mg/dl) do not interfere with the assay.
- Lipaemia, as indicated by triglyceride concentration (up to 1,600 mg/dl), does not interfere with the assay.
- Ascorbic acid (up to 20 mg/dl) does not interfere with the assay.
- Protein, as indicated by human albumin concentration (up to 5 g/dl), does not interfere with the assay.
- EDTA-2K (up to 10 mg/ml) does not interfere with the assay.
- Heparin (up to 100 U/ml) does not interfere with the assay.

#### REFERENCES

1. W. Ganong, Review of medical physiology. 16th ed., Appleton-Lange, East Norwalk, CT; 56-60(1993).
2. P. O. Collinson, F. G. Boa, and D. C. Gaze, Measurement of cardiac troponins. *Ann Clin Biochem* 38; 423-449(2001).
3. J. E. Adams III, G. S. Bordor, V. G. Davila-Roman, et al., Cardiac troponin I: A marker with high specificity for cardiac injury. *Circulation* 88(1); 101-107(1993).
4. K. Nakai, K. Nakai, T. Suzuki, et al., Cardiac troponin I assay evaluation and clinical utility. *JJCLA Vol* 28 (3); 184-192(2003).
5. H. B. Wu, Cardiac troponin T and I in coronary artery diseases. *AACC Endo* 13(3); 79-88(1995).
6. T. H. Lee, G. W. Rouan, M. C. Weisberg, et al. Sensitivity of routine clinical criteria for diagnosing myocardial infarction within 24 h of hospitalization. *Ann Intern Med* 106; 181-186(1987).
7. F. M. Fesmire, D. R. Wharton, and F. B. Calhoun, Instability of ST segments in the early stages of acute myocardial infarction in patients undergoing continuous 12 lead ECG monitoring. *Am Journal Emerg Med* 13; 79-88(1995).
8. Ognibene, et al. Cardiac troponin I in myocardial contusion. *Clin Chem Vol* 44(4); 889-890(1998).
9. E. M. Antman, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. *New England J of Med* 335(18); 1342-1349(1996).
10. J. Mair, D. Morandell, N. Gensar, et al. Equivalent early sensitivity of myoglobin, creatine kinase MB mass, creatine kinase isoform ratios, and cardiac troponins I and T for acute myocardial infarction. *Clin Chem* 41; 1266-1272(1995).
11. J. Mair, I. Wagner, B. Puschendorf, et al. Cardiac troponin I to diagnose myocardial injury. *The Lancet* 341; 838-839 (1993).
12. J. P. Bertinchant, C. Larue, I. Pernel, et al. Release kinetics of serum cardiac troponin in ischemic myocardial injury. *Clin Biochemistry* 29; 587-594(1996).

13. F. S. Apple, R. H. Christenson, A. S. Jaffe, et al., National academy of clinical biochemistry and IFCC committee for standardization of markers of cardiac damage laboratory medicine practice guidelines: Analytical issues for biochemical markers of acute coronary syndromes *Clin Chem* 53:4; 547-551(2007).
14. Morrow, C. P. Cannon, R. L. Jesse, et al. National academy of clinical biochemistry laboratory medicine practice guidelines: Clinical characteristics and utilization of biochemical markers in acute coronary syndromes. *Clin Chem* 53:4; 552-574(2007).
15. Lareu, C. Calzolari, J. P. Bertinchant, et al. Cardiac-specific immunoenzymometric assay of troponin I in the early phase of acute myocardial infarction. *Clin Chem* 39; 972-979(1993).
16. J. Bakker, M. J. W. Koelemey, J. P. M. C. Gorgels, et al. Failure of new biochemical markers to exclude acute myocardial infarction at admission. *The Lancet* 342; 1220-1222(1993).
17. H. B. Wu, Y-J. Feng, R. Moore, et al. Characterization of cardiac troponin subunit release into serum after acute myocardial infarction and comparison of assays for troponin T and I. *Clin Chem* 44:6; 1198-1208(1998).
18. S. Eriksson, H. Halenius, K. Pulkki, et al. Negative interference in cardiac troponin I immunoassays by circulating troponin autoantibodies. *Clin Chem* 51:5; 839-847(2005).

#### TOSOH CORPORATION

3-8-2, Shiba, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

#### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

#### TOSOH BIOSCIENCE, INC.

6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080, USA  
Phone: +1 650 615 49 70

Net Volume (after reconstitution  
for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET

### INTENDED USE

The ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET is intended for in vitro diagnostic use only for the calibration of the ST AIA-PACK cTnI 3rd-Gen assay.

### SUMMARY AND EXPLANATION

The ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET contains buffered bovine serum albumin with assigned levels of cardiac troponin I (cTnI). Calibration should be performed in accordance with the schedule given in the Tosoh AIA System Operators' Manual.

### MATERIAL SUPPLIED

Cat. No.	0025315
2 x 1 ml	ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (1) 0 ng/ml Buffered bovine serum albumin containing no detectable concentration of cTnI, with sodium azide as a preservative. (Liquid)
2 x 1 ml	ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (2) 0.2 ng/ml (approx.) ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (3) 1.0 ng/ml (approx.) ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (4) 5.0 ng/ml (approx.) ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (5) 20.0 ng/ml (approx.) ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (6) 60.0 ng/ml (approx.) Buffered bovine serum albumin containing the assigned concentration of cTnI (described on each vial). (Lyophilised)

### WARNINGS AND PRECAUTIONS

1. The ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. The ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (1) contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always rinse with large volumes of water to prevent any azide accumulation.
4. The material derived from human blood used in the preparations of these calibrators has been tested by FDA-approved methods and found negative for the presence of HBsAg and antibodies to HIV-1 and HCV. Since no testing method can give complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions as those used for patient samples.
5. Do not use beyond the expiry date.
6. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

- Bring the calibrators to 18°-25°C for use.

- Using volumetric pipettes, reconstitute the lyophilised calibrators accurately to the volume of 1 ml with CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline.
- Allow the lyophilised material to dissolve fully, then mix the calibrators gently but thoroughly before performing the calibration.

### STORAGE AND STABILITY

- Always store ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET in an upright position at 2°-8°C when not in use.
- When stored unopened and refrigerated at 2°-8°C, the calibrator set will remain stable until the expiry date on the label.
- The calibrators should be used within 1 day of opening or reconstituting, provided the vials are kept tightly sealed and refrigerated at 2°-8°C.

### PROCEDURE

NOTE: Refer to the CALIBRATION PROCEDURE in the insert sheet of the ST AIA-PACK cTnI 3rd-Gen. For additional procedural instructions regarding calibration, refer to the Tosoh AIA System Operators' Manual.

1. When using new calibrator lots, enter the calibrator concentration values and lot number on the software test file. (Refer to the Tosoh AIA System Operators' Manual for details.)
2. Load the appropriate amount of ST AIA-PACK cTnI 3rd-Gen test cups on the instrument.
3. Add the appropriate amount of each calibrator to sample cups. (Refer to the instrument worksheet for the sample volume.)
4. Print a work list and place the sample cups in the position indicated.
5. Select START. Verify that the positions on the work list match the starting position on the screen.

### ASSIGNMENT OF VALUES

The ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET contains assigned concentrations of cardiac troponin I. The assigned value is determined on a lot-by-lot basis and is designed to provide an assay calibration range of 0.02 to 50 ng/ml cTnI. The calibrators in this set are prepared gravimetrically and compared to internal reference standards.

### RESULTS

- 1) The mean rate for the ST AIA-PACK cTnI 3rd-Gen CALIBRATOR (1) should be  $\leq 0.2$  nmol/(L\*s).
- 2) Since there is a direct relationship between concentration and rate, the rate should rise as the concentration increases.
- 3) Repeat values should be within 10%.

### LIMITATIONS

The ST AIA-PACK cTnI 3rd-Gen CALIBRATOR SET is designed solely for use with ST AIA-PACK cTnI 3rd-Gen assay procedures. Although the approximate value of the highest calibrator is 60 ng/ml, the exact concentration may be slightly different. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of assay range, 50 ng/ml.

**TOSOH CORPORATION**

3-8-2, Shiba, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220


**TOSOH EUROPE N.V.**

Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

**TOSOH BIOSCIENCE, INC.**

6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080, USA  
Phone: +1 650 615-4970

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## AIA-PACK® D-Dimer CONTROL SET

### INTENDED USE

The AIA-PACK D-Dimer CONTROL SET is intended for IN VITRO DIAGNOSTIC USE ONLY, by healthcare professionals, for performing quality control procedures with the ST AIA-PACK D-Dimer assay.

### SUMMARY AND EXPLANATION

The AIA-PACK D-Dimer CONTROL SET contains protein matrix with the assigned levels of D-dimer. The D-dimer controls shall be run at the beginning of each day on which D-dimer assays are scheduled.

### CONTENTS

Cat. No.	0025432
2 x 1 mL	AIA-PACK D-Dimer CONTROL LEVEL 1 Protein matrix containing D-dimer (Lyophilised). See the vial label for the assigned concentration range.
2 x 1 mL	AIA-PACK D-Dimer CONTROL LEVEL 2 Protein matrix containing D-dimer (Lyophilised). See the vial label for the assigned concentration range.

### WARNINGS AND PRECAUTIONS

1. The AIA-PACK D-Dimer CONTROL SET is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local Tosoh sales representative.
3. The control material has been tested by FDA-approved methods and found negative for the presence of HBsAg, antibody to HIV-1/2 and HCV. Because no test method can offer complete assurance that products derived from human origin will not transmit infectious agents, it is recommended that this product be handled with the same precautions as used for patient samples.
4. Do not use beyond the expiry date.
5. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

Bring the controls to 18°-25°C for use. Using volumetric pipettes, reconstitute the lyophilised controls accurately to the volume of 1 mL with CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline. Allow the lyophilised material to dissolve fully.

### STORAGE AND STABILITY

- Always store the AIA-PACK D-Dimer CONTROL SET in an upright position at 2°-8°C when not in use.
- When stored unopened and refrigerated at 2°-8°C, the control set will remain stable until the expiry date on the label.
- The controls should be used within 7 days of opening or reconstituting, provided that the vials are kept tightly sealed and refrigerated at 2°-8°C.

### PROCEDURE

Refer to the Tosoh AIA® System Operators' manual or the Tosoh AIA-CL Analyzer Operator's Manual for additional procedural instructions regarding quality control.

1. Load the appropriate amount of test cups on the instrument.
2. Add the appropriate amount of each control to sample cups. (Refer to the instrument worksheet for the sample volume.)
3. Print a work list and place the sample cups in the position indicated.
4. Select START. Verify that the positions on the work list match the starting position on the screen.

### ASSIGNMENT OF VALUES

The AIA-PACK D-Dimer CONTROL SET contains assigned concentration range of D-dimer. The assigned range is determined on a lot-by-lot basis and is designed to provide target control levels of D-dimer. Since the assay values are dependent upon assay procedures as well as several other factors, each laboratory should establish its own range for the assay procedure being monitored.

### LIMITATIONS

The AIA-PACK D-Dimer CONTROL SET is designed solely for use with ST AIA-PACK D-Dimer and CL AIA-PACK D-Dimer assay procedures.

"AIA" and "AIA-PACK" are the registered trademarks of Tosoh Corporation in the European Union, etc.


### TOSOH CORPORATION

3-8-2, Shiba-Koen, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220


### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO, BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

### TOSOH BIOSCIENCE, INC.

6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080, USA  
Phone: +1 650 615 49 70

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in the European Union

## AIA-PACK DILUENT CONCENTRATE

### INTENDED USE

The AIA-PACK DILUENT CONCENTRATE is intended for in vitro diagnostic use only for the AIA-PACK or the ST AIA-PACK Assay.

### MATERIAL SUPPLIED

Cat. No. 0020956  
4 x 100 ml AIA-PACK DILUENT CONCENTRATE  
Buffer solution with detergent.

### WARNINGS AND PRECAUTIONS

1. The AIA-PACK DILUENT CONCENTRATE is intended for in vitro diagnostic use.
2. Inspect the packaging and the exterior of the bottles for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. This material contains sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always rinse with large volumes of water to prevent any azide accumulation.
4. Do not use beyond the expiry date.
5. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

Add the entire content of the AIA-PACK DILUENT CONCENTRATE (100 ml) to approximately 4.0 l of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline, mix well, and adjust the final volume to 5.0 l.

Alternatively, add 20 ml of the AIA-PACK DILUENT CONCENTRATE (as measured with a graduated cylinder) to 800 ml of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline. Mix well, and adjust the volume to 1000 ml with CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline.

### STORAGE AND STABILITY

- Always store the AIA-PACK DILUENT CONCENTRATE in an upright position at 2°-8°C when not in use.
- When stored unopened and refrigerated at 2°-8°C, the AIA-PACK DILUENT CONCENTRATE will remain stable until the expiry date on the label.
- When diluted the diluent solution is stable for 30 days at 18°-25°C.

### PROCEDURE

Please refer to the insert of the AIA-PACK or the ST AIA-PACK and the Operators' Manual of the Tosoh AIA System analysers for detailed instructions.

## SAFETY INFORMATION

### Hazard statements

H412 Harmful to aquatic life with long lasting effects.

### Precautionary statements

P101 If medical advice is needed, have product container or label at hand.

P102 Keep out of reach of children.

P103 Read label before use.

P273 Avoid release to the environment.

P501 Dispose of contents/container in accordance with local/ regional/ national/ international regulations.

### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

**NET** Net Volume (after reconstitution for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## ST AIA-PACK FT4

### INTENDED USE

ST AIA-PACK FT4 is designed for in vitro diagnostic use only for the quantitative measurement of non-protein-bound (free) thyroxine (FT<sub>4</sub>) in **human serum or heparinised plasma** on Tosoh AIA System Analysers.

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### SUMMARY AND EXPLANATION OF TEST

L-thyroxine (3, 5, 3-, 5-L-tetraiodothyronine (T<sub>4</sub>)) produced by the thyroid gland, circulates in the blood 99.97% bound to plasma proteins including thyroxine-binding globulin (TBG), thyroxine binding prealbumin (TBPA) and albumin (1,2). Approximately 0.03% of total circulating thyroxine is unbound. This free T<sub>4</sub> (FT<sub>4</sub>) is believed to be the physiologically active portion of thyroxine which stimulates the metabolism and controls, via the pituitary, the feedback system involving the release of TSH.(3). Historically, measurement of total serum T<sub>4</sub> (bound + free) has been used to assess the clinical status of the thyroid gland (4,5). However, this analysis is not diagnostically accurate when significant changes occur in the serum binding proteins (6,7). Alterations in TBG concentration, pregnancy, oral contraceptives, estrogen therapy or drugs that alter the binding of thyroxine to the carrier proteins may cause corresponding changes in the total T<sub>4</sub> when unbound free thyroxine levels remain relatively unchanged (8,9). Therefore, measurement of the free T<sub>4</sub> (FT<sub>4</sub>) typically correlates more closely to the patient's actual thyroid status than the total (10).

### PRINCIPLE OF THE ASSAY

The ST AIA-PACK FT4 is a competitive enzyme immunoassay performed entirely in the ST AIA-PACK FT4 test cups. The thyroxine not bound to serum proteins (free T<sub>4</sub>) competes with enzyme-labelled T<sub>4</sub> for a limited number of binding sites on a T<sub>4</sub>-specific antibody immobilised on magnetic beads. After incubation, the beads are washed to remove the unbound enzyme-labelled T<sub>4</sub> and are then incubated with a fluorogenic substrate, 4-methylumbelliferyl phosphate (4MUP). The amount of enzyme-labelled T<sub>4</sub> that binds to the beads is inversely proportional to the free T<sub>4</sub> concentration in the test sample. A standard curve using a range of known standard concentrations is constructed and unknown sample free T<sub>4</sub> concentrations are calculated using this curve.

### MATERIAL SUPPLIED

Cat. No.	0025268
ST AIA-PACK FT4	5 trays x 20 test cups

Plastic test cups containing twelve lyophilised magnetic beads with anti-T<sub>4</sub> rabbit polyclonal antibody and 140 µl of thyroxine (T<sub>4</sub>) conjugated to bovine alkaline phosphatase with sodium azide as a preservative.

### MATERIALS REQUIRED BUT NOT SUPPLIED

The following materials are required to perform Free Thyroxine analysis using the ST AIA-PACK FT4 (Cat. No. 0025268) on the Tosoh AIA system Analysers. The following are available separately from Tosoh

Materials	Cat. No.
AIA Nex•IA or AIA-21	0018539
AIA Nex•IA or AIA-21 LA	0018540
AIA-1800 ST	0019836
AIA-1800 LA	0019837
AIA-2000 ST	0022100
AIA-2000 LA	0022101
AIA-600 II	0019014
AIA-600 II BCR	0019328
AIA-900	0022930
AIA-360	0019945
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK SUBSTRATE REAGENT II/ AIA-PACK SUBSTRATE RECONSTITUENT II	
AIA-PACK FT4 CALIBRATOR SET	0020368
AIA-PACK FT4 CALIBRATOR (1)	0 ng/dl
AIA-PACK FT4 CALIBRATOR (2)	0.4 ng/dl (approx.)
AIA-PACK FT4 CALIBRATOR (3)	1.0 ng/dl (approx.)
AIA-PACK FT4 CALIBRATOR (4)	2.0 ng/dl (approx.)
AIA-PACK FT4 CALIBRATOR (5)	4.0 ng/dl (approx.)
AIA-PACK FT4 CALIBRATOR (6)	9.0 ng/dl (approx.)
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
SAMPLE CUPS	0018581
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970
Additional Requirements for AIA Nex•IA/AIA-21 only:	
PIPETTE TIPS	0018552
PRELOADED PIPETTE TIPS	0018583
Additional Requirements for AIA-600 II, AIA-1800, AIA-900, AIA-2000:	
PIPETTE TIPS	0019215
TIP RACK	0019216
PRELOADED PIPETTE TIPS	0022103

Only materials obtained from Tosoh should be used. Materials obtained elsewhere should not be substituted since assay performance is characterised strictly by Tosoh materials.

### WARNINGS AND PRECAUTIONS

1. The ST AIA-PACK FT4 is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the aluminium pouch for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. Test cups from different lots or different assays shall not be mixed within a tray.
4. The ST AIA-PACK FT4 contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always rinse with large volumes of water to prevent any azide accumulation.
5. Human serum is not used in the preparation of this product; however, since human specimens will be used for samples and other quality control products in the lab may be derived from human serum, please use standard laboratory safety procedures in handling all specimens and controls.
6. Do not use beyond the expiry date.
7. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

8. Serum, dust, metal, or microorganism contamination may cause degradation of reconstituted substrate solution. Store in a clean environment, away from direct sunlight and ultraviolet light.
9. TOSOH recommends that a new pouch of the test cups should be used for calibration.

## STORAGE AND STABILITY

All unopened materials will remain stable until the expiry date on the label when stored at the specified temperature.

Materials	Cat. No.
2°-8°C:	
ST AIA-PACK FT4	0025268
AIA-PACK FT4 CALIBRATOR SET	0020368
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
1°-30°C:	
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970

- After opening the aluminium pouch, ST AIA-PACK FT4 test cups can be left on-board of the TOSOH AIA System Analyzers (18-25°C) for a maximum of 1 day (24 hours). When stored overnight at 2-8°C, the test cups can be used for up to 3 days (3 cycles of 8 hours on board and 16 hours in the refrigerator). Once the aluminium pouch is opened, even if the test cups are stored in the refrigerator, they must be used within 30 days.
- AIA-PACK FT4 CALIBRATOR SET must be kept tightly sealed and refrigerated at 2-8°C. After opening or reconstituting, the calibrators should be used within 1 day.
- Reconstituted substrate solution will remain stable for 3 days at 18°-25°C or 30 days at 2°-8°C.
- Working diluent and wash solutions are stable for 30 days at 18°-25°C.
- Reagents should not be used if they appear cloudy or discoloured.

## SPECIMEN COLLECTION AND HANDLING

- Serum or heparinised plasma is required for the assay. EDTA and citrated plasma SHOULD NOT BE USED.
- When using serum, a venous blood sample is collected aseptically without additives. Store at 18°-25°C until a clot has formed (usually 15-45 minutes), then centrifuge to obtain the serum specimen for assay.
- When using heparinised plasma, a venous blood sample is collected aseptically with the specified additive. Centrifuge and separate plasma from the packed cells as quickly as possible.
- Inadequate centrifugation or the presence of fibrin or particulate matter in the sample may cause an erroneous result.
- Samples containing inhibitors of alkaline phosphatase may cause erroneous results.
- Inspect all samples for air bubbles and foaming. Remove any air bubbles before assay.
- Specimen types should not be used interchangeably during serial monitoring of an individual patient. Measured concentrations may vary slightly between sample types in certain patients.
- **Samples may be stored at 2°-8°C for up to 7 days before analysis. If the analysis cannot be performed within 7 days, the sample should be stored frozen at -20°C or below for up to 60 days.**
- Repeated freeze-thaw cycles should be avoided. Turbid serum samples or samples containing particulate matter should be centrifuged before testing. Before assaying, slowly bring frozen samples to 18°-25°C and mix gently.
- The sample required for analysis is 10 µl.

## PROCEDURE

For the AIA Nex•IA / AIA-21, AIA-600 II, AIA-1800, AIA-900, AIA-2000 and AIA-360, please refer to their Operators' Manual for detailed instructions.

### I. Reagent Preparation

#### A. Substrate Solution

Bring all reagents to 18°-25°C before preparing the working substrate solution. Add one bottle of AIA-PACK SUBSTRATE RECONSTITUENT II to the lyophilised AIA-PACK SUBSTRATE REAGENT II, mix thoroughly to dissolve all solid materials.

#### B. Wash Solution

Add the entire contents of the AIA-PACK WASH CONCENTRATE (100 ml) to approximately 2.0 l of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline, mix well, and adjust the final volume to 2.5 l.

#### C. Diluent

Add the entire contents of the AIA-PACK DILUENT CONCENTRATE (100 ml) to approximately 4.0 l of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline, mix well, and adjust the final volume to 5.0 l.

### II. Calibration Procedure


#### A. Calibration Curve

The calibrators for use with the ST AIA-PACK FT4 are prepared gravimetrically and are compared against internal reference standards.

The calibration curve for the ST AIA-PACK FT4 will remain stable for up to 90 days. Calibration stability is monitored by quality control performance and is dependent on proper reagent handling and Tosoh AIA System maintenance in accordance with the manufacturer's instructions.

Recalibration may be necessary more frequently if controls are out of the established range for this assay or when certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of the wash probe or detector lamp adjustment or change). For further information regarding instrument operation, consult the TOSOH AIA System Operator's Manual.

A sample calibration curve from the AIA-1800 follows and shows the algorithm used to calculate results.

## B. Calibration Procedure

1. Refer to the appropriate Tosoh AIA System Operators' Manual for the procedural instructions.
2. Verify that both the calibrator lot and concentration numbers have been correctly entered on the software.
3. The AIA-PACK FT4 CALIBRATOR SET is lyophilised. All levels should be reconstituted with 1.0 ml of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline.
4. Tosoh recommends that all calibrators should be run in triplicate.

## C. Calibration Acceptability Criteria

1. Since there is an inverse relationship between concentration and rate, the rate should decrease as the concentration increases.
2. The replicate values should be within a 10% range.

## D. Calibration Review and Acceptance

1. Review the calibration curve carefully, applying the criteria listed above.
2. Edit the calibration if necessary, then accept the calibration.

For further information regarding calibration, consult the Tosoh AIA System Operators' Manual.

## III. Quality Control Procedure

### A. Commercially Available Controls

Commercially available controls shall be run at least once per day. It is recommended that at least two (2) levels of controls, normal and abnormal, be used. Laboratory policy for this particular assay specifies the following:

Control Material: \_\_\_\_\_

Frequency: \_\_\_\_\_

Lot number of control material, acceptable limits, and corrective action to be taken if controls do not meet laboratory criteria will be found in a separate quality control document maintained by the laboratory.

### B. Quality Control Procedure

1. Assay quality control specimens as instructed in the specific Operators' Manual for your analyser. In addition, refer to the Tosoh AIA System Operators' Manual for detailed instructions on defining and editing files.
2. Quality control material to be run with this assay is defined by individual laboratory policy.

## IV. Specimen Processing

### A. Preparation

Following the specific instructions in the Operators' Manual for the analyser, place samples on the instrument appropriately. Bar-coded primary tubes as well as sample cups can be run on the AIA Nex•IA / AIA-21, AIA-600 II, AIA-1800, AIA-900, AIA-2000 and AIA-360.

### B. Assay Procedure

1. Provide a sufficient quantity of ST AIA-PACK FT4 test cups for the number of samples to be run.
2. Load patient samples as instructed in the Operators' Manual and proceed with analysis.

## PROCEDURAL NOTES

1. Lyophilised substrate must be completely dissolved.
2. Ligand assays performed by the TOSOH AIA System Analyzers require that the laboratory use water designated by the CAP as Class I or by the CLSI as the clinical laboratory reagent water. Water

should be tested at least once per month and should be free of particulate matter including bacteria. The pH of the water should also be routinely tested. For further information, consult the CLSI document C3-A4, Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline-Fourth Edition.

3. If a specimen free thyroxine concentration is found to be greater than the upper limit of the assay range 8.0 ng/dl, results should be reported as greater than 8.0 ng/dl. Dilution or recovery studies cannot be performed on the Free Thyroxine assay. The equilibrium between bound and free hormone in the Free Thyroxine assay is altered both by dilutions and the addition of measured amount of thyroxine.
4. The Tosoh AIA System Analysers can store two different calibration curves for each analyte at one time. Thus, up to two different lots of ST AIA-PACK FT4 test cups can be used during the same run.
5. If the assay specifications for this test are not ready on the system software, the specifications must be entered under test code **075**.

## CALCULATION OF RESULTS

The Tosoh AIA System Analysers perform all sample and reagent handling operations automatically. The Tosoh AIA System Analysers read the rate of fluorescence produced by the reaction and automatically convert the rate to Free Thyroxine concentration in ng/dl.

## EVALUATION OF RESULTS

### Quality Control

In order to monitor and evaluate the precision of the analytical performance, it is recommended that commercially available control samples shall be assayed in accordance with local regulations.

Minimum recommendations for the frequency of running internal control material are:

After calibration, three levels of the internal control are run in order to accept the calibration curve. The three levels of controls are repeated when certain service procedures are performed (e.g. temperature adjustments, sampling mechanism changes, maintenance of the wash probe or detector lamp adjustments or changes).

After daily maintenance, at least two levels of the control shall be run in order to verify overall performance of the Tosoh AIA System Analysers.

If one or more control sample value(s) is out of the acceptable range, it is necessary to investigate the validity of the calibration curve before reporting patient results.

Standard laboratory procedures should be followed in accordance with the strict regulatory agency under which the laboratory operates.

## LIMITATIONS

- For diagnostic purposes, results obtained from this assay should be used in conjunction with other data (e.g. symptoms, results of other tests, clinical impressions, therapy, etc.).
- Using ST AIA-PACK FT4, the highest measurable concentration of Free Thyroxine in specimens is 8.0 ng/dl, and the lowest measurable concentration in specimens is 0.10 ng/dl (assay sensitivity).
- Although the approximate value of the highest calibrator is 9.0 ng/dl, the exact concentration may be slightly different. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 8.0 ng/dl.
- Although haemolysis has an insignificant effect on the assay, haemolysed samples may indicate mistreatment of a specimen prior to assay and results should be interpreted with caution.
- Lipaemia has an insignificant effect on the assay except in the case of gross lipaemia where spatial interference may occur.
- Circulating autoantibodies to T<sub>4</sub> may interfere with free thyroxine measurements.

- Hormone-binding inhibitors may interfere with the free T<sub>4</sub> assay. Heparin may cause in vivo effects in free T<sub>4</sub> assays. Blood collection for free T<sub>4</sub> assays should not be performed concurrent with administration of heparin therapy.
- Drugs that affect the binding of T<sub>4</sub> to the thyroid hormone carrier proteins or affect the metabolism of T<sub>4</sub> to T<sub>3</sub> may complicate the interpretation of free thyroxine results. In patients with severe non-thyroidal illness (NTI) in whom the free T<sub>4</sub> yields results indicating hypothyroidism, it is suggested that the ST AIA-PACK TSH assay be run to confirm this diagnosis.
- Specimens from patients taking medicines and/or medical treatment may show erroneous results.
- Heterophile antibodies are known to interfere in assays of this type. ST AIA-PACK FT<sub>4</sub> has been designed to minimize the effects of heterophile antibodies, but interference from high titers cannot be ruled out.
- Samples from patients under Asfotase Alfa (Genetical Recombination) treatment may cause falsely elevated/decreased results.
- For a fuller understanding of the limitations of this procedure, please refer to the SPECIMEN COLLECTION AND HANDLING, WARNINGS AND PRECAUTIONS, STORAGE AND STABILITY, and PROCEDURAL NOTES sections in this insert sheet.

## EXPECTED VALUES

Each laboratory should determine a reference interval corresponding to the characteristics of the population being tested. As with all diagnostic procedures, clinical results must be interpreted having due regard to concomitant medications administered to the patient (11)

### I. Reference Ranges

The intervals given here were determined in serum samples from 618 apparently healthy Asian individuals and 130 apparently healthy Caucasian individuals, respectively.

#### • Asian

Number of Samples (n)	618
Reference Interval	0.82 - 1.63 ng/dL (10.6 - 21.0 pmol/L)

#### • Caucasian

Number of Samples (n)	130
Reference Interval	0.61 - 1.57 ng/dL (7.87 - 20.3 pmol/L)

### II. Conversion Factors

FT<sub>4</sub> concentrations in this application are in units of ng/dl. The obtained values in ng/dL must be multiplied by 12.9 when converting to values in pmol/L.

## PERFORMANCE CHARACTERISTICS

### PRECISION

- Within-run precision was determined using four controls in a total of 20 runs. Within each run, one set of duplicates per control was assayed. The mean of each duplicate was used to obtain the pooled standard deviation (SD), which was then used to calculate the coefficient of variation (CV).

Sample	Mean (ng/dl)	Pooled SD (ng/dl)	CV (%)
Serum A	0.599	0.0270	4.5
Serum B	2.00	0.0607	3.0
Serum C	4.61	0.0768	1.7
Plasma A	1.11	0.0348	3.1

- Total precision was determined by the duplicate assay of four controls in 20 separate runs. The means of each run were used to calculate the pooled standard deviation (SD) and coefficient of variation (CV).

Sample	Mean (ng/dl)	Pooled SD (ng/dl)	CV (%)
Serum A	0.599	0.0341	5.7
Serum B	2.00	0.0720	3.6
Serum C	4.61	0.123	2.7
Plasma A	1.11	0.0442	4.0

### CORRELATION

The correlation between serum (x) and heparinised plasma (y) on ST AIA-PACK FT<sub>4</sub> was carried out using 130 patient specimens.

Slope	0.999
y-Intercept	0.004
Correlation Coefficient	0.978
Range of Samples	0.37-2.90
Number of Samples	130

### SPECIFICITY

The following substances were tested for cross-reactivity. Cross-reactivity is expressed in terms of the percentage of concentration of each substance which will be measured as thyroxine.

Compound	Cross-reactivity (%)
Thyroxine (T <sub>4</sub> )	100
Triiodothyronine (T <sub>3</sub> )	2.0
3, 3', 5-triiodothyropropionic acid	1.6
3, 5-diiodothyropropionic acid	< 0.01
3, 5-diiodo-L-tyrosine	< 0.01

### SENSITIVITY

The minimum detectable concentration (MDC) of Free Thyroxine is estimated to be 0.10 ng/dl. MDC is defined as the concentration of FT<sub>4</sub> that corresponds to the rate of fluorescence which is two standard deviations from the mean rate of fluorescence of 5 replicate determinations by AIA-PACK FT<sub>4</sub> CALIBRATOR (1).

### INTERFERENCE

Interference is defined for the purposes of this study with recovery outside of 10% of the known concentration of the specimen after the following substances are added to human specimens.

- Haemoglobin (up to 410 mg/dl), free bilirubin (up to 16 mg/dl) and conjugated bilirubin (up to 17 mg/dl) do not interfere with the assay.
- Lipaemia, as indicated by triglyceride concentration (up to 1,600 mg/dl), does not interfere with the assay.
- Ascorbic acid (up to 20 mg/dl) does not interfere with the assay.
- Heparin (up to 100 U/ml) does not interfere with the assay.

## REFERENCES

1. Oppenheimer, J. H., Role of Plasma Proteins in the Binding, Distribution and Metabolism of the Thyroid Hormones. *New Eng. J. Med.*, 278:1153 (1986).
2. Prince, M. P. and Ramsden, D. B., A New Theoretical Description of the Binding of Thyroid Hormones by Serum Proteins. *Clin. Endocrinol.* 7 :307 (1977).
3. Chopra, I. J., et al, Pathway of Metabolism of Thyroid Hormones. *Rec. Prog. Horm. Res.*, 34: 521 (1978).
4. Ekins, R. P., Free Hormones in Blood: The Concept and the Measurement. *J. Clin. Immunoassay*, 7:163 (1984).
5. Penney, M. D. and O'Sullivan, D. J., Total or Free Thyroxine as a Primary Test of Thyroid Function. *Clin.Chem.* 33:170 (1987).
6. Robbins, J. and Rall, J. E., Hormones in Blood, 1:577 (1967).
7. Ramsden, D. B., et al, Serum Free Thyroxine Concentrations in Normal Euthyroid Subjects and Ones with High Serum Thyroxine Binding Globulin Concentration., *Clin. Chem. Acta*, 130:211 (1983).
8. Wilke, T. J., Diagnostic Value of Three Methods for Assessing Free Thyroxin in Pregnancy. *Ann. Clin. Biochem.* 20:60 (1983).
9. Sparre, L. S., et al, Oestrogen and Thyroxine Binding Globulin Levels in Early Normal Pregnancy. *Acta Endocrinol.* 114:298 (1987).
10. Anderson, B. G., Free Thyroxine in Serum in Relation to Thyroid Function, *JAMA*, 203:577 (1968).
11. Young, D., Effects of Drugs on Clinical Laboratory Tests. 3rd Edition, Washington, DC, American Association for Clinical Chemistry Press (1990).

### TOSOH CORPORATION

3-8-2, Shiba, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220


### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

### TOSOH BIOSCIENCE, INC.

6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080, USA  
Phone: +1 650 615 49 70

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## AIA-PACK FT4 CALIBRATOR SET

### INTENDED USE

The AIA-PACK FT4 CALIBRATOR SET is intended for in vitro diagnostic use only for the calibration of the ST AIA-PACK FT4.

### SUMMARY AND EXPLANATION

The AIA-PACK FT4 CALIBRATOR SET contains human sera with assigned levels of free thyroxine (FT<sub>4</sub>). Calibration should be performed in accordance with the schedule given in the Tosoh AIA System Operators' Manual.

### MATERIAL SUPPLIED

Cat. No.	0020368		
2 x 1 ml	AIA-PACK FT4 CALIBRATOR (1)	0	ng/dl
	Human serum containing no detectable concentration of FT <sub>4</sub> with sodium azide as a preservative. (Lyophilised)		
2 x 1 ml	AIA-PACK FT4 CALIBRATOR (2)	0.4	ng/dl (approx.)
	AIA-PACK FT4 CALIBRATOR (3)	1.0	ng/dl (approx.)
	AIA-PACK FT4 CALIBRATOR (4)	2.5	ng/dl (approx.)
	AIA-PACK FT4 CALIBRATOR (5)	4.0	ng/dl (approx.)
	AIA-PACK FT4 CALIBRATOR (6)	9.0	ng/dl (approx.)
	Human serum containing the assigned concentration of FT <sub>4</sub> (described on each vial) with sodium azide as a preservative. (Lyophilised)		

### WARNINGS AND PRECAUTIONS

1. The AIA-PACK FT4 CALIBRATOR SET is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. This material contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always rinse with large volumes of water to prevent any azide accumulation.
4. The material derived from human blood used in the preparations of these calibrators has been tested by FDA-approved methods and found negative for the presence of HBsAg and antibodies to HIV-1 and HCV. Since no testing method can give complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions as those used for patient samples.
5. Do not use beyond the expiry date.
6. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

Using volumetric pipettes, reconstitute the lyophilised calibrators accurately to the volume of 1 ml with CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI

**C3-A4 guideline.** Allow the lyophilised material to dissolve fully, then mix the calibrators gently but thoroughly before performing the calibration. Bring the calibrators to 18°-25°C for use.

### STORAGE AND STABILITY

- Always store the AIA-PACK FT4 CALIBRATOR SET in an upright position at 2°-8°C when not in use.
- When stored unopened and refrigerated at 2°-8°C, the calibrator set will remain stable until the expiry date on the label.
- The calibrators should be used within 1 day of opening or reconstituting, provided the vials are kept tightly sealed and refrigerated at 2°-8°C.

### PROCEDURE

NOTE: Refer to the CALIBRATION PROCEDURE in the insert sheet of the ST AIA-PACK FT4. For additional procedural instructions regarding calibration, refer to the Tosoh AIA System Operators' Manual.

1. When using new calibrator lots, enter the calibrator concentration values and lot number on the software test file. (Refer to Tosoh AIA System Operators' Manual for details.)
2. Load the appropriate amount of ST AIA-PACK FT4 test cups on the instrument.
3. Add the appropriate amount of each calibrator to sample cups. (Refer to the instrument worksheet for the sample volume.)
4. Print a work list and place the sample cups in the position indicated.
5. Select START. Verify that the positions on the work list match the starting position on the screen.

### ASSIGNMENT OF VALUES

The AIA-PACK FT4 CALIBRATOR SET contains assigned concentrations of FT<sub>4</sub>. The assigned value is determined on a lot-by-lot basis and is designed to provide an assay calibration range of approximately 0.1 to 8.0 ng/dl of FT<sub>4</sub>. The calibrators in this set are prepared gravimetrically and compared to internal reference standards.

### RESULTS

1. Since there is an inverse relationship between concentration and rate, the rate should decrease as the concentration increases.
2. Repeat values should be within 10% range.

### LIMITATIONS

The AIA-PACK FT4 CALIBRATOR SET is designed solely for use with ST AIA-PACK FT4 assay procedures. Although the approximate value of the highest calibrator is 9.0 ng/dl, the exact concentration may vary slightly. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 8.0 ng/dl.

### TOSOH CORPORATION

3-8-2, Shiba, Minato-ku  
TOKYO 105-8623 (JAPAN)

Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

**TOSOH EUROPE N.V.**

Transportstraat 4

B-3980 TESSENDERLO BELGIUM

Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

EC	REP
----	-----

**TOSOH BIOSCIENCE, INC.**

6000 Shoreline Court, Suite 101

South San Francisco, CA 94080, USA

Phone: +1 650 615 49 70

<b>NET</b>	Net Volume (after reconstitution for lyophilised material)
------------	---

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in the European Union

## AIA-PACK MULTI ANALYTE CONTROL MAC

### INTENDED USE

The AIA-PACK MULTI ANALYTE CONTROL is intended for in vitro diagnostic use only for performing quality control procedures, by healthcare professionals, with the Tosoh Automated Immunoassay Analyzers for the analytes listed in the EXTRA INFO AIA-PACK MULTI ANALYTE CONTROL.

### SUMMARY AND EXPLANATION

The AIA-PACK MULTI ANALYTE CONTROL is a stable, lyophilised control serum based on human serum. The performance of the analyser can be checked with the three significant levels for Internal Quality Control (IQC) provided in this product. The AIA-PACK MULTI ANALYTE CONTROL is prepared from pooled human serum with the addition of components of human origin or pure chemicals to achieve levels of clinical significance.

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### WARNINGS AND PRECAUTIONS

- The AIA-PACK MULTI ANALYTE CONTROL is intended for in vitro diagnostic use.
- The control material has been tested and found negative for the presence of HBsAg and antibodies to HIV-1, HIV-2 and HCV. Since no testing method can give complete assurance that products derived from human blood will not transmit infectious agents, this product shall be handled with the same precautions as those used for patient samples.
- Do not use beyond the expiry date.
- For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.
- Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local Tosoh sales representative.

### PREPARATION OF REAGENTS

- Bring the controls to 18° -25 °C for use.
- Remove the screw cap and the rubber stopper carefully in order to prevent the loss of substance when the pressure is equilibrated.
- Add exactly 3.0 ml of reagent water defined by CLSI GP40-A4-AMD guideline using volumetric pipettes, put the stopper back on the vial and allow the reconstituted material to stand for at least 30 minutes at +18°C to +25°C.
- After thoroughly mixing (do not shake) the control serum is ready for use.

### STORAGE AND STABILITY

- Unopened vials of multi analyte control are stable at 2°-8°C up to the expiry date printed on the label.
- The reconstituted control materials remain stable for 7 days, provided that the vials are kept tightly sealed and refrigerated at 2°-8°C. Since some of the analytes have a limited stability one shall check the note of the EXTRA INFO AIA-PACK MULTI ANALYTE CONTROL.
- Freshly reconstituted control sera can be frozen once in aliquots of not less than 0.5 ml at -20°C to -80°C. The frozen aliquots are stable for 4 weeks.

### PROCEDURE

Set up the AIA / AIA-CL analyzer for analysis. Refer to the AIA Analyzer Operators' Manual.

This product should be treated the same as patient specimens and run in accordance with the instructions accompanying the reagent being used.

1. Load the appropriate amount of AIA / CL AIA-PACK test cups on the analyser.
2. Add the appropriate amount of each control to sample cups. (Refer to the analyser work list for the sample volume.)
3. Print a work list and place the sample cups in the position indicated.
4. Select START. Verify that the positions on the work list match the starting position on the screen.

### ASSIGNMENT OF VALUES

1. The confidence limits are the assigned values plus/minus two times the overall standard deviation of the value assignment.
2. Values are provided by the EXTRA INFO AIA-PACK MULTI ANALYTE CONTROL on the webpage: <https://www.diagnostics.eu.tosohbioscience.com/service-support/instructions-for-use>
3. It is recommended that each laboratory establish its own means and acceptable ranges and use the provided values only as guides.
4. Please be aware that the assigned values available in this EXTRA INFO AIA-PACK MULTI ANALYTE CONTROL are likely to get updated during the lifetime of the product. Therefore the assigned values must be regularly consulted on the webpage: <https://www.diagnostics.eu.tosohbioscience.com/service-support/instructions-for-use>

### LIMITATIONS

The AIA-PACK MULTI ANALYTE CONTROL is designed solely for use with the Tosoh Automated Immunoassay Analyzers.

1. This product should not be used past the expiry date.
2. If there is evidence of microbial contamination or excessive turbidity in the reconstituted product, discard the vial.
3. This product is not intended for use as a standard.

0197

### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO (BELGIUM)  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

**NET** Net Volume (after reconstitution for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## ST AIA-PACK BNP CALIBRATOR SET

### INTENDED USE

The ST AIA-PACK BNP CALIBRATOR SET is intended for IN VITRO DIAGNOSTIC USE ONLY for the calibration of the ST AIA-PACK BNP assay.

### SUMMARY AND EXPLANATION

The ST AIA-PACK BNP CALIBRATOR SET contains buffered bovine serum albumin with assigned levels of brain natriuretic peptide (BNP). Calibration should be performed according to the schedule indicated in the TOSOH AIA System Operators' manual.

### DESCRIPTION

Cat. No.	0025328		
2 x 1 ml	ST AIA-PACK BNP CALIBRATOR (1)	0	pg/ml
	Buffered bovine serum albumin containing no detectable concentration of BNP with sodium azide as a preservative (Liquid).		
2 x 1 ml	ST AIA-PACK BNP CALIBRATOR (2)	15	pg/ml (approx.)
	ST AIA-PACK BNP CALIBRATOR (3)	40	pg/ml (approx.)
	ST AIA-PACK BNP CALIBRATOR (4)	150	pg/ml (approx.)
	ST AIA-PACK BNP CALIBRATOR (5)	600	pg/ml (approx.)
	ST AIA-PACK BNP CALIBRATOR (6)	2,300	pg/ml (approx.)
	Buffered bovine serum albumin containing the assigned concentration of BNP (described on each vial) with sodium azide as a preservative (Lyophilised).		

### WARNINGS AND PRECAUTIONS

1. The ST AIA-PACK BNP CALIBRATOR SET is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. This material contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
4. Although material derived from human origin is not used for these calibrators, it is recommended that this product be handled with the same precautions as used for patient samples.
5. Do not use beyond the expiry date.
6. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

Bring the calibrators to 18°-25°C for use. Using volumetric pipettes, reconstitute the lyophilised calibrators accurately to the volume of 1 ml with CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline. Allow the lyophilised material to dissolve fully, then mix the calibrators gently but thoroughly before performing the calibration.

### STORAGE AND STABILITY

- Always store the ST AIA-PACK BNP CALIBRATOR SET in an upright position at 2°-8°C when not in use.
- When stored unopened and refrigerated at 2°-8°C, the calibrator set will remain stable until the expiry date on the label.
- The calibrators should be used within 1 day of opening or reconstituting, provided the vials are kept tightly sealed and refrigerated at 2°-8°C.

### PROCEDURE

NOTE: Refer to the CALIBRATION PROCEDURE in the insert sheet of the ST AIA-PACK BNP. For additional procedural instructions regarding calibration, refer to the TOSOH AIA System Operator's Manual.

1. When using new calibrator lots, enter the calibrator concentration values and lot number into the software test file. (Refer to the TOSOH AIA System Operator's Manual for details.)
2. Load the appropriate amount of ST AIA-PACK BNP test cups on the instrument.
3. Add the appropriate amount of each calibrator to sample cups. (Refer to the instrument worksheet for the sample volume.)
4. Print a work list and place the sample cups in the position indicated.
5. Select START. Verify that the positions on the work list match the starting position on the screen.

### ASSIGNMENT OF VALUES

The ST AIA-PACK BNP CALIBRATOR SET contains assigned concentrations of BNP. The assigned value is determined on a lot-by-lot basis and is designed to provide an assay calibration range of 4.0 to 2 000 pg/ml of BNP. The calibrators in this set are prepared gravimetrically and are compared to internal reference standards.

### RESULTS

1. The mean rate for the ST AIA-PACK BNP CALIBRATOR (1) should be < 2.5 nmol/(L\*s).
2. Since there is a direct relationship between concentration and rate, the rate should increase as the concentration increases.
3. The replicate values should be within a 10% range.

### LIMITATIONS

The ST AIA-PACK BNP CALIBRATOR SET is designed solely for use with ST AIA-PACK BNP assay procedures. Although the approximate value of the highest calibrator is 2 300 pg/ml, the exact concentration may be slightly different. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 2 000 pg/ml.


### TOSOH CORPORATION

3-8-2, Shiba-Koen, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO, BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

NET

Net Volume (after reconstitution for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

**TOSOH BIOSCIENCE, INC.**  
6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080 (USA)  
Phone: +1 650 615-4970

## **AIA-PACK DETECTOR STANDARDIZATION TEST CUP**

---

### **INTENDED USE**

The AIA-PACK DETECTOR STANDARDIZATION TEST CUP is intended for in vitro diagnostic use only. The AIA-PACK DETECTOR STANDARDIZATION TEST CUP is used to measure the substrate background on the TOSOH AIA System Analysers.

### **MATERIAL PROVIDED**

Cat. No. 0020970  
10 Trays x 20 cups AIA-PACK DETECTOR STANDARDIZATION TEST CUP

### **WARNINGS AND PRECAUTIONS**

The AIA-PACK DETECTOR STANDARDIZATION TEST CUP is intended for in vitro diagnostic use only.

### **PREPARATION**

The AIA-PACK DETECTOR STANDARDIZATION TEST CUP is provided ready for use.

### **STORAGE AND STABILITY**

Store the AIA-PACK DETECTOR STANDARDIZATION TEST CUP at 1°– 30°C. The AIA-PACK DETECTOR STANDARDIZATION TEST CUP will remain stable until the expiry date on the label.

### **PROCEDURAL NOTE**

For additional procedural instructions regarding substrate background measurement, refer to the TOSOH AIA System Operators' Manual.

1. Load one (1) AIA-PACK DETECTOR STANDARDIZATION TEST CUP in the appropriate position on the instrument.
2. Place the bottle of the reconstituted AIA-PACK SUBSTRATE REAGENT II in the appropriate position on the instrument.
3. Select the OK button in the material setting dialogue of the daily maintenance to start the automatic measurement of substrate background.

### **TOSOH CORPORATION**

3-8-2, Shiba, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

### **TOSOH EUROPE N.V.**

Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## ST AIA-PACK D-Dimer CALIBRATOR SET

### INTENDED USE

The ST AIA-PACK D-Dimer CALIBRATOR SET is intended for IN VITRO DIAGNOSTIC USE ONLY for the calibration of the ST AIA-PACK D-Dimer assay.

### SUMMARY AND EXPLANATION

The ST AIA-PACK D-Dimer CALIBRATOR SET contains protein matrix with assigned levels of D-dimer. Calibration should be performed according to the schedule indicated in the TOSOH AIA System Operators' manual.

### DESCRIPTION

Cat. No.	0025332		
2 x 1 mL	ST AIA-PACK D-Dimer CALIBRATOR (1)	0	µg/mL
	Protein matrix containing no detectable concentration of D-dimer with sodium azide as a preservative (Liquid).		
2 x 1 mL	ST AIA-PACK D-Dimer CALIBRATOR (2)	0.125	µg/mL (approx.)
	ST AIA-PACK D-Dimer CALIBRATOR (3)	0.5	µg/mL (approx.)
	ST AIA-PACK D-Dimer CALIBRATOR (4)	2.0	µg/mL (approx.)
	ST AIA-PACK D-Dimer CALIBRATOR (5)	8.0	µg/mL (approx.)
	ST AIA-PACK D-Dimer CALIBRATOR (6)	18.4	µg/mL (approx.)
	Protein matrix containing the assigned concentration of D-dimer (described on each vial) (Lyophilised).		

### WARNINGS AND PRECAUTIONS

1. The ST AIA-PACK D-Dimer CALIBRATOR SET is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. The ST AIA-PACK D-Dimer CALIBRATOR (1) contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
4. The calibrator material has been tested by FDA-approved methods and found negative for the presence of HBsAg, antibody to HIV-1/2 and HCV. Because no test method can offer complete assurance that products derived from human origin will not transmit infectious agents, it is recommended that this product be handled with the same precautions as used for patient samples.
5. Do not use beyond the expiry date.
6. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

Bring the calibrators to 18°-25°C for use. **Using volumetric pipettes, reconstitute the lyophilised calibrators accurately to the volume of 1 mL with CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline.** Allow the lyophilised material to dissolve fully, then mix the calibrators gently but thoroughly before performing the calibration.

### STORAGE AND STABILITY

- Always store the ST AIA-PACK D-Dimer CALIBRATOR SET in an upright position at 2°-8°C when not in use.
- When stored unopened and refrigerated at 2°-8°C, the calibrator set will remain stable until the expiry date on the label.
- The calibrators should be used within 1 day of opening or reconstituting, provided the vials are kept tightly sealed and refrigerated at 2°-8°C.

### PROCEDURE

NOTE: Refer to the CALIBRATION PROCEDURE in the insert sheet of the ST AIA-PACK D-Dimer. For additional procedural instructions regarding calibration, refer to the TOSOH AIA System Operators' manual.

1. When using new calibrator lots, enter the calibrator concentration values and lot number into the software test file. (Refer to the TOSOH AIA System Operators' manual for details.)
2. Load the appropriate amount of ST AIA-PACK D-Dimer test cups on the instrument.
3. Add the appropriate amount of each calibrator to sample cups. (Refer to the instrument worksheet for the sample volume.)
4. Print a work list and place the sample cups in the position indicated.
5. Select START. Verify that the positions on the work list match the starting position on the screen.

### ASSIGNMENT OF VALUES

The ST AIA-PACK D-Dimer CALIBRATOR SET contains assigned concentrations of D-dimer. The assigned value is determined on a lot-by-lot basis and is designed to provide an assay calibration range of 0.02 to 16.0 µg/mL of D-dimer. The calibrators for use with ST AIA-PACK D-Dimer are prepared gravimetrically and are compared to internal reference standards.

### RESULTS

1. The mean rate for the ST AIA-PACK D-Dimer CALIBRATOR (1) should be < 0.2 nmol/(L\*s).
2. Since there is a direct relationship between concentration and rate, the rate should increase as the concentration increases.
3. The replicate values should be within a 10% range.

### LIMITATIONS

The ST AIA-PACK D-Dimer CALIBRATOR SET is designed solely for use with ST AIA-PACK D-Dimer assay procedures. Although the approximate value of the highest calibrator is 18.4 µg/mL, the exact concentration may be slightly different. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 16.0 µg/mL.


### TOSOH CORPORATION

3-8-2, Shiba-Koen, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO, BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

**TOSOH BIOSCIENCE, INC.**

6000 Shoreline Court, Suite 101

South San Francisco, CA 94080, USA

Phone: +1 650 615 49 70

**NET**

Net Volume (after reconstitution  
for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## ST AIA-PACK D-Dimer

For Quantitative Measurement of D-dimer in Citrated plasma.

### NAME AND INTENDED USE

ST AIA-PACK D-Dimer is designed for IN VITRO DIAGNOSTIC USE ONLY for the quantitative measurement of D-dimer in human citrated plasma on TOSOH AIA System Analyzers.

### SUMMARY AND EXPLANATION OF TEST

D-dimer is a degradation product of cross-linked fibrin split by plasmin, which is used as a marker reflecting secondary fibrinolysis.

D-dimer is used in the differential diagnosis of disease states of primary hyperfibrinolysis and secondary hyperfibrinolysis. The measurement of D-dimer is especially helpful in the diagnosis and therapeutic follow-up of DIC (disseminated intravascular coagulation) accompanied by secondary hyperfibrinolysis (1–3) and various thrombosis, and helpful in monitoring of fibrinolytic therapy, etc (4, 5).

### PRINCIPLE OF THE ASSAY

The ST AIA-PACK D-Dimer is a two-site immunoenzymometric assay which is performed entirely in the ST AIA-PACK D-Dimer test cups. D-dimer present in the test sample is bound with monoclonal antibody immobilized on a magnetic solid phase and enzyme-labelled monoclonal antibody in the test cups. The magnetic beads are washed to remove unbound enzyme-labelled antibody and are then incubated with a fluorogenic substrate, 4-methylumbelliferyl phosphate (4MUP). The amount of enzyme-labelled monoclonal antibody that binds to the beads is directly proportional to the D-dimer concentration in the test sample. A standard curve is constructed, and unknown sample concentrations are calculated using this curve.

### MATERIAL PROVIDED

Cat. No. 0025232  
ST AIA-PACK D-Dimer 5 trays x 20 test cups

Plastic test cups containing lyophilised twelve magnetic beads coated with anti-D-dimer mouse monoclonal antibody and 100 µL of anti-D-dimer mouse monoclonal antibody conjugated to bovine alkaline phosphatase with sodium azide as a preservative.

### MATERIALS REQUIRED BUT NOT PROVIDED

The following materials are required to perform D-dimer analysis using the ST AIA-PACK D-Dimer (Cat. No. 0025232) on the TOSOH AIA System Analyzers. The following are available separately from TOSOH.

Materials	Cat. No.
AIA-1800 ST	0019836
AIA-1800 LA	0019837
AIA-2000 ST	0022100
AIA-2000 LA	0022101
AIA-600 II	0019014

AIA-600 II BCR			0019328
AIA-900			0022930
AIA-360			0019945
AIA-PACK SUBSTRATE SET II			0020968
AIA-PACK SUBSTRATE REAGENT II / AIA-PACK SUBSTRATE RECONSTITUENT II			
ST AIA-PACK D-Dimer CALIBRATOR SET			0025332
ST AIA-PACK D-Dimer CALIBRATOR (1)	0	µg/mL	
ST AIA-PACK D-Dimer CALIBRATOR (2)	0.125	µg /mL (approx.)	
ST AIA-PACK D-Dimer CALIBRATOR (3)	0.5	µg /mL (approx.)	
ST AIA-PACK D-Dimer CALIBRATOR (4)	2.0	µg /mL (approx.)	
ST AIA-PACK D-Dimer CALIBRATOR (5)	8.0	µg /mL (approx.)	
ST AIA-PACK D-Dimer CALIBRATOR (6)	18.4	µg /mL (approx.)	
AIA-PACK D-Dimer CONTROL SET			0025432
AIA-PACK D-Dimer CONTROL LEVEL 1	0.5	µg /mL (approx.)	
AIA-PACK D-Dimer CONTROL LEVEL 2	8.0	µg /mL (approx.)	
ST AIA-PACK D-Dimer SAMPLE DILUTING SOLUTION			0025532
AIA-PACK WASH CONCENTRATE			0020955
AIA-PACK DILUENT CONCENTRATE			0020956
SAMPLE CUPS			0018581
AIA-PACK DETECTOR STANDARDIZATION TEST CUP			0020970
AIA-PACK SAMPLE TREATMENT CUP			0020971
Additional Requirements for AIA-600 II, AIA-900, AIA-1800 and AIA-2000:			
PIPETTE TIPS			0019215
TIP RACK			0019216
PRELOADED PIPETTE TIPS			0022103

Only materials obtained from TOSOH should be used. Materials obtained elsewhere should not be substituted since assay performance is characterized based strictly on TOSOH materials.

### WARNINGS AND PRECAUTIONS

- The ST AIA-PACK D-Dimer is intended for in vitro diagnostic use only.
- TOSOH AIA Nex·IA / AIA-21 Immunoassay Analyzers can NOT be used to perform the ST AIA-PACK D-Dimer assay.
- Inspect the packaging and the exterior of the aluminium pouch for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
- Test cups from different lots or different assays shall not be mixed within a tray.
- The ST AIA-PACK D-Dimer contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- The material derived from human origin is not used in the preparation of this product; however, since human specimens will be used for samples and other quality control products in the lab may be derived from human origin, please use standard laboratory safety procedures in handling all specimens and controls.
- Do not use beyond the expiry date.
- The ST AIA-PACK D-Dimer has been designed so that the high dose "hook effect" is not a problem for the vast majority of samples. Sample with D-dimer concentration between 16.0 and 400 µg/mL will read > 16.0 µg/mL. The "hook effect" phenomenon may occur only at D-dimer concentrations > 400 µg/mL.
- For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

10. After opening, the vial of ST AIA-PACK D-Dimer SAMPLE DILUTING SOLUTION should be kept tightly sealed with a clean rubber cap. Sealing with dirty material may cause deterioration of the reagent.
11. The remaining sample diluting solution after use should not be mixed with another vial but be discarded to avoid contamination.
12. Serum, dust, metal, or microorganism contamination may cause degradation of reconstituted substrate solution. Store in a clean environment, away from direct sunlight and ultraviolet light.
13. TOSOH recommends that a new pouch of the test cups should be used for calibration.

## STORAGE AND STABILITY

All unopened materials will remain stable until the expiry date on the label when stored at the specified temperature.

Materials	Cat. No.
2°-8°C:	
ST AIA-PACK D-Dimer	0025232
ST AIA-PACK D-Dimer CALIBRATOR SET	0025332
ST AIA-PACK D-Dimer SAMPLE DILUTING SOLUTION	0025532
AIA-PACK D-Dimer CONTROL SET	0025432
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
1°-30°C:	
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970
AIA-PACK SAMPLE TREATMENT CUP	0020971

- After opening the aluminium pouch, ST AIA-PACK D-Dimer test cups can be left on-board of the TOSOH AIA System Analyzers (18°-25°C) for a maximum of 4 days (4 x 24 hours). When stored over night at 2°-8°C, the test cups can be used for up to 12 days (12 cycles of 8 hours on board and 16 hours in the refrigerator). Once the aluminium pouch is opened, even if the test cups are stored in the refrigerator, they must be used within 30 days.
- ST AIA-PACK D-Dimer CALIBRATOR SET must be kept tightly sealed and refrigerated at 2°-8°C. After opening or reconstituting, the calibrators should be used within 1 day.
- After opening, ST AIA-PACK D-Dimer SAMPLE DILUTING SOLUTION can be left on-board of the TOSOH AIA System Analyzers (18°-25°C) for a maximum of 3 days (3 x 24 hours). When stored over night at 2°-8°C, the sample diluting solution can be used for up to 9 days (9 cycles of 8 hours on board and 16 hours in the refrigerator). The sample diluting solution should not be used beyond 90 days after opening, even if it is sealed and stored in the refrigerator.
- AIA-PACK D-Dimer CONTROL SET must be kept tightly sealed and refrigerated at 2°-8°C. After opening or reconstituting, the controls should be used within 7 days.
- Reconstituted substrate solution will remain stable for 3 days at 18°-25°C or 30 days at 2°-8°C.
- Working diluent and wash solutions will remain stable for 30 days at 18°-25°C.
- Reagents should not be used if they appear cloudy or discoloured.

## SPECIMEN COLLECTION AND HANDLING

1. Citrated plasma is required for the assay. Serum, heparinised plasma or EDTA plasma SHOULD NOT BE USED.
2. When using citrated plasma, a venous blood sample is collected aseptically with the designated additive. Centrifuge and separate plasma from the packed cells as soon as possible.

3. Inadequate centrifugation or the presence of fibrin or particulate matter in the sample may cause an erroneous result.
4. Samples containing inhibitors of alkaline phosphatase may cause erroneous results.
5. Inspect all samples for air bubbles and foaming. Remove any air bubbles before assay.
6. Samples may be stored at 2°-8°C for up to 24 hours before analysis. If the analysis cannot be done within 24 hours, the sample should be stored frozen at -20°C or below for up to 60 days.
7. Repeated freeze-thaw cycles should be avoided. Turbid citrated plasma samples or samples containing particulate matter should be centrifuged before testing. Before the assay, bring frozen samples to 18°-25°C slowly and mix gently.
8. The sample required for analysis is 10 µL.

## PROCEDURE

For the AIA-600 II, AIA-900, AIA-1800, AIA-2000, and AIA-360, please refer to their Operators' manual for detailed instructions.

### I. Reagent Preparation

#### A. Substrate Solution

Bring all reagents to 18°-25°C before preparing the working reagent. Add the entire contents of the AIA-PACK SUBSTRATE RECONSTITUENT II (100 mL) to the lyophilised AIA-PACK SUBSTRATE REAGENT II and mix thoroughly to dissolve the solid material.

#### B. Wash Solution

Add the entire contents of the AIA-PACK WASH CONCENTRATE (100 mL) to approximately 2.0 L of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline, mix well, and adjust the final volume to 2.5 L.

#### C. Diluent

Add the entire contents of the AIA-PACK DILUENT CONCENTRATE (100 mL) to approximately 4.0 L of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline, mix well, and adjust the final volume to 5.0 L.

### II. Calibration Procedure


#### A. Calibration Curve

The calibrators for use with the ST AIA-PACK D-Dimer are prepared gravimetrically and are compared to internal reference standards.

The calibration curve for the ST AIA-PACK D-Dimer will remain stable for up to 90 days. Calibration stability is monitored by quality control performance and is dependent on proper reagent handling and TOSOH AIA System maintenance according to the manufacturer's instructions.

Recalibration may be necessary more frequently if controls are out of the established range for this assay or when certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of the wash probe, or detector lamp adjustment or change). For further information regarding instrument operation, consult the TOSOH AIA System Operators' manual.

A sample calibration curve from the AIA-2000 follows and shows the algorithm used for calculating results.

### B. Calibration Procedure

1. Refer to the appropriate TOSOH AIA System Operators' manual for the procedural instructions.
2. Verify that both the calibrator lot and concentration numbers have been correctly entered into the software.
3. The ST AIA-PACK D-Dimer CALIBRATOR (1) is provided ready for use.
4. The ST AIA-PACK D-Dimer CALIBRATOR (2)-(6) are lyophilised. All levels should be reconstituted with 1.0 mL of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline.
5. TOSOH recommends that all calibrators be run in triplicate.

### C. Calibration Acceptability Criteria

1. The mean rate for the ST AIA-PACK D-Dimer CALIBRATOR (1) should be < 0.2 nmol/(L·s).
2. Since there is a direct relationship between concentration and rate, the rate should increase as the concentration increases.
3. The replicate values should be within a 10% range.

### D. Calibration Review and Acceptance

1. Review the calibration curve carefully, using the criteria listed above.
2. Edit the calibration if necessary, then accept the calibration.

For further information regarding calibration, consult the TOSOH AIA System Operators' manual.

## III. Quality Control Procedure

### A. Commercially Available Controls

Commercially available controls shall be run at least once per day. It is recommended that at least two levels of controls, normal and abnormal, be used. Laboratory policy for this particular assay designates the following:

Control Material: \_\_\_\_\_

Frequency: \_\_\_\_\_

Lot number of control material, acceptable limits, and corrective action to be taken if controls do not meet laboratory criteria will be found in a separate quality control document maintained by the laboratory.

### B. Quality Control Procedure

1. Assay quality control specimens as instructed in the specific Operators' manual for your analyzer. In addition, refer to the TOSOH AIA System Operators' manual for detailed instructions on defining and editing the files.

2. Quality control material to be run with this assay is defined by individual laboratory policy.

## IV. Specimen Processing

### A. Preparation

Following the specific instructions in the Operators' manual for the analyzer, place samples on the instrument appropriately. Barcoded primary tubes as well as sample cups can be run on the AIA-600 II, AIA-900, AIA-1800, AIA-2000 and AIA-360.

### B. Assay Procedure

1. Ensure a sufficient quantity of ST AIA-PACK D-Dimer test cups for the number of samples to be run.
2. Load patient samples as instructed in the Operators' manual and proceed with analysis.

Note: The AIA-600 II, AIA-900, AIA-1800 and AIA-2000 will require AIA-PACK SAMPLE TREATMENT CUPS if on-board dilutions are utilized.

## PROCEDURAL NOTES

1. Lyophilised substrate must be completely dissolved.
2. Ligand assays performed by the TOSOH AIA System Analyzers require that the laboratory use water designated by the CAP as Class I or by CLSI as the clinical laboratory reagent water. Water should be tested at least once per month and should be free of particulate matter including bacteria. For further information, consult the CLSI document GP40-A4-AMD, Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline-Fourth Edition.
3. If a specimen D-dimer concentration is found to be greater than the upper limit of the assay range, 16.0 µg/mL, the specimen should be diluted with the ST AIA-PACK D-Dimer SAMPLE DILUTING SOLUTION and reassayed according to the Assay Procedure. The recommended dilution for specimens containing greater than 16.0 µg/mL is 5-fold dilution. It is desirable to dilute the specimen so that the diluted specimen reads between 0.02 and 16.0 µg/mL. The dilution factor should be entered into the software. For further information on the dilution of specimens, refer to the TOSOH AIA System Operators' manual.
4. The TOSOH AIA System Analyzers can store two different calibration curves for each analyte at one time. Therefore, up to two different lots of ST AIA-PACK D-Dimer test cups can be used during the same run.
5. If the assay specifications for this test are not ready in the system software, the specifications must be entered under test code 104.

## CALCULATION OF RESULTS

The TOSOH AIA System Analyzers perform all sample and reagent handling operations automatically. The TOSOH AIA System Analyzers read the rate of fluorescence produced by the reaction and automatically convert the rate to D-dimer concentration in µg/mL.

For samples requiring dilution, the AIA-600 II, AIA-900, AIA-1800 and AIA-2000 will automatically perform dilutions and calculate results if the dilution factors are entered into the software. For detailed information regarding programming dilutions, consult the appropriate TOSOH AIA System Operators' manual.

## EVALUATION OF RESULTS

### Quality Control

In order to monitor and evaluate the precision of the analytical performance, it is recommended that commercially available control samples shall be assayed according to the local regulations.

The minimum recommendations for the frequency of running internal control material are:

After calibration, two levels of the internal control are run in order to accept the calibration curve. The two levels of controls are repeated when certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of the wash probe, or detector lamp adjustment or change).

After daily maintenance, two levels of the control shall be run in order to verify the overall performance of the TOSOH AIA System Analyzers.

If one or more control value(s) is out of the acceptable range, it is necessary to investigate the validity of the calibration curve before reporting patient results. Standard laboratory procedures should be followed in accordance with the strict regulatory agency under which the laboratory operates.

#### LIMITATIONS OF THE PROCEDURE

- For diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g. symptoms, results of other tests, clinical impressions, therapy, etc.).
- Using ST AIA-PACK D-Dimer, the highest measurable concentration of D-dimer in specimens without dilution is 16.0 µg/mL, and the lowest measurable concentration in specimens is 0.02 µg/mL (assay sensitivity).
- Although the approximate value of the highest calibrator is 18.4 µg/mL, the exact concentration may be slightly different. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 16.0 µg/mL.
- Although haemolysis has an insignificant effect on the assay, haemolysed samples may indicate mistreatment of a specimen before assay and results should be interpreted with caution.
- Lipaemia has an insignificant effect on the assay except in the case of gross lipaemia where spatial interference may occur.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show falsely elevated or decreased D-dimer values.
- Samples from patients under Asfotase Alfa (Genetical Recombination) treatment may cause falsely elevated/decreased results.
- For a more complete understanding of the limitations of this procedure, please refer to the SPECIMEN COLLECTION AND HANDLING, WARNINGS AND PRECAUTIONS, STORAGE AND STABILITY, and PROCEDURAL NOTES sections in this insert sheet.

#### EXPECTED VALUES

Each laboratory should determine a reference interval corresponding to the characteristics of the population being tested. As with all diagnostic procedures, clinical results must be interpreted with regard to concomitant medications administered to the patient.

##### I. Reference Ranges

The interval given here was determined in citrated plasma samples from 131 apparently healthy Asian individuals. European studies showed that these reference values are also applicable for a European population.

Reference Interval < 0.5 µg/mL

##### II. Conversion Factors

D-dimer concentrations in this application are in units of µg/mL. The obtained values in µg/mL must be multiplied by 1,000 when converting to values in µg/L.

#### PERFORMANCE CHARACTERISTICS

##### ACCURACY

- Recovery: Three citrated plasma samples were spiked with three different concentration levels of D-dimer and assayed before and after spiking.

Sample	Initial Value (µg/mL)	D-dimer Added (µg/mL)	Expected Value (µg/mL)	Measured Value (µg/mL)	Percent Recovery (%)
Plasma A1	0.0901	0.334	0.424	0.395	93.1
	0.0901	0.667	0.757	0.713	94.1
	0.0901	1.33	1.42	1.36	95.2
Plasma B1	0.101	0.334	0.435	0.420	96.7
	0.101	0.667	0.768	0.741	96.5
	0.101	1.33	1.44	1.44	100.6
Plasma C1	0.0849	0.334	0.419	0.391	93.4
	0.0849	0.667	0.752	0.729	96.9
	0.0849	1.33	1.42	1.42	99.9

- Dilution: Three citrated plasma samples containing high concentrations of D-dimer were serially diluted with the ST AIA-PACK D-Dimer SAMPLE DILUTING SOLUTION and assayed.

Sample	Dilution Factor	Expected Value (µg/mL)	Measured Value (µg/mL)	Percent Recovery (%)
Plasma A2	none	12.7	12.7	100.0
	7.5/10	9.49	9.90	104.4
	5.0/10	6.33	6.84	108.1
	2.5/10	3.16	3.59	113.6
	1.0/10	1.27	1.36	107.6
Plasma B2	none	7.98	7.98	100.0
	7.5/10	5.99	6.31	105.5
	5.0/10	3.99	4.24	106.2
	2.5/10	2.00	2.11	105.8
	1.0/10	0.798	0.862	108.0
Plasma C2	none	2.28	2.28	100.0
	7.5/10	1.71	1.73	101.3
	5.0/10	1.14	1.13	98.7
	2.5/10	0.570	0.537	94.1
	1.0/10	0.228	0.244	106.8

- Linearity: The linearity for ST AIA-PACK D-Dimer was determined, based on guidance from CLSI Protocol EP6-A. The linearity was measured on the AIA-600 II, AIA-900, AIA-1800, AIA-2000 and AIA-360 instrument and has been demonstrated to be linear from 0.02 to 16.0 µg/mL.

##### PRECISION

- Within-run precision was determined using three controls in a total of 20 runs. Within each run, one set of duplicates per control was assayed. The mean of each duplicate was used to

obtain the pooled standard deviation (SD), which was then used to calculate the coefficient of variation (CV).

Sample	Mean	Pooled SD	CV
	(µg/mL)	(µg/mL)	(%)
Plasma A3	0.274	0.0126	4.6
Plasma B3	3.73	0.161	4.3
Plasma C3	8.27	0.398	4.8

- b. Total precision was determined by the duplicate assay of three controls in 20 separate runs. The means of each run were used to calculate the pooled standard deviation (SD) and coefficient of variation (CV).

Sample	Mean	Pooled SD	CV
	(µg/mL)	(µg/mL)	(%)
Plasma A3	0.274	0.0131	4.8
Plasma B3	3.73	0.180	4.8
Plasma C3	8.27	0.428	5.2

#### SPECIFICITY

The following substance was tested for cross-reactivity. The cross-reactivity (%) is the percentage of the compound which will be identified as D-dimer. If this compound is present in the specimen at the same concentration as D-dimer, the final result will be increased by this percentage.

Compound	Concentration added	Cross-reactivity
	(µg/mL)	(%)
Fibrinogen	4,000	0.0047

#### SENSITIVITY

Limit of Detection: The limit of detection for ST AIA-PACK D-Dimer was determined, according to CLSI guideline EP17-A. Blank sample was measured in 60 replicates. Six low level samples were measured in 10 replicates each. As a result, the limit of detection for ST AIA-PACK D-Dimer was estimated to be 0.0132 µg/mL.

#### INTERFERENCE

Interference is defined, for the purposes of this study, with recovery outside of 10% of the known concentration of the specimen after the following substances are added to human specimens.

1. Haemoglobin (up to 440 mg/dL), free bilirubin (up to 17 mg/dL), and conjugated bilirubin (up to 18 mg/dL) do not interfere with the assay.
2. Lipaemia, as indicated by triglyceride concentration (up to 1,600 mg/dL), does not interfere with the assay.
3. Ascorbic acid (up to 20 mg/dL) does not interfere with the assay.
4. Protein, as indicated by human albumin concentration (up to 5 g/dL added to samples from apparently healthy subjects), does not interfere with the assay.
5. Citric acid (up to 20 mg/mL) does not interfere with the assay.

#### REFERENCES

1. Carey MJ and Rodgers GM: Disseminated intravascular coagulation: Clinical and laboratory aspects., Am J Hemat, 59, 65-73 (1998).

2. Mammen EF: Disseminated intravascular coagulation (DIC), Clin Lab Sci, 13, 239-245 (2000).
3. Yu M, Nardella A and Pechet L: Screening tests of disseminated intravascular coagulation: guidelines for rapid and specific laboratory diagnosis., Crit Care Med, 28, 1777-1780 (2000).
4. Anderson DR and Wells PS: D-dimer for the diagnosis of venous thromboembolism., Curr Opin Hematol, 7, 296-301 (2000).
5. Perrier A, Desmarais S, Goehring C, de Moerloose P, Morabia A, Unger PF, et al.: D-dimer testing for suspected pulmonary embolism in outpatients., Am J Respir Crit Care Med, 156, 492-496 (1997).

#### TOSOH CORPORATION

3-8-2, Shiba-Koen, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

#### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO, BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

#### TOSOH BIOSCIENCE, INC.

6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080, USA  
Phone: +1 650 615 49 70

Net Volume (after reconstitution for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## AIA-PACK SUBSTRATE SET II

### INTENDED USE

The AIA-PACK SUBSTRATE SET II is intended for in vitro diagnostic use only for the AIA-PACK or the ST AIA-PACK Assay.

### MATERIAL SUPPLIED

Cat. No.	0020968
2 bottles	AIA-PACK SUBSTRATE REAGENT II 4-methylumbelliferyl phosphate, stabilisers and sodium azide as a preservative (lyophilised).
2 x 100 ml	AIA-PACK SUBSTRATE RECONSTITUENT II Buffer containing sodium azide as a preservative (Liquid).

### WARNINGS AND PRECAUTIONS

1. The AIA-PACK SUBSTRATE SET II is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the bottles for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. These materials contain sodium azide which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always rinse with large volumes of water to prevent any azide accumulation.
4. Do not use beyond the expiry date.
5. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

1. Bring all reagents to 18°-25°C before preparing the working substrate solution.
2. Add one bottle of AIA-PACK SUBSTRATE RECONSTITUENT II to the lyophilised AIA-PACK SUBSTRATE REAGENT II, mix thoroughly to dissolve all solid materials.
3. Each bottle of AIA-PACK SUBSTRATE REAGENT II, when reconstituted, provides sufficient reagent for approximately 400 tests on the TOSOH AIA System Analysers.

### STORAGE AND STABILITY

- Always store the AIA-PACK SUBSTRATE SET II in an upright position at 2°-8°C when not in use.
- Protect it from light.
- When stored unopened and refrigerated at 2°-8°C, the AIA-PACK SUBSTRATE SET II will remain stable until the expiry date on the label.
- The reconstituted substrate solution will remain stable for 3 days at 18°-25°C or 30 days at 2°-8°C.
- Serum, dust, metal, or microorganism contamination may cause degradation of reconstituted substrate solution. Store in a clean environment, away from direct sunlight and ultraviolet light.

### PROCEDURE

Please refer to the insert sheet of the AIA-PACK or the ST AIA-PACK and the TOSOH AIA System Operators' Manual of the Tosoh AIA System Analysers for detailed instructions.

### TOSOH CORPORATION


3-8-2, Shiba, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

### TOSOH BIOSCIENCE, INC.

6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080 (USA)  
Phone: +1 650 615-4970

Net Volume (after reconstitution for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## ST AIA-PACK TSH

### INTENDED USE

ST AIA-PACK TSH is designed for in vitro diagnostic use only for the quantitative measurement of Thyroid Stimulating Hormone (TSH or thyrotropin) in **human serum or heparinised plasma** on Tosoh AIA System Analysers.

### 2. Mëginiq tipas

### SUMMARY AND EXPLANATION

Thyroid Stimulating Hormone is a glycoprotein hormone secreted by the anterior pituitary gland. When feedback suppression of the pituitary is reduced by a reduced production of thyroid hormones (T<sub>4</sub> and T<sub>3</sub>), TSH rises in an attempt to increase thyroid hormone production. This rise occurs while the patient is still asymptomatic and is thus an early and very sensitive indication of hypothyroidism (1-5). TSH is also controlled by the hypothalamic peptide, thyrotropin releasing hormone (TRH).

Accurate determination of serum TSH is the most useful and sensitive test for primary hypothyroidism, where serum thyroid hormone concentrations are depressed and serum TSH concentrations are significantly elevated. Serum TSH determinations may also be used to differentiate between pituitary (secondary) and hypothalamic (tertiary) hypothyroidisms (6-9). Through the use of monoclonal antibody technology which provides the necessary specificity and sensitivity, the usefulness of TSH determination in the diagnosis of hyperthyroidism distinguished from euthyroidism has been well established (10,11).

### PRINCIPLE OF THE ASSAY

The ST AIA-PACK TSH is a two-site immunoenzymometric assay performed entirely in the ST AIA-PACK TSH test cups. TSH present in the test sample is bound with monoclonal antibody immobilised on magnetic beads and monoclonal antibody conjugated with bovine alkaline phosphatase in the test cups. The magnetic beads are washed to remove any non-bound enzyme-labelled monoclonal antibodies and then incubated with a fluorogenic substrate, 4-methylumbelliferyl phosphate (4MUP). The amount of enzyme conjugated with monoclonal antibody that binds to the beads is directly proportional to the TSH concentration in the test sample. A standard curve is constructed, and unknown sample concentrations are calculated using this curve.

### MATERIAL SUPPLIED

Cat. No: 0025294  
ST AIA-PACK TSH 5 trays x 20 test cups

Plastic test cups containing twelve lyophilised magnetic beads coated with anti-TSH mouse monoclonal antibody and 50 µl of anti-TSH mouse monoclonal antibody conjugated to bovine alkaline phosphatase with sodium azide as a preservative.

### MATERIALS REQUIRED BUT NOT SUPPLIED

The following materials are required to perform TSH analysis using the ST AIA-PACK TSH (Cat. No. 0025294) on Tosoh AIA System Analysers. The following are available separately from Tosoh.

Materials	Cat. No.
AIA Nex•IA or AIA-21	0018539

AIA Nex•IA or AIA-21 LA	0018540
AIA-1800 ST	0019836
AIA-1800 LA	0019837
AIA-2000 ST	0022100
AIA-2000 LA	0022101
AIA-600 II	0019014
AIA-600 II BCR	0019328
AIA-900	0022930
AIA-360	0019945
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK SUBSTRATE REAGENT II/ AIA-PACK SUBSTRATE RECONSTITUENT II	
AIA-PACK TSH 3rd-Gen CALIBRATOR SET	0020394
AIA-PACK TSH 3rd-Gen CALIBRATOR (1)	0 µIU/ml
AIA-PACK TSH 3rd-Gen CALIBRATOR (2)	0.2 µIU/ml (approx.)
AIA-PACK TSH 3rd-Gen CALIBRATOR (3)	5.0 µIU/ml (approx.)
AIA-PACK TSH 3rd-Gen CALIBRATOR (4)	25 µIU/ml (approx.)
AIA-PACK TSH 3rd-Gen CALIBRATOR (5)	50 µIU/ml (approx.)
AIA-PACK TSH 3rd-Gen CALIBRATOR (6)	110 µIU/ml (approx.)
AIA-PACK TSH 3rd-Gen SAMPLE DILUTING SOLUTION	0020594
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
SAMPLE CUPS	0018581
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970
AIA-PACK SAMPLE TREATMENT CUP	0020971
Additional Requirements for AIA Nex•IA/AIA-21 only:	
PIPETTE TIPS	0018552
PRELOADED PIPETTE TIPS	0018583
Additional Requirements for AIA-600 II, AIA-900, AIA-1800, AIA-2000:	
PIPETTE TIPS	0019215
TIP RACK	0019216
PRELOADED PIPETTE TIPS	0022103

Only materials obtained from Tosoh should be used. Materials obtained elsewhere should not be substituted since assay performance is characterised strictly by Tosoh materials.

### WARNINGS AND PRECAUTIONS

- The ST AIA-PACK TSH is intended for in vitro diagnostic use only.
- Inspect the packaging and the exterior of the aluminium pouch for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
- Test cups from different lots or different assays shall not be mixed within a tray.
- The ST AIA-PACK TSH contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always rinse with large volumes of water to prevent any azide accumulation.
- Human serum is not used in the preparation of this product; however, since human specimens will be used for samples and other quality control products in the lab may be derived from human serum, please use standard laboratory safety procedures in handling all specimens and controls.
- Do not use beyond the expiry date.
- The ST AIA-PACK TSH has been designed so that the high dose hook effect is not a problem for the vast majority of samples. Samples with TSH concentrations between 100 and 5,000 µIU/ml will read >100 µIU/ml. The "hook effect" phenomenon may occur at TSH concentrations > 5,000 µIU/ml.

8. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.
9. After opening, the vial of the sample diluting solution should be kept tightly sealed with a clean rubber cap. Sealing with dirty material may cause deterioration of the reagent.
10. The remaining sample diluting solution after use should not be mixed with another vial but be discarded to avoid contamination.
11. Serum, dust, metal, or microorganism contamination may cause degradation of reconstituted substrate solution. Store in a clean environment, away from direct sunlight and ultraviolet light.
12. TOSOH recommends that a new pouch of the test cups should be used for calibration.

## STORAGE AND STABILITY

All unopened materials are stable until the expiry date on the label when stored at the specified temperature.

Materials	Cat. No.
2°-8°C:	
ST AIA-PACK TSH	0025294
AIA-PACK TSH 3rd-Gen CALIBRATOR SET	0020394
AIA-PACK TSH 3rd-Gen SAMPLE DILUTING SOLUTION	0020594
AIA-PACK SUBSTRATE SET II	0020968
AIA-PACK WASH CONCENTRATE	0020955
AIA-PACK DILUENT CONCENTRATE	0020956
1°-30°C:	
AIA-PACK DETECTOR STANDARDIZATION TEST CUP	0020970
AIA-PACK SAMPLE TREATMENT CUP	0020971

- After opening the aluminium pouch, ST AIA-PACK TSH test cups can be left on-board of the TOSOH AIA System Analyzers (18-25°C) for a maximum of 10 days (10 x 24 hours). When stored overnight at 2-8°C, the test cups can be used for up to 30 days (30 cycles of 8 hours on board and 16 hours in the refrigerator). Once the aluminium pouch is opened, the test cups must be used within 30 days.
- AIA-PACK TSH 3rd-GEN CALIBRATOR SET must be kept tightly sealed and refrigerated at 2-8°C. After opening, the calibrators should be used within 1 day.
- After opening, AIA-PACK TSH 3rd-GEN SAMPLE DILUTING SOLUTION can be left on-board of the TOSOH AIA System Analyzers (18-25 C) for a maximum of 3 days (3 x 24 hours). When stored overnight at 2-8°C, the sample diluting solution can be used for up to 9 days (9 cycles of 8 hours on board and 16 hours in the refrigerator). The sample diluting solution should not be used beyond 90 days after opening, even if it is sealed and stored in the refrigerator.
- Reconstituted substrate solution will remain stable for 3 days at 18°-25°C or 30 days at 2°-8°C.
- Working diluent and wash solutions will remain stable for 30 days at 18°-25°C.
- Reagents should not be used if they appear cloudy or discoloured.

## SPECIMEN COLLECTION AND HANDLING

- Serum or heparinised plasma is required for the assay. EDTA and citrated plasma SHOULD NOT BE USED.
- When using serum, a venous blood sample is collected aseptically without additives. Store at 18°-25°C until a clot has formed (usually 15-45 minutes), then centrifuge to obtain the serum specimen for assay.
- When using heparinised plasma, a venous blood sample is collected aseptically with the specified additive. Centrifuge and separate plasma from the packed cells as quickly as possible.
- Inadequate centrifugation or the presence of fibrin or particulate matter in the sample may cause an erroneous result.

- Samples containing inhibitors of alkaline phosphatase may cause erroneous results.
- Inspect all samples for air bubbles and foaming. Remove any air bubbles before assay.
- Specimen types should not be used interchangeably during serial monitoring of an individual patient. Measured concentrations may vary slightly between sample types in certain patients.
- Samples may be stored at 2°-8°C for up to 7 days before analysis. If the analysis cannot be performed within 7 days, the sample should be stored frozen at -20°C or below for up to 60 days.
- Repeated freeze-thaw cycles should be avoided. Turbid serum samples or samples containing particulate matter should be centrifuged before testing. Before assaying, slowly bring frozen samples to 18°-25°C and mix gently.
- The sample required for analysis is 100 µl.

## PROCEDURE

For the AIA Nex•IA / AIA-21, AIA-600 II, AIA-900, AIA-1800, AIA-2000 and AIA-360, please refer to their Operators' Manual for detailed instructions.

### I. Reagent Preparation

#### A. Substrate Solution

Bring all reagents to 18°-25°C before preparing the working substrate solution. Add one bottle of AIA-PACK SUBSTRATE RECONSTITUENT II to the lyophilised AIA-PACK SUBSTRATE REAGENT II, mix thoroughly to dissolve all solid materials.

#### B. Wash Solution

Add the entire contents of the AIA-PACK WASH CONCENTRATE (100 ml) to approximately 2.0 L of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline, mix well, and adjust the final volume to 2.5 L.

#### C. Diluent

Add the entire contents of the AIA-PACK DILUENT CONCENTRATE (100 ml) to approximately 4.0 L of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI GP40-A4-AMD guideline, mix well, and adjust the final volume to 5.0 L.

### II. Calibration Procedure


#### A. Calibration Curve

The calibrators for use with the ST AIA-PACK TSH have been standardised on WHO 2nd IRP 80/558 (1983).

The calibration curve for the ST AIA-PACK TSH will remain stable for up to 90 days. Calibration stability is monitored by quality control performance and is dependent on proper reagent handling and Tosoh AIA System maintenance in accordance with the manufacturer's instructions.

Recalibration may be necessary more frequently, if controls are out of the established range for this assay or when certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of wash probe, or detector lamp adjustment or change. For further information regarding instrument operation, consult the TOSOH AIA System Operator's Manual.

A sample calibration curve from the AIA-1800 follows and shows the algorithm used to calculate results.

### B. Calibration Procedure

1. Refer to the appropriate Tosoh AIA System Operators' Manual for procedural instructions.
2. Verify that both the calibrator lot and concentration numbers have been correctly entered on the software.
3. The AIA-PACK TSH 3rd-Gen CALIBRATOR SET is provided ready for use.
4. Tosoh recommends that all calibrators should be run in triplicate.

### C. Calibration Acceptability Criteria

1. The mean rate for the AIA-PACK TSH 3rd-Gen CALIBRATOR (1) should be < 3.0 nmol/(L·s).
2. Since there is a direct relationship between concentration and rate, the rate should rise as concentrations increase.
3. The replicate values should be within a 10% range.

### D. Calibration Review and Acceptance

1. Review the calibration curve carefully, applying the criteria listed above.
2. Edit the calibration if necessary, then accept the calibration.

For further information regarding calibration, consult the Tosoh AIA System Operators' Manual.

## III. Quality Control Procedure

### A. Commercially Available Controls

Commercially available controls shall be run at least once per day. It is recommended that at least two (2) levels of controls, normal and abnormal, be used. Laboratory policy for this particular assay specifies the following:

Control Material: \_\_\_\_\_

Frequency: \_\_\_\_\_

Lot number of control material, acceptable limits, and corrective action to be taken if controls fail to meet laboratory criteria will be found in a separate quality control document kept by the laboratory.

### B. Quality Control Procedure

1. Assay quality control specimens as instructed in the specific Operators' Manual for your analyser. In addition, refer to the Tosoh AIA System Operators' Manual for detailed instructions on defining and editing files.
2. Quality control material to be run with this assay is specified by individual laboratory policy.

## IV. Specimen Processing

### A. Preparation

Following the specific instructions in the Operators' Manual for the analyser, place samples on the instrument appropriately. Barcoded primary tubes as well as sample cups can be run on the AIA Nex•IA / AIA-21, AIA-600 II, AIA-900, AIA-1800, AIA-2000 and AIA-360.

### B. Assay Procedure

1. Provide a sufficient quantity of ST AIA-PACK TSH test cups for the number of samples to be run.
2. Load patient samples as instructed in the Operators' Manual and proceed with analysis.

Note: The AIA Nex•IA / AIA-21, AIA-600 II, AIA-900, AIA-2000 and AIA-1800 will require AIA-PACK SAMPLE TREATMENT CUP if onboard dilutions are used.

## PROCEDURAL NOTES

1. Lyophilised substrate must be completely dissolved.
2. Ligand assays performed by the TOSOH AIA System Analyzers require that the laboratory use water designated by the CAP as Class I or by the CLSI as the clinical laboratory reagent water. Water should be tested at least once per month and should be free of particulate matter including bacteria. For further information, consult the CLSI document GP40-A4-AMD, Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline-Fourth Edition.
3. If a specimen TSH concentration is found to be greater than the upper limit of the assay range, 100 µIU/mL, the specimen should be diluted with the AIA-PACK TSH 3rd-Gen SAMPLE DILUTING SOLUTION and reassayed in accordance with the Assay Procedure. The recommended dilution for specimens containing more than 100 µIU/ml is 10-fold dilution. It is desirable to dilute the specimen so that the diluted specimen reads between 10 and 100 µIU/ml. The dilution factor should be entered on the software. For further information on the dilution of specimens, refer to the Tosoh AIA System Operators' Manual.
4. The Tosoh AIA System Analyzers can store two different calibration curves for each analyte at one time. Thus, up to two different lots of ST AIA-PACK TSH test cups can be used in the same run.
5. If the assay specifications for this test are not ready on the system software, the specifications must be entered under test code 072.

## CALCULATION OF RESULTS

The Tosoh AIA System Analyzers perform all sample and reagent handling operations automatically. The Tosoh AIA System Analyzers read the rate of fluorescence produced by the reaction and automatically convert the rate to TSH concentration in µIU/ml.

For samples requiring dilution, the AIA Nex•IA / AIA-21, AIA-600 II, AIA-900, AIA-2000 and AIA-1800 will automatically perform dilutions and calculate results if the dilution factors are entered on the software. For detailed information regarding programming dilutions, consult the appropriate TOSOH AIA System Operators' Manual.

## EVALUATION OF RESULTS

### Quality Control

In order to monitor and evaluate the precision of the analytical performance, it is recommended that commercially-available control samples shall be assayed in accordance with local regulations.

Minimum recommendations for the frequency of running internal control material are:  
After calibration, three levels of the internal control are run in order to accept the calibration curve.

The three levels of controls are repeated when certain service procedures are performed (e.g. temperature adjustments, sampling mechanism changes, maintenance of the wash probe, or detector lamp adjustments or changes).

After daily maintenance, at least two levels of the control shall be run in order to verify overall performance of the Tosoh AIA System Analysers.

If one or more control sample value(s) is out of the acceptable range, it is necessary to investigate the validity of the calibration curve before reporting patient results.

Standard laboratory procedures should be followed in accordance with the strict regulatory agency under which the laboratory operates.

## LIMITATIONS

- For diagnostic purposes, results obtained from this assay should be used in conjunction with other data (e.g. symptoms, results of other tests, clinical impressions, therapy, etc.).
- Using ST AIA-PACK TSH, the highest measurable concentration of TSH in specimens without dilution is 100 µIU/ml, and the lowest measurable concentration in specimens is 0.01 µIU/ml (assay sensitivity).
- Although the approximate value of the highest calibrator is 110 µIU/ml, the exact concentration may be slightly different. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 100 µIU/ml.
- Although haemolysis has an insignificant effect on the assay, haemolysed samples may indicate mistreatment of a specimen before assaying and results should be interpreted with caution.
- Lipaemia has an insignificant effect on the assay except in the case of gross lipaemia where spatial interference may occur.
- Specimens from patients taking medicines and/or medical treatment may show erroneous results.
- Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show falsely elevated values when tested for TSH.
- Samples from patients under Asfotase Alfa (Genetical Recombination) treatment may cause falsely elevated/decreased results.
- For a fuller understanding of the limitations of this procedure, please refer to the SPECIMEN COLLECTION AND HANDLING, WARNINGS AND PRECAUTIONS, STORAGE AND STABILITY, and PROCEDURAL NOTES sections in this insert sheet.

## EXPECTED VALUES

Each laboratory should determine a reference interval corresponding to the characteristics of the population being tested. As with all diagnostic procedures, clinical results must be interpreted having due regard to concomitant medications administered to the patient (12).

### I. Reference Ranges

The intervals given here were determined in serum samples from 497 apparently healthy Asian individuals and 166 apparently healthy Caucasian individuals, respectively.

- Asian**
  - Number of Samples (n) 497
  - Reference Interval 0.38 - 4.31 µIU/ml (0.38 - 4.31 mIU/l)
- Caucasian**
  - Number of Samples (n) 166
  - Reference Interval 0.40 - 3.60 µIU/mL (0.40 - 3.60 mIU/L)

### II. Conversion Factors

TSH concentrations in this application are in units of µIU/ml. The obtained values in µIU/mL must be multiplied by 1.0 when converting to values in mIU/L.

## PERFORMANCE CHARACTERISTICS

### ACCURACY

- a. Recovery: Three serum samples were spiked with three different levels of TSH and assayed before and after spiking.

Sample	Initial	TSH	Expected	Measured	Percent
	Value (µIU/ml)	Added (µIU/ml)	Value (µIU/ml)	Value (µIU/ml)	Recovery (%)
Serum A1	1.71	20.7	22.4	23.0	102.7
	1.71	41.4	43.1	45.3	105.1
	1.71	82.7	84.4	91.0	107.9
Serum B1	1.37	20.7	22.0	21.3	96.4
	1.37	41.4	42.7	44.8	104.9
	1.37	82.7	84.1	90.7	107.9
Serum C1	1.42	20.7	22.1	24.2	109.4
	1.42	41.4	42.8	48.8	114.1
	1.42	82.7	84.1	89.6	106.5

- b. Dilution: Three serum samples containing high concentrations of TSH were serially diluted with AIA-PACK TSH 3rd-Gen SAMPLE DILUTING SOLUTION and assayed.

Sample	Dilution	Expected	Measured	Percent
	Factor	Value (µIU/ml)	Value (µIU/ml)	Recovery (%)
Serum A2	none		98.6	
	7.5/10	73.9	72.8	98.5
	5.0/10	49.3	49.4	100.2
Serum B2	2.5/10	24.6	24.4	99.2
	none		104	
	7.5/10	77.7	75.3	96.9
Serum C2	5.0/10	51.8	51.6	99.7
	2.5/10	25.9	25.9	100.1
	none		104	
Serum C2	7.5/10	78.0	76.4	97.9
	5.0/10	52.0	51.2	98.5
	2.5/10	26.0	26.3	101.3

### PRECISION

- a. Within run precision was determined using three controls in a total of 20 runs. Within each run, one set of duplicates per control was assayed. The mean of each duplicate was used to obtain the pooled standard deviation (SD), which was then used to calculate the coefficient of variation (CV).

Sample	Mean (µIU/ml)	Pooled SD (µIU/ml)	CV (%)
Serum A3	2.21	0.0555	2.5
Serum B3	9.33	0.283	3.0
Serum C3	80.5	2.13	2.6

- b. Total precision was determined by the duplicate assay of three controls in 20 separate runs. The means of each run were used to calculate the pooled standard deviation (SD) and coefficient of variation (CV).

Sample	Mean ( $\mu\text{IU/ml}$ )	Pooled SD ( $\mu\text{IU/ml}$ )	CV (%)
Serum A3	2.21	0.0699	3.2
Serum B3	9.33	0.307	3.3
Serum C3	80.5	2.39	3.0

#### CORRELATION

The correlation between serum (x) and heparinised plasma (y) on ST AIA-PACK TSH was carried out using 109 patients' specimens.

Slope	0.955
y-Intercept	0.172
Correlation Coefficient	0.999
Range of Samples	0.01 - 94.50
Number of Samples	109

#### SPECIFICITY

The following substances were tested for cross-reactivity. Cross-reactivity is expressed in terms of the percentage of concentration of each substance which will be measured as TSH.

Compound	Cross-reactivity (%)
TSH	100
HCG	< 0.01
FSH	< 0.01
LH	0.17
HGH	< 0.01

#### SENSITIVITY

The minimum detectable concentration (MDC) of TSH is estimated to be 0.01  $\mu\text{IU/ml}$ . The MDC is defined as the concentration of TSH that corresponds to the rate of fluorescence which is two standard deviations from the mean rate of fluorescence of 5 replicate determinations by AIA-PACK TSH 3rd-Gen CALIBRATOR (1).

#### INTERFERENCE

Interference is defined, for the purposes of this study, with recovery outside of 10% of the known concentration of the specimen after the following substances are added to human specimens.

- Haemoglobin (up to 390 mg/dl), free bilirubin (up to 17 mg/dl) and conjugated bilirubin (up to 18 mg/dl) do not interfere with the assay.
- Lipaemia, as indicated by triglyceride concentration (up to 1,600 mg/dl), does not interfere with the assay.
- Ascorbic acid (up to 20 mg/dl) does not interfere with the assay.
- Protein, as indicated by human albumin concentration (up to 5.0 g/dl), does not interfere with the assay.
- Heparin (up to 100 U/ml) does not interfere with the assay.

#### REFERENCES

1. Berger, S. and Quim, J. L., Thyroid Function. In: Fundamentals of Clinical Chemistry. ed. N.W. Tietz. W.B. Saunders Co., Philadelphia, PA. 1976, pp. 824-848.

2. Hoskins, R. G., The Thyroid-Pituitary Apparatus as a Servo (Feed-Back) Mechanism. J.Clin. Endocrinol. Metab. 9, 1429-1431, 1949.
3. Reichlin, S. and Utiger, R. D., Regulation of the Pituitary-Thyroid Axis in Man: Relationship in TSH Concentration to Concentration of Free and Total Thyroxine in Plasma. J. Clin. Endocr. 27, 251-255, 1967.
4. Burger, H. G., and Patel, Y. C., The Value of Serum Thyrotropin Measurement in the Diagnosis and Management of Hypothyroidism. Med. J. Aust. 2,293-297, 1972.
5. Hershman, J. M., and Pittman, J. A., Utility of the Radioimmunoassay of Serum Thyrotropin in Man. Ann. Int. Med. 74, 481-490, 1971.
6. Karlberg, B. E., and Almquist, S., Clinical Experience with Thyrotropin Releasing Hormone (TRH) Stimulation Test in Patients with Thyroid, Pituitary and Hypothalamic Disorders. Acta Endocr. 72, 697-713, 1973.
7. Smith, S. R., Burton, D. W., and Izenstark, J. L., Levels of Serum Thyroid-Stimulating Hormone in Hyperthyroidism. Arch. Intern. Med. 139, 757-760, 1979.
8. Hall, R., Ormston, B. J., Besser, G. M., Cryer, R. J., and McKendrick, M., The Thyrotropin-Releasing Hormone Test in Diseases of the Pituitary and Hypothalamus. Lancet 1 (7754), 759-762, 1972.
9. Lauridsen, U. B., Deckert, T., Friis, T., Kirkegaard, C., Hansen, J. M., and Sierskaek-Nielsen, K., Estimation of Serum Thyrotropin (TSH) and Stimulation with Thyrotropin-Releasing Hormone (TRH) in Thyroid Diseases. Acta Med. Scand. 196, 171-176, 1974.
10. Nicoloff, J. T., Spencer, C. A., The use and misuse of the sensitive thyrotropin assays. J. Clin. Endocrinol. Metab. 71, 553-558, 1990.
11. Taimela, E., Tähtelä, R., Koskinen, P., Nuutila, P., Forsström, J., Taimela, S., et al., Ability of two new thyrotropin (TSH) assays to separate hyperthyroid patients from euthyroid patients with low TSH. Clin. Chem. 40, 101-105, 1994.
12. Young, D., Effects of Drugs on Clinical Laboratory Tests. 3rd Edition, American Association for Clinical Chemistry Press. 1990.


#### TOSOH CORPORATION

3-8-2, Shiba, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220


#### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

#### TOSOH BIOSCIENCE, INC.

6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080, USA  
Phone: +1 650 615 49 70

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in a member state of the European Union

## AIA-PACK TSH 3rd-Gen CALIBRATOR SET

### INTENDED USE

The AIA-PACK TSH 3rd-Gen CALIBRATOR SET is intended for in vitro diagnostic use only for the calibration of the ST AIA-PACK TSH assay.

### SUMMARY AND EXPLANATION

The AIA-PACK TSH 3rd-Gen CALIBRATOR SET contains protein matrix with assigned levels of Thyroid Stimulating Hormone (TSH). Calibration should be performed in accordance with the schedule given in the Tosoh AIA System Operators' Manual.

### MATERIAL SUPPLIED

Cat. No.	0020394		
2 x 1 ml	AIA-PACK TSH 3rd-Gen CALIBRATOR (1)	0	μIU/ml
	Protein matrix containing no detectable concentration of TSH with sodium azide as a preservative.		
2 x 1 ml	AIA-PACK TSH 3rd-Gen CALIBRATOR (2)	0.2	μIU/ml (approx.)
	AIA-PACK TSH 3rd-Gen CALIBRATOR (3)	5.0	μIU/ml (approx.)
	AIA-PACK TSH 3rd-Gen CALIBRATOR (4)	25	μIU/ml (approx.)
	AIA-PACK TSH 3rd-Gen CALIBRATOR (5)	50	μIU/ml (approx.)
	AIA-PACK TSH 3rd-Gen CALIBRATOR (6)	110	μIU/ml (approx.)
	Protein matrix containing the assigned concentration of TSH (described on each vial) with sodium azide as a preservative.		

### WARNINGS AND PRECAUTIONS

1. The AIA-PACK TSH 3rd-Gen CALIBRATOR SET is intended for in vitro diagnostic use only.
2. Inspect the packaging and the exterior of the vials for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. This material contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always rinse with large volumes of water to prevent any azide accumulation
4. Although material derived from human blood is not used for these calibrators, it is recommended that this product be handled with the same precautions as those used for patient samples.
5. Do not use beyond the expiry date.
6. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

The AIA-PACK TSH 3rd-Gen CALIBRATOR SET is supplied ready for use. Bring calibrators to 18°-25°C for use.

### STORAGE AND STABILITY

- Always store the AIA-PACK TSH 3rd-Gen CALIBRATOR SET in an upright position at 2°-8°C when not in use.

- When stored unopened and refrigerated at 2°-8°C, the calibrator set will remain stable until the expiry date on the label.
- The calibrators should be used within 1 day of opening or reconstituting, provided the vials are kept tightly sealed and refrigerated at 2°-8°C.

### PROCEDURE

NOTE: Refer to the CALIBRATION PROCEDURE in the insert sheet of the ST AIA-PACK TSH. For additional procedural instructions regarding calibration, refer to the Tosoh AIA System Operators' Manual.

1. When using new calibrator lots, enter the calibrator concentration values and lot number on the software test file. (Refer to the Tosoh AIA System Operators' Manual for details.)
2. Load the appropriate amount of ST AIA-PACK TSH test cups on the instrument.
3. Add the appropriate amount of each calibrator to sample cups. (Refer to the instrument worksheet for the sample volume.)
4. Print a work list and place the sample cups in the position indicated.
5. Select START. Verify that the positions on the work list match the starting position on the screen.

### ASSIGNMENT OF VALUES

The AIA-PACK TSH 3rd-Gen CALIBRATOR SET contains assigned concentrations of TSH. The assigned value is determined on a lot-by-lot basis and is designed to provide an assay calibration range of 0.01 to 100 μIU/ml of TSH. The calibrators in this set are referred to WHO 2nd IRP 80/558.

### RESULTS

1. The mean rate for the AIA-PACK TSH 3rd-Gen CALIBRATOR (1) should be < 3.0 nmol/(L\*s).
2. Since there is a direct relationship between concentration and rate, the rate should rise as the concentration increases.
3. Repeat values should be within 10% range.

### LIMITATIONS

The AIA-PACK TSH 3rd-Gen CALIBRATOR SET is designed solely for use with ST AIA-PACK TSH assay procedures. Although the approximate value of the highest calibrator is 110 μIU/ml, the exact concentration may vary slightly. The assay specification, ASSAY RANGE HIGH, should be defined as the upper limit of the assay range, 100 μIU/ml.


### TOSOH CORPORATION

3-8-2, Shiba, Minato-ku  
TOKYO 105-8623 (JAPAN)  
Tel.: +81 (0)3 5427 5181 Fax: +81 (0)3 5427 5220

### TOSOH EUROPE N.V.

Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

**TOSOH BIOSCIENCE, INC.**

6000 Shoreline Court, Suite 101

South San Francisco, CA 94080, USA

Phone: +1 650 615 49 70

**NET**

Net Volume (after reconstitution  
for lyophilised material)

**Attention:** This IFU complies with IVD directive 98/79/EC and is intended for use by customers operating in the European Union

## AIA-PACK WASH CONCENTRATE

---

### INTENDED USE

The AIA-PACK WASH CONCENTRATE is intended for in vitro diagnostic use only for the AIA-PACK or the ST AIA-PACK Assay.

### MATERIAL SUPPLIED

Cat. No. 0020955  
4 x 100 ml AIA-PACK WASH CONCENTRATE  
Buffer solution with detergent and bacteriostatic agent.

### WARNINGS AND PRECAUTIONS

1. The AIA-PACK WASH CONCENTRATE is intended for in vitro diagnostic use.
2. Inspect the packaging and the exterior of the bottles for any sign of damage before use. If any damages are visible, contact your local TOSOH sales representative.
3. Do not use beyond the expiry date.
4. For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.

### PREPARATION OF REAGENTS

Add the entire content of the AIA-PACK WASH CONCENTRATE (100 ml) to approximately 2.0 l of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline, mix well, and adjust the final volume to 2.5 l.

Alternatively, add 40 ml of the AIA-PACK WASH CONCENTRATE (as measured with a graduated cylinder) to 800 ml of CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline. Mix well, and adjust the volume to 1000 ml with CAP Class I water or the clinical laboratory reagent water (formerly NCCLS Type I) defined by CLSI C3-A4 guideline.

### STORAGE AND STABILITY

- Always store the AIA-PACK WASH CONCENTRATE in an upright position at 2°-8°C when not in use.
- When stored unopened and refrigerated at 2°-8°C, the AIA-PACK WASH CONCENTRATE will remain stable until the expiry date on the label.
- When diluted the wash solution is stable for 30 days at 18°-25°C.

### PROCEDURE

Please refer to the insert of the AIA-PACK or the ST AIA-PACK and the Operators' Manual of the Tosoh AIA System analysers for detailed instructions.

**SAFETY INFORMATION:** not applicable.

**TOSOH EUROPE N.V.**  
Transportstraat 4  
B-3980 TESSENDERLO BELGIUM  
Tel.: +32 (0)13 66 88 30 Fax: +32 (0)13 66 47 49

**NET** Net Volume (after reconstitution for lyophilised material)

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#### LITERATŪRA

1. Carey MJ and Rodgers GM: Disseminated intravascular coagulation: Clinical and laboratory aspects., Am J Hematol, 59, 65-73 (1998).
2. Mammen EF: Disseminated intravascular coagulation (DIC), Clin Lab Sci, 13, 239-245 (2000).
3. Yu M, Nardella A and Pechet L: Screening tests of disseminated intravascular coagulation: guidelines for rapid and specific laboratory diagnosis., Crit Care Med, 28, 1777-1780 (2000).
4. Anderson DR and Wells PS: D-dimer for the diagnosis of venous thromboembolism., Curr Opin Hematol, 7, 296-301 (2000).

5. Perrier A, Desmarais S, Goehring C, de Moerloose P, Morabia A, Unger PF, et al.: D-dimer testing for suspected pulmonary embolism in outpatients., Am J Respir Crit Care Med, 156, 492-496 (1997).

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

#### LITERATŪRA

- Mukoyama M., et al., Human brain natriuretic peptide, a novel cardiac hormone. *Lancet*, 335, 801-802 (1990).
- Saito Y., et al., Brain natriuretic peptide is a novel cardiac hormone. *Biochem. Biophys. Res. Commun.*, 158, 360-368 (1989).
- Yoshimura M., et al., Hemodynamic, renal, and hormonal responses to brain natriuretic peptide infusion in patients with congestive heart failure. *Circulation*, 84, 1581-1588 (1991).
- Cowie M. R., et al., Clinical applications of B-type natriuretic peptide (BNP) testing. *Eur. Heart J.*, 24, 1710-1718 (2003).
- De Lemos J. A., et al., The prognostic value of B-type natriuretic peptide in patients with acute coronary syndromes. *N. Engl. J. Med.*, 345, 1014-1021 (2001).
- Murdoch D. R., et al., Titration of vasodilator therapy in chronic heart failure according to plasma brain natriuretic peptide concentration: randomized comparison of the hemodynamic and neuroendocrine effects of tailored versus empirical therapy. *Am. Heart J.*, 138, 1126-1132 (1999).
- Tsutamoto T., et al., Attenuation of compensation of endogenous cardiac natriuretic peptide system in chronic heart failure: prognostic role of plasma brain natriuretic peptide concentration in patients with chronic symptomatic left ventricular dysfunction. *Circulation*, 96, 509-516 (1997).
- McCullough P. A., et al., B-Type natriuretic peptide and clinical judgment in emergency diagnosis of heart failure: analysis from breathing not properly (BNP) multinational study. *Circulation*, 106, 416-422 (2002).
- Shimizu H., et al., Degradation of human brain natriuretic peptide (BNP) by contact activation of blood coagulation system. *Clin. Chim. Acta*, 305, 181-186 (2001).

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
- Berger, S. and Quim, J. L., Thyroid Function. In: Fundamentals of Clinical Chemistry, ed. N.W. Tietz. W.B. Saunders Co., Philadelphia, PA. 1976, pp. 824-848.
- Hoskins, R. G., The Thyroid-Pituitary Apparatus as a Servo (Feed-Back) Mechanism. J.Clin. Endocrinol. Metab. 9, 1429-1431, 1949.
- Reichlin, S. and Utiger, R. D., Regulation of the Pituitary-Thyroid Axis in Man: Relationship in TSH Concentration to Concentration of Free and Total Thyroxine in Plasma. J. Clin. Endocr. 27, 251-255, 1967.
- Burger, H. G., and Patel, Y. C., The Value of Serum Thyrotropin Measurement in the Diagnosis and Management of Hypothyroidism. Med. J. Aust. 2,293-297, 1972.
- Hershman, J. M., and Pittman, J. A., Utility of the Radioimmunoassay of Serum Thyrotropin in Man. Ann. Int. Med. 74, 481-490, 1971.
- Karlberg, B. E., and Almquist, S., Clinical Experience with Thyrotropin Releasing Hormone (TRH) Stimulation Test in Patients with Thyroid, Pituitary and Hypothalamic Disorders. Acta Endocr. 72, 697-713, 1973.
- Smith, S. R., Burton, D. W., and Izenstark, J. L., Levels of Serum Thyroid-Stimulating Hormone in Hyperthyroidism. Arch. Intern. Med. 139, 757-760, 1979.
- Hall, R., Ormston, B. J., Besser, G. M., Cryer, R. J., and McKendrick, M., The Thyrotropin-Releasing Hormone Test in Diseases of the Pituitary and Hypothalamus. Lancet 1 (7754), 759-762, 1972.
- Lauridsen, U. B., Deckert, T., Friis, T., Kirkegaard, C., Hansen, J. M., and Sierskaek-Nielsen, K., Estimation of Serum Thyrotropin (TSH) and Stimulation with Thyrotropin-Releasing Hormone (TRH) in Thyroid Diseases. Acta Med. Scand. 196, 171-176, 1974.
- Nicoloff, J. T., Spencer, C. A., The use and misuse of the sensitive thyrotropin assays. J. Clin. Endocrinol. Metab. 71, 553-558, 1990.
- Taimela, E., Tähtelä, R., Koskinen, P., Nuutila, P., Forsström, J., Taimela, S., et al., Ability of two new thyrotropin (TSH) assays to separate hyperthyroid patients from euthyroid patients with low TSH. Clin. Chem. 40, 101-105, 1994.
- Young, D., Effects of Drugs on Clinical Laboratory Tests. 3rd Edition, American Association for Clinical Chemistry Press. 1990.

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
1. Oppenheimer, J. H., Role of Plasma Proteins in the Binding, Distribution and Metabolism of the Thyroid Hormones. *New Eng. J. Med.*, 278:1153 (1986).
2. Prince, M. P. and Ramsden, D. B., A New Theoretical Description of the Binding of Thyroid Hormones by Serum Proteins. *Clin. Endocrinol.* 7:307 (1977).
3. Chopra, I. J., et al, Pathway of Metabolism of Thyroid Hormones. *Rec. Prog. Horm. Res.*, 34: 521 (1978).
4. Ekins, R. P., Free Hormones in Blood: The Concept and the Measurement. *J. Clin. Immunoassay*, 7:163 (1984).
5. Penney, M. D. and O'Sullivan, D. J., Total or Free Thyroxine as a Primary Test of Thyroid Function. *Clin.Chem.*, 33:170 (1987).
6. Robbins J. and Rall, J. E., Hormones in Blood, 1:577 (1967).
7. Ramsden, D. B., et al, Serum Free Thyroxine Concentrations in Normal Euthyroid Subjects and Ones with High Serum Thyroxine Binding Globulin Concentration., *Clin. Chem. Acta*, 130:211 (1983).
8. Wilke, T. J., Diagnostic Value of Three Methods for Assessing Free Thyroxin in Pregnancy. *Ann. Clin. Biochem.*, 20:60 (1983).
9. Sparre, L. S., et al, Oestrogen and Thyroxine Binding Globulin Levels in Early Normal Pregnancy. *Acta Endocrinol.*, 130:298 (1987).
10. Anderson, B. G., Free Thyroxine in Serum in Relation to Thyroid Function, *JAMA*, 203:577 (1968).
11. Young, D., *Effects of Drugs on Clinical Laboratory Tests*. 3rd Edition, Washington, DC, American Association for Clinical Chemistry Press (1990).

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