

# Cytopro® Cytocentrifuge

## Series 2




### **CYTOPRO® ROTOR IS THE WORLD'S MOST ADVANCED CYTOCENTRIFUGE ROTOR**

Automates the cytocentrifugation of cell suspensions to facilitate the sedimentation of cells onto a slide that gets a five-to nine-fold increase in the number of cells on a slide.<sup>1</sup>

# DELIVERING SUPERB CELL RECOVERY WITH TROUBLE-FREE OPERATION



## QUICK AND EASY LOADING & UNLOADING



### PERFECT ALIGNMENT

The pad is always correctly aligned with the tunnel.

### 7 MM DIAMETER CELL AREA


A larger, centered cell deposit area makes screening easier and more sensitive.

### COST-EFFICIENT


Plastic Cytopro chambers are disposable or reusable after decontamination.

### SIMPLE IN SITU FIXATION

Load a few drops of sample into the tunnel port, then load fixative into the sample well. As the rotor accelerates, the cells are fixed after they sediment onto the slide.








**Cytopro® Rotor is the world's most advanced cytocentrifuge rotor with a patented flow control ring that gets a five-to nine -fold increase in the number of cells on a slide.<sup>1</sup> Adaptable with most cytocentrifuge systems and available as a fully integrated system.**

Aerospray® stainers, when used with the Cytopro® Rotor, become dual-purpose stainers equipped with staining AND cytocentrifugation capabilities. One instrument efficiently does the work of two stand-alone instruments - maximizing space and eliminating clutter issues at a fraction of the cost.

## AVAILABLE SAMPLE CHAMBERS


Patented flow control increases cell yield and delivers consistent results.

### Standard Chamber



Combines ease of use with superb cell recovery.

### Dual Chamber



Places two samples on the same slide for increased efficiency and economy.


### Cytopro Magnum®




Accepts sample sizes up to 6 mL.

<sup>1</sup> Huppmann, A. R., Rheingold, S. R., Bailey, L. C., Helfrich, M., & Choi, J. K. (2012). Detection of Leukemic Lymphoblasts in CSF Is Instrument-Dependent. American Journal of Clinical Pathology, 137(5), 795-799. doi:10.1309/ajcpp4hq2fxqezil


## CYTOPRO'S UNIQUE ERGONOMIC DESIGN LEAVES NO CHANCE OF SMEARING CELLS




Urine




Pleural Fluid



Lung Aspirate



Hematology



Buffy Coat

### SIMPLE TO PROGRAM


- Cutting-edge microprocessor control and user programmability provides programming versatility
- In addition to manual programming, you can store up to twelve frequently used combinations of rotation speed (100-2000 rpm), run time (1-99 minutes), and rotor acceleration (Low, Medium or High)
- Programs are easy to retrieve or clear
- The user-friendly display visually confirms inputs, current settings, or operation cycle status

### ROTOR PROTECTS THE SAMPLE AND THE USER

- The Cytopro rotor uses centrifugal force and three proprietary chamber designs to sediment cells onto slides with maximum safety and user convenience
- The shallow profile makes sample ID easy, and loading and unloading samples is convenient in a biological safety hood
- Simply pipette sample fluid directly into the sample well
- Slides and samples load securely into the rotor
- After the run, the chamber release lever cleanly separates the chamber from the slide
- Both the slide and the sample chamber can then be easily removed

### CYTOPRO CHAMBERS

- Loading and pad alignment are simple and foolproof
- The semi-transparent sample chamber allows visual confirmation of correct sample volume
- Patented flow control delivers consistent results. Dual addition ports allow pre-wetting—for excellent cell recovery
- Because damaged spring clips can affect recovery, we use reliable coil springs
- Standard chambers come in single and dual configurations to suit every need



# SPECIFICATIONS


## Cytopro® Cytocentrifuge Model 7622

Chamber Volume	
Dual Chamber	0.3 mL per well
Single Chamber	0.5 mL
Magnum Chamber	6.0 mL
Rotor Slide Capacity	Up to 8 slides and Cytopro chambers
Rotor Speed	100 to 2000 rpm (±5%), user programmable
Run Time	Programmable 1 to 99 mins
Dimensions (W x H x D)	43 cm x 25 cm x 54 cm (17" x 10" x 21")
Clearance Height (lid open)	58 cm (23")
Weight	10.4 kg (22.9 lb.)
Electrical Requirements	100 to 240 V (±10%) @ 50-60 Hz
Fuses	2.0 Amp time-lag (Type T) (5 x 20 mm) (2 required)
Power Consumption	200 Watts maximum
Safety	Lid Interlock: lid must be closed for operation, and is locked down during rotor rotation

# ORDERING INFORMATION

Catalog #	Description
MODEL 7622	Cytopro® Cytocentrifuge
AC-160	Cytopro Rotor
SS-111	Cytopro Cytopad, fast, white, package of 100
SS-112	Cytopro Cytopad, slow, tan, package of 100
SS-113	Sample Chambers w/ fast, white Cytopads, caps, 48pk
<b>SS-114</b>	<b>Sample Chambers w/ slow, tan Cytopads, caps, 48pk</b>
SS-117	Microscope slides, 1/2G, uncoated, for Cytopro
SS-118	Microscope slides, 1/2G, coated, for Cytopro
SS-213	Dual Sample Chambers w/ fast,white pads,caps,48
SS-214	Dual Sample Chambers w/ slow,tan pads,caps,pk/48
SS-217	Microscope Slides, 1/2G, uncoated, Cytopro Dual
SS-218	Microscope Slides, 1/2G, coated, Cytopro Dual
SS-232	Microscope slides, 1/2G, uncoated, for Cytopro Magnum
SS-233	Microscope slides, 1/2G, coated, for Cytopro
SS-234	Cytopro Magnum Sample Chambers, caps, 24pk

*\*Please check additions to catalog numbers made in your country*





# WORLD LEADER IN AEROSPRAY® STAINING TECHNOLOGY AND CYTOCENTRIFUGATION

## 50 YEARS

### “Proven Performance, Time After Time.”

ELITechGroup has been in operation for nearly 50 years and distributes to over 100 countries. Proud to be the World Leader in Aerospray staining technology and Cytocentrifugation, offering the world's most reliable solutions in performance and reliability. ELITechGroup is confident that this automated technology will elevate the performance of any laboratory and facilitate superior results.

Join the thousands of laboratories utilizing the world's best staining technology and take your laboratory to a new level of efficiency and performance.



#### WORLDWIDE OFFICES

Please contact your sales representative for terms, conditions and product availability in your country.

#### Headquarters

Australia  
Belgium & Luxembourg  
Brazil  
France  
Italy  
New Zealand

#### +33 1 41 45 07 10

T: +61 1800 815 098  
T: +32 9 282 05 31  
T: +55 27 3025 1415  
T: +33 4 83 36 10 82  
T: +39 02 48 40 35 42  
T: +64 800 555 611

Serbia  
Switzerland  
The Netherlands  
UK  
United States

T: +381 11 2467119  
T: +41 26 663 86 60  
T: +31 313 430 500  
T: +44 1442 869320  
T: +1 435 752 6011



**ELITechGroup**

EMPOWERING IVD

**AEROSPRAYSTAINING.com**  
**info@elitechgroup.com**