AQT90 FLEX

TnT Test Kit

RADIOMETER 

Contents

en	TnT Test Kit	4
bg	Набор тестове за TnT	8
cs	Testovací souprava TnT.....	14
da	TnT Test Kit	18
de	TnT Testaufrüstungssatz	23
el	ΚΙΤ Εξέτασης TnT	28
es	Kit de Test de TnT.....	33
et	TnT mõõtekomplekt	38
fi	TnT-testaussarja	42
fr	Kit de tests TnT.....	47
hr	Testni komplet TnT.....	52
hu	TnT tesztkészlet.....	57
it	Kit per il test TnT	62
lt	TnT testo rinkinys.....	67
lv	TnT testēšanas komplekts.....	71
nl	TnT Testkit.....	76
no	TnT-testkit	81
pl	TnT Test Kit	85
pt	Kit de Testes TnT.....	90
ro	Kit de testare pentru TnT	95
ru	Набор тестов на TnT.....	100
sk	Súprava testov TnT	105
sv	TnT-testkit	110
tr	TnT Test Kiti	115
zh	TnT 测试卡	120

en

TnT Test Kit

Proprietary name and common name: Cardiac troponin T (TnT)

Catalog no. 942-940

Intended use

For *in vitro* diagnostic use.

The TnT **Test** is an *in vitro* diagnostic assay for the quantitative determination of cardiac TnT in EDTA or lithium-heparin whole-blood or plasma specimens on the AQT90 FLEX analyzer in point of care and laboratory settings. It is indicated for use as an aid in the diagnosis of myocardial infarction and in the risk stratification of patients with acute coronary syndromes with respect to their relative risk of mortality. Current international guidelines should always be followed when interpreting TnT results.

The TnT **CAL** Cartridges are *in vitro* diagnostic reagents intended for the calibration adjustment of the TnT Test on the AQT90 FLEX analyzer by establishing points of reference to estimate TnT values.

NOTICE: The cartridges are ready for use.

Summary

Test – TnT is a protein of the muscular thin filaments and part of the troponin-tropomyosin complex in striated muscle. Troponin complex consists of three components: troponin C (TnC), troponin T (TnT) and troponin I (TnI). This ternary troponin complex is a calcium-sensitive molecular apparatus that regulates the interaction of actin and myosin [1,2]. Both TnI and TnT exist in skeletal and cardiac isoforms. The fact that the cardiac forms are extremely specific to the heart tissue forms the basis for their increased use as specific markers of cardiac muscle damage. The major form of circulating TnT is as free, uncomplexed TnT [3].

Plasma concentrations of TnT in healthy individuals are extremely low or undetectable, which constitutes the basis for the usefulness of TnT in the detection of myocardial infarction (MI). TnT is released from necrotic cardiac myocytes into the bloodstream within hours (~4-8 hours) after the onset of chest pain. The peak TnT concentration is generally reached in 12-48 hours. Plasma TnT levels remain elevated several days after severe MI [4,5,6,7].

The exceptional cardiac specificity and high clinical sensitivity of the troponins shown in clinical studies have markedly increased their role in the diagnostics of MI. TnI and TnT are currently regarded as the most preferred biochemi-

cal markers of myocardial damage and the cornerstones of MI diagnostics [8,9].

Calibration adjustment – Calibration adjustment must be performed every time the Test Cartridge lot changes. The assay has been factory-calibrated using an 9-point calibration and the factory-defined calibration data is provided in the form of a barcode. This data is used as a lot-specific reference curve, which is adjusted to the correct instrument-specific signal level by using a CAL cartridge with the same lot number. The results obtained using the CAL Cartridge are combined with the factory-defined calibration data to produce a lot- and instrument-specific calibration curve for the assay. When a CAL cartridge with a new lot number is loaded, the analyzer requests the factory-defined calibration data so that it can perform the calibration adjustment. If the user does not load a CAL Cartridge when the Test Cartridge lot changes, the analyzer will request it as soon as it detects a Test Cartridge with a new lot number. Test Cartridges with a new lot number can only be used after performing a calibration adjustment with a CAL Cartridge with the same lot number.

NOTICE: Liquid quality controls must be run after every calibration adjustment.

See the *AQT90 FLEX operator's manual* for details.

Principles of the procedure

Test – The TnT assay is an immunoassay based on a concept in which all the specific reagents are provided in a dry stable form within an assay cup. Monoclonal anti-TnT antibodies have been preimmobilized to the surface of the cup, and a separating layer and tracer antibodies have been added on top of the capture antibodies. The separating layer prevents the direct contact of the capture and tracer antibodies in storage. In the assay process, the sample and the assay solution are automatically added to the cup containing the assay-specific reagents. During the incubation period, the tracer and capture antibodies form a sandwich complex with TnT present in the sample. After the incubation, the assay cup is washed with the assay solution and dried. Then the signal from the tracer antibody labeled with europium is measured by means of time-resolved fluorometry (TRF), directly from the dry surface of the assay cup. The concentration of the analyte is directly proportional to the measured europium signal. The measured signal is converted to a concentration using the calibration curve stored in the memory of the instrument [10,11].

The TnT assay is intended for use only with the AQT90 FLEX analyzer. This is a fully automatic continuous-access analyzer utilizing immunoassay technology and time-resolved fluorometric detection.

The analyzer performs sampling from a manually loaded capped primary sample tube, after which the automatic processing takes place. The analyzer performs all assay steps automatically and reports quantitative results. The AQT90 FLEX analyzer includes a touch screen for routine access and advanced user access, such as calibration adjustment and parameter-specific definitions.

See the *AQT90 FLEX operator's manual* for detailed information about the functions of the AQT90 FLEX analyzer.

Calibration adjustment – See the *AQT90 FLEX operator's manual*.

Contents:	<ul style="list-style-type: none"> • Ten TnT Test Cartridges for performing 160 tests • One TnT CAL Cartridge for performing one calibration adjustment • One "TnT Factory-defined Calibration Data" sheet with barcode
Components:	<ul style="list-style-type: none"> • Each Test Cartridge contains 16 analyte-specific assay cups. • Each CAL Cartridge contains eight analyte-specific background cups and eight cups with added antigen.
Storage:	Store unopened at 2-8 °C
Shelf life:	If the protective pouch containing the cartridge remains sealed until just before the cartridge is inserted into the analyzer, it is stable until the end of the month stated in the Expiration date shown on the pouch label.
On-board stability:	See the label on the Test- and CAL Cartridge pouches.
Chemical composition:	
Reactive ingredients:	<ul style="list-style-type: none"> • Mouse monoclonal TnT capturing and tracing antibodies: 850 ng/cup • Recombinant TnT (only in the eight CAL cups with added antigen): approx. 0.12 ng/cup

Other ingredients:	<ul style="list-style-type: none"> • Bovine serum albumin • Bovine γ-globulin • Mouse IgG (blocker substances for heterophilic interferences)
--------------------	--

Main buffer component:	TRIS (Tris Hydroxymethylaminomethane)
------------------------	--

Toxic ingredient:	Sodium azide < 0.1 %
-------------------	----------------------

Data in the cartridge barcodes:

• Checksum	• Cartridge ID
• Barcode type	• Expiration date
• Parameter code	• Barcode version
• Lot no.	• On-board stability

NOTICE: The barcode is lot-specific.

CAUTION – Sensitive to humidity

Do not use any cartridge that has not been kept according to storage specifications or that has been damaged or taken out of its protective pouch.

Product calibrator traceability

Calibration of the AQT90 FLEX TnT assay is comparable to calibration of the TnT assay for the Roche Elecsys 2010 system.

Specimen collection and preparation

Collect blood samples by venipuncture. Use whole-blood or plasma samples with either EDTA or lithium heparin as anticoagulant.

Specimens should be tested as soon as possible after the collection. Whole blood samples must be analyzed within three hours of collection, when kept at ambient temperature (22–25 °C). If this cannot be achieved, it is recommended that plasma is separated and stored cool (2–8 °C for no longer than 24 hours) or frozen. Repeated freezing and thawing should be avoided. **Do not freeze whole-blood specimens.**

Frozen patient specimens (EDTA and lithium-heparin plasma) should be brought to room temperature slowly and gently mixed by hand by inverting the sample tube five times. We recommend that thawed plasma samples are centrifuged to prevent possible micro-clots being aspirated with the sample. The analyzer mixes the sample prior to sample analysis.

The minimum sample volume is 2 mL (independent of the number of AQT90 FLEX tests requested).

Tubes supported by the analyzer: See the *AQT90 FLEX reference manual* for details.

Procedures

NOTICE: To keep the analyzer ready for measurement at all times, the power of the analyzer must always be on.

See the *AQT90 FLEX operator's manual* for handling procedures and detailed descriptions of the following procedures.

Briefly, the procedures consist of the following steps:

Analyzing a patient sample

1. Place the capped sample tube in the sample inlet.
2. Select tube type.
3. Select the required test panel(s) and/or parameter(s) to be tested.
4. Start the analysis.

Calibration adjustment procedure

NOTICE: Calibration adjustment must be performed every time a Test Cartridge with a new lot number is used.

1. Find the CAL Cartridge in the Test Kit box with the new lot number.
2. Load it into the cartridge inlet with the barcode facing downwards.
NOTICE: The analyzer does not accept cartridges after the expiration date printed on its pouch label.

WARNING – Risk of incorrect results

Do not use a cartridge that has already been used on another analyzer, otherwise you risk obtaining incorrect results.

3. Scan the barcode from the "TnT Factory-defined Calibration Data" sheet provided in the Test Kit box.
4. Start the calibration adjustment procedure.
Results are marked on the screen as successful or failed.

Quality control procedure

Liquid quality control measurements should be performed:

- after calibration adjustment has been performed on a new lot of Test Cartridges
 - in accordance with local, state and/or federal regulations
- See the *AQT90 FLEX operator's manual* for details.

Results

Results are displayed on the screen and may be printed as they become available. If the analyzer detects a fault condition, no result will be given and an error message will be displayed.

Performance characteristics

Analytical specificity

The analytical specificity of the AQT90 FLEX TnT assay was determined by studying the cross-reactivity with skeletal muscle troponin T, cardiac troponin I and troponin C at a concentration of 100 µg/L. There was no detectable cross-reaction with skeletal muscle troponin T, cardiac troponin I or troponin C.

Analytical sensitivity and measuring range

The limit of detection has been determined to be 0.01 µg/L.

The reportable range of the assay is 0.010-25 µg/L.

Reference values

Whole blood (lithium-heparin and EDTA) and plasma (lithium-heparin and EDTA) were obtained from 260 apparently healthy individuals (132 women and 128 men) and analyzed using the AQT90 FLEX TnT assay. The 99th percentile was determined to be 0.017 µg/L.

NOTICE: This value should only be used as an example and each laboratory should establish its own reference range.

Imprecision

Within-day and total imprecision was determined by analyzing plasma pools over 20 days, twice a day, 2 replicates per run.

Within run imprecision				Total imprecision			
Sample	Mean µg/L	SD µg/L	% CV	Sample	Mean µg/L	SD µg/L	% CV
Plasma pool 1	0.027	0.0023	8.5	Plasma pool 1	0.027	0.0026	9.6
Plasma pool 2	0.21	0.0073	3.5	Plasma pool 2	0.21	0.0116	5.6
Plasma pool 3	12.0	0.343	2.9	Plasma pool 3	12.0	0.644	5.4

The concentration giving the CV of 10 % of the AQT90 FLEX TnT assay is approximately 0.03 µg/L.

Clinical sensitivity and specificity

Several studies show that cardiac Troponin T concentrations can be used as markers in the diagnosis of myocardial infarction and in the risk stratification of patients with acute coronary syndromes with respect to their relative risk of mortality [12,13,14,15,16,17,18].

The recommended decision threshold value is 0.030 µg/L.

NOTICE: This threshold value should only be used as an example. Each laboratory should establish its own decision threshold level(s).

Hook effect

No hook effect was found when TnT concentrations up to 100 µg/L were measured.

Carry over

Carry over from a sample with high troponin T value (100 µg/L) to an adjacent negative sample was determined to be < 100 ppm.

Interfering substances

Hemolytic, lipemic and icteric samples were shown not to interfere with the TnT assay.

To minimize any potential interferences from HAMA (human anti-mouse antibodies), mouse IgG was added as a blocker substance. Tests showed that the blocker had the desired effect.

WARNING – Risk of incorrect results

Always interpret with caution any assay results that are inconsistent with other available clinical data and a patient's medical history. If samples contain heterophilic antibodies, immunoassays could give falsely elevated or depressed results. Mouse IgG has been added to this assay as a blocker substance to minimize potential interference from HAMA (human anti-mouse antibodies).

The following interfering substances were found to have no notable effect on the AQT90 FLEX TnT assay (interference < 20 %). The interference was tested by using plasma samples with 0.04 µg/L of cardiac troponin T and spiked with the interfering substance at the following concentration (about 5 times the upper therapeutic range):

Substance	Tested concentration mg/L
Acetaminophen	12
Acetylsalicylic acid	200
Allopurinol	240
Amiodarone	20
Ampicillin	200
Ascorbic acid	24
Amlodipine Besylate	4
Atenolol	40
Atorvastatin	32
Biotin	0.05
Caffeine	30
Captopril	40
Chloramphenicol	50
Clopidogrel Bisulfate	30
Creatinine	110
Cyclosporine	40
Diclofenac	60
Digitoxin	60
Digoxin	0.004
Diltiazem	120
Dipyridamole	30
Dopamine	130
Enalapril Maleate	16
Erythromycin	100
Furosemide	16
Heparin	8 IU/mL
Hydralazine	20
Hydrochlorothiazide	20
Indomethacin	16
Isosorbide dinitrate	4
Lisinopril	16
Lovastatin	16
L-thyroxine	46 ng/L
Methyldopa	100
Milrinone lactate	2.4
Nicotine	1.6
Nicotinic acid	40
Nifedipine	36
Nitrofurantoin	40
Nitroglycerin	0.16
Oxytetracycline	100

Substance	Tested concentration mg/L
Phenytoin	40
Probenecid	200
Procainamide	20
Propranolol	64
Quinidine	20
Simvastatin	32
Spironolactone	40
Sulfamethoxazole	320
Theophylline	40
Trimethoprim	64
Verapamil	96
Warfarin	4

Method comparison

The AQT90 FLEX TnT assay (y) was compared to the commercially available cardiac troponin T assay TnT for the Roche Elecsys 2010 system (x) using lithium-heparin plasma samples in the range of 0-25 µg/L (with the AQT90 FLEX TnT assay).

The correlation (Passing-Bablok) was found to be:

$$y = 0.980x - 0.0004; R^2 = 0.9928 (n = 214)$$

Patents

Manufactured under license from Roche Diagnostics. US Patent No. US6,376,206.

In addition, one or more of the following patents and patent applications may apply:

- US Patent No.: US6,429,026
- European Patent No.: EP892927
- Japanese Patent No.: JP3537829
- German Patent No.: DE69729186
- US Patent Application No.: US2008-0226508
- European Patent Application No.: EP2130050
- International Patent Application No.: WO2008/113352

References

See *References* inside the back cover of this insert.

Explanation of graphical symbols

See the *Explanation of graphical symbols* table at the end of this insert.

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References

1. Filatov VL, Katrukha AG, Bulargina TV, Gusev NB. Troponin: Structure, properties, and mechanism of functioning. *Biochemistry (Moscow)* 1999; 64(9): 969-85.
2. Farah CS, Reinach FC. The troponin complex and regulation of muscle contraction. *FASEB J* 1995; 9: 755-67.
3. Gaze DC, Collinson PO. Multiple molecular forms of circulating cardiac troponin: analytical and clinical significance. *Ann Clin Biochem* 2008; 45: 349-55.
4. Wu AHB. Cardiac markers. In: Wild D, ed. *The immunoassay handbook*. London: Nature Publishing Group, 2001: 623-34.
5. Collinson PO, Boa FG, Gaze DC. Measurement of cardiac troponins. *Ann Clin Biochem* 2001; 38(5): 423-49.
6. Apple FS, Falahati A, Paulsen PR, Miller EA, Sharkey SW. Improved detection of minor ischemic myocardial injury with measurement of serum cardiac troponin I. *Clin Chem* 1997; 43: 2047-51.
7. Manuel Del Rey JM, Madrid AH, Valino JM *et al*. Cardiac troponin I and minor cardiac damage: biochemical markers in a clinical model of myocardial lesions. *Clin Chem* 1998; 44: 2270-76.
8. Wu AHB, Apple FS, Gibler WB, Jesse RL, Warshaw MM, Valdes R. National Academy of Clinical Biochemistry Standards of Laboratory Practice: Recommendations for the use of cardiac markers in coronary artery diseases. *Clin Chem* 1999; 45: 1104-21.
9. Thygesen K, Alpert JS, Antman E, Bassand JP. Myocardial infarction redefined – a consensus document of the Joint European Society of Cardiology/ American College of Cardiology Committee for the redefinition of myocardial infarction. *JACC* 2000; 36(3): 959-69.
10. Lövgren T, Pettersson K. Time-resolved fluoroimmunoassay, advantages and limitations. In: Van Dyke K, Van Dyke R, eds. *Luminescence Immunoassay and Molecular Applications*. Boca Raton: CRC Press, 1990; 233-53.
11. Lövgren T, Meriö L, Mitrunen K *et al*. One-step all-in-one dry reagent immunoassays with fluorescent europium chelate label and time-resolved fluorometry. *Clin Chem* 1996; 42: 1196-201.
12. Katus HA, Remppis A, Looser S, *et al*. Enzyme linked immunoassay of cardiac troponin T for the detection of acute myocardial infarction in patients. *Mol Cell Cardiol* 1989; 21(7): 1349-53.
13. Katus HA, Scheffold T, Remppis A, Zehelein J. Proteins of the troponin complex. *Laboratory Medicine* 1992; 23(5): 311-17.
14. Hamm CW, Ravkilde J, Gerhardt W, Jorgensen P, Peheim E, Ljungdahl L, *et al*. The prognostic value of serum troponin T in unstable angina. *N Engl J Med* 1992; 327(3): 146-50.
15. Ohmann EM, *et al*. Risk stratification with admission cardiac troponin T levels in acute myocardial ischemia. *N Engl J Med* 1996; 335: 1333-34.
16. Christenson RH, Duh SH, Newby LK, Ohman EM, Califf RM, Granger CB, *et al*. Cardiac troponin T and cardiac troponin I: relative values in short-term risk stratification of patients with acute coronary syndromes. *Clin Chem* 1998; 44(3): 494-501.
17. Lindahl B, Diderholm E, Lagerqvist B, Venge P, Wallentin L, the FRISC II investigators. Mechanisms behind the prognostic value of Troponin T in unstable coronary artery disease: a FRISC II substudy. *J Am Coll Card* 2001; 38: 979-86.
18. Aviles RJ, Askari AT, Lindahl B, Wallentin L, Jia G, Ohman EM, Mahaffey KW, Newby LK, Califf RM, Simoons ML, Topol EJ, Lauer MS. Troponin T levels in patients with acute coronary syndromes, with or without renal dysfunction. *N Engl J Med* 2002; 346: 2047-52.

Made in Finland for:
Radiometer Medical ApS
Åkandevvej 21
2700 Brønshøj
Denmark
www.radiometer.com

Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, PICO and CLINITUBES are trademarks of Radiometer Medical ApS.

© 2011 Radiometer Medical ApS, 2700 Brønshøj, Denmark.
All rights reserved. 995-639. 201103F.

RADIOMETER 