## AQT90 FLEX

### LQC Multi-CHECK level 2 (MC2)

<b>en</b>
<b>Name</b>
Proprietary name: LQC Multi-CHECK level 2 (MC2)
Catalog no.: 944-232
<b>Intended use</b>
For <i>in vitro</i> diagnostic use.
For use as a liquid quality control serum (LQC) to monitor the precision of laboratory testing procedures for TnI, CKMB, Myo and NT-proBNP on the AQT90 FLEX analyzer.

<b>Contents:</b>	<ul style="list-style-type: none"><li>Six (6) tubes</li> <li>Slip of paper with two barcodes in a plastic sleeve on the outside slip cover</li></ul>
<b>Components:</b>	Each tube contains 0.65 mL of a mixture of Bio-Rad's Liqueich Cardiac Markers Plus Control LT, enough for five (5) tests.
<b>Storage:</b>	Store unopened below –18 °C
<b>Shelf life:</b>	Stable until the expiry date printed on the label when stored unused below –18 °C. Once thawed the product is stable: <ul style="list-style-type: none"><li>for five (5) days if stored unused between 2-8 °C</li> <li>for a maximum of 2 hours if stored unused at room temperature</li></ul>
<b>Limitation of use:</b>	

Product stability is dependent on its storage temperature. "Storage" and "Shelf life" specifications provide details. Do not use the product if it has not been kept according to these specifications, or if it has been damaged

<b>Data in all barcodes:</b>	
<ul style="list-style-type: none"><li>Checksum</li> <li>Barcode type</li> <li>Area code</li></ul>	<ul style="list-style-type: none"><li>Level name</li> <li>Lot no.</li> <li>Expiry date</li></ul>

#### Extra data in Barcodes 1 and 2:

Upper and lower limits of Radiometer's (-R-) default control range for all parameters.

<b>Assigned value:</b>	
This is the center value of a control range. It is a lot-specific and parameter-specific value.	

<b>Quality control procedure</b>	
The enclosed barcodes have to be scanned into the analyzer during this procedure. See the AQT90 FLEX operator's manual for details.	

<b>Reagent preparation</b>	
Remove the LQC solution from the freezer, let it thaw at room temperature for 15 minutes and use immediately.	

<b>WARNING – Risk of infection</b>	
Handle and discard used LQC tubes as biohazardous waste. Failure to do so may result in exposure to potentially infected material. Disposal of all waste should be in accordance with local regulations.	
<b>WARNING – Potentially infectious material</b>	
Each human donor unit used to manufacture this control was tested using FDA-accepted methods and found non-reactive for Hepatitis B Surface Antigen (HBsAg), antibody to Hepatitis C (HCv) and antibody to HIV-1/HIV-2. This product may also contain other human source material for which there are no approved tests. In accordance with good laboratory practice, all human source material should be considered potentially infectious and handled with the same precautions used with patient specimens.	

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