

INSTRUCTIONS FOR USE

HbA1c

VITROS Chemistry Products HbA1c Reagent Kit

Direct % Glycated Hemoglobin

REF 684 2905

Rx ONLY

This device has significant negative interference with fetal hemoglobin (HbF). HbA1c results are invalid for patients with abnormal amounts of HbF including those with known Hereditary Persistence of Fetal Hemoglobin. Refer to the Limitations of the Procedure section of this Instruction For Use for details.

Intended Use

For *in vitro* diagnostic use only.

VITROS Chemistry Products HbA1c reagent is used on VITROS 5,1 FS Chemistry System, VITROS 4600 Chemistry System and the VITROS 5600 Integrated System for the quantitative determination of percent glycated hemoglobin A1c (DCCT/NGSP) and mmol/mol hemoglobin A1c (IFCC) in human whole blood.

The test is to be used as an aid in diagnosis of diabetes, as an aid in identifying patients who may be at risk for developing diabetes mellitus, and for the monitoring of long-term blood glucose control in individuals with diabetes mellitus.

Summary and Explanation of the Test

HbA1c is the major species of glycohemoglobin found in human blood, and is formed by the non-enzymatic glycation of the free amino group at the N-terminal of the hemoglobin A₀ β-chain.¹ The percentage of HbA1c (%A1c) is widely used as a retrospective index of glycemic control in diabetic patients.² Studies have shown that glycemic control results in reduction in diabetic complications.^{3, 4} The glycation of hemoglobin depends on both the lifespan of the red blood cell and the blood glucose concentration. Since the rate of glycated hemoglobin formation is directly proportional to the concentration of glucose in the blood, the level of glycated hemoglobin is an index of the blood glucose concentrations over an extended period of time (~6 to 8 weeks).⁵

Principles of the Procedure

The determination of % glycated hemoglobin (HbA1c) is performed using the VITROS Chemistry Products HbA1c Reagent Kit in conjunction with the VITROS Chemistry Products Calibrator Kit 31 on the VITROS 5,1 FS/4600 Chemistry Systems and the VITROS 5600 Integrated System. The VITROS Chemistry Products HbA1c Reagents are two dual chambered packages containing ready-to-use liquid reagents. Whole blood samples are hemolyzed on the VITROS 5,1 FS/4600 Chemistry Systems and the VITROS 5600 Integrated System. The concentration of HbA1c and total Hb are measured in the hemolyzed samples, controls and calibrators.

Hemoglobin A1c and Hemoglobin

Whole blood samples are hemolyzed on the VITROS 5,1 FS/4600 Chemistry Systems and the VITROS 5600 Integrated System. Calibrators, controls and hemolyzed whole blood samples are mixed with Reagent 1 containing anti-HbA1c antibody to form a soluble antigen-antibody complex. Hemoglobin in the hemolyzed whole blood is converted with Reagent 1 to a hematin derivative that is measured bichromatically at 340 nm and 700 nm. Unbound anti-HbA1c antibody reacts with polyhapten (hexapeptide-glycan, A1c Reagent 2) to form an insoluble antibody-polyhapten immune complex, which is measured turbidimetrically at 340 nm. After a calibration has been performed for each reagent lot, the hemoglobin A1c and Hb concentrations in each unknown sample can be determined using the stored calibration curves and the measured absorbance obtained in the assay of the hemolyzed sample.

%A1c

%A1c is a derived test calculated from the quantitative measurements of hemoglobin and hemoglobin A1c.

HbA1c (SI Units)⁶

HbA1c (mmol of HbA1c/mol of Hb) in SI units is a derived test calculated from the quantitative measurements of hemoglobin and hemoglobin A1c.

HbA1c

Direct % Glycated Hemoglobin

INSTRUCTIONS FOR USE


Warnings and Precautions

Test Type and Conditions

HbA1c and Hb

| Test Type | VITROS System | Approximate Incubation Time | Temperature | Wavelength | Reaction Sample Volume |
|-----------|---------------------|-----------------------------|-----------------|-----------------------------------|------------------------|
| End-Point | 5600, 4600, 5, 1 FS | 7 minutes | 37 °C (98.6 °F) | (HbA1c) 340 nm (Hb) 340/700 nm | 5.0 µL |

Reaction Scheme



Warnings and Precautions

For *in vitro* diagnostic use only.

WARNING:

Take care when handling materials and samples of human origin. Since no test method can offer complete assurance that infectious agents are absent, consider all clinical specimens, controls, and calibrators potentially infectious. Handle specimens, solid and liquid waste, and test components in accordance with local regulations and CLSI Guideline M29⁷ or other published biohazard safety guidelines.

For specific warnings and precautions for calibrators, quality control materials and other components, refer to the Instructions for Use for the appropriate VITROS product or to other manufacturer's product literature.

Reagents

Reactive Ingredients

HbA1c Reagent 1 (R1): HbA1c antibody (ovine serum) ≥ 0.5 mg/mL

HbA1c Reagent 2 (R2): HbA1c Polyhapten ≥ 8 µg/mL

DIL5 (R1): Tetradecyltrimethylammonium bromide (TTAB) $< 1\%$ (w/v).

Other Ingredients

HbA1c Reagent 1 (R1): Buffers, surfactant, stabilizers, and preservatives

HbA1c Reagent 2 (R2): Buffers, surfactant, stabilizers, and preservatives

DIL5 (R1): Surfactants, stabilizers, and preservatives

Reagent Handling

Caution:

Do not use reagent packs with damaged or incompletely sealed packaging.

- Inspect the packaging for signs of damage.
- Be careful when opening the outer packaging with a sharp instrument so as to avoid damage to the individual product packaging.
- The reagents are supplied ready for use.
- Avoid agitation, which may cause foaming or the formation of bubbles.

INSTRUCTIONS FOR USE

Reagent Storage and Stability

HbA1c

Direct % Glycated Hemoglobin

Reagent Preparation

1. Remove from refrigerated storage.
2. Immediately load into Supply 3.

Reagent Storage and Stability

VITROS Chemistry Products HbA1c Reagents are stable until the expiration date on the carton when they are stored and handled as specified. Do not use beyond the expiration date. For additional information, refer to "Calibration".

IMPORTANT: Do Not Freeze.

| Reagent | Storage Condition | | Stability |
|----------|-------------------|-------------------|-----------------------|
| Unopened | Refrigerated | 2–8 °C (36–46 °F) | Until expiration date |
| Opened | On-analyzer | System turned on | ≤ 28 days |
| | On-analyzer | System turned off | ≤ 30 minutes |

Verify performance with quality control materials after reloading reagents that have been removed from Supply 3 and stored for later use.

Specimen Collection, Preparation and Storage

Specimens Recommended

EDTA Whole Blood

- K₂EDTA
- K₃EDTA

IMPORTANT: Certain collection devices have been reported to affect other analytes and assays⁸. Owing to the variety of specimen collection devices available, Ortho-Clinical Diagnostics is unable to provide a definitive statement on the performance of its products with these devices. Confirm that your collection devices are compatible with this assay.

Specimens Not Recommended

- Heparinized Whole Blood
- Sodium Fluoride/Potassium Oxalate Whole Blood

Whole Blood

IMPORTANT: Samples must be transferred to sample cups for analysis. Use of the original blood tube could cause separation of the blood during instrument handling and may lead to invalid results.

Specimen Collection and Preparation

Collect specimens using standard laboratory procedures.^{9, 10}

Patient Preparation

No special patient preparation is necessary.

Special Precautions

None.

Specimen Handling and Storage

- Handle and store specimens in stoppered containers to avoid contamination and evaporation.
- Samples should be fully resuspended prior to use.
- Samples should be free of clots.

Specimen Storage and Stability

| Storage | Temperature | Stability |
|--------------------------|---------------------|-----------|
| Original Specimen | | |
| Room Temperature | 18–28 °C (64–82 °F) | ≤ 3 days |
| Refrigerated | 2–8 °C (36–46 °F) | ≤ 3 days |

HbA1c

Direct % Glycated Hemoglobin

INSTRUCTIONS FOR USE

Testing Procedure

Testing Procedure

Materials Provided

VITROS Chemistry Products HbA1c Reagent Kit

Materials Required but Not Provided

- VITROS Chemistry Products FS Reconstitution Diluent
- VITROS Chemistry Products Calibrator Kit 31
- Quality control materials, such as VITROS Chemistry Products %A1c Performance Verifiers I and II
- Isotonic saline

Operating Instructions

- Check reagent inventories at least daily to ensure that quantities are sufficient for the planned workload.
- For additional information, refer to the operating instructions for your system.

IMPORTANT: *Bring all fluids and samples to room temperature, 18–28 °C (64–82 °F), prior to analysis.*

Specimen Testing

IMPORTANT: *Whole blood specimens must be placed in sample cups for analysis. Do not place an evaporation cap on the sample cup.*

Note: If using VITROS MicroSample Cups for analysis, you must use the VITROS FS Adaptor for MicroSample Cups.

1. Fully resuspend whole blood samples **immediately** prior to transfer to sample cups.
2. Samples should be free of clots.
3. For the system to properly mix the whole blood sample prior to aspiration, pipette 350 µL into a sample cup. Sample volumes less than 300 µL or greater than 400 µL may result in improper mixing of the whole blood sample.

Sample Programming

To report results in %A1c units, select the **%A1c** button in sample programming or via a LIS download.

To report results in HbA1c (mmol/mol) SI units select the **HbA1c** button in sample programming or via a LIS download.

Sample Dilution

On-Analyzer Sample Dilution

On-analyzer dilution is not recommended.

Manual Sample Dilution

If either component concentration (Hb or HbA1c) exceeds the system's measuring (reportable) range:

1. Manually dilute the original, whole blood sample with an equal volume of saline. Ensure original sample is resuspended prior to dilution.
2. Resuspend the diluted whole blood sample by gentle inversion and re-analyze.
3. Do not multiply the results by the dilution factor to obtain an estimate of the original sample's %A1c or HbA1c (mmol/mol) value. Hemoglobin and glycated hemoglobin are both diluted; therefore it is not required to multiply by the dilution factor in order to obtain an estimate of the original sample's value.

Calibration

Required Calibrators

- VITROS Chemistry Products Calibrator Kit 31

Calibrator Preparation, Handling, and Storage

Refer to the Instructions for Use for VITROS Chemistry Products Calibrator Kit 31.

Calibration Procedure

Refer to the operating instructions for your system.

When to Calibrate

Calibrate:

- When the reagent lot changes.

INSTRUCTIONS FOR USE

Calibration

HbA1c

Direct % Glycated Hemoglobin

- When critical system parts are replaced due to service or maintenance.
- When government regulations require.
For example, in the USA, CLIA regulations require calibration or calibration verification at least once every six months.

The VITROS HbA1c assay may also need to be calibrated:

- If quality control results are consistently outside acceptable range.
- After certain service procedures have been performed.

For additional information, refer to the operating instructions for your system.

Calculations

Components

HbA1c

Absorbance is measured at 340 nm after a fixed incubation time. After a calibration has been performed for each reagent lot, HbA1c concentration in the unknown samples can be determined using the stored calibration curve and the measured absorbance obtained in the assay of each sample.

Hb

Absorbance is measured bichromatically at 340 nm and 700 nm after a fixed incubation time. After a calibration has been performed for each reagent lot, Hb concentration in the unknown samples can be determined using the stored calibration curve and the measured absorbance obtained in the assay of each sample.

Derived Tests

Two reporting units are supported by your system.

- HbA1c in SI units
- %A1c in NGSP units

HbA1c

SI Units

$$\text{HbA1c (mmol/mol)} = \frac{\text{HbA1c [g/dL]}}{\text{Hb [g/dL]}} \times 1000$$

%A1c

NGSP units

$$\%A1c = (\text{IFCC}^* \times 0.09148) + 2.152$$

$$^*\text{IFCC} = \text{HbA1c (mmol/mol) SI units}$$

Validity of a Calibration

Calibration parameters are automatically assessed by the system against a set of quality parameters detailed in the Review Assay Data screen (found via Options → Review/Edit Calibrations → Review Assay Data). Failure to meet any of the pre-defined quality parameters results in a failed calibration. The calibration report should be used in conjunction with quality control results to determine the validity of a calibration.

Linearity

Evaluation of the linearity of the VITROS Chemistry Products HbA1c assay was performed based on CLSI EP06-A. ¹¹ For each assay component (Hb, HbA1c, %A1c), a low pool and a high pool were prepared with the assay component outside the measuring range. The low and high concentration pools were mixed to give 16 further pools of intermediate concentrations. The VITROS HbA1c assay was tested using each VITROS System.

The VITROS HbA1c assay is linear through the following ranges:

- Linear range for Hb: 4.972–30.318 g/dL (component)
- Linear range for HbA1c: 0.080–2.530 g/dL (component)
- Linear range for %A1c: 3.034–15.444% (NGSP derived test)
- Linear range for HbA1c: 9.6–145.3 mmol/mol (SI derived test)

HbA1c

Direct % Glycated Hemoglobin

INSTRUCTIONS FOR USE

Quality Control

%A1c, NGSP

| Analyzer | Intercept | Slope | r ² | Concentration Range Tested |
|---------------|-----------|--------|----------------|----------------------------|
| VITROS 5,1 FS | 0.005 | 0.9995 | 0.999 | 2.39–17.35 %A1c |
| VITROS 4600 | 0.030 | 0.9957 | 0.997 | 3.03–15.44 %A1c |
| VITROS 5600 | 0.012 | 0.9988 | 0.999 | 2.89–15.79 %A1c |

mmol/mol, SI

| Analyzer | Intercept | Slope | r ² | Concentration Range Tested |
|---------------|-----------|--------|----------------|----------------------------|
| VITROS 5,1 FS | 0.044 | 0.9995 | 0.999 | 2.6–166.1 mmol/mol |
| VITROS 4600 | 0.223 | 0.9957 | 0.997 | 9.6–145.3 mmol/mol |
| VITROS 5600 | 0.099 | 0.9988 | 0.999 | 8.1–149.1 mmol/mol |

Measuring (Reportable) Range

- Measuring range for Hb: 6.0–22.0 g/dL (component)
- Measuring range for HbA1c: 0.2–2.5 g/dL (component)
- Measuring range for derived tests:

| %A1c NGSP* Units | HbA1c (mmol/mol) SI** Units |
|------------------|-----------------------------|
| 4–14 | 20–130 |

* NGSP – National Glycohemoglobin Standardization Program

** SI – Système International

Traceability of Calibration

The values assigned to the VITROS Chemistry Products Calibrator Kit 31 for %A1c are traceable to the IFCC (International Federation of Clinical Chemistry and Laboratory Medicine) Reference Method.¹ The derived result (%A1c) is calculated from the individual quantitative results for hemoglobin (Hb) and glycated hemoglobin (HbA1c). NGSP results are derived from the International Federation of Clinical Chemistry (IFCC) HbA1c (mmol/mol) SI units using the Master Equation; ^{6, 12, 13}

$$\% \text{NGSP} = (0.09148 \times \text{IFCC}) + 2.152.$$

Quality Control

Quality Control Material Selection

IMPORTANT: VITROS Chemistry Products %A1c Performance Verifiers are recommended for use with the VITROS 5,1 FS/4600 Chemistry and VITROS Integrated Systems. Evaluate the performance of other commercial control fluids for compatibility with this assay before using for quality control.

- Control materials other than VITROS Chemistry Products %A1c Performance Verifiers may show a difference when compared with other %A1c Performance Verifiers methods if they:
 - Depart from a true human matrix.
 - Contain high concentrations of preservatives, stabilizers, or other nonphysiological additives.
- Use of either pre-hemolyzed or whole blood quality control materials is acceptable. Ensure during sample programming that the appropriate configuration is selected.

Quality Control Procedure Recommendations

IMPORTANT: Either pre-hemolyzed or whole blood quality control materials are acceptable for this assay. Ensure that the system is configured correctly for the type of control used. To configure quality control fluids, go to QC→Define Controls→Wh Blood. If using a pre-hemolyzed control, check the box next to "Pretreated." If using a whole blood control, leave the box next to "Pretreated" unchecked.

- Choose control levels that check the clinically relevant range.
- Analyze quality control materials in the same manner as pretreated whole blood, before or during patient sample processing.
- To verify system performance, analyze control materials:
 - After calibration.
 - According to local regulations or at least once each day that the assay is being performed.
 - After specified service procedures are performed. Refer to the operating instructions for your system.
- If control results fall outside your acceptable range, investigate the cause before deciding whether to report patient results.

INSTRUCTIONS FOR USE

Results

HbA1c

Direct % Glycated Hemoglobin

- For general quality control recommendations, refer to *Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline-Third Edition*¹⁴ or other published guidelines.
- For additional information, refer to the operating instructions for your system.

Quality Control Material Preparation, Handling, and Storage

Refer to the Instructions for Use for VITROS Chemistry Products %A1c Performance Verifiers I and II or to other manufacturer's product literature.

Results

Reporting Units

The VITROS 5600 Integrated, VITROS 4600 Chemistry and the VITROS 5,1 FS Chemistry Systems can be programmed to report VITROS Chemistry Products HbA1c assay results in NGSP equivalent %A1c units (by selecting the **%A1c** button) and/or HbA1c (mmol/mol) SI units (by selecting the **HbA1c** button).

Unit Conversion from HbA1c (mmol/mol) SI Units to NGSP %A1c

NGSP %A1c results are derived from the International Federation of Clinical Chemistry (IFCC) HbA1c (mmol/mol) SI units using the Master Equation:^{6, 12, 13} $\%A1c\ (NGSP) = (IFCC^* \times 0.09148) + 2.152$.

where $IFCC^* = HbA1c\ (mmol/mol)\ SI\ Units$

HbA1c in (mmol/mol) SI units can be calculated from the %A1c units using the Master Equation to yield the following equation: $HbA1c\ (mmol/mol) = (\%A1c\ (NGSP) \times 10.93) - 23.50$

Limitations of the Procedure

Known Interferences

Glycated HbF is not detected as it does not contain the glycated β -chain that characterizes %A1c. However, HbF is measured in the total Hb assay and as a consequence, specimens containing high amounts of HbF (>7%) may result in lower than expected mmol/mol HbA1c values (IFCC) and % A1c values (NGSP).

For substances that were tested and did not interfere, refer to "Specificity."

Other Limitations

- The test is intended for the calculation of %A1c or mmol A1c/mol Hb. For patient samples, do not report the individual results for total hemoglobin (Hb) or hemoglobin A1c (HbA1c).
- The result from this or any other diagnostic test should be used and interpreted only in the context of the overall clinical picture.
- The design of VITROS HbA1c assay allows its use in patients with hemoglobin variants S, C, D and E traits. However, in case of disease condition such as HbSS, HbCC, or HbSC, even if the assay does not interfere with their variant, these patients may suffer anemia, increased red blood cell turnover, and transfusion requirements which can adversely affect HbA1c as a marker of long-term glycemic control. Any cause of shortened erythrocyte survival or decrease in mean erythrocyte age (e.g. recovery from acute blood loss, hemolytic anemia) will reduce exposure of erythrocytes to glucose with a consequent decrease in mmol/mol HbA1c values (IFCC) and %A1c values (DCCT/NGSP). Therefore, health care providers should not use the HbA1c test for these patients. Alternative forms of testing such as glycated serum protein or glycated albumin should be considered for these patients.^{15, 16}
- Hemoglobin F does not contain β chains and is not detected by the VITROS HbA1c assay, but it is measured in the total Hb assay. The calculation of mmol A1c/mol Hb values (IFCC) and %A1c values (DCCT/NGSP) relative to total amount of Hb in the sample, produces falsely decreased HbA1c values in specimens with high amounts of HbF.¹⁷
- Hemoglobin A1c should not be used to diagnose gestational diabetes. It reflects the average blood glucose levels over the preceding 3 months (the average life of a red blood cell), and therefore may be falsely low during pregnancy or any other condition associated with recent onset of hyperglycemia.¹⁸
- In very rare cases of rapidly evolving type 1 diabetes, the increase of HbA1c values might be delayed compared to the acute increase in glucose concentrations. In these conditions diabetes mellitus must be diagnosed based on plasma glucose concentrations and/or the typical clinical symptoms.¹⁹
- HbA1c does not provide a measure of glycemic variability or hypoglycaemia. Glycemic control in these type 1 or type 2 diabetic patients prone to glycemic variability is best evaluated by a combination of results from self-monitoring and HbA1c.²⁰
- Hemoglobin A1c should not be used to diagnose diabetes mellitus in patients with malignancies or severe chronic hepatic and renal disease.²¹
- Hemoglobin A1c testing should not replace glucose testing for type 1 diabetes, in pediatric patients and in pregnant women.

HbA1c

Direct % Glycated Hemoglobin

INSTRUCTIONS FOR USE

Expected Values

- Hemoglobin A1c should not be used to diagnose diabetes mellitus in patients with a hemoglobinopathy but normal red cell turnover (e.g. sickle cell trait).
- Hemoglobin A1c should not be used in patients with homozygous sickle cell trait, hemolytic anemia, or other hemolytic diseases and recent significant or chronic blood loss.

Expected Values

Reference Interval

The expected normal HbA1c range in adults is 4.0–6.0% (NGSP) or 20–42 mmol/mol (IFCC units).²² Each laboratory should confirm the validity of these intervals for the population it serves. The Standards of Medical Care in Diabetes - 2015 recommend to diagnose diabetes using a HbA1c method that is NGSP-certified and standardized to the DCCT assay and a cut point of HbA1c $\geq 6.5\%$.^{19, 20}

| | Current* | IFCC traceable methods |
|------------------------------------|----------|------------------------|
| Reference interval (non-diabetics) | 4–6% | 20–42 mmol/mol |

* Refer to methods aligned to the US National Glycohemoglobin Standardization Program.

To provide clinical decision guidance for laboratories reporting results in %A1c (NGSP units), the NGSP interval was published by the National Glycohemoglobin Standardization Program (NGSP).²⁰ The reference interval is applicable to methods traceable to the Diabetes Control and Complications Trial (DCCT).³

| %A1c (NGSP) | Interpretation |
|-------------|------------------|
| ≥ 6.5 | Action Suggested |

To provide clinical decision guidance for laboratories reporting results in HbA1c mmol/mol (SI units), the following values were calculated from the NGSP %A1c cut point using the Master Equation.^{6, 12} The HbA1c cut point is ≥ 48 mmol/mol (SI Units).

| Calculated HbA1c (mmol/mol) SI Units | Interpretation |
|--------------------------------------|------------------|
| ≥ 48 | Action Suggested |

Performance Characteristics

Limit of Quantitation/Detection

The Limit of Detection (LoD) for VITROS Chemistry Products HbA1c Reagent is 2.580% (NGSP Units), determined consistent with CLSI document EP17²³ and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 500 determinations, 5 low-level samples. The Limit of Blank (LoB) is 2.396% (NGSP Units) based on 100 determinations with 5 blank samples. The Limit of Quantitation (LoQ) determined consistent with the guidelines in CLSI document EP17²³, based on 100 determinations with the 5 LoD pools using the classical Westgard model is 2.370% (NGSP Units). The associated bias and precision components were 0.064% NGSP and 0.121% NGSP respectively, which yielded a Westgard model TE of 0.3060 and a %TE of 12.57. Because the determined LoQ is less than the LoD, the final LoQ was set equal to the LoD of 2.580% (NGSP).

The Limit of Detection (LoD) for the Hb component of the VITROS Chemistry Products HbA1c Reagent is 0.312 g/dL, determined consistent with CLSI document EP17²³ and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 500 determinations, 5 low-level samples. The Limit of Blank (LoB) is 0.186 g/dL based on 100 determinations with 5 blank samples. The Limit of Quantitation (LoQ) determined consistent with the guidelines in CLSI document EP17²³, based on 126 determinations with six LoQ pools; and a precision goal of 8.3% using the functional sensitivity method is 2.117 g/dL.

The Limit of Detection (LoD) for the HbA1c component of the VITROS Chemistry Products HbA1c Reagent is 0.072 g/dL, determined consistent with CLSI document EP17²³ and with proportions of false positives (α) less than 5% and false negatives (β) less than 5%; based on 500 determinations, 5 low-level samples. The Limit of Blank (LoB) is 0.042 g/dL based on 100 determinations with 5 blank samples. The Limit of Quantitation (LoQ) determined consistent with the guidelines in CLSI document EP17²³, based on 500 determinations with 5 LoD pools; and a precision goal of 14.6% using the functional sensitivity method is 0.133 g/dL.

INSTRUCTIONS FOR USE

Performance Characteristics

HbA1c

Direct % Glycated Hemoglobin

Limit of Blank, Limit of Detection and Limit of Quantitation – Hb

| LoB | | LoD* | | LoQ | |
|-------|------|-------|------|-------|-------|
| g/dL | g/L | g/dL | g/L | g/dL | g/L |
| 0.186 | 1.86 | 0.312 | 3.12 | 2.117 | 21.17 |

* Proportions of false positives (α) and false negatives (β) were less than 5%; based on 500 determinations, with 5 low-level samples.

Limit of Blank, Limit of Detection and Limit of Quantitation – HbA1c

| LoB | | LoD* | | LoQ | |
|-------|------|-------|------|-------|------|
| g/dL | g/L | g/dL | g/L | g/dL | g/L |
| 0.042 | 0.42 | 0.072 | 0.72 | 0.133 | 1.33 |

* Proportions of false positives (α) and false negatives (β) were less than 5%; based on 500 determinations, with 5 low-level samples.


Limit of Blank, Limit of Detection and Limit of Quantitation – %A1c

| LoB | | LoD* | | LoQ | |
|-------|---------------|-------|---------------|-------|---------------|
| NGSP | SI (mmol/mol) | NGSP | SI (mmol/mol) | NGSP | SI (mmol/mol) |
| 2.396 | 2.67 | 2.580 | 4.68 | 2.580 | 4.68 |

* Proportions of false positives (α) and false negatives (β) were less than 5%; based on 500 determinations, with 5 low-level samples.

Method Comparison

The plots and table show the results of a method comparison study with whole blood samples analyzed on the VITROS 5600 Integrated System with those analyzed using a secondary reference laboratory of the NGSP. The table also shows the results of comparisons with whole blood samples on the VITROS 4600 Chemistry System and the VITROS 5,1 FS Chemistry System with a secondary reference laboratory of the NGSP, and the results of comparisons with whole blood samples on the VITROS 4600 Chemistry System and the VITROS 5,1 FS Chemistry System with the VITROS 5600 Integrated System. The testing followed CLSI Protocol EP9. ²⁴ The VITROS HbA1c assay is certified by the National Glycohemoglobin Standardization Program and is traceable to the Diabetes Control and Complications Trial Reference Method.



HbA1c

Direct % Glycated Hemoglobin

INSTRUCTIONS FOR USE

Performance Characteristics

| | n | Slope | +/- 95% range of slope | Correlation Coefficient | NGSP Units | | | | HbA1c SI Units (mmol/mol) | | | |
|--|-----|-------|------------------------------|----------------------------|---------------------|----------------|----------------------------------|------|---------------------------|----------------|----------------------------------|------|
| | | | | | Range of Samples | Inter- cept | +/- 95% range of intercept | Sy.x | Range of Samples | Inter- cept | +/- 95% range of intercept | Sy.x |
| 5600 vs. comparative method [*] | 357 | 1.00 | 0.99–1.02 | 0.996 | 4.2–13.4 | -0.06 | -0.13–0.02 | 0.18 | 22–123 | -0.4 | -0.9–0.2 | 1.9 |
| 4600 vs. comparative method [*] | 357 | 1.02 | 1.01–1.03 | 0.995 | 4.2–13.4 | -0.14 | -0.22– -0.06 | 0.19 | 22–123 | -0.9 | -1.5– -0.3 | 2.1 |
| 5,1 FS vs. comparative method [*] | 357 | 1.01 | 0.99–1.02 | 0.996 | 4.2–13.4 | -0.08 | -0.16–0.01 | 0.18 | 22–123 | -0.5 | -1.1–0.1 | 1.9 |
| 4600 vs. 5600 | 125 | 1.00 | 0.99–1.01 | 0.998 | 4.3–13.1 | -0.01 | -0.09–0.06 | 0.10 | 23–120 | -0.1 | -0.7–0.4 | 1.1 |
| 5,1 FS vs. 5600 | 126 | 0.99 | 0.97–1.00 | 0.998 | 4.3–13.8 | 0.01 | -0.09–0.12 | 0.12 | 23–127 | -0.2 | -1.0–0.6 | 1.3 |

^{*} Secondary Reference Laboratory of the NGSP

Precision

Precision was evaluated with quality control materials (hemolysate and whole blood-based) for 20-days and whole blood patient samples for 4-days on the VITROS 5,1 FS Chemistry System, VITROS 4600 Chemistry System and VITROS 5600 Integrated System following CLSI Protocol EP5. ²⁵

These results are guidelines. Variables such as instrument maintenance, environment, reagent storage/handling, control material reconstitution, and sample handling can affect the reproducibility of test results.

VITROS 5,1 FS Chemistry System – Combined (%A1c, NGSP)

| Mean HbA1c | Repeatability | | Between Run | | Between Day | | Between Lot | | Between Analyzer | | Total | | N | Days |
|----------------------------------|---------------|------|----------------|------|----------------|------|----------------|------|---------------------|------|-------|------|-----|------|
| | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | | |
| VITROS PV1 [*] 5.65% | 0.033 | 0.58 | 0.037 | 0.65 | 0.039 | 0.69 | 0.045 | 0.80 | 0.085 | 1.50 | 0.115 | 2.03 | 756 | 21 |
| Control [*] 4.96% | 0.026 | 0.52 | 0.032 | 0.64 | 0.029 | 0.58 | 0.043 | 0.87 | 0.063 | 1.27 | 0.091 | 1.83 | 792 | 22 |
| Control [*] 6.38% | 0.048 | 0.75 | 0.052 | 0.82 | 0.050 | 0.78 | 0.057 | 0.89 | 0.084 | 1.32 | 0.133 | 2.09 | 756 | 21 |
| Control [*] 8.17% | 0.065 | 0.80 | 0.041 | 0.50 | 0.065 | 0.80 | 0.081 | 0.99 | 0.178 | 2.18 | 0.220 | 2.69 | 792 | 22 |
| Control [*] 11.85% | 0.147 | 1.24 | 0.092 | 0.78 | 0.119 | 1.00 | 0.102 | 0.86 | 0.058 | 0.49 | 0.241 | 2.03 | 720 | 20 |
| Patient ^{**} 5.02% | 0.023 | 0.46 | 0.031 | 0.62 | 0.032 | 0.64 | 0.033 | 0.66 | 0.094 | 1.87 | 0.111 | 2.21 | 144 | 4 |
| Patient ^{**} 6.51% | 0.038 | 0.58 | 0.038 | 0.58 | 0.042 | 0.64 | 0.053 | 0.81 | 0.099 | 1.52 | 0.131 | 2.01 | 144 | 4 |
| Patient ^{**} 8.16% | 0.054 | 0.66 | 0.071 | 0.87 | 0.031 | 0.38 | 0.064 | 0.78 | 0.107 | 1.31 | 0.157 | 1.92 | 144 | 4 |
| Patient ^{**} 11.93% | 0.171 | 1.43 | 0.173 | 1.45 | 0.090 | 0.75 | 0.114 | 0.96 | 0.100 | 0.84 | 0.301 | 2.52 | 144 | 4 |

Within Lab precision (Total) was determined using three lots of reagents on three analyzers on each platform, using a single calibration per lot.

^{*} Control Fluid

^{**} Whole Blood

INSTRUCTIONS FOR USE

Performance Characteristics

HbA1c

Direct % Glycated Hemoglobin

VITROS 4600 Chemistry System – Combined (%A1c, NGSP)

| Mean HbA1c | Repeatability | | Between Run | | Between Day | | Between Lot | | Between Analyzer | | Total | | N | Days |
|----------------------------------|---------------|------|-------------|------|-------------|------|-------------|------|------------------|------|-------|------|-----|------|
| | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | | |
| VITROS PV1 [*] 5.63% | 0.035 | 0.62 | 0.021 | 0.37 | 0.029 | 0.52 | 0.033 | 0.59 | 0.043 | 0.76 | 0.074 | 1.32 | 756 | 21 |
| Control [*] 4.94% | 0.028 | 0.57 | 0.017 | 0.34 | 0.024 | 0.49 | 0.033 | 0.67 | 0.028 | 0.57 | 0.059 | 1.19 | 792 | 22 |
| Control [*] 6.35% | 0.050 | 0.79 | 0.045 | 0.71 | 0.034 | 0.54 | 0.055 | 0.87 | 0.044 | 0.69 | 0.103 | 1.62 | 756 | 21 |
| Control [*] 8.07% | 0.071 | 0.88 | 0.037 | 0.46 | 0.046 | 0.57 | 0.043 | 0.53 | 0.054 | 0.67 | 0.115 | 1.43 | 792 | 22 |
| Control [*] 11.89% | 0.177 | 1.49 | 0.069 | 0.58 | 0.109 | 0.92 | 0.095 | 0.80 | 0.000 | 0.00 | 0.239 | 2.01 | 720 | 20 |
| Patient ^{**} 4.98% | 0.027 | 0.54 | 0.025 | 0.50 | 0.016 | 0.32 | 0.012 | 0.24 | 0.036 | 0.72 | 0.055 | 1.10 | 144 | 4 |
| Patient ^{**} 6.47% | 0.049 | 0.76 | 0.019 | 0.29 | 0.020 | 0.31 | 0.034 | 0.53 | 0.040 | 0.62 | 0.077 | 1.19 | 144 | 4 |
| Patient ^{**} 8.06% | 0.090 | 1.12 | 0.064 | 0.79 | 0.032 | 0.40 | 0.052 | 0.65 | 0.044 | 0.55 | 0.134 | 1.66 | 144 | 4 |
| Patient ^{**} 11.96% | 0.157 | 1.31 | 0.142 | 1.19 | 0.117 | 0.98 | 0.064 | 0.54 | 0.000 | 0.00 | 0.250 | 2.09 | 144 | 4 |

Within Lab precision (Total) was determined using three lots of reagents on three analyzers on each platform, using a single calibration per lot.

^{*} Control Fluid

^{**} Whole Blood

VITROS 5600 Integrated System – Combined (%A1c, NGSP)

| Mean HbA1c | Repeatability | | Between Run | | Between Day | | Between Lot | | Between Analyzer | | Total | | N | Days |
|----------------------------------|---------------|------|-------------|------|-------------|------|-------------|------|------------------|------|-------|------|-----|------|
| | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | | |
| VITROS PV1 [*] 5.66% | 0.036 | 0.64 | 0.013 | 0.23 | 0.020 | 0.35 | 0.031 | 0.55 | 0.041 | 0.72 | 0.067 | 1.18 | 756 | 21 |
| Control [*] 4.98% | 0.034 | 0.68 | 0.011 | 0.22 | 0.017 | 0.34 | 0.058 | 1.17 | 0.000 | 0.00 | 0.070 | 1.41 | 792 | 22 |
| Control [*] 6.37% | 0.057 | 0.90 | 0.029 | 0.46 | 0.034 | 0.53 | 0.033 | 0.52 | 0.041 | 0.64 | 0.089 | 1.40 | 756 | 21 |
| Control [*] 8.10% | 0.072 | 0.89 | 0.017 | 0.21 | 0.030 | 0.37 | 0.054 | 0.67 | 0.037 | 0.46 | 0.103 | 1.27 | 792 | 22 |
| Control [*] 11.89% | 0.178 | 1.50 | 0.055 | 0.46 | 0.067 | 0.56 | 0.018 | 0.15 | 0.066 | 0.56 | 0.209 | 1.76 | 720 | 20 |
| Patient ^{**} 5.04% | 0.027 | 0.54 | 0.013 | 0.26 | 0.021 | 0.42 | 0.032 | 0.64 | 0.020 | 0.40 | 0.053 | 1.05 | 144 | 4 |
| Patient ^{**} 6.53% | 0.042 | 0.64 | 0.026 | 0.40 | 0.019 | 0.29 | 0.038 | 0.58 | 0.032 | 0.49 | 0.073 | 1.12 | 144 | 4 |
| Patient ^{**} 8.14% | 0.062 | 0.76 | 0.062 | 0.76 | 0.000 | 0.00 | 0.044 | 0.54 | 0.032 | 0.39 | 0.104 | 1.28 | 144 | 4 |
| Patient ^{**} 12.04% | 0.173 | 1.44 | 0.089 | 0.74 | 0.088 | 0.73 | 0.064 | 0.53 | 0.070 | 0.58 | 0.234 | 1.94 | 144 | 4 |

Within Lab precision (Total) was determined using three lots of reagents on three analyzers on each platform, using a single calibration per lot.

^{*} Control Fluid

^{**} Whole Blood

HbA1c

Direct % Glycated Hemoglobin

INSTRUCTIONS FOR USE

Performance Characteristics

VITROS 5,1 Chemistry System – Combined (mmol/mol, SI Units)

| Mean mmol/mol | Repeatability | | Between Run | | Between Day | | Between Lot | | Between Analyzer | | Total | | N | Days |
|---------------------------------|---------------|-----|-------------|-----|-------------|-----|-------------|-----|------------------|-----|-------|-----|-----|------|
| | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | | |
| VITROS PV1 [*] 38.3 | 0.36 | 0.9 | 0.40 | 1.0 | 0.43 | 1.1 | 0.49 | 1.3 | 0.93 | 2.4 | 1.26 | 3.3 | 756 | 21 |
| Control [*] 30.7 | 0.28 | 0.9 | 0.35 | 1.1 | 0.32 | 1.0 | 0.47 | 1.5 | 0.69 | 2.2 | 0.99 | 3.2 | 792 | 22 |
| Control [*] 46.2 | 0.52 | 1.1 | 0.57 | 1.2 | 0.55 | 1.2 | 0.62 | 1.3 | 0.92 | 2.0 | 1.45 | 3.1 | 756 | 21 |
| Control [*] 65.7 | 0.71 | 1.1 | 0.45 | 0.7 | 0.71 | 1.1 | 0.89 | 1.4 | 1.95 | 3.0 | 2.40 | 3.7 | 792 | 22 |
| Control [*] 106.0 | 1.61 | 1.5 | 1.01 | 1.0 | 1.30 | 1.2 | 1.11 | 1.0 | 0.63 | 0.6 | 2.63 | 2.5 | 720 | 20 |
| Patient ^{**} 31.3 | 0.25 | 0.8 | 0.34 | 1.1 | 0.35 | 1.1 | 0.36 | 1.2 | 1.03 | 3.3 | 1.21 | 3.9 | 144 | 4 |
| Patient ^{**} 47.7 | 0.42 | 0.9 | 0.42 | 0.9 | 0.46 | 1.0 | 0.58 | 1.2 | 1.08 | 2.3 | 1.43 | 3.0 | 144 | 4 |
| Patient ^{**} 65.6 | 0.59 | 0.9 | 0.78 | 1.2 | 0.34 | 0.5 | 0.70 | 1.1 | 1.17 | 1.8 | 1.72 | 2.6 | 144 | 4 |
| Patient ^{**} 106.9 | 1.87 | 1.7 | 1.89 | 1.8 | 0.98 | 0.9 | 1.25 | 1.2 | 1.09 | 1.0 | 3.29 | 3.1 | 144 | 4 |

Within Lab precision (Total) was determined using three lots of reagents on three analyzers on each platform, using a single calibration per lot.

^{*} Control Fluid

^{**} Whole Blood

VITROS 4600 Chemistry System – Combined (mmol/mol, SI Units)

| Mean mmol/mol | Repeatability | | Between Run | | Between Day | | Between Lot | | Between Analyzer | | Total | | N | Days |
|---------------------------------|---------------|-----|-------------|-----|-------------|-----|-------------|-----|------------------|-----|-------|-----|-----|------|
| | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | | |
| VITROS PV1 [*] 38.0 | 0.38 | 1.0 | 0.23 | 0.6 | 0.32 | 0.8 | 0.36 | 0.9 | 0.47 | 1.2 | 0.81 | 2.1 | 756 | 21 |
| Control [*] 30.5 | 0.31 | 1.0 | 0.19 | 0.6 | 0.26 | 0.9 | 0.36 | 1.2 | 0.31 | 1.0 | 0.64 | 2.1 | 792 | 22 |
| Control [*] 45.9 | 0.55 | 1.2 | 0.49 | 1.1 | 0.37 | 0.8 | 0.60 | 1.3 | 0.48 | 1.0 | 1.13 | 2.5 | 756 | 21 |
| Control [*] 64.7 | 0.78 | 1.2 | 0.40 | 0.6 | 0.50 | 0.8 | 0.47 | 0.7 | 0.59 | 0.9 | 1.26 | 1.9 | 792 | 22 |
| Control [*] 106.5 | 1.93 | 1.8 | 0.75 | 0.7 | 1.19 | 1.1 | 1.04 | 1.0 | 0.00 | 0.0 | 2.61 | 2.5 | 720 | 20 |
| Patient ^{**} 30.9 | 0.30 | 1.0 | 0.27 | 0.9 | 0.17 | 0.6 | 0.13 | 0.4 | 0.39 | 1.3 | 0.60 | 1.9 | 144 | 4 |
| Patient ^{**} 47.1 | 0.54 | 1.1 | 0.21 | 0.4 | 0.22 | 0.5 | 0.37 | 0.8 | 0.44 | 0.9 | 0.84 | 1.8 | 144 | 4 |
| Patient ^{**} 64.6 | 0.98 | 1.5 | 0.70 | 1.1 | 0.35 | 0.5 | 0.57 | 0.9 | 0.48 | 0.7 | 1.46 | 2.3 | 144 | 4 |
| Patient ^{**} 107.2 | 1.72 | 1.6 | 1.55 | 1.4 | 1.28 | 1.2 | 0.70 | 0.7 | 0.00 | 0.0 | 2.73 | 2.5 | 144 | 4 |

Within Lab precision (Total) was determined using three lots of reagents on three analyzers on each platform, using a single calibration per lot.

^{*} Control Fluid

^{**} Whole Blood

INSTRUCTIONS FOR USE

Performance Characteristics

HbA1c

Direct % Glycated Hemoglobin

VITROS 5600 Integrated System – Combined (mmol/mol, SI Units)

| Mean mmol/mol | Repeatability | | Between Run | | Between Day | | Between Lot | | Between Analyzer | | Total | | N | Days |
|---------------------------------|---------------|-----|-------------|-----|-------------|-----|-------------|-----|------------------|-----|-------|-----|-----|------|
| | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | SD | %CV | | |
| VITROS PV1 [*] 38.4 | 0.39 | 1.0 | 0.14 | 0.4 | 0.22 | 0.6 | 0.34 | 0.9 | 0.45 | 1.2 | 0.73 | 1.9 | 756 | 21 |
| Control [*] 30.9 | 0.37 | 1.2 | 0.12 | 0.4 | 0.19 | 0.6 | 0.63 | 2.0 | 0.00 | 0.0 | 0.77 | 2.5 | 792 | 22 |
| Control [*] 46.1 | 0.62 | 1.3 | 0.32 | 0.7 | 0.37 | 0.8 | 0.36 | 0.8 | 0.45 | 1.0 | 0.97 | 2.1 | 756 | 21 |
| Control [*] 65.0 | 0.79 | 1.2 | 0.19 | 0.3 | 0.33 | 0.5 | 0.59 | 0.9 | 0.40 | 0.6 | 1.13 | 1.7 | 792 | 22 |
| Control [*] 106.4 | 1.95 | 1.8 | 0.60 | 0.6 | 0.73 | 0.7 | 0.20 | 0.2 | 0.72 | 0.7 | 2.28 | 2.1 | 720 | 20 |
| Patient ^{**} 31.5 | 0.30 | 1.0 | 0.14 | 0.4 | 0.23 | 0.7 | 0.35 | 1.1 | 0.22 | 0.7 | 0.58 | 1.8 | 144 | 4 |
| Patient ^{**} 47.8 | 0.46 | 1.0 | 0.28 | 0.6 | 0.21 | 0.4 | 0.42 | 0.9 | 0.35 | 0.7 | 0.80 | 1.7 | 144 | 4 |
| Patient ^{**} 65.4 | 0.68 | 1.0 | 0.68 | 1.0 | 0.00 | 0.0 | 0.48 | 0.7 | 0.35 | 0.5 | 1.14 | 1.7 | 144 | 4 |
| Patient ^{**} 108.1 | 1.89 | 1.7 | 0.97 | 0.9 | 0.96 | 0.9 | 0.70 | 0.6 | 0.77 | 0.7 | 2.56 | 2.4 | 144 | 4 |

Within Lab precision (Total) was determined using three lots of reagents on three analyzers on each platform, using a single calibration per lot.

^{*} Control Fluid

^{**} Whole Blood

Specificity

Substances that Do Not Interfere

- Samples containing hemoglobin variants HbS up to 41%, HbC up to 38%, HbD up to 38% and HbE up to 26% of the total hemoglobin concentration do not interfere.
- The anti-HbA1c antibodies used in this kit do not cross-react with HbA0 up to 90%, HbA1a up to 1.5%, HbA1b up to 4% of the total hemoglobin concentration.
- This method is unaffected by the presence of acetylated hemoglobin, carbamylated hemoglobin, and labile glycated hemoglobin when samples were treated with 50 mg/dL Acetaldehyde, 150 mg/dL Urea, and 1500 mg/dL Glucose for four hours at 37 °C, respectively.
- The substances listed in this table were tested with the VITROS Chemistry Products HbA1c assay at %A1c values of approximately 6.5% and 8.5% (48 to 69 HbA1c (mmol/mol)), using protocols based on CLSI Protocol EP7²⁶, and found not to interfere, bias ≤0.5 %A1c and ≤0.6 %A1c (≤5 and ≤7 HbA1c (mmol/mol)), at the concentration shown.

| Substance [*] | Concentration | |
|------------------------|---------------|----------------|
| Acetaminophen | 20 mg/dL | 1.32 mmol/L |
| Acetylsalicylic acid | 100 mg/dL | 5.55 mmol/L |
| Ampicillin | 100 mg/dL | 2.86 mmol/L |
| Ascorbic Acid | 80 mg/dL | 4.54 mmol/L |
| Bilirubin | 50 mg/dL | 0.86 mmol/L |
| Ca-dobesilate | 20 mg/dL | 0.48 mmol/L |
| Cefoxitin sodium | 250 mg/dL | 5.56 mmol/L |
| Cholesterol | 350 mg/dL | 9.1 mmol/L |
| Cyclosporin | 0.5 mg/dL | 4.16 µmol/L |
| Doxycycline hyclate | 5 mg/dL | 0.10 mmol/L |
| Glucose | 1000 mg/dL | 55.5 mmol/L |
| Glycated Albumin | 500 mg/dL | 0.07 mmol/L |
| Ibuprofen | 50 mg/dL | 2.42 mmol/L |
| Insulin | 592.8 µU/mL | 3926 pmol/L |
| Intralipid | 500 mg/dL | Not Applicable |
| Levodopa | 2 mg/dL | 0.10 mmol/L |

HbA1c

Direct % Glycated Hemoglobin

INSTRUCTIONS FOR USE

Performance Characteristics

| Substance* | Concentration | |
|-----------------------|---------------|-------------|
| Metformin | 4 mg/dL | 0.31 mmol/L |
| Methyl dopa | 2 mg/dL | 94.7 µmol/L |
| Metronidazole | 20 mg/dL | 1.17 mmol/L |
| N-Acetylcysteine | 166.3 mg/dL | 10.2 mmol/L |
| Phenylbutazone | 40 mg/dL | 1.30 mmol/L |
| Rheumatoid Factor | 750 IU/mL | 750 kIU/L |
| Rifampicin | 6 mg/dL | 72.9 µmol/L |
| Rosiglitazone maleate | 0.8 mg/dL | 16.9 µmol/L |
| Theophylline | 10 mg/dL | 0.56 mmol/L |
| Total Protein | 5 g/dL | 50 g/L |
| Total Protein | 9 g/dL | 90 g/L |
| Triglyceride | 1000 mg/dL | 11.3 mmol/L |

* These results are intended to be representative; however, your results may differ somewhat due to reagent lot and specimen-related differences. The degree of interference at concentrations other than those listed for the substance and analyte may not be predictable.

VITROS 5,1 FS Analyzer

| Hemoglobin Variant | # samples tested | Variant Concentration Range (%) | Range of %A1c Concentration | Relative % Difference from Reference Concentration at Low and High HbA1c Concentrations | | | |
|--------------------|------------------|---------------------------------|-----------------------------|---|-------------|----------------|--------------|
| | | | | ~6.0 %HbA1c | | ~9.0 %HbA1c | |
| | | | | Relative %Bias | Range %Bias | Relative %Bias | Range %Bias |
| HbA0 | 70 | 49.3–90.4 | 5.8–13.4 | -1.11 | -6.15–2.91 | -0.90 | -1.93–2.04 |
| HbA1a | | 0.1–1.5 | | | | | |
| HbA1b | | 1.3–4.1 | | | | | |
| HbA2 | 22* | 4.9–6.1 | 5.7–9.0 | -1.80 | -7.25–0.68 | -2.67 | -5.08– -2.13 |
| HbC | 31** | 24.5–38.4 | 5.1–9.8 | -2.56 | -7.13–3.24 | -2.85 | -4.00–0.30 |
| HbD | 21 | 29.0–38.0 | 5.2–11.3 | -1.52 | -7.65–2.86 | -2.95 | -4.34–0.05 |
| HbE | 30*** | 14.3–26.3 | 5.4–9.1 | -0.37 | -6.45–3.95 | -1.00 | -1.26–6.01 |
| HbS | 40 | 28.2–41.6 | 4.6–12.7 | 0.76 | -4.72–5.95 | 0.70 | -1.60–7.86 |
| HbF | 43 | 0.2–34.8 | 5.5–12.8 | -5.31 | -29.17–2.65 | -2.58 | -24.47–1.85 |

* All twenty-two HbA2 patient samples were spiked.

** Nineteen native HbC patient samples and 12 spiked patient samples.

*** Eighteen native HbE patient samples and 12 spiked patient samples.

VITROS 4600 Analyzer

| Hemoglobin Variant | # samples tested | Variant Concentration Range (%) | Range of %A1c Concentration | Relative % Difference from Reference Concentration at Low and High HbA1c Concentrations | | | |
|--------------------|------------------|---------------------------------|-----------------------------|---|-------------|----------------|--------------|
| | | | | ~6.0 %HbA1c | | ~9.0 %HbA1c | |
| | | | | Relative %Bias | Range %Bias | Relative %Bias | Range %Bias |
| HbA0 | 70 | 49.3–90.4 | 5.8–13.4 | -0.91 | -6.91–4.53 | -0.58 | -0.52–3.07 |
| HbA1a | | 0.1–1.5 | | | | | |
| HbA1b | | 1.3–4.1 | | | | | |
| HbA2 | 22* | 4.9–6.1 | 5.7–9.0 | -1.78 | -7.07–1.60 | -1.81 | -4.52– -1.40 |
| HbC | 31** | 24.5–38.4 | 5.1–9.8 | -1.92 | -6.44–2.68 | -2.67 | -2.58– -0.26 |
| HbD | 21 | 29.0–38.0 | 5.2–11.3 | -1.33 | -6.09–1.25 | -3.74 | -4.61– -2.58 |
| HbE | 30*** | 14.3–26.3 | 5.4–9.1 | 2.38 | -3.41–7.92 | -0.69 | -5.13–1.46 |
| HbS | 40 | 28.2–41.6 | 4.6–12.7 | -0.84 | -4.07–4.29 | 0.30 | -1.05–7.85 |
| HbF | 43 | 0.2–34.8 | 5.5–12.8 | -4.88 | -28.65–4.29 | -2.63 | -23.35–2.14 |

* All twenty-two HbA2 patient samples were spiked.

** Nineteen native HbC patient samples and 12 spiked patient samples.

*** Eighteen native HbE patient samples and 12 spiked patient samples.

INSTRUCTIONS FOR USE

HbA1c

References

Direct % Glycated Hemoglobin

VITROS 5600 Analyzer

| Hemoglobin Variant | # samples tested | Variant Concentration Range (%) | Range of %A1c Concentration | Relative % Difference from Reference Concentration at Low and High HbA1c Concentrations | | | |
|--------------------|------------------|---------------------------------|-----------------------------|---|-------------|----------------|--------------|
| | | | | ~6.0 %HbA1c | | ~9.0 %HbA1c | |
| | | | | Relative %Bias | Range %Bias | Relative %Bias | Range %Bias |
| HbA0 | 70 | 49.3–90.4 | 5.8–13.4 | 0.06 | -5.83–4.53 | 0.03 | 0.03–2.79 |
| HbA1a | | 0.1–1.5 | | | | | |
| HbA1b | | 1.3–4.1 | | | | | |
| HbA2 | 22* | 4.9–6.1 | 5.7–9.0 | -3.25 | -8.12–1.33 | -4.12 | -8.03– -1.66 |
| HbC | 31** | 24.5–38.4 | 5.1–9.8 | -2.86 | -6.46–1.18 | -3.52 | -4.05– -1.26 |
| HbD | 21 | 29.0–38.0 | 5.2–11.3 | -2.71 | -7.18–0.45 | -5.37 | -6.65– -4.98 |
| HbE | 30*** | 14.3–26.3 | 5.4–9.1 | -0.49 | -6.15–3.94 | -2.39 | -7.11– -0.62 |
| HbS | 40 | 28.2–41.6 | 4.6–12.7 | -0.74 | -7.40–4.24 | -1.87 | -1.83–4.21 |
| HbF | 43 | 0.2–34.8 | 5.5–12.8 | -4.30 | -30.03–4.24 | -2.03 | -25.67–2.06 |

* All twenty-two HbA2 patient samples were spiked.

** Nineteen native HbC patient samples and 12 spiked patient samples.

*** Eighteen native HbE patient samples and 12 spiked patient samples.

References

1. A. Finke, U. Kobold, W. Hoelzel, C. Weykamp, K. Miedema, J-O. Jeppsson. Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood. *Clin Chem Lab Med*, 40(1): 78-89, 2002.
2. Goldstein DE, Little RR, Wiedmeyer HM, England JD, Rohlfing CL, Wilke AL. Is glycohemoglobin testing useful in diabetes mellitus: Lessons from the diabetes control and complications trial. *Clin Chem*, 40:1637-1640, 1994.
3. The Diabetes Control and Complications Trial Research Group: The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med* 329:977-86; 1993.
4. Intensive blood glucose control with sulphonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33). UK Prospective Diabetes Study (UKPDS) Group. *Lancet* 325:837-53; 1998.
5. Burtis C, Ashwood ER (ed). *Tietz Textbook of Clinical Chemistry*, 3rd ed. Philadelphia, PA; WB Saunders Co; 791, 1999.
6. Hanas R, John G., on behalf of the International HbA1c Consensus Committee. 2010 Consensus Statement on the Worldwide Standardization of the Hemoglobin A1c Measurement. *Clin. Chem.*, 56: 1362 – 1364, 2010.
7. CLSI. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline – Fourth Edition*. CLSI document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
8. Calam RR. Specimen Processing Separator Gels: An Update. *J Clin Immunoassay*. 11:86–90; 1988.
9. CLSI. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Sixth Edition*. CLSI document GP41-A6. Wayne, PA: Clinical and Laboratory Standard Institute; 2007.
10. CLSI. *Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard – Sixth Edition*. CLSI document GP42-A6. Wayne, PA: Clinical and Laboratory Standard Institute; 2008.
11. CLSI. *Evaluation of the Linearity of Quantitative Measurement Procedure: A Statistical Approach; Approved Guideline*. CLSI document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
12. Hoelzel W, Weykamp C, Jeppsson J, Miedema K, Barr JR, Goodall I, Hoshino T, John G, Kobold U, Little R, Mosca A, Mauri P, Paroni R, Susanto F, Takei I, Thienpont L, Umamoto M, Wiedmeyer H. IFCC Working Group on HbA1c Standardization, IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardization Schemes in the United States, Japan and Sweden: A Method Comparison Study, *Clin Chem*, 50 (1): 166-174, 2004.
13. Geistanger A, Arends S, Berding C, Hoshino T, Jeppsson J, Little R, Siebelder C, Weykamp C, on behalf of the IFCC Working Group on Standardization of Hemoglobin A1c. Statistical Methods for Monitoring the Relationship between the IFCC Reference Measurement Procedure for Hemoglobin A1c and the Designated Comparison Methods in the United States, Japan, and Sweden. *Clin Chem*, 54: 1379-1385, 2008.
14. CLSI. *Statistical Quality Control for Quantitative Measurements: Principles and Definitions; Approved Guideline – Third Edition*. NCCLS document C24-A3 [ISBN 1-56238-613-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA; 2006.
15. <http://diabetes.niddk.nih.gov/dm/pubs/hemovari-A1C/index.aspx#hemoglobinopathies>. Accessed on April 28, 2014.
16. <http://www.ngsp.org/factors.asp>. Accessed on April 28, 2014.
17. Bry L, Chen PC, Sacks DB. Effects of hemoglobin variants and chemically modified derivatives on assays for glycated hemoglobin. *Clin Chem*. 47: 153-163; 2001.

HbA1c

Direct % Glycated Hemoglobin








































INSTRUCTIONS FOR USE

Glossary of Symbols

18. International Association of Diabetes and Pregnancy Study Groups Recommendations on the Diagnosis and Classification of Hyperglycemia in Pregnancy. *Diabetes Care*. 33:No 3; 2010.
19. International Expert Committee Report on the Role of the A1c Assay in the Diagnosis of Diabetes. *Diabetes Care*. 32: (7):1327-1334; 2009.
20. Standards of Medical Care in Diabetes. *Diabetes Care*. 38:Supplement 1; S8-S16, 2015.
21. http://www.who.int/diabetes/publications/report-hba1c_2011.pdf. Accessed on April 28, 2014.
22. Implementation of haemoglobin A1c results traceable to the IFCC reference system: the way forward; *Clin Chem Lab Med* 2007; 45(8):942-944.
23. CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition*. CLSI document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
24. CLSI. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline – Second Edition (Interim Revision)*. CLSI document EP9-A2-IR. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
25. CLSI. *Evaluation of Precision Performance of Clinical Chemistry Devices; Approved Guideline*. CLSI document EP05-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2004.
26. CLSI. *Interference Testing In Clinical Chemistry; Approved Guideline – Second Edition*. CLSI document EP7-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.

Glossary of Symbols

The following symbols may have been used in the labeling of this product.

| | | | | | |
|---|---|---|--|---|-----------------------------------|
|  | Do Not Reuse |  | Upper Limit of Temperature |  | Range |
|  | Use by or Expiration Date (Year-Month-Day) |  | Lower Limit of Temperature |  | Range of Means |
|  | Batch Code or Lot Number |  | Temperature Limitation |  | Midpoint |
|  | Serial Number |  | Consult Instructions for Use |  | Revised |
|  | Catalog Number or Product Code |  | Attention: The Instructions for Use (IFU) has been updated |  | Supersedes |
|  | Caution |  | For use in Slide Supply 1 |  | Irritant |
|  | Manufacturer |  | For use in Slide Supply 2 |  | Harmful |
|  | Date of Manufacture |  | SI Units |  | Toxic |
|  | Authorized Representative in the European Community |  | Conventional Units |  | Corrosive |
|  | Contains Sufficient for "n" Tests |  | Value |  | Flammable |
|  | In vitro Diagnostic Medical Device |  | Der Grüne Punkt (the Green Dot). Manufacturer follows certain packaging material waste disposal management regulations |  | Estimated within-lab SD |
|  | Corrosive |  | Flammable |  | Serious Health Hazards |
|  | Health Hazards |  | Acute Toxicity |  | Environmental or Aquatic Toxicity |

INSTRUCTIONS FOR USE

Revision History

HbA1c

Direct % Glycated Hemoglobin

Revision History

| Date of Revision | Version | Description of Technical Changes* |
|------------------|---------|--|
| 2015-09-11 | 3.0 | Substances that Do Not Interfere: Triglyceride changed from 500 mg/dL (5.7 mmol/L) to 1000 mg/dL (11.3 mmol/L) |
| 2015-06-04 | 2.0 | <ul style="list-style-type: none"> • First release of document in US • Prescription Use Statement added • Black Box Warning added • Intended Use: updated wording • Warnings and Precautions: updated to align with the new Safety Data Sheets • Specimens Recommended: added K₂EDTA and K₃EDTA • Calibration: added Linearity section • Results: <ul style="list-style-type: none"> – Reporting Units updated – Unit Conversion updated • Other Limitations: added statements • Expected Values: Reference Interval section updated to include new table and wording • Method Comparison: <ul style="list-style-type: none"> – updated wording – updated plot diagrams – updated table • Precision: <ul style="list-style-type: none"> – updated wording – added new tables • Substances that Do Not Interfere: <ul style="list-style-type: none"> – changed HbD variant concentration level from 37% to 38% – new tables added • References: <ul style="list-style-type: none"> – updated references 7, 19, 20 – added reference 11, 22 • Glossary of Symbols: added Globally Harmonized Symbols to comply with the Classification, Labelling and Packaging (CLP) Regulations |
| 2014-05-19 | 1.0 | First release of document |

* The change bars indicate the position of a technical amendment to the text with respect to the previous version of the document.

When this Instructions For Use is replaced, sign and date below and retain as specified by local regulations or laboratory policies, as appropriate.

Signature


Obsolete Date

HbA1c

Direct % Glycated Hemoglobin

INSTRUCTIONS FOR USE

Revision History



Ortho-Clinical Diagnostics
50-100 Holmers Farm Way
High Wycombe
Buckinghamshire
HP12 4DP
United Kingdom



Ortho-Clinical Diagnostics, Inc.
100 Indigo Creek Drive
Rochester, NY 14626

VITROS is a registered trademark of Ortho-Clinical Diagnostics, Inc.

The third party trademarks used herein are trademarks of their respective owners.

© Ortho-Clinical Diagnostics, Inc., 2014-2015

Ortho Clinical Diagnostics