AQT90 FLEX

Reagent Pack

Content

en	Reagent Pack	4
bg	Касета за реактиви	5
cs	Reagenční souprava	6
da	Reagenspakke	7
de	Reagenzien-Kassette	8
el	Πακέτο αντιδραστηρίων	9
es	Paquete de Reactivos	10
et	Reaktiivipakk	11
fi	Reagenssipakkaus	12
fr	Pack de réactif	13
hu	Reagenscsomag	14
it	Reagent pack	15
lt	Reagentų pakuotė	16
nl	Reagenskit	17
no	Reagens Pakke	18
pl	Pakiet odczynników	19
pt	Cartucho de Reagentes	20
ro	Pachet de reactiv	21
ru	Блок реагентов	22
sv	Reagent Pack	23
tr	Reaktif Paketi	24

Reagent Pack

Proprietary name and common name: Reagent Pack (HEPES buffer)

Catalog no. 903-006

Intended use

For *in vitro* diagnostic use.

The assay buffer in the Reagent Pack is used for the operation of the AQT90 FLEX analyzer. The pack also receives the discarded cups and liquid waste.

The Reagent Pack is used in conjunction with any of the Test Kits and contains sufficient buffer for up to 200 tests.

NOTICE: A number of activities, for example needle clean and startup will reduce the buffer and thus the number of tests.

Summary

The Reagent Pack is used on the AQT90 FLEX analyzer.

See the *AQT90 FLEX operator's manual* for details.

Principles of the procedure

See the *AQT90 FLEX operator's manual*.

Contents

Component:	Reagent Pack
Quantity:	1 pack
Shelf life and storage:	Store at 2-32 °C until the expiration date stated on the label. NOTICE: During transport, the Reagent Pack can be kept at up to 40 °C for a maximum of 2 weeks.
On-board stability:	Up to 25 days

The Reagent Pack contains assay buffer and closed receptacles to contain discarded cups and liquid waste. The same Reagent Pack is used for all assays.

Chemical composition:

Main buffer component:	Hepes
Toxic ingredients:	Germal II (diazolidinylurea - anti-bacterial agent): 1 g/L Propyl p-hydroxybenzoate: 0.2 g/L
pH:	7.5

**WARNING – Risk of infection**

Handle and discard used Reagent Packs as biohazardous waste. They contain the processed cups and the discarded buffer and blood fluids and they should be treated as potentially infectious materials. Failure to follow appropriate caution in handling and disposing these may result in exposure to potentially infectious materials.

Procedure

See the *AQT90 FLEX operator's manual* for details on handling procedures, assay description and a detailed procedure.

Briefly, the procedure consists of the following steps:

1. To remove the used/expired Reagent Pack from the analyzer:
 - a. Press **Analyzer status** on the main screen.
 - b. Press **Consumables > Insert/Replace > Replace Reagent Pack** on the **Analyzer status** screen.
 - c. Pull out the Reagent Pack that has been ejected.
 - d. Dispose the used Reagent Pack as biohazardous waste.
2. Remove the top label on the new Reagent Pack to expose the biohazard label.
3. Insert the pack in its inlet on the analyzer.
4. Press **Accept** to confirm.

Explanations of graphical symbols:

Manufactured by

Reagent Pack

Product name

Contents

Lot no.

Catalog no.

Biohazard

Expiration date

For *in vitro* diagnostic use

Mark of compliance with the applicable EU Directives

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Radiometer, the Radiometer logo, ABL, AQT, TCM, RADIANCE, PICO and CLINITUBES are trademarks of Radiometer Medical ApS.
© 2009 Radiometer Medical ApS. All rights reserved. 990-879. 200909L.