HCG

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000 HCG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of human chorionic gonadotropin (HCG) in serum, and for strictly qualitative determinations in urine, as an aid in the detection of pregnancy.

Catalog Number: **L2KCG2** (200 tests), **L2KCG6** (600 tests)

Test Code: **HCG** Color: **Orange**

Summary and Explanation

Human chorionic gonadotropin (HCG) is a two-chain glycoprotein hormone (MW ~37,000) normally found in blood and urine only during pregnancy. It is secreted by placental tissue, beginning with the primitive trophoblast, almost from the time of implantation, and serves to support the corpus luteum during the early weeks of pregnancy.

According to the literature, circulating HCG typically reaches levels of approximately 2000 mIU/mL one month after conception. A peak level on the order of 100,000 mIU/mL is attained in the third month, after which a gradual decline sets in. Following delivery, the HCG level normally undergoes rapid descent, reaching nonpregnant concentrations (usually less than 5 mIU/mL) some two weeks later.

Ectopic pregnancies and pregnancies terminating in spontaneous abortion tend to have lower than normal circulating HCG levels, while somewhat higher levels are often seen in multiple pregnancies.

Principle of the Procedure

IMMULITE 2000 HCG is a solid-phase, two-site chemiluminescent immunometric assay.

Incubation Cycles: 1 × 30 minutes

Time to First Result: 35 minutes

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 HCG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 5 µL serum or urine

Urine: A freshly voided, first-morning specimen is preferred.²² Centrifuging or filtering (Whatman #1) is recommended to clear cloudy samples.

Storage

Serum: 7 days at 2–8°C, or 2 months at –20°C.¹⁸

Urine: 48 hours at 2–8°C,²² or at –20°C until assayed.¹³

Warnings and Precautions

For *in vitro* diagnostic use.

**CAUTION! POTENTIAL BIOHAZARD**

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.²³⁻²⁵

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations. Contains: sodium azide; HCG Adjustors

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

HCG Bead Pack (L2CG12)

With barcode. 200 beads, coated with monoclonal murine anti-HCG. Stable at 2–8°C until expiration date.

L2KCG2: 1 pack

L2KCG6: 3 packs

HCG Reagent Wedge (L2CGA2)

With barcode. 11.5 mL of alkaline phosphatase (bovine calf intestine) conjugated to polyclonal ovine anti-HCG, in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KCG2: 1 wedge

L2KCG6: 3 wedges

Before use, tear off the top of the label at the perforation, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

HCG Adjustors (LCGL, LCGH)

Two vials (Low and High) 2 mL each, of HCG in an HCG-free human serum matrix, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KCG2: 1 set

L2KCG6: 2 sets

Before running adjustors, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

Kit Components Supplied Separately**HCG Sample Diluent (L2CGZ, L2CGZ4)**

For the on-board dilution of high samples. One vial of concentrated (ready-to-use) HCG-free human serum, with barcode labels. Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C.

Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2CGZ, L2CGZ4: 5 labels

L2CGZ: 50 mL **L2CGZ4:** 100 mL

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes
(16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LCGCM: HCG High Level Control Module

Also Required

Distilled or deionized water; test tubes;
controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
4 weeks

Quality Control Samples: Follow government regulations or accreditation requirements for quality control frequency.

Use controls or serum pools with at least two levels (low and high) of HCG.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

Qualitative Procedure for Urine

The qualitative procedure allows for determining whether a urine sample has an HCG concentration above or below 30 mIU/mL, within the precision of the assay. Concentrations above this cutoff may be considered positive for pregnancy providing other causes of high HCG levels have been ruled out, viz. trophoblastic and nontrophoblastic neoplasia.^{4,9,10,15}

Concentrations below this cutoff are to be expected for samples from individuals who are not pregnant, as well as for women in such an early stage of pregnancy that the

HCG concentration has not yet reached the decision level set for the procedure. Assaying another sample collected two or more days later should settle the matter, since in early pregnancy HCG levels can be expected to double in approximately this time.¹¹

Interpretation of Results

A result of greater than or equal to 30 mIU/mL indicates that the patient sample is **positive** for pregnancy.

A result of less than 30 mIU/mL indicates that the patient sample is **negative** for pregnancy.

Expected Values

Based on its relationship to IMMULITE HCG (see Method Comparison), the assay can be expected to have essentially the same reference ranges.

Males and Nonpregnant Females

In a study performed with the IMMULITE HCG assay and involving serum samples from 797 adult volunteers, including 428 males and 369 nonpregnant females, 95% of the results for males were below that assay's detection limit (1.1 mIU/mL) and all were below 2.5 mIU/mL, while 83% of the results for nonpregnant females were below the detection limit, 95% were below 2.7 mIU/mL, and all were below 5.3 mIU/mL. The subjects were in apparent good health based on a questionnaire.

Pregnant Females

A total of 593 serum samples from apparently healthy pregnant women were processed by the IMMULITE HCG assay. The results are summarized below (in mIU/mL) by gestational age (Gest.) — equivalently, by weeks since last menstrual period (LMP). Also see "Expected Values During Pregnancy" graph.

Weeks of Pregnancy

Gest.	LMP	n	Median	Central 95%
1.3–2	3.3–4	30	71	16–156
2–3	4–5	54	607	101–4870
3–4	5–6	34	5243	1110–31,500
4–5	6–7	34	26,983	2560–82,300
5–6	7–8	36	52,090	23,100–151,000
6–7	8–9	33	93,598	27,300–233,000
7–11	9–13	116	117,678	20,900–291,000
11–16	13–18	72	40,989	6140–103,000
16–21	18–23	80	20,868	4720–80,100
21–39	23–41	104	15,352	2700–78,100

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Pathological Conditions

Ectopic pregnancies or abnormal intrauterine pregnancies have circulating HCG levels below normal for that gestational period; after an initial rise these levels plateau or fall.^{3,17,18}

Limitations

The procedure for urine samples is designed solely as a screen for pregnancy, using a 30 mIU/mL cutoff. Hence urinary results should be reported in a strictly qualitative manner, e.g. as “positive” or “negative” for pregnancy — concentration values should *not* be reported for this sample type.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem* 1988;34: 27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See tables and graphs for data *representative* of the assay's performance. Results are expressed as mIU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Reportable Range: Up to 5000 mIU/mL, in terms of the WHO 3rd IS 75/537 for quantitative determinations in serum.

Analytical Sensitivity: 0.4 mIU/mL

High-Dose Hook Effect: None up to 600,000 mIU/mL

Precision: Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See “Precision” table.)

Linearity: Samples were assayed under various dilutions. (See “Linearity” table for representative data.)

Recovery: Samples spiked 1 to 19 with three HCG solutions (1313, 5720 and 30,340 mIU/mL) were assayed. (See “Recovery” table for representative data.)

Specificity: The assay is highly specific for HCG, with low crossreactivity to other glycoprotein hormones present in patient samples. (See “Specificity” table.)

For concentrations as high as 10,000 ng/mL, the HCG α -subunit (WHO 75/569) showed less than 1% crossreactivity. HCG β -subunit (WHO 75/551), spiked into HCG Sample Diluent, yielded the tabulated results. (See “Specificity” table.)

Bilirubin: Presence of bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 384 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 14 volunteers into plain, heparinized, EDTA vacutainer

tubes, and Becton Dickinson SST[®]. Equal volumes of the matched samples were spiked with various concentrations of HCG, to obtain values throughout the calibration range of the assay, and then assayed by the IMMULITE 2000 HCG procedure.

(Heparin) = 1.06 (Serum) – 7.06 mIU/mL
r = 0.99

(EDTA) = 1.08 (Serum) – 6.7 mIU/mL
r = 0.99

(SST) = 0.98 (PlainTubes) – 3.25 mIU/mL
r = 0.99

Means:

204 mIU/mL (Serum)
211 mIU/mL (Heparin)
197 mIU/mL (EDTA)
214 mIU/mL (SST)

Method Comparison (Serum): The assay was compared to IMMULITE HCG on 115 patient samples. (Concentration range: approximately 1 to over 3000 mIU/mL. See graph.)
By linear regression:

(IML 2000) = 1.0 (IML) – 21 mIU/mL
r = 0.994

Means:

668 mIU/mL (IMMULITE 2000)
689 mIU/mL (IMMULITE)

For the 63 data pairs (out of 116) with results by both methods from 1 to 500 mIU/mL:

(IML 2000) = 0.91 (IML) + 0.9 mIU/mL
r = 0.987

For the 22 data pairs (out of 116) with results by both methods from 1 to 50 mIU/mL:

(IML 2000) = 0.96 (IML) + 0.3 mIU/mL
r = 0.953

Method Comparison (Urine): A total of 123 urine samples from pregnant women were processed by the IMMULITE 2000 HCG assay and by Double Antibody HCG. The samples, some of which were diluted, had HCG values ranging up to approximately 300 mIU/mL. By linear regression:

(IML 2000) = 1.04 (DAb) – 10 mIU/mL
r = 0.988

In another study with the IMMULITE 2000 HCG assay, as a check on the appropriateness of the 30 mIU/mL cutoff used in the qualitative procedure, urinary results for nonpregnant women showed a clean separation from the distribution of

repeated measurements made on a 30 mIU/mL sample. (All results for the nonpregnant women: less than 3 mIU/mL. Mean – 6SD for the 30 mIU/mL sample: 23 mIU/mL.)

References

- 1) Bangham DR, Storrer PL. Standardisation of human chorionic gonadotropin, HCG subunits, and pregnancy tests. *Lancet* 1982;1:390.
- 2) Batzer F. Hormonal evaluation of early pregnancy. *Fertil Steril* 1980;34:1–12.
- 3) Braunstein GD, et al. First-trimester chorionic gonadotropin measurements as an aid in the diagnosis of early pregnancy disorders. *Amer J Obstet Gynecol* 1978;131:25–32.
- 4) Clayton L, et al. Endocrine aspects of trophoblastic neoplasia. *J Reprod Med* 1981;26:192–9.
- 5) Corson SL, Horwitz CA, Lau HL, Soderstrom R. Early diagnosis of pregnancy: a symposium. *J Reprod Med* 1981;26:149–78.
- 6) DeCherney AH, Minkin MJ, Spangler S. Contemporary management of ectopic pregnancy. *J Reprod Med* 1981;26:519–23.
- 7) Goldstein DP, Berkowitz RS, Kosasa TS. Measurement of human chorionic gonadotropin. Major Problems in *Obstet Gynecol* 1982;14:64–97.
- 8) Hager H, Van Weemen BK. Standardisation of human chorionic gonadotropin, HCG subunits, and pregnancy tests. *Lancet* 1982;1:629.
- 9) Kohorn EI. Hydatidiform mole and gestational trophoblastic disease in southern Connecticut. *Obstet and Gynecol* 1982;59:78–84. *Erratum:* The first diagram belongs with the caption to Figure 3, and vice versa.
- 10) Kosasa TS. Measurement of human chorionic gonadotropin. *J Reprod Med* 1981;26:201–6.
- 11) Lenton E, Neal L, Sulaiman R. Plasma concentrations of human gonadotropin from the time of implantation until the second week of pregnancy. *Fertil Steril* 1982;37:773–8.
- 12) Lindstedt G, Janson PO, Thorburn J. Sensitivity of serum chorionic gonadotropin assay for ectopic pregnancy. *Lancet* 1981;1:781–2.
- 13) McCready J, Braunstein GD, et al. Modification of the choriogonadotropin beta-subunit radioimmunoassay for determination of urinary choriogonadotropin. *Clin Chem* 1978;24:1958–61.
- 14) Storrer PL, Gaines-Das R, Bangham DR. International Reference Preparation of Human Chorionic Gonadotrophin for Immunoassay: potency estimates in various bioassay and protein binding assay systems. *J Endocrinol* 1980;84:295–310.
- 15) Vaitukaitis J. Human chorionic gonadotropin – a hormone secreted for many reasons. *N Eng J Med* 1979;301:324–6.
- 16) Danzer H, Braunstein GD, et al. Maternal serum human chorionic gonadotropin concentrations and fetal sex predictions. *Fertil Steril* 1980;34:336–40.
- 17) Hussa RO. *The Clinical Marker HCG*. New York: Praeger Publishers, 1987.
- 18) Greene MF, de Montserrat MF, Tulchinsky D. In: Tietz NW,

editor. Textbook of Clinical Chemistry. Philadelphia: Saunders, 1986: 1745.
 19) Braunstein GD, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. Am J Obstet Gynecol 1976;126:678.
 20) Babson AL The IMMULITE Automated Immunoassay System. J Clin Immunoassay 1991;14:83-8.
 21) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.
 22) Tietz NW, ed. Clinical Guide to Laboratory Tests. 2nd Ed. Philadelphia: W.B. Saunders, 1990.
 23) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8.
 24) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
 25) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

Technical Assistance

For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

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Linearity (mIU/mL)

	Dilution ¹	Observed ²	Expected ³	%O/E ⁴
1	16 in 16 ⁵	204	—	—
	8 in 16	100	102	98%
	4 in 16	50	51	98%
	2 in 16	24	26	92%
	1 in 16	12	13	92%
2	16 in 16	552	—	—
	8 in 16	273	276	99%
	4 in 16	134	138	97%
	2 in 16	68	69	99%
	1 in 16	33	35	94%
3	16 in 16	1144	—	—
	8 in 16	556	572	97%
	4 in 16	267	286	93%
	2 in 16	140	143	98%
	1 in 16	70	72	97%
4	16 in 16	2088	—	—
	8 in 16	1020	1044	98%
	4 in 16	503	522	96%
	2 in 16	240	261	92%
	1 in 16	119	131	91%

Dansk. ¹Fortynding, ²Observeret, ³Forventet, ⁴%Obs./ Forv., ⁵Ufortyndet. **Eesti.** ¹Lahjendus, ²Täheldatud (T), ³Oodatud (O), ⁴% T/O, ⁵16 16-st. **Latviski.** ¹Atšķaidījums, ²Novērots (N), ³Gaidīts (G), ⁴% N/G, ⁵16 : 16. **Lietuviškai.** ¹Skiedimas, ²Nustatyta (N), ³Tikėtasi (T), ⁴% N/T, ⁵16 iš 16. **Norsk.** ¹Fortynning, ²Observeret (O), ³Forventet (E), ⁴% O/E, ⁵Ufortynnet. **Svenska.** ¹Spädning, ²Observerat (O), ³Förväntat (E), ⁴% O/E, ⁵Ospätt.

Recovery (mIU/mL)

	Solution ¹	Observed ²	Expected ³	%O/E ⁴
1	—	270	—	—
	A	339	322	105%
	B	597	543	110%
	C	1770	1774	100%
2	—	554	—	—
	A	617	592	104%
	B	913	812	112%
	C	2114	2043	103%
3	—	949	—	—
	A	984	967	102%
	B	1290	1188	109%
	C	2610	2419	108%

Dansk. ¹Opløsning, ²Observeret, ³Forventet, ⁴%Obs./ Forv. **Eesti.** ¹Lahus, ²Täheledatud (T), ³Oodatud (O), ⁴% T/O. **Latviski.** ¹Šķīdums, ²Novērots (N), ³Gaidīts (G), ⁴% N/G. **Lietuviškai.** ¹Tirpalas, ²Nustatyta (N), ³Tikėtasi (T), ⁴% N/T. **Norsk.** ¹Løsning, ²Observert (O), ³Forventet (E), ⁴% O/E. **Svenska.** ¹Lösning, ²Observerat (O), ³Förväntat (E), ⁴% O/E.

Specificity

	Amount Added (ng/mL) ¹	Apparent Concentration (ng/mL) ²	% Cross-reactivity ³
FSH	26.8	ND	ND
LH	16.5	ND	ND
TSH	860	ND	ND

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

IMMULITE is a trademark of Siemens Healthcare Diagnostics.

© 2015 Siemens Healthcare Diagnostics. All rights reserved.

Made in: UK

Siemens Healthcare
Diagnostics Products Ltd.
Glyn Rhonwy, Llanberis,
Gwynedd LL55 4EL
United Kingdom

0088

when used for Trisomy 21 testing
ved brug til Trisomy 21 testning
kasutatuna trisoomia 21
analüüsimiseks
lietojot 21. hromosomas trisomijas
testçöanai
kai naudojama trisomijos 21 tyrimui
når brukt til Trisomy 21-testing
när använd för trisomi 21-tester

2015-05-21

PINL2KCG – 8 {33}

Changes in this Edition:

cc#EU22534: 1) Removed references to discontinued products. 2) Under Warnings and Precautions, added biohazard symbol and potential biohazard statement, caution statement, required GHS hazard and precautionary codes and statements, offending chemicals and sources. 3) Under References, added citations from the potential biohazard statement. 4) Changed "Origin: UK" to "Made in: UK". 5) Updated the manufacturer address to include Glyn Rhonwy. 6) In Understanding the Symbols, revised the symbol for catalog number, added a GHS pictogram for health hazard, and updated existing symbols with GHS pictograms and terminology for exclamation mark, corrosion, skull and crossbones, and environment.

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