IMMULITE® 2000 Multi-Diluent 2

English

Intended Use: For the on-board dilution of patient samples exceeding the calibration range of the indicated IMMULITE® 2000 assays. See the IMMULITE 2000 Operator's Manual for an explanation of on-board dilution features.

Catalog Number: **L2M2Z, L2M2Z4**

Use with IMMULITE 2000: AFP (L2KAP); Albumin (L2KHA); Beta-2 Microglobulin (L2KBM); BR-MA (L2KBR); Calcitonin (L2KCL); CEA (L2KCE); CK-MB (L2KMB); EPO (L2KEPN); Ferritin (L2KFE); Gastrin (L2KGA); GI-MA (L2KGI); Growth Hormone (L2KGRH); Total IgE (L2KIE); IL2R (L2KIP); IL-6 (L2K6P); Myoglobin (L2KMY); Nicotine Metabolite (L2KNM); PAPP-A (L2KPC); PSA (L2KPS); Free PSA (L2KPF).

Materials Supplied

One vial, concentrated (ready-to-use), nonhuman protein/buffer matrix, with preservative.

L2M2Z: 25 mL **L2M2Z4:** 55 mL

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2M2Z: 3 labels **L2M2Z4:** 5 labels

Storage: 30 days (after opening) at 2–8°C or 6 months (aliquotted) at –20°C. Dispose of in accordance with applicable laws.

Warnings and Precautions

For *in vitro* diagnostic use.

Multi-Diluent 2 is not intended for use in linearity or sensitivity studies, but rather for the dilution of patient samples that exceed the calibration range of the assay.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Technical Assistance

For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2015 Siemens Healthcare Diagnostics. All rights reserved.