# **Syphilis Screen**

**For use on IMMULITE® 2000 systems**

**SIEMENS**

# IMMULITE® 2000 Syphilis Screen

---

## English

---

**Intended Use:** The IMMULITE 2000 Syphilis Screen test is a treponemal testing procedure for the qualitative detection of antibodies to *Treponema pallidum* in human serum or heparinized plasma on the IMMULITE 2000 Systems Analyzers as an aid in the diagnosis of syphilis.

Catalog Numbers: **L2KSY2** (200 tests), **L2KSY6** (600 tests)

Test Code: **SYP** Color: **Dark Pink**

### Summary and Explanation

Syphilis is primarily transmitted via sexual contact, but can also be transmitted from mother to fetus. It is caused by the spirochete *T. pallidum* and has never been successfully cultured in artificial media. Syphilis infections are classified into early (infectious) and late (non-infectious) stages. Early syphilis may be further divided into primary, secondary, and early latent syphilis. The signs and symptoms of syphilis are numerous; before the advent of serological testing, precise diagnosis was very difficult. In fact, the disease was often confused with other diseases, particularly in its tertiary stage. If not treated, syphilis can cause serious effects such as damage to the heart, aorta, brain, eyes, and bones. In some cases these effects can be fatal. Therefore, the serological diagnosis of syphilis is very important.

The serologic diagnosis of syphilis is classified into two groups: nontreponemal tests and treponemal tests. Nontreponemal tests (venereal disease research laboratory [VDRL], rapid plasma reagin [RPR]) detect antibodies formed by the host in response to lipid material released from damaged host cells as well as to lipoprotein-like material released from the spirochete. Treponemal tests detect specific treponemal antibodies, and the techniques used include agglutination (*T. pallidum* haemagglutination [TPHA], *T. pallidum* particle agglutination [TPPA]), immunoassay (enzyme immunoassay

[EIA] or chemiluminescent immunoassay [CLIA], immunofluorescence (fluorescent treponemal antibody absorption [FTA-ABS]), and immunoblotting. Nontreponemal tests have poor sensitivity and specificity, and recombinant antigen-based treponemal tests have higher sensitivity and specificity than native antigen-based treponemal tests.<sup>1,2,3</sup> The IMMULITE 2000 Syphilis Screen assay uses a recombinant treponemal antigen and is a fully automated chemiluminescent immunoassay.

### Principle of the Procedure

IMMULITE 2000 Syphilis Screen assay is a solid phase, one-step chemiluminescent enzyme immunoassay. The solid phase (bead) is coated with purified recombinant *T. pallidum* 17 (Tp17) antigen. The liquid phase consists of alkaline phosphatase (bovine calf intestine) conjugated to purified recombinant *T. pallidum* 17 (Tp17) antigen.

Patient sample and the reagent are incubated together with the coated bead for 30 minutes. During this time, total antibody to *T. pallidum* in the sample forms the antigen sandwich complex with purified recombinant *T. pallidum* 17 (Tp17) antigen on the bead and enzyme conjugated purified recombinant *T. pallidum* 17 (Tp17) antigen in the reagent. Unbound patient sample and enzyme conjugate are then removed by centrifugal washes. Finally, chemiluminescent substrate is added to the reaction tube containing the bead and the signal is generated in proportion to the bound enzyme.

**Incubation Cycles:** 1 × 30 minutes.

**Time to first result:** 35 minutes.

### Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Samples that are cloudy or have particulate material should be clarified by low-speed centrifugation.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 Syphilis Screen has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

**Volume Required:** 100 µL serum or heparinized plasma.

**Storage:** 3 days at 2–8°C; for long-term storage, keep specimens at –20°C.

## Warnings and Precautions

For *in vitro* diagnostic use.

### CAUTION! POTENTIAL BIOHAZARD

Contains human source material. Each donation of human blood or blood component was tested by FDA-approved methods for the presence of antibodies to human immunodeficiency virus type 1 (HIV-1) and type 2 (HIV-2) as well as for hepatitis B surface antigen (HBsAg) and antibody to hepatitis C virus (HCV). The test results were negative (not repeatedly reactive). No test offers complete assurance that these or other infectious agents are absent; this material should be handled using good laboratory practices and universal precautions.<sup>4-6</sup>

**CAUTION:** This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

**H302 + H312,  
H412**

**P280, P273,  
P301 + P312,  
P302 + P312,  
P501**

**Warning!** Harmful if swallowed or in contact with skin. Harmful to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. IF ON SKIN: Call a POISON CENTER or doctor/physician if you feel unwell. Dispose of contents and container in accordance with all local, regional, and national regulations.

**Contains:** sodium azide; Syphilis Screen Adjustor, Syphilis Screen Controls

**Reagents:** Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

EDTA plasma should not be used as a sample type.

Syphilis Positive Control (LSYC2) contains antibodies against *T. pallidum*.

Syphilis Screen Adjustor (LSYR) contains antibodies against *T. pallidum*.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

**Chemiluminescent Substrate:** Avoid contamination and exposure to direct sunlight. (See insert.)

**Water:** Use distilled or deionized water.

## Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

### Syphilis Screen Bead Pack (L2SY12)

With barcode. 200 beads, coated with purified recombinant *T. pallidum* 17 (Tp17) antigen. Stable at 2–8°C until expiration date.

**L2KSY2:** 1 pack. **L2KSY6:** 3 packs.

### Syphilis Screen Reagent Wedge (L2SYA2)

With barcode. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to purified recombinant *T. pallidum* 17 (Tp17) antigen. Stable at 2–8°C until expiration date.

**L2KSY2:** 1 wedge. **L2KSY6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

### Syphilis Screen Adjustor (LSYR)

One vial of lyophilized human serum with antibodies reactive to *T. pallidum*, with preservative. Reconstitute with **4.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 28 days after reconstitution, or for 6 months (aliquotted) at –20°C.

**L2KSY2:** 1 vial. **L2KSY6:** 2 vials.

### Syphilis Screen Controls (LSYC1, LSYC2)

**LSYC1: (Negative Control):** One vial of lyophilized human serum nonreactive to *T. pallidum*, with preservative. **LSYC2: (Positive Control):** One vial of lyophilized human serum with antibodies reactive to *T. pallidum*, with preservative.

Reconstitute each vial with **6.0 mL** distilled or deionized water. Mix by gentle swirling or inversion until the lyophilized material is fully dissolved. Stable at 2–8°C for 28 days after reconstitution, or for 6 months (aliquotted) at –20°C.

**L2KSY2:** 1 set. **L2KSY6:** 2 sets.

Before running adjustors or controls, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

## Kit Components Supplied Separately

**L2SUBM:** Chemiluminescent Substrate

**L2PWSM:** Probe Wash

**L2KPM:** Probe Cleaning Kit

**LRXT:** Reaction Tubes (disposable)

**LSYCM:** Bi-level Syphilis Screen Control Module

Also Required:

Distilled or deionized water; test tubes.

## Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

### Recommended Adjustment Interval:

2 weeks.

**Quality Control:** To monitor system performance and chart trends, as a minimum requirement, quality control materials with at least two levels (low and high) of *T. pallidum* should be assayed on each day that samples are analyzed. Quality control samples should also be assayed when performing adjustment. Treat all quality control samples the same as patient samples.

Siemens Healthcare Diagnostics recommends the use of commercially available quality control materials with at least 2 levels (low and high). A satisfactory level of performance is achieved when the analyte values obtained are within the Acceptable Control Range for the system, or within an established range determined by an appropriate internal laboratory quality control scheme.

### Calculation of Cutoff and S/CO Ratio:

The Master Cutoff of the assay was determined from representative samples to achieve optimal sensitivity and specificity for the assay.

The cutoff is set equal to the average counts per second (mean cps) of the Adjustor (from the most recent adjustment) multiplied by Curve Parameter 1. (See the "Low Adjustor

CPS" and "Curve Parameter 1" fields in the IMMULITE 2000 Kit Information screen, which can be accessed from the menu via Data Entry: Kit Entry.)

Calculation of a signal/cutoff (s/co) ratio is done by using the following formula:

$$\text{S/CO Ratio} = \frac{\text{Sample or Control cps}}{\text{Mean Adjustor cps} \times P1}$$

Calculation and reporting of qualitative results (reactive / nonreactive / indeterminate) and signal/cutoff (s/co) ratio are handled automatically by the IMMULITE 2000.

The result for a sample is reported as "Indeterminate" if the counts per second for that sample fall within  $\pm 10\%$  of the cutoff. The result is reported as "Reactive" if the sample's counts are *above* the indeterminate range, and "Nonreactive" if *below* this range.

## Interpretation of Results

The cutoff of the IMMULITE 2000 Syphilis Screen assay was determined with reactive and nonreactive patient samples by ROC analysis, with a balanced consideration of sensitivity and specificity.

A result of "**Reactive**" (ratio of  $\geq 1.1$ ) indicates that the patient sample is reactive and that *T. pallidum* antibodies were detected in the sample.

A result of "**Nonreactive**" (ratio of  $< 0.9$ ) indicates that the patient sample is nonreactive and that *T. pallidum* antibodies were not detected in the sample.

Any result of "**Indeterminate**" (ratio  $\geq 0.9$  and  $< 1.1$ ) should be retested.

Samples which are reactive at the second test should be considered reactive. Samples which are non-reactive at the second test should be considered non-reactive. A second sample should be collected and tested no less than one week later when the result is repeatedly indeterminate.

The detection of *T. pallidum* total antibodies may indicate recent, past or successfully treated infections. This test cannot distinguish between treated and untreated infections.

The magnitude of the measured results (cps) above the Cutoff is not indicative of the total amount of antibodies detected.

The results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Syphilis Screen chemiluminescent EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

## Expected Values

Serum samples from 157 apparently healthy male and female volunteers were processed by the IMMULITE 2000 Syphilis Screen assay. These samples yielded a mean value of 0.12, a median value of 0.09 and a 95<sup>th</sup> percentile of 1.12. Results are expressed as a signal-to-cutoff ratio.

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

## Limitations

The IMMULITE 2000 Syphilis Screen assay allows testing for the presence of *T. pallidum* total antibodies. It detects both recent and past infections, but it cannot distinguish between different antibody classes.

Detection of *T. pallidum* total antibodies may indicate recent, past or successfully treated syphilis. This test, therefore, cannot discriminate between active and treated disease and, as a consequence, may not be used for determining response to therapy.

Test results are reported qualitatively as reactive or nonreactive for the presence of *T. pallidum* total antibodies. However, diagnosis of syphilis should not be established on the basis of a single test result, but should be determined in conjunction with clinical findings and other diagnostic procedures, as well as in association with medical judgment. Therefore a precise diagnosis should take into consideration clinical history, symptomatology and serological testing.

Assay interference due to circulating antibodies against hepatitis A virus, yaws, pinta and leptospirosis has not been evaluated. The user is responsible for establishing cross-reactivity performance with these infectious agents.

The results in HIV patients, in patients undergoing immunosuppressive therapy, or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

Heterophilic antibodies in human serum and plasma can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

## Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed as a signal-to-cutoff ratio. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

**High-dose Hook Effect:** Whenever samples containing extremely high concentrations are tested in a one-step sandwich assay, the high-dose hook effect can mimic concentrations lower than expected. Data show that the high-dose hook effect did not extend down to the cutoff level of the assay on six high reactive samples, that were serially diluted and tested, indicating no sample misclassification.

**Precision:** Samples were assayed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates. (See "Precision" table.)

**Reproducibility:** Reproducibility was conducted at 3 external clinical laboratories (2 in the U.S., 1 in the Netherlands) and included 2 different kit lots. Five samples targeting the assay range were prepared from pooled serum at Siemens Healthcare Diagnostics. The

same panel was aliquotted, frozen, and provided for assay at all 3 testing sites using at least four replicates per run, 2 runs per day across 5 different days. (See "Reproducibility" tables.)

**Specificity:** A study was conducted to evaluate whether the detection of *T. pallidum* antibody is affected by closely related microorganisms or other substances. Sera containing specific antibodies were tested by the IMMULITE 2000 Syphilis Screen assay with the following results (# reactive / # tested): CMV IgG (1/10), Epstein-Barr virus (1/10), ANA (0/10), HSV IgG (1/10), rubella IgG (1/10), toxoplasma IgG (0/10), rheumatoid factor (0/10), *Borrelia burgdorferi* IgG, US strain (2/10), *Borrelia burgdorferi* IgG, European strain (0/10), *Borrelia burgdorferi* IgM, European strain (0/10), anti-HBs (1/10), HCV (5/20) and HBsAg (0/10). Twelve out of 140 samples yielded reactive results. These 12 samples were further tested by a commercially available blot assay and yielded results positive for *T. pallidum*.

**Bilirubin:** Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

**Hemolysis:** Presence of hemoglobin in concentrations up to 593 mg/dL has no effect on results, within the precision of the assay.

**Lipemia:** Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

**Alternate Sample Type:** To assess the effect of alternate sample types, blood was collected from volunteers into plain, heparinized, and Becton Dickinson SST<sup>®</sup> plastic vacutainer tubes. Forty-eight out of 88 samples were spiked with various concentrations of *T. pallidum* and assayed by the IMMULITE 2000 Syphilis Screen procedure. Results are expressed as a signal-to-cutoff ratio.

(Heparin) = 1.06 (Serum) + 0.172  
 $r = 0.99$   
 $n = 88$

(SST) = 0.91 (Plain Tubes) + 0.105  
 $r = 1.00$   
 $n = 88$

Mean Ratios ( $n = 88$ ):  
3.29 (Serum)

3.65 (Heparin)  
3.09 (SST)

The use of heparinized plasma may increase ratio values at or near the cutoff level.

EDTA plasma is not recommended as a sample type.

**Method Comparison:** The assay was compared to a commercially available treponema assay on 280 samples.

IMMULITE 2000		
Predicate Assay	Reactive	Nonreactive
Positive	54	2
Negative	0	224

Positive Agreement:  $54/56 = 96.4\%$   
95% Confidence Interval: 87.7 – 99.6%  
Negative Agreement:  $224/224 = 100\%$   
Overall Agreement:  $278/280 = 99.3\%$   
95% Confidence Interval: 97.4 – 99.9%

**Clinical Performance:** Three clinical sites (2 U.S. and 1 Netherlands) recruited a total of 1,286 participants to compare performance of the IMMULITE 2000 Syphilis Screen to a commercially available treponema assay. Serum samples were analyzed in parallel with the predicate and test device (split approximately equally between 2 different test lots). Study participants may have qualified for more than one clinically important category (e.g., routine screening, medically confirmed syphilis patients, HIV-positive patients). Conservative estimates of overall agreement are presented that include indeterminate and equivocal results in the denominator.

Test results obtained from both the IMMULITE 2000 and predicate assays for all study participants ( $n = 1,286$ ) are presented below. The test assay demonstrated 99% positive agreement, negative agreement, and overall agreement with the predicate assay.

## Total Study Population

IMMULITE 2000	Predicate Assay			Total
	Positive	Negative	Equivocal	
Reactive	457	10	1	468
Nonreactive	0	814	0	814
Indeterminate	3	1	0	4
Total	460	825	1	1,286

Positive Agreement:  $457/460 = 99.3\%$   
(95% CI = 98.1% – 99.9%)

Negative Agreement:  $814/825 = 98.7\%$   
(95% CI = 97.6% – 99.3%)

Overall Agreement:  $1271/1286 = 98.8\%$   
(95% CI = 98.1% – 99.3%)

Medically diagnosed syphilis samples consisted of 281 patients (100 from the Netherlands and 181 from the U.S.) at various stages of disease ( $n = 78$  (27.8%) primary infection;  $n = 24$  (8.5%) secondary infection;  $n = 47$  (16.7%) latent infection;  $n = 132$  (47.0%) unknown stage of infection). Concordance of test results in this population showed 99% positive agreement, 75% negative agreement, and 98% overall agreement between the IMMULITE 2000 and predicate assays.

## Medically Diagnosed Syphilis Patients

IMMULITE 2000	Predicate Assay			Total
	Positive	Negative	Equivocal	
Reactive	270	2	1	273
Nonreactive	0	6	0	6
Indeterminate	2	0	0	2
Total	272	8	1	281

Positive Agreement:  $270/272 = 99.3\%$   
(95% CI = 97.4% – 99.9%)

Negative Agreement:  $6/8 = 75\%$   
(95% CI = 34.9% – 96.8%)

Overall Agreement:  $276/281 = 98.2\%$   
(95% CI = 95.9% – 99.4%)

924 samples were sent for routine syphilis testing. (243 from the Netherlands and 681 from the U.S.). Concordance of test results in this population showed 99% positive agreement, negative agreement,

and overall agreement between the IMMULITE 2000 and predicate assays.

### Samples Sent for Routine Syphilis Testing

IMMULITE 2000	Predicate Assay			Total
	Positive	Negative	Equivocal	
Reactive	359	5	0	364
Nonreactive	0	558	0	558
Indeterminate	2	0	0	2
Total	361	563	0	924

Positive Agreement: 359/361 = 99.4%  
(95% CI = 98.0% – 99.9%)

Negative Agreement: 558/563 = 99.1%  
(95% CI = 97.9% – 99.7%)

Overall Agreement: 917/924 = 99.2%  
(95% CI = 98.4% – 99.7%)

420 patients with HIV infection participated in the study (94 from the Netherlands and 326 from the U.S.). Concordance of test results in this population showed 99% positive agreement, 96% negative agreement, and 98% overall agreement between the IMMULITE 2000 and predicate assays.

### HIV-Positive Patients

IMMULITE 2000	Predicate Assay			Total
	Positive	Negative	Equivocal	
Reactive	260	7	0	267
Nonreactive	0	152	0	152
Indeterminate	1	0	0	1
Total	261	159	0	420

Positive Agreement: 260/261 = 99.6%  
(95% CI = 97.9% – 100%)

Negative Agreement: 152/159 = 95.6%  
(95% CI = 91.1% – 98.2%)

Overall Agreement: 412/420 = 98.1%  
(95% CI = 96.3% – 99.2%)

306 specimens independently testing positive on treponemal and non-treponemal assays showed 100% positive agreement (302/302) between the IMMULITE 2000 and predicate assays.

Two samples tested negative on both the IMMULITE 2000 and predicate assays, and two samples tested positive on the IMMULITE 2000 assay and negative on the predicate assay (50% negative agreement). Overall agreement was 99% (304/306; 95% Confidence Interval (CI) = 98% – 100%).

## References

- 1) Zrein M, Maure I, Boursier F, Soufflet L. recombinant antigen-based enzyme immunoassay for screening of *Treponema pallidum* antibodies in blood bank routine. J Clin Microbol. 1995 Mar; 33(3) 525-7.
- 2) Young H, Moyes A, Seagar L, McMillan A. Novel recombinant-antigen enzyme immunoassay for serological diagnosis of syphilis. J Clin Microbol. 1998 April; 36(4):913-7.
- 3) Rodriguez I, Alvarez EL, Fernandex C, Miranda A. Comparison of a recombinant-antigen enzyme immunoassay with *Treponema pallidum* hemagglutination test for serological confirmation of syphilis. Mem Inst Oswaldo Cruz. 2002 Apr; 97(3):347-9.
- 4) Centers for Disease Control. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 1988;37:377-82, 387-8.
- 5) Clinical and Laboratory Standards Institute (formerly NCCLS). Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline - Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. NCCLS Document M29-A3.
- 6) Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.

## Technical Assistance

For technical assistance, contact your National Distributor.

[www.siemens.com/diagnostics](http://www.siemens.com/diagnostics)

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

## Tables and Graphs

### Precision (ratio)

	Within-Run <sup>1</sup>			Total <sup>2</sup>	
	Mean <sup>3</sup>	SD <sup>4</sup>	CV <sup>5</sup>	SD	CV
1	0.09	0.004	4.4%	0.004	4.4%
2	0.11	0.006	5.5%	0.006	5.5%
3	0.82	0.051	6.2%	0.055	6.7%
4	1.01	0.065	6.4%	0.071	7.0%
5	1.03	0.063	6.1%	0.075	7.3%
6	1.51	0.084	5.6%	0.109	7.2%
7	2.40	0.138	5.8%	0.159	6.6%
8	6.39	0.488	7.6%	0.493	7.7%

**Dansk.** <sup>1</sup>I samme kørsel, <sup>2</sup>Total, <sup>3</sup>Middelværdi, <sup>4</sup>Standardafvigelse, <sup>5</sup>Variationskoefficient. **Eesti.** <sup>1</sup>Mõõtmisseriese sisene, <sup>2</sup>Seeriaveheline, <sup>3</sup>Keskmine väärtused, <sup>4</sup>Standardhälve, <sup>5</sup>Variatsioonikoeffitsient. **Latviski.** <sup>1</sup>Sērijas robežās, <sup>2</sup>Kopā, <sup>3</sup>Vidējā vērtība, <sup>4</sup>SD, <sup>5</sup>Variāciju koeficients. **Lietuviškai.** <sup>1</sup>Vieno tyrimo metu, <sup>2</sup>Bendras, <sup>3</sup>Vidurkiai, <sup>4</sup>SD, <sup>5</sup>CV. **Norsk.** <sup>1</sup>Innen serien, <sup>2</sup>Total, <sup>3</sup>Middelverdi, <sup>4</sup>SD, <sup>5</sup>CV. **Svenska.** <sup>1</sup>Inom körning, <sup>2</sup>Totalt, <sup>3</sup>Medel, <sup>4</sup>SD, <sup>5</sup>CV.

### Reproducibility (ratio) Lot 1: 3 Sites

Sample <sup>1</sup>	Reps <sup>2</sup>	Mean <sup>3</sup>	Within Day SD <sup>4*</sup>	Within Day %CV <sup>5</sup>	Between Site SD <sup>6**</sup>	Between Site %CV <sup>7</sup>
NC <sup>8</sup>	180	0.11	0.01	4.5	0.01	8.9
PC <sup>9</sup>	180	2.37	0.15	6.3	0.17	7.3
Low Neg <sup>10</sup>	120	0.11	0.01	10.4	0.02	22.6
High Neg <sup>11</sup>	120	0.78	0.04	5.1	0.06	7.7
Cut-off <sup>12</sup>	120	0.98	0.06	6.1	0.07	7.1
Low Pos <sup>13</sup>	120	1.43	0.10	6.6	0.14	9.8
High Pos <sup>14</sup>	120	6.13	0.33	5.4	0.53	8.7

### Reproducibility (ratio) Lot 2: 3 Sites

Sample <sup>1</sup>	Reps <sup>2</sup>	Mean <sup>3</sup>	Within Day SD <sup>4*</sup>	Within Day %CV <sup>5</sup>	Between Site SD <sup>6**</sup>	Between Site %CV <sup>7</sup>
NC <sup>8</sup>	180	0.06	0.00	N/A	0.01	10.9
PC <sup>9</sup>	180	2.63	0.15	5.81	0.16	6.15
Low Neg <sup>10</sup>	120	0.07	0.01	14.9	0.02	29.7
High Neg <sup>11</sup>	120	0.71	0.04	5.2	0.05	6.47
Cut-off <sup>12</sup>	120	0.91	0.05	5.96	0.07	7.51
Low Pos <sup>13</sup>	120	1.35	0.09	6.91	0.12	9.07
High Pos <sup>14</sup>	120	5.50	0.29	5.20	0.49	8.91

\* Multi-site within-run variation (pooled across sites)

\*\* Multi-site total variation (within-run + days + sites)

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.

© 2015 Siemens Healthcare Diagnostics. All rights reserved.

Made in: UK

Siemens Healthcare  
Diagnostics Products Ltd.  
Glyn Rhonwy, Llanberis,  
Gwynedd LL55 4EL  
United Kingdom

2015-05-21

PINL2KSY – 7 {15}

## Changes in this Edition:

cc#EU22534: 1) Under Warnings and Precautions, added biohazard symbol and potential biohazard statement, caution statement, required GHS pictograms/signal words/hazard and precautionary codes and statements, offending chemicals and sources. 2) Under References, added citations from the potential biohazard statement. 3) Changed "Origin: UK" to "Made in: UK". 4) Updated the manufacturer address to include Glyn Rhonwy. 5) In Understanding the Symbols, revised the symbol for catalog number, added a GHS pictogram for health hazard, and updated existing symbols with GHS pictograms and terminology for exclamation mark, corrosion, skull and crossbones, and environment.

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