Toxoplasma Quantitative IgG

**For the Quantitative Measurement of
IgG Antibodies to Toxoplasma gondii
in Human Serum**

For use on IMMULITE® 2000 systems

SIEMENS

IMMULITE® 2000

Toxoplasma Quantitative IgG

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE 2000 Systems Analyzers — for the quantitative measurement of IgG antibodies to *Toxoplasma gondii* in serum, as an aid in the determination of serological status to *Toxoplasma gondii*.

This kit is not FDA-cleared for use in the testing (i.e. screening) of blood or plasma donors.

Catalog Number: **L2KTXP2** (200 tests),
L2KTXP6 (600 tests)
Test Code: **TXP** Color: **Light Gray**

Summary and Explanation

Toxoplasma gondii is an obligate intracellular parasite capable of infecting most mammals, including humans. The organism is transmitted through ingestion of insufficiently cooked meat. Throughout the world, 1 to 90 percent of populations may be infected,¹ with 25 to 30 percent of the adult population infected in the United States.³ While toxoplasmosis may be manifested in several forms, infections are usually clinically inapparent, and latent infections usually persist for life.² Overt clinical symptoms are similar to infectious mononucleosis, with lymphadenopathy, fever, headache, malaise and sometimes pneumonia and myocarditis.⁶

As with other latent infections, acute toxoplasma infection can pose a serious threat to immunocompromised individuals and newborns who acquire the infection *in utero*. Immuno-suppressed patients may develop encephalitis, myocarditis or pneumonitis.³ Congenital infections usually result as a consequence of asymptomatic acute maternal infection. This infection can cause premature delivery, spontaneous abortion or stillbirth.^{4,6} Neonates may manifest chorioretinitis, hydrocephaly, microcephaly, cerebral calcification and psychomotor retardation.¹ The majority of

congenitally infected children will not exhibit any symptoms until later in life.⁴

Management of toxoplasmosis requires serological monitoring of infected individuals,³ as the organism is not readily available for culture. Quantitative testing for the presence of toxoplasma IgG can be useful to determine prior infection and indicate reactivation of the infection. Accurate diagnostic information is important, particularly during pregnancy, as treatment with spiramycin can reduce the risk to the fetus.⁵

Principle of the Procedure

IMMULITE 2000 *Toxoplasma* Quantitative IgG is a solid-phase, enzyme-labeled chemiluminescent sequential immunometric assay. The solid phase (bead) is coated with inactivated, partially purified *Toxoplasma gondii* antigen (RH strain tachyzoites from mouse peritoneum). The liquid phase consists of two reagents: 1) protein based buffer and 2) alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG antibody.

In the first cycle, the on-board diluted patient sample (1-in-20) and the protein based buffer are incubated together with the coated bead for 30 minutes. During this time, IgG antibodies specific to *Toxoplasma gondii* in the sample binds to inactivated *Toxoplasma gondii* antigen on the bead. Unbound sample is then removed by centrifugal washes.

In the second cycle, the enzyme conjugated monoclonal murine anti-human IgG antibody is added to the original reaction tube for additional 30 minutes incubation. The enzyme conjugated monoclonal murine anti-human IgG antibody binds to the immobilized IgG to form antibody sandwich complex. The unbound enzyme conjugate is removed by centrifugal washes. Finally, chemiluminescent substrate is added to the reaction tube containing the bead and the signal is generated in proportion to the bound enzyme.

Incubation Cycles: 2 × 30 minutes.

Time to First Result: 65 minutes.

Specimen Collection

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed samples may indicate mistreatment of a specimen before receipt by the laboratory; hence the results should be interpreted with caution.

Samples which are cloudy or have particulate material should be clarified by low-speed centrifugation.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufacturers may yield differing values, depending on materials and additives, including gel or physical barriers, clot activators and/or anticoagulants.

IMMULITE 2000 Toxoplasma Quantitative IgG has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 10 µL serum. After automatic sample predilution, 10 µL is used for the immunoassay.

Automatic Predilution Factor: 20.

Storage: 3 days at 2–8°C, or 6 months at –20°C.¹²

Warnings and Precautions

For *in vitro* diagnostic use.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative. On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

The bead is coated with *inactivated* toxoplasma antigen. However, caution is advised because of the possible presence of residual organism when working with, or disposing of, the materials supplied.

The anti-toxoplasma IgG results determined for a given specimen with assays from different manufacturers can vary due to differences in assay methods and reagent specificity. Therefore, the results reported by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Toxoplasma Quantitative IgG EIA. Results obtained from other manufacturers' assay methods may not be used interchangeably."

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

Toxoplasma Quantitative IgG Bead Pack (L2TXP12)

With barcode. 200 beads, coated with inactivated, partially purified *Toxoplasma gondii* antigen (RH strain tachyzoites from mouse peritoneum). Stable at 2–8°C until expiration date.

L2KTXP2: 1 pack. **L2KTXP6:** 3 packs.

Toxoplasma Quantitative IgG Reagent Wedge (L2TXPA2)

With barcode. 2 reagents: 11.5 mL of a protein-based buffer, with preservative. 11.5 mL alkaline phosphatase (bovine calf intestine) conjugated to monoclonal murine anti-human IgG antibody in buffer, with preservative. Stable at 2–8°C until expiration date.

L2KTXP2: 1 wedge. **L2KTXP6:** 3 wedges.

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top

of wedge; snap the sliding cover down into the ramps on the reagent lid.

Toxoplasma Quantitative IgG Adjustors (LTXPL, LTXPH)

Two vials (Low and High), 4 mL each, of human serum with IgG reactive to toxoplasma, in buffer, with preservative. Stable at 2–8°C for 30 days after opening, or for 6 months (aliquotted) at –20°C.

L2KTXP2: 1 set. **L2KTXP6:** 2 sets.

Before running adjustors or controls, place the appropriate Aliquot Labels (supplied with the kit) on test tubes, so that the barcodes can be read by the on-board reader.

Toxoplasma IgG Controls (LTXPC1, LTXPC2, LTXPC3)

Three vials (Negative, Low Positive and Positive), 2 mL each. **LTXPC1 (Negative Control):** One vial containing human serum with IgG nonreactive to toxoplasma, with preservative. **LTXPC2, LTXPC3 (Low Positive Control, Positive Control):** Two vials containing human serum with IgG reactive to toxoplasma, with preservative. Stable at 2–8°C for 14 days after opening, or for 6 months (aliquotted) at –20°C.

L2KTXP2: 1 set. **L2KTXP6:** 2 sets.

The IMMULITE 2000 software performs automatic on-board dilution of control samples, and the results will be tracked in the QC database. Enter controls as controls.

Refer to the control insert for concentration levels.

IgG/IgM Sample Diluent (L2IGZ2)

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

L2KTXP2: 1 vial. **L2KTXP6:** 1 vial.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2KTXP2: 3 labels. **L2KTXP6:** 5 labels.

Kit Components Supplied Separately

IgG/IgM Sample Diluent (L2IGZ2)

For the on-board dilution of patient samples and controls. 55 mL concentrated (ready-to-use) nonhuman protein/buffer matrix, with preservative. Stable at 2–8°C for 30 days after opening, or 6 months (aliquotted) at –20°C.

L2IGZ2: 1 vial.

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2IGZ2: 5 labels.

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LTXPCM: Tri-level Toxoplasma IgG Control Module

Also required

Distilled or deionized water; test tubes

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
4 weeks.

Quality Control Samples: Use Toxoplasma IgG Controls supplied with the kit.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control provided will not ensure precision of the cut-off.

In addition to the controls provided, users may wish to run additional controls to cover the upper range of the assay. Additional controls may be tested in accordance with guidelines or

requirements of local, state and/or federal regulations or accrediting organizations.¹¹

It is also recommended that known reactive and nonreactive specimens be run periodically to assure pipetting accuracy for the dilution step.

Interpretation of Results

IMMULITE 2000 Toxoplasma Quantitative IgG assay is standardized in terms of the World Health Organization's Third International Standard for Anti-Toxoplasma Serum. The assay's cutoff is determined from representative reactive and nonreactive samples (confirmed by an ELISA method) using the ROC (Receiver Operating Characteristic) technique.

Reactive: A result greater than or equal to 8 IU/mL is indicative of a past infection.

Nonreactive: A result less than 6.5 IU/mL is considered to be nonreactive. *Results less than 6.5 IU/mL should be reported as Nonreactive only, and no numerical concentration should be reported.*

Indeterminate : A result greater than or equal to 6.5 IU/mL and less than 8 IU/mL is considered indeterminate.

Samples greater than 250 IU/mL should be reported as "greater than 250 IU/mL". Specimens reporting a result of greater than 250 IU/mL, when tested at the required 1-in-21 auto-dilution, may be further diluted with the IMMULITE 2000 IgG/IgM Sample Diluent.

The presence of IgG antibodies to *Toxoplasma gondii* is an indication of previous exposure to the organism. A single specimen can be used to estimate the serological status of the individual.

A result of "indeterminate" indicates an intermediate level of *Toxoplasma gondii* antibodies in the sample: the sample should therefore be retested. Samples which still test as "indeterminate" should be examined by an alternate method, or a second sample should be collected, if possible, within a reasonable period of time (e.g. one week).

Any interpretation of the results is not recommended if the positive or negative control falls outside the range specified in the Quality Control section.

The magnitude of the reported IgG level cannot be correlated to an endpoint titer.

Reports by the laboratory to the physician should include: "The following results were obtained with the IMMULITE 2000 Toxoplasma Quantitative IgG EIA. Values obtained from other manufacturers' assay methods may not be used interchangeably."

Expected Values

Individuals with the toxoplasma organism may not exhibit detectable levels of IgG antibody in the early stages of infection. Levels of IgG antibody to *Toxoplasma gondii* begin to rise one or two weeks after infection. Peak levels are reached in 6 to 8 weeks, then gradually decline over a period of months or even years. Low titers are generally detectable for life. The antibody titer does not correlate with severity of illness.^{2,4,6}

Consider these limits as *guidelines* only. Each laboratory should establish its own reference ranges.

Limitations

The use of IMMULITE 2000 Toxoplasma Quantitative IgG to diagnose recent infection by testing paired sera has not been validated.

The results of the test must be taken within the context of the patient's clinical history, symptomology and other laboratory findings.

The presence of IgG antibodies in a single specimen is not sufficient to distinguish between active and past infection. Patients suspected of having a primary or active toxoplasma infection should be tested for the presence of IgM antibodies to *Toxoplasma gondii*.

For the determination of seroconversion from nonreactive to reactive, a second serum sample should be drawn three to four weeks after the acute, nonreactive sample, during the convalescent stage of the infection. The acute phase sample should be stored and tested in parallel with the convalescent sample.

Individuals with an acute toxoplasma infection may not exhibit any detectable IgG antibodies at the early stage of infection.

The results in HIV patients, in patients undergoing immunosuppressive therapy,

or in patients with other disorders leading to immunosuppression, should be interpreted with caution.

The performance characteristics of this assay have not been established for use with specimens from neonates, cord blood, pretransplant patients, or body fluids other than serum, such as urine, saliva or amniotic fluid.

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34:27-33.] Samples from patients routinely exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in IU/mL. (Unless otherwise noted, all were generated on serum samples collected in tubes without gel barriers or clot-promoting additives.)

Reportable Range: 5.0 – 250 IU/mL. Standardized in terms of WHO 3rd IS for Anti-Toxoplasma Serum, Human.

Analytical Sensitivity: ≤ 5.0 IU/mL.

Precision: Samples were assayed in duplicate in 40 runs for a total of 80 replicates. (See "Precision" table.)

Specificity: The specificity of the IMMULITE 2000 Toxoplasma Quantitative IgG assay was evaluated by testing 73 specimens reactive for *Mycoplasma pneumoniae*, syphilis, varicella-zoster virus, parvovirus B19, cytomegalovirus (CMV), Epstein-Barr virus, HSV 1 and measles. With these specimens, IMMULITE 2000 Toxoplasma Quantitative IgG and IMMULITE Toxoplasma

Quantitative IgG show 99% agreement (72/73).

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay.

Hemolysis: Presence of hemoglobin in concentrations up to 539 mg/dL has no effect on results, within the precision of the assay.

Lipemia: Presence of triglycerides in concentrations up to 3,000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Type: To assess the effect of alternate sample types, blood was collected from 15 volunteers into plain, heparinized, EDTA and Becton Dickinson SST[®] vacutainer tubes. Equal volumes of the matched samples were spiked at various levels with sera containing IgG antibodies to *Toxoplasma*, and then assayed by the IMMULITE 2000 Toxoplasma Quantitative IgG procedure, with the following results.

(Heparin) = 1.1 (Serum) – 0.81 IU/mL
 $r = 0.985$

(EDTA) = 1.0 (Serum) + 1.4 IU/mL
 $r = 0.949$

(SST) = 1.1 (Plain Tubes) – 2.9 IU/mL
 $r = 0.915$

Means:
65 IU/mL (Serum)
68 IU/mL (Heparin)
67 IU/mL (EDTA)
66 IU/mL (SST)

Method Comparison 1: The assay was compared to IMMULITE Toxoplasma Quantitative IgG on 33 samples. (Concentration range: approximately 5 to 112 IU/mL. See "Method Comparison 1" graph.) By linear regression:

(IML 2000) = 1.0 (IML) + 2.0 IU/mL
 $r = 0.994$

Means:
30 IU/mL (IMMULITE 2000)
28 IU/mL (IMMULITE)

Method Comparison 2: The assay was compared to IMMULITE Toxoplasma Quantitative IgG on 40 samples. (Concentration range: approximately 5 to 162 IU/mL. See "Method Comparison 2" graph.) By linear regression:

(IML 2000) = 1.04 (IML) + 1.3 IU/mL
 $r = 0.963$

Means:

35 IU/mL (IMMULITE 2000)

33 IU/mL (IMMULITE)

References

1) Daffos F. Prenatal management of 746 pregnancies at risk for congenital toxoplasmosis. *N Engl J Med* 1988;318: 271-5. 2) Krick JA, Remington JS. Toxoplasmosis in the adult: an overview. *N Engl J Med* 1978;298: 550-3. 3) Krogstad DJ, et al. Blood and tissue protozoa. In: Lennette EH, et al, editors. *Manual of clinical microbiology*. 4th ed. Washington, D.C.: American Society for Microbiology, 1985: 612-30. 4) McCabe RE, Remington JS. *Toxoplasma gondii*. In: Mandell, Douglas, Bennett, editors. *Principles and practice of infectious diseases*. 2nd ed. New York: John Wiley, 1985. 5) Jeannel D, et al. What is known about the prevention of congenital toxoplasmosis? *Lancet* 1990;336:359-61. 6) Remington JS. McLeod R. Toxoplasmosis. In: Braude, editor. *Microbiology and infectious diseases*. Philadelphia: Saunders, 1981. 7) Walls KW. Toxoplasmosis. In: Balows A, et al, editors. *Laboratory diagnosis of infectious diseases*. New York: Springer-Verlag, 1988: 998-1017. 8) Baltz ML, Searcy RL. Clinical significance and advanced serologic diagnosis of ToRCH infections. *Am Clin Lab* 1994 (Mar/Apr): 18-23. 9) Wilson M, et al. Seroepidemiology of toxoplasmosis in the United States. Abstract presented at the American Society for Microbiology, May 1994. 10) Approved standard procedures for the handling and processing of blood specimens, NCCLS Document H18-A, 1991. 11) International Quality Control testing: principles and definitions, NCCLS Document C24-A, 1991. 12) Tietz NW, editor. *Clinical guide to laboratory tests*. 3rd ed. Philadelphia: WB Saunders, 1995:358. 13) National Committee for Clinical Laboratory Standards. *Procedures for the collection of diagnostic blood specimens by venipuncture*; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998.

Additional references : Abdalla KF. Congenital toxoplasmosis among premature infants with different clinical pictures in Saudi Arabia. *J Egypt Soc Parasitol*, 1994 Dec. ♦ Albrecht H. Disseminated toxoplasmosis in AIDS patients – report of 16 cases. *Scand J Infect Dis* 1995. ♦ Attia RA. Toxoplasma IgG & IgM antibodies. A case control study. *J Egypt Soc Parasitol*. 1995 Dec. ♦ Bretagne S. Central nervous system toxoplasmosis in AIDS (letter, comment). *N Engl J Med* 1993 May 6. ♦ Brindle R. Toxoplasma antibodies in HIV-positive patients from Nairobi. *Trans R Soc Trop Med Hyg*, 1991 Nov-Dec. ♦ Camargo ME. Avidity of specific IgG antibodies as markers of recent primary infection caused by *Toxoplasma gondii*. *Rev Inst Med Trop Sao Paulo*, 1991 May-Jun. ♦ Choi WY. Foodborne outbreaks of human toxoplasmosis. *J Infect Dis*,

1997 May. ♦ Derouin F. Predictive value of *Toxoplasma gondii* antibody titres on the occurrence of toxoplasmic encephalitis in HIV-infected patients. *ANRS 005/ACTG 154 Trial Group*. *AIDS* 1996 Nov. ♦ Durmaz R. Seropositivity of toxoplasmosis among reproductive-age women in Malatya, Turkey. *J Egypt Soc Parasitol*, 1995 Dec. ♦ El-Nawawy A. Maternal and neonatal prevalence of toxoplasma and cytomegalovirus (CMV) antibodies and hepatitis-B antigens in an Egyptian rural area. *J Trop Pediatr*, 1996 Jun. ♦ Gilbert RE. Prevalence of toxoplasma IgG among pregnant women in west London according to country of birth and ethnic group (see comments). *BMJ*, 1993 Jan 16. ♦ Gomez-Marin JE. A maternal screening program for congenital toxoplasmosis in Quindio Colombia and application of mathematical models to estimate incidences using age-stratified data. *Am J Trop Med Hyg*, 1997 Aug. ♦ Guerina NG. Neonatal serologic screening and early treatment for congenital *Toxoplasma gondii* infection : The New England Regional Toxoplasma Working Group I see comments. *N Engl J Med*, 1994 Jun 30. ♦ Hajeer AH. *Toxoplasma gondii* : detection of antibodies in human saliva and serum. *Parasite Immunol*, 1994 Jan. ♦ Harning D. Recombinant *Toxoplasma gondii* surface antigen 1 (P30) expressed in *Escherichia coli* is recognized by human Toxoplasma-specific immunoglobulin M (IgM) and IgG antibodies. *Clin Diagn Lab Immunol*, 1996 May. ♦ Hellerbrand C. High predictive value of *Toxoplasma gondii* IgG antibody levels in HIV-infected patients for diagnosis of cerebral toxoplasmosis. *Eur J Clin Microbiol Infect Dis*, 1996 Nov. ♦ Idris MA. Prevalence of IgG antibodies against *Toxoplasma gondii* in human sera from Dhofar, Oman. *Ann Trop Med Parasitol*, 1994 Feb. ♦ Israelski DM. Prevalence of Toxoplasma infection in a cohort of homosexual men at risk of AIDS and toxoplasmic encephalitis (see comments). *J Acquir Immune Defic Syndr*, 1993 Apr. ♦ Jacquier P. Epidemiology of toxoplasmosis in Switzerland : national study of seroprevalence monitored in pregnant women 1990-1991. *Schweiz Med Wochenschr Suppl*, 1995. ♦ Johnson AM. Cloning of *Toxoplasma gondii* gene fragments encoding diagnostic antigens. *Gene*, 1991 Mar 1. ♦ Julvez J. Seroepidemiology of toxoplasmosis in Niamey, Niger. *Med Trop (Mars)*, 1996. ♦ Konishi E. Naturally occurring immunoglobulin M antibodies to *Toxoplasma gondii* in Japanese populations. *Parasitology* 1991 Apr. ♦ Lappalainen M. Outcome of children after maternal primary Toxoplasma infection during pregnancy with emphasis on avidity of specific IgG. *The Study Group. Pediatr Infect Dis J* 1995 May. ♦ Lebech M. Neonatal screening for congenital toxoplasmosis in Denmark : presentation of the design of a prospective study. *Scand J Infect Dis Suppl*, 1992. ♦ Lebech M. Occurrence of

toxoplasmosis in pregnant women in Denmark. A study of 5402 pregnant women. Ugeskr Laeger, 1995 Sep 18. ♦ Low incidence of congenital toxoplasmosis in children born to women infected with human immunodeficiency virus. European Collaborative Study and Research Network on Congenital Toxoplasmosis. 1996 Sep. Eur J Obstet Gynecol Reprod Biol, 1996 Sep. ♦ Petithory JC. Performance of European laboratories testing serum samples for *Toxoplasma gondii*. Eur J Clin Microbiol Infect Dis, 1996 Jan. ♦ Rai SK. Immunoglobulin M and G antibodies in mice in response to *Toxoplasma gondii* (S-273) infection and their antigen recognition patterns in western blotting on various post-infection days. Kansenshogaku Zasshi, 1994 Nov. ♦ Sadaruddin A. Seroepidemiology of *toxoplasma gondii* infection in young school children in Islamabad. JPMA J Pak Med Assoc, 1991 Jun. ♦ Sahwi SY. Toxoplasmosis as a cause of repeated abortion. J Obstet Gynaecol, 1995 Apr. ♦ Vinhal FA. Analysis of aqueous humor in ocular toxoplasmosis : detection of low avidity IgG specific to *Toxoplasma gondii*. Appl Parasitol, 1994 Feb. ♦ Zufferey J. Prevalence of latent toxoplasmosis and serological diagnosis of active infection in HIV-positive patients. Eur J Clin Microbiol Infect Dis, 1993 Aug.

Technical Assistance

For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485:2003.

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IMMULITE® is a trademark of Siemens Healthcare Diagnostics.

©2009 Siemens Healthcare Diagnostics Inc. All rights reserved.

Origin: UK

Siemens Healthcare
Diagnostics Products Ltd.
Llanberis, Gwynedd LL55 4EL
United Kingdom

2009-08-10

PINL2KTXP – 1 {22}

Changes in this Edition:

cc#17730: Updated Company name to Siemens Healthcare Diagnostics. Changed Manufacturer's address to Llanberis, UK. Added country of origin, website address, trademark and copyright statements, symbol glossary section. Standardized Technical Assistance section, changed the phone number. Removed EC Rep (responsibility transferred to Llanberis, UK).

cc#16311: Created a new package insert (PINL2KTXP – 1 {22}) for use in Northern Europe (Danish, Estonian, Latvian, Lithuanian, Swedish).

cc#EU18560: Added references to control module and sample diluent in Kit Components Supplied Separately.

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