12 Suderinama įranga

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Reprocessing of Flexible Endoscopes Olympus Compatibility Statement

1. Reprocessing with the Olympus ETD System

Based on the Olympus Tokyo test reports and our practical experience, we herewith state that below listed automatic processes for endoscope reprocessing have been found compatible with Olympus flexible endoscopes. The processes have been tested for the following machines: ETD4, ETD3, ETD2, miniETD2, miniETD, ETD Double.

The hygienic efficacy of the ETD processes has been thoroughly tested by independent hygienic institutes. For details on test methods and results, please ask for the ETD certificates.

Product Name	Distributor	Major Ingredient
ETD PAA Process with - Olympus EndoDet - Olympus EndoDis - Olympus EndoAct	OLYMPUS	Peracetic Acid
ETD GA Process with (*) - Olympus CLEANER (CETD) - Olympus DISINFECTANT (DETD)	OLYMPUS	Glutaraldehyde

(*) not for ETD Double

REMARKS

It is known from our field experience that applying a combination of different processes for the reprocessing of endoscopes might accelerate damages to materials of the instruments. The above information is only valid for the exclusive influence of the mentioned process using the products named above.

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OLYMPUS ENDODET, ENDOACT AND ENDODIS

PAA process chemistry

Olympus process chemicals for ETD machines provide comprehensive and safe reprocessing of your flexible, heat-sensitive endoscopes. With EndoDet, EndoAct and EndoDis, Olympus has created an optimally adapted, peracetic acid (PAA) based, process chemical system. Fully functional at low temperatures, it is largely effective while being fully compatible with your flexible Olympus endoscopes.

	EndoDet	EndoDis	EndoAct
Product description	EndoDet is a surfactant-based detergent, which allows for optimum cleaning efficacy at lower temperature.	EndoDis is a disinfectant based on peracetic acid (PAA), providing a wide range of efficacy.	EndoAct allows the Olympus PAA process to work in a neutral pH value range and thus efficiently protects sensitive endoscope materials. It is an additive to EndoDis.
			
Composition	Non-ionic surfactants, glycols, solubilizers, cleaning booster	Active ingredient: Peroxyacetic acid Other ingredients: Hydrogen peroxide acetic acid, stabilizer	Phosphates, sodium hydroxide
Spectrum of effect	EndoDet effectively reduces organic debris and protein residues with special cleaning booster.	EndoDis, in combination with EndoAct is largely effective against relevant pathogens in endoscopy and complies with European standards in state-of-the-art endoscope reprocessing. Efficacy (under clean conditions): EN 13727: Bactericidal activity EN 13624: Fungicidal activity EN 14348: Mycobactericidal activity EN 14476: Virucidal activity EN 13704: Sporocidal activity Certificates and documentation available upon request	
Direction for use	EndoDet, EndoDis and EndoAct are all supplied as a concentrate. Place it into your ETD washer-disinfector according to the instruction manual. The concentrate will then be automatically dosed into the machine process by the dosage pumps.		
Dosage (ETD3)	Dosage: 0.6% (6 ml/l) Contact time: 3 min Temperature: 35°C	Dosage: 1.2% (12 ml/l each) Contact time: 5 min Temperature: 35°C	
According to standard ETD process settings			
Presentation	5 l canister (sales unit: 3 canisters)	2.8 l canister (sales unit: 3 canisters)	5 l canister (sales unit: 3 canisters)
Biocompatibility	Biocompatibility has been tested and proven according to ISO 10993. Simulated use tests demonstrated no influence on the viability of cells in cytotoxicity test from residues of the process chemicals		

Note: Protective equipment and clothing should always be worn when handling chemicals. Refer to the material safety data sheet.

Manufactured by: Ecolab Deutschland GmbH, Postfach 130406, D-40554 Düsseldorf, Tel:0211-9893-900, www.ecolab.com

**„OLYMPUS“ CHEMIKALAI „ENDODET“, „ENDOACT“ IR „ENDODIS“
„PAA“ technologiniai chemikalai**

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