# Troponin T hs

cobas®

REF		SYSTEM
05092744 190	200	Elecsys 2010 MODULAR ANALYTICS E170 <b>cobas e 411</b> <b>cobas e 601</b> <b>cobas e 602</b>

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## Nuorodos

- European patent 394816 and US patent 6376206 by Roche Diagnostics GmbH. Specific antibodies to Troponin T, their production and use in a reagent for the determination of myocardial necrosis.
- Katus HA, Remppis A, Looser S, et al. Enzyme linked immunoassay of cardiac troponin T for the detection of acute myocardial infarction in patients. Mol Cell Cardiol 1989;21(7):1349-1353.
- Katus HA, Scheffold T, Remppis A, et al. Proteins of the troponin complex. Laboratory Medicine 1992;23(5):311-317.
- Alpert JS, Thygesen K, Jaffe A, et al. The universal definition of myocardial infarction: a consensus document: Ischaemic heart disease. Heart 2008;94:1335-1341.
- Hamm CW, Ravkilde J, Gerhardt W, et al. The prognostic value of serum troponin T in unstable angina. N Engl J Med 1992;327(3):146-150.
- Ohmann EM, Armstrong PW, Christenson RH, et al. Cardiac troponin T levels for risk stratification with admission cardiac troponin T levels in acute myocardial ischemia. GUSTO IIA Investigators. N Engl J Med 1996;335:1333-1334.
- Christenson RH, Duh SH, Newby LK, et al. Cardiac troponin T and cardiac troponin I: relative values in short-term risk stratification of patients with acute coronary syndromes. Clin Chem 1998;44(3):494-501.
- Lindahl B, Diderholm E, Lagerqvist B, et al. Mechanisms behind the prognostic value of Troponin T in unstable coronary artery disease: a FRISC II substudy. J Am Coll Cardiol 2001;38:979-986.
- Aviles RJ, Askari AT, Lindahl B, et al. Troponin T levels in patients with acute coronary syndromes, with or without renal dysfunction. N Engl J Med 2002;346:2047-2052.
- Lindahl B, Venge P, Wallentin L. Troponin T Identifies Patients With Unstable Coronary Artery Disease Who Benefit From Long-Term Antithrombotic Protection. J Am Coll Cardiol 1997;29(1):43-48.
- Hamm CW, Heeschen C, Goldmann B, et al. Benefit of abciximab in patients with refractory unstable angina in relation to serum troponin T levels. N Engl J Med 1999;340(21):1623-1629.
- Heeschen C, Hamm CW, Goldmann BU, et al. for PRISM study investigators. Troponin concentrations for stratification of patients with acute coronary syndromes in relation to therapeutic efficacy of tirofiban. Lancet 1999;354:1757-1762.
- Lindahl B, Diderholm E, Lagerqvist B, et al. Effects on mortality of long-term treatment with l.m.w. heparin in relation to troponin T level and ECG findings - a FRISC 2 substudy. Eur Heart J 2000;21(Suppl.):521.
- Newby LK, Ohman EM, Christenson RH, et al. Benefit of Glycoprotein IIb/IIIa Inhibition in Patients With Acute Coronary Syndromes and Troponin T-Positive Status: The PARAGON-B Troponin T Substudy. Circulation 2001;103:2891-2896.
- The Task Force for the diagnosis and Treatment of Non-ST-Segment Elevation Acute Coronary Syndromes of the European Society of Cardiology. Guidelines for the diagnosis and treatment of non-ST elevation acute coronary syndromes. European Heart Journal 2007;28:1598-1660.
- Thygesen K, Alpert JS, White HD on behalf of the Joint ESC/ACCF/AHA/WHF Task Force for Redefinition of Myocardial Infarction. JACC 2007;50:2173-95.
- Latini R, Masson S, Anand IS, et al. Prognostic Value of Very Low Plasma Concentrations of Troponin T in Patients with Stable Chronic Heart Failure. Circulation 2007;116:1242-1249.
- Omland T, De Lemos JA, Christophi C, et al. Distribution and determinants of very low levels of cardiac troponin T in patients with stable coronary artery disease: The PEACE trial. Eur Heart J 2008;9(202):1342.
- Omland T, De Lemos JA, Christophi C, et al. Very low cardiac troponin T concentrations and cardiovascular events in patients with stable coronary artery disease: The PEACE trial. Eur Heart J 2008;29(202):1644.
- Allen J, Kleiman NS, Nassif D, et al. Prevalence and prognostic Significance of preprocedural cardiac troponin elevation among patients with stable coronary artery disease undergoing percutaneous coronary intervention. Circulation 2008;118:632-638.

# Troponin T hs

- 21 Sato Y, Kataoka K, Matsumori A, et al. Measuring serum aminoterminal type III procollagen peptide, 7S domain of type IV collagen, and cardiac troponin T in patients with idiopathic dilated cardiomyopathy and secondary cardiomyopathy. *Heart* 1997;78:505-508.
- 22 Dierkes J, Domrose U, Westphal S, et al. Cardiac troponin T predicts mortality in patients with renal disease. *Circulation* 2000;102:1964-1969.
- 23 Ooi DS, Zimmerman D, Graham J, et al. Cardiac troponin T predicts long-term outcomes in hemodialysis patients. *Clin Chem* 2001;47:412-417.
- 24 Apple FS, Wu AHB. Myocardial infarction redefined: Role of cardiac troponin Testing. *Clin Chem* 2001;47:377-379.
- 25 deFilippi C, Wasserman S, Rosanio S, et al. Cardiac Troponin T and C-reactive protein for predicting prognosis, coronary atherosclerosis, and cardiomyopathy in patients undergoing long term hemodialysis. *JAMA* 2003;290:353-359.
- 26 Iliou MC, Fumeron C, Benoit MO, et al. Prognostic value of cardiac Markers in ESRD: Chronic Hemodialysis and New Cardiac markers Evaluation (CHANCE) study. *Am J Kidney Dis* 2003;42:513-523.
- 27 Scott B, Deman A, Peeters P, et al. Cardiac troponin T and malondialdehyde modified plasma lipids in haemodialysis patients. *Nephrol Dial Transplant* 2003;18:737-742.
- 28 Madsen LH, Ladefoged S, Hildebrandt P, et al. Comparison of four different cardiac troponin assays in patients with end-stage renal disease on chronic haemodialysis. *Acute Card Care* 2008;10(3):173-180.
- 29 ver Elst KM, Spapen HD, Nguyen DN, et al. Cardiac troponin I and T are biological markers of left ventricular dysfunction in septic shock. *Clin Chem* 2000;46:650-657.
- 30 Rossing P, Jorsal A, Tarnow L, et al. Plasma hs-Troponin T predicts cardiovascular and all cause mortality as well as deterioration in kidney function in type 1 diabetic patients with nephropathy. Abstract EASD 2008.
- 31 European patent application 890154 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Cardiac Troponin as an indicator of advanced coronary artery disease.
- 32 European patent application 1837659 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Means and methods for the differentiation of acute and chronic myocardial necrosis in symptomatic patients.
- 33 Masson S, Pirelli S, Pulitano G, et al. Circulating Markers of myocyte injury predict first recurrence of atrial fibrillation. *AHA* 2008 November 8-12. New Orleans.
- 34 Lauer B, Niederau C, Kühl U, et al. Cardiac troponin T in patients with clinically suspected myocarditis. *JACC* 1997;30:1354-1359.
- 35 Swaanenburg JCJM, Klaase JM, DeJongste MJL, et al. Troponin I, troponin T, CK-MB-activity and CK-MB mass as markers for the detection of myocardial contusion in patients who experienced blunt trauma. *Clin Chim Acta* 1998;272:171-181.
- 36 Giannitsis E, Müller-Bardorff M, Kurowski V, et al. Independent prognostic value of cardiac troponin T in patients with confirmed pulmonary embolism. *Circulation* 2000;102:211-217.
- 37 Herman EH, Lipshultz SE, Rifai N, et al. Use of cardiac troponin T levels as an indicator of doxorubicin-induced cardiotoxicity. *Cancer Res* 1998;58:195-197.
- 38 European patent 1363128 by F Hoffmann-La Roche AG and Roche Diagnostics GmbH. Diagnosis of myocardial infarction and acute coronary syndrome by combination of markers.
- 39 European patent application 1925943 by F. Hoffmann-La Roche AG and Roche Diagnostics GmbH. Means and methods for optimization of diagnostic and therapeutic approaches in chronic artery disease based on the detection of Troponin T and NT-pro BNP.
- 40 Davis GK, Labugger R, Van Eyk JE, et al. Cardiac troponin T is not detected in western blots of diseased renal tissue. *Clin Chem* 2001;47:782-783.
- 41 Ricchiuti V, Voss EM, Ney A, et al. Cardiac Troponin T isoforms expressed in renal diseased skeletal muscle will not cause false positive results by the second generation cTnT assay by Boehringer Mannheim; *Clin Chem* 1998;44(9):1919-1924.
- 42 Hallermayer K, Klenner D, Vogel R. Use of recombinant human cardiac troponin T for standardization of third generation troponin T methods. *Scand J Clin Invest* 1999;59(Suppl 230):128-131.
- 43 World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria for diagnosis of ischemic heart disease. *Circulation* 1979;59:607-609.
- 44 Müller-Bardorff M, Hallermayer K, Schröder A, et al. Improved troponin T ELISA specific for cardiac troponin T isoform: assay development and analytical validation. *Clin Chem* 1997;43(3):458-466.
- 45 Kampmann M, Rauscher T, Müller-Bardorff M, et al. Clinical Evaluation of Troponin T and CK-MB Mass on the Elecsys 2010 analyzer. Poster presented at Medlab 97, 12th IFCC European Congress of Clinical Chemistry, August 17-22, 1997, Basel, Switzerland.
- 46 Mendis S, Thygesen K, Kuulasmaa K, et al. World Health Organization definition of myocardial infarction: 2008-09 revision. *Int J Epidemiol* 2011;40(1):139-146.
- 47 Antmann E, Bassand J-P, Klein W, et al. Myocardial infarction redefined - a consensus document of the Joint European Society of Cardiology/American College of Cardiology committee for the redefinition of myocardial infarction: The Joint European Society of Cardiology/American College of Cardiology Committee. *Am Coll Cardiol* 2000;36:959-969.
- 48 Kobayashii S, Tanaka M, Tamura N, et al. Serum cardiac troponin T in polymyositis/dermatomyositis. *Lancet* 1992;340(8821):726.
- 49 Mair J, Dienstl F, Puschendorf B. Cardiac Troponin T in the Diagnosis of Myocardial Injury. *Critical Reviews in Clin Lab Sci* 1992;29(1):31-57.
- 50 Carrier M, Solymoss BC, Raymond C, et al. Cardiac Troponin T and Creatine Kinase MB Isoenzyme as Biochemical Markers of Ischemia after Heart Preservation and Transplantation. *J Heart Lung Transplant* 1994;13(4):696-700.
- 51 Löfberg M, Tähtelä R, Härkönen M, et al. Myosin Heavy-Chain Fragments and Cardiac Troponins in the Serum in Rhabdomyolysis. Diagnostic Specificity of New Biochemical Markers. *Arch Neurol* 1995;52:1210-1214.
- 52 Anderson JR, Hossein-Nia M, Brown P, et al. Donor Cardiac Troponin T Predicts Subsequent Inotrope Requirements Following Cardiac Transplantation. *Brief Communications in Transplant* 1994;58(9):1056-1057.
- 53 Franz WM, Remppis A, Kandolf R, et al. Serum Troponin T: Diagnostic Marker for Acute Myocarditis. Letter to the Editor in *Clin Chem* 1996;42(2):340-341.
- 54 Mair P, Mair J, Koller J, et al. Cardiac troponin T in the diagnosis of heart contusion. *Lancet* 1991;338:693.
- 55 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. *J Clin Chem Clin Biochem* 1988 Nov;26(11):783-790.

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© 2013, Roche Diagnostics
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim