cobas b 123 AutoQC Pack Tri-Level

(EN)

REF 05169933 001	8 x 1,0 mL, Level 1 8 x 1,0 mL, Level 2 8 x 1,0 mL, Level 3
-------------------------	---

Intended use

The **cobas b 123** AutoQC pack Tri-Level is a multi-analyte control intended for use as control material to monitor the measurement of pH, PO_2 , PO_3 , SO_2 , Na^+ , K^+ , Cl^- , Ca^{2+} , Hct, Hb and Hb derivatives as well as glucose, lactate and bilirubin on **cobas b 123** POC systems with an AutoQC module.

Product description

The **cobas b 123** AutoQC pack Tri-Level control contains 24 single glass ampoules of multi-analyte controls in 3 levels:

- Level 1**: respiratory acidosis represented by low pH and high PO_2 values, with decreased PO_3 , low values for Na^+ , K^+ , Cl^- , Hb and bilirubin, normal values for glucose and high values for Hct, Ca^{2+} and lactate.
- Level 2**: representing concentrations of normal values for all parameters except glucose, which is low.
- Level 3**: respiratory alkalosis represented by high pH and low PO_2 values, with increased PO_3 , high values for Na^+ , K^+ , Cl^- , Hb, glucose and bilirubin and low values for Hct, Ca^{2+} and lactate.

The **cobas b 123** AutoQC pack Tri-Level control material is formulated using aqueous based, dye containing solution equilibrated with predetermined levels of oxygen, carbon dioxide and nitrogen with salts, organic and carbonate buffers and metabolites.

Precautions and warnings

For in vitro diagnostic use. Exercise the normal precautions required for handling all laboratory reagents. Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request. **Caution:** Handle with care! Store the AutoQC pack upright! Do not insert damaged packs into the AutoQC module. Do not invert the AutoQC pack after partial use in case the pack is reinstalled (used QC ampoules can leak). The reinstallation of a partially-used AutoQC pack which is inverted can damage the AutoQC module.

Storage and stability

- Unopened: up to the stated expiration date at 2-8 °C.
- In use: up to 3 months when installed in the analyzer.
- Do not freeze! Roche Diagnostics cannot guarantee the performance of the quality controls when stored at temperatures outside the recommended range.

Procedure

- Allow a new AutoQC pack to come to equilibrium for at least 24 hours at room temperature prior to use, but not exceeding 7 days. See the information on the instruments "overview" screen.
- Open the analyzer front door and insert the new AutoQC pack into the instrument. Do not invert the AutoQC pack.
- The pack-specific information: lot number, expiration date, target values and ranges are stored on the memory chip on the AutoQC pack. This information is transferred to the analyzer software automatically during installation.
- Close the analyzer front door.

Install the AutoQC pack in the analyzer according to the instructions in the Instructions for use for the respective analyzer. This procedure is guided by the analyzer software. Refer to the latest Instructions for use of the respective analyzer for details on the number of daily QC measurements required for each level.

Limitations

- The figures do not contain red blood cells, and thus may not detect any system malfunctions that can cause lysis of patient red blood cells leading to erroneous results.
- This product is intended for use as quality control, and cannot be used as a calibration standard.
- The listed values only apply to the indicated lot numbers.
- These controls are used to aid in evaluation of the performance of blood gas, electrolyte, metabolite and oximeter instrumentation, and should not be considered as a substitute for other aspects of total quality control.

Target values and ranges

The figures in the value sheet show the mean value and the target range of each level. The values obtained from the analyzer must fall within the defined ranges. Each defined range contains 95 % of all obtained quality control results. Refer to the Quality control section in the Instructions for use for the respective analyzer when the quality control values are outside the expected range. Each laboratory must ensure that corrective measures are taken if values fall outside the limits.

Follow the applicable government regulations and local guidelines for quality control.

Performance characteristics

For information on performance characteristics, refer to the Instructions for use of the respective analyzer. For technical service, contact your local Roche Diagnostics organization.

Last update: 2013-03

FOR US CUSTOMERS ONLY: LIMITED WARRANTY

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects and will be fit for the intended use and labeling. THIS WARRANTY IS LIMITED TO THE EXTENT OF THE SPECIFICATIONS SET FORTH IN THE LABELING. THIS WARRANTY DOES NOT CONSTITUTE AN IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

Significant additions or changes are indicated by a change bar in the margin.

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