

REF		SYSTEM
03051986 190	200	Elecsys 2010 MODULAR ANALYTICS E170 <b>cobas e 411</b> <b>cobas e 601</b> <b>cobas e 602</b>

**English****Intended use**

Immunoassay for the in vitro quantitative determination of free triiodothyronine in human serum and plasma.

The electrochemiluminescence immunoassay "ECLIA" is intended for use on Elecsys and **cobas e** immunoassay analyzers.

**Summary**

Triiodothyronine is one of the thyroid hormones present in serum which regulates metabolism. Determination of this hormone concentration is important for the diagnostic differentiation of euthyroid, hyperthyroid, and hypothyroid states. The major fraction of total triiodothyronine is bound to the transport proteins (TBG, prealbumin, albumin). Free triiodothyronine (FT3) is the physiologically active form of the thyroid hormone triiodothyronine (T3).

The determination of free T3 has the advantage of being independent of changes in the concentrations and binding properties of the binding proteins; additional determination of a binding parameter (T-uptake, TBG) is therefore unnecessary.<sup>1,2,3</sup>

The sequential testing procedure and the use of a labeled antibody reduces the possibility of interference due to altered binding properties of the serum, as can occur with assays employing labeled antigen (analog method).<sup>4,5,6,7</sup>

A variety of methods are available for estimating the free thyroid hormone levels. The direct measurement of fT4 and fT3 via equilibrium dialysis or ultrafiltration is mainly used as a reference method for standardizing the immunological procedures generally used for routine diagnostic purposes.<sup>4</sup>

In the Elecsys FT3 test the determination of free triiodothyronine is made with the aid of a specific anti-T3 antibody labeled with a ruthenium complex<sup>a)</sup>.

a) Tris(2,2'-bipyridyl)ruthenium(II)-complex (Ru(bpy)<sub>3</sub><sup>2+</sup>)

**Test principle**

Competition principle. Total duration of assay: 18 minutes.

- 1st incubation: 15 µL of sample and an anti-T3-specific antibody labeled with a ruthenium complex.
- 2nd incubation: After addition of biotinylated T3 and streptavidin-coated microparticles, the still-free binding sites of the labeled antibody become occupied, with formation of an antibody-hapten complex. The entire complex is bound to the solid phase via interaction of biotin and streptavidin.
- The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell/ProCell M. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier.
- Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent barcode.

**Reagents - working solutions**

The reagent rackpack is labeled as FT3.

- M Streptavidin-coated microparticles (transparent cap), 1 bottle, 12 mL:  
Streptavidin-coated microparticles 0.72 mg/mL; preservative.

- R1 Anti-T3-Ab~Ru(bpy)<sub>3</sub><sup>2+</sup> (gray cap), 1 bottle, 18 mL:  
Monoclonal anti-T3-antibody (sheep) labeled with ruthenium complex 10 ng/mL; phosphate buffer 100 mmol/L, pH 7.0; preservative.
- R2 T3-biotin (black cap), 1 bottle, 18 mL:  
Biotinylated T3 2 ng/mL; phosphate buffer 100 mmol/L, pH 7.0; preservative.

**Precautions and warnings**

For in vitro diagnostic use.

Exercise the normal precautions required for handling all laboratory reagents.

Disposal of all waste material should be in accordance with local guidelines. Safety data sheet available for professional user on request.

Avoid foam formation in all reagents and sample types (specimens, calibrators and controls).

**Reagent handling**

The reagents in the kit have been assembled into a ready-for-use unit that cannot be separated.

All information required for correct operation is read in from the respective reagent barcodes.

**Storage and stability**

Store at 2-8 °C.

Do not freeze.

Store the Elecsys reagent kit **upright** in order to ensure complete availability of the microparticles during automatic mixing prior to use.

Stability:	
unopened at 2-8 °C	up to the stated expiration date
after opening at 2-8 °C	12 weeks
on the analyzers	6 weeks

**Specimen collection and preparation**

Only the specimens listed below were tested and found acceptable.

Undiluted serum collected using standard sampling tubes or tubes containing separating gel.

Undiluted Li-heparin, K<sub>2</sub>-EDTA and K<sub>3</sub>-EDTA plasma.

Criterion: Recovery within 90-110 % of serum value or slope 0.9-1.1 + intercept within < ± 2x analytical sensitivity (LDL) + coefficient of correlation > 0.95.

Stable for 7 days at 2-8 °C, 1 month at -20 °C.<sup>8</sup> Freeze only once.

The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.

Centrifuge samples containing precipitates before performing the assay.

Do not use heat-inactivated samples.

Do not use samples and controls stabilized with azide.

Ensure the samples, calibrators and controls are at 20-25 °C prior to measurement.



Due to possible evaporation effects, samples, calibrators and controls on the analyzers should be analyzed/measured within 2 hours.

#### Materials provided

See "Reagents – working solutions" section for reagents.

#### Materials required (but not provided)

- [REF] 03051994190, FT3 CalSet, for 4 x 1 mL
- [REF] 11731416190, PreciControl Universal, for 2 x 3 mL each of PreciControl Universal 1 and 2
- [REF] 11731416160, PreciControl Universal, for 2 x 3 mL each of PreciControl Universal 1 and 2 (for USA)
- General laboratory equipment
- Elecsys 2010, MODULAR ANALYTICS E170 or **cobas e** analyzer

Accessories for Elecsys 2010 and **cobas e** 411 analyzers:

- [REF] 11662988122, ProCell, 6 x 380 mL system buffer
- [REF] 11662970122, CleanCell, 6 x 380 mL measuring cell cleaning solution
- [REF] 11930346122, Elecsys SysWash, 1 x 500 mL washwater additive
- [REF] 11933159001, Adapter for SysClean
- [REF] 11706802001, Elecsys 2010 AssayCup, 60 x 60 reaction vessels
- [REF] 11706799001, Elecsys 2010 AssayTip, 30 x 120 pipette tips

Accessories for MODULAR ANALYTICS E170, **cobas e** 601 and **cobas e** 602 analyzers:

- [REF] 04880340190, ProCell M, 2 x 2 L system buffer
- [REF] 04880293190, CleanCell M, 2 x 2 L measuring cell cleaning solution
- [REF] 03023141001, PC/CC-Cups, 12 cups to prewarm ProCell M and CleanCell M before use
- [REF] 03005712190, ProbeWash M, 12 x 70 mL cleaning solution for run finalization and rinsing during reagent change
- [REF] 12102137001, AssayTip/AssayCup Combimagazine M, 48 magazines x 84 reaction vessels or pipette tips, waste bags
- [REF] 03023150001, WasteLiner, waste bags
- [REF] 03027651001, SysClean Adapter M

Accessories for all analyzers:

- [REF] 11298500316, Elecsys SysClean, 5 x 100 mL system cleaning solution
- [REF] 11298500160, Elecsys SysClean, 5 x 100 mL system cleaning solution (for USA)

#### Assay

For optimum performance of the assay follow the directions given in this document for the analyzer concerned. Refer to the appropriate operator's manual for analyzer-specific assay instructions.

Resuspension of the microparticles takes place automatically prior to use. Read in the test-specific parameters via the reagent barcode. If in exceptional cases the barcode cannot be read, enter the 15-digit sequence of numbers.

Bring the cooled reagents to approx. 20 °C and place on the reagent disk (20 °C) of the analyzer. Avoid foam formation. The system automatically regulates the temperature of the reagents and the opening/closing of the bottles.

#### Calibration

Traceability: The FT3 assay ([REF] 03051986) has been standardized against the FT3 assay ([REF] 11731386). This in turn was standardized using equilibrium dialysis.<sup>9</sup>

Every Elecsys reagent set has a barcoded label containing specific information for calibration of the particular reagent lot. The predefined master curve is adapted to the analyzer using the relevant CalSet.

**Calibration frequency:** Calibration must be performed once per reagent lot using fresh reagent (i.e. not more than 24 hours since the reagent kit was registered on the analyzer). Renewed calibration is recommended as follows:

- after 1 month (28 days) when using the same reagent lot
- after 7 days (when using the same reagent kit on the analyzer)
- as required: e.g. quality control findings outside the defined limits

#### Quality control

For quality control, use PreciControl Universal.

In addition, other suitable control material can be used.

Controls for the various concentration ranges should be run individually at least once every 24 hours when the test is in use, once per reagent kit, and following each calibration.

The control intervals and limits should be adapted to each laboratory's individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

#### Calculation

The analyzer automatically calculates the analyte concentration of each sample (either in pmol/L, pg/mL or ng/dL).

Conversion factors:	pmol/L x 0.651 = pg/mL
	pg/mL x 1.536 = pmol/L
	pg/mL x 0.1 = ng/dL

#### Limitations - interference

The assay is unaffected by icterus (bilirubin < 564 µmol/L or < 33 mg/dL), hemolysis (Hb < 2.7 mmol/L or < 4.3 g/dL), lipemia (Intralipid < 2000 mg/dL), biotin (< 409 nmol/L or < 100 ng/mL), IgG < 7 g/dL and IgM < 1 g/dL.

Criterion: Recovery within ± 10 % of initial value.

Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

Of 18 commonly used pharmaceuticals tested in vitro, only furosemide caused increased FT3 findings at the daily therapeutic dosage level.

Any influence that might affect the binding behavior of the binding proteins can alter the result of the FT3 tests (e.g. drugs, NTIs (Non-Thyroid-Illness) or patients suffering from FDH (Familial Dysalbuminemic Hyperthyroxinemia)).<sup>10</sup>

Autoantibodies to thyroid hormones can interfere with the assay.

Binding protein anomalies seen with FDH (familial dysalbuminemic hyperthyroxinemia), for example, may cause values which, while characteristic of the condition, deviate from the expected results.<sup>11</sup>

In rare cases, interference due to extremely high titers of antibodies to analyte-specific antibodies, streptavidin or ruthenium can occur. These effects are minimized by suitable test design.

For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.

#### Limits and ranges

##### Measuring range

0.400-50.0 pmol/L or 0.260-32.6 pg/mL (defined by the lower detection limit and the maximum of the master curve). Values below the lower detection limit are reported as < 0.400 pmol/L or < 0.260 pg/mL. Values above the measuring range are reported as > 50.0 pmol/L or > 32.6 pg/mL.



**Lower limits of measurement***Lower detection limit of the test*

Lower detection limit: 0.400 pmol/L or 0.260 pg/mL

The lower detection limit represents the lowest measurable analyte level that can be distinguished from zero. It is calculated as the value lying two standard deviations above that of the lowest standard (master calibrator, standard 1 + 2 SD, repeatability study, n = 21).

**Dilution**

Samples for FT3 determinations cannot be diluted, as T3 in the blood is present in free and protein-bound forms which are in equilibrium. A change in the concentration of the binding proteins alters this equilibrium.

**Expected values**

Establishing values in reference range studies is based primarily on samples obtained from outpatient clinics, hospitals, and commercial laboratories in which TSH and fT4 levels are found to be in the euthyroid range.

These patients often have non-thyroid diseases which might influence the thyroid function in general, and especially the fT3 level.

This may explain the differences observed when comparing the reference ranges based on different population groups using the same fT3 method. Besides local differences in iodine intake the overall health status of the individuals involved is decisive for the outcome of the reference intervals. The Elecsys FT3 assay was used to determine reference ranges in the following groups of individuals from different locations in Germany:

*Adults:*

- 3.1-6.8 pmol/L (2.0-4.4 pg/mL)

From a coastally-situated commercial laboratory in Germany, 5366 routine samples with TSH between 1 and 3 µIU/mL were evaluated by non-parametric calculation of the central 95 % limits and corresponding 95 % confidence intervals (CI) for fT3 concentration:

Median	2.5 <sup>th</sup> percentile	95 % CI of the 2.5 <sup>th</sup> percentile	97.5 <sup>th</sup> percentile	95 % CI of the 97.5 <sup>th</sup> percentile	Unit
4.6	3.1	3.07-3.19	6.8	6.65-6.87	pmol/L
3.0	2.0	2.00-2.08	4.4	4.33-4.47	pg/mL

- 3.9-6.7 pmol/L (2.5-4.3 pg/mL)

870 samples derived from apparently healthy blood donors aged 20 to 69 years from a central German site were evaluated by non-parametric calculation of the central 95 % limits and corresponding 95 % confidence intervals (CI) for fT3 concentration:

Median	2.5 <sup>th</sup> percentile	95 % CI of the 2.5 <sup>th</sup> percentile	97.5 <sup>th</sup> percentile	95 % CI of the 97.5 <sup>th</sup> percentile	Unit
5.1	3.9	3.67-3.99	6.7	6.54-7.00	pmol/L
3.3	2.5	2.39-2.60	4.3	4.26-4.56	pg/mL

The following parameters were recorded in these individuals: the concentration of TSH, fT4, and auto-antibodies to Tg and TPO; the volume and density of the thyroid gland measured by ultrasound; their case history, family, and personal thyroid history; their gender, age, and iodine intake; and whether or not they smoked or were taking oral contraceptives.

For results based on a variety of different inclusion and exclusion criteria, please refer to the separate information given in the brochure "Reference Intervals for Elecsys Thyroid Assays".

- 2.4-6.3 pmol/L (1.5-4.1 pg/mL)

211 samples from dialysis patients were measured with the Elecsys FT3 assay in a multicenter evaluation (pilot study). The data represent the 2.5<sup>th</sup> and the 97.5<sup>th</sup> percentile.

- 1.3-6.3 pmol/L (0.8-4.1 pg/mL), median 3.2 pmol/L or 2.7 pg/mL

94 samples from patient suffering from severe non thyroid illnesses (NTI's) were measured with the Elecsys FT3 assay in a multicenter evaluation (pilot study). The data represent the 2.5<sup>th</sup> and the 97.5<sup>th</sup> percentile.

*Children and adolescence:*

- Samples from newborns, infants, and adolescents up to 19 years of age, characterized as apparently healthy by experts from a medical center in central Germany:

Age	N	Median	2.5 <sup>th</sup> percentile	95 % CI of the 2.5 <sup>th</sup> percentile	Unit
4-30 days	40	5.2 3.4	3.0 2.0	2.6-3.8 1.7-2.5	pmol/L pg/mL
2-12 months	35	5.8 3.8	2.4 1.5	2.4-3.1 1.5-2.0	pmol/L pg/mL
2-6 years	73	6.1 4.0	3.0 2.0	2.9-3.6 1.9-2.4	pmol/L pg/mL
7-11 years	127	6.1 4.0	4.1 2.7	2.5-4.9 1.6-3.2	pmol/L pg/mL
12-19 years	140	5.9 3.9	3.5 2.3	3.1-3.7 2.0-2.4	pmol/L pg/mL

Age	N	Median	97.5 <sup>th</sup> percentile	95 % CI of the 97.5 <sup>th</sup> percentile	Unit
4-30 days	40	5.2 3.4	8.1 5.2	7.3-8.3 4.7-5.4	pmol/L pg/mL
2-12 months	35	5.8 3.8	9.8 6.4	8.8-9.8 5.7-6.4	pmol/L pg/mL
2-6 years	73	6.1 4.0	9.1 6.0	8.2-9.5 5.3-6.2	pmol/L pg/mL
7-11 years	127	6.1 4.0	7.9 5.2	7.6-9.2 5.0-6.0	pmol/L pg/mL
12-19 years	140	5.9 3.9	7.7 5.0	7.3-9.2 4.8-6.0	pmol/L pg/mL

The following exclusion criteria were stipulated in these individuals (both outpatient and hospitalized): no previous or acute thyroid disease, no family history of thyroid disease, no coronary disease, no intensive care, and no postoperative medical care.

Status: Adults: MCE, pilot study and β-site study Elecsys FT3 assay, [REF] 03051986, February 2004, and excerpt "MCE Reference Ranges Thyroid", June 2004. Children and adolescence: Reference range study for thyroid values in children, June 2004.

For detailed information about reference intervals in children, adolescents and pregnant women, refer to the brochure "Reference Intervals for Children and Adults", English: [REF] 04640292, German: [REF] 04625889.

This booklet also contains results of a detailed study about influencing factors on thyroid parameters in a well characterized reference group of adults. Different inclusion and exclusion criteria were applied (e.g. sonographic results (thyroid volume and density) as well as criteria according to the guidelines of the National Academy of Clinical Biochemistry - NACB).

Each laboratory should investigate the transferability of the expected values to its own patient population and if necessary determine its own reference ranges.

**Specific performance data**

Representative performance data on the analyzers are given below. Results obtained in individual laboratories may differ.

**Precision**

Precision was determined using Elecsys reagents, pooled human sera, and controls in a modified protocol (EP5-A) of the CLSI (Clinical and Laboratory



Standards Institute): 6 times daily for 10 days (n = 60). Repeatability on MODULAR ANALYTICS E170 analyzer, n = 21. The following results were obtained:

Elecsys 2010 and cobas e 411 analyzers								
Sample	Repeatability					Intermediate precision		
	Mean		SD		CV	SD		CV
	pmol/L	pg/mL	pmol/L	pg/mL	%	pmol/L	pg/mL	%
HS <sup>b)</sup> 1	2.86	1.86	0.06	0.04	2.1	0.08	0.05	2.8
HS 2	3.85	2.51	0.08	0.05	2.2	0.10	0.07	2.7
HS 3	19.5	12.7	0.29	0.19	1.5	0.36	0.23	1.9
PC U <sup>c)</sup> 1	4.98	3.24	0.08	0.05	1.7	0.11	0.07	2.2
PC U2	22.0	14.3	0.33	0.21	1.5	0.38	0.25	1.7

b) HS = human serum

c) PC U = PreciControl Universal

MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers					
Sample	Repeatability				
	Mean		SD		CV
	pmol/L	pg/mL	pmol/L	pg/mL	%
HS 1	3.06	1.99	0.06	0.04	2.0
HS 2	4.15	2.70	0.08	0.05	2.0
HS 3	20.9	13.6	0.21	0.14	1.0
PC U1	4.92	3.20	0.09	0.06	1.9
PC U2	22.0	14.3	0.29	0.19	1.3

MODULAR ANALYTICS E170, cobas e 601 and cobas e 602 analyzers					
Sample	Intermediate precision				
	Mean		SD		CV
	pmol/L	pg/mL	pmol/L	pg/mL	%
HS 1	3.02	1.97	0.10	0.07	3.4
HS 2	4.04	2.63	0.10	0.07	2.5
HS 3	20.2	13.2	0.40	0.26	2.0
PC U1	5.04	3.28	0.13	0.08	2.6
PC U2	22.2	14.5	0.47	0.31	2.1

The following are the results obtained for the NCCLS (National Committee for Clinical Laboratory Standards) precision measurements in a laboratory in France in the course of the Multicenter Study for the Elecsys FT3 performance evaluation. These data represent the performance from routine measurements on the Elecsys 2010 analyzer in clinical centers.

Elecsys 2010 and cobas e 411 analyzers				
Sample	Mean		Repeatability	Intermediate precision
	CV		CV	
	pmol/L	pg/mL	%	%
HS 1	2.83	1.84	3.3	5.1
HS 2	5.87	3.82	1.3	2.5
HS 3	15.5	10.1	2.2	3.7
PC U1	4.46	2.90	1.0	2.0
PC U2	22.8	14.8	1.5	2.1

### Method comparison

A comparison of the FT3 assay, [REF] 03051986 (y) with the FT3 assay, [REF] 11731386 (x) using clinical samples gave the following correlations (pmol/L):

Number of samples measured: 964

Passing/Bablok<sup>12</sup> Linear regression

y = 0.97x - 0.33

y = 1.0x - 0.44

τ = 0.810

r = 0.990

The sample concentrations were between approx. 0.74 and 48.5 pmol/L (approx. 0.48 and 31.6 pg/mL).

### Analytical specificity

For the antibodies used, the following cross-reactivities were found:

L-T4 0.24 %; D-T4 0.40 %; L-rT3 n.d.<sup>d)</sup>; 3,5-diiodo-L-tyrosine 0.41 %; 3,3',5'-triiodothyroacetic acid 58.2 %; 3,3',5,5'-tetra-iodothyroacetic acid 0.11 %.

d) n.d. = not detectable

### References

- 1 Wheeler MH, Lazarus JH. Diseases of the Thyroid. London, Glasgow, Weinheim, New York, Tokyo, Melbourne, Madras: Chapman and Hall Medical, 1994:107-115.
- 2 Pfannenstiel P, Saller B. Schilddrüsenerkrankungen Diagnose und Therapie. Berliner Medizinische Verlagsanstalt GmbH 1995;2:30-32,60-62.
- 3 Fisher DA. Physiological variations in thyroid hormones; physiological and pathophysiological considerations. Clin Chem 1996;42:135-139.
- 4 Klee GG. Clinical usage recommendations and analytic performance goals for total and free triiodothyronine measurements. Clin Chem 1996;42:155-159.
- 5 Surks MI, Chopra IJ, Mariash CN, et al. American Thyroid Association guidelines for use of laboratory tests in thyroid disorders. JAMA 1990;263:1529-1532.
- 6 Becker DV, Bigos ST, Gaitan E, et al. Optimal use of blood tests for assessment of thyroid function (letter). JAMA 1993;269:273.
- 7 Wild D. The Immunoassay Handbook. Stockton Press, 1994:338.
- 8 Tietz NW. Clinical Guide To Laboratory Tests. 3rd ed. Philadelphia, Pa: WB Saunders Co, 1995:612.
- 9 Method: Nichols Institute, CA, USA.
- 10 Demers LM, Spencer CA. Laboratory support for the diagnosis and monitoring of thyroid disease. National Academy of Clinical Biochemistry, 2002;Section 3B.
- 11 Wada N, Chiba H, Shimizu C, et al. A Novel Missense Mutation in Codon 218 of the Albumin Gene in a Distinct Phenotype of Familial Dysalbuminemic Hyperthyroxinemia in a Japanese Kindred. Journal of Clinical Endocrinology and Metabolism 1997;82(10):3246-3250.
- 12 Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem 1988 Nov;26(11):783-790.

For further information, please refer to the appropriate operator's manual for the analyzer concerned, the respective application sheets, the product information and the Method Sheets of all necessary components (if available in your country).


A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

### Symbols

Roche Diagnostics uses the following symbols and signs in addition to those listed in the ISO 15223-1 standard.



**FT3****FT3 - free triiodothyronine****cobas<sup>®</sup>**

<b>CONTENT</b>	Contents of kit
<b>SYSTEM</b>	Analyzers/Instruments on which reagents can be used
<b>REAGENT</b>	Reagent
<b>CALIBRATOR</b>	Calibrator
	Volume after reconstitution or mixing

**FOR US CUSTOMERS ONLY: LIMITED WARRANTY**

Roche Diagnostics warrants that this product will meet the specifications stated in the labeling when used in accordance with such labeling and will be free from defects in material and workmanship until the expiration date printed on the label. THIS LIMITED WARRANTY IS IN LIEU OF ANY OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE DIAGNOSTICS BE LIABLE FOR INCIDENTAL, INDIRECT, SPECIAL OR CONSEQUENTIAL DAMAGES.

COBAS, COBAS E, ELECSYS, MODULAR and PRECICONTROL are trademarks of Roche. INTRALIPID is a trademark of Fresenius Kabi AB.

All other product names and trademarks are the property of their respective owners.

Significant additions or changes are indicated by a change bar in the margin.

© 2013, Roche Diagnostics



Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim  
[www.roche.com](http://www.roche.com)

Distribution in USA by:  
Roche Diagnostics, Indianapolis, IN  
US Customer Technical Support 1-800-428-2336

