# Troponinas T STAT

**cobas®**
**Troponin T STAT (Short Turn Around Time), CARDIAC T (skubiai atliekamas)**
**04660307 190**
**100 tyrimų**

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cobas®

**Troponin T STAT (Short Turn Around Time), CARDIAC T (skubiai atliekamas)**

4. Ohmann EM, et al. Risk stratification with admission cardiac troponin T levels in acute myocardial ischemia. *N Engl J Med* 1996;335:1333-1334.
5. Christenson RH, Duh SH, Newby LK, Ohman EM, Califf RM, Granger CB, et al. Cardiac troponin T and cardiac troponin I: relative values in short-term risk stratification of patients with acute coronary syndromes. *Clin Chem* 1998;44(3):494-501.
6. Lindahl B, Diderholm E, Lagerqvist B, Venge P, Wallentin L, and the FRISC II investigators. Mechanisms behind the prognostic value of Troponin T in unstable coronary artery disease: a FRISC II substudy. *J Am Coll Card* 2001;38:979-986.
7. Aviles RJ, Askari AT, Lindahl B, Wallentin L, Jia G, Ohman EM, Mahaffey KW, Newby LK, Califf RM, Simoons ML, Topol EJ, Lauer MS. Troponin T levels in patients with acute coronary syndromes, with or without renal dysfunction. *N Engl J Med* 2002;346:2047-52.
8. Lindahl B, Venge P, Wallentin L. Troponin T Identifies Patients With Unstable Coronary Artery Disease Who Benefit From Long-Term Antithrombotic Protection. *J Am Coll Cardiol* 1997;29(1):43-48.
9. Hamm CW, Heeschen C, Goldmann B, Vahanian A, Adgey J, Miguel MC, et al. Benefit of abciximab in patients with refractory unstable angina in relation to serum troponin T levels. *N Engl J Med* 1999;340(21):1623-1629.
10. Heeschen C, Hamm CW, Goldmann BU, et al. for PRISM study investigators. Troponin concentrations for stratification of patients with acute coronary syndromes in relation to therapeutic efficacy of tirofiban. *Lancet* 1999;354:1757-1762.
11. Lindahl B, Diderholm E, Lagerqvist B, Venge P, Wallentin L. Effects on mortality of long-term treatment with I.m.w. heparin in relation to troponin T level and ECG findings - a FRISC 2 substudy. *Eur Heart J* 2000;21(Suppl.):521.
12. Newby LK, Ohman EM, Christenson RH, Moliterno DJ, Harrington RA, White HD, et al. Benefit of Glycoprotein IIb/IIIa Inhibition in Patients With Acute Coronary Syndromes and Troponin T-Positive Status: The PARAGON-B Troponin T Substudy. *Circulation* 2001;103:2891-2896.
13. Alpert JS and Thygesen K. Myocardial Infarction redefined - A consensus document of the joint European Society of Cardiology/American College of Cardiology committee for the redefinition of myocardial infarction. *JACC* 2000;36:959-969.
14. Braunwald E, Antman EM, Beasley JW, Califf RM, Cheitlin MD, Hochman JS, et al. ACC/AHA guidelines for the management of patients with unstable angina and non-ST-elevation myocardial infarction: Executive summary and recommendations. *Circulation* 2000;102:1193-1209.
15. Jaffe AS, Ravkilde J, Roberts R, Naslund U, Apple FS, Galvani M, and Katus H. It's time for a change to a troponin standard. *Circulation* 2000;102:1216-1220.
16. Apple FS, Wu AHB. Myocardial infarction redefined: Role of cardiac troponin Testing. *Clin Chem* 2001;47:377-379.
17. Setsuta K, Seino Y, Takahashi N, Ogawa T, Sasaki K, Harada A, et al. Clinical significance of elevated levels of cardiac troponin T in patients with chronic heart failure. *Am J Cardiol* 1999;84:608-611.
18. Sato Y, Kataoka K, Matsumori A, Sasayama S, Yamada T, Ito H, Takatsu Y. Measuring serum aminoterminal type III procollagen peptide, 7S domain of type IV collagen, and cardiac troponin T in patients with idiopathic dilated cardiomyopathy and secondary cardiomyopathy. *Heart* 1997;78:505-508.
19. Lauer B, Niederau C, Kühl U, Schannwell M, Pauschinger M, Strauer BE, Schultheiss HP. Cardiac troponin T in patients with clinically suspected myocarditis. *JACC* 1997;30:1354-1359.
20. Swaanenburg JCJM, Klaase JM, DeJongste MJL, Zimmerman KW, ten Duis HJ. Troponin I, troponin T, CK-MB-activity and CK-MB mass as markers for the detection of myocardial contusion in patients who experienced blunt trauma. *Clin Chim Acta* 1998;272:171-181.
21. Dierkes J, Domrose U, Westphal S, Ambrosch A, Bosselmann HP, Neumann KH, Luley C. Cardiac troponin T predicts mortality in patients with renal disease. *Circulation* 2000;102:1964-1969.
22. Giannitsis E, Müller-Bardorff M, Kurowski V, Weidtmann B, Wiegand U, Kampmann M, Katus HA. Independent prognostic value of cardiac troponin T in patients with confirmed pulmonary embolism. *Circulation* 2000;102:211-217.
23. James P, Ellis CJ, Whitlock RML, McNeil AR, Henley J, Anderson NE. Relation between troponin T concentration and mortality in patients presenting with an acute stroke: observational study. *BMJ* 320:1502-1504.
24. ver Elst KM, Spapen HD, Nguyen DN, Garbar C, Huyghens LP, Gorus FK. Cardiac troponin I and T are biological markers of left ventricular dysfunction in septic shock. *Clin Chem* 2000;46:650-657.
25. Jenkins DP, Pugsley WB, Alkhulaifi AM, Kemp M, Hooper J, Yellon DM. Ischaemic preconditioning in patients undergoing coronary artery bypass surgery. *Heart* 1997;77:314-318.
26. Lopez-Jimenez F, Goldman L, Sacks DB, Thomas EJ, Johnson PA, Cook EF, Lee TH. Prognostic value of cardiac troponin T after non-cardiac surgery: 6 month follow up data. *JACC* 1997;29:1241-1245.
27. Johansen O, Brekke M, Stromme JH, Valen V, Seljeflot I, Skjaeggstad O, Arensen H. Myocardial damage during percutaneous transluminal coronary angioplasty as evidenced by troponin T measurements. *Eur Heart J* 1998;19:112-117.
28. Herman EH, Lipshultz SE, Rifai N, Zhang J, Papoian T, Yu ZX, et al. Use of cardiac troponin T levels as an indicator of doxorubicin-induced cardiotoxicity. *Cancer Res* 1998;58:195-197.
29. Ooi DS, Zimmerman D, Graham J, Wells GA. Cardiac troponin T predicts long-term outcomes in hemodialysis patients. *Clin Chem* 2001;47:412-417.
30. Apple FS, Murakami MM, Pearce LA, Herzog CA. Predictive value of cardiac troponin I and T for subsequent death in end-stage renal disease. *Circulation* 2002;106:2941-2945.

# Troponinas T STAT

**Troponin T STAT (Short Turn Around Time), CARDIAC T (skubiai atliekamas)**

31. deFilippi C, Wasserman S, Rosanio S, Tiblier W, Sperger H, Tocchi M, Christenson R, Uretsky B, Smiley M, Gold J, Muniz H, Badalamenti J, Herzog C, Henrich W. Cardiac Troponin T and C-reactive protein for predicting prognosis, coronary atherosclerosis, and cardiomyopathy in patients undergoing long term hemodialysis. *JAMA* 2003;290:353-359.
32. Iliou MC, Fumeron C, Benoit MO, Tuppin P, Calonge VM, Moatti N, Buisson C, Jacquot C. Prognostic value of cardiac Markers in ESRD: Chronic Hemodialysis and New Cardiac markers Evaluation (CHANCE) study. *Am J Kidney Dis* 2003;42:513-523.
33. Choy JB, Armstrong PW, Ulan RA, Campbell PM, Gourishankar S, Prosser CI, Tymchak W. Do cardiac troponins provide prognostic insight in hemodialysis patients? *Can J Cardiol* 2003;19:907-911.
34. Scott B, Deman A, Peeters P, Van den Branden C, Stolear JC, Van Camp G, Verbeelen D. Cardiac Troponin T and malondialdehyde modified plasma lipids in haemodialysis patients. *Nephrol Dial Transplant* 2003;18:737-742.
35. Davis GK, Labugger R, Van Eyk JE, Apple FS. Cardiac troponin T is not detected in western blots of diseased renal tissue. *Clin Chem* 2001;47:782-783.
36. Ricchiuti V, Voss EM, Ney A, Odland M, Anderson PAW, Apple FS. Cardiac Troponin T isoforms expressed in renal diseased skeletal muscle will not cause false positive results by the second generation cTnT assay by Boehringer Mannheim; *Clin Chem* 1998;44(9):1919-1924.
37. Wu AHB, Feng Y-J, Moore R, et al. Characterisation of cardiac troponin subunit release into serum after acute myocardial infarction and comparison of assays for troponin T and I. *Clin Chem* 1998;44:1198-1208.
38. Hallermayer K, Klenner D and Vogel R. Use of recombinant human cardiac troponin T for standardization of third generation troponin T methods. *Scand J Clin Invest* 1999;59(Suppl 230):128-131.
39. Müller-Bardorff M, et al. Improved troponin T ELISA specific for cardiac troponin T isoform: assay development and analytical validation. *Clin Chem* 1997;43(3):458-466.
40. Kampmann M, Rauscher T, Müller-Bardorff M, Münch M, Klein G, Katus H. Clinical Evaluation of Troponin T and CK-MB Mass on the Elecsys 2010 analyzer. Poster presented at Medlab 97, 12th IFCC European Congress of Clinical Chemistry, August 17-22, 1997, Basel, Switzerland.
41. Kobayashii S, Tanaka M, Tamura N, Hashimoto H, Hirose S-I. Serum cardiac troponin T in polymyositis/dermatomyositis. *Lancet* 1992;340(8821):726.
42. Mair J, Dienstl F, Puschendorf B. Cardiac Troponin T in the Diagnosis of Myocardial Injury. *Critical Reviews in Clin Lab Sci* 1992;29(1):31-57.
43. Carrier M, Solymoss BC, Raymond C, Leclerc Y, Pelletier LC. Cardiac Troponin T and Creatine Kinase MB Isoenzyme as Biochemical Markers of Ischemia after Heart Preservation and Transplantation. *J Heart Lung Transplant* 1994;13(4):696-700.
44. Löfberg M, Tähtelä R, Härkönen M, Somer H. Myosin Heavy-Chain Fragments and Cardiac Troponins in the Serum in Rhabdomyolysis. Diagnostic Specificity of New Biochemical Markers. *Arch Neurol* 1995;52:1210-1214.
45. Anderson JR, Hossein-Nia M, Brown P, Holt DW, Murday A. Donor Cardiac Troponin T Predicts Subsequent Inotropic Requirements Following Cardiac Transplantation. *Brief Communications in Transplant* 1994;58(9):1056-1057.
46. Franz WM, Remppis A, Kandolf R, Kübler W, Katus HA. Serum Troponin T: Diagnostic Marker for Acute Myocarditis. Letter to the Editor in *Clin Chem* 1996;42(2):340-341.
47. Mair P, Mair J, Koller J, Wieser C, Artner-Dworzak E, Puschendorf B. Cardiac troponin T in the diagnosis of heart contusion. *Lancet* 1991;338:693.
48. Bablok W, et al. A General Regression Procedure for Method Transformation. *J Clin Chem Clin Biochem* 1988;26:783-790.

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©2007 Roche Diagnostics
Roche Diagnostics GmbH, D-68298 Mannheim