

## ACTIVE CLEANING EFFECTIVE DISINFECTION

### Sekusept™ aktiv and Sekusept™ aktiv test strips

#### PERACETIC ACID BASED DISINFECTANT WITH EXCELLENT CLEANING PERFORMANCE

Sekusept aktiv combines active cleaning with effective disinfection and excellent material compatibility in one product to offer the best features for the user, the patient and your instrument.

Sekusept aktiv is also:

- Fast virucidal terminal disinfection
- Actively removes protein stains
- No fixing effect



# Sekusept™ aktiv

Peracetic acid based disinfectant with excellent cleaning performance

## PRODUCT CHARACTERISTICS

- White powder with blue granules
- pH 2% solution tap water at 20°C: approx 8.5 +/- 0.5
- pH 1% solution tap water at 20°C: approx 8.5 +/- 0.5

## MICROBIOLOGICAL PROPERTIES

CLAIM/CONDITIONS	CONCENTRATION	CONTACT TIME
Bactericidal according to EN 13727/EN 14561 Clean & dirty conditions	1%	5 min
Yeasticidal according to EN 13624/EN 14562 Clean & dirty conditions	1%	5 min** 10 min*
Fungicidal according to EN 13624/EN 14562 Clean conditions	2%	30 min
Tuberculocidal according to EN 14348/EN 14563 Dirty conditions	1%	10 min
Mycobactericidal according to EN 14348/EN 14563 Clean conditions	2%	15 min
Limited virucidal according to EN 14476/EN 17111 Clean & dirty conditions	1%	10 min
Virucidal according to EN 14476/EN 17111 Clean conditions	2%	15 min
Sporicidal according to EN 17126 Clean conditions	2%	30 min
Sporicidal against C. difficile according to EN 17126 Clean conditions	1%	15 min
According to VAH (bactericidal, yeastidical and tuberculocidal) Clean & dirty conditions	1%	60 min*
	2%	5 min **
According to VAH (mycobactericidal) Clean & dirty conditions	1%	60 min*
	2%	15 min**

\* = Cleaning-disinfection step (dirty conditions)  
\*\*= Disinfection step (clean conditions)

## PACKAGING AND ORDERING INFORMATION

PRODUCT	PACKAGING UNIT	ORDER CODE
Sekusept aktiv bucket with measuring spoon	4 x 1.5 kg	XXXXXX
Sekusept aktiv bucket with measuring cup	6 kg	XXXXXX
Sekusept aktiv test strips	4 tubes of 50 strips	XXXXXX



## APPLICATION AREA

For the manual cleaning and disinfection as well as the terminal disinfection of medical devices.

## INSTRUCTIONS FOR USE

For the preparation of the application solution use enclosed measuring spoon/ measuring cup.

Add the required powder amount to reach the demanded efficacy into the water (<30°C) during stirring and repeat several times during the first 15 minutes. The active ingredient peracetic acid is formed when the solution is made up. Immerse instruments ensuring all surfaces are completely wet.

Clean items acc. to instrument manufacturer's instruction, including mechanical action (e.g.: brushing, wiping) until visual clean. Remove instruments after the required contact time and rinse thoroughly with water of at least drinking water quality. Dry external surfaces and internal channels (e.g. of endoscopes) afterwards.

For the terminal disinfection, prepare one bath for the cleaning-disinfection step and one bath for the disinfection step.

The residual efficacy of the solution, e.g. in case of increased temperature can be checked using Sekusept™ aktiv test-strips. Not suitable for instruments made of brass or copper or for mechanically predamaged chrome-plated or nickel-plated instruments.

The solution has to be prepared fresh for each working day. In case of visible contamination prepare a fresh solution immediately.

Slight discoloration of powder is possible, this has no influence to the product performance.

## INGREDIENTS

Sodium percarbonate, non-ionic tensides, phosphonates. Active contents of a 2% solution for actual use: > 1,000 ppm peracetic acid.

## STORAGE

Dangerous. Respect the precautions for use (Drawn up according to the European rules in force regarding the classification and labelling of chemical products). Storage: between +0°C and +25°C.

## CERTIFICATES

Medical device. Full antimicrobial activity for Sekusept aktiv is included in the scientific dossier, available on request.

## MEDICAL DEVICE USE SAFELY

Please read label and product information.




FOR PROFESSIONAL USE ONLY.

# Sekusept™ aktiv test strips

## INDICATIONS

Sekusept aktiv 50 test strips.

## INSTRUCTIONS FOR USE



Immerse the reaction zone in the solution for 2 seconds.

Shake off excess liquid from the strip.

Wait 30 seconds, compare the reaction zone with the respective colour field.

Interpretation - compare the colour of the reaction zone with the selected concentration:  
Same or more intense = fully active  
Lighter or less intense = replace solution

## STORAGE CONDITIONS

Close the tube immediately after having removed the required test strips. Store in a cool dry place. Keep away from sunlight. Protect contents against moisture. For single use only.


## PRECAUTIONS FOR USE

For professional use only.


INSTRUMENT  
REPROCESSING

## READING AFTER 30 SECONDS:

Sekusept aktiv 1%



Sekusept aktiv 2%



FOR PROFESSIONAL USE ONLY.

---

**ECOLAB EUROPE**

Richtistrasse 7  
8304 Wallisellen  
Switzerland



**ECOLAB DEUTSCHLAND GmbH**

Ecolab-Allee 1  
40789 Monheim am Rhein  
Germany



[www.ecolab.com](http://www.ecolab.com)