

AMPLATZER™ VALVULAR PLUG III (FORMERLY AVP III)

CLOSES LEAKS OPENS OPPORTUNITIES

ONE-CLASS IMPROVEMENT IN
THE NYHA CLASSIFICATION IN
UP TO 90% OF PATIENTS^{1,2,3,4}

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.



DELIVERING AN EFFECTIVE SOLUTION FOR A WIDE RANGE OF PVL MORPHOLOGIES^{1,5,6}

INTRODUCING THE AMPLATZER VALVULAR PLUG III

Paravalvular leaks can be found in a variety of shapes and sizes. To tackle the leaks and fill the gaps of different morphologies, the Amplatzer Valvular Plug III is designed in a rectangular oval shape that fits a wide range of paravalvular leaks (PVL). As such, the Amplatzer Valvular Plug III can provide effective closure in over 90% of patients at 30 day follow up.^{6,7,8}

A SOLUTION TO A KEY ISSUE

Around the world, paravalvular leaks are a common and challenging problem: paravalvular leaks occur in 7% to 17% of mitral valve replacements (MVRs) and 5% to 10% of aortic valve replacements (AVRs).⁹ By providing an effective solution to this key issue, the Amplatzer Valvular Plug III is improving life quality and longevity to an increasing number of patients.^{1,2,3,4}

BUILT ON THE EXTENSIVE AMPLATZER™ LEGACY OF SAFETY AND EFFICACY

- Pioneered transcatheter cardiovascular and peripheral vascular occlusion
- Over 1.25 million Amplatzer devices implanted worldwide¹¹
- More than 20 years of clinical experience and global leadership

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

CLINICALLY PROVEN OUTCOMES

Clinical studies continue to show that the Amplatzer Valvular Plug III, formerly known as AVP III, is a highly effective solution in closing PVLs near mechanical surgical valves.¹



≥90% OF PVLs EFFECTIVELY CLOSED AT 30 DAY FOLLOW-UP^{6,7,8}



UP TO 90% PATIENTS REPORT ONE-CLASS NYHA CLASSIFICATION IMPROVEMENT^{1,2,3,4}



A SOLUTION THAT HAS RECEIVED CLOSE ATTENTION -
AND STRONGLY POSITIVE REVIEWS.¹¹

NR: Data not reported in the article. NS: Not significant. n/a: Not applicable. Range instead of total was provided for follow-up > 30d as the follow up durations differ greatly between studies.

*In Smolka 2016, small PVLs causing significant hemolysis but no heart failure symptoms was an exclusion criteria; therefore patients with hemolytic anemia only were not included. Total percentages may underrepresent the actual population with hemolytic anemia only and overrepresent the population with heart failure only.

OVERALL – Mitral and Aortic		Cruz-Gonzalez 2014	Smolka 2016		Davidavicius 2014	Swaans 2012	Werner 2018
Number of patients (number of PVLs)		33 (34)	49 (49)		7 (9)	7 (7)	10 (17)
Mitral / Aortic PVLs		27 / 7	29 / 20		9 / 0	6 / 1	12 / 5
Mechanical		32	30		4	4	6
Indication for PVL (% patients (N))	Heart failure	21.2% (7)	89.8% (44)		57.1% (4)	14.2% (1)	50% (5)
	Hemolytic anemia	3% (1)	0%*		0%	42.9% (3)	0%
	Both	75.7% (25)	10.2% (5)		42.9% (3)	42.9% (3)	50% (5)
Access (TA: transapical, TS: transseptal, TF: transfemoral)		Mitral: TF, TS Aortic: TF	Mitral: TS, TA Aortic: TF		Mitral:TA Aortic: n/a	Mitral:TA Aortic: TA	Mitral: TS, TA Aortic: TF
Follow-up		90 d	6 mo, 1 y		40 – 364 d	3 mo	1 y
COMPOSITE ENDPOINT							
Technical success (overall)		90.9%	93.9%		100%	100%	86%
Mitral		92.3%	89.7%		100%	100%	NR
Aortic		100%	100%		n/a	100%	NR
SAFETY ENDPOINTS							
Survival	Intra-procedural	100%	100%		100%	100%	100%
	30 d	100%	98%		100%	100%	80%
	Follow-up	100%	95.9%		85.7%	85.7%	70%
Stroke (30 d)		0	2%		NR	0	0
Conversion to surgery		6%	NR		0	14.2%	0
Bleeding / vascular complications		12%	2%		28.6%	14.2%	20%
EFFECTIVENESS ENDPOINTS							
% patients with absent to moderate PVL regurgitation		100%	93.8%		100%	100%	NR
% patients with reduction in NYHA class		90.3%	90.5%		85.7%	71.4%	NR
Hemolysis		Improved	Improved		Improved	NS	NR

NOTE: Results from clinical studies are not directly comparable.
Information provided for educational purposes only.

EVERY DETAIL IS DESIGNED FOR SUCCESSFUL PVL CLOSURE

ASYMMETRICAL OBLONG SHAPE


In this shape, the disc can fit a large range range of paravalvular leak morphologies.

DENSE NITINOL WIRE LAYER DESIGN

Additional layers of dense Nitinol wire facilitate rapid occlusion


SMALL RETENTION FEATURES

Designed to minimize valve leaflet interference



DESIGN SPECIFICATIONS

The Amplatzer™ Valvular Plug III comes in nine different sizes in order to facilitate optimal PVL closure for a range of gaps.¹⁰



Model Number	Type Amplatzer Valvular Plug III	A Waist Long Axis (mm)	B Waist Short Axis (mm)	C Overhang (mm)	D Plug Length (+1.0 mm / - 0.5 mm)
9-APVL3-042	4x2 mm	4	2	2	6.5
9-APVL3-063	6x3 mm	6	3	2	6.5
9-APVL3-084	8x4 mm	8	4	2	6.5
9-APVL3-103	10x3 mm	10	3	2	6.5
9-APVL3-105	10x5 mm	10	5	2	6.5
9-APVL3-123	12x3 mm	12	3	2	6.5
9-APVL3-125	12x5 mm	12	5	2	6.5
9-APVL3-143	14x3 mm	14	3	2	6.5
9-APVL3-145	14x5 mm	14	5	2	6.5

Expert support at every turn


CLINICAL CASE SUPPORT

- Experienced field personnel
- Over two decades of excellence

CLINICAL TRAINING PROGRAMS

- Training centers and online courses
- Fellows programs

For more information about the Amplatzer™ Valvular Plug III,
contact your Abbott sales representative



DOWNLOAD NOW

Download the free Amplatzer™ Portfolio App from your favorite app store.



REFERENCES:

1. Calvert et al. Percutaneous device closure of paravalvular leak. *Circulation*. 2016;134(13):934-944.
2. Yildirim A, Goktekin O, Gorgulu S, et al. A new specific device in transcatheter prosthetic paravalvular leak closure: a prospective two-center trial. *Catheter Cardiovasc Interv*. 2016;88(4):618-624.
3. Smolka G, Pysz P, Kozłowski M, et al. Transcatheter closure of paravalvular leaks using a paravalvular leak device - a prospective Polish registry. *Postepy Kardiol Interwencyjnej*. 2016;12(2):128-134.
4. Angulo-Llanos R, Sarnago-Cebada F, Rivera AR, et al. Two-Year Follow Up After Surgical Versus Percutaneous Paravalvular Leak Closure: A Non-Randomized Analysis. *Catheterization and cardiovascular interventions : official journal of the Society for Cardiac Angiography & Interventions*. 2016;88(4):626-634.
5. Garcia et al. Outcomes and predictors of success and complications for paravalvular leak closure: an analysis of the SpanisH real-wOrld paravalvular LEaks closure (HOLE) registry. *EuroIntervention*. 2017;12(16):1962-1968.
6. Davidavicius G, Rucinskas K, Drasutiene A, et al. Hybrid approach for transcatheter paravalvular leak closure of mitral prosthesis in high-risk patients through transapical access. *J Thorac Cardiovasc Surg*. 014;148(5):1965-1969.
7. Smolka G, Pysz P, Jasi ski M, et al. Multiplug paravalvular leak closure using Amplatzer Vascular Plugs III: A prospective registry. *Catheter Cardiovasc Interv*. 2016;87(3):478-487.
8. Werner N, Zeymer U, Fraiture B, et al. Interventional treatment of paravalvular regurgitation by plug implantation following prosthetic valve replacement: a single-center experience. *Clin Res Cardiol*. 2018.
9. Cruz-Gonzalez I, Rama-Merchan JC, Rodriguez-Collado J, et al. Transcatheter closure of paravalvular leaks: state of the art. *Neth Heart J*. 2017;25(2):116-124.
10. Nietlispach F, Eckstein F, Seeberger M, Osswald S, Kaufmann BA, Reuthebuch O. Closure of apical access site after transapical, transcatheter paravalvular leak closure. *Can J Cardiol*. 2012;28:516.e5-7.
11. Data on File at Abbott

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION outside of the U.S. ONLY**. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

Abbott

Park Lake, Culliganlaan 2b, 1831 Diegem, Belgium
www.cardiovascular.abbott

™ Indicates a trademark of the Abbott group of companies.

© 2020 Abbott. All Rights Reserved. MAT-2002715 v1.0 | Item approved for OUS use only.

