
NEW AUTOMATED PERIMETERS


Fast and precise perimetry at your fingertips

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PTS 920 | **PTS 925W** | PTS 2000

FAST AND PRECISE PERIMETRY AT YOUR FINGERTIPS

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The PTS automated perimeter series are modern diagnostic instruments for precise and fast testing of visual fields. Beside standard perimetry testing techniques they implement all of the latest visual field testing approaches saving time and deliver cutting edge precision and reliability. Intuitive software with powerful networking and EMR integration, combined with maintenance free and reliable hardware creates the best tools for visual field testing.

EYE TRACKING

Analyses pupil movement and assures supreme reliability of test results.

BLINK CONTROL

No stimuli are omitted due to normal blinking, test is put on hold until the eye is opened again.


VOICE GUIDE

Voice messages assist operator and patient during the examination.

THRESHOLD IN 3 MINUTES WITH ZETA™ **

ZETA™ and ZETA™ Fast strategies - fast and reliable examination with precise threshold estimates.

EYEESEE™ EYE RECORDING

EyeSee™ module records eye preview images during stimulus exposures and displays them when reviewing the test result.


DPA™ ANALYSIS

DPA™ Defect Progression Analysis - inter-test analysis and judgement on defect progression.

HEAD TRACKING

Monitoring and adjusting patient's head position throughout the test.

ALL THAT YOU NEED

- ✓ Complete set of standard perimetry tests.
- ✓ Variety of threshold algorithms and screening techniques.
- ✓ Binocular Driver's Test, Flicker perimetry and Blue on Yellow*.
- ✓ Complete palette of analysis tools for quick and thorough assessment of visual fields.
- ✓ Evaluation of peripheral field loss.
- ✓ Standard test reports and global indices.
- ✓ GHT and sectors analysis.

PTS 925W

Automated Perimeter PTS 925W brings the world's most popular orthogonal fields 30-2 and 24-2 into the compact sized device.

Despite its compact size, the 925W offers cutting-edge design and software capabilities. With the 170° horizontal and 110° vertical testing range (with fixation shifts), there is no scenario in which the PTS 925W could not play a role in glaucoma and other disease management. From precise macula thresholding to the binocular Esterman driving test – everything is within reach.

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- Compact and robust design
- White on white testing
- Goldmann standard bowl and stimulus
- 170° horizontal and 110° vertical testing range
- 30-2, 24-2 and 10-2 testing fields.
- Ptois testing capabilities (Sup 44 test field)
- Thresholding on 24-2 and 30-2 in 3-4 minutes**
- Standard style printouts

PTS 920/BY

Automated Perimeter PTS 920/BY is designed to fulfill all requirements of contemporary perimetry.

Being a continuation of a successful PTS 900 perimeter series, it offers the same ease of use and affordability. With ergonomic design and modern testing and analysis software, the 920 makes a leap forward from the original PTS 900 series. Full sized testing bowl and user-oriented software guarantee the perfect results with no effort.

- 160° horizontal and 100° vertical testing range
- Goldmann standard bowl and stimulus
- Thresholding on 24A and 30A in 3-4 minutes**
- Full featured software in standard
- Statistical package and Dicom interface
- Optional Blue on Yellow testing
- Standard style printouts
- Maintenance free operation

PTS 2000

Automated Perimeter PTS 2000 sets a new standard for full-featured projection perimeters.

This is where the ergonomic design meets the requirements of modern visual field evaluation. Performing everything from a quick screening test to manually controlled kinetic perimetry is simpler than ever before. PTS 2000 brings to you one of the world's largest field testing ranges and options, together with all Goldmann stimuli sizes and stimuli colors used in perimetry. With the modern and intuitive software, all that is available at your fingertips.

- Complete set of static perimetry strategies
- Automatic and manual kinetic perimetry
- 180° horizontal and 130° vertical testing range for peripheral screening
- Aspherical bowl for compact dimensions
- All Goldmann stimuli sizes and colors used in perimetry
- 30-2, 24-2 and 10-2 testing fields
- Ptois testing capabilities (Sup 64 test field)
- Standard style printouts

PTS – AUTOMATED PERIMETER SERIES

THRESHOLD IN 3 MINUTES**

- ZETA™ and ZETA™ Fast strategies - testing based on statistical data for high reliability and precise threshold estimation.
- Advanced Threshold strategy - quicker examination without loss in result resolution.
- 24-2/24A field tested in less than 3 minutes (30-2/30A in less than 4 minutes).**
- Provides more detailed information than screening strategies in comparable duration.

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- High resolution camera for precise monitoring of patients' compliance.
- **Blink control - no stimuli are omitted due to normal blinking, test is put on hold until the eye is opened again.**
- Eye tracking - analyses pupil movement and assures supreme reliability of test results.
- Voice messages assist operator and patient during the examination.
- EyeSee eye recording for even better evaluation of test reliability.
- Head Tracking mechanism to assure appropriate eye position and avoid test artifacts.

STANDARD TEST REPORTS AND GLOBAL INDICES

FLEXIBILITY FOR TIME

- Rich library of predefined test fields and programs.
- Custom test settings can be stored and accessed with one mouse click.
- Follow-up tests can be performed by simple selection of base examination.
- Partial test results of not completed test can be continued at a later date.

DICOM, GDT, EMR, NETWORK INTEGRATION

- **Built-in networking features - data can be accessed from anywhere within your network.**
- Out of the box DICOM support.
- Storing visual field reports on the DICOM image server.
- **DICOM Modality Worklist support.**
- GDT interface and custom EMR's based on text file interface.

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STRUCTURE & FUNCTION ANALYSIS

- Export of visual field data to Optopol SOCT software.
- Invaluable combination of information about the functional quality of vision with comprehensive data on retinal structures.
- S+F provides a comprehensive single page report for glaucoma management.

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
Figure 190. Basic display and the color-coded map of the RAW analysis

The color-coded map is interpolated. The space between the tested points is filled with symbols representing values estimated from the linear interpolation of the tested points. This way it is easier to read the map and to estimate the areas of visual defects.

The interpolated results in the tested field areas are based on mathematical analysis and are not a diagnosis. They represent a probable condition of the test field between the tested points. The actual condition of the visual field can be only determined through testing procedures.

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14.3.5. Interpolation of results to maps 30-2, 24-2 and 10-2

The software allows to present the results of selected examinations in HFA-device style on maps 30-2, 24-2 and 10-2. The function is enabled by using "HFA Interpolation" button (Figure 188. "HFA Interpolation" button).

After enabling, the original examination result will be interpolated. As a result, a continuous retinal sensitivity map is generated, from which the presumptive sensitivity values at points that were not physically tested can be read. Locations belonging to commonly used test areas are used for displaying: 10-2, 24-2 and 30-2; these are popular for HFA devices. Whether the original result of examination is interpolated to map 10-2, 24-2 or 30-2 is dependent on the field range from the original result.

The results interpolation function is only available for tests in which numerical levels of retinal sensitivity may be determined. These include examinations performed with different threshold strategies (Threshold, Fast Threshold, Screening and Advanced). For interpolation to be possible with the selected examination strategy, it must be completed in at least 75%, and the examination must be performed on a sufficiently large field of vision.


Figure 50. Source test preview for follow up

After the settings have been uploaded from the source test, the test parameters can be modified as in any other standard test. You can change the test field value, for example, and leave other parameters unchanged.

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The position of the tested eye should be properly adjusted. It should be placed exactly opposite the fixation point shown in the video image from the eye camera.

To correctly position the eye, adjust the height of the chin rest using the indicators displayed in the camera window.

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Last updated on the 12th of October 2023

Column

CPU Name ▾	CPU Mark (higher is better) ▾	Rank (lower is better) ▾	CPU Value (higher is better) ▾	Price (USD) ▾
Intel Core i5-1335U	17,768	557	52.26	\$340.00*

* - Last price seen from our affiliates.

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Product specifications

HP Data Sheet	HP LaserJet Pro M404 Series printer
Functions	Print
Print speed	Up to 40 ppm (default) ; Up to 42 ppm (HP High Speed)
First page out (ready) black	As fast as 6.1 sec
Resolution (black)	Fine Lines (1200 x 1200 dpi) ← 2.18 p. Skiriamoji geba (raiška) 1200x1200 dpi
Resolution technology	HP FastRes1200, HP ProRes1200, Economode
Monthly duty cycle	Up to 80,000 pages ^[10]
Recommended monthly page volume	750 to 4000 ^[11]
Print Technology	Laser
Display	2-line backlit LCD graphic display
Processor speed	1200 MHz
Number of print cartridges	1 (black) 2.18 p. Nespalvotas

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AC 200 V - 240 V ~50 Hz

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H-K – Heijl–Kraukau blind spot monitoring

SF – short-time fluctuation in the test results

FPOZ – monitoring false positive results (reaction in absence of stimuli) during the test

FNEG – monitoring false negative results (no reaction to stimuli) during the test

Check / uncheck the box next to the test name to enable / disable the selected type of reliability tests to be run during the visual field test.

9.2.5. Setting lens correction for the test

It is important for the patient to clearly see all stimuli during the visual field test. Absence of reaction to stimuli should be attributed solely to eye diseases other than problems with eye refraction. Therefore correction lenses should be used for patients with refractive problems on the test results.

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Test lens of 38 mm in diameter are used for lens correction, placed in a holder in front of the patient. To get reliable test results, make sure that the correction lens option is enabled. The test will be divided into two phases:

- internal diameter of the visual field –within a 30° ring, tested with correction lens
- external diameter of the visual field – outside a 30° ring, tested without correction lens

This division is introduced to eliminate artifacts generated when some of the stimuli is covered by the lens holder and the correction lens frame.

To activate the correction lens option, enter a non-zero value in the correction lens settings:


Figure 40. Settings of correction lens for the test

Correction values can be entered manually or with one of supporting tools. The available lens correction range is between -18D and 18D in 0.25D steps. Cylinder axis can vary within the range of 0° and 180°.

Recalling the recently used correction lens for the tested eye – this function browses all previous tests belonging to current patient. If the selected eye has been already tested, the system will recall the correction lens used in the previous test.

Lens calculator – it calculates the proposed correction lens based on the patient's current clinical refraction and age. The following calculation formulas are used:

- Cylindrical refractive errors $\leq \pm 0.25D$ are ignored.

Background	White (green stimuli) – 10 asb (3.18 cd/m ²) – automatic stabilization (PTS 925 and PTS925 BY only); chromacity (X,Y)=(0.374..0.473;0.366..0.432) Yellow (blue stimuli) – 314 asb (100 cd/m ²) – automatic stabilization (PTS 925 BY only); dominant wavelength (585..595) nm White (white stimuli) – 31.4 asb (10 cd/m ²) – automatic stabilization (PTS 925 W only); chromacity (X,Y)=(0.374..0.473;0.366..0.432)					
Stimuli exposure time	Adjustable within the range of 0.1s to 9.9s or adaptable					
Reaction time	Adjustable within the range of 0.1s to 9.9s or adaptable					
Intervals between exposures	Adjustable within the range of 0.1s to 9.9s or adaptable					
Fixations	Central, single (0°) Central, 4 points Top (15°) Left (35°) Right (35°) Bottom (15°) ¹					
Redefined examination fields	Designation	Nasal range [°]	Temple range [°]	Top range [°]	Bottom range [°]	Number of points
	F50-2 (Full)	50	50	40	40	168
	G50-2 (Glaucoma)	50	27	27	27	104
	24-2	27	21	21	21	54
	30-2	27	27	27	27	76
	5-2	3	3	3	3	16
	10-2	9	9	9	9	68
	P50-2 (Peripheral)	33-50	27-40	27-40	33-50	76
Esterman M	50	75	34	55	100	

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Table 3. Technical data of the PTS 2000 series perimeter

Simulator Bowl Type	Aspherical, according to Goldmann's standard, radius 30 cm
Field range with fixation offset	Left-right: 180° Up 70° Down 70°

²PTS925 Wi only


Figure 53. Test programs

Test Programs includes pre-set values of the following parameters:

- test strategy
- test field
- reliability tests
- central fixation target
- test field modification options
- bracketing of threshold strategies
- determination or pre-set value of the initial stimuli level
- stimulus color
- stimulus time parameters
- BDT sensitivity

After the test program has been loaded and before the test begins, the settings can be modified and adapted on a case-by-case basis – the settings of the stored test program will not be modified.

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