VIASURE

Extraction Kit

Blood Pathogens

CE IVD

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EN For download IFUS from other languages, please enter in **certest.es/viasure/labeling**. Once you be there, please following the instructions for access to the language that you need. If you need additional information, please contact: viasure@certest.es.

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Note: The user should notify the manufacturer and the competent authority of the Member State in which he is established as a user and/or patient of any serious incident related to the product.

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Content

1. Intended purpose	5
2. Safety information	5
3. Summary and Explanation	6
4. Principle of the procedure	7
5. Reagents provided	8
6. Reagents and equipment to be supplied by the user	9
7. Transport and storage conditions	10
8. Warnings and Precautions	11
9. Test procedure	13
9.1. Sample collection, storage and transport	13
9.2. Sample pretreatment	14
9.3. Reagents' resuspension / preparation	14
9.3.1 Buffer's preparation	15
9.3.2 Proteinase K resuspension	15
9.3.3 Carrier resuspension	15
9.3.4 Lyophilized extraction control	15
9.3.5 Magnetics Beads	16
9.4. Extraction procedure	16
9.5 Real-Time PCR	19
10. Result interpretation	19
11. Limitations of the procedure	20
12. Quality control	20
13. Analytical performance characteristics	20
13.1 Analytical linearity	21
13.2 Analytical sensitivity. Limit of Detection (LoD)	21
13.3 Accuracy	22
13.3.1 Trueness (Veracity)	22
13.3.2 Precision	22
13.4 Analytical specificity	25
13.5 Metrological traceability	27
14. Clinical Performance characteristic	27
15. Troubleshooting	29

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ENGLISH

1. Intended purpose

The VIASURE Blood Pathogens Extraction Kit is designed for the isolation and purification of viral, bacterial and parasitic DNA and RNA from different clinical samples (whole blood, plasma, serum, blood culture, and cerebrospinal fluid) for subsequent diagnostic purposes using the automated extraction System. The kit contains different reagents to process clinical samples within automated platforms and facilitate the sensible detection of pathogens with real-time PCR (qPCR and RT-qPCR) assays. This product is not intended for the diagnosis, prevention, or treatment of a disease.

This Kit is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in molecular techniques and diagnostic procedures, and experience with potentially infectious human samples or materials.

2. Safety information

VIASURE Blood Pathogens Extraction Kit contains hazardous contents. For more information, please consult appropriate material safety data sheets (MSDS).

GSH classification

Reagent/Material	GHS symbol	Hazard phrases	Precaution phrases
Lysis Buffer	 Inflammable, corrosive and toxic	Danger H226, H302, H314, H412	P210, P260, P264, P280, P305+P351+P338, P310, P370+P378
Wash 1 Buffer	 Corrosive	Danger H314	P260, P264, P280, P303+P361+P353, P305+P351+P338, P310, P501
Proteinase K	 Danger	Danger H315, H319, H334, H335	P261, P280, P284, P304+P340, P342+P311, P403+P233, P501

Table 1. GHS classification

Hazard phrases

H226; Flammable liquid and vapour.

H302; Harmful if swallowed.

H314; Causes severe skin burns and eye damage.

H315; Causes skin irritation.

H319; Causes serious eye irritation.

H334; May cause allergy or asthma symptoms or breathing difficulties if inhaled.

H335; May cause respiratory irritation.

H412; Harmful to aquatic life with long lasting effects.

Precaution phrases

P210; Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P260; Do not breathe dust/fume/gas/mist/vapours/spray.

P261; Avoid breathing dust/fume/gas/mist/vapours/spray.

P264; Wash hands thoroughly after handling.

P280; Wear protective gloves/protective clothing/eye protection/face protection.

P284; [In case of inadequate ventilation] wear respiratory protection.

P304+340 IF INHALED; Remove person to fresh air and keep comfortable for breathing.

P305+P351+P338; IF IN EYES; Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P303+P361+P353; IF ON SKIN (or hair); Take off immediately all contaminated clothing. Rinse skin with water [or shower].

P310; Immediately call a POISON CENTER or doctor/physician.

P370+P378; In case of fire: Use ... for extinction.

P501; Dispose of contents/container to ...

P342+P311; If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.

P403+P233; Store in a well-ventilated place. Keep container tightly closed.

3. Summary and Explanation

Samples from the bloodstream (serum, plasma and blood) and meninges (cerebrospinal fluid) can be matrices of interest in the diagnosis of infections. The presence of viruses, bacteria or parasites in these matrices may derive from dissemination of a localized infection, microbial growth in immunocompromised patients, contamination during surgical procedures or contact with infected body fluids, through vectors, or by the pathogenesis of the microorganism itself.

The diagnosis of all these microorganisms by conventional methods has great limitations. The detection of bacteremia and fungemia is of great diagnostic importance given the prognosis of mortality between 20% and 50% of cases. Conventional methodologies based on blood culture, identification of microorganisms and subsequent antibiotic sensitivity study are laborious techniques that require a minimum of 48-72 hours. As a therapeutic strategy, empirical treatment with broad-spectrum antibiotics is started until the causative agent of the infection is identified, which is the origin of antimicrobial resistance. In the case of viruses and parasites, culture is more laborious, and serology is usually used to evaluate their presence in blood, which is less sensitive.

Molecular diagnosis reduces time, is compatible with antibiotic treatment unlike cultures, allows the detection of resistance and therefore facilitates the correct selection of the antibiotic to be administered to the patient, and is automatable, easy to implement and has a good cost-effectiveness ratio.

Nevertheless, this technique requires a previous nucleic acid (NA) extraction. The procedure to do so includes three main steps: extracting or releasing the NA from the sample and the cell or organism by lysis; separating or isolating NA from other cellular or sample material; and purification by removal of inhibitory substances. A fourth optional step, concentrating, is important for detection of low-concentration analytes.

An extraction kit for the simultaneous isolation of high quality and quantity DNA and/or RNA, from different clinical samples, and pathogens, is essential for their specific detection through molecular diagnostic assays. Nowadays, the increased number of samples to be processed at clinical laboratories requires the possibility of automating this procedure. The adaptation of the reagents to a robotic platform can improve daily workflow, standardization of sample treatment, decrease variability, reduce biohazards, and avoid contamination by reducing the analyst final hands-on time.

VIASURE Blood Pathogens Extraction Kit is designed for the isolation and purification of viral, bacterial and parasitic DNA and RNA from different clinical samples (whole blood, plasma, serum, blood culture, and cerebrospinal fluid).


4. Principle of the procedure

The VIASURE Blood Pathogens Extraction Kit procedure is based on the spontaneous and reversible adsorption of nucleic acids to superparamagnetic beads with a functionalized surface under appropriate buffer conditions. After some washing steps the nucleic acids are eluted in a purified and concentrated solution. The use of magnetic particles as solid phase for adherence of nucleic acids allows the processing of an elevated number of samples by their simultaneous exposition to a magnetic field, and therefore the automation of the process.

VIASURE Blood Pathogens Extraction Kit has been designed for the automatic platform KingFisher™ Flex Purification System (Thermo Scientific™). This type of platform procedure transfers the magnetic beads between tube/wells from one reagent to the next, by magnetic field attraction. The activation and de-activation of the magnetic field allows beads attraction and beads resuspension within different buffers. Figure 1 shows the basic principle of this type of automated extraction. First, lysis is achieved by incubation with a Lysis Buffer containing chaotropic ions supported by Proteinase K digestion. Then, the Magnetic Beads and Carrier are added to the lysate, and, after the incubation, the adsorption of the nucleic acids occurs (1). The magnetic particles are attracted to the surface of a magnet covered with plastic tips (2), to be transferred to other wells containing Wash Buffer 1 and Wash Buffer 2 (2 times), to eliminate cell debris and salts. After three washing steps (3), beads are dried to eliminate organic solvent residues (4). Finally, the DNA and/or RNA is released from the magnetic particles with an elution buffer (5), which results in a pure sample ready

for quantification and analysis. Purified DNA and RNA can be directly used for downstream applications. It is highly recommended to use extraction controls, positive/negative controls, internal controls to evaluate the complete workflow (See Procedure sections).

Figure 1. Magnetic beads extraction basic principle for a Magnetic beads handling automated instrument.

VIASURE Blood Pathogens Extraction Kit facilitates the simultaneous processing of different samples in the same extraction run, increases daily workflow, robustness and reproducibility, reduces contamination, and contact with biological hazards.

5. Reagents provided

VIASURE Blood Pathogens Extraction Kit includes the following reagents detailed in Table 2.

Reagent/Material	Description	Display	Amount
Lysis Buffer *	Liquid reagents mixture for pathogen lysis and nucleases inactivation	White bottle	1 bottle x 40 ml
Wash 1 Buffer *	Wash solution in which absolute ethanol (not included) must be added	White bottle	1 bottle x 30 ml
Wash 2 Buffer	Wash solution in which absolute ethanol (not included) must be added	White bottle	2 bottles x 20 ml
Elution Buffer	Solution for nucleic acid recovery	White bottle	1 bottle x 20 ml
Magnetic Beads	Magnetized beads suspension for nucleic acid adsorption	Transparent vial with blue cap	1 vial x 1.7 ml
Proteinase K	Enzyme solution for pathogen lysis in stabilized format	Transparent vial with cap	96 preps
Proteinase K Rehydration Buffer	Rehydration buffer for proteinase K vial resuspension	Transparent vial with yellow cap	1 vial x 1.9 ml
Carrier	Carrier for Nucleic Acid Precipitation Product	Transparent vial with orange cap	96 preps
Carrier Rehydration Buffer	Rehydration buffer for carrier vial resuspension	Transparent vial with purple cap	1 vial x 1.9 ml

Extraction control (EC)	Non-infectious nucleic acid lyophilized	Green vial with cap	1 vial
Water RNAse/DNAse Free	Rehydration buffer for Extraction Control vial resuspension	Transparent vial with transparent cap	1 vial x 1 ml

* Contains chaotropic salt. Not compatible with disinfectants containing bleach. See section 8 "Warning and precautions".

Table 2. Reagents and materials provided in VIASURE Blood Pathogens Extraction Kit.

6. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE Blood Pathogens Extraction Kit.

- Sample collection and transport system.
- Absolute ethanol.
- Micropipettes.
- Filter tips.
- Powder free disposable gloves.
- Biosafety cabinet.
- Pipettes and/or graduate test tube.
- Vortex.
- Thermomixer.
- Real-Time PCR instrument (thermocycler).
- Real-Time PCR assay.
- Ribonuclease (RNase)/ Deoxyribonuclease (DNase) free tubes for nucleic acid collection and storage.
- Positive and negative controls for the monitoring of the Extraction process. It is recommended the use of positive and negative extraction controls. (For example, VIASURE Viral Positive Control Kits).
- Sealing film for ELISA and general incubation for covering elution plate and storing nucleic acid. (Optional)

VIASURE Blood Pathogens Extraction Kit was tested and validated with the following workflow:

- Whole blood (diluted 1:2 in water), plasma, serum, blood culture, and cerebrospinal fluid samples + KingFisher™ Flex Purification System (Thermo Scientific™) + VIASURE Real Time PCR Detection Kits (Certest Biotec S.L) + VIASURE V-Lab96 Cyclor (Certest Biotec S.L).

Recommended consumables (Order the parts from your Official Distributor)

The recommended consumables required but not included in the VIASURE Blood Pathogens Extraction Kit when using the for the KingFisher™ Flex Purification System (Thermo Scientific™) system are listed below:

- KingFisher™ deep-well 96 plate - ThermoScientific™
- KingFisher™ 96 plate 200µL - ThermoScientific™
- KingFisher™ 96 tip comb for DW Magnets - ThermoScientific™

Note: In case that the consumable materials used differ from those recommended for this format, they must be validated by user.

7. Transport and storage conditions

- The kit reagents can be shipped at 15-30 °C until the expiration date which is stated on the kit label. Check all components for damages after receiving the kit. Do not use damaged kit components.
- Avoid vibrations during transport to prevent liquid leakage.
- It is recommended to keep Lysis Buffer, Wash 1 Buffer, Wash 2 Buffer, Elution Buffer, Proteinase K Rehydration Buffer, Carrier Rehydration Buffer, Extraction Control and Water RNase/DNase free vial at 15-30 °C away from light and heat sources and store in a vertical position to avoid any possible spillage.
- It is recommended to keep Magnetic Beads vials refrigerated at 2-8°C immediately after receiving the kit and for longer storage. It is recommended not to separate in aliquots.
- It is recommended to keep Proteinase K vials refrigerated at 2-8°C immediately after receiving the kit and until its resuspension. It is recommended not to separate in aliquots. Once re-suspended, keep them at -20°C avoiding freeze-thaw cycles (up to 6 times).
- It is recommended to keep Carrier vials refrigerated at 2-8°C immediately after receiving the kit and until its resuspension. It is recommended not to separate in aliquots. Once re-suspended, keep them at -20°C avoiding freeze-thaw cycles (up to 6 times).
- Once Extraction Control tube has been reconstituted use it immediately or store at -20°C until the expiration date stated on the kit label. It is recommended not to separate in aliquots. Keep it at -20°C avoiding freeze-thaw cycles (up to 6 times).
- Once the extraction has finished to ensure nucleic acid viability, keep eluates refrigerated at 4°C for up 8 hours until its addition to the different downstream applications. For long term storages eluates could be stored in the eluate plate covered with Sealing film for ELISA or transferred to a vial. It is recommended to store the supernatant obtained in the processing at 2-8°C for up to 8 hours and/or -20°C minimizing freeze-thaw cycles (up to 3 times).

The following table summarizes the conditions for transport, storage and use of the overall kit and each component:

Component	Transport Conditions	Storage Conditions	In-use conditions
Entire VIASURE Blood Pathogens Extraction Kit	15-30°C during the shelf life stated in the kit label.	* See storage conditions of each component.	* See in-use conditions of each component.
Lysis Buffer		Store away from light and heat sources. Before use: 15-30°C during the shelf life stated in the kit label. Once opened: 15-30°C during the shelf life stated in the kit label.	Room temperature.
Wash 1 Buffer		Store away from light and heat sources. Before use: 15-30°C during the shelf life stated in the kit label. Once opened: 15-30°C during the shelf life stated in the kit label.	Room temperature
Wash 2 Buffer		Store away from light and heat sources. Before use: 15-30°C during the shelf life stated in the kit label. Once opened: 15-30°C during the shelf life stated in the kit label.	Room temperature
Elution Buffer		Store away from light and heat sources. Before use: 15-30°C during the shelf life stated in the kit label. Once opened: 15-30°C during the shelf life stated in the kit label.	Room temperature
Magnetic Beads		Store away from light and heat sources. Before use: 2-8°C during the shelf life stated in the kit label. Once opened: 2-8°C during the shelf life stated in the kit label.	Immediate use at Room temperature.
Proteinase K		Store away from light and heat sources. Before use: 2-8°C during the shelf life stated in the kit label. Re-suspended component: store at -20°C during the shelf life stated in the kit label avoiding freeze-thaw cycles (up to 6 times).	Immediate use at Room temperature
Proteinase K Rehydration Buffer		Store away from light and heat sources. Before use: 15-30°C during the shelf life stated in the kit label. Once opened: 15-30°C during the shelf life stated in the kit label.	Room temperature.
Carrier		Store away from light and heat sources. Before use: 2-8°C during the shelf life stated in the kit label. Re-suspended component: store at -20°C during the shelf life stated in the kit label avoiding freeze-thaw cycles (up to 6 times).	Immediate use at Room temperature
Carrier Rehydration Buffer		Store away from light and heat sources. Before use: 15-30°C during the shelf life stated in the kit label. Once opened: 15-30°C during the shelf life stated in the kit label.	Room temperature.

Extraction Control		Store away from light and heat sources. Before use: 15-30°C during the shelf life stated in the kit label. Rehydrated control: use it immediately or store at -20°C until the expiration date stated on the kit label avoiding freeze-thaw cycles (up to 6 times).	Immediate use at Room temperature.
Water RNase/DNase free vial		Store away from light and heat sources. Before use: 15-30°C during the shelf life stated in the kit label. Once opened: 15-30°C during the shelf life stated in the kit label	Room temperature.

Table 2. Summary of the conditions for transport, storage, and use of the VIASURE Blood Pathogens Extraction Kit and each component.

8. Warnings and Precautions

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of nucleic acid extraction procedures and real-time PCR and *in vitro* diagnostic procedures, including training on nucleic acid extraction system and the real-time PCR instrument (thermocycler).
- For *in vitro* diagnostic use.
- Instructions for use of the VIASURE Blood Pathogens Extraction Kit must be read carefully before using the product. Do not perform the assay until the information about procedures, safety precautions and limitations described in the Instructions for use have been understood.
- Do not use expired reagents and/or materials.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective external box or the protective pouches are opened or damaged upon arrival.
- Do not use Extraction Control if desiccant is not present or broken inside reagent pouch. Do not remove desiccant from reagent pouch once is open.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Protect reagents against from humidity and sunlight. Prolonged exposure to humidity and sunlight may affect product performance.
- Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase)contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents. In cases where other extraction procedures are performed in the same general area of the laboratory, care must be taken to ensure that the VIASURE Blood Pathogens Extraction Kit and any additional reagents or equipment required for testing are not contaminated.
- Design a unidirectional workflow. Do not reuse reagents and do not recover the remaining liquid volume into the original bottle of the reagent once used.
- Do not reuse any disposable items. Do not refill bottles once empty due to possible contamination between reagents could affect the sensibility of the workflow.
- Use screw-capped tubes and prevent any potential splashing or cross contamination of specimens during preparation.
- After the first use, once the bottles or the tubes have been opened, close each bottle/tube with the bottle caps/tubes caps supplied to store them. Mark each bottle cap/tube cap with an indelible marker. Do not mix the bottle cap from one bottle to the bottle cap from another bottle. Do not mix the tubes caps from one tube to another component tube cap.
- Do not mix tubes from VIASURE Blood Pathogens Extraction Kit with those from VIASURE Real Time PCR Detection Kit. Please consult section 9.3.4 Lyophilized extraction control.
- Avoid possible contamination of reagents with extracted nucleic acids, PCR products, and positive control. To prevent contamination of reagents, use of filter-tips is recommended.

- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles, and mask. Do not eat, drink or smoke or apply cosmetic products in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local, state, and federal regulations.
- All reagents should be handled with care, minimizing the time of exposure to the open, especially those reconstituted with absolute ethanol.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- Reliability of the results depends on adequate specimen collection, storage, transport, and processing procedure.
- Use graduated test tubes to fill the bottles with ethanol.
- In the case of a subsequent qPCR analysis, consult each qPCR kit or real-time PCR instrument's reference manual for additional warnings, precautions, and procedures.
- Lysis Buffer and Wash 1 Buffer contain chaotropic salt. Not compatible with disinfectants containing bleach. Consult the hazards and precautions when handling Lysis Buffer and Wash 1 Buffer in the Material Safety Data Sheets provided with the Device.
- Make sure there is not precipitate in the supplied reagents. Lysis Buffer and Wash 1 Buffer contain guanidine salts that may precipitate. If there is visible precipitate at the bottom of the bottle, resuspend by mild heating and shake until its complete dissolution.
- Make sure that the Magnetic Beads and Carrier are completely resuspended. Shake vigorously the stored vial on a Vortex.
- Proteinase K could precipitate after freeze-thaw cycles. Make sure that Proteinase K is completely resuspended. Shake vigorously the storage vial on a Vortex.
- Before dispensing VIASURE Blood Pathogens Extraction Kit into the platform consumables open carefully each well avoiding reagent spilling and reconstitute the indicated wells as indicated in the procedure section.
- Make sure all tubes or plates have been properly placed following the indicated layout and fit correctly into the specific adapters of the Automated extraction system.
- The waste generated with the VIASURE Blood Pathogens Extraction Kit has not been tested for residual infectious material. A contamination of the liquid waste with residual infectious material is highly unlikely due to strong denaturing Lysis Buffer and Proteinase K treatment but it cannot be excluded completely. Therefore, liquid waste must be considered infectious and should be handled and discarded according to local safety regulations.
- The handling of substances and the disposal of waste may be subject to local, state, or federal law, or to regulations regarding health, environment, or safety. Strictly observe the corresponding provisions. Consult appropriate material safety data sheets (MSDS) for more information.

9. Test procedure

9.1. Sample collection, storage and transport

VIASURE Blood Pathogens Extraction Kit is suitable for various potentially infectious clinical samples: whole blood, plasma, serum, blood culture, and cerebrospinal fluid. For reproducible and high yields appropriate sample storage is essential. Yields may vary from sample to sample depending on factors such as health and age of the donor, kind of sample, collection date, transport, and storage conditions. **Whole blood (diluted 1:2 in water)** and **plasma samples** collected in EDTA (S-Monovette® K3 EDTA, 2.7 ml (Starstedt)), Citrate (S-Monovette® Citrate 3ml 9NC, (Starstedt)), Heparine (S-Monovette® Citrate 3ml 9NC, (Starstedt)), Reveos® LR Set CPD/SAG-M Blood Bag Set with RBC leukoreduction filter (TERUMOBCT), Reveos Automated Blood Processing System, anticoagulant citrate phosphate dextrose (CPD); **serum samples** collected in Clot activator polymer (S-Monovette® Serum Z-Gel, clotting activator/gel, 7.5 ml, cap brown (Starstedt)) and Human serum (Merk); and **blood culture** collected in BD BACTECT™ Plus Aerobic medium in plastic culture vials and BD BACTECT™ Lytic Anaerobic medium in plastic culture vials.

Heparin is a well-known RT-PCR inhibitor, therefore, is recommendable to avoid heparin collection tubes.

As general recommendations, the specimens should be collected as soon as possible, depending on the pathogen during the acute phase of illness or at any time during the clinical course, and before antimicrobial therapy begin, if possible. Specimens should be collected in a sterile container avoiding touching adjacent tissues and organ secretions that are not of interest. All specimens must be labeled with the name and identification number of the person from whom the sample was collected, the source of the specimen, and the date and time it was collected. After collection, specimens should be placed in a biohazard bag and transported to the laboratory as soon as possible.

Depending on the type of specimen, the collection differs. Follow the recommendations by McPherson R.A., Pincus M.R. Henry's Clinical Diagnosis and Management by Laboratory Methods. 23RD edition. Elsevier Health Sciences; 2017.

According to the literature, the conditions of transport, storage and stability depend on the type of the sample matrix and the pathogen to be detected. It is recommended to use fresh specimens for the test. The specimens should be held at room temperature for cerebrospinal fluid (CSF), blood and other body fluids. If a delay is unavoidable (< 2 hours until they can be sent), the general rule is to ship at -20°C or ideally at -70°C or - 80°C.

Optimally, all specimen containers should be opened in a biological safety cabinet. Appropriate barriers should be used to prevent exposure of skin and mucous membranes to the specimen. Gloves and a lab coat must always be worn when handling patient specimens. Masks, goggles (or working behind a plastic shield), and impermeable gowns or aprons must be worn when there is a risk for splashes or droplet formation. Wash your hands once the procedure has finished.

Specimens submitted for molecular testing must be stored in controlled conditions so that nucleic acids do not degrade during storage. Avoid freezing and thawing cycles. Some pathogens from specimens that are frozen and then thawed could be degraded and may result in false-negative test results.

Depending on the clinical specimen and pathogen different sample pretreatments may be needed. Follow the appropriate protocol laboratory guideline in each case.

For Collecting, Handling, and Testing Clinical Specimens consult: the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. Clinical Infectious Diseases, 67(6), e1-e94).

9.2. Sample pretreatment

Depending on the clinical specimen and pathogen different sample pretreatments may be needed. It is important to avoid substances or objects that may interfere in sample pipetting during the automated process, therefore, caps or swabs should be removed before introducing the samples into the sample racks, and solid or highly viscous samples must be pretreated.

- Blood culture: Follow manufacturer instructions for use for microorganism growth. After incubation, transfer supernatant to a vial/tube.
- Whole blood: this matrix must be diluted 1/2 in DNase/RNase free water for performance in molecular biology applications.

The rest of validated samples can be directly introduced into the platform.

9.3. Reagents' resuspension / preparation

Please, carefully read each reagents preparation or resuspension sections to be aware of all the precautions and recommendations it should be taken.

Before starting with the nucleic acid extraction and purification, it is necessary in some reagents to add a certain volume of other reagents (some of them are not included in the kit). The following table indicates which reagent should be added and the required volume:

Reagent/Material	Reagent for resuspension / preparation	Volume
Wash 1 Buffer	Absolute ethanol	70 mL
Wash 2 Buffer (I)	Absolute ethanol	80 mL
Wash 2 Buffer (II)	Absolute ethanol	80 mL
Proteinase K	Proteinase K Rehydration Buffer	1.7 mL
Carrier	Carrier Rehydration Buffer	1.7 mL
Extraction Control	Water RNase/DNase Free	1 mL
Magnetic Beads	-	-

Table 3. Reagents and volumes to be resuspended/prepared in the Reagents.

After the first use, once the bottles or the tubes have been opened, close each bottle/tube with the bottle caps/tubes caps supplied to store them. Mark each bottle cap/tube cap with an indelible marker. Do not mix the bottle cap from one bottle to the bottle cap from another bottle. Do not mix the tubes caps from one tube to another component tube cap. Do not mix tubes from VIASURE Blood Pathogens Extraction Kit with those from VIASURE Real Time PCR Detection Kit.

9.3.1 Buffer's preparation

Lysis Buffer and Wash 1 Buffer contain guanidine salts that may precipitate. If there is visible precipitate at the bottom of the bottle, resuspend by mild heating and shake until its complete dissolution before open the bottles for the first use. In the following uses, ensure the correct closure with the bottle cap before doing this step.

Before starting with the nucleic acid extraction and purification, it is necessary in some reagents to add a certain volume of other reagents (not included in the kit). The following table indicates which reagent should be added and the required volume:

Reagent/Material	Reagent for preparation	Volume
Wash 1 Buffer*	Absolute ethanol	70 mL
Wash 2 Buffer (I)*	Absolute ethanol	80 mL
Wash 2 Buffer (II)*	Absolute ethanol	80 mL

Table 4. Reagents and volumes to be added to the Wash 1 Buffer and Wash 2 Buffer.

**Note: To facilitate the process, a mark has been included on the label of each reagent to tick when the reagents are added.*

Besides, it is recommended to make sure to write with indelible marker the specific automated extraction system bottles used once prepared in the space with the resuspension date and other required information.

Use graduated test tubes to fill the bottles with absolute ethanol.

9.3.2 Proteinase K resuspension

Before starting with the nucleic acid extraction and purification, it is necessary to resuspend Proteinase K (Transparent vial with cap) with 1.7 mL of Proteinase K Rehydration Buffer (Transparent vial with yellow cap) included in the kit.

It is recommended to keep Proteinase K vials refrigerated at 2-8°C immediately after receiving the kit and once resuspended keep at -20°C avoiding freeze-thaw cycles (up to 6 times).

9.3.3 Carrier resuspension

Before starting with the nucleic acid extraction and purification, it is necessary to resuspend Carrier (Transparent vial with orange cap) with 1.7 mL of Carrier Rehydration Buffer (Transparent vial with purple cap) included in the kit. Make sure that the Carrier are completely resuspended. Shake vigorously the storage vial on a Vortex.

It is recommended to keep Carrier vials refrigerated at 2-8°C immediately after receiving the kit and once resuspended, keep at -20°C avoiding freeze-thaw cycles (up to 6 times).

9.3.4 Lyophilized extraction control

The use of the Extraction Control (EC) will be required only for VIASURE Real Time PCR products which do not include an internal control in the lyophilized reaction mix. Please consult Instructions for use of each VIASURE Real Time PCR products to check this information. For these VIASURE Real Time PCR products use the Extraction Control (EC) provided in the VIASURE Blood Pathogens Extraction Kit.

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the PCR-Amplification area. Reconstitute the lyophilized Extraction Control (green vial) adding 1 mL of Water RNase/DNase free (Transparent vial with transparent cap) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C and minimize freeze and thaw cycles (up to 6 cycles).

Note 1: Add 5 µL of the rehydrated EC to the mixture sample and Lysis Buffer to monitor the extraction process and/or discard the inhibition of the polymerase activity.

Note 2: Because of the nature of extraction control synthesis or production do not use the Extraction Control (EC) included in this kit (VIASURE Blood Pathogens Extraction Kit) on samples to be analyzed by qPCR for *E.Coli* and/or *TEM* gene (β-lactamases resistance) detection. In this case, please use the Extraction Control (EC) vial included in the VIASURE Real Time PCR Detection Kit.

9.3.5 Magnetics Beads

Before starting with the nucleic acid extraction and purification, make sure that the magnetic beads are completely resuspended. Shake vigorously the storage vial on a Vortex. Magnetic Beads are included in a 1.7 mL Transparent vial with blue cap.

It is recommended to keep Magnetic Beads vials refrigerated at 2-8°C immediately after receiving the kit.

9.4. Extraction procedure

The procedure below provides instructions on the steps to process samples. Refer to the KingFisher™ Flex Purification System (Thermo Scientific™) user manual for instructions on how to use the platform and the extraction procedure to follow within KingFisher® Flex automatic platforms. *(Note: In case of using VIASURE Blood Pathogens Extraction Kit with another automated extraction system different from the one mentioned in these Instructions for use, the procedure must be validated by the user).*

1. Load with reagents 96 Deep well plates

Load the following five plates with the reagents and volumes in the required wells depending on the number of samples to analyze.

- a. Sample Plate: Lysis Buffer 300 µl, sample 200 µl, rehydrated Extraction Control 5 µl (optional), Proteinase K 15 µl, Magnetic Beads 15 µl and Carrier 15 µl.
- b. Wash 1 Plate: Wash 1 Buffer 800 µl.
- c. Wash 2 (I) Plate: Wash 2 Buffer 800 µl.
- d. Wash 2 (II) Plate: Wash 2 Buffer 800 µl.
- e. Elution Plate: Elution Buffer 100 µl.

Note: It is recommended to identify/mark each plate. Pipette each component separately to avoid possible contamination and to ensure the correct functionality of the reagents.

2. Set Up the instrument.

Start the instrument

Select the Blood_Pathogens Extraction protocol using the cursor keys and press START or use BindIt Software to run the desired protocol via the PC. If not created, please follow the instructions below and create it with BindIt Software.

1. Define a plate Layout (enter name and reagents information as described in 1 Load with reagents 96 Deep well plates).
2. Configure the protocol adding the steps described in the table below:

Step name	Procedure	Setting
Pick Up	96 Tip plate	
Lysis and Binding	Sample plate	
Beginning of step	Precollect	No
	Release beads	Yes
Mixing/Heating parameters	Mixing time (hh:mm:ss), speed	00:10:00, Fast
	Heating during mixing	Yes
	Block temperature (°C)	65
	Preheat	Yes
End of step	Postmix	No
	Collect beads, count	5
	Collect time (s)	0
Collect Beads	Sample plate	
	Collect count	5
	Collect time (s)	4
Wash 1	Wash 1 plate	
Beginning of step	Precollect	No
	Release beads (hh:mm:ss), speed	00:00:20, Bottom mix
Mixing/Heating parameters	First mix time (hh:mm:ss), speed	00:00:20, Bottom mix
	Second mix time (hh:mm:ss), speed	00:00:20, Fast
	Loop count	3
	Heating during mixing	No
End of step	Postmix	No
	Collect beads, count	5
	Collect time (s)	3
Wash 2 (I)	Wash 2 (I) plate	
Beginning of step	Precollect	No
	Release beads (hh:mm:ss), speed	00:00:20, Fast
Mixing/Heating parameters	First mix time (hh:mm:ss), speed	00:00:10, Bottom mix
	Second mix time (hh:mm:ss), speed	00:00:20, Fast
	Loop count	3
	Heating during mixing	No

End of step	Postmix	No
	Collect beads, count	4
	Collect time (s)	2
Wash 2 (II)	Wash 2 (II) plate	
Beginning of step	Precollect	No
	Release beads (hh:mm:ss), speed	00:00:20, Fast
Mixing/heating parameters	First mix time (hh:mm:ss), speed	00:00:10, Bottom mix
	Second mix time (hh:mm:ss), speed	00:00:20, Fast
	Loop count	3
	Heating during mixing	No
End of step	Postmix	No
	Collect beads, count	4
	Collect time (s)	2
Dry	Wash 2 (II) plate	
	Dry time (hh:mm:ss)	00:02:00
	Tip position	Outside well/tube
Elution	Elution Plate	
Beginning of step	Precollect	No
	Release beads (hh:mm:ss)	Yes, 00:00:00
Mixing/Heating parameters	First mix time (hh:mm:ss), speed	00:00:20, Bottom mix
	Second mix time (hh:mm:ss), speed	00:00:45, Medium
	Loop count	6
	Heating during mixing	Yes
	Block temperature (°C)	75
	Preheat	Yes
End of step	Postmix	No
	Collect beads, count	5
	Collect time (s)	2
Collection beads	Elution plate	
Beginning of step	Precollect	No
	Release beads	No
Mixing/Heating parameters	Mixing time (hh:mm:ss), speed	00:02:00, Slow
End of step	Collect beads	No
Leave	Tip plate	

Table 4. Complete settings of the extraction process to be followed within KingFisher® Flex automatic platforms.

3. Plate Set Up into the instrument.

Please follow the instrument graphical interface and fill the platform with the plates following the next order:

Note. Make sure all plates are loaded in the same direction, with the well A1 facing the top right.

- 1st. Tip plate (KingFisher™ 96 plate 200µL - ThermoScientific™ + KingFisher™ 96 tip comb for DW Magnets - ThermoScientific™).
- 2nd. Elution plate.
- 3rd. Wash 2 (II) Buffer plate.
- 4th. Wash 2 (I) Buffer plate.
- 5th. Wash 1 Buffer plate.
- 6th. Sample plate.

Once plates have been loaded click start and the run will be done in 38 minutes.

4. Purified nucleic acids recovery

Once the extraction has finished to ensure nucleic acid viability, keep eluates refrigerated at 4°C for up to 8 hours until its addition to the different downstream applications. For long term storages eluates could be stored in the eluate plate covered with Sealing film for ELISA or transferred to a vial. It is recommended to store the supernatant obtained in the processing at 2-8°C for up to 8 hours and/or -20°C minimizing freeze-thaw cycles (up to 3 times).

Once the extraction and PCR setup transfer the loaded plate to thermocycler for nucleic acid amplification following Manufacturer's indications.

5. Decontamination

Decontaminate properly the platform inner surfaces with appropriate disinfectant solvents and discard used plastic consumables.

For more details consult KingFisher™ Flex Purification System (Thermo Scientific™) user manual.

9.5 Real-Time PCR

The nucleic acids obtained with VIASURE Blood Pathogens Extraction Kit are compatible with any commercially available qPCR or RT-qPCR kit, although it has only been validated with the VIASURE Real Time PCR Detection kits.

Manufacturer's instructions of the VIASURE Real Time PCR Detection Kits must be followed (<https://www.certest.es/es/viasure/>). The positive and negative controls included in the extraction step must not substitute the positive and negative controls provided with the qPCR kits.

Recommendation is to use the thermocycler VIASURE V-Lab96 Cyclor (Certest Biotec S.L) for the analysis and interpretation of the results.

10. Result interpretation

To correctly interpret the results of the controls and clinical samples, please consult the Manufacturer's Instructions of the different downstream applications.

11. Limitations of the procedure

- The efficiency of the performance characteristics of VIASURE Blood Pathogens Extraction Kit for KingFisher™ Flex Purification System (Thermo Scientific™) depends on the specimens' collection, transport, storage, and/or handling methods and the efficiency of the flowchart used by each laboratory: qPCR assay and thermocycler. Each complete procedure, different from the one mentioned in these Instructions for use, must be validated by the end user.
- For appropriate yield and quantity, sample collection, storage, and pre-treatment must be performed following the specific recommendations of official laboratory guidelines.
- If the consumable materials used differ from those recommended for this format, they must be validated by user.
- VIASURE Blood Pathogens Extraction Kit can be automated in the automated platform KingFisher™ Flex Purification System (Thermo Scientific™). Be sure to read this Instruction Manual and the KingFisher™ Flex's User Manual carefully before use. If the automation is performed in a different magnetic bead handling platform with other consumables, it must be validated by user.
- Bottles should be stored at 15-30 °C away from light and heat sources and store in a vertical position to avoid any possible spillage, Magnetic Beads refrigerated at 2-8°C and Proteinase K, Carrier and the Extraction Control stored at -20°C (up to 6 freeze-thaw cycles). It is also important to check the storage conditions of the components once the kit has been opened. Please consult section 7 of the instructions for use of the kit.
- Because of the nature of extraction control synthesis or production do not use the Extraction Control (EC) included in this kit (VIASURE Blood Pathogens Extraction Kit) on samples to be analyzed by qPCR for *E.Coli* and/or *TEM* gene (β -lactamases resistance) detection. In this case, please use the Extraction Control (EC) vial included in the VIASURE Real Time PCR Detection Kit.
- Whole blood must be diluted 1:2 in water due to a correct sample pretreatment. If whole blood is not diluted, it could affect instrument pipetting.
- Heparin is a well-known RT-PCR inhibitor, therefore, is recommendable to avoid heparin collection tubes. (As anticoagulant or collection tubes).

12. Quality control

The expected values of VIASURE Blood Pathogens Extraction Kit depend on the diagnostic procedure (extraction and qPCR assay), so each complete procedure must be validated by the final user.

It is recommended to use the VIASURE Blood Pathogens Extraction Kit in conjunction with appropriate controls to monitor the entire process for different downstream applications. To control the complete process, from sample processing to analysis, perform the following controls: Positive Control (not provided), Negative Control (not provided) and Extraction Control (provided with the kit).

13. Analytical performance characteristics

The target models chosen for the VIASURE Blood Pathogens Extraction Kit design, optimization and validation procedure should generate the maximum number of different specimen types, so that each target model is

associated with a specimen type. These target models have been chosen based on the following criteria: clinical frequency, commercial sales of VIASURE qPCR products, difficulty in lysis, NA target, and NA size.

13.1 Analytical linearity

The linearity of the assay was determined by testing specific dilutions of each validated pathogen in different biological matrix using certified reference strain as shows the table below.

Pathogen	Reference strain	Matrix	Linearity range
Zika Virus	1st WHO International Standard for Zika Virus RNA (PEI code: 11468/16)	Serum	$1 \times 10^4 - 1.1 \text{ IU}/\mu\text{l}$
		Whole blood	$1 \times 10^4 - 1.1 \text{ IU}/\mu\text{l}$
<i>Plasmodium falciparum</i>	1st WHO International Standard for Plasmodium falciparum DNA Nucleic Acid Amplification Techniques (NIBSC code: 04/176)	Whole blood	$5 \times 10^4 - 5 \text{ IU}/\mu\text{l}$
<i>Klebsiella pneumoniae</i>	<i>Klebsiella pneumoniae</i> (ATCC® BAA-2146™)	Aerobic blood culture	$5 \times 10^1 - 5 \times 10^{-3} \text{ CFU}/\mu\text{l}$
		Anaerobic blood culture	$1 \times 10^2 - 1 \times 10^{-2} \text{ CFU}/\mu\text{l}$
<i>Cytomegalovirus</i>	1st World Health Organization Standard for Human Cytomegalovirus for Nucleic acid Amplification techniques (NIBSC code: 09/162)	Plasma	$1 \times 10^2 - 1 \text{ IU}/\mu\text{l}$

Table 5. Linearity.

In conclusion, all real-time PCR assays showed an acceptable efficiency and linearity, (R^2) were >0.9 in all the target reactions tested.

13.2 Analytical sensitivity. Limit of Detection (LoD)

Analytical sensitivity or limit of detection (LoD) was determined based on the LoD obtained in the (RT)-qPCR assay used for the validation procedure. Samples were extracted with VIASURE Blood Pathogens Extraction Kit in the KingFisher™ Flex Purification System (Thermo Scientific™) and analyzed with the thermocycler VIASURE V-Lab96 Cyclor (Certest Biotec S.L).

The results obtained showed that the LoD are the following for the chosen model pathogens, with a positive rate of 95%:

Pathogen	Reference strain	Matrix	LOD
Zika Virus	1st WHO International Standard for Zika Virus RNA (PEI code: 11468/16)	Serum	$0.5 \text{ IU}/\mu\text{l}$
		Plasma	$1 \text{ IU}/\mu\text{l}$
		Whole blood	$0.5 \text{ IU}/\mu\text{l}$
<i>Plasmodium falciparum</i>	1st WHO International Standard for Plasmodium falciparum DNA Nucleic Acid Amplification Techniques (NIBSC code: 04/176)	Whole blood	$5 \text{ IU}/\mu\text{l}$
<i>Klebsiella pneumoniae</i>	<i>Klebsiella pneumoniae</i> (ATCC® BAA-2146™)	Aerobic blood culture	$0.005 \text{ CFU}/\mu\text{l}$
		Anaerobic blood culture	$0.01 \text{ CFU}/\mu\text{l}$
<i>Cytomegalovirus</i>	1st World Health Organization Standard for Human Cytomegalovirus for Nucleic acid Amplification techniques (NIBSC code: 09/162)	Plasma	$1 \text{ IU}/\mu\text{l}$
<i>Neisseria meningitidis</i>	<i>Neisseria meningitidis</i> (ATCC®13091™)	CSF	$2 \text{ CFU}/\mu\text{l}$

Table 6. Limit of Detection (LoD).

13.3 Accuracy

13.3.1 Trueness (Veracity)

The veracity of VIASURE Blood Pathogens Extraction Kit was evaluated with reference material listed below.

Reference strains: the validity of obtained results during the product validation has been verified using different reference strains for each pathogen. Strains used:

- **The American Type Culture Collection ("ATCC®").**
 - *Klebsiella pneumoniae* (ATCC® BAA-2146™)
 - *Neisseria meningitidis* (ATCC®13091™)
 - Mumps Virus (ATCC®VR-106™)
- **National Institute for Biological Standards and Control (NIBSC).**
 - 1st WHO International Standard for Plasmodium falciparum DNA Nucleic Acid Amplification Techniques (NIBSC code: 04/176).
 - 1st WHO International Standard for Human Cytomegalovirus for Nucleic acid Amplification techniques (NIBSC code: 09/162).
 - 1st WHO International Standard for HHV -6B virus DNA Nucleic Acid Amplification Techniques (NIBSC code: 15/266).
 - 1st WHO International Standard for BK virus DNA Nucleic Acid Amplification Techniques (NIBSC code: 14/212).
- **The Paul-Ehrlich-Institut (PEI).**
 - 1st WHO International Standard for Zika Virus RNA (PEI code: 11468/16).
 - 1st WHO International Standard for Chikungunya virus RNA Nucleic Acid Amplification Techniques (PEI code: 11785/16).
 - PEI Reference Preparation WNV Lineage 1 RNA (#13299/19).
- **bei RESOURCES.**
 - Dengue virus type 1 (DEN-1) strain Hawaii (bei RESOURCES Ref: NR-82 (derived from ATCC® VR-1254™).
 - *Leishmania donovani* strain 1S (MHOM/SD/62/1S) (bei RESOURCES Ref: NR-48821).

13.3.2 Precision

To determine the precision of VIASURE Blood Pathogens Extraction Kit, intra-assay (repeatability), inter-assay (reproducibility), inter-batch (intermediate precision) and inter-thermocycler (site-to-site) assays were performed. The panel of samples were analyzed following the workflow VIASURE Blood Pathogens Extraction Kit using the KingFisher™ Flex Purification System (Thermo Scientific™), VIASURE *Malaria differentiation* Real Time PCR Detection Kit, VIASURE Zika Virus Real Time PCR Detection Kit and VIASURE V-Lab96 Cycler (Certest Biotec S.L).

Intra-assay. Repeatability

The repeatability was tested by analyzing replicates all samples will be tested up to twelve replicates in one day using VIASURE Blood Pathogens Extraction Kit and the selected workflow. The following tables shows the results obtained in this assay.

Sample	Target	Viasure channel	\bar{x} (Ct)	σ	CV%
Positive 1 (2-3xLoD)	<i>Plasmodium falciparum</i>	FAM	32.61	0.46	1.41
Positive 2 (10xLoD)	<i>Plasmodium falciparum</i>	FAM	31.14	0.40	1.28
Negative sample	<i>Plasmodium falciparum</i>	FAM	Neg	n.a.	n.a.
	Extraction Control	HEX	21.44	0.34	1.57

Table 7. LoD. Intra-assay repeatability of VIASURE Blood Pathogens Extraction Kit and VIASURE *Malaria differentiation* Real Time PCR Detection Kit. (Ct) = threshold cycle. (\bar{x}) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Sample	Target	Viasure channel	\bar{x} (Ct)	σ	CV%
Positive 1 (2-3xLoD)	Zika Virus	FAM	34.93	1.06	3.04
Positive 2 (10xLoD)	Zika Virus	FAM	33.29	0.43	1.29
Negative sample	Zika Virus	FAM	Neg	n.a.	n.a.
	Extraction Control	HEX	19.70	0.15	0.75

Table 8. Intra-assay repeatability of VIASURE Blood Pathogens Extraction Kit and VIASURE *Zika Virus* Real Time PCR Detection Kit. (Ct) = threshold cycle. (\bar{x}) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Inter-assay. Reproducibility

The reproducibility was tested by analyzing the different samples on three different days by three different operators with the VIASURE Blood Pathogens Extraction Kit and the selected workflow. A summary of results is shown in the tables below.

Sample	Target	Viasure channel	\bar{x} (Ct)	σ	CV%
Positive 1 (2-3xLoD)	<i>Plasmodium falciparum</i>	FAM	33.89	1.09	3.22
Positive 2 (10xLoD)	<i>Plasmodium falciparum</i>	FAM	32.57	1.09	3.34
Negative sample	<i>Plasmodium falciparum</i>	FAM	Neg	n.a.	n.a.
	Extraction Control	HEX	22.54	1.75	7.76

Table 9. Inter-assay reproducibility of VIASURE Blood Pathogens Extraction Kit and VIASURE *Malaria differentiation* Real Time PCR Detection Kit. (Ct) = threshold cycle. (\bar{x}) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Sample	Target	Viasure channel	\bar{x} (Ct)	σ	CV%
Positive 1 (2-3xLoD)	Zika Virus	FAM	34.00	1.07	3.14
Positive 2 (10xLoD)	Zika Virus	FAM	32.49	0.62	1.92
Negative sample	Zika Virus	FAM	Neg	n.a.	n.a.
	Extraction Control	HEX	22.17	1.67	7.55

Table 10. Inter-assay reproducibility of VIASURE Blood Pathogens Extraction Kit and VIASURE Zika Virus Real Time PCR Detection Kit. (Ct) = threshold cycle. (\bar{x}) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Inter-batch. Intermediate precision

The inter-batch values were determined with three replicates of the different samples by using three batches of VIASURE Blood Pathogens Extraction Kit and the selected workflow. The results are described in the following tables.

Sample	Target	Viasure channel	\bar{x} (Ct)	σ	CV%
Positive 1 (2-3xLoD)	<i>Plasmodium falciparum</i>	FAM	34.21	0.48	1.41
Positive 2 (10xLoD)	<i>Plasmodium falciparum</i>	FAM	31.95	0.42	1.32
Negative sample	<i>Plasmodium falciparum</i>	FAM	Neg	n.a.	n.a.
	Extraction Control	HEX	23.67	0.33	1.39

Table 11. Inter-batch Intermediate precision of VIASURE Blood Pathogens Extraction Kit and VIASURE Malaria differentiation Real Time PCR Detection Kit. (Ct) = threshold cycle. (\bar{x}) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Sample	Target	Viasure channel	\bar{x} (Ct)	σ	CV%
Positive 1 (2-3xLoD)	Zika Virus	FAM	34.42	1.50	4.37
Positive 2 (10xLoD)	Zika Virus	FAM	32.48	0.74	2.29
Negative sample	Zika Virus	FAM	Neg	n.a.	n.a.
	Extraction Control	HEX	20.78	1.04	4.99

Table 12. Inter-batch Intermediate precision of VIASURE Blood Pathogens Extraction Kit and VIASURE Zika Virus Real Time PCR Detection Kit. (Ct) = threshold cycle. (\bar{x}) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Site-to-site. Inter-thermocycler

The inter-thermocycler values were determined with three replicates of the same samples used for intra-assay, inter-assay and inter-batch, using the VIASURE Blood Pathogens Extraction Kit, VIASURE *Malaria differentiation* Real Time PCR Detection Kit and VIASURE *Zika Virus* Real Time PCR Detection Kit in the thermocyclers CFX96™ Real-Time PCR Detection System (Bio-Rad), DTprime Real-time Detection Thermal Cycler, AriaMx Real-Time PCR System (Agilent Technologies), VIASURE V-Lab96 Cycler (Certest Biotec S.L.) and QuantGene 9600 qPCR equipment (BIOER). The results are described in the following tables.

Sample	Target	Viasure channel	\bar{x} (Ct)	σ	CV%
Positive 1 (2-3xLoD)	<i>Plasmodium falciparum</i>	FAM	37.06	1.00	2.69
Positive 2 (10xLoD)	<i>Plasmodium falciparum</i>	FAM	34.10	3.85	11.30
Negative sample	<i>Plasmodium falciparum</i>	FAM	Neg	n.a.	n.a.
	Extraction Control	HEX	22.15	1.02	4.63

Table 13. Inter-thermocycler of VIASURE Blood Pathogens Extraction Kit and VIASURE *Malaria differentiation* Real Time PCR Detection Kit. ((Ct) = threshold cycle. (\bar{x}) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

Sample	Target	Viasure channel	\bar{x} (Ct)	σ	CV%
Positive 1 (2-3xLoD)	Zika Virus	FAM	35.71	1.31	3.66
Positive 2 (10xLoD)	Zika Virus	FAM	32.55	0.83	2.56
Negative sample	Zika Virus	FAM	Neg	n.a.	n.a.
	Extraction Control	HEX	22.56	1.89	8.37

Table 14. Inter-thermocycler of VIASURE Blood Pathogens Extraction Kit and VIASURE *Zika Virus* Real Time PCR Detection Kit. ((Ct) = threshold cycle. (\bar{x}) = arithmetic mean Ct value, (σ) = standard deviation, (CV %) = coefficient of variation, Neg = negative, n.a.= not applicable.

13.4 Analytical specificity

Interferences and inhibitors

An interfering substance study was performed to test the possible interfering effect of endogenous and exogenous substances on VIASURE Blood Pathogens Extraction Kit using VIASURE Real Time PCR Detection Kits in the thermocycler VIASURE V-Lab96 Cycler (Certest Biotec S.L.). A total of 27 (for whole blood and serum) and 20 (for plasma) potentially interfering substances were evaluated.

Evaluation of the interference in the whole process (extraction and the (RT)-qPCR)

The following results were obtained in the analysis in whole blood (diluted 1:2 in water).and serum samples:

Substance Name	Concentration tested	Result
Zinc / Zinc sulphate	1.20E-02 mg/ml	N.I
Acetylsalicylic Acid	3.00E-02 mg/ml	N.I
Ibuprofen	2.19E-01 mg/ml	N.I
Chloroquine sulphate/phosphate	6.00E-01 mg/ml	N.I
Acyclovir	6.60E-02 mg/ml	N.I
Foscarnet	7.00E-01 mg/ml	N.I
Oseltamivir	3.99E-04 mg/ml	N.I
Ampicillin	7.50E-02 mg/ml	N.I
Ceftriaxone	8.40E-01 mg/ml	N.I
Doxycycline hyclate	1.80E-02 mg/ml	N.I
Metronidazole	1.23E-01 mg/ml	N.I
Vancomycin	1.20E-01 mg/ml	N.I
EDTA	9.90E-04 mg/ml	N.I
Heparin	3.30E+00 units/ml	N.I
Lorazepam	7.20E-04 mg/ml	N.I
Albumin	1.00E+01 mg/ml	N.I
Bilirubin	4.00E-01 mg/ml	N.I
Genomic DNA	3.50E-03 mg/ml	N.I
Glucose	1.00E+01 mg/ml	N.I
Human Hemoglobine	1.00E+01 mg/ml	N.I
Immunoglobulin G	2.00E+01 mg/ml	N.I
Lactate	9.00E-01 mg/ml	N.I
Leukocytes	6.25E+06 cell/ml	N.I
Triglycerides	1.50E+01 mg/ml	N.I
4-Acetamidophenol	1.56E-01 mg/ml	N.I
Dexamethasone	1.20E-02 mg/ml	N.I
Omeprazole	8.40E-03 mg/ml	N.I

Table 15. Potential interference substances. N.I: No reportable interfere.

The following results were obtained in the analysis in plasma samples:

Substance Name	Concentration tested	Result
Acetylsalicylic Acid	3.00E-02 mg/ml	N.I
Ibuprofen	2.19E-01 mg/ml	N.I
Naproxen	3.60E-01 mg/ml	N.I
Cidofovir	8.10E-02 mg/ml	N.I
Foscarnet	7.00E-01 mg/ml	N.I
Ganciclovir	3.20E-02 mg/ml	N.I
Valganciclovir	2.20E-02 mg/ml	N.I
Ethylenediaminetetraacetic acid (EDTA)	9.90E-04 mg/ml	N.I
Heparin	3.30E+00 units/ml	N.I
	7.50E-09 mg/ml	N.I
Doxycycline hyclate	1.00E+01 mg/ml*	N.I
17- α -Ethinylestradiol	4.00E-01 mg/ml	N.I
Albumin	3.50E-03 mg/ml	N.I

Bilirubin	1.00E+01 mg/ml	N.I
Genomic DNA	1.50E+01 mg/ml	N.I
Human hemoglobine	2.58E-03 mg/ml	N.I
Triglycerides	1.13E+00 mg/ml	N.I
Azathioprine	4.20E-02 mg/ml	N.I
Ciclosporin	9.90E-05 mg/ml	N.I
Mycophenolate mofetil	1.56E-01 mg/ml	N.I
Prednisone	3.00E-02 mg/ml	N.I
4-Acetamidophenol	2.19E-01 mg/ml	N.I

Table 16. Potential interference substances. N.I: No reportable interfere.

In conclusion, different potentially interfering substances have been tested in whole blood (1:2 diluted), serum and plasma. Thus, no interference of any of the evaluated substances was observed at the concentrations indicated for each workflow.

13.5 Metrological traceability

This assay is not designed for measuring purposes.

14. Clinical Performance characteristic

In order to determine the clinical diagnostic accuracy, several evaluations with different clinical sample matrix were conducted at Certest Biotec S.L. (Zaragoza, Spain) through collaboration with national entities.

The clinical performance of VIASURE Blood Pathogens Extraction Kit using the KingFisher™ Flex Purification System (Thermo Scientific™) was tested using clinical spiked samples (79 whole blood, 87 serums, and 82 plasma) and with 26 blood culture of patients with suspected bacteraemia. For the evaluation using contrived samples, the following strains were used to generate at least n= 30 spiked samples at 2x LoD, n=10 at 5x LoD, n=10 at 10x LoD of *Plasmodium falciparum* (NIBSC code: 04/176), Zika Virus (PEI code: 11468/16) and Human Cytomegalovirus (NIBSC code: 09/162), each one in the appropriate clinical matrix (whole blood, serum, and plasma, respectively).

	Site	Workflow	VIASURE assay
1	Certest Biotec S.L. (Zaragoza, Spain)	VIASURE Blood Pathogens Extraction Kit using the KingFisher™ Flex Purification System (Thermo Scientific™) + VIASURE V-Lab96 Cyclers (Certest Biotec S.L.) + VIASURE assays	VIASURE <i>Malaria</i> Real Time PCR Detection Kit
			VIASURE <i>Malaria differentiation</i> Real Time PCR Detection Kit
			VIASURE <i>Zika Virus</i> Real Time PCR Detection Kit
			VIASURE <i>Zika, Dengue & Chikungunya</i> Real Time PCR Detection Kit
			VIASURE CMV q Real Time PCR Detection Kit
			VIASURE <i>Methicillin-resistant Staphylococcus aureus</i> Real Time PCR Detection Kit
			VIASURE <i>P. aeruginosa, K. pneumoniae & P. mirabilis</i> Real Time PCR Detection Kit

Table 17. Site, workflow, and VIASURE assay.

True positive and negative values, false positive and negative values, sensitivity, specificity, PPV, and NPV values for VIASURE Real Time PCR Detection Kit were calculated in relation to initial characterization as shown in the following table:

Site	VIASURE assay	Sample matrix	Target	Overall agreement	TP	TN	FP	FN	Sensitivity	Specificity	PPV	NPV
1	VIASURE Malaria Real Time PCR Detection Kit	Whole blood	<i>Plasmodium falciparum</i>	1 (0.94-1)	62	17	0	0	1 (0.93-1)	1 (0.77-1)	1 (0.93-1)	1 (0.77-1)
	VIASURE Malaria differentiation Real Time PCR Detection Kit	Whole blood	<i>Plasmodium falciparum</i>	1 (0.94-1)	62	17	0	0	1 (0.93-1)	1 (0.77-1)	1 (0.93-1)	1 (0.77-1)
	VIASURE Zika Virus Real Time PCR Detection Kit	serum	ZIKV	1 (0.96-1)	67	20	0	0	1 (0.93-1)	1 (0.80-0.99)	1 (0.93-0.99)	1 (0.80-1)
	VIASURE Zika, Dengue & Chikungunya Real Time PCR Detection Kit	serum	ZIKV	1 (0.96-1)	67	20	0	0	1 (0.93-1)	1 (0.80-1)	1 (0.93-1)	1 (0.80-1)
	VIASURE CMV q Real Time PCR Detection Kit	plasma	CMV	1 (0.96-1)	70	12	0	0	1 (0.94-1)	1 (0.70-1)	1 (0.94-1)	1 (0.70-1)
	VIASURE Methicillin-resistant <i>Staphylococcus aureus</i> Real Time PCR Detection Kit	blood culture	MSSA. (Methicillin-sensitive <i>S. aureus</i>)	1 (0.87-1)	4	22	0	0	1 (0.40-1)	1 (0.82-1)	1 (0.40-1)	1 (0.82-1)
		blood culture	MRCoNS (Methicillin-resistant Coagulase-negative <i>Staphylococci</i>)	1 (0.87-1)	10	16	0	0	1 (0.66-1)	1 (0.76-1)	1 (0.66-1)	1 (0.76-1)
	VIASURE <i>P. aeruginosa</i> , <i>K. pneumoniae</i> & <i>P. mirabilis</i> Real Time PCR Detection Kit	blood culture	<i>P. aeruginosa</i>	1 (0.87-1)	2	24	0	0	1 (0.20-1)	1 (0.83-1)	1 (0.20-1)	1 (0.83-1)
		blood culture	<i>K. pneumoniae</i>	1 (0.87-1)	4	22	0	0	1 (0.40-1)	1 (0.82-1)	1 (0.40-1)	1 (0.82-1)

Table 18. True positive (TP) and negative values (TN), false positive (FP) and false negative (FN) values, sensitivity, specificity, Predictive Positive Values (PPV) and Predictive Negative Values (NPV) of VIASURE Malaria Real Time PCR Detection Kit, VIASURE Malaria differentiation Real Time PCR Detection Kit, VIASURE Zika Virus Real Time PCR Detection Kit, VIASURE Zika, Dengue & Chikungunya Real Time PCR Detection Kit, VIASURE CMV q Real Time PCR Detection Kit, VIASURE Methicillin-resistant *Staphylococcus aureus* Real Time PCR Detection Kit, and VIASURE *P. aeruginosa*, *K. pneumoniae* & *P. mirabilis* Real Time PCR Detection Kit.

The results obtained during the performance evaluation show high agreement using VIASURE Blood Pathogens Extraction Kit using the KingFisher™ Flex Purification System (Thermo Scientific™) in combination with VIASURE Real Time PCR Detection Kits and VIASURE V-Lab96 Cycler.

15. Troubleshooting

Issues/Symptoms/Error Warnings/Observation	Possible Causes and Suggestions
Poor yield or low sensitivity	<p>Degradation of the microorganism during shipping/storage and/or processing.</p> <p>Wrong loading of the reagents into the platform. It is recommended to label the reagents 96 deep Well plates to avoid confusion. Preloaded reagents plates must be introduced by locating the first well (A1) in the up-right corner of the rotatory module positions, and in the required order (follow platform user manual).</p> <p>Incomplete sample lysis</p> <ul style="list-style-type: none"> - Inadequate sample pretreatment. Whole blood must be pretreated to be liquid enough for pipetting. Please consult section 9.2 Sample pretreatment. - No addition or addition of reduced volume of Lysis Buffer, Proteinase K and/or sample. - Not enough heating during the incubation step may reduce the efficiency of lysis and insufficient nucleases inactivation. <p>Incomplete nucleic acid adsorption to beads.</p> <ul style="list-style-type: none"> - Not enough agitation for Magnetic Beads complete resuspension before addition to the plate. - Magnetic Beads degraded because of wrong storage. Do not mix remaining magnetic beads between kits and storage them as recommended. - Incorrect resuspension of the Carrier. <p>Not appropriate washing</p> <ul style="list-style-type: none"> - Not enough volume of ethanol to resuspend Wash 1 Buffer and/or Wash 2 Buffer, may cause the desorption of nucleic acids from magnetic beads or dragging of inhibitory substances to downstream applications. Add the recommended quantity of ethanol to each bottle, and once rehydrated store it as recommended above. <p>Elution</p> <ul style="list-style-type: none"> - Nucleic acid degradation because of contaminated consumables. The elution vials or elution plate should be sterile and free from RNAses. - Wrong storage of the purified nucleic acids could cause their degradation. - Contamination of eluted nucleic acid with ethanol before analysis. Avoid platform decontamination with the opened vials or bottles containing nucleic acid nearby.

Malfunction of the automatic platform	<ul style="list-style-type: none"> - Make sure the labware has been placed in the indicated position before starting an extraction run. - Wrong heating and/or shaking. Perform periodical platform maintenance reviews. Abrasive cleaning agents are not recommended, because they are likely to damage the paint finish. For cleaning use a soft cloth or tissue paper. - Run stop because of wrong fasten tip-comb. Try first to manually stretch the tips both lengthwise and width wise to level the tip comb.
Reagent's degradation	<ul style="list-style-type: none"> - Wrong storage of the Extraction Kit reagents before opening. - Wrong storage after the rehydration of the reagents. - Wrong storage of magnetic beads or insufficient agitation before using.
False positive results in qPCR assay	<p>Cross contamination</p> <ul style="list-style-type: none"> - During sample loading, to avoid caps sample mixing, is recommended to label caps and keep them separately among them. - During Proteinase K or Magnetic Beads on the mixture of Lysis Buffer and sample. Change pipette tip among samples. - Cross-contamination with the qPCR VIASURE Positive Controls (in case a qPCR VIASURE Positive Control is used). - Contamination due to qPCR products or samples from previous processing. It is highly recommended to perform properly the end-of-the-day protocol and clean with ethanol to eliminate microbiological risk and contamination sources. - Contamination with ubiquitous microorganisms. It is recommended to decontaminate the platform surface with ethanol after every run. It is highly recommended to perform properly the end-of-the-day protocol and clean with ethanol to eliminate microbiological risk and contamination sources. - Contamination due to Extraction Control (EC) included in this kit on samples to be analyzed by qPCR for <i>E.Coli</i> and/or <i>TEM</i> gene (β-lactamases resistance).
False negative	<p>Analysis of samples with a microbial concentration under limit of detection (LoD).</p> <p>False negative results because of loss of sensitivity.</p> <p>Incorrect run programming of the extraction by user</p>
Variation of fluorescence values in the qPCR assay	They may vary due to multiple factors such as: PCR equipment, type of sample, previous treatment of the sample, etc... among others.

Table 5. Troubleshooting.

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- **The American Type Culture Collection ("ATCC®").**
 - *Klebsiella pneumoniae* (ATCC® BAA-2146™)
 - *Neisseria meningitidis* (ATCC®13091™)
 - Mumps Virus (ATCC®VR-106™)
- **National Institute for Biological Standards and Control (NIBSC).**
 - 1st WHO International Standard for Plasmodium falciparum DNA Nucleic Acid Amplification Techniques (NIBSC code: 04/176).
 - 1st WHO International Standard for Human Cytomegalovirus for Nucleic acid Amplification techniques (NIBSC code: 09/162).

- 1st WHO International Standard for HHV -6B virus DNA Nucleic Acid Amplification Techniques (NIBSC code: 15/266).
- 1st WHO International Standard for BK virus DNA Nucleic Acid Amplification Techniques (NIBSC code: 14/212).
- **The Paul-Ehrlich-Institut (PEI).**
 - 1st WHO International Standard for Zika Virus RNA (PEI code: 11468/16).
 - 1st WHO International Standard for Chikungunya virus RNA Nucleic Acid Amplification Techniques (PEI code: 11785/16).
 - PEI Reference Preparation WNV Lineage 1 RNA (#13299/19).
- **bei RESOURCES.**
 - Dengue virus type 1 (DEN-1) strain Hawaii (bei RESOURCES Ref: NR-82 (derived from ATCC® VR-1254™).
 - *Leishmania donovani* strain 1S (MHOM/SD/62/1S) (bei RESOURCES Ref: NR-48821).

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

4. Thatcher, S. A. (2015). DNA/RNA preparation for molecular detection. *Clinical Chemistry*, 61(1), 89–99. <https://doi.org/10.1373/clinchem.2014.221374>
5. Vila, J., Gómez, M. D., Salavert, M., & Bosch, J. (2017). [Methods of rapid diagnosis in clinical microbiology: Clinical needs]. *Enfermedades infecciosas y microbiología clínica*, 35(1), 41–46. <https://doi.org/10.1016/j.eimc.2016.11.004>.
6. World Health Organization (WHO). (2018). Antimicrobial. <https://www.who.int/news-room/fact-sheets/detail/antimicrobial-resistance>
7. Abdollahi, H., Savari, M., Zahedi, M. J., Moghadam, S. D., & Hayatbakhsh Abasi, M. (2011). Detection of A2142C, A2142G, and A2143G Mutations in 23s rRNA Gene Conferring Resistance to Clarithromycin among *Helicobacter pylori* Isolates in Kerman, Iran. *Iranian Journal of Medical Sciences*, 36(2), 104–110
8. Ali, N., R. de C. P. Rampazzo, A. D. T. Costa, and M. A. Krieger. 2017. Current Nucleic Acid Extraction Methods and Their Implications to Point-of-Care Diagnostics. *Biomed Res. Int.* 2017:9306564.
9. Abdollahi, H., Savari, M., Zahedi, M. J., Moghadam, S. D., & Hayatbakhsh Abasi, M. (2011). Detection of A2142C, A2142G, and A2143G Mutations in 23s rRNA Gene Conferring Resistance to Clarithromycin among *Helicobacter pylori* Isolates in Kerman, Iran. *Iranian Journal of Medical Sciences*, 36(2), 104–110
10. Di Terlizzi, R., & Platt, S. (2006). The function, composition and analysis of cerebrospinal fluid in companion animals: part I - function and composition. *Veterinary Journal* (London, England: 1997), 172(3), 422–431. <https://doi.org/10.1016/j.tvjl.2005.07.021>
11. Emaus, M. N., M. Varona, D. R. Eitzmann, S.-A. Hsieh, V. R. Zeger, and J. L. Anderson. 2020. Nucleic acid extraction: Fundamentals of sample preparation methodologies, current advancements, and future endeavors. *TrAC Trends Anal. Chem.* 130:115985.
12. Guo, Z., Y. Wang, A. Yang, and G. Yang. 2016. The effect of pH on charge inversion and condensation of DNA. *Soft Matter. The Royal Society of Chemistry* 12:6669–6674.
13. Jurado R, Walker HK. Cerebrospinal Fluid. In: Walker HK, Hall WD, Hurst JW, editors. *Clinical Methods: The History, Physical, and Laboratory Examinations*. 3rd edition. Boston: Butterworths; 1990. Chapter 74. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK398/>
14. Islam, M., A. Aryasomayajula, and P. Selvaganapathy. 2017. A Review on Macroscale and Microscale Cell Lysis Methods. *Micromachines* 8:83.
15. MacVane, S. H., & Nolte, F. S. (2016). Benefits of Adding a Rapid PCR-Based Blood Culture Identification Panel to an Established Antimicrobial Stewardship Program. *Journal of Clinical Microbiology*, 54(10), 2455–2463. <https://doi.org/10.1128/JCM.00996-16>.
16. Marco, F. (2017). Métodos moleculares para el diagnóstico de septicemia. *Enfermedades Infecciosas y Microbiología Clínica*, 35(9), 586–592. <https://doi.org/10.1016/j.eimc.2017.03.002>
17. Niu, M., M. Du, Z. Gao, C. Yang, X. Lu, R. Qiao, and M. Gao. 2010. Monodispersed Magnetic Polystyrene Beads with Excellent Colloidal Stability and Strong Magnetic Response. *Macromol. Rapid Commun.* 31:1805–1810.
18. Rittich, B., and A. Spanová. 2013. SPE and purification of DNA using magnetic particles. *J. Sep. Sci. Germany* 36:2472–2485.
19. Schrader, C., A. Schielke, L. Ellerbroek, and R. Johne. 2012. PCR inhibitors - occurrence, properties and removal. *J. Appl. Microbiol. England* 113:1014–1026.
20. Strich, J. R., Heil, E. L., & Masur, H. (2020). Considerations for Empiric Antimicrobial Therapy in Sepsis and Septic Shock in an Era of Antimicrobial Resistance. *The Journal of Infectious Diseases*, 222(Suppl 2), S119–S131. <https://doi.org/10.1093/infdis/jiaa221>

21. Wang, H.-Y., Kim, S., Kim, J., Park, S.-D., Uh, Y., & Lee, H. (2014). Multiplex real-time PCR assay for rapid detection of methicillin-resistant staphylococci directly from positive blood cultures. *Journal of Clinical Microbiology*, 52(6), 1911–1920. <https://doi.org/10.1128/JCM.00389-14>
22. Wu, D.-W., Li, Y.-M., & Wang, F. (2017). How Long can we Store Blood Samples: A Systematic Review and Meta-Analysis. *EBioMedicine*, 24, 277–285. <https://doi.org/10.1016/j.ebiom.2017.09.024>
23. Zhang, X., M. R. Servos, and J. Liu. 2012. Surface Science of DNA Adsorption onto Citrate-Capped Gold Nanoparticles. *Langmuir* 28:3896–3902.
24. Surveillance Guidelines for Measles, Rubella and Congenital Rubella Syndrome in the WHO European Region. Geneva: World Health Organization; 2012 Dec. Annex 3, Collection, storage and shipment of specimens for laboratory diagnosis and interpretation of results. Available from: <https://www.ncbi.nlm.nih.gov/books/NBK143256/>

Symbols for IVD components and reagents/Símbolos para reactivos y productos para diagnóstico *in vitro*

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Trademarks

Modification rights reserved. All rights reserved. © Certest Biotec, S.L.

All other trademarks that may appear in this package insert are the property of their respective owners.

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