# VIASURE

Extraction Kit

## DNA/RNA Pathogens

CE IVD

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Table A 1. References

## Content

1. Intended use .....	4
2. Safety information .....	4
3. Summary and Explanation .....	5
4. Principle of the procedure .....	6
5. Reagents provided .....	7
6. Reagents and equipment to be supplied by the user .....	8
7. Transport and storage conditions .....	8
8. Warnings and Precautions .....	9
9. Test procedure .....	10
10. Result interpretation .....	15
11. Limitations of the procedure .....	15
12. Quality control .....	15
13. Performance characteristic .....	16
14. Troubleshooting .....	16

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## ENGLISH

### 1. Intended use

The VIASURE DNA/RNA Pathogens Extraction Kit is designed for the isolation and purification of viral, bacterial, parasitic, and fungal DNA and RNA from different clinical samples (swabs, saliva, sputum, feces and urine) for subsequent diagnostic purposes using an automated extraction System. The kit contains different reagents to process clinical samples and facilitate the sensible detection of pathogens with downstream application. This product is not intended for the diagnosis, prevention, or treatment of a disease.

This Kit is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in molecular techniques and diagnostic procedures, and experience with potentially infectious human samples or materials.

### 2. Safety information

VIASURE DNA/RNA Pathogens Extraction Kit contains hazardous contents. For more information, please consult appropriate material safety data sheets (MSDS).

#### GSH classification

Reagent/Material	GHS symbol	Hazard phrases	Precaution phrases
Lysis Buffer	 Inflamable, corrosive y toxic	Danger H226, H302, H314, H412	P210, P260, P264, P280, P305+P351+P338, P310, P370+P378
Wash 1 Buffer	 Corrosive and toxic	Danger H302, H314,	P260, P264, P280, P303+P361+P353, P305+P351+P338, P310, P501
Proteinase K	 Danger	Danger H334	P264, P280, P304+340, P342+311

Table 1. GHS classification

#### Hazard phrases

H226; Flammable liquid and vapour.

H302; Harmful if swallowed.

H314; Causes severe skin burns and eye damage.

H412; Harmful to aquatic life with long lasting effects.

H334; May cause allergy or asthma symptoms or breathing difficulties if inhaled.

### **Precaution phrases**

P210; Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P260; Do not breathe dust/fume/gas/mist/vapours/spray.

P264; Wash hands thoroughly after handling.

P280; Wear protective gloves/protective clothing/eye protection/face protection.

P304+340 IF INHALED; Remove person to fresh air and keep comfortable for breathing.

P305+P351+P338; IF IN EYES; Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P303+P361+P353; IF ON SKIN (or hair); Take off immediately all contaminated clothing. Rinse skin with water [or shower].

P310; Immediately call a POISON CENTER or doctor/physician.

P370+P378; In case of fire: Use ... for extinction.

P501; Dispose of contents/container to ...

P342+P311; If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician.

## **3. Summary and Explanation**

Infectious diseases are communicable disorders caused by pathogenic microorganisms such as bacteria, viruses, parasites, or fungi that can be spread directly or indirectly from one individual to another, with negative economic and social consequences all over the world. The correct and fast diagnosis of infections is crucial to manage the pathology, select the appropriate treatment, reduce morbidity and mortality, as well as prevent or control infectious diseases outbreaks.

The diagnosis of pathogens is traditionally based on cultures and time-consuming methods, which have special requirements. The signs and symptoms of infectious diseases are superimposed; therefore, an etiological diagnosis based only on the clinical symptoms is not reliable. The development of molecular techniques based on nucleic acid amplification (NAAT) are a sensitive and specific diagnostic tool for determining more than one pathogen in one single assay. In this sense, real time PCR is the preferred molecular diagnostic method for infections.

Nevertheless, this technique requires a previous nucleic acid (NA) extraction. The procedure to do so includes three main steps: extracting or releasing the NA from the sample and the cell or organism by lysis; separating or isolating NA from other cellular or sample material; and purification by removal of inhibitory substances. A fourth optional step, concentrating, is important for detection of low-concentration analytes.

An extraction kit for the simultaneous isolation of high quality and quantity DNA and/or RNA, from different clinical samples, and pathogens, is essential for their specific detection through molecular diagnostic assays. Nowadays, the increased number of samples to be processed at clinical laboratories requires the possibility of automating this procedure. The adaptation of the reagents to a robotic platform can improve daily workflow, standardization of sample treatment, decrease variability, reduce biohazards, and avoid contamination by reducing the analyst final hands-on time.


VIASURE DNA/RNA Pathogens Extraction Kit is designed for the isolation and purification of viral, bacterial, parasites and fungal pathogens nucleic acids from different human clinical samples such as swabs, saliva, sputum, urine, and feces.

## 4. Principle of the procedure

The VIASURE DNA/RNA Pathogens Extraction Kit procedure is based on the spontaneous and reversible adsorption of nucleic acids to superparamagnetic beads with a functionalized surface under appropriate buffer conditions. After some washing steps the nucleic acids are eluted in a purified and concentrated solution. The use of magnetic particles as solid phase for adherence of nucleic acids allows the processing of an elevated number of samples by their simultaneous exposition to a magnetic field, and therefore the automation of the process.

VIASURE DNA/RNA Pathogens Extraction Kit has been designed for the automatic platform KingFisher™ Flex (ThermoScientific). This type of platform procedure transfers the magnetic beads between tube/wells from one reagent to the next, by magnetic field attraction. The activation and de-activation of the magnetic field allows beads attraction and beads resuspension within different buffers. Figure 1 shows the basic principle of this type of automated extraction. First, the Magnetic Beads are added to the lysate and after an incubation, the adsorption of the nucleic acids occurs (1). The magnetic particles are attracted to the surface of a magnet covered with plastic tips (2), to be transferred to other wells containing Wash Buffer 1 y Wash Buffer 2, to eliminate cell debris and salts. After some washing steps (3), beads are dried to eliminate organic solvent residues (4). Finally, the DNA or RNA is released from the magnetic particles with an Elution Buffer (5), which results in a pure sample ready for quantification and analysis. Purified DNA and RNA can be directly used for downstream applications. It is highly recommended to use extraction controls, positive/negative controls, internal controls to evaluate the complete workflow (See Procedure sections).

Figure 1. Magnetic beads extraction basic principle for a magnetic beads handling automated instrument.

VIASURE DNA/RNA Pathogens Extraction Kit facilitates the simultaneous processing of different samples in the same extraction run, increases daily workflow, robustness and reproducibility, reduces contamination, and contact with biological hazards.

## 5. Reagents provided

VIASURE DNA/RNA Pathogens Extraction Kit includes the following reagents detailed in Table 2.

Reagent/Material	Comments	Description	Amount
Lysis Buffer *	Liquid reagents mixture for pathogen lysis and nucleases inactivation	White bottle	1 bottle x 65 ml
Wash 1 Buffer	Wash solution in which absolute ethanol (not included) must be added	White bottle	1 bottle x 27 ml
Wash 2 Buffer	Wash solution in which absolute ethanol (not included) must be added	White bottle	1 bottle x 18 ml
Elution Buffer	Solution for nucleic acid recovery	White bottle	1 bottle x 50 ml
Magnetic Beads	Magnetized beads suspension for nucleic acid adsorption	Transparent vial with cap	1 vial x 1.6 ml
Proteinase K	Enzyme solution for pathogen lysis	Transparent vial with cap	1 vial x 1.6 ml
Extraction Control (EC)	Non-infectious nucleic acid lyophilized	Green vial with cap	1 vial
Water RNase/DNase Free	Rehydration buffer for Extraction Control vial resuspension	White vial with cap	1 vial x 1 ml

\* Contains chaotropic salt. Not compatible with disinfectants containing bleach. See section 8 "Warning and precautions".

Table 2. Reagents and materials provided in VIASURE DNA/RNA Pathogens Extraction Kit.

## 6. Reagents and equipment to be supplied by the user

The following list includes the materials that are required for use but not included in the VIASURE DNA/RNA Pathogens Extraction Kit.

- Sample collection and transport system.
- Absolute ethanol.
- Micropipettes.
- Filter tips.
- Powder free disposable gloves.
- Biosafety cabinet.
- Pipettes and/or graduate test tube.
- Vortex.
- Thermomixer.
- Real-Time PCR instrument (thermocycler).
- Real-Time PCR assay.
- Ribonuclease (RNase)/ Deoxyribonuclease (DNase) free tubes for nucleic acid collection and storage.
- Positive and negative controls for the monitoring the Extraction process. It is recommended the use of positive and negative extraction controls. (For example, to test SARS-CoV-2, the kit VIASURE Viral SARS-CoV-2 Positive Control Kit (Ref: VS-VP1NCO or VS-VP101NCO)).

*Note: In case that processes other than those indicated in this Instructions for use are validated, it is recommended to assess whether it is necessary to use materials other than those described in this section.*

The specific materials required when using the for the KingFisher™ Flex (ThermoScientific) system are listed below:

- KingFisher™ deep-well 96 plate - ThermoScientific™ (Number catalogue: 95040450) Units: 4/96 preps
- KingFisher™ 96 plate 200µL - ThermoScientific™ (Number catalogue: 97002540) Units: 1/96 preps
- KingFisher™ 96 tip comb for DW Magnets - ThermoScientific™ (Number catalogue: 97002534) Units: 1/96 preps

## 7. Transport and storage conditions

- The kit reagents can be shipped at 2-40 °C and stored until the expiration date which is stated on the label. Check all components for damages after receiving the kit. Do not use damaged kit components.
- It is recommended to keep components at room temperature (15-30 °C) away from light and store in a vertical position to avoid any possible spillage.
- It is recommended to keep Magnetic Beads vials under refrigeration 4°C-8°C after receiving the kit.
- Lysis Buffer and Wash 1 Buffer contain guanidine salts that may precipitate. If there is visible precipitate at the bottom of the cartridge, resuspend by mild heating and shake.
- Once the Extraction Control has been re-suspended, store it at -20°C.
- It is recommended to store the supernatant obtained in the processing at 4°C for up to 8 hours until its addition to the different downstream applications. This addition should be performed as soon as possible. After use, the remanent of the supernatant can be stored at -20°C minimizing freeze-thaw cycles in the case of retest (up to 2 times). If sample treatment and amplification assay cannot be performed on the same day, the supernatant obtained must be kept at -20°C until use (up to one week).

## 8. Warnings and Precautions

- The product is intended for use by qualified and trained clinical laboratory personnel specifically instructed and trained in the techniques of nucleic acid extraction procedures and real-time PCR and in vitro diagnostic procedures, including training on nucleic acid extraction system and the real-time PCR instrument (thermocycler).
- Do not use past expiration date.
- Do not use the kit if the label that seals the outer box is broken.
- Do not use reagents if the protective external box or the protective pouches are opened or damaged upon arrival.
- Do not use Extraction Control if desiccant is not present or broken inside reagent pouch. Do not remove desiccant from reagent pouch once is open.
- Do not mix reagents from different pouches and / or kits and / or lots and / or another supplier.
- Design a unidirectional workflow. Do not reuse reagents and do not recover the remaining liquid volume into the original bottle of the reagent once used.
- Always avoid microbial and ribonuclease (RNase)/deoxyribonuclease (DNase) contamination of reagents. The use of sterile RNase/DNase-free disposable aerosol resistant or positive displacement pipette tips is recommended. Use a new tip for each specimen. Gloves must be changed before manipulating reagents. In cases where other extraction procedures are performed in the same general area of the laboratory, care must be taken to ensure that the VIASURE DNA/RNA Pathogens Extraction Kit and any additional reagents or equipment required for testing are not contaminated.
- Do not reuse any disposable items.
- Store positive materials used separated from the kit's reagents.
- Use screw-capped tubes and prevent any potential splashing or cross contamination of specimens during preparation.
- Avoid possible contamination of reagents with extracted nucleic acids, PCR products, and positive control. To prevent contamination of reagents, use of filter-tips is recommended.
- Follow Good Laboratory Practices. Wear protective clothing, use disposable gloves, goggles, and mask. Do not eat, drink or smoke or apply cosmetic products in the working area. Once you finish the test wash your hands.
- Specimens must be treated as potentially infectious and/or biohazardous, as well as all the reagents and materials that have been exposed to the samples and they must be handled according to the national safety regulations. Take necessary precautions during the collection, transport, storage, handling, and disposal of samples.
- Samples and reagents must be handled in a biological safety cabinet. Use personal protective equipment (PPE) consistent with current guidelines for the handling of potentially infectious samples. Dispose of waste in compliance with local, state, and federal regulations.
- All reagents should be handled with care, minimizing the time of exposure to the open, especially those reconstituted with absolute ethanol. Consult the hazards and precautions when handling ethanol in the Material Safety Data Sheet.
- Regular decontamination of commonly used equipment is recommended, especially micropipettes and work surfaces.
- When working with the VIASURE DNA/RNA Pathogens Extraction Kit wear suitable protective clothing (e.g., lab coat, disposable gloves, and protective goggles).
- Reliability of the results depends on adequate specimen collection, storage, transport, and processing procedure. Synthetic fibre swabs with plastic shafts are recommended for sample collection. Calcium alginate swabs or

swabs with wooden shafts are not recommended, as they may contain substances that inhibit PCR testing. Make sure that the swab head is completely submerged in the sample collection tube and if necessary, press the swab against the wall of the tube to squeeze out the liquid.

- Use graduated test tubes to fill the bottles with isopropanol or ethanol.
- In the case of a subsequent qPCR analysis, consult each qPCR kit or real-time PCR instrument's reference manual for additional warnings, precautions, and procedures.
- The waste generated with the VIASURE DNA/RNA Pathogens Extraction Kit has not been tested for residual infectious material. A contamination of the liquid waste with residual infectious material is highly unlikely due to strong denaturing Lysis Buffer and Proteinase K treatment but it cannot be excluded completely. Therefore, liquid waste must be considered infectious and should be handled and discarded according to local safety regulations.
- Lysis Buffer contains chaotropic salt. Not compatible with disinfectants containing bleach. Consult the hazards and precautions when handling Lysis Buffer in the Material Safety Data Sheets provided with the Device.
- Make sure there is not precipitate in the supplied reagents. Lysis Buffer and Wash 1 Buffer contain guanidine salts that may precipitate. If there is visible precipitate at the bottom of the bottle, resuspend by mild heating and shake until its complete dissolution.
- If the Proteinase K present a yellowish colour, note that this does not affect the functionality of the kit.
- Make sure that the beads are completely resuspended. Shake vigorously the storage vial on a Vortex.
- Make sure all tubes, plates or cartridges have been properly placed following the indicated layout and fit correctly into the specific adapters of the Automated extraction system.
- Before dispensing VIASURE DNA/RNA Pathogens Extraction Kit into the platform consumables open carefully each reservoir avoiding reagent spilling and reconstitute the indicated wells as indicated in the procedure section.

## 9. Test procedure

### 9.1. Sample collection and transport

VIASURE DNA/RNA Pathogens Extraction Kit is suitable for various potentially infectious human clinical samples: respiratory and sexual swabs samples, saliva, sputum, urine, and feces. For reproducible and high yields appropriate sample storage is essential. Yields may vary from sample to sample depending on factors such as health and age of the donor, kind of sample, collection date, transport, and storage conditions.

As general recommendations, the specimens should be collected as soon as possible, depending on the pathogen during the acute phase of illness or at any time during the clinical course, and before antimicrobial therapy begin, if possible. Specimens should be collected in a sterile container avoiding touching adjacent tissues and organ secretions that are not of interest. All specimens must be labeled with the name and identification number of the person from whom the sample was collected, the source of the specimen, and the date and time it was collected. After collection, specimens should be placed in a biohazard bag and transported to the laboratory as soon as possible.

Depending on the type of specimen, the collection differs. Follow the recommendations by McPherson R.A., Pincus M.R. Henry's Clinical Diagnosis and Management by Laboratory Methods. 23RD edition. Elsevier Health Sciences; 2017.

According to the literature, the conditions of transport, storage and stability depend on the type of the sample matrix and the pathogen to be detected. It is recommended to use fresh specimens for the test. The specimens should be refrigerated for urine, respiratory specimens, and stool; or held at room temperature for cerebrospinal fluid (CSF), blood

and other body fluids. If a delay is unavoidable (< 2 hours until they can be sent), the general rule is to ship at -20°C or ideally at -70°C or 80°C.

Optimally, all specimen containers should be opened in a biological safety cabinet. Appropriate barriers should be used to prevent exposure of skin and mucous membranes to the specimen. Gloves and a lab coat must always be worn when handling clinical specimens. Masks, goggles (or working behind a plastic shield), and impermeable gowns or aprons must be worn when there is a risk for splashes or droplet formation. Wash your hands once the proceeding has finished.

Specimens submitted for molecular testing must be stored in controlled conditions so that nucleic acids do not degrade during storage. Avoid freezing and thawing cycles. Some pathogens from specimens that are frozen and then thawed could be degraded and may result in false-negative test results.

Depending on the clinical specimen and pathogen different sample pretreatments may be needed. Follow the appropriate protocol laboratory guideline in each case. Refer to the CDC guideline (Specimen collection guidelines. Website <https://www.cdc.gov/urdo/downloads/SpecCollectionGuidelines.pdf> and Interim Guidelines.

For Collecting, Handling, and Testing Clinical Specimens for COVID-19 consult: Website <https://www.cdc.gov/coronavirus/2019-ncov/lab/guidelines-clinical-specimens.html>) and the IDSA guideline (Miller, J. M., Binnicker, M. J., Campbell, S., ... & Pritt, B. S. (2018). A guide to utilization of the microbiology laboratory for diagnosis of infectious diseases: 2018 update by the Infectious Diseases Society of America and the American Society for Microbiology. *Clinical Infectious Diseases*, 67(6), e1-e94).

It is recommended to work with inactivated clinical samples. VIASURE DNA/RNA Pathogens Extraction Kit is compatible with inactivating transport media such as guanidinium salts media.

## 9.2. Sample pretreatment

Depending on the clinical specimen and pathogen different sample pretreatments may be needed. It is important to avoid substances or objects that may interfere in sample pipetting during the automated process, therefore, caps or swabs should be removed before introducing the samples into the sample racks, and solid or highly viscous samples must be pretreated.

- Feces: Mix vigorously approximately 100 mg with 900 µL of distilled water, then centrifuge 90 seconds at 1000 g. Recover supernatant with the suspected pathogen for extraction.
- Urine. Centrifuge 10 mL at 5200 g for 15 min, remove supernatant and resuspend pellet in 1 mL of TE 10X Ph
- Viscous samples such as saliva, sputum or bronchoalveolar lavages may require pretreatment with acetyl cysteine 20 mg/mL. Add as many volume as required to make sample is liquid enough to be aspirate with the pipette. If its needed mix and heat during pre-treatment.

Rest of samples can be directly introduced into the platform.

## 9.3. Buffer's addition

Before starting with the nucleic acid extraction and purification, it is necessary in some reagents to add a certain volume of other reagents (not included in the kit). The following table indicates which reagent should be added and the required volume:

Reagent/Material	Reagent for resuspension	Volume
Wash 1 Buffer*	Absolute ethanol	63 mL
Wash 2 Buffer*	Absolute ethanol	72 mL

Table 3. Reagents and volumes to be resuspended in the Wash 1 Buffer and Wash 2 Buffer.

\*Note: To facilitate the process, a mark has been included on the label of each reagent to tick when the reagents are added.

Besides, it is recommended to make sure to write with indelible marker the specific automated extraction system cartridges used once resuspended in the space with the resuspension date and other required information.

Use graduated test tubes to fill the bottles with ethanol.

Lysis Buffer contain guanidine salts that may precipitate. If there is visible precipitate at the bottom of the bottle, resuspend by mild heating and shake.

#### 9.4. Lyophilized extraction control

The use of the Extraction Control (EC) will be required only for VIASURE Real Time PCR products which do not include an internal control in the lyophilized reaction mix. Please consult Instructions for use of each VIASURE product to check this information. For these VIASURE Real Time PCR products use the Extraction Control (EC) provided in the VIASURE DNA/RNA Pathogens Extraction Kit.

Recommendation is to open and manipulate the Extraction Control (EC) in pre-PCR laboratory area away from the PCR-Amplification area. Reconstitute the lyophilized Extraction Control (green vial) adding 550 µL of Water RNase/DNase free (white vial) supplied and vortex thoroughly.

Once the Extraction Control has been re-suspended, store it at -20°C and minimize freeze and thaw cycles (up to 3 cycles).

Add 5 µL of the rehydrated EC to each required well with preloaded Lysis Buffer to monitor the extraction process and/or discard the inhibition of the polymerase activity.

#### 9.5. Extraction procedure

The procedure below provides instructions the steps to process samples. The specific extraction procedure to follow within KingFisher® Flex automatic platforms is listed below. (Note: In case of using VIASURE DNA/RNA Pathogens Extraction Kit with another automated extraction system different from the one mentioned in these Instructions for use, the procedure must be validated by the user).

##### 1. Load with reagents 96 Deep well plates

Load all plates with the following reagents and volumes in the required wells depending on the number of samples to analyze.

- Plate 1: Lysis Buffer 265 µL.
- Plate 2: Wash 1 Buffer 500 µL.
- Plate 3: Wash 2 Buffer 500 µL.
- Plate 4: Elution Buffer 100 µL.

## 2. Sample inactivation

Add 200 µl of sample to the preloaded Lysis Buffer plate making sure that the sample has been properly mixed with the Lysis Buffer, and incubate at room temperature for 15 min. If it is needed to control the complete molecular diagnosis workflow add 5 µl of rehydrated Extraction Control before incubation.

## 3. Magnetic Beads and Proteinase K

Add to each sample 15 µl of Proteinase K and 15 µl of Magnetic Beads.

It is recommended to vortex the Magnetic Beads at least 2 minutes until they are completely resuspended.

## 4. Plate Set Up into the instrument.

Configure the protocol of the Instrument following the table below that describes the complete setup of the extraction process:

Step	Procedure	Setting
<b>Pick Up</b>	<b>96 Tip Comb plate</b>	
<b>Binding</b>	<b>Sample plate</b>	
Beginning of step	Precollect of Magnetic Beads	No
Mixing/Heating	Release beads Mixing time, speed Heating temperature Preheat	5 min, Fast 65°C Yes
End of step	Postmix Collect count Collect time	No 5 0 sec
<b>Wash 1</b>	<b>Wash 1 plate</b>	
Beginning of step	Precollect of Magnetic Beads	No
Mixing/Heating	Release time, speed First shake time, speed Second shake time, speed Loop count Heating during mixing	20 sec, Bottom mix 10 sec, Bottom mix 10 sec, Fast 3 No
End of step	Postmix Collect count Collect time	No 4 1 sec
<b>Wash 2</b>	<b>Wash 2 plate</b>	
Beginning of step	Precollect of Magnetic Beads	No
Mixing/Heating	Release time, speed First shake time, speed Second shake time, speed Loop count Heating during mixing	20 sec, Fast 10 sec, Bottom mix 10 sec, Fast 2 No
End of step	Postmix Collect count Collect time	No 4 1 sec
<b>Dry</b>	<b>Wash 2 plate</b>	
	Dry time Tip position	2 min Outside well/tube
<b>Elution</b>	<b>Elution plate</b>	
Beginning of step	Precollect	No

Mixing/Heating	Release beads First shake time, speed Second shake time, speed Loop count Heating temperature	Yes 15 sec, Bottom mix 45 sec, Medium 6 70 °C
End of step	Preheat Postmix Collect cout Collect time	Yes No 5 0 sec
<b>Collection beads</b>	<b>Elution plate</b>	
Beginning of step	Precollect	No
Mixing/Heating	Release beads Mixing time, speed Heating during mixing	No 1 min 30 sec, Slow No
End of step	Postmix Collect beads	No No
<b>Leave</b>	<b>Tip comb plate</b>	

Table 4. Complete settings of the extraction process to be followed within KingFisher® Flex automatic platforms.

Make sure all plates are loaded in the same direction, with the well A1 facing the top right. Then fill the platform with the plates following the next order:

- 1st. Tip comb.
- 2nd. Elution plate.
- 3rd. Wash 2 Buffer plate.
- 4th. Wash 1 Buffer plate.
- 5th. Lysis Buffer with inactivated sample plate.

Once plates have been loaded click start and the run will be done in 28 minutes.

#### Purified nucleic acids recovery

Eluted samples can be used directly or stored at 4°C if the qPCR is performed immediately after the extraction. For long term storage, keep at -20°C or -80°C minimizing freeze-thaw cycles. Consult Section Transport and storage conditions for more information. Once the extraction and PCR setup transfer the loaded plate to thermocycler for nucleic acid amplification following Manufacturer's indications.

#### 5. Decontamination

Decontaminate properly the platform inner surfaces with appropriate disinfectant solvents and discard used plastic consumables.

For more details consult KingFisher™ Flex (ThermoScientific) user manual.

### 9.6 Real-Time PCR

The nucleic acids obtained with VIASURE DNA/RNA Pathogens Extraction Kit are compatible with any commercially available qPCR or RT-qPCR kit, although it has only been validated with the VIASURE Real Time PCR Detection kits.

Manufacturer's instructions for the VIASURE Real Time PCR Detection Kits must be followed (consult <https://www.certest.es/es/viasure/>). The positive and negative controls included in the extraction step must not substituted the positive and negative controls provided with the qPCR kits.

Recommendation is to use the thermocycler VIASURE V-Lab96 Cyclor (CerTest Biotec S.L) and the decision-support software VIASURE V-Smart (CerTest Biotec S.L) for the analysis and interpretation of the results.

## 10. Result interpretation

To correctly interpret the results of the controls and clinical samples, please consult the Manufacturer's Instructions of the different downstream applications.

If the chosen downstream application is qPCR assay, the recommendation is to use the decision-support software VIASURE V-Smart (CerTest Biotec S.L) for the analysis and interpretation of the results.

## 11. Limitations of the procedure

- The efficiency of the performance characteristics of VIASURE DNA/RNA Pathogens Extraction Kit for KingFisher™ Flex (ThermoScientific) depends on the specimens' collection, transport, storage, and/or handling methods and the efficiency of the flowchart used by each laboratory: qPCR assay and thermocycler. Each complete procedure, different from the one mentioned in these Instructions for use, must be validated by the end user.
- For appropriate yield and quantity, sample collection, storage, and pre-treatment must be performed following the specific recommendations of official laboratory guidelines.
- This sample processing system is compatible with transport media containing guanidium thiocyanate or other guanidium salts. Other transport media different from those indicate in this instruction for use, must be validated by the user.
- If the consumable materials used differ from those recommended for this format, they must be validated by user.
- VIASURE DNA/RNA Pathogens Extraction Kit can be automated in the automated platform KingFisher™ Flex (ThermoScientific). Be sure to read this Instruction Manual and the KingFisher™ Flex's User Manual carefully before use. If the automatation is performed in a different magnetic bead handling platform with other consumables, it must be validated by user.
- Once the extraction has finished to ensure nucleic acid viability, keep under refrigeration until use, and at -20/-80°C for long term storages. Consult Section Transport and storage conditions for more information.

## 12. Quality control

The expected values of VIASURE DNA/RNA Pathogens Extraction Kit depend on the diagnostic procedure (extraction and qPCR assay), so each complete procedure must be validated by the final user.

It is recommended to use the VIASURE DNA/RNA Pathogens Extraction Kit in conjunction with appropriate controls to monitor the entire process for different downstream applications. To control the complete process, from sample processing to analysis, perform the following controls: Positive Control (not provided but recommended), Negative Control (not provided but recommended) and Extraction Control (provided with the kit).

### 13. Performance characteristic

The efficiency of the PCR reaction is not altered with this extraction method.

The clinical performance of VIASURE DNA/RNA Pathogens Extraction Kit was tested using VIASURE Real Time PCR Detection Kits with 750 clinical samples (from the respiratory tract (343 samples): oropharyngeal/nasopharyngeal swabs, sputum, bronchoalveolar lavage (BAL), and bronchial aspirate (BAS); from genital tract (248 samples): vaginal swab and vagino-rectal specimens; and from gastrointestinal tract (159 samples): feces) from patients with suspicion of respiratory tract infections (RTIs), genital tract infections or gastrointestinal infection.

The respiratory pathogen panel tested were SARS-CoV-2, Influenza A, RSV, *Streptococcus pneumoniae*, *Haemophilus influenzae*, *Moraxella catarrhalis* and *Pneumocystis jirovecii*; pathogens which caused genital tract infection: Human Papillomavirus and *Streptococcus B*; and gastrointestinal pathogens: *Salmonella* sp., *Campylobacter* sp., *Clostridium difficile*, *Cryptosporidium* sp., *Giardia lamblia*, Rotavirus, Norovirus GII, and Astrovirus.

VIASURE DNA/RNA Pathogens Extraction Kit was compared with MagMAX™ Viral/Pathogen II (MVP II) Nucleic Acid Isolation Kit for KingFisher™ Flex Purification System (Thermo Scientific™). The VIASURE qPCR tests used were the following: VIASURE SARS-CoV-2 Real Time PCR Detection Kit, VIASURE Flu A, Flu B & RSV Real Time PCR Detection Kit, VIASURE *H. influenzae*, *S. pneumoniae* & *M. catarrhalis* Real Time PCR Detection Kit, VIASURE *Pneumocystis jirovecii* Real Time PCR Detection Kit, VIASURE High-Risk Human Papilloma Virus Real Time PCR Detection Kit and VIASURE *Streptococcus B* Real Time PCR Detection Kit, VIASURE *Salmonella*, *Campylobacter* & *Shigella/EIEC* Real Time PCR Detection Kit, VIASURE *Cryptosporidium*, *Giardia* & *E. histolytica* Real Time PCR Detection Kit, VIASURE Rotavirus Real Time PCR Detection Kit, VIASURE Norovirus GII Real Time PCR Detection Kit, VIASURE Astrovirus Real Time PCR Detection Kit and VIASURE *Clostridium difficile* toxins A+B Real Time PCR Detection Kit.

The results obtained during the performance evaluation show high agreement using VIASURE DNA/RNA Pathogens Extraction Kit in combination with VIASURE Real Time PCR Detection Kits and VIASURE V-Lab96 Cycler.

### 14. Troubleshooting

Issues/Symptoms/Error	Possible Causes and Suggestions
Warnings/Observation	
Poor yield or low sensitivity	<p>Degradation of the microorganism during shipping/storage and/or processing.</p> <p>Wrong loading of the reagents into the platform. It is recommended to label the reagents 96 deep Well plates to avoid confusion. Preloaded reagents plates must be introduced by locating the first well (A1) in the up-right corner of the rotatory module positions, and in the required order (follow platform user manual).</p> <p>Incomplete sample lysis</p> <ul style="list-style-type: none"> <li>- Inadequate sample pretreatment. Sputum and feces must be pretreated to be liquid enough for pipetting.</li> </ul>

	<ul style="list-style-type: none"> <li>- No addition or addition of reduced volume of Lysis Buffer, Proteinase K and/or sample.</li> <li>- Not enough heating during the incubation step may reduce the efficiency of lysis and insufficient nucleases inactivation.</li> </ul> <p>Incomplete nucleic acid adsorption to beads</p> <ul style="list-style-type: none"> <li>- Not enough agitation for Magnetic Beads complete resuspension before addition to the plate.</li> <li>- Magnetic Beads degraded because of wrong storage. Do not mix remaining magnetic beads between kits and storage them as recommended.</li> </ul> <p>Not appropriate washing</p> <ul style="list-style-type: none"> <li>- Not enough volume of ethanol to resuspend Wash 1 Buffer and/or Wash 2 Buffer, may cause the desorption of nucleic acids from magnetic beads or dragging of inhibitory substances to downstream applications. Add the recommended quantity of ethanol to each reservoir, and once rehydrated store it as recommended above.</li> </ul> <p>Elution</p> <ul style="list-style-type: none"> <li>- Nucleic acid degradation because of contaminated consumables. The elution vials or elution plate should be sterile and free from RNAses.</li> <li>- Wrong storage of the purified nucleic acids could cause their degradation.</li> <li>- Contamination of eluted nucleic acid with ethanol before analysis. Avoid platform decontamination with the opened vials or reservoirs containing nucleic acid nearby.</li> </ul>
Malfunction of the automatic platform	<ul style="list-style-type: none"> <li>- Make sure the labware has been placed in the indicated position before starting an extraction run.</li> <li>- Wrong heating and/or shaking. Perform periodical platform maintenance reviews. Abrasive cleaning agents are not recommended, because they are likely to damage the paint finish. For cleaning use a soft cloth or tissue paper.</li> <li>- Run stop because of wrong fasten tip-comb. Try first to manually stretch the tips both lengthwise and width wise to level the tip comb.</li> </ul>
Reagent's degradation	<ul style="list-style-type: none"> <li>- Wrong storage of the Extraction Kit reagents before opening.</li> <li>- Wrong storage after the rehydration of the reagents.</li> </ul>
False positive results in qPCR assay	<p>Cross contamination</p> <ul style="list-style-type: none"> <li>- During sample loading, to avoid caps sample mixing, is recommended to label caps and keep them separately among them.</li> <li>- During Proteinase K or Magnetic Beads on the mixture of Lysis Buffer and sample. Change pipette tip among samples.</li> </ul>

	<ul style="list-style-type: none"> <li>- Cross-contamination with the qPCR VIASURE Positive Controls (in case a qPCR VIASURE Positive Control is used).</li> <li>- Contamination due to qPCR products or samples from previous processing. It is highly recommended to perform properly the end-of-the-day protocol and clean with ethanol to eliminate microbiological risk and contamination sources.</li> <li>- Contamination with ubiquitous microorganisms. It is recommended to decontaminate the platform surface with ethanol after every run. It is highly recommended to perform properly the end-of-the-day protocol and clean with ethanol to eliminate microbiological risk and contamination sources.</li> </ul>
False negative	<p>Analysis of samples with a microbial concentration under limit of detection (LoD).</p> <p>False negative results because of loss of sensitivity.</p> <p>Incorrect run programming of the extraction by user</p>
Variation of fluorescence values in the qPCR assay	<p>They may vary due to multiple factors such as: PCR equipment, treatment procedure, type of sample, previous treatment of the sample, etc... among others.</p>

Table 5. Troubleshooting.

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Modification rights reserved. All rights reserved. © CerTest Biotec, S.L.
All other trademarks that may appear in this package insert are the property of their respective owners.