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**【Storage Conditions & Shelf Life】** Stored at 2°C~30°C in a dry place, sealed and protected from sunlight. Refer to the label for the shelf life. After opened, the shelf life is 30 days.

**【For Instrument】** Urine Sediment Analyzer, FUS Series Urinalysis Hybrid and MUS Series Urinalysis System.

**【Usage】** See user manuals of Urine Sediment Analyzer, FUS Series Urinalysis Hybrid and MUS Series Urinalysis System.

**【Sample Requirement】** Fresh urine sample without centrifugation.

**【Performance Indicator】** pH is about  $6.30 \pm 0.20$  at  $(25 \pm 1) ^\circ\text{C}$

**【Important Notice】**

1. Avoid contact with skin and eyes;
2. Rinse eyes with the flow of warm water for a few minutes if the product splashed into eyes (ensuring eyelids opened). Seek medical care immediately if pain or swelling happened. Rinse with soap water and water once skin contact occurs.

**【Bibliography】**

1. Ferris JA Comparison and standards of urine microscopy Lab Med 14:659-662,1983
2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

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 Please read the User Manual	 Stored at	 Manufacturer
 In vitro diagnostic reagent	 Lot Number	 Valid until
 Conform to EU in vitro diagnostic medical equipment directives 98/79/EC		
 Authorised Representative	 Catalogue number	

**【Manual Approved and Modified Date】** 04/2019

**DIRUI INDUSTRIAL CO., LTD.**

95 Yunhe Street, New & High Tech.

Development Zone Changchun, Jilin 130012 P.R.China

Tel: +86 431 85100409 Fax: +86 431 85172581

E-mail: dirui@dirui.com.cn http://www.dirui.com.cn

EMERGO EUROPE Prinsessegracht 20  
2514 AP The Hague The Netherlands

# Urine Sediment Analyzer Detergent

## Instructions

**[Product Name]** Urine Sediment Analyzer Detergent

**[Basic UDI-DI]** 69257945W0101060299R0QZ

**[Package Specification]** 50mL, 300mL, 500mL, 4x500mL

**[Intended Use]**

It is used for cleaning the pipeline and flow cell of Urine Sediment Analyzer, FUS Series Urinalysis Hybrid and MUS Urinalysis System.

**[Test Principle]**

Sodium hypochlorite has a very high oxidation potential, which can quickly oxidize and decompose grease molecules and macromolecular organics, causing them to lose adhesion and achieve the effect of cleaning.

**[Main Components]** 0.5% sodium hypochlorite

**[Storage Conditions]**

Stored at 2-30°C, protected from sunlight and in seal.

**[Shelf Life]**

The shelf life is 12 months if stored at 2-30°C and in seal. The reagent can be stable for 30 days after being opened if stored at 2-30°C.

**[Date of Manufacture & Expiry Date]** See the label.

**[Applicable Instrument]**

DIRUI FUS Series Urinalysis Hybrid, MUS Series Urinalysis System, FUS Series Urine Sediment Analyzer.

**[Test Method]**

Refer to the applicable user manual of instrument.

**[Product Performance Indices]**

1. Appearance: colorless or light yellow liquid without sediments, particles or floccules
2. Loading volume: no lower than the nominal value
3. pH value: at (25±1)°C, pH is 12.10±0.50
4. Between-batch difference: ΔpH≤0.50

**[Warnings and Precautions]**

1. The product is only for in vitro diagnosis and professional use.
2. Avoid contact with eyes, skin and clothes. If you accidentally come into contact with the skin or clothes, please wash the contact area with water immediately. If it gets into eyes accidentally, rinse immediately with plenty of water and seek medical attention. If the reagent is swallowed accidentally, seek medical attention immediately.
3. Avoid high temperature and direct sunlight.
4. Please affix the hazard labels on the product packaging in accordance with European Community regulations.

R31: contact with acid releases toxic gas

R36/38: irritating to eyes and skin

S50: do not mix with acidic solutions

S24/25: prevent contact with skin and eyes

5. Operators should comply with the laboratory safety operation regulation and dispose of them according to the local medical waste, infectious waste, industrial waste regulations, etc.  
6. If any adverse events occur during use, please report to the competent authority of the Member State where the manufacturer is located.

**Safety precautions:**

Hazard Pictograms: No hazard pictogram is used.

Signal Word(S): No signal word is used.

Hazard Statement: Not applicable.

Precautionary statement: Not applicable.

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**[References]**

CN 111808684 A A Cleaning Solution Suitable for Urine Analyzer or Secretion Analyzer and Its Application

**[Symbol Explanation]**

Please read Instructions	Storage temperature	Manufacturer
Catalogue number	Lot No.	Expiry date
Authorized representative	In vitro diagnostic reagent	County of manufacture
European in Vitro Diagnostic Medical Device Regulation(EU) 2017/746(IVDR)	R31 -36/38 550-24/25	Danger level

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1. Ferris JA urine microscopy comparison & standard: Lab Med 14:659-662, 1983.
2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

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**【Product name】** FE Calibrator II

**【Package specification】** 1×10 mL, 1×12 mL, 1×15 mL, 1×60 mL, 1×125 mL.

**【Intended use】** This product is used to calibrate Urine Sediment Analyzer, FUS Series Urinalysis Hybrid and MUS Series Urinalysis System to ensure test accuracy.

**【Test principle】** Flow cytometry imaging

**【Main component】** QC blood: 900 particles/μL~1,300 particles/μL;

**Note:** fixed values of different batches of formed elements analysis calibrator are different (see the label for details); and the traceability of fixed values is from internal reference methods of the company.

**【Storage conditions and shelf life】**

1. This product shall be sealed and kept in a dark place at the temperature of 2°C~8°C. For shelf life refer to the label.

2. After being unsealed for the first time, it shall be sealed and kept in a dark place at the temperature of 2°C~8°C, and the shelf life is 7 days.

**【Test method】** Refer to the user manual of the applicable instrument.

**【Product performance indices】**

1. Accuracy: total particles: 900 particles/μL~1,300 particles/μL; relative bias shall be within the range of ±10.0%.

2. Homogeneity: CV≤10.0%

**【Cautions】**

1. This product is used for in vitro diagnostics only.

2. Shake the product for 1~3min before use, and do not shake it violently. If it is not mixed well, the number of tested particles may be smaller, thus affecting the calibration result.

3. Freeze is prohibited.

4. To avoid or reduce relevant infection risks, please wear disposable gloves and eyes/face protection articles when using the product.

**【References】**

1. Ferris JA urine microscopy comparison & standard: Lab Med14: 659-662,1983.

2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

**【Instruction Approved and Modified Date】** 05/2023

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Consult instructions for use

Temperature limit

Manufacturer

In vitro diagnostic medical device

Batch code

Use-by date

Keep away from sunlight

Authorised Representative

Catalogue Number

Comply with In Vitro Diagnostic Devices Directive (98/79/EC)

This Way Up

**DIRUI INDUSTRIAL CO., LTD.**

95 Yunhe Street, New & High Tech.

Development Zone Changchun, Jilin 130012 P.R.China

Tel: +86 431 85100409 Fax: +86 431 85172581

E-mail: dirui@dirui.com.cn <http://www.dirui.com.cn>

EC REP

Emergo Europe B.V. Westervoortsedijk 60  
6827 AT Arnhem The Netherlands

**DIRUI**<sup>®</sup>

## FE Control II Negative Control

### Instruction

**[Product name]** FE Control II Negative Control

**[Package specification]** Negative control: 1×20 mL, 1×22 mL, 1×25 mL, 1×60 mL, 1×125 mL.

**[Intended use]** This product is used for QC of test process of Urine Sediment Analyzer, FUS series urinalysis hybrid and MUS series urinalysis system.

**[Test principle]** Flow cytometry imaging

**[Main components]** Negative control: disodium hydrogen phosphate: 6.7 mmol/L.

**Note:** control range varies with the batch (see the label for fixed values of each batch).

**[Storage conditions and shelf life]**

1. This product shall be sealed and kept in a dark place at the temperature of 2°C~8°C. For shelf life refer to the label.
2. After being unsealed for the first time, it shall be sealed and kept in a dark place at the temperature of 2°C~8°C, and the open-vial shelf life is 30 days.

**[Test method]** Refer to the user manual of the applicable instrument.

**[Product performance indices]**

1. Accuracy: total particles number of negative control shall be  $\leq 20/\mu\text{L}$ .
2. Homogeneity: negative control not being tested

**[Cautions]** 1. This product is used for in vitro diagnostics only.

2. The product is prohibited from being frozen.

3. To avoid or reduce relevant infection risks, please wear disposable gloves and eyes/face protection articles when using the product.

**[References]**

1. Ferris JA urine microscopy comparison & standard: Lab Med14: 659-662,1983.

2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

**[Instruction Approved and Modified Date]** 05/2023

**[Symbol Explanation]**

 Consult instructions for use	 Temperature limit	 Manufacturer
 In vitro diagnostic medical device	 Batch code	 Use-by date
 Keep away from sunlight	 Authorised Representative	 Catalogue Number
 Comply with In Vitro Diagnostic Devices Directive (98/79/EC)		 This Way Up

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## FE Control II Positive Control

### Instruction

**[Product name]** FE Control II Positive Control

**[Package specification]** Level 1 of positive control: 1×20 mL, 1×22 mL, 1×25 mL, 1×60 mL, 1×125 mL;  
Level 2 of positive control: 1×20 mL, 1×22 mL, 1×25 mL, 1×60 mL, 1×125 mL;  
Level 3 of positive control: 1×20 mL, 1×22 mL, 1×25 mL, 1×60 mL, 1×125 mL;

**[Intended use]** This product is used for QC of test process of Urine Sediment Analyzer, FUS series urinalysis hybrid and MUS series urinalysis system.

**[Test principle]** Flow cytometry imaging

**[Main components]** Level 1 of positive control: QC blood: 50 particles/μL~300 particles/μL;  
Level 2 of positive control: QC blood: 300 particles/μL~700 particles/μL;  
Level 3 of positive control: QC blood: 900 particles/μL~1,300 particles/μL;

**Note:** control range varies with the batch (see the label for fixed values of each batch).

**[Storage conditions and shelf life]**

1. This product shall be sealed and kept in a dark place at the temperature of 2°C~8°C. For shelf life refer to the label.
2. After being unsealed for the first time, it shall be sealed and kept in a dark place at the temperature of 2°C~8°C, and the open-vial shelf life is 30 days.

**[Test method]** Refer to the user manual of the applicable instrument.

**[Product performance indices]**

1. Accuracy: Level 1 of positive control: relative bias shall be within the range of ±25.0%;  
Level 2 of positive control: relative bias shall be within the range of ±15.0%;  
Level 3 of positive control: relative bias shall be within the range of ±15.0%;
2. Homogeneity: Level 1 of positive control: CV≤15.0%;  
Level 2 of positive control: CV≤10.0%;  
Level 3 of positive control: CV≤10.0%.

**[Cautions]** 1. This product is used for in vitro diagnostics only.

2. Shake the product for 1~3min before use, and do not shake it violently. If it is not mixed well, the number of tested particles may be smaller, thus affecting the QC result.
3. Freeze is prohibited.
4. To avoid or reduce relevant infection risks, please wear disposable gloves and eyes/face protection articles when using the product.

**[References]**

1. Ferris JA urine microscopy comparison & standard: Lab Med14: 659-662,1983.
2. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition.

**[Instruction Approved and Modified Date]** 05/2023

**[Symbol Explanation]**

 Consult instructions for use	 Temperature limit	 Manufacturer	
 In vitro diagnostic medical device	 Batch code	 Use-by date	 Biohazard
 Keep away from sunlight	 Authorised Representative	 Catalogue Number	
 Comply with In Vitro Diagnostic Devices Directive (98/79/EC)		 This Way Up	

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- Ferris JA urine microscopy comparison & standard: Lab Med14:659-662,1983.
- Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition

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**[Package specification]** 1×20 mL, 1×22 mL, 1×25 mL, 1×60 mL, 1×125 mL.

**[Intended use]**

This product is used to determine position of focal plane of microscopic imaging system for Urine Sediment Analyzer, FUS Series Urinalysis Hybrid and MUS Series Urinalysis System.

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**[Main component]** QC blood: 2,500 particles/μL~3,500 particles/μL;

**Note:** fixed values of different batches of formed elements analysis focus are different (see the label for details).

**[Storage conditions and shelf life]**

1.This product shall be sealed and kept in a dark place at the temperature of 2°C~8°C. For shelf life refer to the label.

2.After being unsealed for the first time, it shall be sealed and kept in a dark place at the temperature of 2°C~8°C, and the open-vial shelf life is 30 days.

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**[Product performance indices]**

1. Accuracy: total particles: 2,500 particles/μL~3,500 particles/μL;

2. Homogeneity: CV≤15.0%

**[Cautions]**

1. This product is used for in vitro diagnostics only;

2. Shake the product for 1~3min before use, and do not shake it violently. If it is not being mixed well, the tested number of particles may be smaller, which will affect the focus effect.

3. Freeze is prohibited.

4.To avoid or reduce relevant infection risks, please wear disposable gloves and eyes/face protection articles when using the product.

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**[Symbol Explanation]**

Consult instructions for use	Temperature limit	Manufacturer
In vitro diagnostic medical device	Batch code	Use-by date
Keep away from sunlight	Authorised Representative	Catalogue Number
Comply with In Vitro Diagnostic Devices Directive (98/79/EC)		This Way Up

### DIRUI INDUSTRIAL CO.,LTD.

95 Yunhe Street, New & High Tech.  
Development Zone Changchun, Jilin 130012 P.R.China  
Tel: +86 431 85100409 Fax: +86 431 85172581  
E-mail: dirui@dirui.com.cn <http://www.dirui.com.cn>

Emergo Europe B.V. Westervoortsedijk 60  
6827 AT Arnhem The Netherlands

# Sheath

## Instructions

**[Product Name]** Sheath

**[Basic UDI-DI]** 69257945W010106020302R0FD

**[Package Specification]**

5L/barrel, 10L/ barrel, 15L/ barrel, 20L/ barrel, 28L/ barrel, 30L/ barrel

**[Intended Use]**

It is used for wrapping the sample to form a sheath flow to count formed elements in urine.

**[Test Principle]**

Flat flow cell imaging principle

**[Main Components]**

23mmol/L Tris buffer, 0.6% sodium chloride

**[Storage Conditions]**

Stored at 2-30°C in a dry place, protected from sunlight and in seal.

**[Shelf Life]**

1. The shelf life is 18 months if stored at 2-30°C and in seal.

2. The shelf life is 60 days after being opened if stored at 2-30°C.

**[Date of Manufacture & Expiry Date]** See the label

**[Applicable Instrument]**

DIRUI FUS Series Urinalysis Hybrid (Except FUS-2000), MUS Series Urinalysis System, Urine Sediment Analyzer (Model:FUS-360) .

**[Test Method]**

Use and install the reagent according to user manual of applicable instrument.

**[Sample Requirements]**

The sample type is fresh urine sample without centrifugation, placed at room temperature, and the urine should be tested within 2 hours.

**[Sample Processing]**

The sample does not require special treatment and can be used directly for testing.

**[Product Performance Indices]**

1. Appearance: colorless transparent liquid without sediments, particles or floccules

2. Net weight: no lower than the nominal value

3. pH value: at (25±1)°C, pH is 7.55±0.20

4. Osmotic concentration: (275±10)mmol/L, i.e. (275±10)mOsm/kg

5. Between-batch difference:

Test item	Between-batch difference
pH value	≤0.40
Osmotic concentration	≤20mmol/L (20mOsm/kg)

**[Warnings and Precautions]**

1. The product is only for in vitro diagnosis and professional use.

2. Avoid contact with eyes and skin;

3. Rinse eyes with the flow of warm water for a few minutes if the product splashes into eyes, ensure that the eyes are open. Seek medical care immediately if pain or swelling happens.

Once in contact with the skin, rinse with soap water and then rinse with water.

4. The test waste liquid should be treated in accordance with the medical waste management regulations.

5. If any adverse events occur during use, please report to the competent authority of the Member State where the manufacturer is located.

**Safety precautions:**

Hazard pictograms: No hazard pictogram is used

Signal word(S): No signal word is used.

**Preventive measure:**

Not applicable

**Emergency measure:**

Not applicable

For information about the safety technical manual, contact your local supplier or customer service representative.

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**[Symbol Explanation]**

 Please read Instructions	 Storage temperature	 Manufacturer
 Catalogue number	 Lot No.	 Expiry date
 Authorized representative	 In vitro diagnostic reagent	 County of manufacture
 European in Vitro Diagnostic Medical Device Regulation(EU) 2017/746(IVDR)		

**[Instructions Approved & Modified Date]** 2023-07

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B.V. Westervortsedijk 60