Isotonac-3, Isotonac-4, Hemolynac-3, Hemolynac-3N, Hemolynac-5, Hemolynac-310, Hemolynac-510 Cleanac-3, Cleanac, Cleanac-710

ITALIANO

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Nome: Isotonac-3, Isotonac-4, Hemolynac-3, Hemolynac-3N, Hemolynac-5, Hemolynac-310, Hemolynac-510, Cleanac-3, Cleanac, Cleanac-710

DESTINAZIONE D'USO

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BIBLIOGRAFIA

1. Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition, CLSI Document H26-A2, Vol. 30, No. 14, 2010.
2. Health Care Financing Administration. Clinical Laboratory Improvement Amendments of 1998; Final Rule. Federal Register, Vol. 57, [42 CFR § 493.1217, 493.1218], 1992.

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NIHON KOHDEN FIRENZE S.r.l

Via Torta, 72/74 - 50019 Sesto Fiorentino-Firenze, Italy
Tel +39 055 3045 1 - Fax +39 055 308548

Isotonac·3, Isotonac·4, Hemolynac·3, Hemolynac·3N, Hemolynac·5, Hemolynac·310, Hemolynac·510 Cleanac·3, Cleanac, Cleanac·710

ENGLISH

PRODUCT IDENTIFICATION

name: Isotonac·3, Isotonac·4, Hemolynac·3, Hemolynac·3N, Hemolynac·5, Hemolynac·310, Hemolynac·510, Cleanac·3, Cleanac, Cleanac·710

INTENDED USE

For in vitro diagnostic use only. Reagents to be used on NIHON KOHDEN analyzers in combination as indicated in Table 1.

Isotonac·4, Hemolynac·310, Hemolynac·510 and Cleanac·710 are only intended for human blood. In regard to the intended use, please refer to Table 2.

STORAGE AND SHELF LIFE AFTER FIRST OPENING

Store between 1 and 30 °C. Check reagent integrity and expiry date, before installing the device. Once opened, the reagent must be used within the expiry date indicated in Table 4. When not in use, use the original cap to close the container except Hemolynac·310, Hemolynac·510 and Cleanac·710. Do not force.

COMPOSITION

Active ingredients

Ref. Table 3

WARNING AND PRECAUTIONS

Labelling information

Ref. Table 5.

Disposal information

Unused product and emptied containers should be disposed in accordance with local environmental regulations. Instrument wastes, as potential biohazardous, must be handled with care and disposed in accordance with local environmental regulations.

SPECIMEN

The devices have been designed to be used with whole blood samples. Refer to Operator's Manual for sample collection and storage conditions.

PROCEDURE

Refer to Operator's Manual for reagents installation and priming, and for instrument operation. When installing a new batch of diluent or lysing reagent, check the background according to Operator's Manual. Refer to Operator's Manual for working temperature.

Use each reagent in its own packaging (do not transfer). Follow GLP recommendations. If the instrument connections were not compatible, work according to Good Laboratory Practice (GLP), in order to avoid any kind of contamination, either from external environment to the device, or from the device to the external environment.

METHODOLOGY AND CALCULATION OF ANALYTICAL RESULTS

Refer to Operator's Manual for information on methodology and calculation of analytical results.

INTERNAL QUALITY CONTROL

Background verification

Before running testing samples check background count according to Operator's Manual, and verify that all the values for each parameter are within Operator's Manual limits. If not, follow Operator's Manual prescriptions. In case of problems attributed to one of the reagents, install a new unit of reagent and verify the background again. Do not run testing samples until the background is within the specified limits.

Calibration and calibration verification procedures

Instrument calibration should be verified in the following conditions:

1. after first installation
 2. at least every six months (or according to manufacturer's requirements for calibration frequency, if more stringent)
 3. after a major repair or a major change (replacement of critical parts that may influence test performance, software upgrade, complete change of reagents)
 4. when quality control data reflect a significant change of baseline accuracy (drift), or an out-of-control condition
 5. in case of replacement of the lysing reagent Hemolynac·3 with the lysing reagent Hemolynac·3N and vice versa
- Calibration procedures should be preceded by a preventive maintenance (verification of baseline accuracy and precision: ref. CLSI document H26-A2). Calibration and calibration verification must be performed using appropriate calibration material (fresh whole blood samples which have been analyzed by reference methods, or appropriate calibrator), following Operator's Manual prescriptions. When calibration is required, verify instrument calibration immediately after calibrating, according to Operator's Manual instructions (refer also to CLSI document H26-A2).

Daily Quality Control Procedures

In order to assure the accuracy and reliability of test results, it is recommended to perform quality control procedures, according to Operator's Manual instructions, using appropriate control material (retained patient specimens or stabilized control blood), immediately before, during, or immediately after each analysis run (intended as an analysis sequence). It is recommended to run two different levels of stabilized control blood (refer to CLSI document H26-A2).

BIBLIOGRAPHY

1. Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition, CLSI Document H26-A2, Vol. 30, No. 14, 2010.
2. Health Care Financing Administration. Clinical Laboratory Improvement Amendments of 1998; Final Rule. Federal Register, Vol. 57, [42 CFR § 493.1217, 493.1218], 1992.

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Name: CLEANAC-810

INTENDED USE

For in vitro diagnostic use only. Reagents to be used on a Nihon Kohden MEK-9100 automated hematology analyzer.

CLEANAC-810 is only intended for human blood.

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Store between 2 and 8 °C. Check reagent integrity and expiry date before installing the device. Once opened, the reagent must be used within the day.

COMPOSITIONActive ingredients

Sodium hypochlorite (% active chlorine) 6-7 % w/v

WARNING AND PRECAUTIONSLabelling information

No hazard label required according to 1272/2008 regulation (CLP) art. 1,5,d). Active ingredients: Sodium hypochlorite (%Cl active) 6-7% w/v. Do NOT mix with acid. Shipment labelling: UN 1791 PG III – Hypochlorite solution in limited quantity.

Disposal information

Unused product and emptied containers should be disposed of in accordance with local environmental regulations. Instrument waste, as a potential biohazard, must be handled with care and disposed of in accordance with local environmental regulations.

PROCEDURE

Refer to the Operator's Manual for reagent installation and priming, and for instrument operation. Refer to the Operator's Manual for working temperature.

Use the reagent in its own packaging (do not transfer).

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Some specifications of this product may differ depending on the destination country or region. Therefore descriptions in the Japanese manual and English and other language manuals may also differ.

The CE mark is a protected conformity mark of the European Community. Products marked with this symbol comply with the requirements of the In vitro Diagnostic Medical Device Directive 98/79/EC.

General Handling Precautions

WARNING

Dispose of the hematology calibrator by following your local laws for disposing of infectious medical waste. Otherwise, it may affect the environment.

WARNING

Do not swallow the hematology calibrator. It may be infectious.

WARNING

When handling the hematology calibrator, wear rubber gloves to protect yourself from infection.

NOTE

- This hematology calibrator is potential biohazardous material. MEK-CAL hematology calibrator is intended solely for in vitro diagnostic use by trained personnel. This product contains human-sourced and/or potentially infectious components. Components from human donors used in preparation of this product were tested by FDA approved methods for the presence of the antibodies to human immunodeficiency virus (HIV-1/HIV-2) and Hepatitis C virus (HCV) as well as for Hepatitis B virus surface antigen and found to be negative.
- Do not use the hematology calibrator after the expiration date. The measurement data may be inaccurate.
- Store the hematology calibrator in a medical refrigerator.
- Store the hematology calibrator at temperatures between 2 and 8°C (36 and 46°F).
- Use the hematology calibrator at temperature from 18 to 30°C (64 to 86°F). Otherwise, the measurement data may be inaccurate.
- Mix the hematology calibrator thoroughly before use.
- Platelet cannot be measured by precipitation or centrifuge.

General

The MEK-CAL hematology calibrator is for calibration of Nihon Kohden hematology analyzers. For calibration, refer to the provided assay table for specific instrument models and the Operator's Manual of the hematology analyzer.

NOTE

MEK-CAL hematology calibrator is designed for Nihon Kohden hematology analyzers. Do not use with any other instrument.

 **NIHON KOHDEN**

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku, Tokyo 161-8560, Japan
Phone +81 (3) 5996-8036 Fax +81 (3) 5996-8100

 European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003 827-0 Fax +49 6003 827-599

Reagent Production
NIHON KOHDEN FIRENZE S.r.l.
Via Torna 72/74, 50019 Sesto Fiorentino (FI), Italy
Phone +39 055 3045 1 Fax +39 055 308548

SHANGHAI KOHDEN MEDICAL ELECTRONIC
INSTRUMENT CORPORATION
7th Floor, Dawning Centre Tower A No.500, Hongbaoshi Road
Changning District, Shanghai 201103, China
Phone +86 (21) 6270-0909 Fax +86 (21) 6270-9700

MEK-CAL hematology calibrator contains human erythrocytes, simulated leukocytes, and mammalian platelets in a plasma-like fluid. MEK-CAL hematology calibrator is for in vitro diagnostic use.


To achieve optimum performance and safety, thoroughly read this instruction and the hematology analyzer operator's manual before using this hematology calibrator.

Instruction for Use

NOTE

- Dark red supernatant fluid, discoloration of the product or unacceptable results may indicate deterioration. Do not use the hematology calibrator if deterioration is suspected.
- Confirm that the lot number of the product matches the lot number on the assay sheet.

- Take the hematology calibrator out of the medical refrigerator and allow it to warm at a temperature of 18 to 30°C for 15 minutes before mixing.
- Mix the hematology calibrator by holding it horizontally between the palms of the hands. Roll it back and forth for 20 to 30 seconds and occasionally turn it upside down.

NOTE

- Do not mix the hematology calibrator roughly.
- Do not mix the hematology analyzer using a mixer.
- Mix the hematology calibrator until the red cells are completely suspended. Mix thoroughly but do not shake it.
- Immediately before measuring, gently turn the hematology calibrator upside down 8 to 10 times.

NOTE

Incomplete mixing of the hematology calibrator affects the measurement result.

- Calibrate the hematology analyzer using the hematology calibrator as instructed in the hematology analyzer Operator's Manual.
- Return the hematology calibrator to the medical refrigerator within 30 minutes of use.

Expected Results

For calibration values, refer to the assay sheet.

Storage

Store at 2 to 8°C.

NOTE

Do not freeze the hematology calibrator.

Expiration Date

The expiration date is shown on the container label.

Expiration Date After Opening the Package

Opened hematology calibrator should be used within 7 days.

Copyright Notice

The entire contents of this manual are copyrighted by Nihon Kohden. All rights are reserved.

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