# AFIAS-10

## Quick Operation Guide


# AFIAS-10 Start-up

- 1)  Turn on the switch 'Power On/OFF'.  
(The power switch is located at the bottom left of AFIAS-10 back)
- 2)  Press the power button on the front of the equipment to turn it on.
- 3) When turned on, AFIAS-10 briefly displays the Boditech Med logo and enters the login window.
  - ① For the first log in, please insert the password **0000**.  
The instrument will request to change the password after the first login.\*

- 4) When the login is completed, after going through the device initialisation process the main screen is displayed.

- 5) On the main screen, check the temperature of the equipment at the bottom.


If it is confirmed that the equipment temperature is green (**35 ± 2 °C**), the equipment is ready for testing.

- ▶ To conduct the AFIAS test, prepare the following items:
- ▶ AFIAS ID Chip.  
: If the LOT to be tested has not been registered before, proceed with registering the ID Chip.  
Registration is needed only once per LOT.
- ▶ AFIAS cartridge.
- ▶ Patient sample or Boditech control.

\* refer to separate manual for more detailed instructions.


# Main Screen Features


## ID Chip Registration

① Insert the ID Chip you want to register into the ID Chip port and touch the 'ID Chip' button on the main screen. When the ID Chip List window appears, the ID Chip information is automatically saved.

\* Insert the ID Chip into the ID Chip port (side with the item name visible).


② When registering more than one item, insert the second ID chip while the ID chip List window is open and press the Load button to register any additional item.

\* AFIAS-10 can register 100 ID chip information.

\* When changing Lot, the Lot specific ID Chip needs to be registered.

# Cartridge Structure


## Pre-inspection Precautions

- ▶ Verify that the lot numbers printed on the test cartridge and ID Chip match.
- ▶ Verify the expiry date.
- ▶ Use cartridges immediately after opening packaging.
- ▶ If the cartridge was refrigerated, leave it at room temperature for at least 15 minutes before proceeding with the test.
- ▶ Use a separate tip for dispensing samples onto the cartridge well (inside the instrument only insert tips from the box).
- ▶ Do not use tips if they are curved or damaged.


## Pre-inspection Preparations

[Tip Fill]


### 1) Tip Fill & Cartridge placing on Tip Station

- ① Insert the tip in the tip stations on the side of the holder.

- ② Touch the tip station part on the main screen and correctly insert the number of tips installed.

## [Cartridge Inserting]


\* When a cartridge is installed, the slot changes to activated.

① Insert cartridges.

② Insert the tips into the tip holder of the inserted cartridge.



③ Click the 'Load' button on the holder where the cartridges are inserted. (AFAS-10 will read the cartridge barcode and display the cartridge item name).

After checking the barcode, the corresponding item will be shown.

## 2) Specimen and Rack Preparation

### [Sample preparation]

\* Start from the right, where the rack barcode is located.

① Blood tube rack.

(A) Prepare the blood tube rack and specimen.

\* Blood tube available size: 13mm x 75~100 mm

\* Minimum sample volume in the tube: 1mL (Dead volume) + a

(B) Open the lid of the blood tube and place the tubes on the rack.

\* In the case of a whole blood sample, before taking off the lid slowly invert twice up and down then proceed. (Be careful not to create bubbles in the tube when inverting)

\* When attaching a patient barcode, position it vertically in a similar way as the container's own label.

\* Sample barcode recognition area: 1cm lower than the smaller tube edge to 1cm higher than the smaller tube bottom.  
\* smaller tube size: 75mm.

## ② Eppendorf tube rack

(A) Prepare the Eppendorf tube rack and the Eppendorf tube containing the specimen.

Minimum sample volume in the tube: 40 uL (Dead volume) + a (refer to the product insert for the volume)


(B) After fully opening the lid of the Eppendorf tube, place it in the rack following the instruction below.

\* Start from the right, where the rack barcode is located.


## [Installing the Specimen Tube Rack]

① Press the rack 'unload' button on the bottom right corner, in order for the rack holder to be ejected.

② Insert the tube rack.

\*When loading the rack, it automatically detects the presence of the tubes and displays them.

# Test Method 1) Automatic order designation (bidirectional LIS)

[Place the cartridges first]

① After placing the cartridge in the holder, select "unloaded". **Unloaded**

② Place the tip into the tip position.

③ Open the lid of the sample tube and place in the sample tube rack.

→ The sample rack enters the equipment, automatically recognises the tube barcodes, select the order, and matches them directly to the cartridges.

[place the sample tube first]

① Open the cap of the sample tube and install the sample tube rack.


→ The sample rack enters the equipment, automatically recognises the tube barcodes, select the order, and the required quantity of cartridge per item appears in the pop-up window.

② After placing the requested cartridges, press "unloaded". **Unloaded**

→ The cartridges and the sample tubes will match automatically.

③ Insert the tip into the tip position.

[In case of changing the sample type is needed]


(A) Touch the sample rack area to open the sample window.  
(B) Select sample type.


(C) Select the sample tube you want to change the sample type.


D) Touch the 'x' button to save the entered information.


④ Touch the 'Start' button.

⑤ Select 'OK'.

# Test Method 2) Manual Ordering


If not all cartridges are of the same item, it is required to manually match sample tubes and cartridges.

- ① When the cartridges have been loaded, select. Unloaded
- ② Insert tips in the tip positions.
- ③ Load the sample tube rack.
- ④ Touch the sample rack area on the screen to designate sample tubes to the cartridges.


Sample Window

Sample type	BloodTube	WB	S/P							X					
Sample tubes	1	2	3	4	5	6	7	8	9	10					
Patient ID (name)															
Patient ID (barcode)															
item list	<table style="width: 100%; border-collapse: collapse;"> <tr> <td style="border: 1px solid gray; padding: 2px;">CK-MB (2)</td> <td style="border: 1px solid gray; padding: 2px;">PCT Plus (1)</td> <td style="border: 1px solid gray; padding: 2px;">Tn-I Plus (2)</td> <td style="border: 1px solid gray; padding: 2px;">D-Dimer (2)</td> <td style="border: 1px solid gray; padding: 2px;">NT-proBNP (2)</td> </tr> </table>										CK-MB (2)	PCT Plus (1)	Tn-I Plus (2)	D-Dimer (2)	NT-proBNP (2)
CK-MB (2)	PCT Plus (1)	Tn-I Plus (2)	D-Dimer (2)	NT-proBNP (2)											
Add/Remove option	<input type="radio"/> Add <input type="radio"/> Remove														

\* The item list shows the item name of the inserted cartridge (with available quantity).


- ⑤ Select a sample type. WB S/P
- ⑥ Touch the item from the list (select add).
- ⑦ Select a sample (tube) which corresponds to the selected parameter that needs to be tested.


\* It is possible to match up to 10 items with one sample.


e.g.1) Match S/P (serum/plasma) and Tn-I to samples 1 and 2.

e.g.2) Match whole blood and CRP for samples 1 to 3 and whole blood for samples 4 to 6 and Tn-I Plus.


e.g.3) Matching S/P and NT-proBNP to CK-MB in samples 3 and 4.

\* Remove designated items.


- ① Select the corresponding item from the item list part.
- ② Change to 'Remove'.
- ③ Touch and remove the sample (if you add an item again later, change to 'Add' and proceed).

- ⑧ If you wish to add a patient ID, it is possible to type it using the keyboard on the screen. When finished, touch the X button to save the information.

- ⑨ Select the button 'Start'.

- ⑩ Press 'OK'.


# Test Method 3) Emergency Test


Emergency test allows to manually insert the sample, this can be done in case of small sample volume or emergency examination.

- ① Insert the cartridge of the test needed for an emergency check.
- ② Put the required sample amount into the cartridge sample well.



- ③ Insert the tip in the tip position.

- ④ Click the 'Load' button on the corresponding bay.

- ⑤ Click the keyboard button on the bottom of the cartridge.

- ⑥ Enter Patient ID (optional).

- ⑦ Select Emergency Test.

- ⑧ Select General Tip.

- ⑨ Select the correct sample type.


- ⑩ Press Set.


- ⑪ The emergency test will then acquire the priority and will be processed immediately.

# Delete the Test

## 1) How to cancel while a test is processing


- ① Touch the icon 'Stop' → Select OK.


- ② AFIAS-10 will cancel all tests, with the exception of those that are already at the reaction time step.

## 2) Canceling an individual test


- ① Touch the center of the cartridge slot for at least 3 seconds → Select OK.

- ② Only the selected cartridge will be cancelled.

## 3) Canceling all bay tests

- ① Touch for at least 3 seconds the button under the Bay that you want to cancel → press OK.

- ② All tests will be cancelled.

# Maintenance

## 1) Clean Tip removal area

\* Slide on the left side of the instrument.

## 2) Cleaning and emptying the tip waste box

## 3) System check

- ① Remove any remains from the tip removal slide on the left side of the equipment using an alcohol solution.

- ① Remove the waste box attached to the side of the equipment, discard the tip, and clean the inside of the tip box.
- ② Attach the waste box in place using the magnet.

- ① Prepare the System check cartridge, System check ID chip  
When to use it: When AFIAS-10 is first installed, when AFIAS-10 is moved to another location, when an error occurs, when the test result is considered abnormal or different from the reference value.


- ② How to use.

(A) Register the System Check ID chip.

(B) Insert the System check cartridge into the cartridge slot and press the 'Unloaded' button.


C) Press the 'Start' button to proceed with the inspection and check the result.

# Register and Test Quality Control (QC)

- ① Press the **Analyzer** button at the top of the main screen to switch to **QC** mode.
- ② Press the button **Control Material** to open the control material window.
- ③ If the control substance is not registered or the control Lot has changed, press **Add** to add controls (repeat for each level) and press the 'Register' button.


Item	Lot	Exp.	Sub Item	Mean	Low	High
HbA1c	AACQK16_LV.2	2021.11.30		9.9	8.42	11.39
HbA1c	AACQK16_LV.1	2021.11.30		5.3	4.51	6.1

- ④ Insert the cartridge to be used and press the **Unloaded** button to recognise the cartridge.
- ⑤ If 2 or more cartridges have been inserted, the values Lv1 and Lv2 will be assigned automatically. If manual application of the control material is required, press,  to set the control material.
- ① Control Lot (Lv1 / Lv2) ② Tip type ③ Sample type.


- ⑥ Add the control material to be inspected into the sample well.
- ⑦ Add the tip into the tip hole.
- ⑧ Press the 'Start' button to proceed with the QC.
- ⑨ QC results can be found both in History and in ID-Chip information.

## - History

✓	No	Patient ID	Item	Sample Type	Result	Date	Lot	Print
	4	21225 506	HbA1c	WB	4.9	2021.02.28	AAQGB18F	Search
		21225 506	AAQGB18F	OC	%	11.15.28		
✓	5	21225 505	HbA1c	WB	4.2	2021.02.28	AAQGB18F	LIS
		21225 504	AAQGB18F	OC	%	11.13.42		
✓	4	21225 504	HbA1c	WB	Invalid	2021.02.28	AAQGB18F	
		21225 504	AAQGB18F	OC	%	10.29.52		
✓	3	21225 003	HbA1c	WB	Invalid	2021.02.28	AAQGB18F	
		21225 002	AAQGB18F	OC	%	10.28.16		
✓	2	21225 002	HbA1c	WB	6.2	2021.02.28	AAQGB18F	
		21225 002	AAQGB18F	OC	%	09.50.03		
✓	1	21225 001	HbA1c	WB	6.3	2021.02.28	AAQGB18F	
		21225 001	AAQGB18F	OC	%	09.49.18		

## - ID Chip

Item	Lot	Exp	Rev	QC
HbA1c	AAQGB18F	12.28.2022	0	Pass
Tn-I Plus	TNRDB10G	12.28.2022	0	Required
hsCRP	CAPLR01A	12.28.2022	0	Required
PSA	PSRDA39F	12.02.2022	0	Required
Tn-I Plus	TNRCB05G	11.29.2022	0	Required

Item	Lot	L1	L2	Status	Date
HbA1c	AAQGB18F	AAQGB18F_LV1 (4.31 - 6.17) 5.0	AAQGB18F_LV2 (8.76 - 11.85) 10.0	Pass	05.26.2020
HbA1c	AAQGB18F	AAQGB18F_LV1 (4.31 - 6.17) 10.0	AAQGB18F_LV2 (8.76 - 11.85) 5.0	Fail	05.26.2020

# AFIAS-10 Shut Down

## 1) Sleep mode

- ① Press briefly the power button  to enter sleep mode.

## 2) Reboot

- ① Press and hold the button  for at least 3 seconds and touch 'Reboot' in the active window.
- ② The equipment reboots again.

## 3) Shut down

- ① -1  Press and hold the power button for at least 4 seconds to enter the 'Power Off' window.
- ① -2 Press and hold the button  button or at least 3 seconds and touch 'Power off' in the active window.
- ② When you touch the 'OK' button, the instrument will proceed with the shutting down procedure.
- ③ Turn off the main power on/off switch on the back of the instrument.


# Frequently Asked Questions

- The below issues can appear in the pop-up window or on the screen.

Issue	Cause	Solution
<b>No ID Chip</b> 	If the inserted cartridge and the lot on the ID chip do not match.	Please register the ID Chip and proceed with the test again.
	If the inserted cartridge and the lot on the ID chip do not match.	Insert the cartridges of Lot that match the ID Chip and proceed with the test.
<b>Unknown</b> 	Barcode printed on the cartridge: If the barcode is obscured by other material or the printing is damaged and the equipment does not recognise the cartridge barcode.	<ol style="list-style-type: none"> <li>1. Remove any other material that obscures the barcode and reinsert the cartridge.</li> <li>2. Recognition failure: Press the center of the cartridge slot picture on the screen to select manually the correct ID chip information (ID chip and cartridge LOT must match)</li> <li>3. If the same error code continues to occur, test using a new cartridge.</li> </ol>
	Barcode module failure.	Contact Boditech Med or the local representative for technical support.
<b>Tip required</b> 	If the test was carried out without a tip.	Check whether the tip has been inserted, add tip(s) to the cartridge or Extra Tip station, and touch the center of the slot in question to re-examine.
<b>Sample Empty</b> 	Sample detection error.	Check the sample and proceed with the examination again.
	If sample amount is small.	
<b>Sample Error</b> 	If the sample tube lid is not removed when inserted in the sample rack.	Check the tube lid removal status and proceed with the inspection again.
	Sample Error. (Bubble Detection or Blood Clot)	<p>If there is too much foam in the sample tube, remove the bubbles in the tube and re-examine the sample.</p> <p>If there are clots, take the sample again and re-examine.</p>
<b>WB (B)</b> <b>S/P (B)</b>	Sampling attempted 2-3 times before processing (possible problem in the sample, such as bubble).	Please consider examining again and report the final result.
<b>Invalid TF</b> <b>/ Invalid 25~ 34</b>	In case of abnormal graph, the instrument shows one of these errors.	Please test with the new cartridge. If the same error occurs again, call Boditech Med or the local representative for technical support.
<b>Extra Tip 1~X</b>	If the Extra Tip quantity is insufficient to process specific items. (T3, Testosterone, Total Anti-infliximab ...)	After replenishing the Extra Tip station, proceed with the inspection.
<b>Abnormal Temperature,</b> <b>Test banned</b>	If the test is started with the temperature inside the device not reaching the appropriate temperature.	When the temperature inside the instrument reaches the appropriate temperature, the test will resume.

# Trouble shooting

Issue	Cause	Solution
Tip is not removed	After tip removing step, the device still recognises that there is a tip.	After the technician cleans the tip sensor contaminated area, proceed with the test again.
	The instrument attempted to remove the Tip, but it could not actually remove it.	Clean the contaminated area of the tip removal slide.
If the sample rack goes in but comes out immediately	Sample rack type not recognised.	Attach a new sample rack barcode and proceed with the inspection again.
If the sample rack does not enter automatically	If the sample rack itself is not recognised.	Contact Boditech Med or the local representative for technical support.
If the equipment does not recognize the sample tube barcode	The sample barcode location is outside the barcode recognition range.	If the barcode is outside the recognition area, reattach it to the correct position and proceed with the inspection again. -13 mm x 75 mm tube : attach it at least 1 cm above the bottom of the tube -13 mm x 100 mm Tube : 1cm above the bottom of the tube
	If the sample barcode is damaged and cannot be read from the equipment.	Re-print the barcode, attach it, and proceed with the testing it again.
	If the equipment fails to rotate the sample barcode.	1. If the sample barcodes are pasted in layers, remove the existing barcodes, attach a single barcodes, and proceed with the re-examination.  2. If the sample barcode gets detached, secure it and conduct again the inspection. If the problem rises when the sample barcode is not loosen or detached, contact Boditech Med or the local representative for technical support.
If the cartridge is not present but the instrument recognizes it as present	Possible cartridge sensor malfunctions due to a light source entering the equipment.	Please change the position of the equipment in a direction where no light source enters directly.
If the cartridge itself is not recognized (inactivated)	Cartridge sensor failure.	Contact Boditech Med or the local representative for technical support.
Initialization fail	If the equipment moving parts are not in the correct position.	Please proceed with the equipment reboot.

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**European Representative:** **Obelis s.a.**  
 Bd. Général Wahis, 53, 1030 Brussels, Belgium  
 (SRN: BE-AR-000000106)

**Product:** **AFIAS-10**  
 - Catalog No.: FPRR038  
 - Basic UDI-DI: 880613301031BM

**Intended use:** AFIAS-10 is an analyzer intended for use in conjunction with fluorescence immunoassay (FIA) kits for quantitative, semi-quantitative and qualitative measurements of various analytes.  
 For *in vitro* diagnostic use only.

**Classification:** Class A (Rule 5)  
 - According to the Annex VIII of the REGULATION (EU) 2017/746

**Conformity Assessment Route:** According to the Annex II and Annex III of the REGULATION (EU) 2017/746

We, Boditech Med Inc., herewith declare under our sole responsibility that the above-mentioned product is in conformity with the following European Union harmonisation legislation.

- REGULATION (EU) 2017/746 OF EUROPEAN PARLIAMENT AND OF THE COUNCIL OF 5 April 2017 on *in vitro* diagnostic medical devices and repealing Directive 98/79/EC and Commission Decision 2010/227/EU
- DIRECTIVE 2011/65/EU OF EUROPEAN PARLIAMENT AND OF THE COUNCIL OF 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment
- COMMISSION DELEGATED DIRECTIVE (EU) 2015/863 of 31 March 2015 amending Annex II to Directive 2011/65/EU of the European Parliament and of the Council as regards the list of restricted substances

All supporting documentation is retained under the premises of the manufacturer

**Standards applied:** EN ISO 13485:2016, EN 13612:2002, EN ISO 14971:2019,  
 EN ISO 15223-1:2021, EN ISO 18113-1:2011, EN ISO 18113-3:2011,  
 EN 61010-1:2010, EN 61010-2-101:2017, EN IEC 61326-1:2021,  
 EN IEC 61326-2-6:2021, EN 62304:2006, EN 62366-1:2015,  
 ISO/TR 20416:2020, IEC 62321-1:2013

**Place, Date of Issue:** Chuncheon, Republic of Korea, February 27, 2023

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6 November 2023

To whom it may concern,

We hereby authorize INTERLUX Ltd. With its principal place of business at **Aviečių str. 16, LT-08418, Vilnius, Lithuania** to act as the exclusive distributor for the below-stated products manufactured by Boditech Med Inc. with its principal place of business at 43 Geodudanji 1-gil, Dongnae-myeon, Chuncheon, Gangwon-do, Republic of Korea, Postal Code 24398 (the "Company").

INTERLUX Ltd. Is authorized to register, sell, market, provide full guarantees and technical service in Lithuania, Latvia, and Estonia market.

- **ichroma™ analyzers and reagents**
- **AFIAS analyzers and reagents**

This Letter of Authorization shall become effective as of January 1, 2024 and expires on December 31, 2024. The terms of this Letter of Authorization may be renewed for an additional year period, subject to the terms of any definitive agreement which shall be entered into by and between the Distributor and the Company. No change or amendment to this Letter of Authorization shall be valid unless agreed in writing by the Company.

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Vaidas Vilmantas (MB „Beikeris“, įm .k. 304539005)

# Certificate Of Training

**This is to certify that Mr. Rokas Arlauskas at Medita Co. Ltd. – Interlux Group has successfully completed the basic training requirement in accordance with the provisions of Boditech Med, Inc. to use and service ichroma™ and AFIAS platforms.**

**This entitles the recipient to perform limited services for the above products.**

**Certificate No. BT231005-01**

**Date of Issuance : October 05<sup>th</sup>, 2023**

**Authorized by Mr. Kun Heo**

---

**Boditech Med, INC.**

Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon,  
Chuncheon-si, Gangwon-do 24398, Korea  
T. +82 33 243 1400 / F. +82 33 243 9373

Wednesday, March 6, 2024

## STATEMENT

To whom it may concern,

We **Boditech Med Inc.**, a corporation duly organized and existing under the laws of the Republic of Korea, having its principal office at 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si Gangwon-do, 24398 Korea, hereby certify that:

- The product **SMFP-43 AFIAS Rota/Adeno** is a qualitative product (results are positive, negative, or indeterminate). For qualitative products we do not provide calibrator materials.

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## STATEMENT

To whom it may concern,

We **Boditech Med Inc.**, a corporation duly organized and existing under the laws of the Republic of Korea, having its principal office at 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si Gangwon-do, 24398 Korea, hereby certify that:

- The **AFIAS-10** analyzer incorporates a feature allowing users to fill in the User ID field during login through either manual input or barcode scanning. Patient ID field can also be populated by either manual input or barcode scanning.

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# AFIAS CRP

## INTENDED USE

**AFIAS CRP** is a fluorescence immunoassay (FIA) for the quantitative determination of C-Reactive Protein (CRP) in human whole blood/serum/plasma. It is useful as an aid in management and monitoring of autoimmune diseases and infectious processes, such as rheumatoid arthritis.

For *in vitro* diagnostic use only.

## INTRODUCTION

C-reactive protein (CRP) is a protein found in blood; the level of which rises in response to inflammation. CRP is first acute-phase protein to be described and is an exquisitely sensitive systemic marker of inflammation and tissue damage. The serum CRP level may rise from a normal level of < 5 mg/L to 500 mg/L during the body's general, nonspecific response to infections and other acute inflammatory conditions. Measurement of CRP concentration has been widely used as a clinical tool for monitoring the status of inflammation, effectiveness of treatment of various infections and autoimmune diseases such as rheumatoid arthritis. And high sensitive CRP (hsCRP) has been suggested that along with serum cholesterol, may be marker of diagnosis in cardiovascular diseases (CVD). And hsCRP is emerging as the strongest and most independent predictive risk factor for CVD.

## PRINCIPLE

The test uses a sandwich immunodetection method. The detector antibodies in detector bind to antigens in the sample, forming antigen-antibody complexes and migrate onto nitrocellulose matrix to be captured by the other immobilized-antibodies on a test strip.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by the instrument for AFIAS tests to show CRP concentration in the sample.

## COMPONENTS

**AFIAS CRP** consists of 'cartridges'.

- Each sealed aluminum pouch contains two cartridges.
- Each sealed cartridge packaged in an aluminum pouch has three components including a cartridge part, a detector part and detector diluent part.
- The cartridge part contains the membrane called a test strip which has mouse monoclonal anti-human CRP at the test line, human C-reactive protein at the antigen line, and streptavidin at the control line.
- The detector part has a granule containing mouse monoclonal anti-human CRP-fluorescence conjugator, mouse monoclonal anti-human CRP, Biotin-BSA fluorescence conjugator, bovine serum albumin (BSA) and sucrose as a stabilizer and sodium azide as a preservative in phosphate buffered saline (PBS).
- The detector diluent part contains antifoam and tween 20 as a detergent, sodium azide as a preservative in phosphate buffered saline (PBS).

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this 'Instructions for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).
- Do not reuse cartridges, pipette tips, and C-tips. It should be used for testing one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use the cartridge, if pouch is damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- If test components and/or sample are stored in refrigerator, then allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used cartridges and pipette tips, and C-tip should be handled carefully and discarded by an appropriate method in accordance with relevant local regulations.
- The cartridge contains sodium azide (NaN<sub>3</sub>), and they may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.
- Biotin interference was not observed in **AFIAS CRP** when biotin concentration in the sample was up to 3,500 ng/mL. In reference, if biotin is taken at a level of 300 mg/day, the biotin concentration in the body is about 1,200 ng/mL. However, if a patient has been taking biotin at dosage of more than 300mg a day, it is recommended to collect blood again 24 hours after discontinuation of biotin intake.
- AFIAS CRP** will provide accurate and reliable results subject to the below conditions.
  - AFIAS CRP** should be used only in conjunction with the instrument for AFIAS tests.
  - Have to use recommended anticoagulant.
 

Recommended anticoagulant
K <sub>2</sub> EDTA, K <sub>3</sub> EDTA, Sodium heparin, Lithium heparin, Sodium citrate
- C-tip should be used when the following conditions are met.**
  - C-tip provided with the kit is recommended to obtain correct test result.
  - Whole blood should be immediately tested after collection.
  - Do not perform a test with C-tip on General Mode. It might cause an erroneous result.
  - Excess whole blood around the C-tip should be wiped off.
  - In order to avoid cross-contamination, please do not re-use C-tip for multiple samples.
  - AFIAS cartridge should be inserted and positioned in the cartridge holder prior to the blood sample collection.
  - While collecting blood, be careful not to create air bubbles in the C-tip.

## LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is the most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician in conjunction with clinical symptoms and other relevant test results.

## STORAGE AND STABILITY

Component	Storage condition		
	Storage Temperature	Shelf life	Note
Cartridge	2 - 30 °C	20 months	Unopened
		1 month	Resealed

- Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal.

## MATERIALS SUPPLIED

REF	SMFP-2	
Components of <b>AFIAS CRP</b>		
Cartridge box :		
	- Cartridge	24
	- Pipette tip (zipper bag)	24
	- C-tip (10 µL)	24
	- Spare cartridge zipper bag	1
	- ID chip	1
	- Instructions for use	1

## MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from **AFIAS CRP**. Please contact our sales division for more information.

Instrument for AFIAS tests	
- <b>AFIAS-1</b>	REF FPRR019
- <b>AFIAS-3</b>	REF FPRR040
- <b>AFIAS-6</b>	REF FPRR020
- <b>AFIAS-10</b>	REF FPRR038
■ <b>Boditech CRP Control</b>	REF CFPO-100
■ <b>Boditech CRP Calibrator</b>	REF CFPO-112

## SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS CRP** is human whole blood/serum/plasma.

- It is recommended to test the sample within 24 hours after collection.
- The samples (serum, plasma) should be separated from the clot by centrifugation within 3 hours after the collection of whole blood.

- The samples (whole blood, serum, plasma) may be stored for a week at 2-8 °C prior to being tested. If testing will be delayed more than a week, samples should be frozen at below -20 °C.
- The samples stored frozen at below -20 °C for 3 months showed no performance difference.
- However, the whole blood sample should not be kept in a freezer in any case.
- As a repeated freeze-thaw cycle may affect the test result, do not refreeze previously frozen samples.
- Collection of whole blood sample using C-tip
  - Hold the C-tip horizontally and touch the surface of the blood with the tip of the C-tip.
  - Capillary action will automatically draw the blood sample to C-tip and stop.
  - Wipe off any excess blood around the tip.
  - Double-check if whole blood is filled accurately in the C-tip and AFIAS reader is ready for a test on the 'C-tip mode'.

## TEST SETUP

- Check the components of the **AFIAS CRP** as described below.: Cartridges, pipette tips, C-tips, an ID chip, a spare cartridge zipper bag and an instructions for use.
- Ensure that the lot number of the cartridge matches that of the ID chip.
- If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for AFIAS tests.
- Empty the tip box.
- Insert the ID chip into the 'ID chip port'.
- ✳ **Please refer to the instrument for AFIAS tests operation manual for complete information and operating instructions.**

## TEST PROCEDURE

- **AFIAS-1, AFIAS-3, AFIAS-6,**
- General mode**
- Insert the cartridge into the cartridge holder.
  - Insert a tip into the tip hole of the cartridge.
  - Select the 'General mode' in the instrument for AFIAS tests.
  - Take 100 µL of the sample (whole blood/serum/plasma/control) using a pipette and dispense it into the sample well of the cartridge.
  - Tap the 'Start' button on the screen.
  - The test result will be displayed on the screen after 3 minutes.
- C-tip Mode**
- Insert a cartridge into the cartridge holder.
  - Take 10 µL of whole blood using a C-tip.
  - Insert the whole blood-filled C-tip into the tip hole of the cartridge.
  - Select the 'C-tip mode' in the instrument for AFIAS tests.
  - Tap the 'Start' button on the screen.
  - The test result will be displayed on the screen after 3 minutes.
- **AFIAS-10**
- Normal mode**
- Insert a cartridge into the cartridge holder.

- 2) Insert a tip into the tip hole of the cartridge.
- 3) Tap the 'Load' button of the bay that holds the cartridge with the tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- 4) Insert the sample tube into the tube rack.
- 5) Insert the tube rack into the loading part of the sampling station.
- 6) Tap the 'Start' button on the screen.
- 7) The test result will be displayed on the screen after 3 minutes.

**Emergency mode (General tip)**

- 1) The test procedure is same with the 'Normal mode 1) - 3)'.
  - 2) Convert the 'Emergency mode' in AFIAS-10.
  - 3) Select the tip type (General tip) on the screen.
  - 4) Select the sample type (whole blood/serum/plasma) on the screen.
  - 5) Take 100 µL of the sample using a pipette and dispense it into the sample well of the cartridge.
  - 6) Tap the 'Start' button on the screen.
  - 7) The test result will be displayed on the screen after 3 minutes.

**Emergency mode (C-tip)**

- 1) Insert a cartridge into the cartridge holder.
- 2) Take 10 µL of whole blood using a C-tip.
- 3) Insert the C-tip with sample into the tip hole of the cartridge.
- 4) Tap the 'Load' button of the bay that holds the cartridge with a tip to read the barcode of the cartridge and please confirm the item name written on the cartridge.
- 5) Convert the 'Emergency mode' in AFIAS-10.
- 6) Select the tip type (C-tip) on the screen.
- 7) Tap the 'Start' button on the screen.
- 8) The test result will be displayed on the screen after 3 minutes.

**INTERPRETATION OF TEST RESULT**

- The instrument for AFIAS tests calculates the test result automatically and displays CRP concentration of the test sample in terms of mg/L.
- Cut-off : 10 mg/L
- Working range : 0.5-200 mg/L
- Effect of Hematocrit

The CRP Whole Blood of the instrument for AFIAS tests is calibrated to read the CRP serum concentration of a blood sample with a hematocrit of 40 %. If the actual hematocrit value deviates from 40 %, the result should be corrected by multiplying with the respective factor in the table: deviates from 40 %, the result should be corrected by multiplying with the respective factor in the table:

Hct (%)	Factor	Hct (%)	Factor
20-29	0.8	56-58	1.4
30-36	0.9	59-61	1.5
37-42	1.0	62-63	1.6
43-47	1.1	64-65	1.7
48-51	1.2	66-67	1.8
52-55	1.3	68-69	1.9

**Reference range, HCT:**

- Woman: 35-44 %
- Men: 39-48 %

**QUALITY CONTROL**

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand with **AFIAS CRP**. For more information regarding obtaining the control materials, contact **Boditech Med Inc.'s Sales Division for assistance**. (Please refer to the instructions for use of control material.)

**PERFORMANCE CHARACTERISTICS**

- **Analytical sensitivity**
  - Limit of Blank (LoB) 0.19 mg/L
  - Limit of Detection (LoD) 0.32 mg/L
  - Limit of Quantitation (LoQ) 0.50 mg/L
- **Analytical Specificity**
  - **Cross-reactivity**  
 Biomolecules such as below the ones in the table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS CRP** test results did not show any significant cross-reactivity with these biomolecules.

Cross reactants	Concentration
IL-6	50 µg/mL
PCT	50 µg/mL
Serum amyloid P	500 µg/mL
Serum amyloid A	500 µg/mL
Ferritin	500 µg/mL

- **Interference**  
 Interferents listed in the following table were added to the test sample at the concentration mentioned below. **AFIAS CRP** test results did not show any significant interference with these materials.

Interferents	Concentration
Ascorbic acid	350 µmol/L
Bilirubin (conjugated)	475 µmol/L
Albumin	60 g/L
Glucose	1,000 mg/dL
Triglyceride mixture (Triglyceride, total)	1,500 mg/dL
Hemoglobin	10 g/L
Biotin	3,500 ng/ml

- **Precision**
  - **Single site study**  
 3 Lots of **AFIAS CRP** were tested for 20 days. Each standard material was tested 2 times per day. For each test, each material was duplicated.  
 Repeatability (within-run precision)  
 Within-laboratory precision (Total precision)  
 Lot to lot precision
  - **Between persons**  
 Three different persons tested one lot of **AFIAS CRP**, ten times at each concentration of the control standard.
  - **Between sites**  
 One person tested **AFIAS CRP** at three different sites, ten times at each concentration of the control standard.
  - **Between readers**  
 One person tested **AFIAS CRP** with three different readers, ten times at each concentration of the control standard.

Expected Value [mg/L]	Repeatability			Within-laboratory precision		
	Mean [mg/L]	SD	CV (%)	Mean [mg/L]	SD	CV (%)
2.0	1.97	0.16	8.0	2.02	0.17	8.5
10.0	9.96	0.63	6.3	10.01	0.61	6.1
150.0	148.76	13.95	9.4	149.33	13.79	9.2

Expected Value [mg/L]	Lot to lot precision			Between site		
	Mean [mg/L]	SD	CV (%)	Mean [mg/L]	SD	CV (%)
2.0	2.00	0.18	8.9	2.04	0.16	7.6
10.0	9.99	0.6	6.0	10.02	0.71	7.1
150.0	149.04	14.39	9.65	152.32	9.09	6.0


Expected Value [mg/L]	Between person			Between reader		
	Mean [mg/L]	SD	CV (%)	Mean [mg/L]	SD	CV (%)
2.0	2.02	0.13	6.5	2.00	0.18	9.0
10.0	10.03	0.81	8.1	9.93	0.70	7.0
150.0	155.47	11.92	7.7	152.49	11.79	7.7

- **Accuracy**  
 The accuracy was confirmed by testing 3 different lots, three times with international standard material (ERM-DA474/IFCC, human serum).

	Lot 1		Lot 2		Lot 3	
	Measured Value [mg/L]	Bias (%)	Measured Value [mg/L]	Bias (%)	Measured Value [mg/L]	Bias (%)
Run 1	41.82	1.5	39.88	-3.2	42.68	3.6
Run 2	40.09	-2.7	41.42	0.5	40.52	-1.7
Run 3	42.06	2.1	41.96	1.8	40.95	-0.6

- **Comparability**  
 CRP concentration of 100 clinical samples were quantified independently with **AFIAS CRP (AFIAS-6)** and **comparator A** as per prescribed test procedures. Test results were compared, and their comparability was investigated with linear regression and correlation coefficient (R). The regression equation and correlation coefficient are as follows.

**REFERENCES**

1. Pepys MB and Hirschfield GM. C-reactive protein: a critical update. J Clin. Invest 2003; 111:1805-1812.
2. Volanakis JE. Human C-reactive protein: expression, structure, and function. Mol Immunol 2001;38:189-197.
3. Koenig W, Sund M, Frohlich M, et al. C-reactive protein, a sensitive marker of inflammation, predicts future risk of coronary heart disease in initially healthy middle-aged men. Circulation 1999; 99:237-242.
3. Rifai N, Ridker PM. Proposed Cardiovascular Risk Assessment Algorithm Using High-Sensitivity C-reactive protein and Lipid Screening. Clin. Chem. 2001; 47:28-30
4. Rifai N and Ridker PM. High-Sensitivity C-Reactive Protein: A novel and Promising Marker of Coronary Heart Disease. Clin. Chem.

- 2001; 47(3): 403-411.
5. Biasucci LM, Liuzzo G, Grillo RL, et al. Elevated levels of C-reactive protein at discharge in patients with unstable angina predict recurrent instability. Circulation 1999; 99:855-860
6. Taubes G. Does inflammation cut to the heart of the matter Science 2002; 296:242-245
7. Ridker PM, Hennekens CH, Buring JE, and Rifai N. C-reactive protein and other markers of inflammation in the prediction of cardiovascular disease in women. N Engl J Med 2000;342(12): 836-843
8. Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point-of-care testing. Clin Chem 1999; 45:1676-1678
9. Oh SW, Moon JD, Park SY, et al. Evaluation of fluorescence hs-CRP immunoassay for point-of-care testing. Clin Chim Acta 2005; 356:172-177
10. Claus DR, Osmond AP, Gewurz H. Radioimmunoassay of human C-reactive protein and levels in normal sera. J. Lab. Clin Med 1976;87:120-128
11. Kindmark CO. The concentration of C-reactive protein in sera from healthy individuals. Scand J Clin Lab Invest 1972;29:407-411

**Note:** Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For technical assistance, please contact:  
**Boditech Med Inc.'s Technical Services**  
 Tel: +(82) -33-243-1400  
 E-mail: sales@boditech.co.kr

**Boditech Med Inc.**  
 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, Republic of Korea  
 Tel: +(82) -33-243-1400  
 Fax: +(82) -33-243-9373  
 www.boditech.co.kr

**EC REP Obelis s.a**  
 Bd. Général Wahis 53,1030 Brussels, Belgium  
 Tel: +(32) -2-732-59-54  
 Fax: +(32) -2-732-60-03  
 E-Mail: mail@obelis.net

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
**LITERATŪROS ŠARŠAS**

1. Pepys MB and Hirschfield GM. C-reactive protein: a critical update. J Clin. Invest 2003; 111:1805-1812.
2. Volanakis JE. Human C-reactive protein: expression, structure, and function. Mol Immunol 2001;38:189-197.
3. Koenig W, Sund M, Frohlich M, et al. C-reactive protein, a sensitive marker of inflammation, predicts future risk of coronary heart disease in initially healthy middle-aged men. Circulation 1999; 99:237-242.
3. Rifai N, Ridker PM. Proposed Cardiovascular Risk Assessment Algorithm Using High-Sensitivity C reactive protein and Lipid Screening. Clin. Chem. 2001; 47:28-30
4. Rifai N and Ridker PM. High-Sensitivity C-Reactive Protein: A novel and Promising Marker of Coronary Heart Disease. Clin. Chem. 2001; 47(3): 403-411.
5. Biasucci LM, Liuzzo G, Grillo RL, et al. Elevated levels of C-reactive protein at discharge in patients with unstable angina predict recurrent instability. Circulation 1999; 99:855-860
6. Taubes G. Does inflammation cut to the heart of the matter Science 2002; 296:242-245
7. Ridker PM, Hennekens CH, Buring JE, and Rifai N. C-reactive protein and other markers of inflammation in the prediction of cardiovascular disease in women. N Engl J Med 2000;342(12): 836-843
8. Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point-of-care testing. Clin Chem 1999; 45:1676-1678
9. Oh SW, Moon JD, Park SY, et al. Evaluation of fluorescence hs-CRP immunoassay for point-of-care testing. Clin Chim Acta 2005; 356:172-177
10. Claus DR, Osmond AP, Gewurz H. Radioimmunoassay of human C-reactive protein and levels in normal sera. J. Lab. Clin Med 1976;87:120-128
11. Kindmark CO. The concentration of C-reactive protein in sera from healthy individuals. Scand J Clin Lab Invest 1972;29:407-411

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# Boditech CRP Control

## INTENDED USE

Boditech CRP Control is intended for the quality control of CRP assay kits manufactured/provided by Boditech Med Inc. For *in vitro* diagnostic use only.

## COMPONENTS

Boditech CRP Control consists of 'Boditech CRP Control Level 1', 'Boditech CRP Control Level 2', 'Instructions for use' and 'Control value & Barcode sheet'.

- **Boditech CRP Control is provided in liquid form.** 4 p.

- The control contains CRP antigen and horse serum.
- **The control materials are contained in vials, and the vials are further packaged in a box.**

## SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech CRP Control should not be used past the expiration date.
- Boditech CRP Control is solely designed for the quality control of CRP cartridges manufactured/provided by Boditech Med Inc.
- Human-derived materials in Boditech CRP Control were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be non-reactive. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human-derived materials and clinical samples should be handled as though they are capable of transmitting infectious diseases and should be disposed of as hazardous wastes.
- All waste materials should be disposed of in accordance with the requirements of your local waste management authorities.
- Boditech CRP Control contains sodium azide (NaN<sub>3</sub>), and they may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure. Avoid contact with skin, eyes, and clothing. In case of contact, rinse immediately with running water.

## STORAGE AND STABILITY

- Storage and stability condition of Boditech CRP Control.

	Unopened	Opened (After reconstitution)
<b>Temperature</b>	2 ~ 8 °C	2 ~ 8 °C
<b>Expiration date</b>	Until expiration date on the label	28 days

- Close the opened vial tightly after use.
- After use, any residual substance should not be returned to the original vial.
- Bacterial contamination of opened Boditech CRP Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be used.

## TEST PROCEDURE

Boditech CRP Control is supplied in liquid form.

1. Thoroughly mix the contents of the vial before each use by gently inverting for several times.
2. When removed from the refrigerator, allow the vial to stand for about 30 minutes before using the control.

**Please refer to the instructions for use of the test cartridges for detailed test procedure.**

In the event of damage to the package, contact the **Boditech Med Inc.'s Technical Sales.**

## MATERIALS SUPPLIED

**REF** CFPO-100

Boditech CRP Control Box (2 vials)

- Boditech CRP Control Level 1 1
- Boditech CRP Control Level 2 1
- Instructions for use 1
- Control value & Barcode sheet 1

## INTERPRETATION OF THE RESULT

The test result of the 'Boditech CRP Control' should be consistent with the expected result of control value sheet.

If the test results fall outside the expected result, check the following for potential sources of error, and retest after resolving these matters. If the error still persists, contact the **Boditech Med Inc.'s Technical Sales.**

- ※ Potential sources of error
  - Errors in a testing process
  - Incorrect storage condition of Boditech CRP Control
  - Use of expired or contaminated Boditech CRP Control.
  - Faulty Boditech's CRP assay kits.
  - Faulty Boditech's instruments.

## QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.

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	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices

For technical assistance, please contact:

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# Boditech CRP Calibrator

## INTENDED USE

Boditech CRP Calibrator is intended for *in vitro* diagnostic use in the calibration of CRP Assay Kit.

**For *in vitro* diagnostic use only.**

## INTRODUCTION

The use of Boditech CRP Calibrator may be considered as an objective assessment of the precision of CRP Assay Kits and is an integral part of Good Laboratory Practices. Boditech CRP Calibrator is provided in liquid form.

## COMPONENTS

Boditech CRP Calibrator consists of 'Boditech CRP Calibrator level 1', 'Boditech CRP Calibrator level 2', 'Instruction for Use' and 'Barcode Sheet'.

- The calibrator contains CRP antigen and Horse serum
- Each control vial packed in a box.

## SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech CRP Calibrator should not be used past the expiration date.
- Boditech CRP Calibrator is solely designed to be provided instrument-specific calibration curves of Boditech Readers and CRP Assay Kits.
- When a certain parameter of the instrument for AFIAS test is calibrated from Cal/QC mode, the system displays fluorescence intensity of each test instead of its actual result.
- Human source materials from which Boditech CRP Calibrator is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

## STORAGE AND STABILITY

- Opened: Store refrigerated (+2 to +8 °C). Once reconstituted, stable for 28 days at (+2 to +8 °C) if kept capped in original container and free from contamination. After use, any residual product should not be returned to the original tube.
- Unopened: Store refrigerated (+2 to +8°C).
- Boditech CRP Calibrator is stable to expiration date on the label.
- Bacterial contamination of reconstituted Boditech CRP Calibrator will cause reductions in the stability of many components. If bacterial contamination is suspected, the tube should be discarded and a fresh tube needs to be reconstituted.

## INSTRUCTIONS FOR USE

Boditech CRP Calibrator is supplied in liquid form.

1. Thoroughly mix the contents of the tube before each use by gently inverting for several times.
2. Follow the procedure according to the instruction for use provided with the kit.

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For Technical Assistance  
**Boditech Med Inc.'s Technical Services** at  
Tel: +82 (33) 243-1400  
E-mail: sales@boditech.co.kr

Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon,  
Chuncheon-si, Gang-won-do, 24398  
Republic of Korea  
Tel: +82 -33-243-1400 / Fax: +82 -33-243-9373  
www.boditech.co.kr

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
**LITERATŪROS ŠARŠAS**

- Goldstein DE, Little RR, Lorenz RA, Malone JI, Nathan D, Peterson CM. Tests of glycemia in diabetes. Diabetes Care 1995; 18:896-909.
- Bunn HF. Nonglycosylated glycosylation of protein: relevance to diabetes. Am J Med 1981; 70:325-30.
- Jovanovic L, Peterson CM. The clinical utility of glycosylated hemoglobin. Am J Med 1981; 70:331-8.
- Nathan DM, Singer DE, Hurxthal K, Goodson JD. The clinical information value of the glycosylated hemoglobin assay. N Engl J Med 1984; 310:341-6.
- Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. Clin Chem 1986; 32:B64-70.
- Goldstein DE, Little RR, England JD, Wiedmeyer H-M, McKenzie E. Methods of glycosylated hemoglobins: high performance liquid chromatography and thiobarbituric acid colorimetric methods. In: Clarke WL, Lerner J, Pohl SL, eds. Methods in diabetes research, Vol. 2. New York: John Wiley, 1986:475-504.
- Tahara Y, Shima K. The response of GHb to stepwise plasma glucose change over time in diabetic patients. Diabetes Care 1993; 16:1313-4.
- Svendsen PA, Lauritzen T, Soegaard U, Nerup J. Glycosylated haemoglobin and steady-state mean blood glucose concentration in type 1 (insulin-dependent) diabetes. Diabetologia 1982; 23:403-5.
- Cefalu WT, Wang ZQ, Bell-Farrow A, Kiger FD, Izlar C. Glycohemoglobin measured by automated affinity HPLC correlates with both short-term and long-term antecedent glycemia. Clin Chem 1994; 40:1317-21.

- Singer DE, Coley CM, Samet JH, Nathan DM. Tests of glycemia in diabetes mellitus. Their use in establishing a diagnosis and in treatment. Ann Intern Med 1989; 110:125-37.
- Molnar GD. Clinical evaluation of metabolic control in diabetes. Diabetes 1978; 27:216-25.
- UK Prospective Diabetes Study. Reduction in HbA1c with basal insulin supplement, sulfonylurea or biguanide therapy in maturity-onset diabetes. Diabetes 1985; 34:793-8.
- Baker JR, Johnson RN, Scott DJ. Serum fructosamine concentrations in patients with type II (non-insulin-dependent) diabetes mellitus during changes in management. BMJ (Clin Res Ed) 1984; 288:1484-6.
- Tahara Y, Shima K. Kinetics of HbA1c, glycated albumin, and fructosamine and analysis of their weight functions against preceding plasma glucose level. Diabetes Care 1995; 18:440-7.
- Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point of care testing. Clin Chem 1999; 45:1676-1678.

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2.3.2.2 p.

## INTENDED USE

**AFIAS HbA1c is a fluorescence Immunoassay (FIA) for the quantitative determination of HbA1c (Hemoglobin A1c) in human whole blood. It is useful as an aid in management and monitoring of the long-term glycemic status in patients with diabetes mellitus.**

For *in vitro* diagnostic use only.

## INTRODUCTION

Glycated protein is formed post-translationally through the slow, nonenzymatic reaction between glucose and amino groups on proteins. HbA1c is a clinically useful index of mean glycemia during the preceding 120 days, the average life span of erythrocytes. Carefully controlled studies have documented a close relationship between the concentrations of HbA1c and mean glycemia. HbA1c is considered as a more reliable parameter in monitoring glycemia over the glycemic reading with the conventional glucometer.

## PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibody in buffer binds to antigen in sample, forming antigen-antibody complexes, and migrates onto nitrocellulose matrix to be captured by the other immobilized-antibody on test strip.

The more antigen in sample forms the more antigen-antibody complex and leads to stronger intensity of fluorescence signal on detector antibody, which is Instrument for the instrument for AFIAS tests displays the content of glycated hemoglobin in terms of percent of the total hemoglobin in blood.

## COMPONENTS

**AFIAS HbA1c** consists of 'Cartridge', 'C- tip', 'ID chip' and 'Instruction for use'.

- Each sealed aluminum pouch contains two cartridges.
- Each cartridge packaged in an aluminum pouch has three components, detection buffer part, hemolysis buffer part and cartridge part.
- Cartridge part contains a test strip, the membrane which has anti human HbA1c at the test line, while rabbit IgG at the control line.
- Detection buffer part contains anti human HbA1c-fluorescence conjugate, anti rabbit IgG-fluorescence conjugate, bovine serum albumin (BSA) as a stabilizer and sodium azide in phosphate buffered saline (PBS) as a preservative.
- Hemolysis buffer part contains non-ionic detergent and sodium azide in phosphate buffered saline (PBS) as a preservative.
- C-tip (Capillary tip) is a useful tool for point of care testing which requires a small volume of capillary blood from

fingertip, heel site (in infants) or ear-lobe.

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this 'Instruction for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect of test result(s).
- Do not reuse cartridges. A cartridge should be used for testing one sample only.
- The cartridge should remain sealed in its aluminum pouch until just before use. Do not use the cartridge, if pouch is damaged.
- For shipping, samples must be packed in accordance with local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used and should be recollected.
- Allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used C-tips, and cartridges should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
- AFIAS HbA1c** will provide accurate and reliable results subject to the below conditions.
  - AFIAS HbA1c** should be used only in conjunction with the instrument for AFIAS tests.
  - Have to use recommended anticoagulant sample.

Sample type	Recommended anticoagulant
Whole blood	Na <sub>2</sub> EDTA, K <sub>2</sub> EDTA, K <sub>3</sub> EDTA, Sodium citrate, Lithium heparin

- C-tip should be used when the following conditions are met.**

- C-tip provided with the kit is recommended to obtain correct test result.
- Capillary blood should be immediately tested after collection.
- Do not leave C-tip after collection of blood, test immediately.
- Excess capillary blood around the C-tip should be wiped off.
- In order to avoid cross-contamination, please do not re-use C-tip for multiple samples.
- AFIAS cartridge should be inserted and positioned in the cartridge holder prior to the blood sample collection.

## LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample

components to the capture/detector antibodies.

- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

## STORAGE AND STABILITY

- The cartridge is stable for 20 months (while sealed in an aluminum pouch) if stored at 2-8 °C.
- Return an unused cartridge to the spare cartridge zipper bag containing the desiccant pack. Reseal along entire edge of zip-seal. May be stored for up to 1 month at 2-8°C.

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Following items can be purchased separately from **AFIAS HbA1c**

Please contact our sales division for more information.

- AFIAS-1** **REF** FPRR019
- AFIAS-6** **REF** FPRR020
- Boditech HbA1c Control** **REF** CFPO-96
- Boditech HbA1c Calibrator** **REF** CFPO-108

## SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS HbA1c** is human whole blood.

- It is recommended to test the sample within 12 hours after collection.
- Samples may be stored for up to a week at 2-8 °C prior to being tested.
- If testing will be delayed more than a week, samples should be frozen at -70 °C or below. Samples stored frozen at -70 °C or below for 3 months showed no performance difference.
- Once the sample was frozen, it should be used one time only for test, because repeated freezing and thawing can

result in the change of test values.

- Collection of capillary blood sample using C-tip (Fingertip blood)
  - Clean the area with a pre-injection swab.
  - Pierce with a sterile lancet.
  - Wipe away first drop of blood.
  - Gently massage the surrounding area towards a C-tip for a second drop.
  - Hold a C-tip horizontally and touch the tip of C-tip to the blood drop.
  - Capillary action will automatically draw the blood sample to C-tip and stop
  - Wipe off any excess blood around the tip.
  - Double-check if capillary blood is fully filled in the C-tip and AFIAS reader is ready for a test on the 'C-tip mode'.

## TEST SETUP

- Check the components of the **AFIAS HbA1c** as described below.: Cartridge, C-tip, ID chip and instruction for use.
- Keep the sealed cartridge (if stored in refrigerator) at room temperature for at least 30 minutes just prior to the test. Place the cartridge on a clean, dust-free and flat surface.
- Turn on the instrument for AFIAS tests.
- Empty the tip box.
- Insert the ID chip into the "ID chip port".
- Please refer to the instrument for AFIAS tests 'Operation Manual' for complete information and operating instructions.

## TEST PROCEDURE

- Select "Test Mode" in the instrument for AFIAS tests.
  - General Mode : Tube blood (with anticoagulant)
  - C-tip Mode : Fingertips blood (non-anticoagulant blood)
  - You must use C-tip only to collect test samples in both test mode (General / C-tip mode).**
  - Do not use a general pipette tip to collect test samples.**
- Insert the cartridge into the cartridge holder.
- Take 10 μL of whole blood (tube blood or fingertip blood) with a C-tip. (Please refer to "SAMPLE COLLECTION AND PROCESSING")
- Insert the whole blood-filled C-tip into the tip hole of the cartridge.
- Tap the 'START' icon on the screen.
- The test result will be displayed on the screen after 10 minutes.

※ Note: Refer to the instrument for AFIAS tests Operation Manual to select sample type.

## INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculate the test result automatically and displays HbA1c concentration of the test sample in terms of percentage (%), NGSP) and mmol/mol (IFCC).
- The reference value
  - NGSP (%): 4.5-6.5 %
  - IFCC (mmol/mol): 26-48 mmol/mol

- The working range
  - NGSP (%): 4-15 %
  - IFCC (mmol/mol): 20.2-140.4 mmol/mol
  - eAG (mg/dL): 68.1-383.8 mg/dL

**QUALITY CONTROL**

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are not provided with **AFIAS HbA1c**. For more information regarding obtaining the control materials, contact **Boditech Med Inc.'s Sales Division** for assistance. (Please refer to the instruction for use of control material.)

**PERFORMANCE CHARACTERISTICS**

- Analytical sensitivity**
  - Limit of Blank (LoB) 2.26%
  - Limit of Detection (LoD) 2.39%
  - Limit of Quantification (LoQ) 3.0 %

- Analytical Specificity**
  - Cross-reactivity
  - There was no significant cross-reactivity from these materials with the **AFIAS HbA1c** test measurements.

cross-reactivity material	Standard material conc. (%)		
	5.3	8.1	13.5
	Recovery (%)		
HbA0 (2 mg/mL)	98.9	102.0	101.0
HbA1a&A1b mixture (5 mg/mL)	100.9	101.9	100.1
Acetylated hemoglobin (100 mg/mL)	100.7	100.2	100.5
Carbamylated hemoglobin (100 mg/mL)	100.2	101.2	100.6
Glycated h-Albumin (2.2 mg/mL)	100.1	100.6	102.0
HbA1d (100 mg/mL)	101.6	101.3	100.5
Acetylaldehyde hemoglobin (100 mg/mL)	99.8	101.4	100.7

- Cross-reactivity
- There was no significant interference from these materials with the **AFIAS HbA1c** test measurements.

Interference material	Standard material conc. (%)		
	5.3	8.1	13.5
	Recovery (%)		
Non-interference	100.0	99.8	99.5
Acetaminophen (1.5 mmol/L)	100.3	100.6	99.6
L-Ascorbic acid (5 mg/dL)	100.1	100.6	100.8
Bilirubin (20 mg/mL)	100.3	100.3	100.0
D-Glucose (500 mg/dL)	99.8	100.6	102.0
Intralipid (1 mg/mL)	101.0	99.7	100.5
Triglyceride (3000 mg/mL)	100.4	99.6	99.9

Urea (100 mmol/L)	100.3	100.0	100.8
-------------------	-------	-------	-------

- Precision**
  - Between lot  
One person tested three different lots of **AFIAS HbA1c**, five times at each concentration of the control standard.
  - Between reader  
Three different persons tested **AFIAS HbA1c** at three different reader; five times at each concentration of the control standard.
  - Between day  
One person tested **AFIAS HbA1c** for five days; five times at each concentration of the control standard.
  - Between site  
One person tested **AFIAS HbA1c** at three different sites; five times at each concentration of the control standard.


HbA1c (%)	Between lot		Between reader		Between day		Between site	
	Avg.	CV(%)	Avg.	CV(%)	Avg.	CV(%)	Avg.	CV(%)
5.3	5.21	2.7	5.25	3.6	5.31	2.5	5.27	2.3
8.1	8.15	3.1	8.12	2.7	8.20	2.9	8.23	3.0
13.5	13.48	2.8	13.38	3.3	13.60	2.8	13.57	2.8

- Accuracy**  
The accuracy was confirmed each different concentration.

HbA1c (%)	Between lot		Between reader	
	Avg.	Bias(%)	Avg.	Bias(%)
5.3	5.21	-0.02	5.25	-0.01
8.1	8.15	0.01	8.12	0
13.5	13.48	0	13.38	-0.01

HbA1c (%)	Between day		Between site	
	Avg.	Bias(%)	Avg.	Bias(%)
5.3	5.31	0	5.27	-0.01
8.1	8.2	0.01	8.23	0.02
13.5	13.6	0.01	13.57	0.01

- Comparability**  
HbA1c concentrations of 165 clinical samples were quantified independently with **AFIAS HbA1c** and Tosoh G7 HPLC (Tosoh Bioscience Inc. Japan) as per prescribed test procedures. Test results were compared and their comparability was investigated with linear regression and coefficient of correlation (R). Linear regression and coefficient of correlation between the two tests were  $Y = 0.997X + 0.0309$  and  $R = 0.9982$  respectively.

**REFERENCES**

- Goldstein DE, Little RR, Lorenz RA, Malone JI, Nathan D, Peterson CM. Tests of glycemia in diabetes. Diabetes Care 1995; 18:896-909.
- Bunn HF. Nonenzymatic glycosylation of protein: relevance to diabetes. Am J Med 1981; 70:325-30.
- Jovanovic L, Peterson CM. The clinical utility of glycosylated hemoglobin. Am J Med 1981; 70:331-8.
- Nathan DM, Singer DE, Hurxthal K, Goodson JD. The clinical information value of the glycosylated hemoglobin assay. N Engl J Med 1984; 310:341-6.
- Goldstein DE, Little RR, Wiedmeyer HM, England JD, McKenzie EM. Glycated hemoglobin: methodologies and clinical applications. Clin Chem 1986; 32:B64-70.
- Goldstein DE, Little RR, England JD, Wiedemeyer H-M, McKenzie E. Methods of glycosylated hemoglobins: high performance liquid chromatography and thiobarbituric acid colorimetric methods. In: Clarke WL, Larner J, Pohl SL, eds. Methods in diabetes research, Vol. 2. New York: John Wiley, 1986:475-504.
- Tahara Y, Shima K. The response of GHb to stepwise plasma glucose change over time in diabetic patients. Diabetes Care 1993; 16:1313-4.
- Svendsen PA, Lauritzen T, Soegaard U, Nerup J. Glycosylated haemoglobin and steady-state mean blood glucose concentration in type 1 (insulin-dependent) diabetes. Diabetologia 1982; 23:403-5.
- Cefalu WT, Wang ZQ, Bell-Farrow A, Kiger FD, Izlar C. Glycohemoglobin measured by automated affinity HPLC correlates with both short-term and long-term antecedent glycemia. Clin Chem 1994; 40:1317-21.
- Singer DE, Coley CM, Samet JH, Nathan DM. Tests of glycemia in diabetes mellitus. Their use in establishing a diagnosis and in treatment. Ann Intern Med 1989; 110:125-37.
- Molnar GD. Clinical evaluation of metabolic control in diabetes. Diabetes 1978; 27:216-25.
- UK Prospective Diabetes Study. Reduction in HbA1c with basal insulin supplement, sulfonylurea or biguanide therapy in maturity-onset diabetes. Diabetes 1985; 34:793-8.
- Baker JR, Johnson RN, Scott DJ. Serum fructosamine concentrations in patients with type II (non-insulin-dependent) diabetes mellitus during changes in management. BMJ (Clin Res Ed) 1984; 288:1484-6.
- Tahara Y, Shima K. Kinetics of HbA1c, glycated albumin, and fructosamine and analysis of their weight functions against preceding plasma glucose level. Diabetes Care 1995; 18:440-7.
- Brooks DE, Devine DV, Harris PC, et al. RAMP(TM): A rapid, quantitative whole blood immunochromatographic platform for point of care testing. Clin Chem 1999; 45:1676-1678.

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For technical assistance; please contact:  
**Boditech Med Inc.'s Technical Services**  
 Tel: +82 33 243-1400  
 E-mail: sales@boditech.co.kr

**Boditech Med Incorporated**  
 43, Geodudanji 1-gil, Dongnae-myeon,  
 Chuncheon-si, Gang-won-do, 24398  
 Republic of Korea  
 Tel: +(82) -33-243-1400  
 Fax: +(82) -33-243-9373  
 www.boditech.co.kr

**Obelis s.a**  
 Bd. Général Wahis 53,  
 1030 Brussels, BELGIUM  
 Tel: +(32) -2-732-59-54  
 Fax: +(32) -2-732-60-03  
 E-Mail: mail@obelis.net

# Boditech HbA1c Control

## INTENDED USE

Boditech HbA1c Control is intended for *in vitro* diagnostic use in the quality control of HbA1c Assay Kit. **For *in vitro* diagnostic use only.**

## INTRODUCTION

The use of Boditech HbA1c Control may be considered as an objective assessment of the precision of HbA1c Assay Kits and is an integral part of Good Laboratory Practices. Boditech HbA1c Control is provided in lyophilized form.

## COMPONENTS

Boditech HbA1c Control consists of 'Boditech HbA1c Control level 1', 'Boditech HbA1c Control level 2', 'Instruction for Use' and 'Barcode Sheet'.

- The control contains HbA1c stock and Horse serum
- Each control vial packed in a box. **4 p.**

## SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech HbA1c Control should not be used past the expiration date.
- Boditech HbA1c Control is solely designed to be provided instrument-specific calibration curves of Boditech Readers and HbA1c Assay Kits.
- Human source materials from which Boditech HbA1c Control is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

## STORAGE AND STABILITY

- Storage and stability condition of Boditech HbA1c Control.

	Unopened	Opened (after reconstitution)	
Temperature	+2 to +8 °C	+2 to +8 °C	-20 to -80 °C
Expiration date	Until expiration date on the label.	7 days	20 days

- Close the opened Boditech HbA1c Control bottle tightly after use.
- Once the Boditech HbA1c Control was frozen, it should be used ONE TIME ONLY for test, because repeated freezing and thawing can result in the change of test values.
- After use, any residual product should NOT BE RETURNED to the original vial.
- Bacterial contamination of reconstituted Boditech HbA1c Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

## INSTRUCTIONS FOR USE

Boditech HbA1c Control is supplied in lyophilized form.

1. Carefully reconstitute each vial of lyophilized with exactly 0.5 mL of sterilized distilled water.
2. Close the bottle and allow to stand for 30 minutes before use. Ensure contents are completely dissolved by swirling gently. Avoid formation of foam. Do not shake.

Please refer to package inserts of the test cartridges for detailed test procedure.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

In the event of damage to the package, contact the **Boditech Med Inc.'s Technical Services.**

## MATERIALS SUPPLIED

### REF CFPO-96

Boditech HbA1c Control Box (2 vials)		
Boditech HbA1c Control level 1 (0.5 mL)	1	
Boditech HbA1c Control level 2 (0.5 mL)	1	
Instruction For Use		1
Control value & Barcode Sheet		1

## QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.

For Technical Assistance

**Boditech Med Inc.'s Technical Services** at

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# Boditech HbA1c Calibrator

## INTENDED USE

Boditech HbA1c Calibrator is intended for *in vitro* diagnostic use in the calibration of HbA1c Assay Kit.

**For *in vitro* diagnostic use only.**

## INTRODUCTION

The use of Boditech HbA1c Calibrator may be considered as an objective assessment of the precision of HbA1c Assay Kits and is an integral part of Good Laboratory Practices. Boditech HbA1c Calibrator is provided in lyophilized form.

## COMPONENTS

Boditech HbA1c Calibrator consists of 'Boditech HbA1c Calibrator level 1', 'Boditech HbA1c Calibrator level 2', 'Instruction for Use' and 'Barcode Sheet'.

- The Calibrator contains HbA1c stock and Horse serum
- Each calibrator vial packed in a box.

## SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech HbA1c Calibrator should not be used past the expiration date.
- Boditech HbA1c Calibrator is solely designed to be provided instrument-specific calibration curves of Boditech Readers and HbA1c Assay Kits.
- When a certain parameter of the instrument for AFIAS test is calibrated from Cal/QC mode, the system displays fluorescence intensity of each test instead of its actual result.
- Human source materials from which Boditech HbA1c Calibrator is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

## STORAGE AND STABILITY

- Storage and stability condition of Boditech HbA1c Calibrator.

	Unopened	Opened (After reconstitution)	
<b>Temperature</b>	+2 to +8 °C	+2 to +8 °C	-20 to -80 °C
<b>Expiration date</b>	Until expiration date on the label.	7 days	20 days

- Close the opened Boditech HbA1c Calibrator bottle tightly after use.
- Once the Boditech HbA1c Calibrator was frozen, it should be used ONE TIME ONLY for test, because repeated freezing and thawing can result in the change of test values.
- After use, any residual product should NOT BE RETURNED to the original vial.
- Bacterial contamination of reconstituted Boditech HbA1c calibrator will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

## INSTRUCTIONS FOR USE

Boditech HbA1c Calibrator is supplied in lyophilized form.

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## MATERIALS SUPPLIED

### REF CFPO-108

Boditech HbA1c Calibrator Box (2 vials)

Boditech HbA1c Calibrator level 1 (0.5 mL)	1
Boditech HbA1c Calibrator level 2 (0.5 mL)	1
Instruction for Use	1
Calibrator value & Barcode Sheet	1

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2.3.2.3 p.

### INTENDED USE

**AFIAS Tn-I Plus is a fluorescence Immunoassay (FIA) for the quantitative determination of cardiac troponin-I (Tn-I) in human whole blood/serum/plasma. It is useful as an aid in management and monitoring of acute myocardial infarction (AMI).  
 For *in vitro* diagnostic use only.**

### INTRODUCTION

Cardiac troponins are currently the most sensitive and specific biochemical markers of myocardial necrosis. There are three types of troponin in heart muscle fibers. Those are troponin-C, -I, and -T. Together they contribute to make cardiac muscle fibers contract. The clinical measurement of serum Tn-I has become an important tool in the diagnosis of acute myocardial infarction. Serum Tn-I is a more reliable than creatine kinase as a prognostic marker in people with ischemic chest pain. National and international scientific organizations have suggested the use of troponins, Tn-I and Tn-T, when implementing new diagnostic strategies in patients with acute coronary syndrome.

### PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibodies in buffer bind to antigens in sample, forming antigen-antibody complexes and migrate onto nitrocellulose matrix to be captured by the other immobilized antibodies on the strip.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by instrument for AFIAS tests to show Tn-I concentration in the sample.

### COMPONENTS

**AFIAS Tn-I Plus** consists of a 'cartridges'.

- Each cartridge packaged in an aluminum pouch has three components including cartridge part, detector part and diluent part.
- Cartridge part contains the membrane called a test strip, which has streptavidin at the test line, and chicken IgY at the control line.
- Detector part contains anti Tn-I-fluorescence conjugate, anti Tn-I-biotin conjugate, anti chicken IgY-fluorescence conjugate, bovine serum albumin (BSA), sucrose, bromophenol blue, and sodium azide as a preservative.
- Diluent part contains tween 20, CA-630, and sodium azide as a preservative.

### WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow the instructions and procedures described in this 'Instruction for use'.
- Use fresh samples only and avoid direct sunlight.
- It is possible to use frozen samples. Please refer to "SAMPLE COLLECTION AND PROCESSING."

- Lot numbers of all the test components (cartridge and ID chip) must match each other.
- Do not interchange test components between different lots or use test components after the expiration date, either of which might yield incorrect test result(s).
- The cartridge should remain sealed in its aluminum pouch until just before use. Do not use the cartridge if the pouch found damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with local requirement. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- Allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used cartridges and pipette tips should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
- No Biotin interference was observed in **AFIAS Tn-I Plus** when biotin concentration in the sample was below 5 ng/mL. If a patient has been taking biotin at dosage of more than 0.03 mg a day, it is recommended to test again 24 hours after discontinuation of biotin intake.

- AFIAS Tn-I Plus** will provide accurate and reliable results subject to the below conditions.
  - AFIAS Tn-I Plus** should be used only in conjunction with the instrument for AFIAS tests.
  - Have to use recommended anticoagulant sample.

Recommended anticoagulant
Sodium heparin, Lithium heparin, Sodium citrate

### STORAGE AND STABILITY

Component	Storage condition		
	Storage Temperature	Shelf life	Note
Cartridge	2 - 8 °C	20 months	Unopened
		1 month	Resealed

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### LIMITATION OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause the false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors,

degradation of the test components/reagents or presence of interfering substances in the test samples.

- Any clinical diagnosis based on the test result should be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

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Following items can be purchased separately from **AFIAS Tn-I Plus**.

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### SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS Tn-I Plus** is human whole blood/serum/plasma.

- It is recommended to test the sample within 24 hours after collection.
- The serum or plasma should be separated from the clot by centrifugation within 3 hours after the collection of whole blood. If longer storage is required, e.g. if the test could not be performed within 24 hours, serum or plasma should be immediately frozen below -20 °C. The freezing storage of sample up to 3 months does not affect the quality of results.
- However, the whole blood sample should not be kept in a freezer in any case.
- Once the sample was frozen, it should be thawed one time and only for test, because repeated freezing and thawing can result in the changed test values.

### TEST SETUP

- Check the components of the **AFIAS Tn-I Plus** as described below: Cartridges, Pipette tips, ID chip, Cartridge spare zipper bag and Instruction for use
- If the sealed cartridge has been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Turn on the instrument for AFIAS tests
- Empty the tip box.
- Insert the ID chip into the "ID chip port".  
 (Please refer to the instrument for AFIAS tests 'Operation Manual' for complete information and operating instructions.

### TEST PROCEDURE

- Take 100 µL of sample with a pipette and dispense it into the sample well in the cartridge.
- Insert the cartridge into the cartridge holder
- Insert a tip into the tip hole of the cartridge.
- Tap the 'START' icon on the screen.
- The test result will be displayed on the screen after 12 minutes.

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### INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculates the test result automatically and displays Tn-I concentration of the test sample in terms of ng/mL.
- The working range of the **AFIAS Tn-I Plus** is 0.01-15.00 ng/mL.
- Expected Values**
  - In studies performed with the **AFIAS Tn-I Plus** assay involving 70 healthy volunteers in Korea, the upper reference limit (99th percentile) for Tn-I was 0.04 ng/mL. The lowest concentration with a CV less than or equal to 10 % with the **AFIAS Tn-I Plus** assay was 0.04 ng/mL.
  - Due to the release kinetics of Tn-I, a result below the decision limit within the first hours of the onset of symptoms does not rule out myocardial infarction with certainty. If myocardial infarction is still suspected, repeat the test at appropriate intervals.
  - A cut-off of 0.3 ng/mL Tn-I is recommended for diagnosis of AMI, as this yields optimal performance of 91 % of sensitivity and 92.1 % of specificity. However, laboratories should establish their own diagnostic cut-off concentration based on the clinical practice at their perspective institutions.

### QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are not provided with AFIAS Tn-I Plus. For more information regarding obtaining the control materials, contact [Boditech Med Inc.'s Sales Division for assistance](#). (Please refer to the instruction for use of control material.)

### PERFORMANCE CHARACTERISTICS

- Analytical sensitivity**
  - Limit of Blank (LoB): 0.004 ng/mL
  - Limit of Detection (LoD): 0.01 ng/mL
  - Limit of Quantitation (LoQ): 0.03 ng/mL

**Analytical specificity**

- Cross-reactivity

Biomolecules such as below the ones in the table were added to the test sample(s) at concentrations much higher than their normal physiological levels in the blood. **AFIAS Tn-I Plus** test results did not show any significant cross-reactivity with these biomolecules.

Cross-reactivity material	Conc. (ng/mL)
CK-MB	60
NT-proBNP	1,000
Myoglobin	1,000
D-Dimer	1,000

- Interference

Interference materials such as below the ones in the table were added to the test sample(s) the same as the below concentrations. **AFIAS Tn-I Plus** test results did not show any significant interference with these materials except for EDTA.

Interference material	Conc.
Bilirubin	350 μmol/L
Cholesterol	13 mmol/L
D-glucose	1,000 mg/dL
Hemoglobin	2 g/L
L-Ascorbic acid	350 μmol/L
Triglyceride mixture	500 mg/dL
Heparin	3000 U/L
Sodium citrate	2 mg/mL
EDTA	3.4 μmol/L
Biotin	5 ng/mL

**Precision**

- Intra-assay (Between lot)

One person tested three different lots of **AFIAS Tn-I Plus**, ten times at each concentration of the control standard.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.94	0.95	0.92
7.50	7.65	7.32	7.51
Conc. (ng/mL)	Mean	SD	CV (%)
0.23	0.23	0.013	5.5
0.94	0.93	0.057	6.1
7.50	7.50	0.444	5.9

- Inter-assay (Between person)

Three different persons tested **AFIAS Tn-I Plus**; five times at each concentration of the control standard.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.95	0.96	0.94
7.50	7.75	7.50	7.48
Conc. (ng/mL)	Mean	SD	CV (%)
0.23	0.23	0.013	5.4
0.94	0.95	0.056	5.9
7.50	7.58	0.400	5.3



**Accuracy**

The accuracy was confirmed by 3 different lots testing ten times each different concentration.

Mean of 3 lots			
Conc. (ng/mL)	Lot 1	Lot 2	Lot 3
0.23	0.23	0.23	0.23
0.94	0.95	0.91	0.93
7.50	7.56	7.35	7.37
Conc. (ng/mL)	Mean	CV (%)	Recovery(%)
0.23	0.23	5.8	101.0
0.94	0.93	5.3	99.0
7.50	7.43	5.7	99.0

**Comparability**

Using Access2 (Beckman Coulter) as a comparison reference device for **AFIAS Tn-I Plus**, 100 clinical samples were independently tested for its Tn-I concentration following each instrument's procedure. Results of both the test methods were analyzed and their comparability was investigated with linear regression and coefficient of correlation (R).

**REFERENCES**

1. Mauro Panteghini, Franca Pacani, Kiang-Teck J.Yeo, Fred S. Apple, Robert H. Christenson, Francesco Dati, Johannes Mair, Jan Ravkilde, and Alan H.B. We. Evaluation of Imprecision for Cardiac Troponin Assays at Low-Range Concentrations. 2004;50:2:327-332.
2. Alan McNeil, PhD, FRACP, FRCPA. The Trouble with Troponin. Heart, Lung and Circulation 2007;16:S13-S16.
3. David M. Bunk and Micahel J. Welch. Characterization of a New Certified Reference Material for Human Cardiac Troponin I. Clinical Chemistry 2002;52:2:212-219
4. Jaffe AS, Ravkilde J, Roberts R, Naslund U, Apple FS, Galvani M, Katus H. It's time for a change to a troponin standard. Circulation 2000;102:1216 -1220.
5. Jillan R. Tate, David Heathcote, Gus Koerbin, Gary Thean,

David Andriske, Jone Bonar, Janice Gill. The harmonization of cardiac troponin I measurement is independent of sample time collection but is dependent on the source of calibrator. Clinica Chimica Acta 324:2002:13-23.

6. Ohman EM, Armstrong PW, Christenson RH, et al. Cardiac troponin T levels for risk stratification in acute myocardial ischemia. N Engl J Med 1996;335:1333-41.
7. Antman EM, Tanasijevic MJ, Thompson B, et al. Cardiac-specific troponin I levels to predict the risk of mortality in patients with acute coronary syndromes. N Engl J Med 1996;335:1342-9.

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# Boditech Tn-I Plus Control

## INTENDED USE

Boditech Tn-I Plus Control is intended for *in vitro* diagnostic use in the quality control of Tn-I Assay Kit. **For *in vitro* diagnostic use only.**

## SUMMARY AND PRINCIPLE

The use of Boditech Tn-I Plus Control may be considered as an objective assessment of the precision of Tn-I Assay Kits and is an integral part of Good Laboratory Practices. Boditech Tn-I Plus Control is provided in lyophilized form.

## COMPONENTS

Boditech Tn-I Plus Control consists of 'Boditech Tn-I Plus Control level 1', 'Boditech Tn-I Plus Control level 2', 'Instruction for Use' and 'Control value & Barcode Sheet'.

- The control contains Human Cardiac Troponin I/C Complex Antigen and horse serum

- Each control vial packed in a box. 4 p.

## SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech Tn-I Plus Control should not be used past the expiration date.
- Boditech Tn-I Plus Control is solely designed to be provided instrument-specific calibration curves of Boditech Readers and Tn-I Assay Kits.
- Human source materials from which Boditech Tn-I Plus Control is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

## STORAGE AND STABILITY

- Storage and stability condition of Boditech Tn-I Plus Control.

Unopened	
Temperature	Expiration date
+2 to +8 °C	Until expiration date on the label.

Opened (After reconstitution)	
Temperature	Expiration date
+2 to +8 °C	1 day
-20 to -80 °C	7 days

- Close the opened Boditech Tn-I Plus Control bottle tightly after use.
- Once the Boditech Tn-I Plus Control was frozen, it should be used ONE TIME ONLY for test, because repeated freezing and thawing can result in the change of test values.
- After use, any residual product should NOT BE RETURNED to the original vial.

- Bacterial contamination of reconstituted Boditech Tn-I Plus Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

## INSTRUCTIONS FOR USE

Boditech Tn-I Plus Control is supplied in lyophilized form.

1. Carefully reconstitute each vial of lyophilized with exactly 1 mL of sterilized distilled water.
2. Close the bottle and allow to stand for 30 minutes before use. Ensure contents are completely dissolved by swirling gently. Avoid formation of foam. Do not shake.

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For Technical Assistance  
**Boditech Med Inc.'s Technical Services at**

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# Boditech Tn-I Plus Calibrator

## INTENDED USE

Boditech Tn-I Plus Calibrator is intended for *in vitro* diagnostic use in the calibration of Tn-I Assay Kit.

**For *in vitro* diagnostic use only.**

## INTRODUCTION

The use of Boditech Tn-I Plus Calibrator may be considered as an objective assessment of the precision of Tn-I Assay Kits and is an integral part of Good Laboratory Practices. Boditech Tn-I Plus Calibrator is provided in lyophilized form.

## COMPONENTS

Boditech Tn-I Plus Calibrator consists of 'Boditech Tn-I Plus Calibrator level 1', 'Boditech Tn-I Plus Calibrator level 2', 'Instruction for Use' and 'Control value & Barcode Sheet'.

- The calibrator contains Human Cardiac Troponin I/C Complex Antigen in Human Plasma.
- Each calibrator vial packed in a box.

## SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech Tn-I Plus Calibrator should not be used past the expiration date.
- Boditech Tn-I Plus Calibrator is solely designed to be provided instrument-specific calibration curves of Boditech Readers and Tn-I Assay Kits.
- When a certain parameter of the instrument for AFIAS test is calibrated from Cal/QC mode, the system displays fluorescence intensity of each test instead of its actual result.
- Human source materials from which Boditech Tn-I Plus Calibrator is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HBsAg) and Hepatitis C Virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

## STORAGE AND STABILITY

- Storage and stability condition of Boditech Tn-I Plus Calibrator.

Temperature	Unopened	Opened (After reconstitution)	
	2 ~ 8 °C	2 ~ 8 °C	-20 ~ -80 °C
Expiration date	Until expiration date on the label.	1 day	7 days

- Close the opened Boditech Tn-I Plus Calibrator bottle tightly after use.
- Once the Boditech Tn-I Plus Calibrator was frozen, it should be used ONE TIME ONLY for test, because repeated freezing and thawing can result in the change of test values.
- After use, any residual product should NOT BE RETURNED to the original vial.
- Bacterial contamination of reconstituted Boditech Tn-I Plus Calibrator will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

## INSTRUCTIONS FOR USE

Boditech Tn-I Plus Calibrator is supplied in lyophilized form.

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2.3.2.4 p.

### INTENDED USE

**AFIAS Rota/Adeno** is a fluorescence Immunoassay (FIA) for the qualitative determination of both rotavirus and adenovirus in human feces. It is useful as an aid in management and monitoring of viral gastroenteritis.

For *in vitro* diagnostic use only.

### INTRODUCTION

The major symptoms of viral enteritis are diarrhea and vomiting. Viral enteritis is infectious disease caused by many viruses such as Rotavirus, Adenovirus.

Rotaviruses (RVs) are the main etiologic agents of serious diarrheal disease in infants and young children under while 2 years of age throughout the world. Group A RVs are the major cause of human infections. Outbreaks with a strict seasonal winter pattern occur in tropic climates. Infections are spread more evenly throughout the year. After a short incubation period of 24–48 h, the onset of illness is sudden, with watery diarrhea, vomiting, and rapid dehydration. Untreated RV infection is a major cause of infantile death in developing countries. <sup>1)</sup>

Adenoviruses may cause epidemics, endemics and sporadic infections in all geographical regions of the world. They do not show seasonal outbreak and can be seen throughout entire year. Adenovirus infection cause variety symptom in several region. Especially, the adenovirus type 40 and 41 cause acute gastroenteritis primarily in children like as Rotavirus. <sup>2)</sup>

Also, it is known that Rotavirus-Adenovirus co-infection rate is above 5% in gastroenteritis patients. Diagnosing of the viral disease is important in reducing unnecessary use of antibiotics. <sup>2,3)</sup>

**AFIAS Rota/Adeno** is an immunoassay for the detection of both Rotavirus and Adenovirus in human stool sample.

### PRINCIPLE

The test uses a sandwich immunodetection method; the detector antibodies in buffer bind to antigens in the sample, forming antigen-antibody complexes and migrate onto nitrocellulose matrix to be captured by the other immobilized- antibodies on the test line.

More antigens in the sample will form more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by instrument for AFIAS test to display 'rotavirus positive' and 'adenovirus positive' in the sample.

### COMPONENTS

**AFIAS Rota/Adeno** consists of 'cartridge', 'pipette tip', 'sample collection tube', 'ID chip' and 'Instruction for use'.

- Each cartridge is packaged in an aluminum pouch, has three components; a cartridge part, detector part and diluent part.
- The cartridge part contains the membrane called a test strip which has rotavirus antibody at the each test line, adenovirus antibody at the each test line and Chicken IgY at the control line.
- The detector part contains rotavirus & adenovirus antibody fluorescence conjugate, anti-chicken IgY fluorescence conjugate, bovine serum albumin (BSA) as a stabilizer and sodium azide as a preservative in phosphate buffered saline (PBS).

- The diluent part contains detergent, bovine serum albumin (BSA), sodium azide as a preservative in phosphate buffered saline (PBS).
- The sample collection tube contains bovine serum albumin (BSA), detergent, sodium chloride and sodium azide in Tris-HCl and DDW.
- The sample dilution buffer is pre-dispensed in a sample collection tube. 25 sample collection tubes are packaged in a box and enclosed in a Styrofoam box with an ice-pack for the shipment.

### WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this 'Instruction for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (cartridge, ID chip and sample collection tube) should agree.
- Do not interchange test components between different lots or use test components after the expiration date, either of which might yield incorrect of test result(s).
- Do not reuse cartridge or sample collection tube. A sample collection tube should be used for processing one sample only. A cartridge should be used for testing one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use the cartridge if pouch is damaged or has already been opened.
- For shipping, samples must be packed in accordance with local regulations.
- Allow cartridge and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used pipette tips, sample collection tubes, and cartridges should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
- AFIAS Rota/Adeno will provide accurate and reliable results subject to the below conditions.
  - **AFIAS Rota/Adeno** should be used only in conjunction with the instrument for AFIAS test.

### LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause false negative result as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

### STORAGE AND STABILITY

Component	Storage condition		Note
	Storage Temperature	Shelf life	
Cartridge	2 - 8 °C	20 months	Unopened
Sample collection tube	2 - 8 °C	20 months	Unopened

### MATERIALS SUPPLIED

REF SMFP-43

Components of **AFIAS Rota/Adeno**

- Cartridge Box Contains
  - Cartridge 24
  - Pipette Tip (Zipper bag) 24
  - Spare Cartridge Zipper bag 1
  - ID Chip 1
  - Instruction for Use 1
- Sample Collection Tube Box Contains
  - Sample collection tube 25

### MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from **AFIAS Rota/Adeno**.

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
### SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS Rota/Adeno** is human feces.

- Fecal samples must be taken as soon as the symptoms (diarrhea) appear.
- Collect random samples of feces in a clean, dry container or a receptacle it, making sure to exclude large solid lumps.


- Fill up the groove of a sampling stick with fecal samples and please check whether the quantity is too much or not.
- It is recommended to test the sample after collection for immediately.
- Samples may be stored for up to a week at 2-8 °C prior to being tested.
- If testing will be delayed more than a week, samples should be frozen at -20 °C or below. Samples stored frozen at -20 °C or below for 2 months showed no performance difference.
- Once the sample was frozen, it should be used one time only for test, because repeated freezing and thawing can result in erroneous results.

### TEST SETUP

- Check the components of the **AFIAS Rota/Adeno** as described below: cartridge, pipette tip, sample collection tube, ID chip and instruction for use
- Keep the sealed cartridge (if stored in the refrigerator) at room temperature for at least 30 minutes just prior to the test. Place the cartridge on a clean, dust-free and flat surface.
- Turn on the instrument for AFIAS test.

- Empty the tip box.
- Insert the ID chip into the "ID chip port".
- Please refer to the instrument for AFIAS test 'Operation Manual' for complete information and operating instructions.

### TEST PROCEDURE

- Collect random samples of feces in a clean, dry container or a receptacle, making sure to exclude large solid lumps.
- Invert a sample collection tube and loosen the cap which is attached a sampling stick (yellow color).
- Introduce the sampling stick into the fecal sample about 5~6 times at different sites. In order to get sampling even in the spirals of the stick and to ensure appropriate specimen to buffer ratio, try to avoid obtaining clumps of fecal matter.
- Return the stick to the sample collection tube. Tighten the cap thoroughly and shake the tube vigorously so as to disperse the specimen throughout the buffer in the tube.
- Break off the black tip on the outside of the black cap.
- Discard 3 drops of reagent onto the paper towel before applying to the cartridge.
- Select 'General Mode' in the instrument for AFIAS test.
- Load 3~4 drops (about 50~100 µL) of the mixture onto the sample well of the cartridge.
- Insert the cartridge into the cartridge holder.
- Insert a tip into the tip hole of the cartridge.
- Input 'Patient ID' and select the 'Sample type' into the instrument for AFIAS test.
- Tap the 'START' icon on the screen.
- Tap the 'Confirm start' icon on the screen after reconfirmation of 'Patient ID' and 'Sample type'.
- The test result will be displayed on the screen after 12 minutes.

※ Note: Refer to the instrument for AFIAS test Operation Manual to select a sample type.

### INTERPRETATION OF TEST RESULT

- The instrument for AFIAS test calculates the test result automatically and displays 'Positive'/'Negative'/'Indeterminate'.
- Ancillary value is served in the form of a cut-off index (COI) value.

- Rotavirus line		
Cut-off Index (COI)	Result	Note
≤ 0.90	Negative	No need for additional test.
> 0.90, < 1.0	Indeterminate	Need to retest.
≥ 1.0	Positive	Need confirmation test.

- Adenovirus line		
Cut-off Index (COI)	Result	Note
≤ 0.90	Negative	No need for additional test.
> 0.90, < 1.0	Indeterminate	Need to retest.
≥ 1.0	Positive	Need confirmation test.

**QUALITY CONTROL**

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are not provided with AFIAS Rota/Adeno. For more information regarding obtaining the control materials, contact Boditech Med Inc.'s Sales Division for assistance. (Please refer to the instruction for use of control material.)

**PERFORMANCE CHARACTERISTICS**

- Analytical sensitivity**
  - Cut-Off Index (COI) ≥ 1.0
- Hook effect**
  - AFIAS Rota/Adeno was not observed hook effect for these virus types..

Virus type	Concentration
#1 Adenovirus type 40	2.5 X 10 <sup>6</sup> PFU/mL
#2 Adenovirus type 41	1.6 X 10 <sup>6</sup> PFU/mL
#3 Rotavirus WA strain	1.2 X 10 <sup>6</sup> PFU/mL
#4 Rotavirus DS-1 strain	1.4 X 10 <sup>6</sup> PFU/mL

- Analytical inclusivity**
  - Rotavirus
    - Two Rotavirus genotyped samples (genotype 40,41) were tested with the AFIAS Rota. All of genotyped samples were shown positive.

Virus type	Concentration
#1 Rotavirus WA strain	1.2 X 10 <sup>6</sup> PFU/mL
#2 Rotavirus DS-1 strain	1.4 X 10 <sup>6</sup> PFU/mL

- Adenovirus
  - Two Adenovirus type samples (40 & 41) were tested with the AFIAS Rota/Adeno. The Enteric Adenovirus (40&41) samples were shown positive.

Virus type	Concentration
#1 Adenovirus type 40	2.5 X 10 <sup>6</sup> PFU/mL
#2 Adenovirus type 41	1.6 X 10 <sup>6</sup> PFU/mL

- Analytical specificity**
  - Cross-reactivity**
    - There was no false positive result from 19 species virus samples and 24 species bacteria samples containing potentially cross-reactive substances with the AFIAS Rota/Adeno test.

Bacteria	OD280	Conc.
#1 Staphylococcus aureus (ATCC 29213)	1.35	≥10 <sup>6</sup> CPE/ml
#2 Enterococcus faecalis(ATCC 29212)	0.94	
#3 Escherichia coli (ATCC 25922)	0.76	
#4 Kleb-siella oxytoca (ATCC 700432)	0.66	
#5 Pseudomonas aeruginosa(ATCC 27853)	0.72	
#6 Neisseria gonorrhoeae (ATCC 27853)	0.82	
#7 Aeromonas hydrophila (KCCM 32586)	0.73	
#8 Enterobacter cloacae (KCCM 32586)	0.92	
#9 Vibrio parahaemolyticus (KCCM11965)	0.56	
#10 Salmonella group B (Clinical isolate from patient)	0.81	
#11 Salmonella group C (Clinical isolate from patient)	0.89	
#12 Salmonella group D (Clinical isolate from patient)	0.63	
#13 Salmonella group E (Clinical isolate from patient)	0.72	
#14 Shigella group D	0.69	

(Clinical isolate from patient)	
#15 Staphylococcus epidermidis (Clinical isolate from patient)	0.65
#16 Serratia marcescens (Clinical isolate from patient)	0.72
#17 Yersinia enterocolitica (Clinical isolate from patient)	0.91
#18 Salmonella typhi (Clinical isolate from patient)	1.02
#19 Clostridium difficile (Clinical isolate from patient)	1.63
#20 Candida albicans (Clinical isolate from patient)	0.64
#21 Candida parapsilosis (Clinical isolate from patient)	0.78
#22 Campylobacter spp	0.99
#23 Proteus vulgaris	0.67
#24 Proteus mirabilis	0.83

Virus	Conc. PFU/mL
#1 ADV 1	1.7X10 <sup>6</sup> PFU/mL
#2 ADV2 (KBPV-VR-58)	8X10 <sup>6</sup> PFU/mL
#3 ADV3 (KBPV-VR-2)	1.43X10 <sup>6</sup> PFU/mL
#4 ADV4 (KBPV-VR-60)	4X10 <sup>6</sup> PFU/mL
#5 ADV5	7.9X10 <sup>7</sup> PFU/mL
#6 ADV31 (ATCC VR-1109)	3.2X10 <sup>6</sup> PFU/mL
#7 ADV40 (ATCC VR-931)	2.5X10 <sup>6</sup> PFU/mL
#8 ADV41 (ATCC VR-930)	1.6X10 <sup>6</sup> PFU/mL
#9 enterovirus type 71 (ATCC VR-784)	3.4X10 <sup>6</sup> PFU/mL
#10 Cytomegalovirus (ATCC-VR-538)	2.0X10 <sup>6</sup> PFU/mL
#11 poliovirus type 1 (ATCC VR-58)	3.5X10 <sup>6</sup> PFU/mL
#12 coxsackievirus B type 5 (ATCC VR-1036)	1.7X10 <sup>6</sup> PFU/mL
#13 coxsackievirus B type 6 (ATCC VR-1037)	1.3X10 <sup>6</sup> PFU/mL
#14 herpes simplex virus type 1 (ATCC-VR-733)	4.4X10 <sup>7</sup> PFU/mL
#15 herpes simplex virus type2	8.0X10 <sup>5</sup> PFU/mL
#16 Norovirus VLP GI	2.74mg/ml
#17 Norovirus VLP GII	3.48mg/ml
#18 Rotavirus WA strain	1.2X10 <sup>6</sup> PFU/mL
#19 Rotavirus DS-1 strain	1.4X10 <sup>6</sup> PFU/mL

- Interference**
  - Interference materials such as below the ones in the table were added to the test samples at concentrations much higher than their normal physiological levels in human feces. AFIAS Rota/Adeno test results did not show any significant interference with these biomolecules and chemical drugs.

Biomolecules	Conc.
#1 Bilirubin	500 μM
#2 Hemoglobin	2 g/L
#3 Triglycerides	1.5 mg/mL
#4 Cholesterol	20 mM
#5 BSA	60 g/L
Chemicals	Conc.
#1 cephradine	0.25 mg/mL
#2 cefuroxime axetil	
#3 Cefpodoxime proxetil	
#4 cefixime	
#5 tetracycline hcl	
#6 levofloxacin	
#7 amoxicillin	
#8 loperamide oxide	
#9 metronidazole	
#10 ibuprofen	
#11 acetaminophen	
#12 barium sulfate	
#13 DMSO	
#14 DMF	

#15	DDW
#16	PBS

- Precision**
  - Between lots
    - One person tested three different lots of AFIAS Rota/Adeno, ten times at each concentration of the control standard.
  - Between persons
    - Three different persons tested one lot of AFIAS Rota/Adeno, Five times at each concentration of the control standard.

Sample	Cal.	Between lots		Between persons	
		Positive/ Sample No.	Positive	Positive/ Sample No.	Positive
Rotavirus	Negative	0/10	0 %	0/5	0 %
	Low	10/10	100 %	5/5	100 %
	Mid	10/10	100 %	5/5	100 %
Adenovirus	Negative	0/10	0 %	0/5	0 %
	Low	10/10	100 %	5/5	100 %
	Mid	10/10	100 %	5/5	100 %

- Between days
  - One person tested one lot of AFIAS Rota/Adeno during three days, Five times at each concentration of the control standard.
- Between sites
  - One person tested AFIAS Rota/Adeno at three different site, five times at each concentration of the control standard.

Sample	Cal.	Between days		Between sites	
		Positive/ Sample No.	Positive	Positive/ Sample No.	Positive
Rotavirus	Negative	0/5	0 %	0/5	0 %
	Low	5/5	100 %	5/5	100 %
	Mid	5/5	100 %	5/5	100 %
Adenovirus	Negative	0/5	0 %	0/5	0 %
	Low	5/5	100 %	5/5	100 %
	Mid	5/5	100 %	5/5	100 %

Rotavirus	Result	95 % CI	Sample no.	
			Positive	Negative
Clinical sensitivity	99.00%	94.6-100%	100	100
Clinical specificity	99.00%	94.6-100%		
Adenovirus	Result	95 % CI	Sample no.	
			Positive	Negative
Clinical sensitivity	100%	96.2-100%	95	100
Clinical specificity	97.00%	91.5-99.4%		

AFIAS Rota/Adeno	Reference analyzer assay result		
	Rotavirus		
	Positive	Negative	Total
Positive	98	2	100
Negative	0	100	100
Total	98	102	200
Percent positive agreement (%)		100.0	
Percent negative agreement (%)		98.0	
Overall percent agreement (%)		99.0	
AFIAS Rota/Adeno	Reference analyzer assay result		
	Adenovirus		
	Positive	Negative	Total
Positive	91	7	98
Negative	0	97	97
Total	91	104	195
Percent positive agreement (%)		100.0	
Percent negative agreement (%)		93.3	
Overall percent agreement (%)		96.4	

**REFERENCES**

- Rotavirus Methods and Protocols. James Gray et al., Methods in Molecular Medicine., 2000, 6-7 pp.
- Diarrheagenic pathogens in polymicrobial infections, Brianna Lindsay et al., Emerging infectious disease, 2011, 17:4 606-611
- Rotavirus and adenovirus frequency among patients with acute gastroenteritis and their relationship to clinical parameters: a retrospective study in Turkey, Asia Pacific Family Medicine, 2009, 8:8 1-8.

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For technical assistance; please contact:  
**Boditech Med Inc.'s Technical Services**  
 Tel: +(82) 33 243-1400  
 E-mail: sales@boditech.co.kr

**Boditech Med Incorporated**  
 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398 Republic of Korea  
 Tel: +(82) -33-243-1400  
 Fax: +(82) -33-243-9373  
 www.boditech.co.kr

**Obelis s.a**  
 Bd. Général Wahis 53, 1030 Brussels, BELGIUM  
 Tel: +(32) -2-732-59-54  
 Fax: +(32) -2-732-60-03  
 E-Mail: mail@obelis.net

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# Boditech Rota/Adeno Control

## INTENDED USE

Boditech Rota/Adeno Control is intended for *in vitro* diagnostic use in the quality control of Rota/Adenovirus Assay Kit.  
**For *in vitro* diagnostic use only.**

## INTRODUCTION

The use of Boditech Rota/Adeno Control may be considered as an objective assessment of the precision of Rota/Adenovirus Assay Kits and is an integral part of Good Laboratory Practices. Boditech Rota/Adeno Control is provided in lyophilized form.

## COMPONENTS

Boditech Rota/Adeno Control consists of 'Boditech Rota/Adeno Control level 1', 'Boditech Rota/Adeno Control level 2', 'Instruction for Use' and 'Barcode Sheet'.

- The level 1 control contains bovine serum albumin (BSA), Tris, Tween 20 and sodium azide as a preservative in DDW.
- The level 2 control contains Rotavirus WA strain inactivated virus, Adeno 41 type inactivated virus bovine serum albumin (BSA), Tris, Tween20 and sodium azide as a preservative in DDW.

## SAFETY PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech Rota/Adeno Control should not be used past the expiration date.
- Boditech Rota/Adeno Control is solely designed to be provided instrument-specific calibration curves of Boditech Readers and Rota/Adenovirus Assay Kits.
- Human source materials from which Boditech Rota/Adeno Control is derived were tested at a donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody and Hepatitis B Surface Antigen (HBsAg) and Hepatitis C virus (HCV) antibody, and found to be NON-REACTIVE. FDA-approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

## STORAGE AND STABILITY

- Storage and stability condition of Boditech Rota/Adeno Control.

	Unopened	Opened (After reconstitution)
Temperature	2 ~ 8 °C	2 ~ 8 °C
Expiration date	Until expiration date on the label.	14 days

- Close the opened Boditech Rota/Adeno Control bottle tightly after use.
- After use, any residual product should NOT BE RETURNED to the original vial.
- Bacterial contamination of reconstituted Boditech Rota/Adeno Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

## INSTRUCTIONS FOR USE

Boditech Rota/Adeno Control is supplied in lyophilized form.

1. Carefully reconstitute each vial of lyophilized with exactly 1 mL of sterilized distilled water.
2. Close the bottle and allow to stand for 30 minutes before use. Ensure contents are completely dissolved by swirling gently.  
 Avoid formation of foam. Do not shake.

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In the event of damage to the package, contact the **Boditech Med Inc.'s Technical Services**.

## MATERIALS SUPPLIED

REF **CFPO-164**

Boditech Rota/Adeno Control Box (2 vials) 4 p.

Boditech Rota/Adeno Control level 1 (1 mL)	1
Boditech Rota/Adeno Control level 2 (1 mL)	1
Instruction for Use	1
Barcode sheet	1

## CONTROL VALUE

Boditech Rota/Adeno Control composed of negative control and positive control.

Control Level	Result
Level 1	Negative
Level 2	Positive

## INTERPRETATION OF THE RESULT

The testing result of the 'Boditech Rota/Adeno Control' should be in agreement with expected result. If the test results fall outside the expected result, repeat the test.

Indication for the wrong test result

- Errors in a manner that testing is performed
- Use of too cold or too warm Boditech Rota/Adeno control.
- Use of expired or contaminated Boditech Rota/Adeno Control.
- Errors in AFIAS Rota/Adeno or ichroma™ Rota/Adeno.
- Errors of AFIAS reader (AFIAS-1, AFIAS-6) or ichroma™ II.

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2.3.25 p.

## INTENDED USE

**AFIAS COVID-19/Flu Ag Combo** is a fluorescence Immunoassay (FIA) for simultaneous detection of Influenza A virus, Influenza B virus, and/or SARS-CoV-2 infection or respective detection of each Influenza A virus, Influenza B virus or SARS-CoV-2 in human nasopharyngeal swab. It is helpful as an aid in the distinguish novel corona virus and influenza infection in flu-like symptom.

For *in vitro* diagnostic use only.

## INTRODUCTION

The third zoonotic human coronavirus (CoV) of the century emerged in December 2019. This virus, the newly identified coronavirus 2019-nCoV, could cause risky pneumonia so that prevention and control of the infection has become highly required. The 2019-nCoV is a member of the Betacoronavirus Genus, that also includes Severe Acute Respiratory Syndrome coronavirus (SARS-CoV) and Middle East Respiratory Syndrome coronavirus (MERS-CoV). Since it is identified that symptoms become rapidly severe without a proper treatment after onset of illness, early diagnosis of the virus infection is quite crucial. Currently, the spread of the viral transmission become fast so that the prevention of local transmission requires a point-of care test (POCT).

Influenza, or flu, known as 'febrile respiratory illness' can cause mild to severe symptoms, such as a high fever, chills, headache, muscle pains coughing and even death. These illness typically begins after exposure to the influenza virus in the respiratory epithelial cell from person-to-person through sneezing, coughing, or touching contaminated surfaces. Within 48 hours after the onset of symptoms, the patient is strongly recommended to visit the nearest medical facility for the diagnosis of Influenza A or B and to take the antiviral medication. The preventive measure is highly required for those at increased risk for severe illness, so that early and differential diagnosis between influenza types A or B is very essential.

**AFIAS COVID-19/Flu Ag Combo** is an *in vitro* diagnostic medical device that helps you to diagnose and the distinguish novel corona virus and influenza infection by detecting the specific antigen of SARS-CoV-2, influenza A, and influenza B.

## PRINCIPLE

AFIAS-1/6/10 is automated device that uses microfluidic to do dispensing and sample management. Once sample is loaded into sample well, automated microfluidic liquid handler transfers sample and process the subsequent immunoassay (FIA).

The principle of immunoassay uses a sandwich immunodetection method; the detector antibodies in buffer bind to antigens in the sample, forming antigen-antibody complexes, and migrates onto nitrocellulose matrix to be captured by the other immobilized antibodies on test strip.

More antigens in the sample will form the more antigen-antibody complexes which lead to stronger fluorescence signal by detector antibodies, which is processed by instrument for AFIAS tests to show concentration of SARS-CoV-2, influenza A, influenza B antigens in sample respectively. This signal then is interpreted by the reader to display the 'Positive' / 'Negative' in the sample.

## COMPONENTS

**AFIAS COVID-19/Flu Ag Combo** consists of 'Cartridges'.

- Each cartridge packaged in an aluminum pouch has two components, a detector part and cartridge part.
- Cartridge part contains the membrane called a test strip which has mouse monoclonal anti-nCoV, anti-influenza A, and anti-influenza B antibodies at the test line, and chicken IgY at the control line.
- Detector part contains the mouse monoclonal anti-nCoV, anti-influenza A, and anti-influenza B fluorescence conjugate, anti-chicken IgY-fluorescence conjugate as a lyophilized granule type, bovine serum albumin (BSA) as a stabilizer and sodium azide as a preservative.
- The extraction tube contains sodium chloride, sodium azide as a preservative in Tris-HCl buffer.

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow instructions and procedures described in this 'Instruction for use'.
- Do not reuse cartridge and extraction set.
- Lot numbers of all the test components (Cartridge, ID chip, Tip) must match each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield incorrect test result(s).
- The cartridge should remain sealed in its aluminum pouch until just before use. Do not use the cartridge if pouch is damaged or has already been opened.
- Allow the cartridge and sample to be at room temperature for approximately 30 minutes.
- The instrument for AFIAS tests may generate slight vibration during use.
- Used cartridges, extraction tubes, nozzles, pipette tips and swabs should be handled carefully and disposed by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.

## WARNINGS AND PRECAUTIONS FOR SAMPLE

- Use the fresh samples.
- It is recommended to test the sample immediately after sample collection.
- Refrain from smoking or eating, while sample is collected.
- Do not collect samples outside of the nasopharynx / oropharynx. In any cases, pre-education for user is required for the proper sample collection.
- Please use fresh swab to avoid the cross-reactivity between samples. Never reuse the sterile swab.
- The improper samples such as those from an individual who has recently taken any interfering medicine or samples mistakenly mixed up with different patients shall cause inaccurate test results.

## LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result(s) due to the non-responsiveness of the antigen to the antibodies which is most common if the epitope is masked by some unknown components, so therefore not being able to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may also cause false negative as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the

- test samples.
- If the test result is "Negative" even though the patient has significant infectious symptoms, it should be recommended to conduct additional test including PCR or culture test.
- The accurate determination of test result as "Positive" should be confirmed by additional clinical evaluation.
- "Negative" result should be considered with possibilities of other infections. Positive result should be considered with additional infections by another pathogenic bacterium.
- If the product has positive results, any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.
- In case of antigen concentration is low, the test may yield false negative results. Therefore, the negative results cannot exclude the possibility of infection completely.
- This product is only to detect the presence of a SARS-CoV-2, influenza A, and influenza B antigen.

## STORAGE AND STABILITY

Storage condition			
Component	Storage Temperature	Shelf life	Note
Cartridge	4 - 30 °C	20 months	Unopened
		1 month	Resealed
Extraction buffer	4 - 30 °C	20 months	Disposable

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## MATERIALS SUPPLIED

REF	SMFP-77	Components of AFIAS COVID-19/Flu Ag Combo	
		Cartridge Box Contains	
		- Cartridge	24
		- Pipette tip (Zipper bag)	24
		- Extraction set	
		Extraction buffer	24
		Nozzle	24
		- Spare Cartridge Zipperbag	1
		- ID chip	1
		- Instruction for use	1

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Following items can be purchased separately from **AFIAS COVID-19/Flu Ag Combo**.

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## SAMPLE COLLECTION AND PROCESSING

The sample type for **AFIAS COVID-19/Flu Ag Combo** is human nasopharyngeal.

- Collection method for sample

Nasopharyngeal swab

### Nasopharyngeal swab sample collection


To collect samples, insert a sterile swab in the nasal cavity and spin it smoothly in the nasopharynx.

- It is recommended to test the sample immediately after collection. If do not use sample immediately, should store at 2-8 °C.
- Samples stored at 2-8 °C for 2 days showed no performance difference.
- It is highly recommended that **AFIAS COVID-19/Flu Ag combo** test would be performed on the nasopharyngeal swab specimens directly collected from patients with provided extraction buffer.

## TEST SETUP

- Check the components of the **AFIAS COVID-19/Flu Ag Combo** as described below: Cartridge, Extraction tube set, ID chip, Spare cartridge zipper bag and instruction for use
- If the sealed cartridge and extraction buffer have been stored in a refrigerator, place them on a clean and flat surface at room temperature for at least 30 minutes before testing.
- Empty the tip box.
- Insert the ID chip into the "ID chip port". (Please refer to the instrument for AFIAS tests 'operation manual' for complete information and operating instructions.)

## TEST PROCEDURE

### [General mode: AFIAS-6, AFIAS-1]

- Open the lid of the aluminum foil of the extraction buffer tube.
- Collect samples with a sterile swab and then put it into the extraction buffer tube. Spin the sterile swab 5 times and squeeze the sterile swab to extract the sample into the buffer.
- Pull the sterile swab out of the extraction buffer tube.
- Continue squeezing and pushing the swab to the top of extraction buffer tube to pull it out of tube.
- Assemble a nozzle to the top of extraction tube.
- Squeeze the extraction buffer tube gently to completely fill the sample well.
- Select 'General Mode' in the instrument for AFIAS tests
- Insert the cartridge into the cartridge holder.
- Insert a tip into the tip hole of the cartridge.
- Tap the 'START' icon on the screen.
- The test results will be displayed on the screen after 20 minutes.

### [Emergency mode: AFIAS-10]

- The test procedure is same with "General mode: AFIAS-6/AFIAS-1". (General mode: AFIAS-6/AFIAS-1 ① - ⑥)
- Select 'Emergency Mode' in the instrument for AFIAS tests.
- Select the tip type(Normal) on the screen.
- Insert the cartridge into the cartridge holder.
- Insert a tip into the tip hole of the cartridge.
- Tap the 'START' icon on the screen.
- The test result will be displayed on the screen after 20 minutes.

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## INTERPRETATION OF TEST RESULT

- The instrument for AFIAS tests calculates the test result automatically and displays as below.

COI (Cut-off index)	Display	
	< 1	≥ 1
COVID-19	COVID-19 Negative	COVID-19 Positive
Influenza A	Flu A Negative	Flu A Positive
Influenza B	Flu B Negative	Flu B Positive

- If test result is Invalid, you need to perform a new test on a new test cartridge with a new test sample.

**QUALITY CONTROL**

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are provided on demand. For more information regarding obtaining the control materials, contact [Boditech Med Inc.'s Sales Division for assistance.](#) (Please refer to the instruction for use of control material.)

**PERFORMANCE CHARACTERISTICS**

**Analytical sensitivity**

- Cut-off**  
The cut-off value is 1 as COI (Cut off index) that is obtained from algorithm of the instrument.

LOD		
Virus type	Conc.	
COVID-19 SARS-CoV-2 (USA/WA1/2020)	0.55x10 <sup>2</sup> TCID50/mL	
Influenza A A/Brisbane/02/2018 (H1N1)pdm09-like virus	4.4 x 10 <sup>2</sup> pfu/mL	
Influenza A A/Kansas/14/2017 (H3N2)-like virus	1.5 x 10 <sup>3</sup> pfu/mL	
Influenza B B/Colorado/06/2017-like virus (B/Victoria/2/87 lineage)	2.0 x 10 <sup>4</sup> pfu/m	
Influenza B B/Phuket/3073/2013-like virus (B/Yamagata/16/88 lineage)	2.9 x 10 <sup>3</sup> pfu/mL	

**Analytical specificity**

- Cross-reactivity**  
There was no significant cross-reactivity on 30 various other viruses and 36 various bacteria with the **AFIAS COVID-19/Flu Ag Combo** test.

Virus		
1 Adenovirus type1	16	Echovirus 6
2 Adenovirus type2	17	Echovirus 9
3 Adenovirus type3	18	Enterovirus 71
4 Adenovirus type4	19	HCMV-AD-169
5 Adenovirus type6	20	HSV-1 - F(3A20)
6 Adenovirus type7	21	HSV-2 - MS(4A6)
7 Corona virus - FCV(3A2)	22	Meales virus
8 Corona virus - FIP(2A4)	23	Mumps virus
9 Coxsackievirus A2	24	Polio virus - sin(3A4)
10 Coxsackievirus A4	25	Respiratory Syncytial virus A
11 Coxsackievirus B1 - conn5	26	Rhinovirus - RV21
12 Coxsackievirus B3-nancy (5A1)	27	Rhinovirus -RV14
13 Dengue virus	28	Rhinovirus -RV71
14 Echovirus 25	29	Rubella virus
15 Echovirus 3	30	Zika virus
Bacteria		
1 <i>Candida albicans</i>	19	<i>Neisseria gonorrhoeae</i>
2 <i>Candida glabrata</i>	20	<i>Neisseria meningitidis</i>
3 <i>Candida tropicalis</i>	21	<i>Neisseria sicca</i>
4 <i>Citrobacter freundii</i>	22	<i>Proteus mirillis</i>
5 <i>Corynebacterium sp.</i>	23	<i>Proteus vulgaris</i>
6 <i>Corynebacterium diphtheriae</i>	24	<i>Pseudomonas aeruginosa</i>
7 <i>Enterococcus faecalis</i>	25	<i>Serratia marcescens</i>
8 <i>Enterococcus gallinarum</i>	26	<i>Staphylococcus aureus</i>
9 <i>Escherichia coli</i>	27	<i>Staphylococcus epidermidis</i>

10 <i>Hemophilus influenzae</i>	28	<i>Stenotrophomonas maltophilia</i>
11 <i>Hemophilus parainfluenzae</i>	29	<i>Streptococcus sp. (Group D)</i>
12 <i>Klebsiella oxytoca</i>	30	<i>Streptococcus agalactiae (Group B)</i>
13 <i>Klebsiella pneumoniae</i>	31	<i>Streptococcus anginosus (Group F)</i>
14 <i>Lactobacillus sp.</i>	32	<i>Streptococcus dysgalactiae (Group C)</i>
15 <i>Legionella spp</i>	33	<i>Streptococcus dysgalactiae (Group G)</i>
16 <i>Listeria monocytogenes</i>	34	<i>Streptococcus mutans</i>
17 <i>Moraxella catarrhalis</i>	35	<i>Streptococcus pneumoniae</i>
18 <i>Mycobacterium tuberculosis</i>	36	<i>Streptococcus pyogenes</i>

- Interference**  
There was no significant interference effect on from these substances.

Interference materials	Conc.
1 Nasal sprays drop	20%
2 Nasal corticosteroids	20%
3 Homeopathic allergy relief medicine	20%
4 Mouth wash (Listerin)	5 mg/mL
5 Throat lozenges, oral anesthetic & analgesic	5 mg/mL
6 Antiviral drugs (Tamiflu; Oseltamivir)	5 mg/mL
7 Antibiotic nasal ointment (Bactroban; mupirocin)	5 mg/mL
8 Whole blood	1%
9 Analgesic (Acetaminophen)	10 mg/mL
10 Analgesic (Ibuprofen)	10 mg/mL
11 Povidone-iodine	1%
12 Acetylsalicylic acid (Aspirin)	20 mg/mL
13 Antibacterial (cefadroxil)	5 mg/mL
14 Mucin (Porcine stomach)	0.50%
15 Throat lozenge (VICKS; cetylpyridinium chloride)	20 mg/mL
16 Throat lozenge (dipotassium glycyrrhizinate)	20 mg/mL
17 Throat lozenge (Nandina extraction)	20 mg/mL

**Precision**

- Between lots**  
One person tested three different lots of **AFIAS COVID-19/Flu Ag Combo**, ten times at each concentration of the control standard.
- Between persons**  
Three different persons tested one lot of **AFIAS COVID-19/Flu Ag Combo**, ten times at each concentration of the control standard.
- Between days**  
One person tested one lot of **AFIAS COVID-19/Flu Ag Combo** during three days, ten times at each concentration of the control standard.
- Between sites**  
One person tested **AFIAS COVID-19/Flu Ag Combo** at three different site, ten times at each concentration of the control standard.

	Conc.	Between Lot		Between person	
		Positive /No	Positive rate	Positive /No	Positive rate
COVID-19	Neg.	0/30	0%	0/30	0%
	Low pos.	30/30	100%	30/30	100%
	Mid pos.	30/30	100%	30/30	100%
Flu A	Neg.	0/30	0%	0/30	0%
	Low pos.	30/30	100%	30/30	100%
	Mid pos.	30/30	100%	30/30	100%
Flu B	Neg.	0/30	0%	0/30	0%
	Low pos.	30/30	100%	30/30	100%
	Mid pos.	30/30	100%	30/30	100%

COVID-19	Neg.	0/30	0%	0/30	0%
	Low pos.	30/30	100%	30/30	100%
	Mid pos.	30/30	100%	30/30	100%
Flu A	Neg.	0/30	0%	0/30	0%
	Low pos.	30/30	100%	30/30	100%
	Mid pos.	30/30	100%	30/30	100%
Flu B	Neg.	0/30	0%	0/30	0%
	Low pos.	30/30	100%	30/30	100%
	Mid pos.	30/30	100%	30/30	100%

**Clinical performance evaluation**

**AFIAS COVID-19/Flu Ag Combo** has demonstrated the following clinical performance results. (predicate: RT-PCR)

	COVID-19	Influenza A	Influenza B
Percent Positive Agreement	96.8 % (30/31)	96.9 % (31/32)	93.8 % (30/32)
Percent Negative Agreement	100 % (117/117)	100 % (74/74)	100 % (74/74)

**\* COVID-19 Ag**

- Percent Positive Agreement(%)=96.8% (95%CI:81.5% - 99.8%)
- Percent Negative Agreement(%)=100% (95%CI: 96.0% - 100%)

**\* Influenza A**

- Percent Positive Agreement(%)=96.9% (95%CI:84.3% - 99.5%)
- Percent Negative Agreement(%)=100% (95%CI: 95.1% - 100%)

**\* Influenza B**

- Percent Positive Agreement(%)=93.8% (95%CI:79.9% - 98.3%)
- Percent Negative Agreement(%)=100% (95%CI: 95.1% - 100%)

**REFERENCES**

- Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases Interim guidance (2020) 17 Jan, WHO.
- Wölfel et al. Virological assessment of hospitalized cases of coronavirus disease 2019 (2020) Nature. [Epub ahead of print]
- Trivedi SU et al. Development and Evaluation of a Multiplexed Immunoassay for Simultaneous Detection of Serum IgG Antibodies to Six Human Coronaviruses (2019) Sci Rep. 9: 1390
- Yongchen et al. Different longitudinal patterns of nucleic acid and serology testing results based on disease severity of COVID-19 patients (2020) Emerg Microbes Infect 20: 1
- Patric J Gavin, Richard B Thomson. Review of rapid diagnostic tests for influenza. Clinical and applied Immunology Reviews 4 (2003) 151-172
- Suzanne E. Dale, Christine Mayer, Marie C. Mayer and Marilyn A. Menegus. Analytical and clinical sensitivity of the 3M rapid detection influenza A+B assay. Journal of clinical microbiology, Nov. 2008, p. 3904-3807
- Christine C. Ginocchio, Frank Zhang, Ryhana Manji et al., Evaluation of multiple test methods for the detection of the novel 2009 influenza A (H1N10 during the New York City outbreak. Journal of clinical virology 45 (2009) 191-195
- Chang Kye Lee, Chi Hyun Cho, et al., Evaluation of Sofia fluorescent immunoassay analyzer for influenza A/B virus. Journal of Clinical Virology 55 (2012) 239-243
- Michael A. Di Maio, Malaya K. Sahoo, Jesse Waggoner, Benjamin A. Pinsky. Comparison of Xpert Flu rapid nucleic acid testing with rapid antigen testing for the diagnosis of influenza A and B. Journal of virological Methods 186 (2012) 137-140
- Gary P. Leonardi, Adele M Wilson, Alejandro R. Zuretti. Comparison of conventional lateral-flow assays and a new fluorescent immunoassay to detect influenza viruses. Journal of virological methods 189 (2013) 379-382
- Comparison of SD BIOLINE rapid influenza antigen test using two different specimens Nasopharyngeal swabs and nasopharyngeal aspirates. Korean J Clin Microbio. Vol. 13(4):147-150(2010)

**Note:** Please refer to the table below to identify various

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For technical assistance; please contact:

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#### LITERATŪRA

1. Laboratory testing for 2019 novel coronavirus (2019-nCoV) in suspected human cases Interim guidance (2020) 17 Jan, WHO.
2. Wölfel et al. Virological assessment of hospitalized cases of coronavirus disease 2019 (2020) Nature. [Epub ahead of print]
3. Trivedi SU et al. Development and Evaluation of a Multiplexed Immunoassay for Simultaneous Detection of Serum IgG Antibodies to Six Human Coronaviruses (2019) Sci Rep. 9: 1390
4. Yongchen et al. Different longitudinal patterns of nucleic acid and serology testing results based on disease severity of COVID-19 patients (2020) Emerg Microbes Infect 20: 1
5. Patric J Gavin, Richard B Thomson. Review of rapid diagnostic tests for influenza. Clinical and applied Immunology Reviews 4 (2003) 151-172

6. Suzanne E. Dale, Christine Mayer, Marie C. Mayer and Marilyn A. Menegus. Analytical and clinical sensitivity of the 3M rapid detection influenza A+B assay. Journal of clinical microbiology, Nov. 2008, p. 3904-3807
7. Christine C. Ginocchio, Frank Zhang, Ryhana Manji et al., Evaluation of multiple test methods for the detection of the novel 2009 influenza A (H1N1) during the New York City outbreak. Journal of clinical virology 45 (2009) 191-195
8. Chang Kye Lee, Chi Hyun Cho, et al., Evaluation of Sofia fluorescent immunoassay analyzer for influenza A/B virus. Journal of Clinical Virology 55 (2012) 239-243
9. Michael A. Di Maio, Malaya K. Sahoo, Jesse Waggoner, Benjamin A. Pinsky. Comparison of Xpert Flu rapid nucleic acid testing with rapid antigen testing for the diagnosis of influenza A and B. Journal of virological Methods 186 (2012) 137-140
10. Gary P. leonardi, Adele M wilson, Alejandro R. zuretti. Comparison of conventional lateral-flow assays and a new fluorescent immunoassay to detect influenza viruses. Journal of virological methods 189 (2013) 379-382
11. Comparison of SD BIOLINE rapid influenza antigen test using two different specimens Nasopharyngeal swabs and nasopharyngeal aspirates. Korean J Clin Microbio. Vol. 13(4):147-150(2010)

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# Boditech COVID-19/Flu Ag Control

## INTENDED USE

Boditech COVID-19/Flu Ag Control is intended for in vitro diagnostic use in the quality control of SARS-CoV-2/ influenza A/ influenza B antigen assay kit. **For in vitro diagnostic use only.**

## INTRODUCTION

The use of Boditech COVID-19/Flu Ag Control may be considered as an objective assessment of the precision of SARS-CoV-2/influenza A/influenza B antigen assay kits and is an integral part of Good Laboratory Practices. Boditech COVID-19/Flu Ag Control is provided in lyophilized form.

## COMPONENTS

Boditech COVID-19/Flu Ag Control consists of 'Boditech COVID-19/Flu Ag Control level 1', 'Boditech COVID-19/Flu Ag Control level 2', 'Instruction for Use' and 'Barcode Sheet'.

- The level 1 control contains bovine serum albumin (BSA) and sodium azide as a preservative in PBS.
- The level 2 control contains recombinant SARS-CoV-2 NP protein, recombinant influenza A NP protein, recombinant influenza B NP protein, bovine serum albumin (BSA), and sodium azide as a preservative in PBS.
- Each control vial packed in a box.

## SAFETY PRECAUTIONS AND WARNINGS

- For in vitro diagnostic use only.
- Do not pipette by mouth.
- Exercise proper precautions that would be normally required for handling laboratory reagents.
- Boditech COVID-19/Flu Ag Control should not be used past the expiration date.
- Boditech COVID-19/Flu Ag Control is not derived from human-derived substances. However, since no method can offer complete assurance as to the absence of infectious agents, these human source materials and all patient samples should be handled as though capable of transmitting infectious diseases and should be disposed of as hazardous wastes.

## STORAGE AND STABILITY

- Storage and stability condition of Boditech COVID-19/Flu Ag Control.

	Unopened	Opened
Temperature	2 ~ 8 °C	2 ~ 8 °C
Expiration date	Until expiration date on the label.	1 day

- Close the opened Boditech COVID-19/Flu Ag Control bottle tightly after use.
- After use, any residual product should NOT BE RETURNED to the original vial.
- Bacterial contamination of reconstituted Boditech COVID-19/Flu Ag Control will cause reductions in the stability of many components. If bacterial contamination is suspected, the vial should be discarded and a fresh vial needs to be reconstituted.

## INSTRUCTIONS FOR USE

Boditech COVID-19/Flu Ag Control is supplied in lyophilized form.

1. Carefully reconstitute each vial of lyophilized with exactly 0.5 mL of sterilized distilled water.
2. Close the bottle and allow to stand for 30 minutes before use. Ensure contents are completely dissolved by swirling gently. Avoid formation of foam. Do not shake.
3. Open the lid of the aluminum foil of the extraction buffer tube.
4. Transfer 50 µL of control using a transfer pipette to the extraction buffer tube (450 µL).
5. Mix well by pipetting 10 times.
6. Assemble a nozzle to the top of extraction buffer tube.

\*Next step is please follow to the instruction for use of ichroma™ COVID-19/Flu Ag Combo & AFIAS COVID-19/Flu Ag Combo.

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## MATERIALS SUPPLIED

**REF** CFPO-298

Boditech COVID-19/Flu Ag Control Box (2 vials)	
Boditech COVID-19/Flu Ag Control level 1 (0.5 mL)	1
Boditech COVID-19/Flu Ag Control level 2 (0.5 mL)	1
Instruction for Use	1
Barcode Sheet	1

## CONTROL VALUE

Boditech COVID-19/Flu Ag Control composed of negative control and positive control.

Boditech COVID-19/Flu Ag Control	Result		
	COVID-19	Flu A	Flu B
Level 1	Negative	Negative	Negative
Level 2	Positive	Positive	Positive

## INTERPRETATION OF THE RESULT

The testing result of the 'Boditech COVID-19/Flu Ag Control' should be in agreement with expected result. If the test results fall outside the expected result, repeat the test.

Indication for the wrong test result.

- Errors in a manner that testing is performed.
- Use of too cold or too warm Boditech COVID-19/Flu Ag Control.
- Use of expired or contaminated Boditech COVID-19/Flu Ag Control.
- Errors in ichroma™ COVID-19/Flu Ag Combo or AFIAS COVID-19/Flu Ag Combo.
- Errors of Boditech's readers.

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For Technical Assistance

**Boditech Med Inc.'s Technical Services at**

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
# EU DECLARATION OF CONFORMITY

- Manufacturer:** **Boditech Med Inc.**  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea  
(SRN: KR-MF-000011028)
- European Representative:** **Obelis s.a.**  
Bd. Général Wahis, 53, 1030 Brussels, Belgium  
(SRN: BE-AR-000000106)
- Product:** **AFIAS-10**  
- Catalog No.: FPRR038  
- Basic UDI-DI: 880613301031BM
- Intended use:** AFIAS-10 is an analyzer intended for use in conjunction with fluorescence immunoassay (FIA) kits for quantitative, semi-quantitative and qualitative measurements of various analytes.  
For *in vitro* diagnostic use only.
- Classification:** Class A (Rule 5)  
- According to the Annex VIII of the REGULATION (EU) 2017/746
- Conformity Assessment Route:** According to the Annex II and Annex III of the REGULATION (EU) 2017/746

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European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech CRP Control  
Cat. No. : CFPO-100

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing  
device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

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European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech CRP Calibrator  
Cat. No. : CFPO-112

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing  
device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: AFIAS HbA1c  
Cat. No. : SMFP-28

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech HbA1c Control  
Cat. No. : CFPO-96

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing  
device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

**Boditech Med Inc.** [www.boditech.co.kr](http://www.boditech.co.kr)

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, KOREA  
바디텍메드(주) 강원도 춘천시 동내면 거두단지 1길 43 Tel +82-33-243-1400 Fax +82-33-243-9373

RA-DOC-III-22 (Rev. 03)

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech HbA1c Calibrator  
Cat. No. : CFPO-108

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing  
device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: AFIAS Tn-I Plus  
Cat. No. : SMFP-35

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

**Dr. Eui Yul Choi / CEO**

**Boditech Med Inc.** [www.boditech.co.kr](http://www.boditech.co.kr)

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, KOREA  
바디텍메드(주) 강원도 춘천시 동내면 거두단지 1길 43 Tel +82-33-243-1400 Fax +82-33-243-9373

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech Tn-I Plus Control  
Cat. No. : CFPO-212

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

  
Dr. Eui Yul Choi / CEO

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech Tn-I Plus Calibrator  
Cat. No. : CFPO-213

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing  
device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

\_\_\_\_\_

# DECLARATION OF CONFORMITY

Manufacturer:	Boditech Med Inc. 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, Republic of Korea
European Representative:	OBELIS S.A Bd. Général Wahis 53, 1030 Brussels, Belgium
Product:	AFIAS Rota/Adeno Cat. No. : SMFP-43
Classification:	Others (Neither listed in the annex II of the IVDD, Non-self-testing device)
Conformity Assessment Route:	Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

**Boditech Med Inc.** [www.boditech.co.kr](http://www.boditech.co.kr)

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, KOREA  
바디텍메드(주) 강원도 춘천시 동내면 거두단지 1길 43 Tel +82-33-243-1400 Fax +82-33-243-9373

3BD-IRex.1

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech Rota/Adeno Control  
Cat. No. : CFPO-164

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing  
device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

---

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: AFIAS COVID-19/Flu Ag Combo  
Cat. No. : SMFP-77

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

—

**Boditech Med Inc.** [www.boditech.co.kr](http://www.boditech.co.kr)

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, KOREA  
바디텍메드(주) 강원도 춘천시 동내면 거두단지 1길 43 Tel +82-33-243-1400 Fax +82-33-243-9373

RA-DOC-II-174 (Rev. 01)

# DECLARATION OF CONFORMITY

Manufacturer: Boditech Med Inc.  
43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si,  
Gang-won-do, 24398, Republic of Korea

European Representative: OBELIS S.A.  
Bd. Général Wahis 53, 1030 Brussels, Belgium

Product: Boditech COVID-19/Flu Ag Control  
Cat. No. : CFP0-298

Classification: Others  
(Neither listed in the annex II of the IVDD, Non-self-testing  
device)

Conformity Assessment Route: Self-Declaration Route According to the Annex III of the IVDD

We herewith declare that the above mentioned products meet the provisions of the Council Directive 98/79/EC for In Vitro Diagnostic medical devices. All supporting documentation is retained under the premises of the manufacturer and the manufacturer is exclusively responsible for the declaration of conformity.

Standards applied: EN ISO 15223-1:2021, EN ISO 13485:2016, EN 13612:2002,  
EN ISO 23640:2015, EN 13641:2002, EN ISO 14971:2019,  
EN 13975:2003, EN ISO 17511:2021, EN ISO 18113-1:2011,  
EN ISO 18113-2:2011

Place, Date of Issue: Chuncheon, Korea, May 20, 2022

Signature:

\_\_\_\_\_

**Boditech Med Inc.** [www.boditech.co.kr](http://www.boditech.co.kr)

43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gang-won-do, 24398, KOREA  
바디텍메드(주) 강원도 춘천시 동내면 거두단지 1길 43 Tel +82-33-243-1400 Fax +82-33-243-9373

RA-DOC-III-118 (Rev. 01)

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS401
	Rev. No.	04
	Rev. Date	2019.01.17

**I. General Information**

- A. Product Name/Catalogue Number  
: AFIAS CRP / SMFP-2
- B. Recommended use of the chemical and restriction on use
  - Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

**II. Hazard Identification**

GHS classification of the product and national or regional

	Cartridge	Detection buffer
GHS classification	Not applicable	Not applicable
IMDG Code	Not applicable	Not applicable

GHS label elements including precautionary

	Cartridge	Detection buffer
Symbols	No data available	No data available
Signal word	No data available	No data available
Hazard statements	No data available	No data available

Precautionary statements

	Cartridge	Detection buffer
Prevention	No data available	No data available
Reaction	No data available	No data available
Storage	No data available	No data available
Disposal	No data available	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS401
	Rev. No.	04
	Rev. Date	2019.01.17

Other hazards

	Cartridge		Detection buffer
	Test strip	Housing	
NEPA	No data available	Health: 1 Flammability: 1 Physical hazard: 0	No data available

**III. Composition/Information on Ingredients**

[Cartridge]

	Component	Synonyms	Classification No.	Amount (%)
Test strip	Test line : Mouse monoclonal Anti-CRP	No data available	No data available	0.0000097
	Control line : Streptavidin	No data available	No data available	0.0000097
Housing	Poly(methyl methacrylate)	PMMA	9011-14-7	55 ~ 70
	Acrylonitrile butadiene styrene	ABS	9003-56-9	28 ~ 45

[Buffer]

	Component	Synonyms	Classification No.	Amount (%)
	Water	Dihydrogen Oxide	7732-18-5	99.77
	Mouse monoclonal anti-CRP	Fluorescence conjugate	No data available	0.22

**IV. First Aid Measures**

- In case of eye contact  
: Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
  - In case of skin contact  
: Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
  - If inhaled  
: Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring
  - If swallowed  
: If irritation or signs of toxicity occur, seek medical attention.
- \* The following symptoms may occur: Irreversible eye damage.  
\* Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS401
	Rev. No.	04
	Rev. Date	2019.01.17

\* This material should be considered as being potentially infectious.

**V. Fire Fighting Measures**

- Recommended extinguishing media  
: Suitable extinguishing agents including CO2, WATER SPRAY or regular form.
- Specific hazard from the chemical  
: Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters  
: Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters  
: Chemical resisted protective equipment for fire-fighter.

**VI. Accidental Release Measures**

- Personal precautions
  - \* If not danger, stop to leak
  - \* Do not contact with exposed material without protective equipment
  - \* Note the avoid conditions and materials
  - \* Do not breathe dust / fume / gas / mist / vapours / spray
- Environmental precautions  
: No data available
- Methods and materials for purification and cleaning up  
: After absorbed exposed materials, clean up with cleaner and water

**VII. Handling and Storage**

- Precaution for safe handling
  - \* Avoid skin friction.
  - \* Note the avoid conditions and materials from fire or flame.
- \* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
  - \* Note the avoid conditions and materials.
  - \* Store container tightly closed in a well-ventilated area.
- \* Storage temperature: low temperature

**VIII. Exposure controls and Personal protection**

- Engineering management  
: Good insulation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

---

Respiratory protection	Chemical resisted respiratory protection
------------------------	--

# MATERIAL SAFETY DATA SHEET

Document No.	BT-MSDS401
Rev. No.	04
Rev. Date	2019.01.17

Hand protection	Chemical resisted protective gloves
Eye protection	Chemical resisted eye protection
Skin and body protection	Chemical resisted protective clothes
- Standard of exposure	
Chemical	No data available
Biological	No data available

## IX. Physical and Chemical properties

	Cartridge		Detection buffer
	Test strip	Housing	
State	No data available	Natural color or pigment processed color	Liquid at 20 °C
Odor	Odorless	Little odorless	No data available
Odor threshold value	No data available	No data available	No data available
pH	No data available	Natural	7.83
Melting/ Freezing point	No data available	180 ~ 200°C/ No data available	No data available
Initial boiling point or boiling range	No data available	No data available	No data available
Flash point	No flash occurred under 90 °C (Closed-cup)	No data available	No flash occurred under 90 °C (Closed-cup)
Evaporating rate	No data available	No data available	No data available
Lower explosion limit	No data available	No data available	No data available
Upper explosion limit	No data available	No data available	No data available
Vapor pressure	No data available	No data available	No data available
Water solubility	No data available	Not soluble in water, soluble in acetone	No data available
Density	No data available	1.11	No data available
Vapor density	No data available	No data available	No data available
Specific gravity	No data available (water=1)	No data available (water=1)	1.01 (water=1)
N-octan/water partition coefficient	No data available	No data available	No data available
Auto-ignition temperature	No data available	No data available	No data available
Decomposition	No data available	No data available	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS401
	Rev. No.	04
	Rev. Date	2019.01.17

temperature			
Viscosity	No data available	No data available	No data available
Molecular weight	No data available	50,000 ~ 250,000	No data available

### X. Stability and Reactivity

	Cartridge		Detection buffer
	Test strip	Housing	
Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful	- Stable at high temperature and pressure	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful
Conditions to avoid	Heat, high temperature	Heat, high temperature above 200°C	Heat, high temperature
Materials to avoid	No data available	Oxidizing agent: fire or explosive risk	No data available
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency conditions		

### XI. Toxicological Information

	Cartridge	Detection buffer
Information on the likely routes of exposure	No data available	No data available
Acute toxicity	No data available	No data available
Skin corrosion/Irritation	No data available	No data available
Serious eye damage/Eye irritation	No data available	No data available
Respiratory/Skin sensitization	No data available	No data available
Germ cell mutagenicity	No data available	No data available
Carcinogenicity	No data available	No data available
Reproductive toxicity	No data available	No data available
Specific target organ toxicity(single exposure)	No data available	No data available
Specific target organ toxicity(repeated exposure)	No data available	No data available
Aspiration hazard	No data available	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS401
	Rev. No.	04
	Rev. Date	2019.01.17

**XII. Ecological Information**

	Cartridge	Detection buffer
Toxicity	No data available	No data available
Persistence and degradability	No data available	No data available
Bio-accumulative potential	No data available	No data available
Mobility in soil	No data available	No data available
Other adverse effects	No data available	No data available

**XIII. Disposal Considerations**

Cartridge	Test strip	Not available
	Housing	Consider the notice on Waste management act
Detection buffer		Not available

**XIV. Transport Information**

IMDG Code	No data available
DRF	No data available
RID	No data available
ADR	No data available
AND	No data available
IATA	Not regulated as dangerous goods

- Not a dangerous good under DOT, IATA, ADR, IMDG or RID.

**XV. Regulatory Information**

Korea Industrial Safety and Health Act	No data available
Korea Hazardous Materials Safety Control Act	No data available
Korea Toxic Chemical Control Act	No data available
Korea Wastes Control Act	No data available
Other internal and foreign acts	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS401
	Rev. No.	04
	Rev. Date	2019.01.17

#### **XVI. Other Notes**

Limitations: The information and recommendations set forth in this MSDS are believed to be correct as of this date.

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## I. General Information

- A. Product Name/Catalogue Number  
: Boditech CRP Control / CFPO-100
- B. Recommended use of the chemical and restriction on use
  - Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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### IV. First Aid Measures

- In case of eye contact  
: Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact  
: Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled  
: Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring
- If swallowed  
: If irritation or signs of toxicity occur, seek medical attention.
- \* The following symptoms may occur: Irreversible eye damage.
- \* Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- \* This material should be considered as being potentially infectious.

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**VI. Accidental Release Measures**

- Personal precautions
  - \* If not danger, stop to leak
  - \* Do not contact with exposed material without protective equipment
  - \* Note the avoid conditions and materials
  - \* Do not breathe dust / fume / gas / mist / vapours / spray
- Environmental precautions
  - : No data available
- Methods and materials for purification and cleaning up
  - : After absorbed exposed materials, clean up with cleaner and water

**VII. Handling and Storage**

- Precaution for safe handling
  - \* Avoid skin friction.
  - \* Note the avoid conditions and materials from fire or flame.
  - \* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
  - \* Note the avoid conditions and materials.
  - \* Store container tightly closed in a well-ventilated area.
- \* Storage temperature: low temperature

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- Engineering management
  - : Good insolation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

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**XI. Toxicological Information**

	Detection buffer
Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available
Aspiration hazard	No data available

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### I. General Information

- A. Product Name/Catalogue Number  
: Boditech CRP Calibrator / CFPO-112
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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GHS classification	Not applicable
IMDG Code	Not applicable

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Component	Synonyms	Classification No.	Concentration
Glass	Silicon dioxide	7631-86-9	> 73%

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#### V. Fire Fighting Measures

- Recommended extinguishing media
  - : Suitable extinguishing agents including CO<sub>2</sub>, WATER SPRAY or regular form.
- Specific hazard from the chemical
  - : Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters
  - : Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters
  - : Chemical resisted protective equipment for fire-fighter.

#### VI. Accidental Release Measures

- Personal precautions
  - \* If not danger, stop to leak
  - \* Do not contact with exposed material without protective equipment
  - \* Note the avoid conditions and materials
  - \* Do not breathe dust / fume / gas / mist / vapours / spray
- Environmental precautions
  - : No data available
- Methods and materials for purification and cleaning up
  - : After absorbed exposed materials, clean up with cleaner and water

#### VII. Handling and Storage

- Precaution for safe handling
  - \* Avoid skin friction.

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS519
	Rev. No.	02
	Rev. Date	2020. 09. 14

- \* Note the avoid conditions and materials from fire or flame.
- \* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
  - \* Note the avoid conditions and materials.
  - \* Store container tightly closed in a well-ventilated area.
- \* Storage temperature: low temperature

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- Engineering management
  - : Good insulation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

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## XII. Ecological Information

Toxicity	No data available
Persistence and degradability	No data available
Bio-accumulative potential	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS519
	Rev. No.	02
	Rev. Date	2020. 09. 14

Mobility in soil	No data available
Other adverse effects	No data available

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**I. General Information**

- A. Product Name/Catalogue Number  
: AFIAS HbA1c / SMFP-28
- B. Recommended use of the chemical and restriction on use
  - Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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[Buffer]

	Component	Synonyms	Classification No.	Amount (%)
	Water	Dihydrogen Oxide	7732-18-5	99.65
	HbA1c peptide	Fluorescence conjugate	No data available	0.34

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Bio-accumulative potential	No data available	No data available
Mobility in soil	No data available	No data available
Other adverse effects	No data available	No data available

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- Not a dangerous good under DOT, IATA, ADR, IMDG, or RID.

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### I. General Information

- A. Product Name/Catalogue Number  
: Boditech HbA1c Calibrator / CFPO-108
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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- \* Note the avoid conditions and materials from fire or flame.
- \* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
- \* Note the avoid conditions and materials.
- \* Store container tightly closed in a well-ventilated area.
- \* Storage temperature: low temperature

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- Engineering management
- : Good insulation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

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## XI. Toxicological Information

Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available
Aspiration hazard	No data available

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## I. General Information

- A. Product Name/Catalogue Number  
: AFIAS Tn-I Plus / SMFP-35
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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Other hazards

	Cartridge	
	Test strip	Housing
		Health: 1
NEPA	No data available	Flammability: 1
		Physical hazard: 0

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### IV. First Aid Measures

- In case of eye contact  
: Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact  
: Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled  
: Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring.
- If swallowed  
: If irritation or signs of toxicity occur, seek medical attention.
- \* The following symptoms may occur: Irreversible eye damage.
- \* Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- \* This material should be considered as being potentially infectious.

### V. Fire Fighting Measures

- Recommended extinguishing media  
: Suitable extinguishing agents including CO<sub>2</sub>, WATER SPRAY or regular form.
- Specific hazard from the chemical  
: Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS433
	Rev. No.	04
	Rev. Date	2020. 09. 14

- Special measure for fire-fighters

: Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.

- Protective Equipment and precautions for Fire fighters

: Chemical resisted protective equipment for fire-fighter.

## VI. Accidental Release Measures

- Personal precautions

\* If not danger, stop to leak

\* Do not contact with exposed material without protective equipment

\* Note the avoid conditions and materials

\* Do not breathe dust / fume / gas / mist / vapours / spray

- Environmental precautions

: No data available

- Methods and materials for purification and cleaning up

: After absorbed exposed materials, clean up with cleaner and water

## VII. Handling and Storage

- Precaution for safe handling

\* Avoid skin friction.

\* Note the avoid conditions and materials from fire or flame.

\* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask

- Method for safe handling

\* Note the avoid conditions and materials.

\* Store container tightly closed in a well-ventilated area.

\* Storage temperature: low temperature

## VIII. Exposure controls and Personal protection

- Engineering management

: Good insulation and ventilation should be sufficient to control airborne levels.

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## XI. Toxicological Information

Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available
Aspiration hazard	No data available

## XII. Ecological Information

Toxicity	No data available
----------	-------------------

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS433
	Rev. No.	04
	Rev. Date	2020. 09. 14

Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	No data available

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### I. General Information

- A. Product Name/Catalogue Number  
: Boditech Tn-I Plus Control / CFPO-212
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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### VI. Accidental Release Measures

- Personal precautions
- \* If not danger, stop to leak
- \* Do not contact with exposed material without protective equipment
- \* Note the avoid conditions and materials
- \* Do not breathe dust / fume / gas / mist / vapours / spray

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDSC42
	Rev. No.	00
	Rev. Date	2018. 01. 26

- Environmental precautions
  - : No data available
- Methods and materials for purification and cleaning up
  - : After absorbed exposed materials, clean up with cleaner and water

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### XI. Toxicological Information

	Detection buffer
Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDSC42
	Rev. No.	00
	Rev. Date	2018. 01. 26

Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available
Aspiration hazard	No data available

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## I. General Information

- A. Product Name/Catalogue Number  
: Boditech Tn-I Plus Calibrator / CFPO-213
- B. Recommended use of the chemical and restriction on use
  - Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

## II. Hazard Identification

GHS classification of the product and national or regional	
GHS classification	Not applicable
IMDG Code	Not applicable
GHS label elements including precautionary	
Symbols	No data available
Signal word	No data available
Hazard statements	No data available
Precautionary statements	
Prevention	No data available
Reaction	No data available
Storage	No data available
Disposal	No data available
Other hazards	
NEPA	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS525
	Rev. No.	01
	Rev. Date	2020. 03. 05

### III. Composition/Information on Ingredients

Component	CAS Number	Weight %
Bovine Serum Albumin	N/A	95%
Sucrose	57-50-1	4.9%
Sodium Azide	26628-22-8	0.1%

We recommend handling all chemicals with caution.

### IV. First Aid Measures

- In case of eye contact
  - : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
  - : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
  - : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring
- If swallowed
  - : If irritation or signs of toxicity occur, seek medical attention.
- \* The following symptoms may occur: Irreversible eye damage.
- \* Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- \* This material should be considered as being potentially infectious.

### V. Fire Fighting Measures

- Recommended extinguishing media
  - : Suitable extinguishing agents including CO2, WATER SPRAY or regular form.
- Specific hazard from the chemical
  - : Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters
  - : Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters
  - : Chemical resisted protective equipment for fire-fighter.

### VI. Accidental Release Measures

- Personal precautions
  - \* If not danger, stop to leak
  - \* Do not contact with exposed material without protective equipment
  - \* Note the avoid conditions and materials
  - \* Do not breathe dust / fume / gas / mist / vapours / spray

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS525
	Rev. No.	01
	Rev. Date	2020. 03. 05

- Environmental precautions  
: No data available
- Methods and materials for purification and cleaning up  
: After absorbed exposed materials, clean up with cleaner and water

### VII. Handling and Storage

- Precaution for safe handling
  - \* Avoid skin friction.
  - \* Note the avoid conditions and materials from fire or flame.
  - \* Protective equipment: Chemical resisted protective goggles, gloves, clothes and mask
- Method for safe handling
  - \* Note the avoid conditions and materials.
  - \* Store container tightly closed in a well-ventilated area.
  - \* Storage temperature: low temperature

### VIII. Exposure controls and Personal protection

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Respiratory protection	Chemical resisted respiratory protection
Hand protection	Chemical resisted protective gloves
Eye protection	Chemical resisted eye protection
Skin and body protection	Chemical resisted protective clothes

- Standard of exposure

Chemical	No data available
Biological	No data available

### IX. Physical and Chemical properties

State	Liquid at 20 °C
Odor	No data available
Odor threshold value	No data available
pH	7.83
Melting/Freezing point	No data available
Initial boiling point or boiling range	No data available
Flash point	No flash occurred under 90 °C (Closed-cup)
Evaporating rate	No data available
Lower explosion limit	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS525
	Rev. No.	01
	Rev. Date	2020. 03. 05

Upper explosion limit	No data available
Vapor pressure	No data available
Water solubility	No data available
Density	No data available
Vapor density	No data available
Specific gravity	1.01 (water=1)
N-octan/water partition coefficient	No data available
Auto-ignition temperature	No data available
Decomposition temperature	No data available
Viscosity	No data available
Molecular weight	No data available

### X. Stability and Reactivity

Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful
Conditions to avoid	Heat, high temperature
Materials to avoid	No data available
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency

### XI. Toxicological Information

	Detection buffer
Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
Specific target organ toxicity(single exposure)	No data available
Specific target organ toxicity(repeated exposure)	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS525
	Rev. No.	01
	Rev. Date	2020. 03. 05

Aspiration hazard	No data available
-------------------	-------------------

**XII. Ecological Information**

	Detection buffer
Toxicity	No data available
Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	No data available

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## I. General Information

- A. Product Name/Catalogue Number  
: AFIAS Rota/Adeno / SMFP-43
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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	Cartridge
GHS classification	Not applicable
IMDG Code	Not applicable

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Other hazards

		Cartridge	
		Test strip	Housing
NEPA		No data available	Health: 1 Flammability: 1 Physical hazard: 0

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**IV. First Aid Measures**

- In case of eye contact
  - : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact
  - : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled
  - : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring.
- If swallowed
  - : If irritation or signs of toxicity occur, seek medical attention.
- \* The following symptoms may occur: Irreversible eye damage.
- \* Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- \* This material should be considered as being potentially infectious.

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS438
	Rev. No.	02
	Rev. Date	2020. 09. 14

### V. Fire Fighting Measures

- Recommended extinguishing media  
: Suitable extinguishing agents including CO<sub>2</sub>, WATER SPRAY or regular form.
- Specific hazard from the chemical  
: Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters  
: Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters  
: Chemical resisted protective equipment for fire-fighter.

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### VIII. Exposure controls and Personal protection

- Engineering management  
: Good insolation and ventilation should be sufficient to control airborne levels.
- Personal protective equipment

Respiratory protection	Chemical resisted respiratory protection
Hand protection	Chemical resisted protective gloves

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS438
	Rev. No.	02
	Rev. Date	2020. 09. 14

Eye protection	Chemical resisted eye protection
Skin and body protection	Chemical resisted protective clothes
- Standard of exposure	
Chemical	No data available
Biological	No data available

**IX. Physical and Chemical properties**

	Cartridge	
	Test strip	Housing
State	No data available	Natural color or pigment processed color
Odor	Odorless	Little odorless
Odor threshold value	No data available	No data available
pH	No data available	Natural
Melting/ Freezing point	No data available	180 ~ 200°C/ No data available
Initial boiling point or boiling range	No data available	No data available
Flash point	No flash occurred under 90 °C (Closed-cup)	No data available
Evaporating rate	No data available	No data available
Lower explosion limit	No data available	No data available
Upper explosion limit	No data available	No data available
Vapor pressure	No data available	No data available
Water solubility	No data available	Not soluble in water, soluble in acetone
Density	No data available	1.11
Vapor density	No data available	No data available
Specific gravity	No data available (water=1)	No data available (water=1)
N-octan/water partition coefficient	No data available	No data available
Auto-ignition temperature	No data available	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS438
	Rev. No.	02
	Rev. Date	2020. 09. 14

Decomposition temperature	No data available	No data available
Viscosity	No data available	No data available
Molecular weight	No data available	50,000 ~ 250,000

## X. Stability and Reactivity

	Cartridge	
	Test strip	Housing
Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful	- Stable at high temperature and pressure
Conditions to avoid	Heat, high temperature	Heat, high temperature above 200°C
Materials to avoid	No data available	Oxidizing agent: fire or explosive risk
Hazardous decomposition products formed under fire conditions	Irritation or toxic gases and fume may be occurred in fire-emergency conditions	

## XI. Toxicological Information

Information on the likely routes of exposure	No data available
Acute toxicity	No data available
Skin corrosion/Irritation	No data available
Serious eye damage/Eye irritation	No data available
Respiratory/Skin sensitization	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
Specific target organ toxicity (single exposure)	No data available
Specific target organ toxicity (repeated exposure)	No data available
Aspiration hazard	No data available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS438
	Rev. No.	02
	Rev. Date	2020. 09. 14

## XII. Ecological Information

Toxicity	No data available
Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	No data available

## XIII. Disposal Considerations

Cartridge	Test strip	Not available
	Housing	Consider the notice on Waste management act

## XIV. Transport Information

IMDG Code	No data available
DRF	No data available
RID	No data available
ADR	No data available
AND	No data available
IATA	Not regulated as dangerous goods
- Not a dangerous good under DOT, IATA, ADR, IMDG or RID	

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## I. General Information

- A. Product Name/Catalogue Number  
: Boditech Rota/Adeno Control / CFPO-164
- B. Recommended use of the chemical and restriction on use
  - Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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GHS classification	Chronic aquatic toxicity (Category 3), H412
IMDG Code	Not applicable

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Symbols	No data available
Signal word	No data available
Hazard statements	H412 Harmful to aquatic life with long lasting effects.

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**IV. First Aid Measures**

- In case of eye contact  
 : Remove from source of exposure. Wash with copious amounts of water for at least 15 minutes. If irritation or signs of toxicity occur, seek medical attention.
- In case of skin contact  
 : Remove from source of exposure. Wash affected area with soap and water. If irritation or signs of toxicity occur, seek medical attention.
- If inhaled  
 : Get medical attention immediately. Remove to fresh air. If not breathe, give cardiopulmonary resuscitation. If breathing is difficult, give oxygen and continue to monitoring
- If swallowed  
 : If irritation or signs of toxicity occur, seek medical attention.
- \* The following symptoms may occur: Irreversible eye damage.
- \* Medical conditions aggravated by exposure: Pre-existing eye ailments. Hypersensitivity.
- \* This material should be considered as being potentially infectious.

**V. Fire Fighting Measures**

- Recommended extinguishing media  
 : Suitable extinguishing agents including CO2, WATER SPRAY or regular form.
- Specific hazard from the chemical  
 : Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Special measure for fire-fighters  
 : Rescuers need to note the personal protective equipment. If not danger, remove containers. Non-flammable, Corrosive or toxic gases and fume may be occurred in fire-emergency. Inhalation may be harmful.
- Protective Equipment and precautions for Fire fighters  
 : Chemical resisted protective equipment for fire-fighter.

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDSC58
	Rev. No.	01
	Rev. Date	2020. 03. 05

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**XII. Ecological Information**

Toxicity	No data available
Persistence and degradability	No data available
Bio-accumulative potential	No data available
Mobility in soil	No data available
Other adverse effects	Harmful to aquatic life with long lasting effects.

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Not a dangerous good under DOT, IATA, ADR, IMDG, or RID.

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### . General Information

- A. Product Name/Catalogue Number  
: AFIAS COVID-19/Flu Ag Combo / SMFP-77
- B. Recommended use of the chemical and restriction on use
- Recommended use: In vitro diagnostic Medical Device
  - Restriction on use: For in vitro diagnostic use only.
- C. Manufacturer  
: Boditech Med Inc.
- D. Address  
: 43, Geodudanji 1-gil, Dongnae-myeon, Chuncheon-si, Gangwon-do, Korea
- E. Emergency Phone No.  
: +82-33-243-1400

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	Cartridge	Extraction buffer
GHS classification	Not applicable	Not applicable
IMDG Code	Not applicable	Not applicable

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Other hazards

	Cartridge		Extraction buffer
	Test strip	Housing	
NEPA	No data available	Health: 1 Flammability: 1 Physical hazard: 0	No data available

# MATERIAL SAFETY DATA SHEET

Document No.	BT-MSDS459
Rev. No.	00
Rev. Date	2020. 10. 16

## III. Composition/Information on Ingredients

[Cartridge]

	Component	Synonyms	Classification No.	Amount
Housing	Poly(methyl methacrylate)	PMMA	9011-14-7	55 ~ 70 %
	Acrylonitrile butadiene styrene	ABS	9003-56-9	28 ~ 45 %

[Extraction Buffer]

	Component	Synonyms	Classification No.	Amount (%)
	Water	Dihydrogen Oxide	7732-18-5	> 96

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- \* If not danger, stop to leak
- \* Do not contact with exposed material without protective equipment
- \* Note the avoid conditions and materials
- \* Do not breathe dust / fume / gas / mist / vapours / spray
- Environmental precautions
  - : No data available
- Methods and materials for purification and cleaning up
  - : After absorbed exposed materials, clean up with cleaner and water

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## X. Stability and Reactivity

	Cartridge		Extraction buffer
	Test strip	Housing	
Chemical stability	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful	- Stable at high temperature and pressure	- Irritation or toxic gases and fume may be occurred in fire-emergency - Inhalation may harmful
Conditions to avoid	Heat, high temperature	Heat, high temperature above 200°C	Heat, high temperature
Materials to avoid	No data available	Oxidizing agent: fire or explosive risk	No data available
Hazardous decomposition products formed under fire	Irritation or toxic gases and fume may be occurred in fire-emergency		

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS459
	Rev. No.	00
	Rev. Date	2020. 10. 16

conditions

### XI. Toxicological Information

	Cartridge	Extraction buffer
Information on the likely routes of exposure	No data available	No data available
Acute toxicity	No data available	No data available
Skin corrosion/Irritation	No data available	No data available
Serious eye damage/Eye irritation	No data available	No data available
Respiratory/Skin sensitization	No data available	No data available
Germ cell mutagenicity	No data available	No data available
Carcinogenicity	No data available	No data available
Reproductive toxicity	No data available	No data available
Specific target organ toxicity (single exposure)	No data available	No data available
Specific target organ toxicity (repeated exposure)	No data available	No data available
Aspiration hazard	No data available	No data available

### XII. Ecological Information

	Cartridge	Extraction buffer
Toxicity	No data available	No data available
Persistence and degradability	No data available	No data available
Bio-accumulative potential	No data available	No data available
Mobility in soil	No data available	No data available
Other adverse effects	No data available	No data available

### XIII. Disposal Considerations

Cartridge	Test strip	Not available
	Granule	Not available
	Housing	Consider the notice on Waste management act
Extraction buffer		Not available

<b>MATERIAL SAFETY DATA SHEET</b>	Document No.	BT-MSDS459
	Rev. No.	00
	Rev. Date	2020. 10. 16

#### XIV. Transport Information

IMDG Code	No data available
DRF	No data available
RID	No data available
ADR	No data available
AND	No data available
IATA	Not regulated as dangerous goods
- Not a dangerous good under DOT, IATA, ADR, IMDG, or RID	

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